



Deluxe Automatic Blood Pressure Monitor, Wrist

HCPCS: A4670



RETAIL BOX

- WHO Classification
- Irregular Heart Beat Detection
- One Touch Operation
- Large Easy-to-Read Screen
- Date & Time Display
- Pulse Rate Indicator
- Digital Error Messages
- Travel Bag
- Auto-Off Function
- Memory Stores 120 Readings (2 x 60)
- Average Last 3 Results

PRODUCTS

Description

BP3200	Deluxe Automatic Blood Pressure Monitor, Wrist, 1/ea, 24/carton
--------	---

SPECIFICATIONS

Description

Unit Dimensions	3.31" x 2.52" x 1.14"
Carton Size	14.4" x 13.8" x 8.2"
Measurement Method	Oscillometric Method
Measurement Range	Pressure: 0mmHg~300mmHg Pressure: +/- 3mmHg
Power Supply	2x AAA Batteries
Auto Power Off	3 Minutes After Last Reading
Cuff Size	5"-8.5"
Memory Capacity	2 x 60 Reading with Date & Time
Operation Environment	50°F to 104°F, <85% RH
LCD Size	1.93" x 1.50"
Warranty	2 Years