



G30270A

### Locking Raised Toilet Seats

- » Front clamping mechanism locks the seat onto the toilet, so it won't come loose or shift during use or transfer
- » Wide, contoured surface raises seat level by 5"
- » Fits most toilet bowls from 11" to 14"L
- » Optional padded armrests assist patient on and off the toilet

Item No.	Armrests	Weight Capacity	Seat Elevation	Overall W x L	Opening W x L	Pkg.
G30260	No	350 lbs	5"	20.5 x 18"	8 x 9.5"	3/cs
G30270A	Yes	350 lbs	5"	20.5 x 18"	8 x 9.5"	3/cs



G30300

### Toilet Safety Rail

- » Assist patients with joint pain or at risk of falling
- » Handles adjust and rotate back for a wide range of positions
- » Foam armrests provide a comfortable, sure grip
- » Easy-to-clean aluminum frame mounts securely onto bowl
- » Adjustable bracket range of 18 to 24 inches

Item No.	Description	Weight Capacity	Height Range	Width Range	Overall Depth Arms	Overall Depth Legs	Pkg.
G30300	Aluminum	300 lbs	26 x 31"	18 x 24"	16.75"	12"	2/cs



G30250

### Economy Raised Toilet Seat

- » Makes using the toilet easier for those with trouble bending
- » Simple design made of polypropylene resin is easy to clean
- » Non-slip pads and lip inside toilet offer additional stability

Item No.	Weight Capacity	Overall W x L	Opening W x L	Pkg.
G30250	250 lbs	14.5 x 15.5"	8 x 9.5"	3/cs



MDT9130TL1

Compatible with Alarms sold on page 25.

### Toilet Sensor Pads

- » Compact toilet sensor pad allows caregivers to monitor patients while maintaining their dignity and privacy
- » Position under the toilet seat
- » 30-day warranty
- » Alarm sold separately

Item No.	Pad Size	Pkg.
MDT9130TL1	8 x 6 x 4"	ea
MDT9130TL5	8 x 6 x 4"	5/cs