

Inflatable Rubber Cushion



RETAIL BOX

- Durable, heavy-gauge cushion conforms to the natural contours of the body for comfort
- Cleans with just a damp cloth
- Easy to inflate and deflate
- Product contains natural rubber
- Retail Packaged
- Dimensions: 14" diameter (35.56 cm) when inflated, 3.75" (8.89 cm) high when inflated
- Weight Capacity 300 lbs

PRODUCTS**Description**

RTLPC23346

Inflatable Cushion Rubber Red, 6bx/cs