



**KNOW YOU'RE
COVERED.**

Personal Protective Apparel
Head-to-Toe Protection for Healthcare









PROTECT YOUR PEOPLE AND YOUR BOTTOM LINE.

You want to know you're protected. Our protective apparel goes through rigorous testing and meets strict FDA guidelines. We follow the highest quality standards to make sure you get the protection you need.

As a leader in protective apparel, we also have the most comprehensive product line and low acquisition costs. You can find the exact products you need—and maintain cost efficiencies.

AAMI Level isolation gowns are color-coded by neck binding

CONTENTS

	Recommended Gowns	4-5
	AAMI Isolation Gowns	6-7
	Cover Gowns	8-9
	Protective Gowns	10-11
	Footwear	12-13
	Headwear and Beard Cover	14-15
	Scrubs	16
	Coveralls	17
	Lab Coats	18-19

CHOOSE THE RIGHT GOWN.

Our recommendations for your workplace environment.



A Housewide, Med/Surg and ICU

Our QuickComply™ AAMI Level 2 gowns feature an over-the-head style which improves compliance while reducing donning and removing time. Side ties eliminate the difficulties of tying behind the back, while thumb loop wrists allow for a good gown and glove fit.

NONLV240-Regular/Large
NONLV240XL-X-Large
Additional items on page 7

B Endoscopy and Emergency

AAMI Level 3 gowns are made from comfortable and breathable premium multi-ply material. Features a classic tie waist, hook and loop neck and knit cuffs.

NONLV325-Regular/Large
NONLV325XL-X-Large
Additional items on page 7



A Housewide, Med/Surg and ICU

B Endoscopy and Emergency



C Decontamination

Made from premium polyethylene coated materials that have been tested against ASTM F1670 and ASTM F1671. An over-the-head, open back style reduces donning and removing time, while thumb loop wrists allow for a good gown and glove fit.

NON26276-Regular/Large
 NON26276XL-X-Large
 NON26276XXL-XX-Large
Additional items on page 11

D Chemotherapy

Protective polyethylene coated material tested against ASTM F1670, F1671 and the top 11 chemotherapy drugs per ASTM F739-12. Features hook and loop neck, waist ties and knit cuffs.

NON25457M-Medium
 NON25457-Regular/Large
 NON25457XL-X-Large
 NON25457XXL-XX-Large
Additional items on page 11

C Decontamination

D Chemotherapy



NONLV100

NONLV200

NONLV315

Levels of Barrier Protection – AAMI PB70:2012¹

Choose the right gown by matching the color on the neck binding to the AAMI level on the chart.

	🔴 LEVEL 1	🟢 LEVEL 2		🟡 LEVEL 3	
Test	AATCC 42:2000	AATCC 42:2000	AATCC 127:1998	AATCC 42:2000	AATCC 127:1998
Requirements at 4% AQL	Water impact ≤ 4.5 g	Spray impact ≤ 1.0 g	Hydrostatic Pressure ≥ 20 cm	Spray Impact ≤ 1.0 g	Hydrostatic pressure ≥ 50 cm
Anticipated Risk of Exposure	Low	Moderate		Moderate to High	

AAMI Level 1 Isolation Gown

Recommended for: Med/Surg, Visitors, Contact Precautions

- » Made from multi-ply material
- » Side ties eliminate the difficulties of tying behind the back
- » Choose between over-the-head or neck tie styles

Item No.	Description	Pkg.
Medium Weight, Neck Ties, Elastic Wrists		
NONLV100	Regular/Large, Yellow	100/cs
NONLV100XL	X-Large, Yellow	100/cs
Medium Weight, Over-the-Head, Thumb Loop Wrists		
NONLV140	Regular/Large, Yellow	100/cs
NONLV140XL	X-Large, Yellow	100/cs

AAMI Level 2 Isolation Gown

Recommended for: Med/Surg, ICU, Emergency, Airborne and Contact Precautions

- » Made from multi-ply material with side ties to eliminate the difficulties of tying behind the back
- » Choose between over-the-head or tape tab neck styles
- » Available with elastic wrists or thumb loops

Item No.	Description	Pkg.
Medium Weight, Tape Tab Neck, Elastic Wrists		
NONLV200M	Medium, Yellow	100/cs
NONLV200	Regular/Large, Yellow	100/cs
NONLV200XL	X-Large, Yellow	100/cs
NONLV200XXL	XX-Large, Yellow	100/cs
NONLV200XXXL	XXX-Large, Yellow	100/cs
Medium Weight, Tape Tab Neck, Thumb Loops, Individually Folded, Dispensing Box		
NONLV200TL	Regular/Large, Yellow	40/cs
NONLV200TLXL	X-Large, Yellow	40/cs
Medium Weight, Over-the-Head, Thumb Loops		
NONLV240	Regular/Large, Yellow	100/cs
NONLV240XL	X-Large, Yellow	100/cs
Lightweight, Over-the-Head, Thumb Loops		
NONLV242	Regular/Large, Yellow	100/cs
NONLV242XL	X-Large, Yellow	100/cs

AAMI Level 3 Isolation Gown

Recommended for: ICU, Emergency/Trauma, Airborne, Contact and Droplet Precautions

- » Comfortable and breathable materials with classic tie waist
- » Choose between tape tab or hook and loop neck styles
- » Available with elastic wrists or knit cuffs

Item No.	Description	Pkg.
Heavyweight Multi-Ply, Tape Tab Neck, Elastic Wrists		
NONLV315	Regular/Large, Yellow	100/cs
NONLV315XL	X-Large, Yellow	100/cs
Premium Multi-Ply, Hook and Loop Neck, Knit Cuffs		
NONLV325	Regular/Large, Blue	50/cs
NONLV325XL	X-Large, Blue	50/cs
Microporous, Tape Tab Neck, Knit Cuffs		
NONLV350	Regular/Large, White	50/cs
NONLV350XL	X-Large, White	50/cs

Gown Features



Side ties
Available on AAMI Level 1 and 2 gowns



Neck tie style
Available on AAMI Level 1 gowns



Thumb loop wrists
Available on AAMI Level 1 and 2 gowns



Tape tab neck style
Available on AAMI Level 2 and 3 gowns



Neck ties help ensure the gown remains securely in place.

Classic Cover Gown

Recommended for: Visitors, Food Services

- » Made from lightweight polypropylene and multi-layer materials
- » Generously sized for full coverage and flexibility
- » Side ties eliminate difficulties of tying behind the back
- » Neck ties and elastic wrists

Item No.	Description	Pkg.
----------	-------------	------

Lightweight Polypropylene

CRI4000	Regular/Large, Yellow	50/cs
CRI4000B	Regular/Large, Blue	50/cs
CRI4001	X-Large, Yellow	50/cs
CRI4001G	X-Large, Green	50/cs
CRI4003	XXX-Large, Yellow	50/cs

Lightweight Multi-Layer

NON27236	Regular/Large, Yellow	50/cs
NON27236XL	X-Large, Yellow	50/cs

Medium Weight Multi-Layer

NON27SMS2	Regular/Large, Yellow	100/cs
NON27SMS2XL	X-Large, Yellow	100/cs
NON27SMS2XXXL	XXX-Large, Yellow	100/cs



NON27SMS5



Convenient over-the-head, open back style reduces donning and removing time.

Over-the-Head, Open Back Gown

Recommended for: Visitors, Food Services, Airborne and Contact Precautions

- » Available in lightweight polypropylene, multi-layer and medium weight multi-layer materials
- » Generously sized for coverage and flexibility
- » Over-the-head style reduces donning and removing time
- » Thumb loop wrists

Item No.	Description	Pkg.
Lightweight Polypropylene		
CRI4200	Regular/Large	50/cs
CRI4201	X-Large	50/cs
Lightweight Multi-Layer		
NON27SMS25	Regular/Large	100/cs
NON27SMS25XL	X-Large	100/cs
Medium Weight Multi-Layer		
NON27SMS5	Regular/Large	100/cs
NON27SMS5XL	X-Large	100/cs

Polyethylene Gowns

Recommended for: Decontamination, Emergency, ICU, Food Services, Airborne and Contact Precautions and Droplet Precautions

- » Material tested against ASTM F1670 for blood and bodily fluid penetration
- » Over-the-head style reduces donning and removing time
- » Thumb loop wrists help to prevent sleeves from sliding
- » Available in bulk or individually folded in a dispenser box

Item No.	Description	Pkg.
Polyethylene, Blue		
CRI5000	Regular/Large	75/cs
CRI5010	Universal	75/cs
CRI5001	X-Large	75/cs
CRI5003	XXX-Large	75/cs
Premium Polyethylene, Blue		
NONTH150	Regular/Large	75/cs
NONTH180	Universal	75/cs
NONTH200	X-Large	75/cs
NONTH100B	Regular/Large, Bulk Packaged	100/cs
NONTH200B	X-Large, Bulk Packaged	75/cs

Over-the-Head Microporous Gown

Recommended for: Endoscopy, GI Lab, ICU, Emergency, Maternity

- » Made of microporous breathable material tested against ASTM F1670
- » Over-the-head style reduces donning and removing time
- » Thumb loop wrists prevent sleeve slide
- » Individually folded and in a dispenser box

Item No.	Description	Pkg.
NONTH350	Regular/Large	40/cs
NONTH400	X-Large	40/cs



CRI5010

NONTH350

**LEARN MORE ABOUT
BOUFFANT CAPS
AND SHOE COVERS.**

SEE PAGES 12-14

Polyethylene-Coated Gowns

Recommended for: Decontamination, ICU, Endoscopy, GI Lab

- » Material tested against ASTM F1670 for blood and bodily fluid penetration
- » Neck and side ties
- » Available with elastic wrists or knit cuffs

Item No.	Description	Pkg.
Polyethylene Coated, Yellow		
CRI4010	Elastic Wrists, Regular/Large	50/cs
CRI4011	Elastic Wrists, X-Large	50/cs
CRI4010CF	Knit Cuffs, Regular/Large	50/cs
CRI4011CF	Knit Cuffs, X-Large	50/cs
CRI4012CF	Knit Cuffs, XX-Large	50/cs

Premium Polyethylene Coated, Blue

NON27114	Elastic Wrists, Regular/Large	50/cs
NON27114XL	Elastic Wrists, X-Large	50/cs
NON27116	Knit Cuffs, Regular/Large	50/cs
NON27116XL	Knit Cuffs, X-Large	50/cs
NON27116XXL	Knit Cuffs, XX-Large	50/cs

Chemotherapy Gowns

Recommended for: Oncology, Radiology, Lab, Pharmacy

- » Protective polyethylene coated material tested against ASTM F1670, F1671, and the top 11 chemotherapy drugs per ASTM F739-12
- » Premium Breathable Material tested against ASTM F1671, and the top 11 chemotherapy drugs per ASTM F739-99a
- » Hook and loop neck, waist ties and knit cuffs

Item No.	Description	Pkg.
Protective Polyethylene Coated		
NON25457M	Medium	30/cs
NON25457	Regular/Large	30/cs
NON25457XL	X-Large	30/cs
NON25457XXL	XX-Large	30/cs
Premium Breathable Material		
NON27457M	Medium	30/cs
NON27457	Regular/Large	30/cs
NON27457XL	X-Large	30/cs
NON27457XXL	XX-Large	30/cs

Chemotherapy gowns tested against the following chemotherapy drugs per ASTM F739-99a/ASTM F739-12:

- | | |
|------------------------------|-----------------------|
| » Dacarbazine | » Cisplatin |
| » Carmustine (BiCNU) | » Etoposide (Toposar) |
| » Cyclophosphamide (Cytoxan) | » Paclitaxel (Taxol) |
| » Doxorubicin Hydrochloride | » Thio-Tepa |
| » 5-Fluorouracil | » Mitomycin C |
| | » Vincristine Sulfate |

Over-the-Head Polyethylene-Coated Gowns

Recommended for: Decontamination, Emergency, ICU, Maternity

- » Polyethylene Coated materials tested against ASTM F1670
- » Premium Polyethylene Coated materials tested against ASTM F1670, ASTM F1671
- » Over-the-head style reduces donning and removing time
- » Available with thumb loop wrists or knit cuffs

Item No.	Description	Pkg.
Polyethylene Coated		
NON272790H	Knit Cuffs, Regular/Large, Blue	50/cs
NON272790HXL	Knit Cuffs, X-Large, Blue	50/cs
NON272790HTL	Thumb Loops, Regular/Large, Blue	50/cs
NON272790HTLX	Thumb Loops, X-Large, Blue	50/cs
CRI4500	Thumb Loops, Regular/Large, Yellow	100/cs

Premium Polyethylene Coated

NON26278	Knit Cuffs, Regular/Large, Blue	50/cs
NON26278XL	Knit Cuffs, X-Large, Blue	50/cs
NON26278XXL	Knit Cuffs, XX-Large, Blue	50/cs
NON26276	Thumb Loops, Regular/Large, Blue	50/cs
NON26276XL	Thumb Loops, X-Large, Blue	50/cs
NON26276XXL	Thumb Loops, XX-Large, Blue	50/cs



NON25457



NON29859



CRI2010

Multi-Layer Shoe Covers

- » Made of multi-layer material
- » Heavyweight material features stripe non-skid pattern for increased traction
- » Range of sizes to fit all shoe types

Item No.	Description	Pkg.
NON29858	Heavyweight, Regular/Large	100/bx, 300/cs
NON29859	Heavyweight, X-Large	100/bx, 200/cs
NON28858	Medium Weight, Regular/Large	100/bx, 300/cs
NON28859	Medium Weight, X-Large	100/bx, 200/cs
NON28852	Medium Weight, XX-Large	100/bx, 200/cs
CRI2006	Lightweight, Regular/Large	100/bx, 300/cs
CRI2007	Lightweight, X-Large	100/bx, 200/cs

Polyethylene Shoe Cover

- » Smooth bottom for situations where traction is not essential
- » Made from polyethylene material designed to withstand high levels of fluids

Item No.	Description	Pkg.
CRI2010	One Size Fits Most	1,000/cs



CRI2002



CRI2000



NON27145



NON27145

Spunbond Polypropylene Shoe Covers

- » Made of spunbond polypropylene
- » Choice of non-skid pattern or smooth bottom

Item No.	Description	Pkg.
NON28758	Non-Skid, Regular/Large	100/bx, 300/cs
NON28759	Non-Skid, X-Large	100/bx, 200/cs
NON28752	Non-Skid, XX-Large	100/bx, 200/cs
CRI2002	Non-Skid, Regular/Large	100/bx, 300/cs
CRI2003	Non-Skid, X-Large	100/bx, 200/cs
CRI2000	Smooth Bottom, Regular/Large	100/bx, 300/cs
CRI2001	Smooth Bottom, X-Large	100/bx, 200/cs

Ankle Covers

- » Poly-coated bottom with material tested against blood and bodily fluids per ASTM F1670
- » Multi-layer material above ankle for breathability
- » Non-skid tread for increased traction

Item No.	Description	Pkg.
NON27145	Regular/Large	100/bx, 300/cs
NON27145XL	X-Large	100/bx, 300/cs

Knee-High Boot Covers

- » Premium breathable material tested against viral penetration per ASTM F1671
- » Blue laminate material tested against blood and bodily fluids per ASTM F1670 and viral penetration per ASTM F1671
- » Poly-coated material tested against blood and bodily fluids per ASTM F1670
- » Choose between non-skid bottom or foam strips for increased traction
- » Available with hook and loop closure or elastic ankle

Item No.	Description	Pkg.
NON27348P	Premium Breathable, Foam Strip Bottom, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27348PXL	Premium Breathable, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27200	Blue Laminate, Foam Strip Bottom, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27200XL	Blue Laminate, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27144	Poly-coated, Non-skid, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27144XL	Poly-coated, Non-skid, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27150	Poly-coated, Foam Strip Bottom, Hook and Loop Closure, Regular/Large	50/bx, 150/cs
NON27150XL	Poly-coated, Foam Strip Bottom, Hook and Loop Closure, X-Large	50/bx, 150/cs
NON27143	Poly-coated to the Ankle, SMS Material to the Knee, Non-skid, Elastic, Regular/Large	50/bx, 150/cs
NON27143XL	Poly-coated to the Ankle, SMS Material to the Knee, Non-skid, Elastic, X-Large	50/bx, 150/cs



NON27348P



NON27144



NON27200



NON27143



NON28229



CR11008



NON30233B



NON28226



CR11004



CR11005



NON28233F



NON28230R



CR11009



CR11002

Multi-Layer Bouffant Caps

- » Variety of prints and sizes; all come 100/bx, 500/cs
- » Choice of heavy or medium weight multi-ply or heavyweight single-ply

Item No.	Description	Pkg.
NON30232B	21", Blue	500/cs
NON30233B	24", Blue	500/cs
NON30223B	28", Blue	500/cs
NON28226	24", US Flag	500/cs
NON28228	21", Brush Stroke	500/cs
NON28229	24", Brush Stroke	500/cs
NON28230R	24", Breast Cancer Ribbon	500/cs
NON28233F	24", Butterfly	500/cs
NON28533B	24", Blue, Heavyweight	500/cs
NON28234P	24", Multi-color, Heavyweight Single-ply	500/cs

Polypropylene Bouffant Caps

- » Variety of colors and sizes, all come 100/bx, 500/cs
- » Choice of medium weight polypropylene or lightweight polypropylene

Item No.	Description	Pkg.
CR11000	21", Green	500/cs
CR11001	21", Blue	500/cs
CR11002	21", White	500/cs
CR11003	24", White	500/cs
CR11004	24", Blue	500/cs
CR11005	24", Green	500/cs
CR11008	24", Pink	500/cs
CR11009	24", Red	500/cs
CR11202	24", Blue, Medium Weight	500/cs
CR11203	24", White, Heavyweight	500/cs

NEW!

Surgical headwear that meets AORN Guidelines for Surgical Headwear

Surgical Headwear

- » Better coverage from the back of the neck to sideburn area
- » One size fits most—provides all-day breathable comfort and a secure fit
- » No elastic band in front means cap won't leave a mark
- » Quality, longer-lasting design creates less waste than bouffant caps

Item No.	Description	Pkg.
NON61950	Premium Blue Material	600/cs
NON61990	Scrim Material	500/cs
NON61980	Multi-layer Material	500/cs



NON61980



NON61950



NON61990

Surgeon Caps

- » Available in scrim, multi-layer or spunlace materials
- » One size fits most
- » Ties in the back
- » Polypropylene crown for breathability

Item No.	Description	Pkg.
NON28625	Scrim Material for Absorbency	500/cs
NON28626	Multi-layer Material for Fluid-resistance	500/cs
NON28628	Spunlace Material for Comfort	500/cs



NON28625

Surgeon Hood and Beard Covers

- » Made of multi-layer or polypropylene material
- » One size fits most

Item No.	Description	Pkg.
NONSH100C	Surgeons Hood with Multi-Layer Material, Polypropylene Crown	300/cs
NONSH200	Multi-layer Hood	300/cs
NONSH600	Polypropylene Hood	300/cs
NONSH400	Polypropylene Beard Cover	1,000/cs



NONSH100C



NONSH200



NONSH600



NONSH400

Scrubwear

- » Choice of v-neck or round neck top with chest pocket
- » Choice of elastic waist or drawstring pant with two pockets
- » Made from multi-ply material

Item No.	Description	Pkg.
----------	-------------	------

Scrub Top, Round Neck

NON27212S	Small	30/cs
NON27212M	Medium	30/cs
NON27212L	Large	30/cs
NON27212XL	X-Large	30/cs
NON27212XXL	XX-Large	30/cs
NON27212XXXL	XXX-Large	30/cs

Scrub Pant, Elastic

NON27213S	Small	30/cs
NON27213M	Medium	30/cs
NON27213L	Large	30/cs
NON27213XL	X-Large	30/cs
NON27213XXL	XX-Large	30/cs
NON27213XXXL	XXX-Large	30/cs

Scrub Top, V-Neck

NON27202S	Small	30/cs
NON27202M	Medium	30/cs
NON27202L	Large	30/cs
NON27202XL	X-Large	30/cs
NON27202XXL	XX-Large	30/cs

Scrub Pant, Drawstring

NON27203S	Small	30/cs
NON27203M	Medium	30/cs
NON27203L	Large	30/cs
NON27203XL	X-Large	30/cs
NON27203XXL	XX-Large	30/cs



NON27212M
NON27213M

NON27202L
NON27203L

Microporous Coveralls

Recommended for: Lab, O.R., Maintenance, Maternity

- » Choose between elastic wrists and elastic ankles or straight wrists and straight ankles

Item No.	Description	Pkg.
----------	-------------	------

Elastic Wrists and Ankles

NONCV700L	Large	25/cs
NONCV700XL	X-Large	25/cs
NONCV700XXL	XX-Large	25/cs
NONCV700XXXL	XXX-Large	25/cs

Straight Wrists and Ankles

NONCV300XL	X-Large	25/cs
NONCV300XXL	XX-Large	25/cs

Multi-Layer Coveralls

Recommended for: Lab, O.R., Maintenance

- » Made from multi-layer material
- » Elastic wrist and straight ankle

Item No.	Description	Pkg.
NONCV250XL	X-Large	25/cs
NONCV250XXL	XX-Large	25/cs

Spunbond Coveralls

Recommended for: O.R., Maintenance, Visitors

- » Available in medium weight or heavyweight polypropylene material
- » Choose between elastic wrist and straight ankle or straight wrist and straight ankle

Item No.	Description	Pkg.
----------	-------------	------

Medium Weight, Straight Wrists, Straight Ankles

CR13000	X-Large	25/cs
CR13001	XX-Large	25/cs

Heavyweight, Elastic Wrists, Straight Ankles

NONCV150L	Large	25/cs
NONCV150XL	X-Large	25/cs
NONCV150XXL	XX-Large	25/cs
NONCV150XXXL	XXX-Large	25/cs
NONCV150XXXXL	XXXX-Large	25/cs



NONCV150XL



NONSW500M



NONSW100L



NONRP600M

Classic Lab Coats with Knit Collar and Cuffs

- » Multi-layer material provides comfort and protection
- » Three reinforced pockets
- » Convenient loop for hanging

Item No.	Description	Pkg.
White		
NONSW500S	Small	30/cs
NONSW500M	Medium	30/cs
NONSW500L	Large	30/cs
NONSW500XL	X-Large	30/cs
NONSW500XXL	XX-Large	30/cs
NONSW500XXXL	XXX-Large	30/cs

Blue		
NONSW600S	Small	30/cs
NONSW600M	Medium	30/cs
NONSW600L	Large	30/cs
NONSW600XL	X-Large	30/cs
NONSW600XXL	XX-Large	30/cs

Classic Lab Jackets with Knit Collar and Cuffs

- » Multi-layer material provides comfort and protection
- » Three reinforced pockets
- » Convenient loop for hanging

Item No.	Description	Pkg.
White		
NONRP500S	Small	30/cs
NONRP500M	Medium	30/cs
NONRP500L	Large	30/cs
NONRP500XL	X-Large	30/cs
NONRP500XXL	XX-Large	30/cs

Blue		
NONRP600S	Small	30/cs
NONRP600M	Medium	30/cs
NONRP600L	Large	30/cs
NONRP600XL	X-Large	30/cs
NONRP600XXL	XX-Large	30/cs

Classic Lab Coats with Traditional Collar and Knit Cuffs

- » Multi-layer material provides comfort and protection
- » Three reinforced pockets
- » Convenient loop for hanging

Item No.	Description	Pkg.
White, Snap Front		
NONSW100S	Small	30/cs
NONSW100M	Medium	30/cs
NONSW100L	Large	30/cs
NONSW100XL	X-Large	30/cs
NONSW100XXL	XX-Large	30/cs
NONSW100XXXL	XXX-Large	30/cs

Blue, Snap Front		
NONSW400S	Small	30/cs
NONSW400M	Medium	30/cs
NONSW400L	Large	30/cs
NONSW400XL	X-Large	30/cs
NONSW400XXL	XX-Large	30/cs

Blue, Zipper Front		
NONSW480S	Small	30/cs
NONSW480M	Medium	30/cs
NONSW480L	Large	30/cs
NONSW480XL	X-Large	30/cs
NONSW480XXL	XX-Large	30/cs
NONSW480XXXL	XXX-Large	30/cs

MEDLINE PPE COMPANION™

Protection That's Always By Your Side

Streamline caregiving. Move PPE closer to the point of care. This unique isolation cart has all you need to protect yourself—and your patients.

Isolation Cart

A Precaution Identification

Multilingual instructions show the clinically appropriate level of personal protection.

B Organized Storage

Multilingual labels on each drawer.

C Hand Protection

Three glove sizes. More storage behind.

D Face Protection

At-a-glance supply of face masks with shields and earloop masks, ample storage for N95 respirators and more.

E Body Protection

Gowns dispense one at a time with a touch and a pull.

F Gown Storage

Additional storage for a minimum of 50 gowns.

G Quick Disinfection

Sanitize hands. Disinfect surfaces.



Know you're protected. Contact your Medline Representative for more information.







Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US    

REFERENCE:

1. Considerations for selecting protective clothing used in healthcare for protection against microorganisms in blood and body fluids. Centers for Disease Control and Prevention website. Table 3: ANSI/AAMI PB 70:12 Classification of barrier performance of surgical gowns, other protective apparel, surgical drapes and drape accessories. Available at: <http://www.cdc.gov/niosh/nppt/topics/protectiveclothing/>. Accessed May 3, 2016.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.

© 2016 Medline Industries, Inc. Companion and QuickComply are trademarks and Medline is a registered trademark of Medline Industries, Inc. MKT1662404 / LIT171 / 2.5M / JPG / 3