

Low Bed Overbed Table



- Designed specifically for use with low beds with 3.25" frame to floor clearance
- Walnut woodgrain melamine top
- Top can be raised or lowered in infinite settings between 19.5"- 27.5"
- Tabletop can be raised with the slightest upward pressure
- The frame is constructed with tubular steel and durable, attractive chrome finish
- Swivel casters allow for easy maneuverability

PRODUCTS	Description	
13081	1/cs	
SPECIFICATIONS	Description	
Height (adjustable)	19.5" - 27.5"	
Weight	25 lbs.	
Carton Shipping Weight	28 lbs.	
Weight Limit	100 lbs.	
Base Dimensions	15" × 28"	
Tabletop Dimensions	15" × 30"	
Warranty	Limited Lifetime	