

MOBILITY SCOOTER CATALOG



Get **Out**
There!

PRIDE® **MOBILITY** SCOOTERS
LIVE YOUR BEST®

[pridemobility.com](https://www.pridemobility.com)



Don't let a change in mobility affect your lifestyle. Live your best with Pride® Mobility Scooters and enjoy sleek, sporty style in a high-performance package.

With an array of models and colors to choose from, Pride Mobility Scooters make it easier for you to get out there and enjoy life.



What's Inside

Travel Mobility	3-7	Full-Size Luxury	11-12
Jazzy® Zero Turn	3	Pursuit	11
Jazzy® Passport	3	Pursuit XL.....	11
Go-Chair®	4	Raptor	12
Revo™ 2.0	4	Heavy Duty	12
Go-Go® ES 2	5	Maxima	12
Go-Go® Ultra X	5	Accessories	13
Go-Go Elite Traveller®.....	5	Innovative Technology	14
Go-Go Elite Traveller® Plus.....	6	Specifications	15-18
Go-Go® LX with CTS.....	6	Travel Mobility.....	15-16
Go-Go® Sport.....	6	Full/Mid-Size and Heavy Duty....	17-18
Go-Go® Folding Scooter.....	7		
Mid-Size Luxury	8-10		
Victory® 9.2.....	8		
Victory® 10.2.....	8		
Victory® 9.....	9		
Victory® 10.....	9		
Victory® 10 LX w/CTS	10		
Victory® Sport.....	10		

Choosing the right Pride® Mobility Scooter

When choosing a Pride® Mobility Scooter, specific requirements and conditions can affect your process.

Indoor maneuverability

- Travel mobility to mid-size luxury scooters offer enhanced maneuverability

Stability

- 4-wheel scooters provide increased stability

Indoor and outdoor versatility

- Mid-size and full-size luxury scooters are great for both indoor and outdoor use

Outdoor capabilities

- Mid-size and full-size luxury scooters offer excellent outdoor capabilities

Increased weight capacity

- Mid-size and full-size luxury scooters provide increased weight capacity

Easy portability

- Travel mobility scooters are easy to disassemble, transport and store

Speed and performance

- Mid-size and full-size luxury scooters offer increased speed and performance

TRAVEL MOBILITY
*Performance and
 Convenience on the Go*



Jazzy® ZERO TURN

4-WHEEL SCOOTER

- 300 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- iTurn Technology™ provides 4-wheel stability with 3-wheel maneuverability
- 38" turning radius navigates tight corners and small spaces effortlessly
- Dual motors for added traction in rough terrain
- Standard under seat storage
- Delta tiller with wraparound handles
- Charger port in tiller
- LED lighting package
- USB charger built into the tiller for convenient charging of smart phone and portable devices


Color Choices

INTERCHANGEABLE
SHROUDS*

Jazzy® PASSPORT

- 250 lb. weight capacity
- Maximum speed up to 3.6 mph
- Up to 9.6 miles per charge range
- A compact and lightweight design makes traveling a breeze
- Back seat pocket for extra storage while on-the-go
- Folds in just a few simple steps for effortless transport
- Includes convenient, under seat mesh storage bag, mesh cup holder and 40" lap belt
- Large foot platform to accommodate the user's comfort
- Specialty foam seat design for extreme comfort

Color Choices



Feather-touch disassembly

GOCHAIR®

Features

- #1 for portability
- 300 lb. weight capacity
- One hand feather-touch disassembly
- 18" x 17", 5-point swivel seat with latch release
- Faster speed up to 3.7 mph
- Compact turning radius 25.5"
- Dual swing-away, under seat storage bins
- Height and width adjustable armrests
- Heaviest piece is only 36 lbs.

Color Choices

INTERCHANGEABLE
SHROUDS*

REVO™ 2.0

3 AND 4-WHEEL SCOOTER

Features

- 375 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- CTS Suspension (Comfort-Trac front and rear suspension)
- Up to 12 miles per charge range (3W), 13 miles per charge range (4W)
- Standard under seat storage
- Delta tiller with wraparound handles
- Charger port in tiller
- LED headlight
- Easily charge the batteries on or off the unit
- USB mobile device charging port

Color Choices





GOGO[®] es2

3-WHEEL SCOOTER

Features

- 250 lb. weight capacity
- Two speed settings high and low (indoor/ outdoor), maximum speed up to 4 mph
- Up to 6.8 miles per charge range
- Lightweight, one-piece frame for added stability and easy portability

Color Choices



Black



3-wheel

GOGO[®] ULTRAX

3 & 4-WHEEL SCOOTER

Features

- 260 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- Maximum speed up to 4 mph
- Up to 6.9 miles per charge range (3W), 7.2 per charge range (4W)

Color Choices



Red



4-wheel



3-wheel

GOGO[®] ELITE TRAVELLER[®]

3 & 4-WHEEL SCOOTER

Features

- 300 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- Maximum speed up to 4 mph
- Up to 10.9 miles per charge range (3W), 12 per charge range (4W)
- Charger port in tiller

Color Choices

INTERCHANGEABLE SHROUDS*



Red



Blue



4-wheel



3-wheel

GOGO[®] ELITE TRAVELLER[®] PLUS

3 AND 4-WHEEL SCOOTER

Features

- 300 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- Maximum speed up to 4.5 mph
- Up to 12.7 miles per charge range (3W), 12.7 per charge range (4W)
- Delta tiller with wraparound handles
- Charger port in tiller
- LED headlight

Color Choices

INTERCHANGEABLE SHROUDS*



Red



Blue



4-wheel



3-wheel

GOGO[®] LX WITH CTS SUSPENSION

3 AND 4-WHEEL SCOOTER

Features

- 300 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- CTS Suspension (Comfort-Trac front and rear suspension)
- Maximum speed up to 4.4 mph
- Up to 11.7 miles per charge range (3W), 12.6 miles per charge range (4W)
- Delta tiller with wraparound handles
- Charger port in tiller
- LED headlight

Color Choices

INTERCHANGEABLE SHROUDS*



Red



Blue

Optional



Orion Yellow



Onyx Black



Pearl Pink



Sunburst Orange



4-wheel



3-wheel

GOGO[®] SPORT

3 AND 4-WHEEL SCOOTER

Features

- 325 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into five super lightweight pieces for convenient transport
- Maximum speed up to 4.7 mph
- Up to 15.9 miles per charge range (3W), 14 miles per charge range (4W)
- Delta tiller with wraparound handles
- Charger port in tiller
- LED headlight

Color Choices

INTERCHANGEABLE SHROUDS*



Red



Blue



4-wheel

GO-GO® FOLDING SCOOTER

4-WHEEL SCOOTER

Features

- 250 lb. weight capacity
- Maximum speed up to 4 mph
- Up to 13 miles per charge range (lithium ion), 9.3 miles per charge range (sealed lead acid)
- Height-adjustable tiller
- Lighting package
- Optional lithium ion battery
- Optional armrests

Color Choices



White



Note: Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go-Go® Folding Scooter and its lithium-ion battery is permitted.



VICTORY® SERIES

*The Ultimate in Style
and Performance®*

NEW



3-wheel



4-wheel

VICTORY® 9.2

3 AND 4-WHEEL SCOOTER

Features

- 300 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speeds up to 5.6 mph
- Up to 17.5 miles per charge
- Delta tiller with wraparound handles
- Frosted LED headlights and low-mounted LED pathway light brightens your route
- Under-tiller puddle light improves visibility to USB charging port and tiller-mounted storage
- Integrated USB mobile device charging port in tiller
- Integrated storage/cup holder in tiller

Color Choices

Candy
Apple RedOcean
Blue

NEW



3-wheel



4-wheel

VICTORY® 10.2

3 AND 4-WHEEL SCOOTER

Features

- 400 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speeds up to 5.6 mph
- Up to 13.5 miles per charge
- Delta tiller with wraparound handles
- Frosted LED headlights and low-mounted LED pathway light brightens your route
- Under-tiller puddle light improves visibility to USB charging port and tiller-mounted storage
- Integrated USB mobile device charging port in tiller
- Integrated storage/cup holder in tiller

Color Choices

Candy
Apple RedOcean
Blue

VICTORY[®] 9

3 AND 4-WHEEL SCOOTER



3-wheel



4-wheel

Features

- 300 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speed up to 5.3 mph
- Up to 13 miles per charge range
- Delta tiller with wraparound handles
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Easy access tiller-mounted charger port

Color Choices



Candy Apple Red



Viper Blue

VICTORY[®] 10

3 AND 4-WHEEL SCOOTER



3-wheel



4-wheel

Features

- 400 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speed up to 5.3 mph
- Up to 16 miles per charge range
- Delta tiller with wraparound handles
- Long lasting, high intensity LED headlight offers optimal pathway illumination
- Easy access tiller-mounted charger port

Color Choices



Candy Apple Red



Viper Blue

VICTORY[®] 10

LX WITH CTS SUSPENSION

4-WHEEL SCOOTER



Features

- 400 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- CTS Suspension (Comfort-Trac front and rear suspension)
- Maximum speed up to 5.6 mph
- Up to 16 miles per charge range
- Delta tiller with wraparound handles
- Black vinyl reclining high-back seat with sliders, height adjustable headrest, flip-up, height and width adjustable armrests
- Full LED lighting package that includes curb and rear lights and directional signals
- Easy access tiller-mounted charger port

Color Choices



Candy Apple Red



True Blue

VICTORY[®] SPORT

4-WHEEL SCOOTER



Features

- 350 lb. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speed up to 8 mph
- Up to 11 miles per charge range
- Delta tiller with wraparound handles
- Black vinyl reclining high-back seat with sliders, height adjustable headrest, flip-up, height and width adjustable armrests
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination with directional signals
- Front and rear suspension
- Easy access tiller-mounted charger port

Color Choices



Candy Apple Red



Viper Blue



FULL-SIZE LUXURY

*Outstanding Outdoor
Performance with
Luxurious Features*



4-WHEEL SCOOTER

Features

- 400 lb. weight capacity
- Infinitely-adjustable tiller angle with conveniently located handle
- Maximum speed up to 8 mph
- Front and rear suspension
- Up to 14 miles per charge range (NF-22), 19 miles per charge range (G24)
- Feather-touch disassembly
- Delta tiller with wraparound handles
- Full lighting package that includes directional signals and lower headlight
- Reclining high-back seat with sliders and headrest

Color Choices



Candy
Apple Red



Silver



4-WHEEL SCOOTER

Features

- 400 lb. weight capacity
- Infinitely-adjustable tiller angle with conveniently located handle
- Maximum speed up to 7.5 mph
- Up to 30 miles per charge range (G27), 22 miles per charge range (G24)
- Front and rear suspension
- Delta tiller with wraparound handles
- Deluxe reclining high-back seat with sliders and headrest
- 13" low-profile, pneumatic tires for excellent outdoor performance

Color Choices



Candy
Apple Red



Silver



FULL-SIZE LUXURY



3-WHEEL SCOOTER

Features

- 400 lb. weight capacity
- Front and rear suspension
- Maximum speed up to 14 mph
- Full LED lighting package with directional signals and brake lights that offer full visibility for increased safety
- Hand brake with regenerative braking system
- Secure, lockable under seat storage compartment with seat support rod
- Easy to read LCD console which shows estimated speed, mileage and battery display
- Limited recline fold-flat seat with manual adjustment and angle/width adjustable armrests
- Sleek and solid design for enhanced stability
- Front storage compartment located on the tiller
- Tubeless pneumatic tires for a smooth ride

Color Choices



Silver



3 AND 4-WHEEL SCOOTER

HEAVY-DUTY

*Durable Performance at
an Exceptional Value*

Features

- 500 lb. weight capacity
- Ultra heavy-duty drivetrain with 110 amp controller
- Maximum speed up to 5.3 mph
- Easily adjustable tiller
- Up to 18 miles per charge range (3W, NF-22), 9.6 miles per charge range (power elevating seat, U-1), 16 miles per charge range (4W)
- Durable and stylish rear polyurethane bumper
- User-friendly console

Color Choices



Candy
Apple Red



Viper
Blue



ACCESSORIES

A wide selection of Pride Mobility Scooter accessories are available to personalize your Pride Mobility Scooter to better meet your needs and make it uniquely yours.



Cane clip attachment



Cup holder



Single cane holder



Safety flag



Double cane holder



Under seat storage bins



Rear basket



XLR USB charger



Oxygen holder



RAM® X-Grip®
cell phone holder



Wishbone crutch holder



Additional battery box
(12 amp)



Lapbelt



Saddlebag



Walker holder



Additional battery box
(18 amp)



Rearview mirror



Weather breaker



Weather cover

Military Patches Available



DISCLAIMER: Not all accessories are available for all products. Please speak with your Provider for details and availability.

INNOVATIVE TECHNOLOGY



Patented, intelligent-turning iTurn Technology™ offers a 38" turning radius to ensure effortless navigation within tight corners and small spaces.



CTS SUSPENSION



Front Suspension



Rear Suspension

Comfort-Trac Suspension (CTS) is the next generation in mobility scooter ride comfort. Pride's **industry-first** technology incorporates advanced suspension components for improved absorption over rough surfaces, and a smoother ride experience over various terrains.

FEATHER-TOUCH DISASSEMBLY



Scooter/travel mobility
disassembly



Go Chair® disassembly

Most Go-Go® Travel Mobility and Pride® Mid-Size Scooters feature **Feather-Touch** — one-hand disassembly for hassle-free transport and storage.

SPECIFICATIONS

TRAVEL MOBILITY

Get Out There!

	Jazzy® Zero Turn	Revo™ 2.0	Jazzy® Passport	Go Chair®	Go-Go® ES 2	Go-Go® Ultra X	Go-Go Elite Traveller®	Go-Go Elite Traveller® Plus	Go-Go® Sport	Go-Go® LX w/ CTS Suspension	Go-Go® Folding Scooter
FRONT TIRES	2" x 7" SOLID	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	8" SOLID (CASTER)	6" SOLID (CASTER)	2" x 7.5" SOLID	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2" x 8" SOLID	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	1.5" x 7" SOLID
REAR TIRES	2.5" x 9" SOLID	3" x 9" SOLID	9" SOLID (PNEUMATIC DRIVE)	9" SOLID (DRIVE)	3" x 7.5" SOLID	2.5" x 8" SOLID	2.5" x 8" SOLID	3" x 9" SOLID	3" x 9" SOLID	2.5" x 8" SOLID	2" x 8" SOLID
MAXIMUM SPEED²	UP TO 4 MPH	UP TO 5 MPH	UP TO 3.6 MPH	UP TO 3.7 MPH	UP TO 4 MPH	UP TO 4 MPH	UP TO 4 MPH	UP TO 4.5 MPH	UP TO 4.7 MPH	UP TO 4.4 MPH	UP TO 4 MPH
GROUND CLEARANCE^{4,11,14}	1.4" at motor 1.5" at anti-tip 3.5" at side deck	1.25" at transaxle 3.25" at mid-frame (3W) 3.5" at mid-frame (4W)	3" at motor	1.8" at motor 1.7" at anti-tip 2.7" at foot plate bracket	1.5" at motor 0.75" at anti-tip	1.3" at transaxle 0.67" at anti-tip 3.1" at deck (3W) 3" at deck (4W)	1.38" 0.75" at anti-tip 2.25" at mid-frame	1.7" at transaxle (3W) 0.83" at anti-tip (3W) 2.6" at side deck (3W)	1.5" at motor 0.75" at anti-tip 1.25" at rear frame lock (3W) 3" at controller (4W)	1.5" at motor 0.75" at anti-tip 2.5" at frame	2.3"
3-W TURNING RADIUS⁴	N/A	40"	23.5"	N/A	31"	33"	33"	34.5"	36"	33"	32"
4-W TURNING RADIUS⁴	38.25"	52"	N/A	25.5" with foot platform	N/A	44"	44"	45.75"	50.5"	44.75"	N/A
3-W LENGTH/WIDTH^{4,10}	N/A	46"/20.75"	N/A	N/A	36.75"/19.5"	37"/19.5"	37"/19.5"	40"/21.25"	42"/21.25"	39.75"/20.5"	37"/19"
4-W LENGTH/WIDTH^{4,10}	45.75"/20"	47"/20.75"	40"/24.25"	28.5" without footplate 33"/25.5"	N/A	39.5"/19.5"	39.5"/19.5"	40.25"/21.25"	43"/21.25"	41"/20.5"	N/A
SEAT-TO-GROUND HEIGHT RANGE	21.5"-23.5"	21"-23.5"	19"	22"-23"	19"-20"	18.25"-20.25"	18.25"-20.25"	20.5"-22"	20.5"-22"	20.75"-22.25"	20.5"
SEAT TYPE (STANDARD)	COMPACT, FOLDABLE (17"X17")	FOLDABLE, MOLDED PLASTIC (20"X18")	COMPACT, FOLDABLE (17.5"X16.5")	COMPACT, FOLDABLE (18"X17")	COMPACT, FOLDABLE (17"X15")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (16"X13")
(3/4-WHEEL) WEIGHT^{3,4,12}	120 LBS.	120 LBS./132 LBS.	54 LBS.	67 LBS.	58 LBS.	76 LBS./80 LBS.	76 LBS./80 LBS.	86 LBS./87 LBS.	93 LBS./99 LBS.	85 LBS./87 LBS.	51 LBS.
SEAT-TO-DECK HEIGHT RANGE	16"-18"	15"-18"	14.5"	16.25"-17.25"	16.25"-17.25"	15.75"-17.25"	15.75"-17.25"	16.25"-17.75"	17"-18"	16.25"-17.25"	13"
BATTERY WEIGHT⁹	13.5 LBS. EA.	25 LBS. EA.	7 LBS. EA.	14 LBS. EA.	9 LBS. EA.	9 LBS. EA. ⁵	9 LBS. EA. ⁵ / 14 LBS. EA. ⁷	14 LBS. EA. ⁷	14 LBS. EA. ⁷	9 LBS. EA. ⁵ / 14 LBS. EA. ⁷	9 LBS. EA. (SEALED LEAD ACID / 5 LBS. EA. (LITHIUM-ION)
BATTERY REQUIREMENTS⁶	(2) 22AH	(2) U1	(1) 24V, 18AH (LITHIUM ION)	(2) 18AH	(2) 12AH	(2) 12AH	(2) 12AH OR (2) 18AH	(2) 18AH	(2) 18AH	(2) 12AH OR (2) 18AH	(2) 12AH OR (1) 12AH (LITHIUM-ION)
PER CHARGE RANGE (MAX.)^{2,8}	UP TO 10 MILES	3W: UP TO 12 MILES 4W: UP TO 13 MILES	UP TO 9.6 MILES	UP TO 8.7 MILES (STANDARD) 10.9 MILES with extended range	UP TO 6.8 MILES	3W: UP TO 6.9 MILES 4W: UP TO 7.2 MILES	3W: UP TO 8.5 MILES (12AH) UP TO 10.9 MILES (18AH) 4W: UP TO 8.2 MILES (12AH) UP TO 12 MILES (18AH)	3W: UP TO 12.7 MILES 4W: UP TO 12.7 MILES	3W: UP TO 15.9 MILES 4W: UP TO 14 MILES	3W: UP TO 8 MILES (12AH) UP TO 11.7 MILES (18AH) 4W: UP TO 8.6 MILES (12AH) UP TO 12.6 MILES (18AH)	UP TO 13 MILES (LITHIUM-ION) UP TO 9.3 MILES (SEALED LEAD ACID)
BATTERY CHARGER	OFF-BOARD, 3.5A	OFF-BOARD, 3.5A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 3.5A	OFF-BOARD, 3.5A	OFF-BOARD, 2A	OFF-BOARD, 2A (SEALED LEAD ACID) OFF-BOARD, 2.5A (LITHIUM-ION)
SUSPENSION	FRONT & REAR	FRONT & REAR	LIMITED	LIMITED	N/A	N/A	N/A	N/A	N/A	FRONT & REAR	N/A
MAXIMUM WEIGHT CAPACITY	300 LBS.	375 LBS.	250 LBS.	300 LBS.	250 LBS.	260 LBS.	300 LBS.	300 LBS.	325 LBS.	300 LBS.	250 LBS.
STANDARD BODY COLORS	FIRE OPAL, GREEN MOONSTONE, ONYX BLACK, PEWTER, PEARL WHITE, SAPPHIRE BLUE	GREY STREET, TRUE BLUE	PEWTER	RUBY RED, SAPPHIRE BLUE, CITRINE YELLOW, AMBER ORANGE, PEARL WHITE, ROSE QUARTZ	BLACK	RED	RED, BLUE ¹	RED, BLUE ¹	RED, BLUE ¹	RED, BLUE ¹	WHITE
WEIGHT OF HEAVIEST PIECE⁴	48 LBS. (REAR SECTION)	48 LBS. (3-WHEEL) (REAR SECTION) 57 LBS. (4-WHEEL) (FRONT SECTION)	N/A	36 LBS. (REAR SECTION)	41 LBS. (3-WHEEL) (SCOOTER BASE)	29 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	29 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	32 LBS. (3-WHEEL) (REAR SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	35 LBS. (3-WHEEL) (FRONT SECTION) 44 LBS. (4-WHEEL) (FRONT SECTION)	34 LBS. (3-WHEEL) (FRONT SECTION) 37 LBS. (4-WHEEL) (FRONT SECTION)	51 LBS. (UNIT DOES NOT DISASSEMBLE - UNIT WEIGHT WITHOUT BATTERY PACK)
HCPCS	K0899	K0899	K0899	K0899	K0800 - GROUP 1 STANDARD (3W)	K0899	K0800 - GROUP 1 STANDARD	K0806 - GROUP 2 STANDARD ¹³	K0801 - GROUP 1 HEAVY DUTY	K0899	K0899

1. All shroud panels included.
2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition.
This specification can be subject to a variance of (+ or -) 10%.
3. Battery weight not included.
4. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
5. Battery pack is 18 lbs. with 12AH batteries.
6. Deep cycle (AGM or gel cell type recommended).
7. 18 AH battery pack weighs 29.5 lbs.
8. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance.
Test conducted at maximum weight capacity.

9. Battery weight may vary based on manufacturer.
10. Width of base (not including seat).
11. Does not account for anti-tip clearance.
12. Includes standard seat weight.
13. Must be billed to Medicare as either a beneficiary requested upgrade with a properly executed ABN, or as an upgrade at no additional charge to the beneficiary.
14. At midframe.
15. Light package.

SPECIFICATIONS

MID/FULL-SIZE, HEAVY DUTY

Get Out There!

	Victory® 9	Victory® 10	Victory® 9.2	Victory® 10.2	Victory 10® LX w/ CTS Suspension ¹⁵	Victory® Sport ¹⁵	Maxima	Pursuit® ¹⁵	Pursuit® XL ¹⁵	Raptor ¹⁵
FRONT TIRES	3W: 3" x 8" SOLID 4W: 3" X 9" SOLID	3.5" x 10" SOLID	3W: 3" x 8" SOLID 4W: 3" X 9" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID	3" x 10" SOLID	3" x 13" SOLID	13" PNEUMATIC	3" x 16" PNEUMATIC
REAR TIRES	3" x 9" SOLID	3.5" x 10" SOLID	3" x 9" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID	4" x 10" SOLID	13" SOLID	13" PNEUMATIC	3" x 16" PNEUMATIC
MAXIMUM SPEED²	UP TO 5.3 MPH	UP TO 5.3 MPH	UP TO 5.6 MPH	UP TO 5.2 MPH	UP TO 5.6 MPH	UP TO 8 MPH	UP TO 5.3 MPH	UP TO 8 MPH	UP TO 7.5 MPH	UP TO 14 MPH
GROUND CLEARANCE^{4,11,14}	1.8" at transaxle 1.4" at anti-tip 2.7" at side deck (3W) 2.1" at motor bracket bolt (4W)	2.2" at motor (3W) 1.26" at anti-tip (3W) 3.4" at side deck (3W) 2" at motor (4W) 1.5" at anti-tip (4W) 3.1 at side deck (4W)	1.8" at transaxle 1.4" at anti-tip 2.7" at side deck	2.4" at transaxle 2.1" at anti-tip 3.0" at side deck	2.4" at transaxle 0.9" at anti-tip 2.9" at front frame	2.25" at motor 1" at anti-tip 3" at rear lock up	2.5" at motor 2" at anti-tip (3W) 4.5" at side deck (3W) 1.5" at anti-tip (4W) 3.25" (4W)	3.7" at motor 1.8" at anti-tip 4.5" at side deck	2.75" at motor 5.25" at mid-frame	3.5" at center of frame 1.37" at anti-tip 6.18"
3-W TURNING RADIUS⁴	40.25"	42"	40.25"	42"	N/A	N/A	50.5"	N/A	N/A	60"
4-W TURNING RADIUS⁴	51.5"	54"	52"	52.25"	60"	51.75"	54"	67.75" LEFT 66.75" RIGHT	82.5"	N/A
3-W LENGTH/ WIDTH^{4,10}	43.5"/22.25"	45.25"/22.25"	43.5"/22.25"	45.25"/22.25"	N/A	N/A	46"/25.5"	N/A	N/A	59"/24.75"
4-W LENGTH/ WIDTH^{4,10}	45.5"/22.5"	47"/22.25"	46"/22.25"	47"/22.25"	47"/22.25"	47"/26"	47"/25.5"	53.25"/25.25"	56"/27.25"	N/A
SEAT-TO-GROUND HEIGHT RANGE	19.25"-20.75"	20.5"-22"	21.5"-24.5"	22"-25"	23.5"- 25.5"	21.5"-24.25"	22.5"-24.5"	24.75"-26.25"	26.25"-27.75"	28"
SEAT TYPE (STANDARD)	FOLDABLE, MOLDED PLASTIC (18"X17")	FOLDABLE, MOLDED PLASTIC (18"X17")	FOLDABLE, MOLDED PLASTIC (18"X17")	FOLDABLE, MOLDED PLASTIC (18"X17")	RECLINING H/B W/SLIDERS & HEADREST (18"X18")	RECLINING H/B W/SLIDERS & HEADREST (18"X18")	DELUXE CONTOURED MEDIUM-BACK (22"X18")	RECLINING H/B W/SLIDERS & HEADREST (18"X18")	HIGH-BACK PILLOW TOP W/SLIDERS & HEADREST (19.5"X16")	HIGH BACK (20" X 18") (100" LIMITED RECLINE)
(3/4-WHEEL) WEIGHT^{3,4,12}	117 LBS./131 LBS.	133 LBS./137 LBS.	125.75LBS./136 LBS.	130.25 LBS./139 LBS.	179 LBS.	173 LBS.	171 LBS./176 LBS.	225 LBS.	254 LBS.	208 LBS.
SEAT-TO-DECK HEIGHT RANGE	17.25"-18.75"	17.5"-19"	15.25"-18.25"	15.25" - 18.25"	16.5"-18.5"	16"-18.5"	3W: 16"-18" 4W: 16"-18"	17"-18.75"	21"-22.5"	19"
BATTERY WEIGHT⁹	24 LBS. EA	24 LBS. EA/ 32 LBS. EA.	24.5 LBS. EA	24 LBS. EA/ 32 LBS. EA.	32 LBS. EA.	32 LBS. EA.	38 LBS. EA. (NF-22 OR 55 AH) 24 LBS. EA. (U-1)	38 LBS. EA. (NF-22 OR 55 AH) 53 LBS. EA. (75 AH)	53 LBS. EA. (75 AH) 60 LBS. EA. (100 AH)	53 LBS. EA.
BATTERY REQUIREMENTS⁶	(2) U1	(2) U1 OR (2) 40AH	(2) U1	(2) U1 OR (2) 40AH	(2) 40AH	(2) 40AH	(2) NF-22 (2) U-1 (POWER ELEVATING SEAT)	(2) NF-22 OR (2) GROUP 24	(2) GROUP 27 OR (2) GROUP 24	(2) 75AH
PER CHARGE RANGE (MAX.)^{2,8}	UP TO 13 MILES	UP TO 16 MILES	3W: UP TO 17.5 MILES 4W: UP TO 16.5 MILES	3W4W: UP TO 11 MILES (U-1) 3W: UP TO 13 MILES (40AH) 4W: UP TO 13.5 MILES (40AH)	UP TO 16 MILES	UP TO 11 MILES	3W: UP TO 18 MILES, U-1 (POWER ELEVATING SEAT) UP TO 9.6 MILES 4W: UP TO 16 MILES	G24: UP TO 19 MILES NF-22: UP TO 14 MILES	G27: UP TO 30 MILES G24: UP TO 22 MILES	19 MILES (UP TO 400 LBS. AND FULL SPEED) 31 MILES (UP TO 220 LBS. AND AT LOW SPEED)
BATTERY CHARGER	OFF-BOARD, 3.5A (U-1)	OFF-BOARD, 3.5A (U-1/ OFF-BOARD, 5A (40 AH)	OFF-BOARD, 3.5A (U-1/ OFF-BOARD, 5A (40 AH)	OFF-BOARD, 3.5A (U-1/ OFF-BOARD, 5A (40 AH)	OFF-BOARD, 5A (40 AH)	OFF-BOARD, 5A	OFF-BOARD, 5A	OFF-BOARD, 5A (NF-22) OFF-BOARD, 8A (GROUP 24)	OFF-BOARD, 8A	OFF-BOARD, 8A
SUSPENSION	N/A	N/A	N/A	N/A	FRONT & REAR	FRONT & REAR	3-W REAR/ 4-W FRONT & REAR	FRONT & REAR	FRONT & REAR	FRONT & REAR
MAXIMUM WEIGHT CAPACITY	300 LBS.	400 LBS.	300 LBS.	400 LBS.	400 LBS.	350 LBS.	500 LBS.	400 LBS.	400 LBS.	400 LBS.
STANDARD BODY COLORS	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, TRUE BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, SILVER	CANDY APPLE RED, SILVER	SILVER
WEIGHT OF HEAVIEST PIECE⁴	48 LBS. (3-WHEEL) (FRONT SECTION) 59 LBS. (4-WHEEL) (FRONT SECTION)	54 LBS. (3-WHEEL) (FRONT SECTION) 61 LBS. (4-WHEEL) (FRONT SECTION)	51.75 LBS. (3-WHEEL) (FRONT SECTION) 61.5 LBS. (4-WHEEL) (FRONT SECTION)	56.5 LBS. (3-WHEEL) (FRONT SECTION) 64.5 LBS. (4-WHEEL) (FRONT SECTION)	78 LBS. (4-WHEEL) (FRONT SECTION)	68 LBS. (4-WHEEL) (FRONT SECTION)	76 LBS. (3-WHEEL) (REAR SECTION) 76 LBS. (4-WHEEL) (REAR SECTION)	95 LBS. (FRONT SECTION)	205 LBS. (MAIN FRAME)	N/A (UNIT DOES NOT DISASSEMBLE)
HCPCS	K0800 - GROUP 1 STANDARD	K0801 - GROUP 1 HEAVY DUTY	K0800 - GROUP 1 STANDARD	K0801 - GROUP 1 HEAVY DUTY	K0899	K0899	K0802 - GROUP 1 VERY HEAVY DUTY (3W) K0808 - GROUP 2 ¹⁴ VERY HEAVY DUTY (4W)	K0899	K0899	K0899

- All shroud panels included.
- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- Battery weight not included.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
- Battery pack is 18 lbs. with 12AH batteries.
- Deep cycle (AGM or gel cell type recommended).
- 18 AH battery pack weighs 29.5 lbs.
- Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

- Battery weight may vary based on manufacturer.
- Width of base (not including seat).
- Does not account for anti-tip clearance.
- Includes standard seat weight.
- Must be billed to Medicare as either a beneficiary requested upgrade with a properly executed ABN, or as an upgrade at no additional charge to the beneficiary.
- At midframe.
- Light package.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement.

To find a compatible lift for your vehicle/mobility product, use our new

LIFT FINDER

pridemobility.com/lift-finder



PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

AUTHORIZED PRIDE® PROVIDER

PROVIDER INFO