

## Med-Aire Plus 8" Alternating Pressure and Low Air Loss Mattress System with 10" Defined Perimeter



INSET A

- Active combination therapy system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control
- Designed to prevent, treat and heal pressure ulcers in the home or long term care setting
- Digital control unit offers 4 cycle times and adjustable comfort weight settings for optimal immersion and patient satisfaction
- 3 diagnostic audio and visual alarms alert for issues with air pressure, the system, or power
- 10" deep static perimeter tubing surrounds the length of the mattress preventing patient falls or entrapment
- Pillow feature provides head stability & comfort
- Advanced Clinical and Care Giver Features: Auto-Firm, Seat Inflate, microprocessing feedback, lock-out and more!
- Cell-on-cell designed air cells prevent "bottoming out" and provide up to 24 hours of power outage protection
- 8" deep air cells are easily removed with quick disconnects for convenient servicing
- Removable 4 way stretch PU cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design

PRODUCTS	Description
14029DP	Alternating Pressure/Low Air Loss Mattress System, E0277, 1/ea

SPECIFICATIONS	Description
Control Unit Airflow	8 LPM
Control Unit Cycle Time	10, 15, 20, 25 Minutes
Control Unit Type	Compressor
Control Unit Alarms	Visual / Audible
Mattress Cover Material	Fluid-Resistant Stretch
Weight Capacity	450 lbs.
Warranty	18 Month Limited

ACCESSORIES	Description
14029DPC	Cover, 1/ea
14029XP	Pump Only for 14029/14029DP, 1/ea
14029DP-CELL	Replacement Bladder for 14029DP, 1/ea
14029DP-CELL 1	Replacement Bladder with Laser Holes for 14029DP, 1/ea
14427	Electric Pump (Sold Separately) (Replaces Mason # 8151-EP), 1/ea
14029DPMN	Mattress Only, 14029DP, 3 Tube System, 1/ea
14030XP	Replacement pump for 14030, 14048, 14054 & 14060, 1/ea
14029C	Cover, 1/ea
14029B	Mattress Carry Bag, 1/ea