



Platform Walker/Crutch Attachment



10105CB SOLD

- Designed for those who are unable to grip the walker
- Cradles forearm with soft vinyl padding for comfort
- Large top bracket allows the platform attachment to be placed over the hand grip of the walker
- Adjustable strap securely holds arm in place
- 10105CB (Crutch Attachment) sold separately

PRODUCTS

Description

10105-1	Walker Platform Attmnt, 1/ea
10105-2	alker Platform Attatcmt, 2/cs
10116	Universal Platform Walker/Crutch Attachment, 1/pr

SPECIFICATIONS

Description

Overall Width (10105 10116)	5.75" 6"
Overall Depth	18"
Weight (10105 10116)	2 7 lbs.
Carton Shipping Weight (1010501 10105-2 10116)	3 6 10 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

ACCESSORIES

Description

10116BRACKET	Brackets for 10116, 1/ea
10105CB	CrutchBracket 10105, 1/ea