


ESSENTIAL ACCESSORIES FOR PATIENT CARE AND WEANING

Shiley™ Tracheostomy Tube Accessories



Effective and easy to use, accessories include:

- Disposable inner cannula (DIC) with integral 15 mm snap-lock connector for use with DCT, DFEN, PERC, DCFS, and DCFN
- Spare inner cannula (SIC) with red, integral 15 mm twist-lock connector — for temporary use during cleaning or removal of Shiley™ reusable inner cannulae, for use with LPC, FEN, CFS, and CFN
- Fenestrated disposable inner cannula (DICFEN) for use with DFEN and DCFN
- White 15 mm cap (CAP) attaches to the inner cannula and used during weaning to help facilitate breathing through the upper airway; use with Shiley™ fenestrated inner cannulae, or Shiley™ flexible tracheostomy tubes
- Decannulation plug (DCP), closes the proximal end of the FEN and CFN tubes to permit breathing through fenestration and the upper airway
- Disposable decannulation plug (DDCP), closes the proximal end of DCFS, DCFN, and DFEN to permit breathing through fenestration and upper airway

DIC and DICFEN available in Jackson sizes 4, 6, 8, and 10. DDCP and CAP fit all sizes of Shiley™ disposable inner cannula tracheostomy tubes.

Latex free. Sold in boxes of 10. Packaged individually, sterile.

