

For older children and/or larger pouch sizes

SenSura® Mio Flat MIDI Drainable Pouch

- Standard wear
- EasiClose™ WIDE Outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window for easy placement
- **Not made with natural rubber latex**



Opaque 10¼" (26 cm) / 485 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
10461	Cut-to-fit	¾ - 1¾" (10-45 mm)	A4424

SenSura® Mio Flat MIDI Closed Pouch

- Standard wear
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window for easy placement
- **Not made with natural rubber latex**



Opaque 7¾" (19 cm) / 340 mL Units: 30

Code	Cut	Stoma Size	HCPCS*
10871	Cut-to-fit	⅝ - 1¾" (15-45 mm)	A4416

SenSura® Mio Flat MINI Closed Pouch

- Standard wear
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window for easy placement
- Comfortable and discreet stoma cap alternative
- **Not made with natural rubber latex**



Opaque 5¾" (13 cm) / 110 mL Units: 30

Code	Cut	Stoma Size	HCPCS*
10870	Cut-to-fit	⅝ - 1¾" (15-45 mm)	A4416

SenSura® Mio Flat MAXI Urostomy Pouch

- Standard wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for people in a wheelchair
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



Transparent 10½" (27 cm) / 400 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
10586	Cut-to-fit	¾ - 1¾" (10-45 mm)	A5071