



PHYSICIAN OFFICE



Comprehensive Practice Essentials

LET US HELP YOU FOCUS ON PATIENT CARE.

Low unit-of-measure quantities. Limited storage space. Declining reimbursements. The demands of running a busy physician practice never stop. And neither do we.

Our dedicated representatives understand your challenges first hand and have the solutions to help you focus more on your patients and less time running your practice.

- » **Manufacturer-direct benefits**, including guaranteed savings
- » **Specialists for lab testing** from point-of-care testing, CLIA-waived to complex testing for better visibility and patient outcomes
- » **Latest products and best practice ideas** to keep your office in top form
- » **Unmatched logistical expertise** and world-class service
- » **Help with integration** into a health network and support to optimize individual practice expertise

Find everything you need from general med/surg supplies and pharmaceuticals to lab equipment and much more.

Contact your dedicated Medline physician office representative • Call 1-855-294-9618 • Visit Medline.com/physician

DEDICATED FIELD SALES REPS

physician office specialists,
lab equipment &
supply consultants









EASY ORDERING

online, phone, fax




PRODUCTS YOU NEED

- » Point-of-care testing 
- » Diagnostic supplies 
- » Exam gloves 
- » Needles & syringes 
- » Capital equipment 
- » Injectables and vaccines 
- » Instruments 
- » Physical therapy 

NEXT-DAY & SAME-DAY DELIVERY





■ Point-of-Care Testing	5-22
■ Lab Supplies	23-26
■ Diagnostic Equipment.....	27-44
■ Exam Room Supplies	45-68
■ Gloves.....	69-72
■ Hand Hygiene.....	73-78
■ Instruments	79-92
■ Patient Care.....	93-118
■ Pharmaceuticals.....	119-138
■ Syringes & Blood Draw.....	139-152
■ Wound Care	153-172
■ Capital Equipment	173-198
■ Office Supplies.....	199-208

ORDER ONLINE AND GET WHAT YOU NEED — WHERE AND WHEN YOU NEED IT.

Online ordering offers you the convenience of 24/7 shopping. You can easily search for products, get the latest availability and pricing, save shopping lists, place orders, reviews invoices and credit memos, set up email reminders, and create and manage repeating orders.


Register now by following these six simple steps:

1 Open your web browser and go to www.medline.com/physician.


2 In the upper right hand corner of the home page, click “Login/Register”




3 A popup will appear, click on “Register”.



4 This will bring you to the Registration page. Under User Login Request, click on “Create User Login”. If you need your Medline account number, please contact us at 1-855-294-9618 or ask your Medline Physician Office Representative.



5 On the “User Login Request” page, complete the form and enter your Medline account number.



6 When you finish, click “Submit Request.” You will receive an email with your User ID and password within one business day.

Questions?

Contact our Help Desk at helpdesk@medline.com or your Medline Physician Office Representative at: www.medline.com/physician.

POINT-OF-CARE TESTING

Medline Point-of-Care Tests	6-7
Urinalysis	7-8
Pregnancy Tests	9-10
Fecal Occult Blood Tests	10
Strep A Tests/Influenza Tests	11-12
RSV/Mono Tests	14
H. Pylori and Women's Health Tests	14-15
HIV/Drug Tests	15-16
Cholesterol Tests	16
Lead Tests	16
Hemoglobin Tests	17
Hemoglobin A1c Tests	18
Blood Glucose Tests	18-19
PT/INR Systems	19
Point-of-Care Chemistry	20
Lancets	22
Monofilament	22



Medline point-of-care tests provide the accurate, rapid results necessary to evaluate and deliver value-based care and a better patient experience. Our point-of-care tests can also help improve staff efficiencies and enhance patient satisfaction, which may result in increased referrals and revenue.



MPH22025HCG



MPH11025HCG

hCg Pregnancy Tests

- » Highly accurate rapid test that detects human chorionic gonadotropin in multiple formats to fit the needs of your staff and facility
- » Results in three minutes using urine and five minutes using serum
- » Sensitivity: >99%
- » Urine Strips and Cassettes: 25 mIU/ml
- » Combo Cassettes: 20 mIU/mL (urine); 10 mIU/ml (serum)
- » Built in controls ensure accuracy
- » CLIA waived using urine; Moderately Complex using serum
- » Compatible external control: ALR92218

Item No.	Description	Pkg
MPH22025HCG	Combo+ Test Cassette	25/bx
MPH12025HCG	Test Cassette	25/bx
MPH11025HCG	Test Strip	25/bx



MPH30025SA



Strep A Dipstick Test

- » Rapid test that helps in early detection of Strep A pharyngitis using a throat swab
- » Results in just 5 minutes
- » Built-in procedural control
- » Direct detection of bacterial antigen not dependent on bacterial growth in culture
- » Sensitivity: 97%, Specificity: 95%
- » CLIA waived

Item No.	Description	Pkg
MPH30025SA	Strep A Test Strips	25/bx



MPH40020HP



H Pylori Cassette Test

- » For in vitro diagnostic use only
- » Results in 10 minutes
- » Rapid test for the qualitative detection of IgG antibodies to *Helicobacter pylori* (*H. pylori*) in whole blood, serum and plasma

Item No.	Description	Pkg
MPH40020HP	H. Pylori Rapid Test	20/bx



MPH50025M



Mono Test

- » Rapid detection for Infectious Mononucleosis (IM)
- » Results in 8 minutes
- » One-step antibody test for IM that uses direct solid-phase immunoassay technology to detect IM heterophile antibodies in human serum, plasma or whole blood
- » Relative sensitivity: 99.9%
- » Relative specificity: 98.8%
- » Moderate using serum or plasma
- » CLIA waived using whole blood

Item No.	Description	Pkg
MPH50025M	Mono Test Cassette	25/bx

Our dedicated PO lab specialists can help you with understanding the best testing options for your lab. For a free consultation, contact your Medline PO sales specialist.



iFOB Test Kits

MPH70025CTM

- » Rapid detection for GI disorders and colorectal cancer screening
- » Patient-friendly, and non-invasive, requires only one sample per test
- » Results in five minutes
- » No dietary or medication restrictions
- » Built-in control ensures accuracy
- » 50ng hHB/ml analytical sensitivity
- » 99.7% relative accuracy
- » CLIA waived

Item No.	Description	Pkg
MPH70025CTM	Test Kit: Test Cassettes, Collection Tubes and Mailers	25/bx
MPH70025CT	Test Kit: Test Cassettes and Collection Tubes	25/bx
MPH70025M	Mailers	25/bx
MPH70025T	Collection Tubes	25/bx



ACNUC071

120 Mini Urine Analyzer

- » Runs 60 test/hour in single test mode or 120 test/hour in continuous test mode
- » User and QC lockout
- » Automatic calibration for accurate results and easy operation
- » Saves and recalls the last 2,000 results
- » Reads Medline Urinalysis ALB/CRE, 4OB, 5OB, URO7 and 10SG parameter strips
- » Accessories available: optional barcode reader, hermal printer paper and sticky paper

Item No.	Description	Pkg
MPHUA120MP	120 Mini Urine Analyzer Promo (Includes Analyzer and 5 Bottles of 10SG Strips)	1/kt
ACNUC071	Control, Urine, Dip Tub, 2 x 12 ml	1/kt
ACNUC011	Control, Urine Liquid, Combo, 2 x 5 ml	1/kt



Urinalysis Reagent Strips

- » Multiple strip configurations to fit your testing needs
- » Compatible for visual and analyzer reading
- » Analytes appear in the same order as leading brands
- » CLIA waived

Item No.	Description	Pkg
MPHUA10SG	10SG Urinalysis Strips	100 ea/bt
MPHUAURO7	URO7 Urinalysis Strips	100 ea/bt
MPHUA5OB	5OB Urinalysis Strips	100 ea/bt
MPHUA4OB	4OB Urinalysis Strips	100 ea/bt
MPHUA2GP	2GP Urinalysis Strips (Visual Read Only)	100 ea/bt
MPHUAMAC	ALB/CRE Urinalysis Strips (Analyzer Read Only)	25 ea/bt



MPHUA10SG

Analytes

Item No.	Leukocytes	Nitrite	Urobilinogen	Protein	pH	Blood	Specific Gravity	Ketone	Bilirubin	Glucose	Albumin	Creatinine
MPHUA10SG	•	•	•	•	•	•	•	•	•	•		
MPHUAURO7	•	•		•	•	•		•		•		
MPHUA5OB	•	•		•		•				•		
MPHUA4OB	•	•		•						•		
MPHUA2GP				•						•		
MPHUAMAC											•	•



BAH1780

SIEMENS

CLINITEK® Status+ Analyzer Urinalysis

- » Offers a broad menu for routine point-of-care urinalysis testing on one analyzer
- » Most results are complete in 60 seconds
- » Performs Multistix® 10 SG and family of test strips, CLINITEK
- » Microalbumin, and Multistix PRO® detection of early kidney disease and CLINITEST hCG pregnancy tests
- » Printed results
- » Connectivity option to EMR/Data Management Systems
- » Patient and operator ID, operator management, and automatic calibration

Item No.	Description	Pkg
BAH1780	Analyzer	1/ea
BAH1759	CLINITEK Label Paper	5/pk



AMV2161Z



BAH1760



AMV2184Z

SIEMENS

CLINITEK Reagent Strips and Accessories

- » CLIA waived

Item No.	Description	Pkg
AMV2161	Multistix 10SG Reagent Strips	100 ea/btl, 24 btl/cs
AMV2161Z	Multistix 10SG Reagent Strips	100 ea/btl
BAH2083Z	CLINITEK Reagent Strips	25 ea/btl
BAH1554Z	Multistix PRO 10 LS Reagent Strips	100 ea/btl
BAH1760	CLINITEST hCG Pregnancy Test	25 ea/bx
AMV2166Z	Uristix® 4 Reagent Strips	100 ea/btl
AMV2184Z	Uristix Reagent Strips	100 ea/btl



SIEMENS

Chek-Stix® Combo Pack Strips

- » Combination kit (positive and negative control strips) for urinalysis testing
- » Positive, negative, or defined results for bilirubin, blood, glucose, ketones, leukocytes, nitrites, pH, protein, specific gravity, and urobilinogen
- » Control strips are reconstituted with distilled water and remain stable for eight hours
- » CLIA waived

Item No.	Description	Pkg
AMV1364Z	Chek-Stix Combo Pack	50 ea/pk



Dipper Urinalysis Dipstick Control

- » Simulates patient testing
- » Liquid, made with human urine
- » 18-month refrigerated shelf life from date of manufacture
- » 3 months (or 20 dips) open vial stability when stored at 2 - 8°C
- » Designed for use with most urinalysis reagent strips



Item No.	Description	Pkg
QMX144001	Urine Dipstick Control, 6 x 15 ml	1 pk/kt



Urisys 1100 Urine Analyzer & Accessories

- » Fast and automatic
- » Easy-to-read color changes
- » Prints patient report
- » Can be stored for one year at room temperature after vial has been opened
- » CLIA waived



BID0326076316

Item No.	Description	Pkg
BID260763602	Starter Kit; Meter and 6 Vials of Chemstrip 10	1/ea
BID3617556001	Urisys 1100® Urine Meter	1/ea
BID1008552	Chemstrip 7	600/cs
BID1895427160	Chemstrip 9	6/cs
BID1895362160	Chemstrip 10	600/cs
BID0326076316	Chemstrip 10MD	100/bx



QDL20109



QDL20110A



QuickVue® One-Step hCG Urine Test

- » Results in three minutes for urine, five minutes for serum in combo cassette
- » Greater than 99% sensitivity, specificity, and accuracy
- » 25 mIU/ml sensitivity in serum and urine
- » CLIA waived for urine

Item No.	Description	Pkg
QDL20109	QuickVue hCG Urine Test Cassette	25/bx
QDL0272	Urine hCG Control Set	1 kt
QDL20110A	QuickVue One-Step Pregnancy Test Cassettes	50/bx
QDL0272	Urine hCG Control Set	1 kt
QDL281	Serum hCG Control Set	1 kt



ALR4581235020



hCG Pregnancy Test

- » Suitable for laboratory and point-of-care testing
- » Delivers fast, accurate results in as little as three minutes
- » Greater than 99% sensitivity
- » hCG detection:
 - Urine strip and cassette tests detect 25 mIU/ml
 - Combo+ test detects 20 mIU/ml for urine and 10 mIU/ml for serum
- » CLIA complexity:
 - Waived using urine
 - Moderate complexity using serum

Item No.	Description	Pkg
ALR4581235020	hCG Dipstick (92211)	50/bx
ALR4581255020	hCG Cassette (92217)	40/bx
ALR4581245020	hCG Combo Cassette (92215)	40/bx



MPH22025HCG

MPH11025HCG



Pregnancy Tests

- » Multiple formats to fit the needs of your staff and facility
- » Results in three minutes using urine and five minutes using serum
- » Sensitivity: >99%
- » Urine strips and cassettes: 25 mIU/ml
- » Combo Cassettes: 20 mIU/ml (urine); 10 mIU/ml (serum)
- » Built in controls ensure accuracy
- » CLIA waived using urine; Moderately complex using serum
- » Compatible external control: ALR92218

Item No.	Description	Pkg
MPH22025HCG	Combo+ Test Cassette	25/bx
MPH12025HCG	Test Cassette	25/bx
MPH11025HCG	Test Strip	25/bx



OSOM101

OSOM102



OSOM® hCG Pregnancy Tests

- » Urine Dipsticks
 - Results in five minutes or less
 - 25 mIU/mL sensitivity
- » Urine Card
 - Results in three minutes
 - 20mIU/ml sensitivity
- » CLIA waived

Item No.	Description	Pkg
OSOM101	OSOM hCG Urine Dipstick Test	50/bx
OSOM102	OSOM Card Pregnancy Test	25/bx
OSOM134	OSOM hCG Urine Control Set	3/bx



MPH70025CTM

iFOB Test

- » Rapid fecal occult blood detection for GI disorders and colorectal cancer screening. Patient-friendly, non-invasive test, requires only one sample per test
- » Results in five minutes
- » No dietary or medication restrictions
- » Built-in control ensures accuracy
- » Specific to human hemoglobin
- » 50ng hHB/ml analytical sensitivity
- » 99.7% relative accuracy
- » CLIA waived

Item No.	Description	Pkg
MPH70025CTM	Test Kit: Test Cassettes, Collection Tubes and Mailers	25/bx
MPH70025CT	Test Kit: Test Cassettes and Collection Tubes	25/bx
MPH70025M	Mailers	25/bx
MPH70025T	iFOB Collection Tubes	25/bx



HLS5072



ColoScreen Fecal Occult Blood Test Kits

- » Test for occult blood in fecal matter in batches or as single specimen.
- » 100 specimen applicators
- » 100 Single ColoScreen Slides with On-Slide
- » CLIA waived

Item No.	Description	Pkg
HLS5072	Single Slides, Developer, and Applicators	100/bx
HLS5073	Take-Home Kits Including Return Envelope, Triple Slide, Applicators, and Developer	80/bx



HMET1CK30



Fecal Occult Blood Test

- » Rapid fecal occult blood test (FIT)
- » Results available in five minutes
- » No dietary or medication restrictions
- » Specific to human hemoglobin
- » 50ng hHB/ml analytical sensitivity
- » CLIA waived

Item No.	Description	Pkg
HMET1CK30T	30 Cassettes, 30 Tubes	30/bx
HMET1CK30	30 Cassettes, 50 Tubes, 50 Mailers	30/bx
HMET1CK50	50 Cassettes, 50 Tubes	30/bx
HMET1CM50	Patient Specimen Collection and Return Kits*	30/bx
HMET1CT30	30 Test Cassettes	30/bx
HMET1TB50	50 Sample Collection Tubes Only	30/bx
HMET1TC01	Negative and Positive Controls, 5 ml/ea One Negative, One Positive Control Per Box	1 kt

*Each kit consists of: One patient instruction card, two collection papers, one specimen pouch, and one return mailer.



C64151AZ



Hemocult® Fecal Occult Blood Test

- » Rapid FOBT for colorectal cancer screening
- » Results available in 60 seconds
- » Built-in performance monitors
- » 14-day on-slide sample stability
- » Comprehensive patient educational support
- » CLIA waived

Item No.	Description	Pkg
B-C60151AZ	Hemocult Single Slides Kit	100/bx
B-C61130AZ	Hemocult II DispensaPak Patient Screening Kit	40/bx
B-C64151AZ	Hemocult SENA Single Slides Kit	40/bx
B-C64130AZ	Hemocult II SENA Dispensapak Plus Patient Screening Kit	100/bx



MPH30025SA



Strep A Dipstick Test

- » Simple and rapid test that helps in early detection of Strep A pharyngitis
- » Accurate and reliable results in just five minutes
- » Built-in procedural control
- » Direct detection of bacterial antigen not dependent on bacterial growth in culture
- » Sensitivity: 97%, Specificity: 95%
- » CLIA waived

Item No.	Description	Pkg
MPH30025SA	Strep A Test Strips	25/bx



OTC92001A



Acceava® Strep A Dipstick

- » Easy-to-read rapid test helps with early detection of Strep A pharyngitis. Accurate and reliable results in just five minutes
- » Built-in procedural control
- » Sensitivity: 97%, Specificity: 95%
- » CLIA waived

Item No.	Description	Pkg
OTC92001A	Acceava Strep A Dipsticks	50/bx



OSOM147

OSOM

OSOM Ultra Strep A Test

- » Two color immunochromatographic assay test to detect group A streptococcus antigen
- » Results in five minutes
- » Easy-to-read two color results
- » Specimen Type: Throat swab
- » Controls included in kit
- » CLIA waived

Item No.	Description	Pkg
OSOM147	OSOM Ultra Strep A Dipstick Test	25/bx
OSOM149	OSOM Ultra Strep A Dipstick Test	50/bx



QDL20108



QuickVue Strep A Tests

- » Rapid Strep A test to detect group a streptococcal antigen directly from patient throat swab specimen
- » Results in five minutes
- » CLIA waived

Item No.	Description	Pkg
QDL0343	QuickVue In-Line Strep A Test	25/bx
QDL0345	Strep A Liquid Control Set	1/bx
QDL20108	QuickVue Dipstick Strep A Test	50/bx
QDL20227	Strep A Swabs	75/bx
QDL0354	Strep A Liquid Control Set	1/bx

Antigens can be taken directly from throat swabs and beta-hemolytic colonies on blood agar plates.



PBOBSP510



BioSign Rapid Flu A + B Antigen Panel Test

- » Provides quick diagnosis of influenza A and B
- » Easy-to-use cassette format and rapid results allow for point-of-care testing
- » Uses nasal swab, nasopharyngeal swab and nasal aspirate/wash specimens
- » Positive results are available in 10 minutes
- » CLIA waived for nasal swab and nasopharyngeal swab samples

Item No.	Description	Pkg
PBOBSP510	BioSign A + B Influenza Test	25/bx



QDL20183



QuickVue Influenza A+B Test

- » Detects and differentiates influenza type A and type B antigens directly from nasal swab, nasopharyngeal swab, nasal wash and/or nasal aspirate specimens
- » Results in 10 minutes
- » Sensitivity: Flu A 94%; Flu B 70%
- » Specificity: Flu A 90%; Flu B 97%
- » CLIA waived

Item No.	Description	Pkg
QDL20183	QuickVue Influenza A+B Test	25/bx



QDL20299

Sofia 2

- » Results in three minutes
- » Advanced result technology
- » Wireless connectivity
- » CLIA waived

Item No.	Description	Pkg
QDL20299	Sofia 2 Reader	1/ea
QDL20218	Influenza A & B Test	25/bx
QDL20274	Strep A+ Test	25/bx
QDL20260	RSV Test	25/bx



BID7341920190



Cobas Liat

- » Results in 20 minutes
- » Single assay runs
- » Closed system with multiple controls
- » CLIA waived

Item No.	Description	Pkg
BID7341920190	Cobas Liat Analyzer	1/ea
BID7341890190	Influenza A & B Test	20/bx
BID7341911190	Strep A Test	20/bx
BID8160104190	Influenza A & B + RSV Test	20/bx



BDD256055



BDD256045

Veritor® System

- » Incorporates the unique technologies that provide an objective and lab-quality test result for Flu A + B, RSV, and Strep A
- » Testing requires minimal hands-on time
- » Unique Adaptive Read Technology helps reduce false-positive results
- » Easy-to-read digital display
- » CLIA waived

Item No.	Description	Pkg
BDD256055	BD Veritor System Reader	1/ea
BDD256045	Influenza A & B Test	30/bx
BDD256038	RSV Test	30/bx
BDD256040	Strep A Test	30/bx



ALRNAT024

Alere™ i

- » Alere i is a rapid, instrument-based, isothermal system for the qualitative detection of infectious diseases
- » Provides molecular results in just minutes to make decisions faster
- » Intuitive and user friendly system
- » Lab accurate results
- » CLIA waived for influenza A + B and Strep A

Item No.	Description	Pkg
ALRNAT024	System	1/ea
ALR427000	Influenza A & B Test	24/bx
ALRI734000	Strep A Test	24/bx



Design. Deliver. Install.

We do it all. World-class furnishings and capital solutions for your physician office. We can be your partner to help provide quality products and navigate logistical obstacles to optimize patient care.

- » Premium products
- » Planning and strategy
- » Support and design
- » Fast quotes
- » Delivery and install services
- » White glove service options

Contact your Medline representative to get started.
1-855-294-9618 | www.medlinecapitalquote.com



OTC430122

BinaxNOW RSV Test

- » Intended for in vitro diagnosis of neonatal and pediatric patients under the age of five
- » Results available in 15 minutes
- » Critical care instructions available to patient
- » Sensitivity 92%, Specificity 92%
- » CLIA waived

Item No.	Description	Pkg
OTC430122	BinaxNOW RSV Test	22/bx



QDL20193

QuickVue® RSV Test

- » Detects respiratory syncytial virus (RSV) antigen (viral fusion protein) directly from nasopharyngeal swab, nasopharyngeal aspirate, or nasal/nasopharyngeal wash specimens for symptomatic pediatric patients (eighteen years of age and younger)
- » A single sample can be used to run both the QuickVue RSV Test and QuickVue brand Influenza tests

Item No.	Description	Pkg
QDL20193	QuickVue RSV Test	20/bx



MPH50025M



Mono Test

- » Rapid detection for Infectious Mononucleosis (IM)
- » One-step antibody test for IM that uses direct solid-phase immunoassay technology to detect IM heterophile antibodies in human serum, plasma or whole blood
- » Results in eight minutes
- » Relative sensitivity: 99.9%
- » Relative specificity: 98.8%
- » Moderate using serum or plasma
- » CLIA waived using whole blood

Item No.	Description	Pkg
MPH50025M	Mono Test Cassette	25/bx



MPH40020HP

H. Pylori Cassette Test

- » For in vitro diagnostic use only
- » Results in 10 minutes
- » Rapid test for the qualitative detection of IgG antibodies to *Helicobacter pylori* (*H. pylori*) in whole blood, serum and plasma

Item No.	Description	Pkg
MPH40020HP	<i>H. Pylori</i> Rapid Test	20/bx



OSOM145

OSOM Mono Test

- » OSOM® Rapid Mononucleosis (Mono) test in serum, plasma or whole blood
- » Result in five minutes
- » Easy-to-read two-color results
- » Includes two extra test sticks in each kit
- » Sensitivity: 100%
- » Specificity: 99%
- » CLIA waived (whole blood)

Item No.	Description	Pkg
OSOM145	OSOM Mono Test	25/bx



ALR4581505020



Clearview® H. Pylori Test

- » *H. pylori* rapid chromatographic immunoassay test for use with whole blood, serum or whole blood
- » Results in 10 minutes
- » Sensitivity: 93%
- » Specificity: 89.2%
- » CLIA waived (whole blood)

Item No.	Description	Pkg
ALR4581505020	Clearview <i>H. Pylori</i> Test (92405)	30/bx



OSOM181

SEKISUI

OSOM® Trichomonas Test

- » Results in 10 minutes
- » Sensitivity: 83%
- » Culture: 95%
- » CLIA waived

Item No.	Description	Pkg
OSOM181	Trichomonas Rapid Test	25/bx



OSOM183

SEKISUI

OSOM BVBLUE® Test

- » Detects elevated vaginal fluid sialidase activity, an enzyme produced by bacterial pathogens associated with Bacterial Vaginosis
- » Results in 10 minutes
- » Easy-to-read results
- » Sensitivity: 92.8%
- » Specificity: 98%
- » Detects vaginitis caused by Gardnerella, Bacteroides, Prevotella and Mobiluncus
- » CLIA waived

Item No.	Description	Pkg
OSOM183	OSOM BVBLUE Test	25/bx



Toxicology Drug Screening Devices

- » Fast results
- » CLIA waived

Item No.	Description	Pkg
DRUG10DIPDXC	iScreen® 10-Panel Dip Card (AMP, BAR, BUP, BZO, COC, THC, MAMP, MTD, MOP, OXY)	25/bx
DRUG1127023	iCup® Dx 12-Panel All Inclusive Cup with Adulteration (AMP, BAR, BZO, COC, MAMP, MDMA, MTD, OPI, OXY, PCP, TCA, THC (CR, NI, PH, BL, SG))	25/bx
DRUG1147022	iCup Dx 14-Panel All Inclusive Cup (AMP, BAR, BUP, BZO, COC, MAMP, MDMA, MTD, OPI, OXY, PCP, PPX, THC, TCA)	25 ea/bx
DRUG8CUPDX	iCup Dx 8-Panel All Inclusive Cup with Adulteration (AMP, BZO, COC, MAMP, MTD, OPI, OXY, THC (CR, NI, PH, BL, SG))	25/bx



DRUG8CUPDX



ALR7D2648



Determine™ HIV-1/2 Ag/Ab Combo

- » Rapid test that detects both HIV-1/2 antibodies and free HIV-1 p24 antigen
- » Results in 20 minutes
- » Multiple sample types; fingerstick and venous whole blood, serum or plasma
- » CLIA waived (whole blood)

Item No.	Description	Pkg
ALR7D2648	Alere Determine HIV-1/2Ag/Ab Combo	25/bx



HMESCCUP13124



Drug Screening Devices

- » Results in five minutes
- » CLIA waived

Item No.	Description	Pkg
HMESCCUP13124	12 Panel Cup (AMP, BARB, BUP, BZO, COC, mAMP, MDMA, MOR, MTD, OXY, PCP, THC)	25/bx
HMEFCUP4144A2	14 Panel Cup with Adulteration (AMP, BARB, BUP, BZO, COC, EDDP, mAMP, MDMA, MTD, OPI, OXY, PCP, TCA, THC, PH-SG-CR)	25/bx
HMESCDIP12104	10 Panel Dip Card (AMP, BARB, BUP, BZO, COC, mAMP, MTD, OPI, OXY, THC)	25/bx



CHLS10989



CHLS10959



CHLS88773



CHLS11010



CHLS52193



CHLS10311



Cholestech LDX System and Accessories

- » Testing for cholesterol and related lipids, glucose and glycemic control, and liver enzymes
- » Results available in five minutes
- » Fingerstick sampling and 40 µl sample
- » CLIA waived

Item No.	Description	Pkg
CHLS10959	Cholestech LDX System	1/ea
CHLS13454	Cholestech LDX Analyzer Only	1/ea
CHLS10989	Lipid Profile	10/bx
CHLS10991	Lipid Profile & Glucose	10/bx
CHLS10986	Total Cholesterol	10/bx
CHLS10987	Total Cholesterol and HDL	10/bx
CHLS10988	Total Cholesterol and Glucose	10/bx
CHLS10990	Total Cholesterol, HDL, Glucose	10/bx
CHLS88773	Multi-Analyte Controls	1 kt

Item No.	Description	Pkg
Accessories		
CHLS26333	Thermal Printer Labels	1/r1, 400/r1
CHLS10228	Optics Cassette with Case	1/ea
CHLS11010	Pipette Tips for use with 35µL, 40µL, and 50µL MiniPet Pipettes	50/bx
CHLS11846	35µL MiniPet Pipette	1/ea
CHLS10311	Capillary Plungers, to be used with Cholestech Capillary Tube	50/bx
CHLS52193	40µL Capillary Tubes (Lipids/Glucose)	50/bx



CardioChek® Plus Family of Products

- » CardioChek Plus analyzer provides rapid on-site lipid panel screening
- » Results in 90 seconds
- » Portable, handheld device
- » No refrigeration of strips required
- » Wireless communication
- » CLIA waived



PTH2700

Item No.	Description	Pkg
PTH2700	CardioChek Plus Analyzer	1/ea
PTH2729	CardioChek® Plus Smart Bundle™ 2-Pack: 1 Vial of Lipid + 1 vial of eGLU® Test Strips	15/bx
PTH1710	Lipid Panel Test Strips	15/bx
PTH2713	eGLU Test Strips – Electrochemical Glucose	50/bx
PTH2750	CardioChek Printer	1/ea



M-B706760



LeadCare® II Lead Analyzer

- » A rapid test for lead blood testing
- » Results in three minutes
- » CLIA waived

Item No.	Description	Pkg
M-B706760	LeadCare II Analyzer Kit Includes: Portable Blood-Lead Testing instrument, AC Power Cord, AA Batteries, User's Guide (English Quick Reference Guide, Instructional Video, CDC Blood-Collection Video)	1/ea
M-B706762	LeadCare II Test Kit Includes: 48 Sensors, 48 Tubes of Treatment Reagent, Capillary Tubes and Plungers, Transfer Droppers, 2 Levels of Quality Control, Calibration Button	48/kt
M-B703447	LeadCare II Printer	1/ea
M-B703443	LeadCare II Printer Labels	2/bx

HMQ121721



Hb 201+ Systems

- » Provides rapid and reliable hematology tests
- » Results in less than one minute
- » EMR connectivity
- » Built in self test
- » CLIA waived

Item No.	Description	Pkg
HMQ121721	System	1/ea
HMQ111716	Microcuvettes (Vial)	200/bx
HMQ111715	Microcuvettes (Individual)	100/bx



OG3000CONTC

HemoPoint® H2 Meter

- » Use to manage anemia and red blood cell disorders with this hemoglobin A hematocrit meter
- » Results in 30 seconds
- » Easy-to-read touchscreen
- » 4,000 test recall
- » Small sample size: 8 µl
- » CLIA waived

Item No.	Description	Pkg
S-OG3000CONTC	Meter	1/ea
S-03015100	Microcuvettes	100/bx
S-03015200	Microcuvettes	200/bx
S-03025050	Microcuvettes (Individual)	50/bx



EVENCARE® G3 Professional

- » Results in six seconds
- » Test strip ejector
- » No coding required
- » Small blood sample size: 0.8 µl
- » Stores 300 test times



MPH3560Z

Item No.	Description	Pkg
MPH3560Z	Control Solution High/Low	1/bx
RAYCR2032	Replacement Battery, 3V	10/bx
MPH3540	Glucose Meter	1/ea
MPH3550	Test Strips	50/bx, 12 bx/cs
MPH3550Z	Test Strips	50/bx
MPH3550FL	Foil Wrapped Test Strips	50/bx, 12 bx/cs
MPH3550FLZ	Foil Wrapped Test Strips	50/bx
MPH3560	Control Solution High/Low	1/bx, 6 bx/cs



ALR1116554

Afinion™ HbA1c

- » Designed for quantitative determination of glycated hemoglobin (HbA1c) in human whole blood used to monitor metabolic control in patients with diabetes
- » Results in three minutes
- » Small blood sample: 1 µl
- » CLIA waived

Item No.	Description	Pkg
ALR1116554	Afinion System	1/ea
ALR1115015	Afinion HbA1c	15/test
ALR1115178	Afinion HbA1c Control	1 kt



BAH5075US

DCA Vantage® Analyzer

- » Point-of-care solution for diabetes management
- » Results in six minutes
- » Small blood sample: 1 µl
- » CLIA waived

Item No.	Description	Pkg
BAH5075US	DCA Vantage Analyzer	1/ea
AMV5035	HbA1c Reagent Test Kits	4 bx/cs, 10/bx
AMV5035Z	HbA1c Reagent Test Kits	10/bx
AMV5068AZ	Control Solution Set for DCA 2000 Plus and Vantage	1/bx

PTH3021



A1CNow+® System

- » Designed to provide a fast, easy and accurate way to obtain HbA1c results to support timely decisions at the point-of-care
- » Results in five minutes
- » Small (5µL) blood sample
- » Portable—battery-operated allows use in multiple exam rooms
- » No refrigeration required (four month room temperature)

Item No.	Description	Pkg
PTH3021	A1CNow+ 20-Test Kit	20/bx
PTH3024	A1CNow+ 10-Test Kit	10/bx

HMQ120706



Glucose 201 System

- » Enables the highest accuracy in just three simple steps
- » Reduces serious contamination risks
- » CLIA waived

Item No.	Description	Pkg
HMQ120706	Glucose 201 System	1/ea
HMQ110706	Glucose 201 Microcuvettes (Vial)	100/bx
HMQ110723	Glucose 201 Microcuvettes (Individual)	50/bx

MPH6550



MPH6540



Harmony™ Blood Glucose Monitoring System

- » Utilizes the newest technology in blood glucose testing in order to deliver more accurate results
- » Test strips are designed to meet greater accuracy standards and limit interference with therapies such as oxygen
- » Wider hematocrit range allows for testing a wider patient population
- » Test strip ejector: Eliminates contact with used strip
- » Large backlit display: can be clearly read in low light

Item No.	Description	Pkg
MPH6540	Glucose Meter	1/ea
MPH6550	Test Strips	50/bx, 12 bx/cs
MPH6550Z	Test Strips	50/bx
MPH6560	Control Solution High/Low	1/bx, 6 bx/cs
MPH6560Z	Control Solution High/Low	1/bx
RAYCR2032	Replacement Battery, 3V	10/bx

MPH1540



MPH1550



EVENCARE® G2® Blood Glucose Monitoring System

- » Advanced technology makes diabetic care easier
- » Fast test time: Results in six seconds
- » Easy-to-use: Voice assistance guides you through the testing process and provides results in either English or Spanish
- » No coding required: Minimizes errors and inaccurate results
- » Small blood sample size: 1 µL
- » Alternate site testing: Option to test on the forearm, palm or fingertips
- » Memory option: Stores 300 test times and can provide 7-, 14- and 30-day averaging
- » Alarm option: Set reminders for up to three different times

Item No.	Description	Pkg
MPH1540	Glucose Meter	1/ea
MPH1550	Test Strips	50/bx, 12 bx/cs
MPH1550Z	Test Strips	50/bx
MPH1560	Control Solution High/Low	1/bx, 6 bx/cs
MPH1560Z	Control Solution High/Low	1/bx
MPHBAAZ	Replacement Battery, AAA	24/bx



EVENCARE® G3 Professional Blood Glucose Monitoring System

- » Quality control lock mode: Ensures quality control testing is performed
- » Fast test time: Results in six seconds
- » Test strip ejector: Helps reduce cross contamination by eliminating contact with used strip
- » No coding required: Minimizes errors and inaccurate results
- » Small blood sample size: 0.8µl
- » Alternate site testing: Option to test on the forearm, palm or fingertips

MPH3540



MPH3550



Item No.	Description	Pkg.
MPH3540	Glucose Meter	1/ea
MPH3550	Test Strips	50/bx, 12 bx/cs
MPH3550Z	Test Strips	50/bx
MPH3550FL	Foil Wrapped Test Strips	50/bx, 12 bx/cs
MPH3550FLZ	Foil Wrapped Test Strips	50/bx
MPH3560	Control Solution High/Low	1/bx, 6 bx/cs
MPH3560Z	Control Solution High/Low	1/bx
RAYCR2032	Replacement Battery, 3V	10/bx



CoaguChek® XS System

- » Uses exclusive technology to help ensure accurate results, and offers the flexibility and control that healthcare professionals need to confidently manage a wide range of anticoagulated patients
- » Accurate, precise results in about one minute
- » Neutralizes therapeutic levels of heparin and LMWH
- » INR corrected for hematocrit within specified range—manage a broad array of patients with one device
- » Onboard controls—no external QC necessary
- » Easy and convenient to use
- » Small sample size: 8µl for easy dosing
- » 300 test memory
- » CLIA waived



BID4837975001

Item No.	Description	Pkg
Meters		
BID4837975001	Meter	1/bx
Accessories		
BID4625315160	XS Strips	48/bx
BID4348150001	XS Lancets	200/bx
BID1621173001	Capillary Tubes	100/bx



PT/INR System

- » Directly detects clotting endpoint providing greater confidence in patient results
- » Results in one minute or less
- » Small sample size: 10µl
- » 100 test memory
- » CLIA waived



ABO03P6001

Item No.	Description	Pkg
ABO03P6001	PT/INR Professional System Includes: 1 Meter 1 Carrying Case 1 Professional User Manual 1 Professional Quick Reference Guide Auto Single-Usage Lancets – Bag of 12 4 AA Alkaline Batteries	1/ea
ABO03PTESTKIT	Coag-Sense Test Kit Includes: Coag-Sense Test Strips (50/bx) Control Test Strips (4/bx) 1 Vial Control Solution Transfer Tubes (54/bx)	50/bx

medline.com/physician is your source for product information, including instructions and videos.



AB004J6020



i-STAT System

- » Performs a full menu of the most commonly performed diagnostic tests
- » Results in two minutes
- » Handheld, portable device
- » Small sample size of 2 or 3 drops of blood
- » Flexible connectivity options

Item No.	Description	Pkg
AB004J6020	System Waived Kit	1/ea
AB009P3125*	Test Cartridge Chem 8+	25/bx
AB003P8425*	Test Cartridge Creatinine	25/bx
AB003P8325*	Test Cartridge Glucose	25/bx
AB003P8225*	Test Cartridge E3+	25/bx
AB003P8125*	Test Cartridge EC4+	25/bx
AB003P8025*	Test Cartridge 6+	25/bx

*Indicates CLIA waived cartridge.
For more cartridge options, please contact your Medline representative.



AB007P0501



Piccolo Xpress Chemistry Monitor

- » Provides a wide range of clinical chemistry tests
- » Results in 12 minutes
- » Portable instrument
- » Requires 100ul of blood, plasma or serum
- » Stores up to 5,000 patient records
- » CLIA waived

Item No.	Description	Pkg
AB007P0501	Piccolo Xpress Chemistry Monitor	1/ea
AB007P0208*	Piccolo Reagent Disc Comprehensive Metabolic Panel	10/bx
AB007P0204*	Piccolo Reagent Disc Basic Metabolic Panel	10/bx
AB007P0205*	Piccolo Reagent Disc Lipid Panel	10/bx
AB007P0210*	Piccolo Reagent Disc Liver Panel Plus	10/bx
AB007P0201*	Piccolo Reagent Disc General Chem 6	10/bx
AB007P0207*	Piccolo Reagent Disc Renal Panel	10/bx

*Indicates CLIA waived disc.
For more cartridge options, please contact your Medline representative.

Don't forget about infection prevention with point-of-care testing

- » Hand sanitizers on pages 74-75
- » Gloves on pages 69-72
- » Disinfectants on pages 62-67



 MEDLINE

ENHANCE PATIENT CARE, GROW YOUR BUSINESS.

With Your Own Lab, Receive Immediate Results and Increase Revenue.

Whether you're building a new lab, expanding your current practice or updating your equipment, our team of laboratory specialists can help.

Find the quality brands of lab diagnostic equipment you trust to improve patient satisfaction and your bottom line.



Ask about our complimentary laboratory analysis and flexible financing options

Contact your Medline representative · Call 1-855-294-9618 · Visit Medline.com/physician



Heel Incision Devices

- » High-quality steel blades make small, precise incisions for reduced bruising
- » Disposable, one-time use devices
- » Easy trigger activates blade quickly, safely and precisely
- » Color-coded for easy identification of devices

Item No.	Description	Pkg
MPHMPREEMIE2	Micro-Preemie Heel Lancet, Light Blue	50/bx, 4 bx/cs
MPHMPREEMIE2Z	Micro-Preemie Heel Lancet, Light Blue	50/bx
MPHPREEMIE2	Preemie Heel Lancet, White	50/bx, 4 bx/cs
MPHPREEMIE2Z	Preemie Heel Lancet, White	50/bx
MPHNEWBORN2	Newborn Heel Lancet, Pink	50/bx, 4 bx/cs
MPHNEWBORN2Z	Newborn Heel Lancet, Pink	50/bx
MPHTODDLER2	Toddler Heel Lancet, Purple	50/bx, 4 bx/cs
MPHTODDLER2Z	Toddler Heel Lancet, Purple	50/bx



MPHADJDEPTH



Adjustable-Depth Safety Lancets

- » Innovative rotating dial allows for three easy depth adjustments: 1.3mm, 1.8mm, 2.3mm
- » Pressure activated
- » Compares to the Accu-Chek Safe-T-Pro Plus

Item No.	Description	Pkg
MPHADJDEPTH	Adjustable-Depth Safety Lancet	1,000/cs
MPHADJDEPTHZ	Adjustable-Depth Safety Lancet	200/bx



MSC0002

Monofilament

- » Reliable, cost effective method of testing for the presence or absence of sensation
- » Helps identify patients at risk of developing neuropathy or diabetic foot ulcers
- » 5.07/10 g monofilament

Item No.	Description	Pkg
MSC0002	Foot Monofilament, Diabetic Neuropathy Test	25/bx



MPHSAFETY28



MPHST28



MPHSTSIDE28



MPHPRESS28

Safety Lancets

- » Simple and safe for healthcare professionals and patients
- » Needle retracts automatically and cannot be re-used
- » Different sizes allow healthcare professionals to obtain the proper amount of blood for capillary testing
- » Double spring action and controlled needle depth provides greater patient comfort
- » Quick and easy-to-use

Item No.	Description	Pkg
----------	-------------	-----

Push Button

MPHSAFETY28	28 Gauge, 1.6 mm	200/bx, 18 bx/cs
MPHSAFETY28Z	28 Gauge, 1.6 mm	200/bx
MPHSAFETY285	28 Gauge, 1.6 mm	50/bx, 24 bx/cs
MPHSAFETY285Z	28 Gauge, 1.6 mm	50/bx
MPHSAFETY281	28 Gauge, 1.6 mm	100/bx, 20 bx/cs
MPHSAFETY281Z	28 Gauge, 1.6 mm	100/bx
MPHSAFETY23	23 Gauge, 1.8 mm	200/bx, 10 bx/cs
MPHSAFETY23Z	23 Gauge, 1.8 mm	200/bx
MPHSAFETY21	21 Gauge, 1.8 mm	18 bx/cs
MPHSAFETY21Z	21 Gauge, 1.8 mm	200/bx
MPHSAFETYBL	Blade, 1.6 mm	200/bx, 18 bx/cs
MPHSAFETYBLZ	Blade, 1.6 mm	200/bx

Pressure Activated

MPHST28	28 Gauge, 1.8 mm	10 bx/cs
MPHST28Z	28 Gauge, 1.8 mm	200/bx
MPHST21	21 Gauge, 2.2 mm	10 bx/cs
MPHST21Z	21 Gauge, 2.2 mm	200/bx

Side Button

MPHSTSIDE28	28 Gauge, 1.8 mm	200/bx, 10 bx/cs
MPHSTSIDE28Z	28 Gauge, 1.8 mm	200/bx
MPHSTSIDE21	21 Gauge, 2.4 mm	200/bx, 10 bx/cs
MPHSTSIDE21Z	21 Gauge, 2.4 mm	200/bx

Contact Activated

MPHPRESS30	30 Gauge, 1.5 mm	100/bx, 20 bx/cs
MPHPRESS30Z	30 Gauge, 1.5 mm	100/bx
MPHPRESS28	28 Gauge, 1.8 mm	100/bx, 20 bx/cs
MPHPRESS28Z	28 Gauge, 1.8 mm	100/bx
MPHPRESS26	26 Gauge, 1.8 mm	100/bx, 20 bx/cs
MPHPRESS26Z	26 Gauge, 1.8 mm	100/bx
MPHPRESS23	23 Gauge, 2.2 mm	100/bx, 20 bx/cs
MPHPRESS23Z	23 Gauge, 2.2 mm	100/bx
MPHPRESS21	21 Gauge, 2.2 mm	100/bx, 20 bx/cs
MPHPRESS21Z	21 Gauge, 2.2 mm	100/bx



LAB SUPPLIES

Bags/Containers.....	24
Specimen Containers.....	24-25
Transfer Pipettes.....	25
Culture Tubes.....	25
Petri Dishes.....	25
Microscope Slides.....	26
Lab Accessories.....	26
Cover Glass.....	26
Chemicals.....	26



Biohazard Bags

- » Unique two-pouch design isolates patient specimen and documents
- » Zip closure for containment
- » Biohazard symbol printed on bag

Item No.	Description	Pkg
DYND30260	6 x 6" Bag w/Document Pouch	1,000/cs
DYND30261	6 x 9" Bag w/Document Pouch	1,000/cs
DYND30271	8 x 8" Bag w/Document Pouch	1,000/cs
DYND30281	12 x 15" Bag w/Document Pouch	500/cs



DYND69STAT

Red STAT Biohazard Bag

- » Zip closure
- » Document pouch
- » Biohazard logo printed on bag
- » 2 mil thickness

Item No.	Description	Pkg
DYND69STAT	6 x 9" Red, STAT Biohazard Bag	1,000 ea/cs



DYNDPVE912

Patient Belonging Valuables Bag

- » Tamper evident adhesive closure
- » Writable chart to list personal property
- » Numbered bag with patient receipt
- » Adhesive seal closure

Item No.	Description	Pkg
DYNDPVE912	9.25 x 14.5" Valuables Bag	1,000 ea/cs



Pathology Containers with Lids

- » High-grade polypropylene
- » May be frozen without cracking
- » Secure, snap-lock caps

Item No.	Description	Pkg
DYND34255	8-oz./236.6 ml	250/cs
DYND34260	16-oz./473.2 ml	100/cs
DYND34265	32-oz./946.4 ml	100 ea/cs
DYND34270	64-oz./1.9 L	50/cs
DYND34275	165-oz./4.9 L	25/cs



Small and Medium Pre-filled Formalin Containers

- » Leak resistant and 95 kPa-compliant
- » Filled halfway with 10% neutral buffered formalin fixative
- » Packaged in transport packs for easy distribution

Item No.	Description	Pkg
MCHEMPF20	0.7-oz./20 ml	256 ea/cs, 8 bx/cs
MCHEMPF40	1.4-oz./40 ml	256 ea/cs, 8 bx/cs
MCHEMPF60	2-oz./60 ml	256 ea/cs, 8 bx/cs
MCHEMPF90	3-oz./90 ml	192 ea/cs, 8 bx/cs
MCHEMPF120	4-oz./120 ml	96 ea/cs, 6 bx/cs
MCHEMPF180	6.1-oz./180 ml	64 ea/cs
MCHEMPF250	8.5-oz./250 ml	48 ea/cs
MCHEMPF500	16.9-oz./500 ml	27 ea/cs
MCHEMPF940	31.8-oz./940 ml	15 ea/cs



DYND30330



DYND30382



Specimen Containers

Item No.	Description	Pkg
DYND30382	3-oz./90 ml	400/cs
DYND30330	4-oz./120 ml	100/cs
DYND30385	4-oz./120 ml	300/cs
DYND30350	4.5-oz./130 ml	100/cs

Sterile Fluid Pathway

DYND30382	3-oz./90 ml Pneumatic-tube Safe, Click-n-Close	400/cs
DYND30331	4-oz./120 ml	100/cs
DYND30389	4-oz./120 ml Pneumatic-tube Safe, Click-n-Close	100/cs
DYND30351	4.5-oz./130 ml	100/cs

Packaged Clean, Nested in Sleeve, 25 Containers and Lids per Sleeve

DYND30335	4-oz./120 ml	500/cs
-----------	--------------	--------

Clear Plastic Urinalysis Container

DYND30100	6-oz./177.4 ml, No Lid	500/cs
DYND30102	6-oz./177.4 ml, with Pre-Printed Lids for Specimen Identification	500/cs
DYND30103	6-oz./177.4 ml, with Pre-Printed Lids for Specimen Identification	500/cs



DYND30510



Stool Specimen Container with Snap Lid

- » Red plastic container comes with a printed white lid

Item No.	Description	Pkg
DYND30510	Container with Snap Lid	250 ea/cs



MLABTP5NGRAD

Transfer Pipettes

- » Produced from a custom grade of low-density polyethylene (LDPE)
- » Unbreakable one-piece pipettes draw and dispense liquids in controlled, reproducible drop sizes
- » Packaged in a dispenser box, ideal for bench top use

Item No.	Description	Pkg
MLABTP15PED	1.5 ml, Pediatric,	5,000/cs, 10 bx/cs
MLABTP5GRAD	5 ml, Graduated, 2 mil Bulb Draw	5,000/cs, 10 bx/cs
MLABTP5LGBLB	5 ml, Graduated, 2 mil Bulb Draw, Large Bulb	5,000/cs, 10 bx/cs
MLABTP7GRAD	7 ml, Graduated, 3.2 mil Bulb Draw	5,000/cs, 10 bx/cs
MLABTP5NGRAD	5 ml, Non-Graduated, 2 mil Bulb Draw	5,000/cs, 10 bx/cs
DLS137135	3 ml, Sterile, Individually Wrapped	400/cs
DLS137038	5 ml, 3.4 ml Bulb Draw, Sterile, Individually Wrapped	400/cs



Culture Test Glass Tube

- » Made of high quality borosilicate glass
- » Well-rounded bottoms with smooth, fire polished rims
- » Packaged in boxes of 250 tubes
- » Boxes are shrink-wrapped to keep tubes particle-free and fit into most test tube dispensers

Item No.	Description	Pkg
MLAB10X75B	10 x 75 mm, 3 ml	1,000 ea/cs
MLAB12X75B	12 x 75 mm, 5 ml	1,000 ea/cs
MLAB13X100B	13 x 100 mm, 7 ml	1,000 ea/cs




Petri Dishes

- » Manufactured from pure virgin polystyrene
- » Sterilized by gamma irradiation
- » Engineered for optimum flatness
- » Highly transparent

Item No.	Description	Pkg
MLAB3501	100 x 15 mm, Sterile, Slippable	500 ea/cs
MLAB3561	100 x 15 mm, Sterile, Semi-Stackable	500 ea/cs
MLAB3571	100 x 15 mm, Sterile, Stackable	500 ea/cs
MLAB3525	150 x 15 mm, Sterile, Large Dish	500 ea/cs
MLAB3515	60 x 15 mm, Sterile, Small Dish	600 ea/cs

**Additional transfer pipettes and lab supplies available.
Contact your Medline Sales Representative today!**



MLABSLIDEP

Microscope Slides

- » Produced from premium glass sheets; ground 90° edges
- » Pre-cleaned
- » Dimensions: 25 x 75 x 1 mm

Item No.	Description	Pkg
MLABSLIDEFW	White Frosted	1,440/cs, 20 bx/cs
MLABSLIDEF	Frosted	1,440/cs, 20 bx/cs
MLABSLIDEP	Plain	1,440/cs, 20 bx/cs



MLAB140020

Black Lab Marker, Fine Tip

- » Used for writing on plastic, metal, or porcelain products
- » Autoclavable

Item No.	Description	Pkg
MLAB140020	Fine Tip, Non-Sterile	10 ea/bx

For more lab supplies, visit medline.com/physician.



Cover Glass

- » Three thicknesses:
No. 1: 0.14 mm +/- 0.015 mm
No. 1.5: 0.175 mm +/- 0.015 mm
No. 2: 0.21 mm +/- 0.02 mm
- » Each box contains a desiccant; vacuum-sealed in aluminum
- » Packaged 1-oz. per box, 10 boxes per case



Item No.	Description	Pkg
MLAB22X221	22 x 22 mm, #1 Thickness	10 ea/cs
MLAB22X2215	22 x 22 mm, #1.5 Thickness	10 ea/cs
MLAB22X222	22 x 22 mm, #2 Thickness	10 ea/cs
MLAB24X401	24 x 40 mm, #1 Thickness	10 ea/cs
MLAB24X4015	24 x 40 mm, #1.5 Thickness	10 ea/cs
MLAB24X402	24 x 40 mm, #2 Thickness	10 ea/cs
MLAB24X501	24 x 50 mm, #1 Thickness	10 ea/cs
MLAB24X5015	24 x 50 mm, #1.5 Thickness	10 ea/cs
MLAB24X502	24 x 50 mm, #2 Thickness	10 ea/cs
MLAB24X601	24 x 60 mm, #1 Thickness	10 ea/cs
MLAB24X6015	24 x 60 mm, #1.5 Thickness	10 ea/cs
MLAB24X602	24 x 60 mm, #2 Thickness	10 ea/cs



MCHM110



MCHM130

Chemicals

- » Solutions produce repeatable and accurate results
- » For in vitro diagnostic and laboratory use only

Item No.	Description	Pkg
AZRES600	Acetone, 1-gal./3.8 L	4/cs
AZRES609	Xylene, 4 x 1-gal./3.8 L/cs	4/cs
AZRES921	Wright Giemsa Stain, 32-oz./946.4 ml Bottle	1/ea
MCHM100	Deionized Water, Reagent-Grade, 5-gal./18.9 L Cube	1/ea
MCHM101	Deionized Water, Reagent Grade, 1-gal.	4/cs
MCHM110	Formalin, 10% Neutral Buffered, 4 x 1-gal./3.8 L	4/cs
MCHM111	Formalin, 10% Neutral Buffered, 5-gal./18.9 L Cube	1/ea
MCHM130	100% Reagent-Grade Alcohol, 4 x 1-gal./3.8 L	4/cs
MCHM131	95% Reagent-Grade Alcohol, 4 x 1-gal./3.8 L	4/cs
MCHM140	Alcohol, Isopropyl, 99%, 4 x 1-gal./3.8 L	4/cs
MCHM145	Alcohol, Isopropyl, 70%, 4 x 1-gal./3.8 L	4/cs
MCHM150	Methanol, 100%, 4 x 1-gal./3.8 L	4/cs



DIAGNOSTIC EQUIPMENT

Blood Pressure	28-31
Stethoscopes	32
Thermometry	33-34
Otosopes/Ophthalmoscopes.....	34-38
Patient Monitoring/Diagnostic Equipment	38-40
Defibrillators	40-41
Bladder Scanners	41-42
Diagnostic Equipment & Supplies	42-43
Spirometry.....	44



MDS5001



MedlineElite Automatic Digital Blood Pressure Unit

- » Quick, automatic systolic/diastolic blood pressure and pulse rate customizing inflation levels for each patient
- » Universal cuff to fit a range of patients
- » Easy-to-use, one-push-button design
- » Stores the 90 most recent measurements
- » Large digital readout
- » Cuff and four AA batteries included

Item No.	Description	Pkg
MDS1001	Monitor with Adult Cuff, Small & Compact	1/ea
MDS3001	Automatic Inflation with Adult Cuff, Fits Arms 8.7–11.8"/22–30 cm	1/ea
MDS3001LA	Automatic Inflation with Large Adult Cuff, Fits Arms 11.8–16.5"/30–42 cm	1/ea
MDS3001PLUS	Automatic Inflation with Adult and Large Adult Cuffs	1/ea
MDS3001U	Monitor with Universal Sized Cuff	1/ea
MDS4001	Automatic with Backlight Display, Includes Adult Cuff, Size 8.7–11.8"/22–30 cm	1/ea
MDS4001LAT	Automatic with Talking Feature and Backlight Display, Includes Large Adult Cuff, Size 11.8–16.5"/30–42 cm	1/ea
MDS5001	Bedside Monitor with Universal Sized Cuff	1/ea
MDS9970	Small Adult Cuff, Fits Arms 6.7–8.7"/17–22 cm	1/ea
MDS9971	Adult Cuff, Fits Arms 8.7–11.8"/22–30 cm	1/ea
MDS9972	Large Adult Cuff, Fits Arms 11.8–16.5"/30–42 cm	1/ea
MDS9973	Extra Large Adult Cuff, Fits Arms 16.5–18.9"/42–48 cm	1/ea
MDS9974	Universal Sized Cuff, Fits Arms 22–42 cm	1/ea



MDS3002



MedlinePro Semi-Automatic Digital Blood Pressure Monitor

- » Quick, manual inflated systolic/diastolic blood pressure unit
- » Easy-to-use, one-push-button design
- » Stores the 90 most recent measurements
- » Large digital readout
- » Cuff and four AA batteries included

Item No.	Description	Pkg
MDS3002	Monitor, Adult	1/ea
MDS3002LA	Monitor with Cuff, Large Adult	1/ea



MDS9380



MedlinePlus PVC Handheld Aneroids

- » Great quality. Great price.
- » PVC tubing and bladder
- » Nylon range finder cuff for quick, convenient patient fitting
- » Carrying case included
- » Not made with natural rubber latex

Item No.	Description	Pkg
MDS9386	Handheld Aneroid, Infant	1/ea
MDS9387	Handheld Aneroid, Child	1/ea
MDS9380	Handheld Aneroid, Adult	1/ea
MDS9388	Handheld Aneroid, Large Adult	1/ea



MDS9410



MedlineElite Handheld Aneroids

- » Our finest handheld aneroids deliver exceptional accuracy with a top-quality, zinc-alloy precision gauge
- » Easy-to-read, even in dim light
- » Shatter-resistant
- » Convenient patient fitting
- » No-stop pin action for visual calibration check and added accuracy
- » Carrying case included
- » Made with natural rubber latex

Item No.	Description	Pkg
MDS9411	Child	1/ea
MDS9410	Adult	1/ea
MDS9413	Large Adult	1/ea



MDS9380LF



MedlinePro Neoprene Handheld Aneroids

- » High-quality handheld aneroids.
- » Strong neoprene tubing and bladder
- » Nylon range finder cuff for quick, convenient patient fitting
- » Carrying case included
- » Not made with natural rubber latex

Item No.	Description	Pkg
MDS9386LF	Handheld Aneroid, Infant	1/ea
MDS9387LF	Handheld Aneroid, Child	1/ea
MDS9380LF	Handheld Aneroid, Adult	1/ea
MDS9388LF	Handheld Aneroid, Large Adult	1/ea
MDS9389LF	Handheld Aneroid, Thigh	1/ea



MDS4003



MedlinePlus Digital Wrist Blood Pressure Monitor

- » Easy-to-use automatic wrist blood pressure monitor
- » Automatically inflates and deflates
- » Provides quick digital readings of systolic/diastolic pressure and pulse rate
- » Memory with date and time
- » Color gauge for quick detection
- » Includes two AAA batteries

Item No.	Description	Pkg
MDS4003	Digital Wrist Blood Pressure Monitor	1/ea



MDSP3043



Pedometer and Heart Rate Watch

- » Takes pulse measurements
- » Offers three pedometer display options:
 - Number of steps and current time
 - Number of steps during exercise
 - Calories burned with distance walked
- » Automatically stores last 50 days of steps
- » LCD display and four-second backlight

Item No.	Description	Pkg
MDSP3043	Fitness Watch	1/ea



Wall-Mount and Mobile Aneroids

- » Large, raised numbers and a luminescent coating make these wall-mount and mobile aneroids easy to read in dim light
- » Pivoting bracket allows easy viewing
- » Nylon range finder cuff, bladder bulb, valve and 8-ft. coiled tubing included
- » Inflation system and bulb store conveniently behind gauge
- » Stand included with mobile models



MDS9400LF

Item No.	Description	Pkg
MedlinePro Neoprene Aneroids		
MDS9400LF	Wall-Mount Aneroid with Adult Cuff	1/bx
MDS9407LF	Mobile Aneroid with Adult Cuff	1/ea
MDS9418	Mobile Stand	1/ea

Item No.	Description	Pkg
MedlinePlus PVC Aneroids		
MDS9400	Wall-Mount Aneroid with Adult Latex Cuff	1/bx
MDS9407	Mobile Aneroid with Adult Latex Cuff	1/ea
MDS9418	Mobile Stand	1/ea
MDS9472	Basket for Cuff	1/ea



MDS9140



MedlineElite Compli-Mates®

- » Kit includes:
- Aneroid sphygmomanometer
 - Matching 22"/56 cm stainless steel dual-head stethoscope
 - Matching nylon carrying case

Item No.	Description	Pkg
MDS9140	Steel/Black (Stainless Steel)	1/ea



MDS9124



MedlineDeluxe Compli-Mates

- » Kit includes:
- Aneroid sphygmomanometer
 - Matching 22"/56 cm Sprague-Rappaport stethoscope
 - Matching nylon carrying case

Item No.	Description	Pkg
MDS9121	Magenta	1/ea
MDS9123	Purple	1/ea
MDS9124	Royal Blue	1/ea
MDS9125	Black	1/ea
MDS9128	Hunter Green	1/ea



MDS9106



MedlinePlus Compli-Mates

- » Kit includes:
- Aneroid sphygmomanometer
 - Matching nylon carrying case

Item No.	Description	Pkg
MDS9104	Red	1/ea
MDS9106	Navy Blue	1/ea
MDS9108	Hunter Green	1/ea
MDS9109	Orchid	1/ea
MDS9111	Magenta	1/ea
MDS9112	Teal	1/ea
MDS9113	Purple	1/ea
MDS9114	Royal Blue	1/ea



MDS9133



MedlinePro Compli-Mates

- » Kit includes:
- Aneroid sphygmomanometer
 - Matching dual-head stethoscope
 - Matching nylon carrying case

Item No.	Description	Pkg
MDS9116	Magenta	1/ea
MDS9117	Teal	1/ea
MDS9118	Purple	1/ea
MDS9119	Royal Blue	1/ea
MDS9120	Black	1/ea
MDS9133	Hunter Green	1/ea
MDS9136	Red	1/ea



Aneroid Accessories



MDS9414



MDS9630

Handheld Aneroids			Wall-Mount Aneroids	
MedlineElite	MedlinePro	MedlinePlus	MedlinePro	MedlinePlus
MDS9410	MDS9386LF	MDS9386	MDS9407LF MDS9400LF	MDS9400 MDS9407
MDS9411	MDS9387LF	MDS9387		
MDS9413	MDS9380LF	MDS9380		
	MDS9388LF	MDS9388		

Premier Gauge and Parts

Gauge (MDS9414)	•			
Glass Face (MDS9415)	•			

Standard Gauge and Parts

Aneroid Gauge (MDS9381)		•	•	
-------------------------	--	---	---	--

Bulbs and Valves

PVC Bulb (MDS91481)			•	•
Neoprene Bulb (MDS91481LF)		•	•	
Premier Bulb and Valve Assembly (MDS9630)	•			
PVC Bulb and Valve Assembly (MDS91482)		•	•	•

Two-Tube Inflation Bag and Nylon Range Finder Cuff

Neoprene, Infant (MDS91422LF)		•		•
Neoprene, Child (MDS91421LF)		•		•
Neoprene, Adult (MDS91420LF)		•		•
Neoprene, Large Adult (MDS91423LF)		•		•
PVC, Infant (MDS91422)			•	•
PVC, Child (MDS91421)			•	•
PVC, Adult (MDS91420)			•	•
PVC, Large Adult (MDS91423)			•	•
PVC, Thigh (MDS91424)			•	•

Two-Tube Seamless Inflation Bag

Neoprene, Adult (MDS9505LF)		•		•
Neoprene, Large Adult (MDS9506LF)		•		•
PVC, Adult (MDS9505)			•	•
PVC, Large Adult (MDS9506)			•	•

Nylon Range Finder Hook-and-Loop Cuffs

Child (MDS9421)		•	•	•	•
Adult (MDS9420)		•	•	•	•
Large Adult (MDS9423)		•	•	•	•

8' Coiled Tubing with Connector

PVC (MDS9475)					•
Neoprene (MDS9475LF)				•	

Miscellaneous Parts

Metal Connector/Coupler (MDS9481)				•	•
Basket for Wall-Mounted or Mobile Aneroids (MDS9472)				•	•
Mobile Stand (MDS9418)				•	•



MDS926501



Synergy Dual-Frequency Stethoscope

- » High end, stainless steel stethoscope that provides excellent auscultation and acoustics
- » Binaurals and chest piece help trap sound
- » Durable, thick-walled tubing to eliminate sound interference
- » Soft, narrow eartips for everyday comfort

Item No.	Description	Pkg
MDS926501	Synergy, Dual-Head, Black	1/ea
MDS926502	Synergy, Dual-Head, Blue	1/ea
MDS926503	Synergy, Dual-Head, Hunter Green	1/ea
MDS926504	Synergy, Dual-Head, Burgundy	1/ea
MDS926601	Synergy, Single-Head, Black	1/ea
MDS926801	Synergy, Cardiology, Black	1/ea



MDS926201



Dual-Head Stethoscopes

- » High-quality dual-head stethoscopes feature a non-chill ring on the bell for increased patient comfort
- » Thick-walled 22-inch/56 cm tubing minimizes background interference

Item No.	Description	Pkg
MDS926201	Black	1/ea
MDS926202	Gray	1/ea
MDS926203	Blue	1/ea
MDS926205	Lavender	1/ea
MDS926206	Red	1/ea
MDS926207	Hunter Green	1/ea
MDS926208	Yellow	1/ea
MDS926209	Pink	1/ea
MDS926210	Orange	1/ea



MDS926301



Sprague-Rappaport Stethoscopes

- » Component kit includes two sets of eartips, large and small diaphragm, and adult, pediatric, and infant bells
- » 22-inch/56 cm dual-tubing transmits even the faintest sounds

Item No.	Description	Pkg
MDS926301	Black	1/ea
MDS926302	Gray	1/ea
MDS926303	Blue	1/ea
MDS926305	Lavender	1/ea
MDS926306	Red	1/ea
MDS926307	Hunter Green	1/ea
MDS926311	Yellow	1/ea
MDS926312	Orange	1/ea



MDS92250



Stainless Steel Stethoscopes

- » Quality, handcrafted stainless steel stethoscopes
- » Color coordinated non chill rings
- » Durable 22-inch/56 cm tubing
- » Units come with two extra sets of eartips and one extra diaphragm

Item No.	Description	Pkg
MDS92250	Adult, Gray	1/ea
MDS92260	Adult, Black	1/ea
MDS92270	Adult, Blue	1/ea
MDS92280	Adult, Burgundy	1/ea
MDS92290	Adult, Pink	1/ea
MDS92500	Cardiology, Black	1/ea
MDS92150	Neonatal, Gray	1/ea
MDS92160	Neonatal, Pink	1/ea
MDS92200	Pediatric, Gray	1/ea



MDS926101



Single-Head Stethoscopes

- » Single-head stethoscopes have a lightweight design and thick-walled 22-inch/56 cm tubing to reduce background interference

Item No.	Description	Pkg
MDS926101	Black	1/ea
MDS926102	Gray	1/ea
MDS926103	Blue	1/ea
MDS926105	Lavender	1/ea
MDS926106	Red	1/ea
MDS926107	Hunter Green	1/ea
MDS926108	Yellow	1/ea
MDS926109	Pink	1/ea
MDS926110	Orange	1/ea



MDS131040



Disposable Pen Lights

- » High-intensity exam light
- » Lightweight and disposable
- » Can be clipped to pocket

Item No.	Description	Pkg
MDS131040	Disposable Pen Light, White	6/pk
MDS131040H	Disposable Penlight, White	1/ea



A06000200



Thermoscan Pro Thermometer

- » Results in 2–3 seconds
- » Position the probe tip in the ear canal and push a button for a virtually instantaneous and accurate reading
- » Allows a wide range of patient temperatures in a variety of ambient temperatures
- » Ergonomic design is comfortable and durable
- » Wall-mount base and probe covers are available
- » Designed for patient comfort, these disposable covers help protect the unit, user and patient from the potential risk of cross contamination
- » Three-year warranty

Item No.	Description	Pkg
W-A06000200	Thermometer w/Small Cradle	1/ea
W-A06000300	Thermometer w/Large Cradle	1/ea

Disposable Probe Covers

A06000005	Thermoscan 6000 Probe Covers	200/bx, 25bx/cs, 5,000/cs
W-A06000005ZZ	Thermoscan 6000 Probe Covers	200/bx



MDS9950



30-Second Digital Thermometers

- » Accurate digital thermometers
- » Results in 30 seconds
- » LCD screen
- » Simple push-button start
- » Disposable sheaths

Item No.	Description	Pkg
Oral		
MDS9950	Oral °F & °C	144/cs
MDS9950Z	Oral °F & °C	12/bx
MDS9950H	Oral °F & °C	1/ea
MDS9950B	Oral °F & °C, Bulk	200/cs
MDS9950BZ	Oral °F & °C, Bulk	20/bx
MDS9607	Oral Sheaths	100/bx
MDS9928	Oral Kit w/20 Sheaths, °F & °C	1/ea

Rectal

MDS9952	Rectal °F & °C	144/cs
MDS9952Z	Rectal °F & °C	12/bx
MDS9952H	Rectal °F & °C	1/ea
MDS9952B	Rectal °F & °C Bulk	200/cs
MDS9952BZ	Rectal °F & °C Bulk	20/bx
MDS9608	Pre-lubricated Rectal Sheaths	100/bx
MDS9929	Rectal Kit w/20 Sheaths, °F & °C	1/ea

Flex-tip

MDS99901	Flex-tip, Water-Resistant, °F & °C	240/cs
MDS99901Z	Flex-tip, Water-Resistant, °F & °C	12/bx
MDS99901H	Flex-tip, Water-Resistant, °F & °C	1/ea
MDS99902	Flex-tip, Water-Resistant, °F & °C	240/cs
MDS99902Z	Flex-tip, Water-Resistant, °F & °C	12/bx
MDS99902H	Flex-tip, Water-Resistant, °F & °C	1/ea



MDS8700



RightTemp™ Tympanic Thermometer

- » Clinically proven accuracy
- » Non-invasive thermometer takes temperature from the tympanic membrane
- » All-in-one cradle probe cover dispensing system
- » Probe cover eject button
- » Readings in °F/°C

Item No.	Description	Pkg
MDS8700	Ear Thermometer and Cradle Dispenser	1/ea
MDS8701	Probe Covers f/ RightTemp Thermometer	5,000/cs
MDS8701Z	Probe Covers f/ RightTemp Thermometer	200/bx
MDS8701H	Probe Covers f/ RightTemp Thermometer	20/pk



MDSNOTOUCH



No Touch Forehead Thermometer

- » Takes measurements when holding the probe parallel to the forehead and approximately 1-inch away
- » Helps lower risk of cross contamination by eliminating patient contact
- » No probe covers needed
- » Displays green light for healthy temperature and red light for fever

Item No.	Description	Pkg
MDSNOTOUCH	Forehead Thermometer	1/ea



Tympanic Thermometer

- » Innovative, quick-response ear thermometer with quick probe release
- » One-touch operation
- » Mechanical probe cover quick-release allows easy probe cover change
- » Safe to use on patients of any age
- » Readings in Fahrenheit and Celsius
- » Takes an accurate, non-invasive temperature from the tympanic membrane



MDS9700

Item No.	Description	Pkg
MDS9700	Thermometer	1/ea
MDS9701	Probe covers	100/bx



Digital Temple Thermometer

- » Instant read digital temple thermometer accurately measures body heat using R.A.T.E.™ technology (Rapid, Accurate, Temperature Establishment)
- » Sensors rapidly track heat flow generated from blood vessels to the skin's surface and convert measurements to body temperature
- » Easy-to-use on asleep or alert patients
- » Meets accuracy standards
- » Large digital display
- » Uses two AAA batteries
- » Memory stores previous reading



MDS9698L

Item No.	Description	Pkg
MDS9698L	Temple Thermometer	1/ea
MDS9699	Probe Covers	250/bx



SureTemp® Plus 690 Thermometer

- » Captures accurate, repeatable temperatures in approximately four to six seconds for oral, 10 to 13 seconds for rectal, and 12 to 15 seconds for adult axillary
- » Optional removable/interchangeable probes and probe wells
- » Takes approximately 6,000 readings on three AA batteries (not included)



W-A01690200

Item No.	Description	Pkg
W-A01690200	Oral/Axillary Thermometer	1/ea
W-A02892000	Rectal Probe	1/ea
W-A02893000	Oral Probe	1/ea
W-A05031750	Disposable Probe Covers	250/bx, 30 bx/cs, 7,500/cs
W-A05031750Z	Disposable Probe Covers	250/bx



MDS3703

ri-scope L2 Otoscopes

- » Fiber optic illumination
- » Bi-directional swiveling, sealing optical glass lens with 3-fold magnification

Item No.	Description	Pkg
MDS3716	3.5V LED for Batteries, with AA-Handle	1/ea
MDS3703	3.5V LED for Batteries, with C-Handle	1/ea



MDS3704

ri-scope L3 Otoscope

- » Mini operation lens included, increasing magnification from 3x to 4x while allowing for easy insertion of instruments
- » Easy-to-operate SpecEjec specula ejection device
- » Enlarged conduit for instruments

Item No.	Description	Pkg
MDS3704	3.5V, LED, for 2 Lithium Batteries, C-handle	1/ea
MDS3701	2.5V, Xenon, for 2 Alkaline Batteries, C-handle	1/ea



MDS3725



EliteVue Otoscope

- » 5x magnification allows for significantly improved view of the auditory canal and full view of the tympanic membrane
- » Reflection-optimized high performance glass optics
- » Wheel to adjust focal depth
- » Suitable for pneumatic otoscopy (supplied without ball)
- » Must use SpecEjec speculas MDS10803533 (3mm) or MDS10803534 (4mm)

Item No.	Description	Pkg
MDS2200202E	2.5V, XL, for 2 Alkaline Batteries, C-handle	1/ea
MDS2200204E	2.5V, LED, for 2 Alkaline Batteries, C-handle	1/ea



MDS3725



ri-scope L2 Ophthalmoscope

- » Diopter wheel with 29 corrective lenses (plus 1-10, 12, 15, 20, 40, minus 1-10, 15, 20, 25, 30, 35)
- » Easy-to-operate aperture hand-wheel with semi-circle, small/medium/large circle, fixation star, and slit
- » Includes filter wheel engageable for all apertures with symbol display, red-free filter, blue filter and polarization filter
- » Optimized high-performance optics with aspherical condenser lens

Item No.	Description	Pkg
MDS3725	3.5V XL Batteries, with C-Handle	1/ea



MDS3726



ri-scope® L3 Ophthalmoscope

- » Corrective lenses for 89 diopter values (plus 1-45 in single steps, minus 1-44 in single steps)
- » Easy-to-operate aperture hand-wheel with grid, semi-circle, small/medium/large circle, fixation star, and slit
- » Includes filter wheel engageable for all apertures with symbol display, red-free filter, blue filter and polarization filter
- » Optimized high-performance optics with aspherical condenser lens
- » Focusing device with zero retention force for the quick correction of ametropias

Item No.	Description	Pkg
MDS3726	3.5V XL Batteries, with C-Handle	1/ea



MDSK10321



MDSK10331



Disposable Universal Ear Specula

- » Black plastic and single use accessory
- » Fits all popular brands (such as Welch Allyn)
- » Packed in individual sleeves of 25 for convenient restocking of the speculum dispenser


Item No.	Description	Pkg
MDSK10321*	2.5 mm	8,500/cs
MDSK10331*	4.25 mm	8,500/cs
MDS10803533	3 mm SpecEjec Speculas for EliteVue and L3 Heads Only	1,000/bg
MDS10803534	4 mm SpecEjec Speculas for EliteVue and L3 Heads Only	1,000/bg

*Available in child items H (34/pk) and Z (850/bg)



Diagnostic Center

- » Customize your center for what you want
- » Fully modular system, customized to the end user needs; can be expanded/upgraded at any point in time
- » Non-slip, ergonomically designed handles with patented technology for 100% light output with variable light adjustment and automatic safety shut-down after two minutes
- » Easy cleaning of the instruments
- » Illuminated diode signals whether handle is on or off
- » All SKUs ending in "B" come with grey mounting board



Model Featured:
MDS36533018B

Medline Integrated Diagnostic System and Transformer Sets Components Chart

Item No.	Wall Transformer	Otosopes				Ophthalmoscopes		Thermometer		Specula Dispenser	Blood Pressure	
	MDS3653300	MDS10565301	MDS10580301	MDS10513301E	MDS10571301	MDS10573301	MDS3656301	MDS3658	MDS3654	MDSL3655109	MDS1784003	
		L2 LED	L2 XENON	Xenon EliteVue Otoscope	L2 XENON	L3 XENON	Tympanic Therm	Covidien RPT		Aneroid	Automatic	
MDS36533001	•	•			•							
MDS36533002	•		•		•							
MDS36533003	•	•			•				•			
MDS36533004	•		•		•				•			
MDS36533005	•	•			•				•	•		
MDS36533006	•		•		•				•	•		
MDS36533007	•	•				•			•	•		
MDS36533008	•	•			•		•		•	•		
MDS36533009	•	•			•		•		•	•		
MDS36533010	•	•				•	•		•	•		
MDS36533011	•	•			•				•		•	
MDS36533012	•		•		•				•	•	•	
MDS36533013	•	•				•			•		•	
MDS36533014	•	•			•		•		•	•		
MDS36533015	•	•			•		•		•		•	
MDS36533016	•	•				•	•		•	•		
MDS36533017	•	•				•	•		•		•	
MDS36533018	•			•	•		•		•	•		
MDS36533018B	•			•	•		•		•	•		
MDS36533019	•			•	•		•		•	•		
MDS36533019B	•			•	•		•		•	•		

If you don't see what you're looking for, please contact your Medline rep or diagnostics.medline.com to create you own wall system.



WelchAllyn®
Green Series 777 Integrated Diagnostic System

- » Convenient, modular package for fast, accurate diagnoses
- » Includes wall aneroid, ophthalmoscope, otoscope, speculum dispenser and thermometer
- » Customized to meet your needs

Item No.	Description	Pkg
W-A791MXPROBP	Includes: GS 777 Wall Transformer, Coaxial Ophthalmoscope, MacroView Otoscope, Wall Aneroid, KleenSpec Dispenser, SureTemp Plus Thermometer, and Connex ProBP Device	1/ea
W-A77791MX	Includes: GS 777 Wall Transformer, Coaxial Ophthalmoscope, MacroView Otoscope, Wall Aneroid, KleenSpec Dispenser, and SureTemp Plus Thermometer	1/ea
W-A777912MPX	Includes: GS 777 Wall Transformer, PanOptic Ophthalmoscope, MacroView Otoscope, Wall Aneroid, KleenSpec Dispenser, and SureTemp Plus Thermometer	1/ea
W-A77792M	Includes: GS 777 Wall Transformer, Coaxial Ophthalmoscope, MacroView Otoscope, Wall Aneroid, and KleenSpec Dispenser	1/ea
W-A777921MP	Includes: GS 777 Wall Transformer, PanOptic Ophthalmoscope, MacroView Otoscope, Wall Aneroid, and KleenSpec Dispenser	1/ea
W-A77710	Wall Transformer	1/ea

WelchAllyn®
3.5V Standard Ophthalmoscope

- » Halogen light provides true tissue color and consistent, long-lasting illumination, and six apertures for generalist use: micro, small, and large spot sizes, fixation target, slit aperture, and red-free filter



W-A11710

Item No.	Description	Pkg
W-A11710	3.5V Standard Ophthalmoscope (Head Only)	1/ea



W-A11720

WelchAllyn®
3.5V Coaxial Ophthalmoscope

- » Coaxial Vision System facilitates ophthalmoscopy by enabling easier entry into the eye, a larger field of view, and reduced glare compared to standard ophthalmoscopes

Item No.	Description	Pkg
W-A11720	3.5V Coaxial Ophthalmoscope (Head Only)	1/ea



W-A20200

WelchAllyn®
Pneumatic Otoscope

- » Halogen HPX® lamp provides light output for true tissue color and consistent, long-lasting illumination
- » Rotating lens with ultraseal enables easy pneumatic otoscopy

Item No.	Description	Pkg
W-A20200	3.5V Pneumatic Otoscope (Head Only)	1/ea
W-A21501	Insufflation Bulb Without Tip	1/ea



W-A71900

WelchAllyn®
3.5V Rechargeable Handles

Item No.	Description	Pkg
W-A71900	3.5V, Lithium Ion	1/ea
W-A71000A	3.5V, NiCad	1/ea

For other Welch Allyn diagnostic equipment and parts, please contact your Medline Sales Representative



Welch Allyn

3.5V Diagnostic Otoscope

- » Halogen HPX lamp provides light output for true tissue color and long-lasting performance
- » Sealed system for pneumatic otoscopy and accurate diagnosis of otitis media

Item No.	Description	Pkg
W-A25020	3.5V Diagnostic Otoscope (Head Only)	1/ea
W-A23820	MacroView Otoscope with Throat Illuminator (Head Only)	1/ea
W-A21504	Insufflator Bulb with Tip	1/ea

Accessories

W-A52401H	KleenSpec® Dispenser	1/ea
W-A52434UZ	0.17"/4.25 mm)	850/bg
W-A52432UZ	0.11"/2.75 mm) Pediatric	850/bg
W-A72200	Replacement NiCad Rechargeable Battery (Black)	1/ea
W-A72300	Orange Replacement NiCad Battery	1/ea

Replacement Bulbs

W-A03000U	Replacement Bulb, Halogen 3.5V, Ophth.	1/ea
W-A03100U	Replacement Bulb, 3.5V for Otoscope	1/ea
W-A06500U	Replacement Bulb for Macroview Otoscope	1/ea
BHM03300	Replacement Bulb, #33	6/bx
BHM03100	Replacement Bulb, #31	6/bx
BHM03000	Replacement Bulb, #30	6/bx
BHM03000H	Replacement Bulb, #30	1/ea



Masimo

Pronto Spot-Check Hemoglobin (SpHb) Sensors

- » Use these sensors with the Pronto MM09167 and starter bundles

Item No.	Description	Pkg
MM03419	Pediatric, 3/7.6 cm Rainbow DCI, 200 Spot Checks	1/ea
MM03418	Adult, 3/7.6 cm Rainbow DCI, 200 Spot Checks	1/ea
MM03421	Pediatric, 3/7.6 cm Rainbow DCI, 400 Spot Checks	1/ea
MM03420	Adult, 3/7.6 cm Rainbow DCI, 400 Spot Checks	1/ea



Masimo

Masimo® Pronto® Spot-Check Tester

- » Provides noninvasive and quick spot-check testing of total hemoglobin (SpHb), SpO₂, pulse rate and perfusion index
- » Sensors included in the starter bundle or purchased separately

Item No.	Description	Pkg
MM09167	Pulse CO-Oximeter Device	1/ea

Pronto with 200 Tests

MM03581	Pediatric	1/ea
MM03580	Adult	1/ea

Pronto with DCIP, 400 Tests

MM03583	Pediatric DCIP	1/ea
---------	----------------	------



Masimo

Handheld Pulse Oximeters

- » Masimo SET Measure-through Motion and Low Perfusion™ pulse oximetry for oxygen saturation (SpO₂), pulse rate (PR), and perfusion index (PI)
- » Lightweight handheld device with FastStart™, to allow for rapid measurement at startup
- » Perfusion Index (PI) is an assessment of the pulsatile strength at a specific monitoring site (e.g. the hand, finger, or foot), and, as such, PI is an indirect and noninvasive measure of peripheral perfusion
- » Signal I.Q. (SIQ) provides an assessment of the confidence in the displayed SpO₂ value
- » Audible and visual alarms for Sensor Off and Low Battery

Item No.	Description	Pkg
----------	-------------	-----

Rad-5V

MM09540	Pediatric	1/ea
MM09711	Adult	1/ea



Vital Signs Monitor

- » One monitor that provides spot and continuous measurement for NIBP, SpO₂ and temperature
- » Portable and lightweight
- » Easy-to-read and user-friendly monitor
- » Available in different configurations to meet your clinical needs

Choose Your Configuration

Item No.	BP	SpO ₂	Oral Temp	Ear Temp	Printer	Stand	Pkg
MDSM3ASN1	•	•					1/ea
MDSM3ASNTC	•	•	•				1/ea
MDSM3ASNTCP	•	•	•		•		1/ea
MDSM3ASNTTP	•	•		•	•		1/ea
MDSM3ASNTT	•	•		•			1/ea
MDSM3ASN	•	•				•	1/ea
MDSM3ASNT	•	•	•			•	1/ea
MDSM3ASNTP	•	•	•		•	•	1/ea
MDSM3ASNTPIR	•	•		•	•	•	1/ea
MDSM3ASNTTR	•	•		•		•	1/ea

Choose Your Accessories

Item No.	Description	Pkg
MDSM3ASTAND01	Rolling Stand w/Storage Basket	1/ea
MDSM3AWALL01	Wall Mount w/Storage Basket	1/ea
MDSSPO2ADL	Adult Pulse Oximeter Finger Sensor, Reusable	1/ea
MDSSPO2DISP**	Adult Pulse Oximeter Finger Sensor, Disposable	1/ea
MDSSPO2PED**	Pediatric Pulse Oximeter Finger Sensor, Reusable	1/ea
MDSSPO2PEDDSP**	Pediatric Pulse Oximeter Finger Sensor, Disposable	1/ea
MDSSPO2EAR**	Edan Pulse Oximeter Finger Sensor	1/ea
MDSM3ACUFFI*	Edan Blood Pressure Cuff, Infant (14–21 cm)	1/ea
MDSM3ACUFFC*	Edan Blood Pressure Cuff, Child (20.5–28 cm)	1/ea
MDSM3ACUFFA*	Edan Blood Pressure Cuff, Adult (25–35 cm)	1/ea
MDSM3ACUFFL*	Edan Blood Pressure Cuff, Large Adult (33–47 cm)	1/ea
MDSM3ACUFFT*	Edan Blood Pressure Cuff, Thigh (42–54 cm)	1/ea
MDSM3AWRRNTY	Additional 1-Year Extended Warranty	1/ea
KDL500500	Filac FasTemp Thermometer Probe Covers	500/cs
KDL500500Z	Filac FasTemp Thermometer Probe Covers	20/bx
MDSM3A7208058	Ear Probe Covers	200/bx
MDSM3A78035	Printer Paper	1/ea


*Medline cuffs also available.

**Will need to purchase MDSSPO2EXT

Additional Parts

Item No.	Description	Pkg
MDSBPTUBE	NIBP Extension Tube, 3 meters	1/ea
MDSSPO2EXT	SpO ₂ Extension Cable, 2 meters	1/ea
MDSM3A101927	Replacement Single Caster Wheel for Vital Signs Stand	1/ea
SNT91000300	Standard Power Cord for Vital Signs Monitor	1/ea
MDSBATTERY42	Rechargeable Lithium Ion Battery, 14.8V, 2.5AH	1/ea
MDSM3ABASKET1	Primary Basket for Vital Signs Stand	1/ea
MDSM3ABASKET	Secondary Basket for Vital Signs Stand	1/ea





MDSM611903



MedGel® Stress/Diaphoretic Foam Electrode

- » Designed for situations where a strong adhesive is required
- » Clear tracings with our high performance, specially formulated gel and extended surface area
- » Fluid-impervious
- » Easy removal with full-edge lift tab

Item No.	Description	Pkg
MDSM611903	ECG Stress Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611905	ECG Stress Foam	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611950	ECG Stress Foam	50/pk, 1,000/ct, 6,000/cs



MDSM616101



MedGel® Diagnostic Resting Tab Electrodes

- » Disposable tab electrode manufactured with a silver/silver chloride gel for reliable ECG signals
- » Large, solid conductive area offers maximum adhesion to maintain contact with skin
- » Specially designed for use with resting ECGs

Item No.	Description	Pkg
MDSM616101	ECG Resting Tab, Small	100/pk, 500/bx, 5000/cs
MDSM616101Z	ECG Resting Tab, Small	500/bx
MDSM616201	ECG Resting Tab, Adult	100/pk, 500/bx, 5000/cs
MDSM616201Z	ECG Resting Tab, Adult	500/bx



MDS0157040172



ECG Adapter Clips

- » Wide contact area for exceptional trace quality

Item No.	Description	Pkg
MDS0157040173	Replacement Alligator Clips for Leadwire Cable (Resting)	10/pk
MDS0157040172	Replacement Banana Clips for Leadwire Cable (Resting/Snap)	10/pk



CRU911447



Ultrasound Probe Covers

- » Helps protect equipment, prevent cross-contamination and save disinfection time
- » Wide open end for inserting the transducer
- » Pre-gelled on the inside for faster setup
- » Disposable for convenience and savings

Item No.	Description	Pkg
CRU911447	3D/4D, Endocavity	100/bx
PLI3801Z	2.5"/1.75 x 9.5"/6.4 cm/4.4 x 24.1 cm, Latex Free	100/bx
PLI3803Z	3D/4D, Latex Free	100/bx



DEFM5066HOME

PHILIPS

OnSite AEDs

- » Designed for the ordinary person in the extraordinary moment—easy-to-use device that arrives ready-to-rescue
- » Built with on-demand CPR coaching to provide guidance during resuscitation

Item No.	Description	Pkg
DEFM5066HOME	Philips OnSite AED	1/ea
DEFM5066PACK	Philips OnSite AED Package with Wall Cabinet and Carrying Case	1/ea

Philips Replacement Accessories

DEFM5071A	Adult Pads	1/ea
DEFM5072A	Pediatric Pads	1/ea
DEF5070A	OnSite Battery Pack	1/ea



ZOL2110201CC

ZOLL

AED Plus Defibrillators

- » Equipped with Real CPR Help, the AED Plus provides real-time feedback on the quality of CPR

Item No.	Description	Pkg
ZOL21110201CC	AED Plus Carrying Case Package (Semi-Automatic)	1/ea
ZOL21110201PK	AED Plus Wall-Cabinet Package (Semi-Automatic)	1/ea

Zoll® AED Replacement Accessories

ZOL8900080001	Adult Pads	1/ea
ZOL8900081001	Pediatric Pads	1/ea
ZOL8000080701	Battery	1/ea



PHY0403000149



LIFEPAK CR® Plus AED

- » Advanced technology trusted by emergency medical professionals
- » Designed specifically for the first person to respond to a victim of sudden cardiac arrest
- » Two steps required turn the device on and apply the electrode pads
- » Includes carry case and two sets of adult electrodes
- » Best for multi-unit AED programs
- » Eight-year warranty

Item No.	Description	Pkg
PHY0403000149	Fully Automatic	1/ea
PHY80403148	Semi-Automatic	1/ea



PHY0427000134



LIFEPAK Express®

- » Includes one set of adult electrodes
- » Carry case not included
- » Ideal for any physician office or LTC facility looking to protect patients and visitors
- » Five-year warranty

Item No.	Description	Pkg
PHY0427000134	Semi-Automatic	1/ea



PHY9425000023



LIFEPAK 1000

- » Rugged, most durable LIFEPAK AED primarily used by first responders in loud or harsh environments
- » Five-year warranty
- » Graphical display model or available with 3-wire ECG
- » Programs under the direction of on-site clinical professionals or users with higher training levels

Item No.	Description	Pkg
PHY9425000023	Graphical Display	1/ea



DFTDCF1310EN



Lifeline Series

- » Shows you what to do in an emergency with a full-motion color video
- » Easy to maintain and upgradeable for future upgrades and enhancements
- » Meet U.S. military “drop and shock” specifications
- » Durable and rugged with an ergonomic design

Item No.	Description	Pkg
DFTDCF1310EN	Lifeline View AED® (Semi-Automatic)	1/ea
DFTDCFA100EN	Lifeline AED (Semi-Automatic)	1/ea
DFTDCFA120	Lifeline AUTO AED (Fully-Automatic)	1/ea

Exclusively from Medline

In our commitment to bringing you high-quality products at a good value, Medline exclusively offers the BioCon bladder scanners. These portable, 3-D ultrasound devices provide quick, safe and non-invasive measurements of bladder volume and post-void urinary retention (PVR) to help in the daily care of patients.

CAUTI Prevention

Research has shown that 80% of urinary tract infections acquired at healthcare facilities are associated with an indwelling urethral catheter* What’s more, Medicare no longer reimburses for treatment of CAUTI if it happens while a patient is hospitalized, giving facilities a major incentive to prevent it.

*Lo E, Nicolle L, Classen D, Arias A, Podgorny K, Anderson DJ, et al. SHEA/IDSA practice recommendation: strategies to prevent catheter-associated urinary tract infections in acute care hospitals. Infect Control Hosp Epidemiol. 2008;29:S41–S50. Available at <http://www.jstor.org/stable/10.1086/675718>. Accessed August 8, 2016.



MDSMRC900



BioCon 900 Ultrasonic Bladder Scanner with Printer

- » Built-in versatility with ability to measure bladder volume in three seconds
- » Bladderpoint™ targeting, self-diagnosis without need for calibration, and five-year warranty with drop protection coverage
- » Comes with training materials, competency testing

Item No.	Description	Pkg
MDSB10900	Bladder Scanner with Printer	1/ea
MDSMRC900	Stand with Basket	1/ea



MDSB10750

BioCon 750 Bladder Scanner

- » Easy-to-use and durable, allowing you to get quick and accurate bladder volume readings
- » Bladderpoint™ patent pending bladder targeting technology
- » Self diagnostic feature that checks the hardware of the device
- » Industry Leading Support with on-going education, training and service through the device's lifespan
- » 5-year warranty that covers drop protection

Item No.	Description	Pkg
MDSB10750	Scanner w/Printer	1/ea
MDSB10700STND	Stand	1/ea



QRSZ70000301



Universal ECG™

- » Connect Universal ECG to a PC via USB to acquire, store and analyze up to 12 channels of high-resolution ECG data
- » Interfaces with Office Medic, offering eight report options, which are transferable to PDF, JPEG or TIFF formats
- » Uses the Louvain algorithm to obtain accurate and reliable ECG data
- » Weighs less than 12-oz./354.9 ml

Item No.	Description	Pkg
QRSZ70000301	12-Channel ECG	1/ea



All-in-One Holter Software Kit

- » Contains all essentials needed to install, connect and analyze data from the Holter monitor with Holter LX Analysis
- » Includes: Holter software CD software key (dongle), media reader and extension cable
- » Three-year warranty on hardware and software

Item No.	Description	Pkg
QRS70100000	Holter Software Kit	1/ea



Medline 12 Channel Interpretive ECG Machine

- » Interpretive, 12 channel ECG unit streamlines patient monitoring with accurate recordings
- » High resolution color display screen
- » One-touch operation and shortcut keys
- » Portable
- » 2-year warranty



MDSECGSE1201 with MDSECGSTAND01

Item No.	Description	Pkg
MDSECGSE1201	Monitor, ECG, 12 Lead, SE 1201	1/ea
MDSECGSTAND01	Stand, ECG, F/SE 1201	1/ea
GCO32017349	Z-Fold ECG Replacement Paper	20 pd/cs
MDS0157040172	Replacement Banana Clips for Leadwire Cable (Resting/Snap)	10/pk
MDS0157040173	Replacement Alligator Clips for Leadwire Cable	1/pk
MDS0157107048	Replacement Leadwire Cable	1/ea



QRS70000500



Q200/HE™ Combination Holter/Event Recorder

- » 14-day Holter recorder plus 30-day Event monitor can be deployed to a patient in either Holter or Event mode
- » Large, easy view LCD screen, soft-keys and a very simple user interface
- » Low Cost, Shielded Cables
- » Patient data is encoded on a memory card
- » Reliable digital technology is used throughout recording to enhance service life and lower costs

Item No.	Description	Pkg
QRS70000500	Q200/HE Combination Holter/Event Recorder	1/ea



WelchAllyn

PC-Based Resting ECG

- » Large display saves time and eliminates repeat ECGs
- » Automatic ECG measurements speed physician review and interpretation
- » Programmable report formats print on plain paper
- » On-screen ECG report editor makes changes and additions to diagnosis easy
- » Serial comparison allows previous ECGs to be overlaid and compared for changes
- » Optional ECG interpretation available
- » Can be upgraded with stress, spirometry and/or ABPM modules for a complete cardiopulmonary suite

Item No.	Description	Pkg
W-ACPRUIUBD	Digital ECG with Interpretation	1/ea
W-ACPRUNUBD	Digital ECG Without Interpretation	1/ea
W-ACP1501ENB	CP150 Box ECG, No Interpretation	1/ea
W-ACP150A1ENB	CP150 ECG, Includes Interpretation	1/ea
W-A150AS1ENB	CP150 Box ECG, Includes Interpretation, with Spirometry	1/ea
W-A105660	Spirometry Upgrade f/ CP150	1/ea

EDAN



PC-Based ECG with Software and Sampling Box

- » Digital 12- or 15-lead resting ECG solution with one-button operation
- » Includes reliable SEMIP® algorithm
- » Color-coded lead map and ECG waveforms facilitate instant signal quality indications
- » The built-in ECG glossary consists of rich diagnostics templates with detailed feature descriptions and integrated results

Item No.	Description	Pkg
MDSECGDP12	12-Lead ECG Sampling Box	1/ea
MDSECGDX12	12-Lead ECG Sampling Box, Wireless	1/ea
MDSECGDE15	15-Lead ECG Sampling Box	1/ea

EDAN



SE-2003/SE-2012 Holter Monitor

- » Has advanced software analysis including pacemaker detection in a lightweight design
- » IP27 rated to help prevent unexpected damage from water or dust
- » OLED screen displays waveforms in real time for up to eight days on one battery

Item No.	Description	Pkg
MDSECG2003	3 Channel Monitor w/software	1/ea
MDSECG2012	12 Channel Monitor w/software	1/ea



HCS93250



Incentive Spirometers

- » Helps patients perform and monitor breathing exercises
- » One-way valve ensures patients inhale rather than exhale
- » Highly visible pistons and universal graphics
- » Two volume capacities available

Item No.	Description	Pkg
HCS93250	2,500 ml	12/cs
HCS93250H	2,500 ml	1/ea
HCS93500	5,000 ml	12/cs
HCS93500H	5,000 ml	1/ea



QRSZ70000101



Orbit™ Digital Spirometer

- » Improved accuracy, portability and USB 2.0 Connectivity
- » Easily connects with a Windows-based PC and rests on the tabletop

Item No.	Description	Pkg
QRSZ70000101	Portable Orbit Spirometer	1/ea
QRSZ50002608	Mouthpiece, Pre-Calibrated	100/bx



MIM40000026

IQspiro® Digital Spirometer

- » Cost-effective device for detecting and managing pulmonary disease with speed, accuracy and ease-of-use
- » Runs on Windows® PC or laptop computer
- » Connects to today's top EHR/EMR systems or with IQmanager®

Item No.	Description	Pkg
MIM40000026	IQspiro Digital Spirometer	1/ea
MIM21001206Z	Mouthpiece, IQspiro Digital Spirometer	25/bx



Replacement Spirotubes

Item No.	Description	Pkg
S-D297000Z	Spirometer, Type A	100/bg
S-D298020	Pneumotach, FloSense Bidirectional	25/bx
S-D298040050	Pneumotach, FloSense II, Disp., for Spirometer	50/bx



HCSPF100

Peak Flow Meters

- » Help patients monitor asthma
- » Individually calibrated to ensure accuracy
- » Not made with natural rubber latex

Item No.	Description	Pkg
HCSPF100	Flow Meter, Airzone	12/cs
HCSPF100H	Flow Meter, Airzone	1/ea



HUD1801



POCKETPEAK® Peak Flow Meter

- » Compact and easy-to-use and comes complete with durable, sliding zone indicators
- » Colored markers allow patients to easily monitor changes in airflow, and healthcare professionals to recommend a treatment plan of action
- » Meets National Asthma Education Program and National Heart, Lung, and Blood Institute standards
- » ABS plastic housing and tempered stainless steel vane are completely dishwasher-safe
- » Ergonomic design reduces user confusion and misuse
- » Patient education video available

Item No.	Description	Pkg
HUD1801	Universal Range	12/cs
HUD1801H	Universal Range	1/ea
HCSPF100	Meter, Airzone	12/cs



EXAM ROOM SUPPLIES

Paper Products & Supplies.....	46-49
Trash Liners.....	49-50
Table Paper & Linens.....	50-52
Disposable Protective Apparel.....	53-55
Patient Apparel.....	56-58
Lab Coats.....	59-60
Scrubs.....	61
Cleaners & Disinfectants.....	62-67
Batteries.....	67
Eyewash.....	68
Fluid Clean-Up Kit.....	68



NON245288



NON245277

NON243276



Premium Facial Tissue

- » Medline's softest, strongest and most absorbent 2-ply tissue
- » Boxes come in stylish designs
- » Boxes have a shield to repel dust particles

Item No.	Description	Pkg
NON245277	8 x 8.3"	100 sheets/bx, 36 bx/cs
NON245277Z	8 x 8.4"	100 sheets/bx
NON245278	8 x 8.3"	85 sheets/bx, 36 bx/cs
NON245288	Animal Pals, Cube Box	75 sheets/bx, 36 bx/cs



NON243277



Standard Facial Tissue

- » Two-ply tissue from 100% recycled paper

Item No.	Description	Pkg
NON243277	8 x 8.3"	90 sheets/bx, 30 bx/cs
NON243276	6 x 8"	60 sheets/bx, 72 bx/cs
NON243275	5.7 x 7"	40 sheets/bx, 200 bx/cs



NON26800

NON28800



NON26805



Green Tree® Basics

- » Made from recycled materials or 100% renewable sources
- » Soft and very affordable

Item No.	Description	Pkg
1-Ply		
NON26801	Standard Roll, 1,000 Sheets/Roll	96 rls/cs
2-Ply		
NON27800	Standard Roll, 100% Recycled, 500 Sheets/Roll	96 rls/cs
NON26800	Standard Roll, 4.5 x 3.8", 500 Sheets/Roll	96 rls/cs
NON26800H	Standard Roll, 4.5 x 3.8", 500 Sheets/Roll	1/ea
NON28800	Deluxe Roll, 4.5 x 3.25"	96 rls/cs
NON28805	Deluxe Roll, 9" Roll	12 rls/cs
NON26805	Jumbo Roll, 3.5" x 1,000'	8 rls/cs
NON28830	Small Core Toilet Paper, Also Fits Coreless Dispensers	36 rls/cs



Roll Towels

- » Strong enough to stand up to tough cleaning needs
- » Made of recycled pulp and rapidly renewable materials
- » Disposable, single-use towels

Item No.	Description	Pkg
NON25825	8" x 350', Natural	12 rl/cs
NON26830	8.25 x 12", Center-Pull	600 sht/rl, 6 rl/cs
NON26871	10" x 425', White	6 rl/cs
NONPBM350	8" x 350', White	12 rl/cs
NONPBM800B	8" x 800', White	6 rl/cs
NONPBM800N	8" x 800', Natural	6 rl/cs



When you choose our GreenSmart® paper line, you're helping the environment—one roll at a time! Even though they are gentle on the environment, based on their recycled content, Green Tree products are strong enough to stand up to tough cleaning needs.



NON26835



greensmart

Perforated Roll Towels

- » Made from rapidly renewable sources
- » Perforated

Item No.	Description	Pkg
NON26835	2-Ply, 8.8 x 11", White	85 sht/rl, 30 rl/cs



NON26820



greensmart

C-Fold Towels

- » Both sides of the towel are folded inwards
- » Stacked individually
- » 240 per pack, 10 packs per case

Item No.	Description	Pkg
NON26820	White	2,400/cs

*All greensmart facial tissues and towels have a minimum of 70% recycled content.



NON26810



greensmart

Multi-Fold Towels

- » Inter-folded for uniform dispensing
- » 9.125 x 9.5-inch sheet size
- » 250 per pack, 16 packs per case

Item No.	Description	Pkg
NON26810	White, Standard	4,000/cs
NON26810H	White, Standard	250/pk
NON26813	White, Deluxe	4,000/cs
NON26813H	White, Deluxe	250/pk
NON26818*	White, Premium	2,400/cs
NON25815	Natural, Deluxe	4,000/cs

*All greensmart facial tissues and towels have a minimum of 70% recycled content. NON26818 is not a greensmart product.



EVSDSMFB



Manual C-Fold/Multifold Towel Dispenser

- » Low-profile design
- » Touchless, smooth, one at a time presentation of C-fold or Multifold towels
- » Designed to eliminate over-stuffing of towels
- » Optional adhesive back is great for mounting on the mirror to replace towel baskets. Uses NON26810, NON26820, NON25815 or NON26813

Item No.	Description	Pkg
EVSDSMFB	Manual C-Fold/Multifold Towel Dispenser	1/ea



EVSDSTP92



Jumbo Bath Tissue Dispenser

- » Slim and compact
- » For use with NON26805

Item No.	Description	Pkg
EVSDSTP92	Jumbo Jr., Holds 2 Rolls of 9" Toilet Paper	1/ea
EVSDSTP91	Jumbo Jr. Single, Holds 1 Roll of 9" Toilet Paper	1/ea



GPC59447A



Georgia-Pacific

enMotion Impulse® 10 Automated Towel Dispenser

- » Designed to fit in small spaces
- » Helps reduce waste and the risk of cross contamination
- » Adjustable settings: sheet length, time delay and dispense mode

Item No.	Description	Pkg
GPC59447A	Automated Towel Dispenser	1/ea
GPC89460	High-Capacity Roll Towel	6/cs

Don't forget Hand Hygiene supplies.
Available on pages 74-76



GPC59488A



Georgia-Pacific

enMotion Automated Touchless Towel Dispenser

- » Highly reliable, affordable solution for automated touchless towel dispensing
- » Delivers paper towel in response to hand wave

Item No.	Description	Pkg
GPC59488A	Automated Touchless Towel Dispenser	1/ea



GPC89480



Georgia-Pacific

enMotion High Capacity Roll Towels

- » Specifically engineered to be used our enMotion wall-mount automated towel dispenser or enMotion Impulse® 10 towel dispenser

Item No.	Description	Pkg
GPC89480	High Capacity Roll Towels	6/cs



GPCGP27700



Georgia-Pacific

Preference® Jumbo Perforated Roll Towel

- » Ideal for cleaning and wiping in healthcare settings
- » Jumbo-sized 250-count roll offers the same full-weight, full-strength paper as standard, 2-ply preference

Item No.	Description	Pkg
GPCGP27700	Jumbo Perforated Roll Towel	12/cs



GPC20887



Georgia-Pacific

BigFold Z® Premium C-Fold Replacement Paper Towels

- » Premium towels fit into existing C-fold and multifold dispensers, with no adaptor required
- » Cost-saving C-fold substitute
- » Unique fold allows for reliable dispensing that reduces excess and the risk of cross contamination

Item No.	Description	Pkg
GPC20887	Replacement Paper Towels	10/cs



GPC20885



Georgia-Pacific

BigFold Z C-Fold Replacement Paper Towels

- » Unique fold allows for reliable, touchless dispensing of one fully-opened, absorbent towel
- » Interfold design cannot be loaded incorrectly
- » Fits into existing C-fold and multifold dispensers with the use of dispenser assistance device

Item No.	Description	Pkg
GPC20885	C-fold Paper Towels	10/cs



GPC33587



Georgia-Pacific

BigFold C-Fold Replacement Paper Towels

- » Meets EPA's federal procurement guidelines for post-consumer fiber and total recycled content
- » Fits into existing C-fold and combination C-fold multifold dispensers
- » No adaptor required

Item No.	Description	Pkg
GPC33587	C-fold Replacement Paper Towels	10/cs



GPC56797A



GPC19375



Compact® Side-By-Side Double Roll Bathroom Tissue Dispenser

- » Side-by-side coreless bathroom tissue dispensing system provides up to six times the capacity of single standard 2-ply rolls
- » Delivers 3,000 sheets of quality 2-ply tissue or 6,000 sheets of 1-ply tissue
- » Features a transfer paddle to prevent access to a new roll before the current roll is completely used up

Item No.	Description	Pkg
GPC56797A	Bathroom Tissue Dispenser	1/ea
GPC19375	Coreless 2-Ply Bathroom Tissue	36/cs
GPC1988001H	Envision® 2-Ply Embossed Bathroom Tissue	1/ea



GPC19378



Compact Coreless 2-Ply Bathroom Tissue

- » Reduces service calls to refill dispensers
- » High capacity

Item No.	Description	Pkg
GPC19378	1,000 Sheets 2-Ply Tissue Per Roll	18/cs



GPC16840



Angel Soft ps® 2-Ply Premium Embossed Bathroom Tissue

- » EPA compliant for recycled fiber content

Item No.	Description	Pkg
GPC16840	2-Ply Premium Embossed Bathroom Tissue	40/cs



GPC16620



Angel Soft ps Ultra 2-Ply Premium Bathroom Tissue

- » Perforated EZ Access cases to simplify storage and handling
- » Quality softness, two-ply
- » Designer embossing

Item No.	Description	Pkg
GPC16620	Ultra 2-Ply Premium Bathroom Tissue	20/cs



GPC48580HH



Angel Soft ps Premium Facial Tissue

- » Delivers the appearance and softness for today's discerning business customers and employees

Item No.	Description	Pkg
GPC48580HH	Flat Box, 100 Sheets/Box	1/bx



NON122424



Red Liners

- » Printed with the universal biohazard symbol and "Infectious Waste" in English and Spanish
- » Meet the 165-gram Dart Drop Test (ASTM-D1709-85) and 480-gram Elmendorf Tear Test (ASTM-D1922-67)

Item No.	Description	Pkg
NON151717H	17 x 17", 4-gal.	500/cs
NON122433	24 x 33", 15-gal.	250/cs
NON122424	24 x 24", 10-gal.	500/cs
NON123143	31 x 43", 33-gal.	200/cs



NONCXB36H



Black Low-Density Liners

- » Exceptional stretch- and tear-resistance
- » Ideal for heavier and sharper waste
- » Roll packs are easy to store and distribute, and are ideal for use on carts

Item No.	Description	Pkg
NONCXB36H	30 x 36", 20-30-gal.	250/ea
NONTGG58X	38 x 58", 60-gal.	100/ea
NONTGG39X	33 x 40", 33-gal.	250/ea



NONHCR24MC



High-Density Liners

- » Made of high-density polyethylene resins
- » Available in assorted sizes for step-on cans, linen hampers or garbage drums

Item No.	Description	Pkg
Black		
NONHCR24LB	24 x 24", 10-gal.	1,000/ea
Clear		
NONHCR24MC	24 x 24", 10-gal.	1,000/ea
NONHCR33MC	24 x 33", 15-gal.	1,000/ea
NONHCR40HC	33 x 40", 33-gal.	500/ea



Exam Table Paper

- » Strong, absorbent table paper protects exam tables from dirt and moisture
- » Disposable, single-use paper
- » All rolls are mechanically metered for exact length

Item No.	Description	Pkg
Deluxe Rolls		
NON24322	Smooth, 18" x 225'	12/cs
NON24323	Smooth, 20" x 225'	12/cs
NON24326	Smooth, 21" x 225'	12/cs
NON23328	Smooth, 8.5 x 225'	25/cs
NON24329	Crepe, 14.5" x 125'	12/cs
NON24324	Crepe, 18" x 125'	12/cs
NON24325	Crepe, 21" x 125'	12/cs



Item No.	Description	Pkg
Standard Rolls		
NON23319	Smooth, 14.5" x 225'	12/cs
NON23322	Smooth, 18" x 225'	12/cs
NON23323	Smooth, 20" x 225'	12/cs
NON23326	Smooth, 21" x 225'	12/cs
NON23324	Crepe, 18" x 125'	12/cs
NON23321	Crepe, 20" x 125'	12/cs
NON23325	Crepe, 21" x 125'	12/cs

Infant Scale Paper Sheets

NON253200CS	Crepe, White, 18 x 24"	1,000/cs
NON243200	Crepe, White, 20 x 30"	1000/cs




MDT219683



Medsoft Pillows

- » High-quality, low-cost pillows are soft, plump and fluid-proof
- » Soft vinyl ticking is flame-retardant, antistatic and fluid-proof
- » Vacuum packed
- » Please allow four hours for pillow to fully expand

Item No.	Description	Pkg
MDT219683	White, 20 x 26"	12/cs
MDT219683H	White, 20 x 26"	1/ea



PM2127-22



ComfortMed® Disposable Pillows

- » Single-use disposable pillow
- » Made of flame-resistant, garnetted polyester fill with white ticking and a blue hem thread
- » Filling is made from partially recycled material

Item No.	Description	Pkg
PM2127-22	Light-Weight, 16 x 22"	12/cs
PM1824-15	Medium-Weight, 18 x 24"	12/cs
PM1824-15H	Medium-Weight, 18 x 24"	1/ea
PM2127-18	Medium-Weight, 21 x 27"	12/cs
PM1622-22	Heavy-Weight, 21 x 27"	12/cs
PM2127-22H	Heavy-Weight, 21 x 27"	1/ea



NON24356B

Professional Disposable Towels

» Tissue/poly towels catch fluid and hold it, making them ideal for washing, wiping or draping

Item No.	Description	Pkg
2-Ply Tissue/Poly		
NON24356W	13 x 18", White	500/cs
NON24356B	13 x 18", Blue	500/cs
NON24361	17 x 19", White	500/cs

3-Ply Tissue and Tissue/Poly

NON24357W	Tissue, 13" x 18", White	500/cs
NON24357B	Tissue, 13" x 18", Blue	500/cs
NON24359	Tissue, 17" x 19", White	500/cs
NON24358W	Tissue/Poly, 13" x 18", White	500/cs
NON24358B	Tissue/Poly, 13" x 18", Blue	500/cs



NON24345

Pillowcases

Item No.	Description	Pkg
NON24345	Tissue/Poly, 21 x 30", White	100/cs
NON24346	Tissue/Poly, 21 x 30", Blue	100/cs
NON32300	Multilayer, 20 x 29", White	100/cs
NON32500	Multilayer, 20 x 29", Blue	100/cs
NON25300	Ultracel, Tissue/Poly, 21 x 30", White	100/cs



Interblend Percale Linens

- » Hydrogen peroxide bleached for whiter color and softer feel
- » Long wearing, wrinkle-resistant
- » Intimate blend of 55% cotton/ polyester percale fabric

Item No.	Description	Pkg
Flat Sheets		
MDTFS4P04N	60 x 104"	5 dz/cs
MDTFS4P04NZ	60 x 104"	2 dz/cs
MDTFS4P08N	60 x 108"	5 dz/cs
MDTFS4P15N	60 x 115"	5 dz/cs
MDTFS4P04	66 x 104"	5 dz/cs
MDTFS4P04Z	66 x 104"	2 dz/cs
MDTFS4P08	66 x 108"	5 dz/cs
MDTFS4P15	66 x 115"	5 dz/cs



Item No.	Description	Pkg
Contour Sheets		
MDTCS4P80	36 x 80 x 9"	5 dz/cs
MDTCS4P80Z	36 x 80 x 9"	2 dz/cs
MDTCS4P84	36 x 84 x 9"	5 dz/cs

Draw Sheets

MDTDS4P72	54 x 72"	5 dz/cs
MDTDS4P72Z	54 x 72"	2 dz/cs
MDTDS4P81	54 x 81"	5 dz/cs

Pillowcases

MDTPC4P34	42 x 34"	12 dz/cs
MDTPC4P34Z	42 x 34"	2 dz/cs



Classic Muslin Linens

- » Hydrogen peroxide bleached for whiter color and softer feel
- » Durable blend of 55% cotton/ 45% polyester muslin fabric

Item No.	Description	Pkg
Flat Sheets		
MDTFS4M04N	60 x 104"	5 dz/cs
MDTFS4M04NZ	60 x 104"	2 dz/cs
MDTFS4M04	66 x 104"	5 dz/cs
MDTFS4M04Z	66 x 104"	2 dz/cs
MDTFS4M08	66 x 108"	5 dz/cs
MDTFS4M15	66 x 115"	5 dz/cs

Item No.	Description	Pkg
Contour Sheets		
MDTCS4M80	36 x 80 x 9"	5 dz/cs
MDTCS4M80Z	36 x 80 x 9"	2 dz/cs
MDTCS4M84	36 x 84 x 9"	5 dz/cs

Draw Sheets

MDTDS4M72	54 x 72"	5 dz/cs
MDTDS4M72Z	54 x 72"	2 dz/cs
MDTDS4M81	54 x 81"	5 dz/cs

Pillowcases

MDTPC4M34	42 x 34"	12 dz/cs
MDTPC4M34Z	42 x 34"	2 dz/cs



Disposable Blanket

Item No.	Description	Pkg
Polyester/Cellulose		
NONDB4080	40 x 80", Gray	10/cs



NON24339C



Drape Sheets

Item No.	Description	Pkg
2-Ply Tissue		
NON23339	40 x 48", Deluxe, White	100/cs
NON24339M	40 x 48", Mauve	100/cs
NON24339A	40 x 60", White	100/cs
NON24339B	40 x 72", White	50/cs
3-Ply Tissue		
NON24336	40 x 48", White	100/cs
NON24336A	40 x 60", White	100/cs
NON24339C	40 x 72", White	50/cs
NON24336C	40 x 90", White	50/cs
Tissue/Poly		
NON24340	40 x 48", Blue	100/cs
NON24338	40 x 60", Blue	100/cs



NON24330BLU



Flat Sheets and Draw Sheets

Item No.	Description	Pkg
Flat Sheet, Tissue/Poly		
NON24330A	60 x 96", White	25/cs
NON24330B	58 x 102", White	25/cs
NON24330BLU	58 x 96", Blue	25/cs
Disposable Draw Sheets, Rayon/Poly		
NON32100	32 x 48", White	50/cs
NON32200	40 x 60", White	50/cs



Tissue/Poly Sheet,
various sizes available



Stretcher Sheets

Item No.	Description	Pkg
Flat, Tissue/Poly		
NON24333	40 x 72", Blue	50/cs
NON24335	40 x 90", Blue	50/cs
NON24337	46 x 96", White	50/cs
NON34700	40 x 90", White	50/cs
Fitted, Polypropylene		
NON34000	32 x 72", Blue	50/cs
NON34600	40 x 80", Blue	50/cs
Fitted, Reusable		
MDTNC4J10BLU	Jersey, Blue, 55% Cotton/45% Polyester, 24 x 72 x 12", 10-oz.	5 dz/cs
MDTNC4J10WHT	Jersey, White, 55% Cotton/45% Polyester, 24 x 72 x 12", 10-oz.	5 dz/cs
Flat, Multilayer		
NON33025	36 x 60", Dark Blue	50/cs
NON33100	40 x 84", Dark Blue	50/cs
NON33100H	40 x 84", Dark Blue	1/ea
Set		
NON37250	Top Sheet (NON33100), Fitted Sheet (NON34000), Pillowcase (NON32300)	24/cs



NON24398



Equipment Covers

Item No.	Description	Pkg
Polypropylene		
NON24398	Armboard Cover with Strap, 13 x 28.5", Blue	100/cs
NON24398H	Armboard Cover with Strap, 13 x 28.5", Blue	1/ea
NON24398NV	Armboard Cover, No Strap, 13 x 28.5", Blue	100/cs
NON24399	Headboard Cover, 14.5 x 15.25", Blue	100/cs

CRI4000



Classic Cover Gowns

- » Breathable polypropylene, lightweight, gowns are generously sized for optimum comfort, flexibility and mobility
- » Classic neck and side ties
- » Elastic wrists for secure fit

Item No.	Description	Pkg
----------	-------------	-----

Neck and Side Tie, Elastic Wrists

CRI4000	Regular/Large, Yellow	10/bg, 50/cs
CRI4000Z	Regular/Large, Yellow	10/pkg
CRI4000B	Regular/Large, Blue	10/bg, 50/cs
CRI4001	X-Large, Yellow	10/bg, 50/cs
CRI4001G	X-Large, Green	10/bg, 50/cs
CRI4003	3X-Large, Yellow	10/bg, 50/cs

Over-the-Head Neck, Dual-Side Tie, Thumb Loop Wrists

CRI4200	Regular/Large, Yellow	50/cs
CRI4201	X-Large, Yellow	50/cs

CRI5000



Standard Polyethylene Gowns

- » Blue polyethylene material
- » Over-the-head, open back and thumb loops

Item No.	Description	Pkg
----------	-------------	-----

CRI5000	Regular/Large, Blue	15/bx, 75/cs
CRI5000Z	Regular/Large, Blue	15/bx
CRI5001	X-Large, Blue	15/bx, 75/cs
NONTH150	Regular/Large, Blue	75/cs
NONTH180	Universal, Blue	75/cs
NONTH200	X-Large, Blue	75/cs

NON27SMS2



Medium Weight Multi-Ply Fluid Resistant Isolation Gowns

- » Side ties help eliminate the difficulties of tying behind the back
- » Neck tie style with elastic wrists

Item No.	Description	Pkg
----------	-------------	-----

Yellow

NON27SMS2	Regular/Large	10/bg, 100/cs
NON27SMS2XL	X-Large	10/bg, 100/cs
NON27SMS2XXXL	3X-Large	10/bg, 100/cs



NON30233B



Multilayer Bouffant Caps

Item No.	Description	Pkg
----------	-------------	-----

Medium-Weight Multi-Ply

NON30232B	21", Blue	100/bx, 500/cs
NON30233B	24", Blue	100/bx, 500/cs
NON30223B	28", Blue	100/bx, 500/cs

Heavy-Weight Multi-Ply

NON28533B	24", Blue	100/bx, 500/cs
NON28234P	24", Multicolor Print	100/bx, 500/cs

CRI1000



CRI1002



Boundary Bouffant Caps

Item No.	Description	Pkg
----------	-------------	-----

Lightweight Polypropylene

CRI1000	21", Green	100/bx, 500/cs
CRI1001	21", Blue	100/bx, 500/cs
CRI1002	21", White	100/bx, 500/cs
CRI1003	24", White	100/bx, 500/cs
CRI1004	24", Blue	100/bx, 500/cs

Polypropylene

CRI1202	24", White	500/cs
CRI1203	24", White	100/bx, 500/cs



NON28625Z



Disposable Surgeon's Cap

- » Feature a lightweight single layer crown material allowing for breathability
- » Less linting
- » Not made with natural rubber latex
- » One size fits most

Item No.	Description	Pkg
NON28625Z	Multilayer, Tie Back, Blue	100/bx
NON28626Z	Multilayer, Tie Back, Blue	100/bx
NON28628Z	Spunlace, Tie Back, Blue	144/cs



NON28858



Multilayer Shoe Covers

- » Made of multilayer material
- » Features non-skid for increased traction
- » Range of sizes to fit all shoe types

Item No.	Description	Pkg
----------	-------------	-----

Heavyweight

NON29858	Regular/Large	300/cs
NON29859	X-Large	200/cs

Mediumweight

NON28858	Regular/Large	300/cs
NON28859	X-Large	200/cs
NON28852	Sports Size	200/cs

Lightweight

CRI2006	Regular/Large	300/cs
CRI2007	X-Large	200/cs



CRI2002

CRI2000



Spunbond Polypropylene Shoe Covers

- » Choice of non-skid pattern or smooth bottom

Item No.	Description	Pkg
----------	-------------	-----

Premium Non-Skid Tread

NON28758	Regular, Up to Men's Size 12	100/bx, 300/cs
NON28759	X-Large, Up to Men's Size 15	100/bx, 200/cs
NON28752	Sport, Fits Athletic Shoes	100/bx, 200/cs

Non-Skid Tread

CRI2002	Regular/Large	300/cs
CRI2003	X-Large	200/cs

Spunbond Polypropylene

CRI2000	Regular/Large	200/cs
CRI2001	X-Large	200/cs



NON24775



Safety Glasses & Fluid Protective Goggles

- » Clear lenses available for optical clarity
- » Fit comfortably and securely
- » Larger size glasses can be worn over prescription eyewear

Item No.	Description	Pkg
----------	-------------	-----

Goggles

NON24775	Deluxe Fluid Protection Lab Goggles, Fog-Free	36/bx
NON24776	Standard Fluid Protection Lab Goggles	36/bx

Safety Glasses

NON24777V	Visitor Safety Glasses	144/cs
NON24774	Safety Glasses with Black Frame	12/bx
NON24770	Safety Glasses, Clear Frame	144/cs



Polyethylene Shoe Covers

- » Smooth bottom

Item No.	Description	Pkg
----------	-------------	-----

CRI2010	One Size Fits Most	1,000/cs
---------	--------------------	----------



NONFS300



Disposable Eyewear and Face Shields

- » Assorted color reusable frames to go with our clear disposable lenses
- » Resistant to fogging
- » Optically clear polyester for maximum visibility

Item No.	Description	Pkg
----------	-------------	-----

NONLENS	Disposable Lens for Use with NONFRAME	250/bx
NONFRAME	Frame for Use with NONLENS, Assorted Colors	100/bx
NONASMBLD	Assembled Frames and Lenses, Assorted Colors	50/cs
NONFS100	Extra-Long Face Shield, 9.5" Length	40/cs
NONFS300	Full-Length Face Shield	96/cs
NONFS400	¾-Length Face Shield	96/cs



Fluid-Protection Face Masks

- » Choose from three levels of fluid resistance: 80 mmHg (level 1), 120 mmHg (level 2) and 160 mmHg (level 3)
- » Anti-fog strip prevents eyeglasses from fogging
- » Facemasks with an eye shield feature an anti-glare strip
- » Not made with natural rubber latex

Item No.	Fluid Resistance Level	Shield	Ties/Earloops	Pkg
----------	------------------------	--------	---------------	-----

Breathable Film Mask

NON27810	3	Y	Ties	25/bx, 250/cs
NON27810EL	3	N	Earloops	25/bx, 250/cs
NON27812EL	3	N	Earloops	50/bx, 500/cs

Surgical Face Mask

NON27710	3	Y	Ties	25/bx, 100/cs
NON27712	3	N	Ties	50/bx, 300/cs
NON27410	2	Y	Ties	25/bx, 100/cs
NON27412	2	N	Ties	50/bx, 300/cs
NON27405	1	Y	Ties	25/bx, 100/cs
NON27408	1	N	Ties	50/bx, 300/cs

Procedure Face Mask

NON27710EL	3	Y	Earloops	25/bx, 100/cs
NON27712EL	3	N	Earloops	50/bx, 300/cs
NON27410EL	2	Y	Earloops	25/bx, 100/cs
NON27412EL	2	N	Earloops	50/bx, 300/cs
NON27408EL	1	N	Earloops	50/bx, 300/cs

Chamber-Style Face Mask

NON27411AG	2	Y	Ties	25/bx, 100/cs
NON27382FR	2	N	Ties	50/bx, 300/cs



NON27812EL



Face Masks

- » Light and comfortable
- » A full-width nosepiece guarantees proper fit
- » Elastic earloops for donning and positioning
- » Not made with natural rubber latex

Item No.	Description	Pkg
----------	-------------	-----

Procedure Face Mask, Blue

NON27375	Cellulose Inner, Spunbond Polypropylene Outer Facing	50/bx, 300/cs
NON27378	Spunbond Polypropylene Inner and Outer Facings	50/bx, 300/cs

Isolation Face Mask, Yellow

NON27120	Cellulose Inner, Spunbond Polypropylene Outer Facing	50/bx, 300/cs
NON27122	Spunbond Polypropylene Inner and Outer Facings	50/bx, 300/cs



NON27110



NON24244



NON27046



NON27346SL



NON24247



NON27148



Item No.	Description	Color	Pkg
----------	-------------	-------	-----

Patient Gowns, Waist Tie, Tissue/Poly/Tissue

NON24356	Deluxe Gown, 30 x 42"	Blue	50/cs
NON24355	Deluxe Gown, 30 x 42"	White	50/cs
NON24354	Deluxe X-Ray Gown, T/P, 30 x 42"	Blue	50/cs

Patient Gowns, Multilayer Material

NON27046	Sleeveless, Waist Tie, Regular/Large	Blue	50/cs
NON27048	Sleeveless, Waist Tie, Regular/Large	Pink	50/cs
NON27146	Sleeveless, Neck and Waist Tie, Regular/Large	Blue	50/cs
NON27146SL	Short Sleeves, Neck and Waist Tie, Regular/Large	Blue	50/cs
NON27146SLXL	Short Sleeves, Neck and Waist Tie, X-Large	Blue	50/cs
NON27346SL	Three-Armhole Patient Gown, Sleeveless, Regular/Large	Blue	50/cs

Scrim Patient Gowns, Waist Tie

NON254572	Bariatric Patient Gown, X-Large, 72 x 45"	Blue	25/cs
NON24247	Patient Gown, 30 x 42"	Blue	50/cs

Item No.	Description	Color	Pkg
----------	-------------	-------	-----

Patient Capes, Tissue/Poly/Tissue

NON25249	Deluxe Cape, 30 x 21"	Blue	100/cs
NON24248	Deluxe Cape, 30 x 21"	White	100/cs

Patient Capes, Sleeveless, Tape Tab Closure, Multilayer Material

NON27036	30 x 21"	Blue	100/cs
NON27038	30 x 21"	Pink	100/cs

Disposable Exam Shorts, Elastic Waist

NON27209*	Available Small-XXL	Blue	30/cs
-----------	---------------------	------	-------

Disposable Robes and Telemetry Pouch, Multilayer Material

NON27148	Patient Robe, Regular/Large	Blue	25/cs
NON27148XL	Patient Robe, X-Large	Blue	25/cs
NON249005	Telemetry Pouch, Disposable	Blue	100/cs

*Add size code S (Small), M (Medium), L (Large), XL (X-Large), XXL (2X-Large).

Looking for additional apparel options?
Visit us on medline.com/physician.



Adult IV Gowns

- » Comfort, style and modesty
- » Contemporary-designed unisex IV gowns
- » Soft 100% polyester fabric
- » MR-safe plastic snaps
- » One-size-fits-most; 64" sweep, 47" length, 5" in-sleeve

Item No.	Description	Pkg
MDTPG6IVRSB	Sage and Beige	6 dz/cs
MDTPG6IVRRB	Royal Blue	6 dz/cs
MDTPG6IVCAG	Cascade, Green	6 dz/cs
MDTPG6IVCAB	Cascade, Blue	6 dz/cs


MDTPG6IVCAB



Sage and Beige



Royal Blue



Cascade, Green



Cascade, Blue



Unisex Teen Gowns

- » Let your patient's personality shine through
- » Designed in partnership with a young patients in our fresh, fun patterns are age appropriate and help make your patients more comfortable and relaxed
- » One-size-fits-most designed for smaller patients; 60" sweep, 44" length, 4.5" in-sleeve)
- » MR-safe plastic snaps means no need to change gowns for x-rays

Item No.	Description	Pkg
MDTPG6IVCHEV	Chevron, Blue	6 dz/cs
MDTPG6IVTIDI	Tie-dye, Blue/green	6 dz/cs
MDTPG6IVPOP	Pop, Multi-Color	6 dz/cs

MDTPG6IVPOP



Chevron, Blue



Tie-dye, Blue/Green



Pop, Multi-Color



IV Gown

- » The sleeves have snaps for easy access to IV lines
- » Telemetry pocket placed slightly lower for modesty
- » Side ties and a generous sweep provide full back side coverage
- » MRI and laundry safe plastic snaps
- » Durable plastic snaps
- » 64" sweep; 47" length; 4.5" in-sleeve

Item No.	Description	Pkg
MDTPG7PTSPKT	Honeycomb Pattern	6 dz/cs



Modest Criss-Cross Nursing Gown

- » This nursing gown features a stretch knit criss-cross insert that provides modesty
- » Features durable, high quality plastic snaps that resist breakage when laundered and true tie side coverage
- » 68" sweep; 45" length; 7" in-sleeve

Item No.	Description	Pkg
MDTSG7ITSFLB	Blue Floral	6 dz/cs



Women's Front Open Gown

- » Ideal for mammography or physician offices
- » It ties in front, making it easy to dress and providing full back side coverage
- » 70" sweep; 44" length; 11.5" in-sleeve

Item No.	Description	Pkg
MDTPG7RF0FLG	Green Floral Print	6 dz/cs



MDTHR9R04SAG



Women's Wellness Gowns

- » Maternity services. Mammograms. Breast cancer treatment. Our new fashion forward gowns are designed specifically with women's comfort and confidence in mind
- » Modest cuts and extra coverage for enhanced dignity
- » Ties at side or the front for full back and side coverage
- » Softer colors and attractive styles increase patient satisfaction
- » Durable 85% poly/15% nylon fabric stands up to repeated washings

Item No.	Description	Pkg
MDTHR9R04SAG	Green Suede	6 dz/cs
MDTHR9R04PINK	Pink Suede	6 dz/cs
MDTHR9R04LIL	Lilac Suede	6 dz/cs



MDTHR9R04PINK



MDTHR9R04LIL



One-piece, loose back belt



Unisex Knee-Length Lab Coat


- » Performance poplin blend, 80% polyester/20% cotton (white only)
- » Navy and Ciel Blue colors are made of 65% polyester/35% cotton
- » Five-button front; left chest pocket with pencil slot; two lower pockets
- » Side vent openings
- » 41" length

Women's Classic Staff-Length Lab Coat

- » Made of 80% polyester/20% cotton poplin for easy care, resists wrinkles
- » White button-front with pencil divided left breast pocket, two lower pockets and princess-seamed back with one-piece backbelt
- » Side slashes for access to trouser pockets

Item No.	Size	Item No.	Size
White			
83044QHWWXS	XX-Small	83044RNNS	Small
83044QHWXS	X-Small	83044RNNM	Medium
83044QHWS	Small	83044RNNL	Large
83044QHWM	Medium	83044RNNXL	X-Large
83044QHWL	Large	83044RNNXXL	2X-Large
83044QHWXL	X-Large	83044RNNXXXL	3X-Large
83044QHWXXL	2X-Large		
83044QHWXXXL	3X-Large		
Gray			
83044GRYXS	X-Small	83044RCWXS	X-Small
83044GRYS	Small	83044RCWS	Small
83044GRYM	Medium	83044RCWM	Medium
83044GRYL	Large	83044RCWL	Large
83044GRYXL	X-Large	83044RCWXL	X-Large
83044GRYXXL	2X-Large	83044RCWXXL	2X-Large
83044GRYXXXL	3X-Large	83044RCWXXXL	3X-Large
Navy			
Light Blue			

Item No.	Dress Size	Recommended Categorical Size
MDT11WHT2E	2	X-Small
MDT11WHT4E	4	Small
MDT11WHT6E	6	
MDT11WHT8E	8	Medium
MDT11WHT10E	10	
MDT11WHT12E	12	
MDT11WHT14E	14	Large
MDT11WHT16E	16	
MDT11WHT18E	18	X-Large
MDT11WHT20E	20	
MDT11WHT22E	22	2X-Large
MDT11WHT24E	24	
MDT11WHT26E	26	3X-Large
MDT11WHT28E	28	
MDT11WHT30E	30	4X-Large
MDT11WHT32E	32	
MDT11WHT34E	34	5X-Large
MDT11WHT36E	36	



Looking for a coat with an iPad pocket?
Ask your Medline Representative for details and ordering information



MDT12WHT44E



Men's Classic Staff-Length Lab Coat

- » Made of 80% polyester/20% cotton poplin for easy care, resists wrinkles
- » White button-front with pencil divided left breast pocket, two lower pockets and princess-seamed back with one-piece backbelt
- » Side slashes for access to trouser pockets

Item No.	Chest Size	Recommended Categorical Size
MDT12WHT32E	32	X-Small
MDT12WHT34E	34	Small
MDT12WHT36E	36	
MDT12WHT38E	38	Medium
MDT12WHT40E	40	
MDT12WHT42E	42	Large
MDT12WHT44E	44	
MDT12WHT46E	46	X-Large
MDT12WHT48E	48	
MDT12WHT50E	50	2X-Large
MDT12WHT52E	52	
MDT12WHT54E	54	3X-Large
MDT12WHT56E	56	
MDT12WHT58E	58	4X-Large
MDT12WHT60E	60	



MDT10WHT44E



Unisex Consultation Coats

- » Three-button front coat with left chest and two lower pockets
- » Performance poplin blend, 80% polyester/20% cotton for easy care
- » 30" length

Item No.	Chest Size	Recommended Categorical Size
MDT10WHT32E	32	X-Small
MDT10WHT34E	34	
MDT10WHT36E	36	Small
MDT10WHT38E	38	
MDT10WHT40E	40	Medium
MDT10WHT42E	42	
MDT10WHT44E	44	Large
MDT10WHT46E	46	
MDT10WHT48E	48	X-Large
MDT10WHT50E	50	
MDT10WHT52E	52	2X-Large
MDT10WHT54E	54	
MDT10WHT56E	56	3X-Large
MDT10WHT58E	58	



88018QHW8



Women's Consultation Coats

- » Performance poplin blend, 80% polyester/20% cotton for easy care
- » Button-front with lapel collar
- » Left breast pocket and two lower pockets for ample, easy storage
- » Princess seams front and back for a tailored look

Item No.	Dress Size	Recommended Categorical Size
88018QHW4	4	X-Small
88018QHW6	6	
88018QHW8	8	Small
88018QHW10	10	
88018QHW12	12	Medium
88018QHW14	14	
88018QHW16	16	Large
88018QHW18	18	
88018QHW20	20	X-Large
88018QHW22	22	
88018QHW24	24	2X-Large
88018QHW26	26	
88018QHW28	28	3X-Large
88018QHW30	30	



ComfortEase™ Unisex Two Pocket Scrub Top

- » Slightly longer v-neck pullover tunic
- » Two lower pockets for easy storage and side vents for easy movement
- » Consider ordering one size up for comfort and fit



Item No.	Description	Pkg
8800 STYLE*	Scrub Top	1/ea

*See below to select your size and color.



ComfortEase™ Unisex Drawstring Pants

- » Braided drawstring makes it easy to tie and untie
- » Useful back pocket and deluxe four-panel construction with sewn in gussets on sizes L and up for greater comfort
- » Stain-resistant finish
- » Regular 31-1/2" inseam



Item No.	Description	Pkg
400 STYLE*	Drawstring Pants	1/ea

*See below to select your size and color.



ComfortEase™ Unisex Crew Neck Warm-Up Jacket

- » Classic styling with self-colored snaps, knit cuffs and two convenient lower pockets
- » Soil release finish



Item No.	Description	Pkg
8832 STYLE*	Warm-up Jacket	1/ea

*See below to select your size and color.

ComfortEase

Sizes: XXS-7XL

WNE	Wine	JRL	Royal
JEG	Evergreen	JNT	Midnight Blue
JCB	Caribbean Blue	XTQ	White
JTH	Ceil Blue	DKW	Black

For our complete color selection, contact your Medline Representative.

How to order

- » Pick your style: 8800, 400, 8832
- » Choose your color
- » Select your size - XXS - 7XL

Style Number + Color + Size



Women's Berkeley Tops

- » Designed to keep shape and color over time
- » Roomy two lower welt pockets and chest pocket with hook-and-loop closures



5582 STYLE

Item No.	Description	Pkg
5582 STYLE*	Berkeley Tops	1/ea

*See below to select your size and color.



Women's Pacific Pants

- » Flattering design curved through the hip and calf to create a long leg look
- » Wide elastic waistband with inner drawstring for fit and comfort
- » Cargo pockets provide ample room



5570 STYLE

Item No.	Description	Pkg
5570 STYLE*	Pacific Pants	1/ea

*See below to select your size and color.



Unisex Lincoln Jacket

- » Versatile style with knit cuffs and hidden snaps
- » Three patch pockets for convenient storage



5542 STYLE

Item No.	Description	Pkg
5542 STYLE*	Lincoln Jacket	1/ea

*See below to select your size and color.

AVE

Sizes: XXS-3XL

RED	Red	GRY	Light Gray
WNE	Wine	CHR	Charcoal
HTR	Hunter	OLV	Olive
CRB	Caribbean	KHK	Khaki
CBL	Ceil	CHC	Chocolate
RYL	Royal	BLK	Black
NVY	Navy		

For our complete color selection, contact your Medline Representative.



MRXXM1405H

MAP101400H



Metrex®

» These sterilizing and disinfecting solutions and supplies can be used to disinfect surfaces across all the critical touch points in your facility

Item No.	Description	Pkg
----------	-------------	-----

MetriCide® OPA Plus

MAP106000	1-gal.	4/cs
MAP106000H	1-gal.	1/ea

MetriCide 28, 2.5%, 28 Day

MAP102800	1-gal.	4/cs
MAP102800H	1-gal.	1/ea
MAP102805	1-qt./0.9 L	16/cs
MAP102805H	1-qt./0.9 L	1/ea

MetriCide 14, 2.6%, 14 Day

MAP101400H	1-gal.	1/ea
MRXXM1405H	1-qt./0.9 L	1/ea

CaviWipes®, Classic, Disinfectant

MAP131100Z	Small, 6 x 6.75"	160/bx
MAP131155Z	XL, 9 x 12"	50/bx
MAP131150	XL, 9 x 12"	12/cs

Strip Glutaraldehyde, MetriTest™

MAP10303H	1.5%	60/ea
MAP10304H	1.8%	60/ea

Disinfectant, Enviroid®

MAP133300H	1-gal.	1/btl
------------	--------	-------

CaviCide®

MAP131000H	Disinfectant, 1-gal.	1/ea
MAP131024CS	Spray Bottle: 24-oz.	12/cs



J-J20390H



Sterilants and Test Strips

- » Provides exceptional sanitation in medical environments
- » Cidex® OPA Solution is glutaraldehyde-free and effective at room temperature (68°F)
- » Enzol® is a mild, low-foaming enzymatic detergent
- » Cidex Plus has a built-in corrosion/rust inhibitor

Item No.	Description	Pkg
----------	-------------	-----

J-J20390H	Cidex OPA, 1-gal.	1/ea
J-J20392H	Cidex OPA Test Strips	60/btl
J-J2252HH	Enzol Enzyme Detergent, 1-gal.	1/ea
J-J2920CS	Cidex Test Strips, Active Dialdehyde, 60 Bottle	2/btl
MDS200360	Cidex/Metricide OPA Test Strips	100/cs



SDCML020137



SDCML020127

Rapicide® OPA/28 High-Level Disinfectant

- » 28 day reuse period and ten minute manual disinfection time

Item No.	Description	Pkg
----------	-------------	-----

SDCML020127	Rapicide OPA/28 HLD Solution	4 gal./cs
SDCML020137	Rapicide OPA/28 Test Strips	100/cs



Micro-Kill™ Bleach Germicidal Bleach Wipes

- » These hard surface disinfectant wipes have a clean scent
- » For use on countertops, tables and other equipment that requires the use of bleach for disinfection
- » Meets CMS Tag F441 guidelines for a 1:10 bleach dilution
- » EPA Reg. No 37549-1

Item No.	Description	Pkg
MSC351400ANZ	Wipes, 150-Count Canister: 7 x 8"	1/ea
MSC351410AN	Wipes, 150-Count Canister: 5 x 6"	6/cs
MSC351420ANZ	Wipes, Individually Wrapped: 3 x 3"	100/box
MSC351440AN	Wipes, 55-Count Canister: 12 x 12"	8/cs

Microorganism	Kill Time
<i>Acinetobacter baumannii</i> -MDR	30 seconds
Adenovirus Type 5	30 seconds
Avian Influenza A (H3N2)	30 seconds
Avian Influenza A (H5N1)	30 seconds
Avian Influenza A (H7N9)	30 seconds
<i>Bordetella pertussis</i>	30 seconds
CA-Methicillin Resistant <i>Staphylococcus aureus</i>	30 seconds
<i>Campylobacter jejuni</i>	30 seconds
Canine Parvovirus	30 seconds
<i>Enterobacter aerogenes</i>	30 seconds
<i>Enterobacter cloacae</i>	30 seconds
<i>Enterococcus faecalis</i> -VRE	30 seconds
<i>Enterococcus faecium</i> -MDR	30 seconds
<i>Escherichia coli</i>	30 seconds
<i>Escherichia coli</i> O157:H7	30 seconds
<i>Escherichia coli</i> -ESBL	30 seconds
<i>Escherichia coli</i> -NDM-1	30 seconds
Feline Panleukopenia Virus	30 seconds
Hepatitis A Virus	30 seconds
Hepatitis B Virus	30 seconds
Hepatitis C Virus	30 seconds
Herpes Simplex Virus Type 2	30 seconds
Human Coronavirus	30 seconds
Human Immunodeficiency Virus Type 1	30 seconds
Influenza A Virus (Hong Kong)	30 seconds
Influenza A Virus (Novel H1N1)	30 seconds

Microorganism	Kill Time
<i>Klebsiella pneumoniae</i>	30 seconds
<i>Klebsiella pneumoniae</i> -CRKP (CRE)	30 seconds
<i>Klebsiella pneumoniae</i> -ESBL	30 seconds
<i>Klebsiella pneumoniae</i> -NDM-1	30 seconds
<i>Legionella pneumophila</i>	30 seconds
Linezolid Resistant <i>Staphylococcus aureus</i>	30 seconds
<i>Listeria monocytogenes</i>	30 seconds
Methicillin Resistant <i>Staphylococcus aureus</i> (MRSA)	30 seconds
Norovirus (Feline Calicivirus Virus)	30 seconds
Polio Virus Type 1	30 seconds
<i>Proteus mirabilis</i>	30 seconds
<i>Pseudomonas aeruginosa</i>	30 seconds
Respiratory Syncytial Virus	30 seconds
Rhinovirus Type 37	30 seconds
Rotavirus	30 seconds
<i>Salmonella enterica</i>	30 seconds
<i>Serratia marcescens</i>	30 seconds
<i>Shigella dysenteriae</i>	30 seconds
<i>Staphylococcus aureus</i>	30 seconds
<i>Staphylococcus aureus</i> -VISA	30 seconds
<i>Staphylococcus aureus</i> -VRSA	30 seconds
<i>Streptococcus pneumoniae</i>	30 seconds
<i>Streptococcus pyogenes</i>	30 seconds
<i>Candida albicans</i>	1 min
<i>Trichophyton mentagrophytes</i>	1 min
<i>Clostridium difficile</i> spores	3 min

Kill (i.e., contact) time for a disinfectant is the amount of time a surface must remain wet with the product to achieve disinfection.



Micro-Kill Bleach Germicidal Bleach Solution

- » Kills C. diff Spores in Five Minutes*
- » Micro-Kill Bleach is a disinfecting sporicidal solution that meets CDC, CMS Tag F441, SHEA, IDSA, APIC and Joint Commission guidelines for killing *Clostridium difficile* spores²⁻⁵
- » Common uses of Micro-Kill Bleach solution disinfects hard, non-porous surfaces, including countertops, tables, carts and equipment
- » Bleach solution equivalent to 1:10 bleach dilution
- » Kills 52 microorganisms
- » EPA Reg. No. 37549-2

Item No.	Description	Pkg
EVSCH900**	Liquid Bleach, 1-gal.	4/cs
EVSCH910**	Liquid Bleach, 32-oz.	12/cs

*Ensure that the surface remains wet for the entire contact time.

**May not be available in all 50 states.

Microorganism	Kill Time
<i>Acinetobacter baumannii</i> -MDR	30 seconds
Adenovirus Type 5	30 seconds
Avian Influenza A (H3N2)	30 seconds
Avian Influenza A (H5N1)	30 seconds
Avian Influenza A (H7N9)	30 seconds
<i>Bordetella pertussis</i>	30 seconds
CA-Methicillin Resistant <i>Staphylococcus aureus</i>	30 seconds
<i>Campylobacter jejuni</i>	30 seconds
Canine Parvovirus	30 seconds
<i>Enterobacter aerogenes</i>	30 seconds
<i>Enterobacter cloacae</i>	30 seconds
<i>Enterococcus faecalis</i> -VRE	30 seconds
<i>Enterococcus faecium</i> -MDR	30 seconds
<i>Escherichia coli</i>	30 seconds
<i>Escherichia coli</i> O157:H7	30 seconds
<i>Escherichia coli</i> -ESBL	30 seconds
<i>Escherichia coli</i> -NDM-1	30 seconds
Feline Panleukopenia Virus	30 seconds
Hepatitis A Virus	30 seconds
Hepatitis B Virus	30 seconds
Hepatitis C Virus	30 seconds
Herpes Simplex Virus Type 2	30 seconds
Human Coronavirus	30 seconds
Human Immunodeficiency Virus Type 1	30 seconds
Influenza A Virus (Hong Kong)	30 seconds
Influenza A Virus (Novel H1N1)	30 seconds

Microorganism	Kill Time
<i>Klebsiella pneumoniae</i>	30 seconds
<i>Klebsiella pneumoniae</i> -CRKP (CRE)	30 seconds
<i>Klebsiella pneumoniae</i> -ESBL	30 seconds
<i>Klebsiella pneumoniae</i> -NDM-1	30 seconds
<i>Legionella pneumophila</i>	30 seconds
Linezolid Resistant <i>Staphylococcus aureus</i>	30 seconds
<i>Listeria monocytogenes</i>	30 seconds
Methicillin Resistant <i>Staphylococcus aureus</i> (MRSA)	30 seconds
Norovirus (Feline Calicivirus Virus)	30 seconds
Polio Virus Type 1	30 seconds
<i>Proteus mirabilis</i>	30 seconds
<i>Pseudomonas aeruginosa</i>	30 seconds
Respiratory Syncytial Virus	30 seconds
Rhinovirus Type 37	30 seconds
Rotavirus	30 seconds
<i>Salmonella enterica</i>	30 seconds
<i>Serratia marcescens</i>	30 seconds
<i>Shigella dysenteriae</i>	30 seconds
<i>Staphylococcus aureus</i>	30 seconds
<i>Staphylococcus aureus</i> -VISA	30 seconds
<i>Staphylococcus aureus</i> -VRSA	30 seconds
<i>Streptococcus pneumoniae</i>	30 seconds
<i>Streptococcus pyogenes</i>	30 seconds
<i>Candida albicans</i>	1 min
<i>Trichophyton mentagrophytes</i>	1 min
<i>Clostridium difficile</i> spores	5 min

Kill (i.e., contact) time for a disinfectant is the amount of time a surface must remain wet with the product to achieve disinfection.



Micro-Kill AF²

- » No bleach or alcohol added
- » 2-minute kill time*
- » Kills 47 microorganisms
- » Quat-based disinfectant for hard non-porous surfaces
- » EPA Reg. No. 6836-372-37549

Item No.	Wipe Dimensions	Pkg
MSC351000**	6 x 6.75" Wipes, 160 Count Canister	12 ea/cs
MSC351010**	7 x 15" Wipes, 65 Count Canister	6 ea/cs

*Ensure that the surface remains wet for the entire contact time.
 **May not be available in all 50 states.



MSC351000

MSC351010



Microorganism	Kill Time
<i>Acinetobacter baumannii</i>	2 minutes
<i>Acinetobacter baumannii</i> Multi Drug Resistant (MDR)	2 minutes
Avian Influenza A Virus H3N2	2 minutes
Avian Influenza A Virus H5N1	2 minutes
<i>Burkholderia cepacia</i>	2 minutes
<i>Campylobacter jejuni</i>	2 minutes
<i>Candida albicans</i>	2 minutes
Community Associated, Methicillin-Resistant <i>Staphylococcus aureus</i> (CA-MRSA) (USA 300)	2 minutes
Community Associated, Methicillin-Resistant <i>Staphylococcus aureus</i> (CA-MRSA) (USA 400)	2 minutes
<i>Enterobacter aerogenes</i>	2 minutes
<i>Enterococcus faecalis</i>	2 minutes
<i>Escherichia coli</i> (E.coli)	2 minutes
<i>Escherichia coli</i> O157:H7	2 minutes
<i>Escherichia coli</i> NDM-1 Positive	2 minutes
ESBL <i>Escherichia coli</i>	2 minutes
ESBL <i>Klebsiella pneumoniae</i>	2 minutes
Feline calicivirus	2 minutes
Hepatitis B Virus (HBV)	2 minutes
Hepatitis C Virus (HCV)	2 minutes
Herpes Simplex Virus Type 1	2 minutes
Herpes Simplex Virus Type 2	2 minutes
HIV-1	2 minutes
Human Coronavirus	2 minutes
Human Influenza A Virus (H7N9)	2 minutes

Microorganism	Kill Time
Influenza A Virus/Hong Kong (H3N2)	2 minutes
<i>Klebsiella pneumoniae</i>	2 minutes
<i>Klebsiella pneumoniae</i> NDM-1 Positive	2 minutes
<i>Listeria monocytogenes</i>	2 minutes
Methicillin-Resistant <i>Staphylococcus aureus</i> (MRSA)	2 minutes
Multi Drug-Resistant <i>Staphylococcus aureus</i>	2 minutes
Norwalk Virus (Norovirus)	2 minutes
Penicillin Resistant <i>Streptococcus pneumoniae</i>	2 minutes
<i>Proteus vulgaris</i>	2 minutes
<i>Pseudomonas aeruginosa</i>	2 minutes
Respiratory Syncytial Virus	2 minutes
Rotavirus	2 minutes
<i>Salmonella enterica</i>	2 minutes
<i>Salmonella enterica</i> serovar Typhimurium	2 minutes
SARS Associated Coronavirus (SARS)	2 minutes
<i>Serratia marcescens</i>	2 minutes
<i>Shigella dysenteriae</i>	2 minutes
<i>Staphylococcus aureus</i>	2 minutes
<i>Streptococcus pyogenes</i> associated with scarlet fever	2 minutes
Vancomycin-Resistant <i>Streptococcus pyogenes</i> (VRSA)	2 minutes
Vancomycin Intermediate-Resistant <i>Staphylococcus aureus</i> (VISA)	2 minutes
Vancomycin-Resistant <i>Enterococcus faecalis</i> (VRE)	2 minutes
Vancomycin-Resistant <i>Staphylococcus aureus</i> (VRSA)	2 minutes
<i>Vibrio cholerae</i>	2 minutes

Kill (i.e., contact) time for a disinfectant is the amount of time a surface must remain wet with the product to achieve disinfection.



MSC35120

MSC351200

Micro-Kill+[†] Disinfectant, Deodorizing Cleaning Wipes with Alcohol

- » This one-step disinfectant cleaner for hard non-porous surfaces kills 12 different infectious microorganisms, within two minutes**
- » Heavy-duty fabric withstands hard-surface scrubbing and contains enough liquid to ensure maximum surface coverage
- » EPA Reg. No. 59894-10-37549

Item No.	Description	Pkg
MSC351200	6 x 6.75" Wipes, 160 Count Canister	12 ea/cs
MSC351230	7 x 10" Wipes, 50 Count Narrow Canister	12 ea/cs
MSC351210	10 x 10" Wipes, 65 Count Canister	12 ea/cs
MSC351240	10 x 14" Wipes, 42 Count Canister	12 ea/cs

*Medline Micro-Kill+ Disinfectant, Deodorizing Cleaning Wipes with Alcohol is not recommended for use on porous surfaces or lacquered surfaces as it might remove the finish. Spot-testing on an inconspicuous location is recommended before use on plastic or painted surfaces **Ensure that the surface remains wet for the entire contact time. †ADA compliant when wall mounted.

Microorganism	Kill Time
HIV-1 (AIDS virus)	1 minute
<i>Escherichia coli</i> (<i>E. coli</i>) O157:H7	2 minutes
Hepatitis B Virus (HBV)	2 minutes
Hepatitis C Virus (HCV)	2 minutes
Herpes Simplex Virus Type 2	2 minutes
Influenza A (A2 Japan)	2 minutes
Methicillin-Resistant <i>Staphylococcus aureus</i> (MRSA)	2 minutes
<i>Mycobacterium Tuberculosis</i> (TB)	2 minutes
<i>Pseudomonas aeruginosa</i>	2 minutes
<i>Salmonella enterica</i>	2 minutes
<i>Staphylococcus aureus</i>	2 minutes
Vancomycin-Resistant <i>Enterococcus faecalis</i> (VRE)	2 minutes
<i>Trichophyton mentagrophytes</i> (Athlete's foot fungus)	5 minutes

Kill (i.e., contact) time for a disinfectant is the amount of time a surface must remain wet with the product to achieve disinfection.



MSC351300

MSC351310

Micro-Kill One

- » The quaternary ammonium and alcohol solution dries evenly on hard surfaces with no beading or residue left behind
- » EPA-registered to kill 25 microorganisms in a one-minute* contact time
- » EPA Reg. No. 88494-2-37549

Item No.	Description	Pkg
MSC351300	6 x 6.75" Wipes, 160 Count Canister	12 ea/cs
MSC351310	7.5 x 15" Wipes, 65 Count Canister	12 ea/cs

*Ensure that the surface remains wet for the entire contact time.

Microorganism	Kill Time
<i>Acinetobacter baumannii</i>	1 minute
<i>Burkholderia cepacia</i>	1 minute
<i>Campylobacter jejuni</i>	1 minute
<i>Escherichia coli</i> O157:H7 (<i>E. coli</i>)	1 minute
Hepatitis B Virus (HBV)	1 minute
Hepatitis C Virus (HCV)	1 minute
Herpes Simplex Virus Type 1	1 minute
Human Immunodeficiency Virus (HIV-1) (AIDS)	1 minute
Human Coronavirus	1 minute
Influenza virus type A2 (Flu)	1 minute
<i>Kiebsiella pneumoniae</i>	1 minute
<i>Listeria monocytogenes</i>	1 minute
Methicillin Resistant <i>Staphylococcus aureus</i> (MRSA)	1 minute
<i>Mycobacterium bovis</i> (<i>Tuberculosis</i> (TB))	1 minute
Norovirus	1 minute
Poliovirus	1 minute
<i>Pseudomonas aeruginosa</i>	1 minute
Respiratory Syncytial Virus (RSV)	1 minute
Rhinovirus	1 minute
Rotavirus Strain WA	1 minute
<i>Salmonella enterica</i>	1 minute
<i>Staphylococcus aureus</i> (Staph)	1 minute
<i>Streptococcus pyogenes</i> (Strep)	1 minute
<i>Trichophyton mentagrophytes</i> (Athlete's Foot fungus)	1 minute
Vancomycin Resistant <i>Enterococcus faecium</i> (VRE)	1 minute

Kill (i.e., contact) time for a disinfectant is the amount of time a surface must remain wet with the product to achieve disinfection.

MSC351128



Micro-Kill+ Wipes Canister Wall-Mount Brackets

» Conveniently holds Micro-Kill and Micro-Kill+ canisters

Item No.	Description	Pkg
MSC351128	1.5 L Canister Bracket	12/cs

DURACELL

Duracell® Procell® Alkaline

» Long lasting and dependable
 » Date coded to ensure freshness



Item No.	Description	Pkg
DRCPC1500Z	AA Battery	24/bx
DRCPC2400Z	AAA Battery	24/bx
DRCPC1604Z	9V Battery	12/bx
DRCPC1400Z	C Battery	12/bx
DRCPC1300Z	D Battery	12/bx



MedCell Advantage Alkaline

» Mercury-free alkaline batteries

Item No.	Description	Pkg
MPHBAAAZ	AAA Battery	24/bx
MPHBAAZ	AA Battery	24/bx
MPHB9VZ	9V Battery	12/bx
MPHBCZ	C Battery	12/bx
MPHBDZ	D Battery	12/bx



DRC4133301201

DURACELL

Hearing Aid Batteries

» Duracell EasyTab® is the first easy-to-replace hearing aid battery
 » 3 Volt, Lithium, Coin Battery #2032, used in glucose monitors, remotes, calculators, and much more
 » Easy to remove, hold and insert

Item No.	Size	Pkg
DRC333661186	H/A, 1.4V, ZA, DA10B8ZM10	48/cs
DRC333661247Z	H/A, 1.4V, ZA, DA312B8ZM09	6/bx
DRC333661261	H/A, 1.4V, ZA, DA675B6ZM10	36/cs
DRC4133301201	Lithium, 3V	6/bx

AHS2131090H



Acme United Corporation

Eyewash

- » Sterile isotonic solution
- » Tamper-evident cap

Item No.	Description	Pkg
AHS2131090H	Squeeze Bottle: 4-oz.	1/ea



FAO90499001

FIRST AID ONLY

Eye/Face Wash Stations

- » Stations include bottles of saline solution with break-away cap
- » Wide-mouth bottles provide initial flushing for alkali and acid splash
- » Safety sealed and dated to guard against contamination
- » Complies with ANSI Z358 1-1990 safety standard for personal eyewash stations

Item No.	Description	Pkg
FAO24000001	Eye Wash Station, One 16-oz. Bottle	6 ea/cs
FAO2402001	Eye Wash Station, One 32-oz. Bottle	6 ea/cs
FAO24102001	Eye Wash Station, Two 16-oz. Bottles	6 ea/cs
FAO24300001	Eye Wash Station, Two 32-oz. Bottles	6 ea/cs
FAO24101001	Eye Wash, 16-oz. Bottle, Breakaway Top	12 ea/cs
FAO24201001	Eye Wash, 32-oz. Bottle, Re-closable Breakaway Bottle	12 ea/cs



MPH17CE210H



Fluid Clean-Up Kit

- » For safe containment and clean-up of spills and the decontamination of exposed hard surfaces

Item No.	Description	Pkg
MPH17CE210H	Kit, Spill, Biohazard, Standard	1/ea
MPH17CD410H	Kit, Spill, Biohazard, Deluxe	1/ea
MPHCP610H	Kit, Spill, Biohazard, Premium	1/ea
MPH17CE210*	Standard Spill Kit—Quat Pair Non-Latex Gloves, Bio-Hazard Bag, 2 Twist Ties, Scoop with Scraper, 2g Absorbent Powder, Hand Wipe, Quat Disinfectant Wipe	36/cs
MPH17CD410	Deluxe Spill Kit—Quat Quat Disinfectant Wipe, Pair Non-Latex Gloves, Bio-Hazard Bag, 2 Twist Ties, Scoop with Scraper, 8g Absorbent Powder, Mask with Face Shield, Gown/Apron	6/cs
MPH17CP610	Premium Spill Kit—2 pair Non-Latex gloves, 2 Red Biohazard Bags, 2 Twist Ties, 1 Scoop with Scraper, 10 grams of Absorbent Powder, 2 Hand Sanitizing Wipes, 2 Disinfectant Wipes, 1 Mask with Face Shield, 1 Gown/Apron and 1 Instructions Page	6/cs

*Not made with natural rubber latex.



MPH17CD410

EXAM GLOVES

Restore Gloves	70
Aloetouch Gloves	70
Nitrile Exam Gloves	70
Generation Pink Exam Gloves	71
Vinyl Exam Gloves	71
Latex Exam Gloves	71

SURGICAL GLOVES

Latex & Synthetic Surgical Gloves	72
---	----





OAT2585



OAT6801

Restore® Nitrile Exam Gloves

- » Contains a layer of maxOat+®, a proprietary blend of colloidal oatmeal that helps soothe and hydrate skin
- » Tested for use with chemotherapy drugs**
- » SmartGuard® film helps keep gloves in the box
- » Color-coded sizing quickly identifies the right size glove every time

Item No.	Description	Pkg
----------	-------------	-----

Restore Sense® Nitrile Exam Gloves

OAT2583*	X-Small	250/bx, 2,500/cs
OAT2584*	Small	250/bx, 2,500/cs
OAT2585*	Medium	250/bx, 2,500/cs
OAT2586*	Large	250/bx, 2,500/cs
OAT2587*	X-Large	230/bx, 2,300/cs

Restore™ Nitrile Exam Gloves

OAT6800*	X-Small	250/bx, 2,500/cs
OAT6801*	Small	250/bx, 2,500/cs
OAT6802*	Medium	250/bx, 2,500/cs
OAT6803*	Large	250/bx, 2,500/cs
OAT6804*	X-Large	230/bx, 2,300/cs

*Add "H" to be purchased by the box

**Data on file. Not for use with Carmustine or Thiotepa.



MDS195283



AloeTouch® Ice Nitrile Exam Gloves

- » Combines durable nitrile protection with a soothing inner coating of organic freeze-dried aloe vera
- » Fingertip textured for advanced grip

Item No.	Description	Pkg
MDS195283*	X-Small	200/bx, 2,000/cs
MDS195284*	Small	200/bx, 2,000/cs
MDS195285*	Medium	200/bx, 2,000/cs
MDS195286*	Large	200/bx, 2,000/cs
MDS195287*	X-Large	180/bx, 1,800/cs

*Add "H" to be purchased by the box.



MDS195175



AloeTouch® 3G Vinyl Exam Gloves

- » Third generation synthetic vinyl exam gloves that offer exceptional tactile sensitivity
- » Lined with organic freeze-dried aloe vera

Item No.	Description	Pkg
MDS195173	X-Small	100/bx, 1,000/cs
MDS195174*	Small	100/bx, 1,000/cs
MDS195175*	Medium	100/bx, 1,000/cs
MDS195176*	Large	100/bx, 1,000/cs
MDS195177*	X-Large	100/bx, 1,000/cs

*Add "H" to be purchased by the box.



MDS2585



SensiCare® Silk Nitrile Exam Gloves with SmartGuard

- » A thinner nitrile glove that offers durable protection
- » Tested for use with chemotherapy drugs**
- » Textured fingertips
- » SmartGuard film over the box opening and color-coded sizing around the rim

Item No.	Description	Pkg
MDS2583*	X-Small	250/bx, 2,500/cs
MDS2584*	Small	250/bx, 2,500/cs
MDS2585*	Medium	250/bx, 2,500/cs
MDS2586*	Large	250/bx, 2,500/cs
MDS2587*	X-Large	230/bx, 2,300/cs

*Add "H" to be purchased by the box.

**Data on file. Not for use with Carmustine or Thiotepa.



MDS6802



SensiCare® Ice Nitrile Exam Gloves with SmartGuard

- » Offers durable nitrile protection with a slightly thicker glove
- » Tested for use with chemotherapy drugs**
- » Fingertip textured for advanced grip
- » SmartGuard film over the box opening and color-coded sizing around the rim

Item No.	Description	Pkg
MDS6800*	X-Small	250/bx, 2,500/cs
MDS6801*	Small	250/bx, 2,500/cs
MDS6802*	Medium	250/bx, 2,500/cs
MDS6803*	Large	250/bx, 2,500/cs
MDS6804*	X-Large	230/bx, 2,300/cs

*Add "H" to be purchased by the box.

**Data on file. Not for use with Carmustine or Thiotepa.



FG3002



FitGuard® Touch Nitrile Exam Gloves

- » Textured fingertips and advanced tactile sensitivity
- » Powder-free
- » Tested for use with chemotherapy drugs**

Item No.	Description	Pkg
FG3000*	X-Small	300/bx, 3,000/cs
FG3001*	Small	300/bx, 3,000/cs
FG3002*	Medium	300/bx, 3,000/cs
FG3003*	Large	300/bx, 3,000/cs
FG3004*	X-Large	250/bx, 2,500/cs

*Add "H" to be purchased by the box.

**Please contact Medline for a complete list of the drugs and concentrations tested and the testing method used. Not for use with Carmustine and Thiotepa.



MDS193094



Wire Glove Dispensers

- » Universal design can mount horizontally
- » White-coated wire

Item No.	Description	Pkg
MDS192094B	Double	6/cs
MDS193094	Triple	10/cs



PINK2585

Generation Pink® Sense Nitrile Exam Gloves

- » Enhanced tactile sensitivity for greater control and touch
- » Textured fingertips for a better grip
- » Tested for use with select chemotherapy drugs**
- » Color-coded sizing helps quickly identify the right size glove

Item No.	Description	Pkg
PINK2583	X-Small	250/bx, 2,500/cs
PINK2584	Small	250/bx, 2,500/cs
PINK2585	Medium	250/bx, 2,500/cs
PINK2586	Large	250/bx, 2,500/cs
PINK2587	X-Large	230/bx, 2,300/cs

**Data on file. Not for use with Carmustine or Thiotepa.



PINK5085

Generation Pink® Pearl Nitrile Exam Gloves

- » Extra protection of nitrile
- » A layer of aloe vera helps soothe hands
- » Full texture for better grip

Item No.	Description	Pkg
PINK5083*	X-Small	100/bx, 1,000/cs
PINK5084*	Small	100/bx, 1,000/cs
PINK5085*	Medium	100/bx, 1,000/cs
PINK5086*	Large	100/bx, 1,000/cs
PINK5087*	X-Large	90/bx, 900/cs

*Add "H" to be purchased by the box.

For our complete selection of gloves, visit www.medline.com/physician.



PINK6802

Generation Pink Nitrile Exam Gloves

- » Excellent tactile sensitivity and barrier protection
- » Textured fingertips for better grip
- » Holds twice the amount of gloves and reduces trips to the supply closet

Item No.	Description	Pkg
PINK6800*	X-Small	250/bx, 2,500/cs
PINK6801*	Small	250/bx, 2,500/cs
PINK6802*	Medium	250/bx, 2,500/cs
PINK6803*	Large	250/bx, 2,500/cs
PINK6804*	X-Large	230/bx, 2,300/cs

*Add "H" to be purchased by the box.



PINK6075

Generation Pink 3G Vinyl Exam Gloves

- » Excellent tactile sensitivity
- » Stretches for flexibility and comfort
- » Not made with natural rubber

Item No.	Description	Pkg
PINK6073*	X-Small	100/bx, 1,000/cs
PINK6074*	Small	100/bx, 1,000/cs
PINK6075*	Medium	100/bx, 1,000/cs
PINK6076*	Large	100/bx, 1,000/cs
PINK6077*	X-Large	90/bx, 900/cs

*Add "H" to be purchased by the box.



CUR8105

Latex Exam Gloves

- » Have excellent barrier properties combined with great fit and feel
- » Powder-free
- » Fully textured for superior grip



MDS192075

Accutouch™ Vinyl Exam Gloves

- » Perfect cost-effective alternative is ideal for general patient care
- » Powder-free

Item No.	Description	Pkg
MDS192074*	Small	100/bx, 1,000/cs
MDS192075*	Medium	100/bx, 1,000/cs
MDS192076*	Large	100/bx, 1,000/cs
MDS192077*	X-Large	90/bx, 900/cs

*Add "H" to be purchased by the box.



CUR9225

Stretch Vinyl Exam Gloves


- » Advanced formulation stretch vinyl that conforms to the hands for excellent sensitivity
- » Economy vinyl option

Item No.	Description	Pkg
CUR9224*	Small	150/bx, 1,500/cs
CUR9225*	Medium	150/bx, 1,500/cs
CUR9226*	Large	150/bx, 1,500/cs
CUR9227*	X-Large	130/bx, 1,300/cs

*Add "H" to be purchased by the box.

Item No.	Description	Pkg
CUR8103*	X-Small	100/bx, 1,000/cs
CUR8104*	Small	100/bx, 1,000/cs
CUR8105*	Medium	100/bx, 1,000/cs
CUR8106*	Large	100/bx, 1,000/cs
CUR8107*	X-Large	90/bx, 900/cs

*Add "H" to be purchased by the box.



MSG2270




Triumph® Gloves

- » Sterile latex and Powder-free
- » Smooth for enhanced tactile sensitivity
- » Synthetic Polymer coating to improve donning

Item No.	Description	Pkg
MSG2255*	Size 5.5	50/bx, 200/cs
MSG2260*	Size 6	50/bx, 200/cs
MSG2265*	Size 6.5	50/bx, 200/cs
MSG2270*	Size 7	50/bx, 200/cs
MSG2275*	Size 7.5	50/bx, 200/cs
MSG2280*	Size 8	50/bx, 200/cs
MSG2285*	Size 8.5	50/bx, 200/cs
MSG2290*	Size 9	50/bx, 200/cs

*Add "Z" to be purchased by the box.



MSG6070



Neolon® 2G Gloves

- » Sterile, Powder-free, Synthetic Neoprene
- » Lightly textured, standard thickness
- » Made with 2nd generation polychloroprene material

Item No.	Description	Pkg
MSG6055*	Size 5.5	50/bx, 200/cs
MSG6060*	Size 6	50/bx, 200/cs
MSG6065*	Size 6.5	50/bx, 200/cs
MSG6070*	Size 7	50/bx, 200/cs
MSG6075*	Size 7.5	50/bx, 200/cs
MSG6080*	Size 8	50/bx, 200/cs
MSG6085*	Size 8.5	50/bx, 200/cs
MSG6090*	Size 9	50/bx, 200/cs

*Add "Z" to be purchased by the box.



MSG9070



SensiCare® PI Gloves

- » Sterile, Powder-free, Synthetic Polyisoprene
- » Smooth standard thickness
- » Manufactured using Medline's revolutionary polyisoprene material, Isolex™, for an optimal synthetic option

Item No.	Description	Pkg
MSG9055*	Size 5.5	50/bx, 200/cs
MSG9060*	Size 6	50/bx, 200/cs
MSG9065*	Size 6.5	50/bx, 200/cs
MSG9070*	Size 7	50/bx, 200/cs
MSG9075*	Size 7.5	50/bx, 200/cs
MSG9080*	Size 8	50/bx, 200/cs
MSG9085*	Size 8.5	50/bx, 200/cs
MSG9090*	Size 9	50/bx, 200/cs

*Add "Z" to be purchased by the box.



MDS108070LT



Triumph LT® Gloves

- » Sterile latex and Powder-free
- » Lightly textured, standard thickness

Item No.	Description	Pkg
MDS108055LT*	Size 5.5	50/bx, 200/cs
MDS108060LT*	Size 6	50/bx, 200/cs
MDS108065LT*	Size 6.5	50/bx, 200/cs
MDS108070LT*	Size 7	50/bx, 200/cs
MDS108075LT*	Size 7.5	50/bx, 200/cs
MDS108080LT*	Size 8	50/bx, 200/cs
MDS108085LT*	Size 8.5	50/bx, 200/cs
MDS108090LT*	Size 9	50/bx, 200/cs

*Add "Z" to be purchased by the box.



MSG6570



DermAssure™ Green Gloves

- » Not made with chemical accelerators or natural rubber latex
- » Dark green color and enhanced sensitivity make it an ideal underglove option
- » Smooth micro microthin surgical glove

Item No.	Description	Pkg
MSG6555*	Size 5.5	50/bx, 200/cs
MSG6560*	Size 6.0	50/bx, 200/cs
MSG6565*	Size 6.5	50/bx, 200/cs
MSG6570*	Size 7.0	50/bx, 200/cs
MSG6575*	Size 7.5	50/bx, 200/cs
MSG6580*	Size 8.0	50/bx, 200/cs
MSG6585*	Size 8.5	50/bx, 200/cs
MSG6590*	Size 9.0	50/bx, 200/cs

*Add "Z" to be purchased by the box.



MSG1070



SensiCare® with Aloe Gloves

- » Sterile, Powder-free, Synthetic Polyisoprene
- » Smooth standard thickness
- » Coated with 100% pure aloe vera to help soothe and moisturize dry skin

Item No.	Description	Pkg
MSG1055*	Size 5.5	25/bx, 100/cs
MSG1060*	Size 6	25/bx, 100/cs
MSG1065*	Size 6.5	25/bx, 100/cs
MSG1070*	Size 7	25/bx, 100/cs
MSG1075*	Size 7.5	25/bx, 100/cs
MSG1080*	Size 8	25/bx, 100/cs
MSG1085*	Size 8.5	25/bx, 100/cs
MSG1090*	Size 9	25/bx, 100/cs

*Add "Z" to be purchased by the box.

HAND HYGIENE

Hand Sanitizers.....	74-75
Antibacterial Soaps.....	76
Lotion Soaps.....	76-77
Foaming Soaps.....	77
Dispensers.....	77





HH85G02



Spectrum™ Clinical Gel Hand Sanitizer

- » Contains 85% ethyl alcohol
- » Contains aloe and vitamin E
- » Gel formula evaporates quickly on hands
- » >99.9999 log reduction on average

Item No.	Description	Pkg
HH85G02	2-oz.	48/cs
HH85G15	15-oz.	6/cs
HH85G1000	1,000 ml	4/cs



Spectrum Advance Gel Hand Sanitizer

- » Contains 70% ethyl alcohol along with Aloe and Vitamin E
- » Moisturizing formula for soft hands
- » Kills 99.99% of common germs on hands
- » Made with added emollients to help soothe skin

Item No.	Description	Pkg
HH70G02*	2-oz./59.1 ml	48/cs
HH70G04*	4-oz./118.3 ml	24/cs
HH70G08*	8-oz./236.6 ml	12/cs
HH70G15*	15-oz./443.6 ml	6/cs
HH70G1000*	33.8-oz./1,000 ml	4/cs
HH70GPACKET*	Packet	2,000/cs

*Add H to end of item number to order by the each.



HH62G02



Spectrum Gel Hand Sanitizer

- » Contains 62% ethyl alcohol along with aloe vera and vitamin E
- » Gel formula evaporates quickly on hands
- » Fragrance-free

Item No.	Description	Pkg
HH62G02	2-oz.	48/cs
HH62G15	15-oz.	6/cs
HH62G1000	1,000 ml	4/cs



Sterillium Comfort Gel Hand Sanitizers

- » Kills 99.999% of germs in 15 seconds without water
- » Helps increase skin moisture level
- » Leaves hands feeling soft and healthy
- » Dermatologist tested
- » Product is intended for use on the body only

Item No.	Description	Pkg
MSC097065	50 ml	42/cs
MSC097063	475 ml	20/cs
MSC097062	1,000 ml	8/cs



Spectrum Advance Foam Hand Sanitizer

- » Contains 70% ethyl alcohol
- » Added moisturizers help reduce dry skin
- » Foam formula evaporates quickly on hands

Item No.	Description	Pkg
HH70F75	7.5-oz.	12/cs
HH70F18	18-oz.	12/cs
HH70F1000	1,000 ml	4/cs



MSC097039A



HealthGuard Instant Hand Sanitizers

- » Contains ethyl alcohol and moisturizers to clean hands as it kills most common germs
- » Aloe vera and vitamin E to soothe the skin

Item No.	Description	Pkg
MSC097034NA	1.5-oz.	48/cs
MSC097039A	1,000 ml	12/cs



MSC097047



Aerosol Foam Hand Sanitizer

- » Leaves hands feeling soft and healthy
- » Contains 62% v/v ethyl alcohol
- » Drip-free nozzle design
- » Fragrance-free
- » Product is intended for use on the body only

Item No.	Description	Pkg
MSC097047	1.5-oz.	48/cs
MSC097048	8-oz.	24/cs
MSC097049	16-oz.	12/cs



KUT68241MEDA

HealthGuard Alcohol-Free Foaming Hand Sanitizer

- » Nonflammable
- » Benzethonium chloride formula with aloe vera and vitamin E
- » Has a light linen scent

Item No.	Description	Pkg
KUT68241MEDA	1 Liter	6/cs



MSC68841A

HealthGuard Foaming Hand Sanitizers

- » Nonaerosol foam contains 62% v/v ethyl alcohol
- » Foam formula quickly disperses on hand
- » With aloe vera and vitamin D to help soothe the skin

Item No.	Description	Pkg
MSC68817A	1.7-oz.	24/cs
MSC097618A	18-oz., Pump	6/cs
MSC68841A	1,000 ml	6/cs
MSC097688A	1 Liter, For Manual Dispenser	6/cs
MSC68824BA	1 Liter, For Auto Dispenser	2/cs



MSC350320

Spectrum Sanitizing Hand Wipes

- » Sanitize while providing friction needed to remove visible soils
- » Contains 65% v/v ethyl alcohol
- » Enriched with soothing moisturizers, aloe and vitamin E
- » Product is intended for use on the body only

Item No.	Description	Pkg
MSC350320	20 Count	48/cs
HH70W50	50 Count	6/cs
HH70W160	160 Count	6/cs



MSC097070



Sterillium Med Liquid Hand Sanitizer

- » Kills germs while providing exceptional dermal tolerance
- » 85% w/w ethyl alcohol
- » Evaporates quickly without floor-damaging thickening agents

Item No.	Description	Pkg
MSC097072	50 ml	50/cs
MSC097071	500 ml, Pump	20/cs
MSC097070	1,000 ml	8/cs

See pages 74-75 for hand sanitizers and page 77 for soap dispensers.



Spectrum Antibacterial Soap

- » Spectrum antibacterial soap helps kill harmful germs
- » Product is intended for use on the body only
- » Active ingredient: 0.13% benzalkonium chloride

Item No.	Description	Pkg
HHABSP04	4-oz./118.3 ml	12/cs
HHABSP75	7.5-oz./221.8 ml	12/cs
HHABSP16	16-oz./473.2 ml	12/cs
HHABSP1000	33.8-oz./1,000 ml	4/cs
HHABSPG	1-gal./3.8 L	4/cs



KUT65036

Skintegritty Antimicrobial Soap

- » Antibacterial
- » Enhanced with moisturizing aloe vera and nourishing Vitamin E
- » Product is intended for use on the body only

Item No.	Description	Pkg
KUT65036	1,000 ml	8/cs



MSC098200A

HealthGuard Antibacterial Soap

- » Contains 0.13% benzalkonium chloride to kill 99.9% of the most common micro-organisms that can cause illness
- » Moisturizes and deodorizes while leaving hands clean
- » Product is intended for use on the body only

Item No.	Description	Pkg
MSC098204A	4-oz.	24/cs
MSC098200A	800 ml	12/cs
MSC098203A	1,000 ml	12/cs



HHABSP1000FM

Spectrum Foaming Antibacterial Soap

- » Spectrum antibacterial soap helps kill harmful germs
- » Product is intended for use on the body only
- » Active ingredient: 0.10% benzalkonium chloride

Item No.	Description	Pkg
HHABSP1000FM	33.8-oz./1,000 ml	4/cs



KUT64041MEDA

HealthGuard Antibacterial Foaming Soap

- » Spring citrus fragrance
- » Formula contains 0.13% benzalkonium chloride that kills 99.9% of common bacteria that can cause illness
- » Product is intended for use on the body only

Item No.	Description	Pkg
KUT64078MEDA	950 ml	4/cs
KUT64041MEDA	1,000 ml	6/cs
KUT64024MEDA	1,000 ml, For Auto Dispenser	2/cs



Spectrum Lotion Soap

- » Aloe vera to condition skin
- » Gently moisturizes and cleans hands
- » Product is intended for use on the body only

Item No.	Description	Pkg
HHSP04	4-oz./118.3 ml	12/cs
HHSP75	7.5-oz./221.8 ml	12/cs
HHSP16	16-oz./473.2 ml	12/cs
HHSP1000	33.8-oz./1,000 ml	4/cs
HHSPG	1 gal./3.8 L	4/cs



MSC098103



HealthGuard Enriched Lotion Soap

- » Provides moisture to frequently washed hands
- » Gentle formula and soothing fragrance nurture both the skin and the senses
- » pH-balanced to help maintain proper skin and moisture ratio

Item No.	Description	Pkg
MSC098104	4-oz.	24/cs
MSC098114	16-oz.	12/cs
MSC098103	1,000 ml	12/cs



HHSP1000FM



General Foaming Soap

- » Easy foaming
- » Gentle enough for everyday use
- » Product is intended for use on the body only

Item No.	Description	Pkg
HHSP1000FM	33.8-oz./1,000 ml	4/cs



MSC69024MED



HealthGuard Foaming Soap for Auto Dispenser

- » For use with dispenser MSC9962GRA or MSC9962BLK dispensers
- » Product is intended for use on the body only

Item No.	Description	Pkg
MSC69024MED	1 Liter	2/cs



KUT69041MED



HealthGuard Luxury Foaming Soap

- » Offers a clean contemporary dispenser with a unique push pad that dispenses 0.75 milliliters of rich foam per push
- » Provides twice the number of hand washings per 1,000 milliliters bag
- » Product is intended for use on the body only

Item No.	Description	Pkg
KUT69041MED	1,000 ml	6/cs
KUT64341MED	1,000 ml, Fragrance Free	6/cs



AUTODISPW



MANDISPB

Spectrum Dispensers

- » Easy installation

Item No.	Description	Pkg
----------	-------------	-----

Manual Dispensers

MANDISPW	White	6/cs
MANDISPB	Black	6/cs

Automatic Dispensers

AUTODISPW	White	1/cs
AUTODISPB	Black	1/cs



MSC9941GRA



MSC9962BLK

Dispensers

- » Can be used with a variety of our hand soaps and hand antiseptics

Item No.	Description	Pkg
----------	-------------	-----

Bag-In-Box

MSC9950	White	6/cs
MSC9951	Black	6/cs

Foaming Manual Push

MSC9941BLK	Black	1/ea
MSC9941GRA	Gray	1/ea

Foaming Automatic

MSC9962BLK	Black	1/ea
MSC9962GRA	Gray	1/ea



LET'S PARTNER TOGETHER TO HELP STOP THE FLU.



One stop for all your flu vaccine needs. Keep your patients and staff healthy with everything you need to fight the flu.

- » Pre-book – 100% of pre-books filled and shipped in last four years
- » Leading brands to meet all your needs
- » Flexible returnability on all prebooks
- » Custom shipping options
- » Open stock availability during flu season
- » Three easy ways to order – online, phone or fax

Contact your Medline Representative, Customer Service at 1-855-294-9618 or visit medline.com/flu for more information.



INSTRUMENTS

Family & General Practice	80-83
Obstetrics & Gynecology	84-87
Ear, Nose & Throat	87-88
Vein & Vascular	88-89
Blades and Disposables	89-90
Instrument Accessories	90
Procedure Trays	90-91
Cleaning & Sterilization	91-92

König instruments are handcrafted in Germany, offering durability and precision.

- » We warranty every König surgical instrument to be free of functional defects in workmanship and materials
- » Instruments may be purchased individually or in custom sets to meet any need or preference
- » The König line includes more than 18,000 surgical instruments. Contact us to learn more, and request the König catalog for more options.

INSTRUMENTS



MDS0228820

Percussion Hammer

Item No.	Description	Pkg
MDS0228820	Hammer, Percussion, Taylor	1/ea



MDS0610301

Knife Handle

Item No.	Description	Pkg
MDS0610301	#3, Graduated	1/ea



MDS0800411

Scissors

Item No.	Description	Pkg
MDS0839011	Iris, Knapp, Straight, Sharp/Sharp, 4"/10.2 cm	1/ea
MDS0800411	Fine Operating, Straight, 4.75"/12 cm	1/ea
MDS0838510	Iris, Knapp, Curved, Sharp/Sharp, 4"/10.2 cm	1/ea
MDS0834011	Iris, Straight, Sharp/Sharp, 4.5"/11 cm	1/ea
MDS0816014	Mayo, Straight, 5.5"/14 cm	1/ea
MDS0816114	Mayo, Curved, 5.5"/14 cm	1/ea
MDS0816017	Mayo, Straight, 6.75"/ 17 cm	1/ea
MDS0816117	Mayo, Curved, 6.75"/ 17 cm	1/ea
MDS0880209	Stitch, Spencer, 3.5"/9 cm	1/ea
MDS0880111	Stitch, Spencer-Delicate, 4.75"/12 cm	1/ea
MDS0880414	Stitch, Littauer, 5.5"/14 cm	1/ea
MDS0728981	Iris, Knapp, Straight, Light, 4"/10.2 cm	1/ea
MDS0728986	Iris, Knapp, Curved, Light, 4"/10.2 cm	1/ea
MDS0828115	Metzenbaum, Blunt/Blunt, Curved, 5.75"/14.6 cm	1/ea
MDS0820909	Gradle, 3.75"/9.5 cm	1/ea

Bandage Scissors

Item No.	Description	Pkg
MDS0890114	Lister, 5.5"/14 cm	1/ea
MDS0890118	Lister, Angled, 7.25"/18.4 cm	1/ea
MDS0895014	Universal, 5.75"/14.5 cm	1/ea
MDS0895018	Universal, 7.5"/19 cm	1/ea
MDS0895016	Universal, 6.25"/16 cm	1/ea
MDS0895015	Universal, 6.5"/16.5 cm	1/ea



MDS1017012

Forceps

Item No.	Description	Pkg
MDS1017012	Adson, Dressing, 4.75"/12 cm	1/ea
MDS1018612	Adson, Tissue, Micro, 1 x 2 Teeth, 4.75"/12 cm	1/ea
MDS1018012	Adson, 1 x 2, 4.75"/12 cm	1/ea
MDS1080209	Splinter, Serrated, Straight, 3.5"/9 cm	1/ea
MDS1080211	Splinter, Serrated, Straight, 4.5"/11 cm	1/ea
MDS1010214	Dressing, Serrated, 5.5"/14 cm	1/ea
MDS1012014	Tissue, Standard, 1 x 2 Teeth, 5.5"/14 cm	1/ea
MDS0725964	Dressing/Tissue, Bishop-Harmon, 3.375"/8.6 cm	1/ea
MDS0746773	Swiss, Cilia, 3 mm, Diag. Jar, 3.75"/9.5 cm	1/ea

Hemostats

Item No.	Description	Pkg
MDS1222012	Halsted, Mosquito, Straight, 5"/13 cm	1/ea
MDS1222112	Halsted, Mosquito, Curved, 5"/13 cm	1/ea
MDS1224014	Kelly, Straight, 5.5"/14 cm	1/ea
MDS1224114	Kelly, Curved, 5.5"/14 cm	1/ea
MDS1225014	Crile, Straight, 5.5"/14 cm	1/ea
MDS1225114	Crile, Curved, 5.5"/14 cm	1/ea



MDS2410013

Needle Holders

Item No.	Description	Pkg
MDS2410013	Webster, Smooth, 5"/13 cm	1/ea
MDS2418014	Mayo-Hegar, Serrated, 5.5"/14 cm	1/ea
MDS2441414	Olsen Hegar, TC, 5.5"/14 cm	1/ea
MDS2413013	Halsey, Smooth, 5"/13 cm	1/ea
MDS2410413	Webster, Smooth, TC, 4.75"/12.1 cm	1/ea
MDS2414413	Halsey, Serrated, TC, 5"/12.2 cm	1/ea
MDS2418414	Mayo-Hegar, TC, 5"/12.2 cm	1/ea



MDS4641518

Elevators


Item No.	Description	Pkg
MDS4641518	Freer, Septum, Double-ended, SH/BL, 7.75"/19.7 cm	1/ea
MDS1211951	Freer, Double-ended, Sharp, Slight Curved, 8.5"/21.5 cm	1/ea
MDS1211856	Periosteal, No. 856, 7"/17.8 cm	1/ea



MDS7433614

Nail Nippers, Double Spring

Item No.	Description	Pkg
MDS7433614	Concave Jaws, 5.5"/14 cm	1/ea
MDS7436013	Straight Jaws, 5"/12.7 cm	1/ea
MDS7406013	Concave Jaws, 4.75"/12 cm	1/ea
MDS7436213	English Anvil Pattern, 5"/12.7 cm	1/ea



MDS7434113

Nail Splitters

Item No.	Description	Pkg
MDS7434113	Concave Jaws, 5"/12.7 cm	1/ea
MDS7411661	#60, Angled, 2 mm, 5"/12.7 cm	1/ea



MDS7436511

Cuticle Nippers

Item No.	Description	Pkg
MDS7436511	Convex Jaws, 4.25"/10.8 cm	1/ea



MDS7432711

Platypus Nail-Pulling Forceps


Item No.	Description	Pkg
MDS7432711	Standard Dusted Jaws, 5.5"/14 cm	1/ea



MDS7411110

Extractor, Comedone, Schamberg

Item No.	Description	Pkg
MDS7411110	Double, 3.75"/9.5 cm	1/ea



MDS205777

Disposable Ear Curettes

- » Seven color-coded styles included
- » Variety pack has 15 white curettes and 10 of each of the remaining six colors

Item No.	Description	Pkg
MDS209555	Loop, 4 mm, White	50/bx
MDS204888	Small Spoon, 2 mm, Blue	50/bx
MDS201222	Circle Loop, 3 mm, Green	50/bx
MDS204111	Flex Loop, 3 mm, Purple	50/bx
MDS206333	Spoon, 4 mm, Yellow	50/bx
MDS203444	Angled Loop, 4 mm, Red	50/bx
MDS202999	Bent Loop, 4 mm, Orange	50/bx
MDS205777	Variety Pack	75/bx



MDS0033320

Disposable Sterile Biopsy Punches

- » Razor-sharp stainless steel blade for precise excisions

Item No.	Description	Pkg
MDS0033310	1 mm	50/bx
MDS0033315	1.5 mm	50/bx
MDS0033320	2 mm	50/bx
MDS0033325	2.5 mm	50/bx
MDS0033330	3 mm	50/bx
MDS0033335	3.5 mm	50/bx
MDS0033340	4 mm	50/bx
MDS0033350	5 mm	50/bx
MDS0033360	6 mm	50/bx
MDS0033380	8 mm	50/bx



Disposable Dermal Curettes

Item No.	Description	Pkg
MDS0033520	2 mm	50/bx
MDS0033530	3 mm	50/bx
MDS0033540	4 mm	50/bx
MDS0033550	5 mm	50/bx
MDS0033570	7 mm	50/bx



Curettes-Excavators

Item No.	Description	Pkg
Single-Ended Without hole		
MDS0740581	1 mm Dia, 5"/13 cm, Without Hole	1/ea
MDS2591211	1 mm Dia, 5"/13 cm, with Hole	1/ea
MDS7440572	2 mm Dia, 5"/13 cm, Without Hole	1/ea
MDS0740852	2 mm Dia, 5"/13 cm, with Hole	1/ea
MDS7440573	2.5 mm Dia, 5"/13 cm, Without Hole	1/ea
MDS0740503	2.5 mm Dia, 5"/13 cm, with Hole	1/ea
MDS7440574	4 mm Dia, 5"/13 cm, Without Hole	1/ea
MDS0740504	4 mm Dia, 5"/13 cm, with Hole	1/ea

Double End

MDS7405812	with Hole, 1.5-2 mm	1/ea
MDS7405813	with Hole, 1.5-2.55 mm	1/ea
MDS74113420	Double End, 1.5-2 mm	1/ea



Nucleus Knives

Item No.	Description	Pkg
MDS0391657	S/E, 2 mm	1/ea
MDS0391656	S/E, 3 mm	1/ea
MDS0391650	S/E, 4 mm	1/ea



Excavators

Item No.	Description	Pkg
MDS7444089	Probe, Double End, 6"/15 cm	1/ea
MDS7444090	Packer, Double End, 5.75"/14.6 cm	1/ea
MDS0744091	Elevator, Spatula, Packer, Straight, 5.75"/14.6 cm	1/ea
MDS0740975	Probe, Packer, Double End, 6"/15 cm	1/ea



Integra® Miltex® Products

Item No.	Description	Pkg
Reflex Hammer		
MIL1202	Taylor, 7.5"/19 cm	1/ea
Disposable Scalpels		
MIL4410	#10	10/bx
MIL4411	#11	10/bx
MIL4415	#15	10/bx
Knife Handle		
MIL47	#3 with Scale	1/ea
Scissors		
MIL516	OR, Straight, Sharp/Blunt, 5.5"/14 cm	1/ea
MIL5304	Iris, Straight, 4.5"/11.4 cm	1/ea
MIL5306	Iris, Curved, 4.5"/11.4 cm	1/ea
MIL5514	Lister Bandage, 5.5"/14 cm	1/ea
MIL5516	Lister Bandage, 7.25"/18.4 cm	1/ea
MIL51000	Bandage and Utility, 7.5"/19 cm	1/ea
MIL9100	Spencer Stitch, 3.5"/8.9 cm	1/ea
MIL9104	Littauer Stitch, 5.5"/14 cm	1/ea



Item No.	Description	Pkg
EKG Caliper		
MIL13015	Economy	1/ea
Instrument Brush		
MIL31000	Nylon, ST/3	3/set
Forceps		
MIL6118	Adson Dressing, 4.75"/12.1 cm	1/ea
MIL6120	Adson Tissue, 4.75"/12.1 cm, 1 x 2	1/ea
Needle Holders		
MIL86	Webster, Smooth Jaw, 5"/12.7 cm	1/ea
MIL842	Mayo-Hegar, 5"/12.7 cm	1/ea
Nail Nippers		
MILV940200	4.5"/11.4 cm	1/ea
MILV940205	4.5"/11.4 cm BRL SPR	1/ea
MILV940207	5.5"/14 cm BRL SPR	1/ea
Disposable Safety Scalpels		
MIL4510	#10	10/bx
MIL4511	#11	10/bx
MIL4515	#15	10/bx



Integra® Miltex® Products

Item No.	Description	Pkg
Carbon Surgical Blades		
MIL4110	#10	100/bx
MIL4111	#11	100/bx
MIL4115	#15	100/bx
Scissors		
MIL5300	Iris, Straight, 4"/10.2 cm	1/ea
MIL5302	Iris, Curved, 4"/10.2 cm	1/ea
MIL5180	Metzenbaum, Curved, 5.5"/14 cm	1/ea
MIL181652	Eye Suture, Curved, 3.75"/9.5 cm	1/ea
MIL18SC1652	S-C Gradle, 3.75"/9.5 cm	1/ea
Forceps		
MIL18868	Bishop-Harmon, 3.375", 1 x 2 Teeth	1/ea
MIL181107	Swiss Cilia, 3.75"/9.5 cm	1/ea
Needle Holders		
MIL86TC	Webster, Smooth, 5"/12.7 cm	1/ea
MIL88ATC	Halsey, Serrated, 5"/12.7 cm	1/ea
MIL842TC	Mayo-Hegar, 5"/12.7 cm	1/ea



Item No.	Description	Pkg
Schamberg Comedone Extractors		
MIL33200	Square	1/ea
MIL33201	Crimp	1/ea
MIL33202	X-Fine	1/ea
Disposable Biopsy Punches		
MIL3331	2 mm	50/bx
MIL3331HH	2 mm	1/ea
MIL3332	3 mm	50/bx
MIL3332HH	3 mm	1/ea
MIL3333	3.5 mm	50/bx
MIL3333HH	3.5 mm	1/ea
MIL3334	4 mm	50/bx
MIL3334HH	4 mm	1/ea
MIL3335	5 mm	50/bx
MIL3335H	5 mm	1/ea
MIL3336	6 mm	50/bx
MIL3336HH	6 mm	1/ea
Disposable Dermal Curettes		
MIL3352	2 mm	50/bx
MIL3353	3 mm	50/bx
MIL3354	4 mm	50/bx
MIL3355	5 mm	50/bx
MIL3357	7 mm	50/bx



Single-Use Sterile Disposable Instruments

Item No.	Description	Pkg
Forceps		
TRI66590	Alligator Forceps, Serrated Jaw, 3 1/2"/8.9 cm	20/cs
TRI66790	Magill Forceps, Serrated, Curved, 9 3/4"/24.8 cm	20/cs
TRI66195	Adson Forceps, Serrated, 4 3/4"/12.1 cm	20/cs
Hemostats		
TRI66140	Mosquito Forceps, Straight, 5"/12.7 cm	20/cs
TRI66145	Mosquito Forceps, Curved, 5"/12.7 cm	20/cs
Scissors		
TRI66155	Iris Scissors, Straight, 4 1/2"/11.4 cm	20/cs
TRI66855	Lister Bandage Scissors, Angled, 7 1/4"/18.4 cm	20/cs
Misc		
MDS10749	Taylor Percussion Hammer	12/cs
MDS10755	Probe with Eye, 5"/12.7cm	12/cs
Instrument Trays		
DYND71010	Suture Removal Tray	40/cs
DYNJ03014	Laceration Tray	16/cs

For more information about our full line of single sterile disposable instruments, please contact your Medline representative.



MDS7011001



MDS1617025



MDS7036024



MDS7058132



MDS7059440



MDS2010616



Instruments

Item No.	Description	Pkg
----------	-------------	-----

Vaginal Specula

MDS7011001	Graves, Small, 29.9 x 7.9"/76 x 20 mm	1/ea
MDS7011002	Graves, Medium, 3.7 x 1.4"/95 x 35 mm	1/ea
MDS7011003	Graves, Large, 45.3 x 13.8"/115 x 35 mm	1/ea
MDS7012001	Pederson, Small, 29.5 x 5.1"/75 x 13 mm	1/ea
MDS7012002	Pederson, Medium, 40.2 x 8.7"/102 x 22 mm	1/ea
MDS7012003	Pederson, Large, 47.2 x 9.8"/120 x 25 mm	1/ea

Forceps

MDS1617025	Sponge, Foerster, Straight, 9.75"/24.8 cm	1/ea
MDS1623126	Packing, Bozeman, Curved, 10.25"/26 cm	1/ea
MDS7050425	Tenaculum, Braun, 10"/25 cm	1/ea
MDS7067201	Biopsy, Kevorkian-Younger, Alligator	1/ea
MDS7067211	Tischler, Cervical, Oblong Spec., 8.75"/22.2 cm	1/ea

Item No.	Description	Pkg
----------	-------------	-----

Uterine Sound

MDS7058132	Sims, Malleable, Graduated, 12.5"/31.8 cm	1/ea
------------	---	------

Endospeculum

MDS7036024	Gynecology, Kogan, 11"/27.9 cm, Screw	1/ea
------------	---------------------------------------	------

Curette

MDS7059440	Kevorkian, Biopsy, Curved, 11.75"/29.8 cm	1/ea
------------	---	------

Cervical Biopsy Punch

MDS7063101	Tischler, 3 x 7.9"/7.6 x 20 cm	1/ea
------------	--------------------------------	------

Probe

MDS2010616	Buttoned, Double-ended 6.25"/2 mm, 16cm	1/ea
------------	---	------



MDS6300300

MDS6300500

MDS6300400



Cube Pessaries

Item No.	Description	Pkg
----------	-------------	-----

Without Drain

MDS6300300	Size 0, 25 mm	1/ea
MDS6300301	Size 1, 29 mm	1/ea
MDS6300302	Size 2, 33 mm	1/ea
MDS6300303	Size 3, 37 mm	1/ea
MDS6300304	Size 4, 41 mm	1/ea
MDS6300305	Size 5, 45 mm	1/ea
MDS6300306	Size 6, 50 mm	1/ea
MDS6300307	Size 7, 56 mm	1/ea
MDS6300308	Size 8, 63 mm	1/ea
MDS6300309	Size 9, 70 mm	1/ea
MDS6300310	Size 10, 75 mm	1/ea

With Drain

MDS6300400	Size 0, 25 mm	1/ea
MDS6300401	Size 1, 29 mm	1/ea
MDS6300402	Size 2, 33 mm	1/ea
MDS6300403	Size 3, 37 mm	1/ea
MDS6300404	Size 4, 41 mm	1/ea
MDS6300405	Size 5, 45 mm	1/ea
MDS6300406	Size 6, 50 mm	1/ea
MDS6300407	Size 7, 56 mm	1/ea
MDS6300408	Size 8, 63 mm	1/ea
MDS6300409	Size 9, 70 mm	1/ea
MDS6300410	Size 10, 75 mm	1/ea



Donut Pessaries

Item No.	Description	Pkg
----------	-------------	-----

MDS6300500	Size 0, 2"/51 mm	1/ea
MDS6300501	Size 1, 2.25"/57 mm	1/ea
MDS6300502	Size 2, 2.5"/64 mm	1/ea
MDS6300503	Size 3, 2.75"/70 mm	1/ea
MDS6300504	Size 4, 3"/76 mm	1/ea
MDS6300505	Size 5, 3.25"/83 mm	1/ea
MDS6300506	Size 6, 3.5"/89 mm	1/ea
MDS6300507	Size 7, 3.75"/95 mm	1/ea



The König line includes more than 18,000 surgical instruments.

Contact us to learn more, and to receive the König catalog for more options.



Dish Pessaries

Item No.	Description	Pkg
Without Support		
MDS6300600	Size 0, 50 mm	1/ea
MDS6300601	Size 1, 55 mm	1/ea
MDS6300602	Size 2, 60 mm	1/ea
MDS6300603	Size 3, 65 mm	1/ea
MDS6300604	Size 4, 70 mm	1/ea
MDS6300605	Size 5, 75 mm	1/ea
MDS6300606	Size 6, 80 mm	1/ea
MDS6300607	Size 7, 85 mm	1/ea

With Support

MDS6300700	Size 0, 50 mm	1/ea
MDS6300701	Size 1, 55 mm	1/ea
MDS6300702	Size 2, 60 mm	1/ea
MDS6300703	Size 3, 65 mm	1/ea
MDS6300704	Size 4, 70 mm	1/ea
MDS6300705	Size 5, 75 mm	1/ea
MDS6300706	Size 6, 80 mm	1/ea
MDS6300707	Size 7, 85 mm	1/ea



Gehring Pessaries

Item No.	Description	Pkg
MDS6300900	Size 0, 50 mm	1/ea
MDS6300901	Size 1, 55 mm	1/ea
MDS6300902	Size 2, 60 mm	1/ea
MDS6300903	Size 3, 65 mm	1/ea
MDS6300904	Size 4, 70 mm	1/ea
MDS6300905	Size 5, 75 mm	1/ea
MDS6300906	Size 6, 80 mm	1/ea
MDS6300907	Size 7, 85 mm	1/ea
MDS6300908	Size 8, 90 mm	1/ea
MDS6300909	Size 9, 95 mm	1/ea

Gellhorn Pessaries

Item No.	Description	Pkg
With Stem		
MDS6300800	Size 0, 1.5"/38 mm	1/ea
MDS6300801	Size 1, 1.75"/44 mm	1/ea
MDS6300802	Size 2, 2"/51 mm	1/ea
MDS6300803	Size 3, 2.25"/57 mm	1/ea
MDS6300804	Size 4, 2.5"/64 mm	1/ea
MDS6300805	Size 5, 2.75"/70 mm	1/ea
MDS6300806	Size 6, 3"/76 mm	1/ea
MDS6300807	Size 7, 3.25"/83 mm	1/ea

With Short Stems

MDS6301000	Size 0, 1.5"/38 mm	1/ea
MDS6301001	Size 1, 1.75"/44 mm	1/ea
MDS6301002	Size 2, 2"/51 mm	1/ea
MDS6301003	Size 3, 2.25"/57 mm	1/ea
MDS6301004	Size 4, 2.5"/64 mm	1/ea
MDS6301005	Size 5, 2.75"/70 mm	1/ea
MDS6301006	Size 6, 3"/76 mm	1/ea
MDS6301007	Size 7, 3.25"/83 mm	1/ea
MDS6301008	Size 8, 3.5"/89 mm	1/ea
MDS6300808	Size 8, 3.5"/89 mm	1/ea



Oval Pessaries with Support

Item No.	Description	Pkg
MDS6301501	Size 1, 2"/51 mm	1/ea
MDS6301502	Size 2, 2.25"/57 mm	1/ea
MDS6301503	Size 3, 2.5"/64 mm	1/ea
MDS6301504	Size 4, 2.75"/70 mm	1/ea
MDS6301505	Size 5, 3"/76 mm	1/ea
MDS6301506	Size 6, 3.25"/83 mm	1/ea
MDS6301507	Size 7, 3.5"/89 mm	1/ea
MDS6301508	Size 8, 3.75"/95 mm	1/ea
MDS6301509	Size 9, 4"/102 mm	1/ea

Hodge Pessaries

Item No.	Description	Pkg
With Support		
MDS6301200	Size 0, 65 mm	1/ea
MDS6301201	Size 1, 70 mm	1/ea
MDS6301202	Size 2, 75 mm	1/ea
MDS6301203	Size 3, 80 mm	1/ea
MDS6301204	Size 4, 85 mm	1/ea
MDS6301205	Size 5, 90 mm	1/ea
MDS6301206	Size 6, 95 mm	1/ea
MDS6301207	Size 7, 100 mm	1/ea
MDS6301208	Size 8, 105 mm	1/ea
MDS6301209	Size 9, 110 mm	1/ea

Without Support

MDS6301100	Size 0, 65 mm	1/ea
MDS6301101	Size 1, 70 mm	1/ea
MDS6301102	Size 2, 75 mm	1/ea
MDS6301103	Size 3, 80 mm	1/ea
MDS6301104	Size 4, 85 mm	1/ea
MDS6301105	Size 5, 90 mm	1/ea
MDS6301106	Size 6, 95 mm	1/ea
MDS6301107	Size 7, 100 mm	1/ea
MDS6301108	Size 8, 105 mm	1/ea
MDS6301109	Size 9, 110 mm	1/ea



Shatz Pessaries

Item No.	Description	Pkg
MDS6302000	Size 0, 1.5"/38 mm	1/ea
MDS6302001	Size 1, 1.75"/44 mm	1/ea
MDS6302002	Size 2, 2"/51 mm	1/ea
MDS6302003	Size 3, 2.25"/57 mm	1/ea
MDS6302004	Size 4, 2.5"/64 mm	1/ea
MDS6302005	Size 5, 2.75"/70 mm	1/ea
MDS6302006	Size 6, 3"/76 mm	1/ea
MDS6302007	Size 7, 3.25"/83 mm	1/ea
MDS6302008	Size 8, 3.5"/89 mm	1/ea



Marland Pessaries

Item No.	Description	Pkg
Without Support		
MDS6301302	Size 2, 2.25"/57 mm	1/ea
MDS6301303	Size 3, 2.5"/64 mm	1/ea
MDS6301304	Size 4, 2.75"/70 mm	1/ea
MDS6301305	Size 5, 3"/76 mm	1/ea
MDS6301306	Size 6, 3.25"/83 mm	1/ea
MDS6301307	Size 7, 3.5"/89 mm	1/ea
MDS6301308	Size 8, 3.75"/95 mm	1/ea

With Support

MDS6301402	Size 2, 2.25"/57 mm	1/ea
MDS6301403	Size 3, 2.5"/64 mm	1/ea
MDS6301404	Size 4, 2.75"/70 mm	1/ea
MDS6301405	Size 5, 3"/76 mm	1/ea
MDS6301406	Size 6, 3.25"/83 mm	1/ea
MDS6301407	Size 7, 3.5"/89 mm	1/ea
MDS6301408	Size 8, 3.75"/95 mm	1/ea



Knobbed Ring Pessaries

Item No.	Description	Pkg
Without Support		
MDS6301700	Size 0, 1.75"/44 mm	1/ea
MDS6301701	Size 1, 2"/51 mm	1/ea
MDS6301702	Size 2, 2.25"/57 mm	1/ea
MDS6301703	Size 3, 2.5"/64 mm	1/ea
MDS6301704	Size 4, 2.75"/70 mm	1/ea
MDS6301705	Size 5, 3"/76 mm	1/ea
MDS6301706	Size 6, 3.25"/83 mm	1/ea
MDS6301707	Size 7, 3.5"/89 mm	1/ea

With Support

MDS6301800	Size 0, 1.75"/44 mm	1/ea
MDS6301801	Size 1, 2"/51 mm	1/ea
MDS6301802	Size 2, 2.25"/57 mm	1/ea
MDS6301803	Size 3, 2.5"/64 mm	1/ea
MDS6301804	Size 4, 2.75"/70 mm	1/ea
MDS6301805	Size 5, 3"/76 mm	1/ea
MDS6301806	Size 6, 3.25"/83 mm	1/ea
MDS6301807	Size 7, 3.5"/89 mm	1/ea



Ring Pessaries

Item No.	Description	Pkg
Without Support		
MDS6301600	Size 0, 1.75"/44 mm	1/ea
MDS6301601	Size 1, 2"/51 mm	1/ea
MDS6301602	Size 2, 2.25"/57 mm	1/ea
MDS6301603	Size 3, 2.5"/64 mm	1/ea
MDS6301604	Size 4, 2.75"/70 mm	1/ea
MDS6301605	Size 5, 3"/76 mm	1/ea
MDS6301606	Size 6, 3.25"/83 mm	1/ea
MDS6301607	Size 7, 3.5"/89 mm	1/ea
MDS6301608	Size 8, 3.75"/95 mm	1/ea
MDS6301609	Size 9, 4"/102 mm	1/ea

With Support

MDS6301900	Size 0, 1.75"/44 mm	1/ea
MDS6301901	Size 1, 2"/51 mm	1/ea
MDS6301902	Size 2, 2.25"/57 mm	1/ea
MDS6301903	Size 3, 2.5"/64 mm	1/ea
MDS6301904	Size 4, 2.75"/70 mm	1/ea
MDS6301905	Size 5, 3"/76 mm	1/ea
MDS6301906	Size 6, 3.25"/83 mm	1/ea
MDS6301907	Size 7, 3.5"/89 mm	1/ea
MDS6301908	Size 8, 3.75"/95 mm	1/ea
MDS6301909	Size 9, 4"/102 mm	1/ea




Cup Pessaries without Support

Item No.	Description	Pkg
Without Support		
MDS6300100	Size 0, 50 mm	1/ea
MDS6300101	Size 1, 55 mm	1/ea
MDS6300102	Size 2, 60 mm	1/ea
MDS6300103	Size 3, 65 mm	1/ea
MDS6300104	Size 4, 70 mm	1/ea
MDS6300105	Size 5, 75 mm	1/ea
MDS6300106	Size 6, 80 mm	1/ea
MDS6300107	Size 7, 85 mm	1/ea


With Support

MDS6300200	Size 0, 50 mm	1/ea
MDS6300201	Size 1, 55 mm	1/ea
MDS6300202	Size 2, 60 mm	1/ea
MDS6300203	Size 3, 65 mm	1/ea
MDS6300204	Size 4, 70 mm	1/ea
MDS6300205	Size 5, 75 mm	1/ea
MDS6300206	Size 6, 80 mm	1/ea
MDS6300207	Size 7, 85 mm	1/ea



Pessary Accessories

Item No.	Description	Pkg
MDS6303000	OS Cervical Dilator Set, Reusable	1/ea
MDS6303001	Large Dilator Set	1/ea
MDS6303002	Medium Dilator Set	1/ea
MDS6303003	Small Dilator Set	1/ea
MDS6303010	Endometrial Biopsy Cannula, Easy Flex	25/bx
MDS6303012	Endometrial Biopsy Cannula, Hard Flex	25/bx
MDS6303015	Endometrial Sampling Kit, with Curette and Syringe	10/ea
MDS6305001	OS Locator	25/bx
MDS6305002	Canal Dilator	25/bx
MDS6305003	Fundus Sound™	25/bx
MDS6306000	Uterine Sound	25/bx
MDS6309001	Pessary Remover	1/ea
MDS6310000	Pessary Fitting Set	1/ea



MIL303015

Integra Miltex Products


Item No.	Description	Pkg
Speculums		
MIL3010	Graves Vaginal, Small	1/ea
MIL3015	Graves Vaginal, Medium	1/ea
MIL3020	Graves Vaginal, Large	1/ea
MIL3050	Pederson Vaginal, Small	1/ea
MIL3055	Pederson Vaginal, Medium	1/ea
MIL3060	Pederson Vaginal, Large	1/ea
MIL3062	Pederson Vaginal, X-Narrow	1/ea
MILV93060	Vantage Pederson, Large	1/ea
MILV93061	Vantage Pederson, X-Large	1/ea
MIL301351	MX Kogan Impr Endo	1/ea



MIL3010



MIL301382



MIL30650

MIL303010




MIL30965

Item No.	Description	Pkg
Forceps		
MIL7600	Foerster, 9.5"/24.1 cm, STR SER	1/ea
MIL7622	Bozeman, 10.25"/26 cm CVD	1/ea
MIL30965	Schroeder-Braun Tenaculum, 9.75"/24.8 cm	1/ea
MIL301482	K-Y Biopsy, 9.5"/24.1 cm, 3.5 x 8	1/ea
MIL301442WL	Tischler Biopsy	1/ea
MIL301443WL	Burke Biopsy	1/ea

Item No.	Description	Pkg
Curette		
MIL301382	K-Y with Basket	1/ea

Item No.	Description	Pkg
Endo Biopsy Cannulas		
MIL303010	SoftFlex™	25/bx
MIL303012	FirmFlex™	25/bx
Endometrial Sampling Set		
MIL303015	Endometrial Sampling Set	10/pk
Uterine Sounds		
MIL306000	Disposable, 9.75"/24.8 cm	25/pk
MIL30650	SIMS CM 13	1/ea



MIL19102



MIL19284



MIL19314



Integra Miltex Products

Item No.	Description	Pkg
Tuning Forks		
MIL19102	C-128	1/ea
MIL19104	C-256	1/ea
MIL19106	C-512	1/ea
Buck Curettes		
MIL19282	6.5"/16.5 cm ST BL 0	1/ea
MIL19284	6.5"/16.5 cm ST BL 1	1/ea
MIL19286	6.5"/16.5 cm ST BL 2	1/ea
Billeau Ear Loops		
MIL19314	Small Flex	1/ea
MIL19316	Medium Flex	1/ea



MIL19385



MIL19370



MIL19322



MIL19582

Item No.	Description	Pkg
Forceps		
MIL6197	Gruenwald/Jansen, 6.5"/16.5 cm	1/ea
MIL19370	Lucae Ear, 5.5"/14 cm Bayo	1/ea
MIL19371	Jansen, 6.25"/15.9 cm Bayo	1/ea
MIL19434	Hartmann-Noyes, 3.25"/8.3 cm	1/ea
MIL192100	Micro Alligator Ear, Serrated, 3.25"/8.3 cm	1/ea
Ear Syringe		
MIL19385	4-oz./113.4 g with Shield	1/ea

Item No.	Description	Pkg
Baron Suction Tubes		
MIL19580	3 Fr	1/ea
MIL19582	5 Fr	1/ea
MIL19584	7 Fr	1/ea
Disposable Ear Curettes		
MIL19320	Ring Tip	50/bx
MIL19321BX	Loop Tip	50/bx
MIL19322	Small Spoon Tip	50/bx
MIL19323	Large Spoon Tip	50/bx



MDS4406001



MDS4418001



MDS4420614



MDS4426108



Instruments

Item No.	Description	Pkg
Tuning Forks		
MDS4406001	Aluminum, C 128	1/ea
MDS4407202	Hartmann, C1, 256	1/ea
MDS4406002	Aluminum, C, 256	1/ea

Item No.	Description	Pkg
Ear Curettes		
MDS7412803	Fox, 3 mm, 5-5/8"/14.3 cm	1/ea
MDS7412804	Fox, 4 mm, 5-5/8"/14.3 cm	1/ea
MDS4418002	Loop, Billeau, Medium, 6.5"/16.5 cm	1/ea
MDS4418001	Loop, Billeau, Small, 6.5"/16.5 cm	1/ea

Item No.	Description	Pkg
Forceps		
MDS4615720	Dressing, Gruenwald (Jansen), 8"/20.3 cm	1/ea
MDS4420614	Dressing, ENT, Lucae, Serrated, 5.5"/14 cm	1/ea
MDS4426108	Ear, Hartmann, Alligator, 3.5"/8.9 cm	1/ea

Item No.	Description	Pkg
Suction Tubes		
MDS4410203	Baron, 3 Fr, 4"/10.2 cm	1/ea
MDS4410205	Baron, 5 Fr, 4"/10.2 cm	1/ea
MDS0440317	Baron, 7 Fr, 4"/10.2 cm	1/ea



MDS9284010



Crochet Phlebectomy Hooks

Item No.	Description	Pkg
MDS9284001	Style 00, 3.00 mm, 6"/15 cm	1/ea
MDS9284002	Style 0, 2.75 mm, 6"/15 cm	1/ea
MDS9284003	Style 1, 2.60 mm, 6"/15 cm	1/ea
MDS9284004	Style 2, 2.45 mm, 6"/15 cm	1/ea
MDS9284005	Style 3, 2.30 mm, 6"/15 cm	1/ea
MDS9284006	Style 4, 2.15 mm, 6"/15 cm	1/ea
MDS9284007	Style 5, 2.00 mm, 6"/15 cm	1/ea
MDS9284008	Style 6, 1.85 mm, 6"/15 cm	1/ea
MDS9284009	Style 7, 1.70 mm, 6"/15 cm	1/ea
MDS9284010	Style 8, 1.55 mm, 6"/15 cm	1/ea



MDS9283101



Mueller-Style Phlebectomy Hooks

Item No.	Description	Pkg
Right		
MDS9283101	#1, 5"/12.5 cm	1/ea
MDS9283102	#2, 5"/12.5 cm	1/ea
MDS9283103	#3, 5"/12.5 cm	1/ea
MDS9283104	#4, 5"/12.5 cm	1/ea
Left		
MDS9283201	#1, 5"/12.5 cm	1/ea
MDS9283202	#2, 5"/12.5 cm	1/ea
MDS9283203	#3, 5"/12.5 cm	1/ea
MDS9283204	#4, 5"/12.5 cm	1/ea



MDS1821413



Frazier Dura Hook

Item No.	Description	Pkg
MDS1821413	Hook, Sharp, 3 mm, 5"/13 cm	1/ea



MDS9283501



Ramelet Phlebectomy Hooks

Item No.	Description	Pkg
Right		
MDS9283501	#1, 4"/10.2 cm	1/ea
MDS9283502	#2, 4"/10.2 cm	1/ea
MDS9283301	#1, 5"/12.2 cm	1/ea
MDS9283302	#2, 5"/12.2 cm	1/ea
MDS9283401	#3, 5"/12.2 cm	1/ea
MDS9283402	#4, 5"/12.2 cm	1/ea

Item No.	Description	Pkg
Left		
MDS9283601	#1, 4"/10.2 cm	1/ea
MDS9283602	#2, 4"/10.2 cm	1/ea

Item No.	Description	Pkg
Ramelet Phlebectomy Clamp		
MDS92842	5"/12.2 cm	1/ea



MDS9283801



Oesch-Style Phlebectomy Hooks

Item No.	Description	Pkg
Right		
MDS9283701	#1, 6.5"/16.5 cm	1/ea
MDS9283702	#2, 6.5"/16.5 cm	1/ea
MDS9283703	#3, 6.5"/16.5 cm	1/ea
MDS9283704	#4, 6.5"/16.5 cm	1/ea

Item No.	Description	Pkg
Left		
MDS9283801	#1, 6.5"/16.5 cm	1/ea
MDS9283802	#2, 6.5"/16.5 cm	1/ea
MDS9283803	#3, 6.5"/16.5 cm	1/ea
MDS9283804	#4, 6.5"/16.5 cm	1/ea



MDS1821601



Nasal Tenacula Joseph Hooks

Item No.	Description	Pkg
MDS1821601	Delicate, 6.25"/15.9 cm	1/ea
MDS1821602	Delicate, 2-Prong, 2 mm, 6.25"/15.9 cm	1/ea

MDS9284901



Oesch Pin-Strippers

Item No.	Description	Pkg
MDS9284801	Parva, 1.8 mm, 300 mm	1/ea
MDS9284802	Magna, 2.0 mm, 520 mm	1/ea
MDS9284803	Retriever, 2.0 mm, 520 mm	1/ea
MDS9284804	Parva, 2.0 mm, 520 mm	1/ea
MDS9284901	With Rinse Catheter, 60 mm	1/ea

MDS9283004



Varady Phlebectomy Extractors

Item No.	Description	Pkg
MDS9283001	Dissecting Spatula, 1.9 and 2.9 mm, 7"/17.8 cm	1/ea
MDS9283002	Large Hook with Ball Tip and Spatula, 1.9 and 2.8 mm, 7"/17.8 cm	1/ea
MDS9283003	Medium Crochet-Style Hook and Spatula, 1.7 and 1.8 mm, 7"/17.8 cm	1/ea
MDS9283004	Micro Spatula, 0.7 and 1.2 mm, 6.75"/17.1 cm	1/ea
MDS9283005	Micro Hook and Spatula, 0.7 and 1.3 mm, 6.75"/17.1 cm	1/ea
MDS9283006	Small Hook and Spatula, 1.1 and 1.5 mm, 6.75"/17.1 cm	1/ea

MDS1821302



Kleinert-Kutz Skin Hooks

Item No.	Description	Pkg
MDS1821302	Hook, Size 2, 6.25"/15.9 cm	1/ea



MDS15010 Size 10 MDS15210 Size 10 MDS15110L Size 11



Carbon and Stainless Steel Blades

- » Carbon steel blades made from premium Swedish steel
- » Engineered to fit all handles
- » Color coded packaging
- » ISO 7740 certified

Item No.	Description	Pkg
----------	-------------	-----

Carbon Steel Blades

MDS15010	Size 10	100/bx
MDS15011	Size 11	100/bx
MDS15012	Size 12	100/bx
MDS15015	Size 15	100/bx
MDS15020	Size 20	100/bx
MDS15021	Size 21	100/bx
MDS15022	Size 22	100/bx
MDS15023	Size 23	100/bx

Carbon Steel Twin-Back Blades

MDS15010	Size 10	150/bx
MDS15011	Size 11	150/bx
MDS15015	Size 15	150/bx
MDS15020L	Size 20	150/bx



MDS0610301



Scalpel Knife Handle

Item No.	Description	Pkg
MDS0610301	#3, Graduated	1/ea
MDS0610300	#3, Non Graduated	1/ea
MDS0610400	#4	1/ea
MDS0610700	#7	1/ea
MDS0610116	Corwin	1/ea



MDS15210 MDS15215 MDS15410 BLDGLV3



Surgical Blades and Disposable Scalpels

- » Premium Swedish steel blades ergonomically designed for better feel and precision.
- » ABS plastic handle to ensure greater strength
- » Graduated scale on handle

Item No.	Description	Pkg
----------	-------------	-----

Stainless Steel Blades

MDS15110	Size 10	100/bx
MDS15111	Size 11	100/bx
MDS15112	Size 12	100/bx
MDS15115	Size 15	100/bx
MDS15120	Size 20	100/bx
MDS15121	Size 21	100/bx
MDS15122	Size 22	100/bx

Stainless Steel Wide Foil Blades

MDS15110L	Size 10	150/bx
MDS15111L	Size 11	150/bx
MDS15115L	Size 15	150/bx

Disposable Scalpels

MDS15210	Size 10	10/bx
MDS15211	Size 11	10/bx
MDS15212	Size 12	10/bx
MDS15215	Size 15	10/bx
MDS15220	Size 20	10/bx
MDS15221	Size 21	10/bx

Safety Scalpels


MDS15310	Size 10	10/bx
MDS15311	Size 11	10/bx
MDS15315	Size 15	10/bx

Blade Glove Cartridges

MDS15410	Size 10	50/bx
MDS15411	Size 11	50/bx
MDS15415	Size 15	50/bx

Blade Glove Metal Handles

BLDGLV3	Size 3	5/bx
---------	--------	------



ASC720001Z



DermaBlade® Shave Biopsy

- » Disposable blade to easily remove skin protuberances and deep lesions while reducing ragged skin edges
- » MicroCoat® surgical blade
- » Three digit sure-grip teeth

Item No.	Description	Pkg
ASC720001	Shave Biopsy	400 ea/cs
ASC720001Z	Shave Biopsy	50/bx



DYNJSM01



Sterile Surgical Skin Marker

- » Medical-grade, non-irritating violet ink
- » Large ink reservoir for longer product life
- » Constant ink flow for clear marking
- » Rule included

Item No.	Description	Pkg
DYNJSM01	6"/15 cm	50/cs
DYNJSM01H	6"/15 cm	1/ea
DYNJSM11	Prep-Resistant Skin Marker	100/cs




STAPLER35W



Disposable Surgical Skin Stapler

- » Includes a counter to indicate the exact number of staples used
- » Contains 35 regular or wide medical grade stainless steel staples

Item No.	Description	Pkg
STAPLER35RB	Regular, Sterile	22/cs
STAPLER35W	Wide, with Counter, Sterile	6/cs
STAPLER35WA	Wide	6/cs
STAPLER35WPB	Wide, Sterile	22/cs



DYNJ07900H



Incision and Drainage Tray

- » Use these trays to standardize product usage and control supply costs
- » Tray includes:
 - 8-ply 4 x 4-inches/10 x 10 cm Gauze Dressings, 4/ea
 - Fenestrated Drape, 1/ea
 - 13 x 19-inches/33 x 48 cm White Polybacked Towel, 1/ea
 - Kelly Forceps Straight, 1/ea
 - Adson Thumb Forceps Serrated, 1/ea
 - Scalpel No. 11, 1/ea
 - Alcohol Prep Pad, 1/ea
 - PVP Swabsticks and Polybag with Red Twist-Tie in Ridged Tray, 3/pk

Item No.	Description	Pkg
DYNJ07900H	Incision and Drainage Tray	1/ea



MDS708555H

DYNJ07254AH



Suture Removal Tray

- » Suture Removal Trays help you run your procedure more efficiently
- » Helps reduce the risk of cross contamination

Item No.	Description	Pkg
MDS708555H	Economy	1/ea
DYNJ07254AH	Floor Grade	1/ea



Medline is the No. 1 provider of several major products to the U.S. healthcare market, including exam room kits.

Contact us to learn more, and to receive the Solutions for Minor Procedures catalog for more options.



Laceration Tray

- » All the necessary components conveniently assembled in a compact, sterile package to help you maximize efficiency in the exam room
- » Tray includes:
 - 2-Compartment Tray, 1/ea
 - 8-ply, 2 x 2"/5 x 5 cm Gauze, 5/ea
 - 8-ply, 4 x 4"/10 x 10 cm Gauze, 10/ea
 - Blue Blotting Towels 19.5 x 21"/49.5 x 53.3 cm, 2/ea
 - Fenestrated Drape, 1/ea
 - Adson Thumb Forceps 1 x 2 Teeth, 1/ea
 - Comfort Loop Straight Iris Scissors, 1/ea

- Comfort Loop Curved Mosquito Forceps, 1/ea
- Comfort Loop Webster Needle Holder, 1/ea
- 10 ml Syringe, 1/ea
- 18 g x 1.5"/3.8 cm Needle, 1/ea
- 25 g x 0.625"/1.5 cm Needle, 1/ea
- 27 g x 0.5"/1.2 cm Needle, 1/ea
- 2-oz./59 ml Medicine Cups, 2/ea

Item No.	Description	Pkg
DYNJ03000	Laceration Tray	16/cs
DYNJ03000H	Laceration Tray	1/ea



SDCPP012



PassPort® Plus Sterilizer Monitoring Service

- » Spore testing with one test and one control strip
- » Steam test results documented after 24 hours
- » Test reports sent on a quarterly basis
- » Certified for use in Dry Heat, Chemical Vapor, and EO Gas Processes

Item No.	Description	Pkg
SDCPP012	Sterilizer Monitoring Service	12/cs



MDS200590



Load Monitoring

Item No.	Description	Pkg
MDS200565	Biological Indicator, Steam	1/ea
MDS200595	Biological Test Pack, Steam	1/ea
MDS200590	Biological Monitoring, Smart-Read	1/ea



Decontamination

- » Quality solutions and detergents recommended for use with our instruments
- » Fast, high level disinfection for use with reusable medical devices



MDS880PA28

Item No.	Description	Pkg
MDS8800FOAMH	Foam, 24-oz.	1/ea
MDS8800B101H	Bio-Zolve, 32-oz.	1/ea

Detergents

MDS88000B91H	Enzymatic, 1 Gallon	1/ea
MDS88000B1H	Neutral, 1 Gallon	1/ea

Lubricants

MDS8800T32H	For Instruments, 32-oz.	1/ea
MDS8800T2	For Instruments, 1 Gallon	1/ea

Disinfectants

MDS88GLUT28MH	Gluteraldehyde	1/ea
MDS8800PA28H	OPA, 4 x 1 Gallon	1/ea



MDS88000B1H



MDS200615



Surgical Instrument Sterilization Steam Indicator Strips

- » Measure the time, temperature and presence of steam in steam sterilizers

Item No.	Description	Pkg
MDS200100	Indicator Strip, 4"/10 cm	250/bx
MDS200615	Class 5 Integrator, Migrating	1,000/bg



MDS88INSTWP

Stain-Away Instrument Wipes

- » Remove rust and stains to restore instrument luster
- » Safe for all stainless steel or tungsten carbide instruments

Item No.	Description	Pkg
MDS88INSTWP	Stain Away	6 ea/cs



MDS88BRUSH



Surgical Instrument Cleaning Brushes

- » During manual instrument cleaning effectively remove debris from hard-to-reach areas

Item No.	Description	Pkg
MDS88BRUSHSS	7"/17.8 cm Cleaning Brush	3/pk
MDS88BRUSH	7"/17.8 cm Nylon Bristles	3/pk



MDS200836H



Surgical Instrument Sterilization Steam Autoclave Tape

- » Dual-purpose adhesive tape seals the sterilization pack, then visually indicates that pack has been exposed to the steam sterilization cycle
- » Distinctive color change from beige to black shows that it is not reversible

Item No.	Description	Pkg
MDS200848H	0.75" x 60-yds./ 1.9 cm x 54.86 m, Roll	1/ea
MDS200836H	1" x 60-yds./ 2.5 cm x 54.86 m, Roll	1/ea
MDS200846H	1" x 60-yds./ 2.5 cm x 54.86 m, Blue	1/ea



Pouches

- » Pouches are transparent and have easy-to-read process indicators, with distinct color changes
- » Constructed from strong, medical grade materials that easily withstand heavy instrumentation
- » Pre-fold feature on all self-seal pouches means a quick, easy and accurate fold every time



Sterilization tubing/Rolls and Pouches



MPP200215

Item No.	Description	Pkg
Self-Seal Pouch		
MPP100520GS	3.5 x 9"/8.9 x 22.9 cm	2,000/cs
MPP100525GS	5.25 x 10"/13.3 x 25.4 cm	2,000/cs
MPP100545GS	7.5 x 13"/19 x 33 cm	1,000/cs
MPP100546GS	8 x 16"/20.3 x 40.6 cm	1,000/cs
MPP100550GS	5 x 15"/12.7 x 38.1 cm	1,000/cs
MPP100555GS	12 x 15"/30.5 x 38.1 cm	1,000/cs
MPP100558GS	12 x 18"/30.5 x 45.7 cm	1,000/cs
MPP100565GS	3.5 x 22"/8.9 x 55.9 cm	1,000/cs

Heat Seal Pouch

MPP200215	4 x 8"/10.2 x 20.3 cm	2,000/cs
MPP200225	6 x 10"/15.2 x 25.4 cm	2,000/cs
MPP200245	7.5 x 13"/19 x 33 cm	1,000/cs
MPP200246	8 x 16"/20.3 x 40.6 cm	1,000/cs
MPP200250	5 x 15"/12.7 x 38.1 cm	1,000/cs
MPP200255	12 x 15"/30.5 x 38.1 cm	1,000/cs
MPP200258	12 x 18"/30.5 x 45.7 cm	1,000/cs
MPP200265	3.5 x 22"/8.9 x 55.9 cm	1,000/cs
MPP200270	6 x 22"/15.2 x 55.9 cm	2,000/cs
MPP200275	18 x 22"/45.7 x 55.9 cm	1,000/cs

Item No.	Description	Pkg
Sterilization Tubing/Roll		
MPP200410S	3" x 100'/7.6 cm x 30.48 m	8/cs
MPP200410R	3" x 656'/7.6 cm x 200 m	1/cs
MPP200415S	4" x 100'/10.2 cm x 30.48 m	16/cs
MPP200415R	4" x 656'/10.2 cm x 200 m	1/cs
MPP200420S	6" x 100'/15.2 cm x 30.48 m	8/cs
MPP200420R	6" x 656'/15.2 cm x 200 m	1/cs
MPP200422S	8" x 100'/20.3 cm x 30.48 m	4/cs
MPP200422R	8" x 656'/20.3 cm x 200 m	1/cs
MPP200425S	10" x 100'/25.4 cm x 30.48 m	4/cs
MPP200425R	10" x 656'/25.4 cm x 200 m	1/cs

Tyvek Pouch, Self-Seal

MPP100650N	3.5 x 9"/8.9 x 22.9 cm	2,000/cs
MPP100655N	5.25 x 10.25"/13.3 x 26 cm	2,000/cs
MPP100660N	7.5 x 13"/19 x 33 cm	1,000/cs
MPP100670N	12 x 18"/30.5 x 45.7 cm	1,000/cs

Tyvek Tubing

MPP100680N	6" x 100'/15.2 cm x 30.48 m	8/cs
MPP100685N	10" x 100'/25.4 cm x 30.48 m	4/cs



PATIENT CARE

Liquids and Wet Pack Products	94-96
Applicators and Tongue Blades	96-97
Hot/Cold Therapy.....	97-98
Vaginal Speculums.....	98
Oxygen Delivery	99
Resuscitation	99
Nebulizers	100
Pulse Oximetry	101
Handheld Oximeters	101
Catheters and Drain Bags.....	102-105
Specimen Containers.....	105-106
Disposable Drapes	106
Personal Care	106-107
Sutures	107-108
Pain Management	108
Regional Anesthesia and Pain Management	109-111
Cryosurgical Systems	111
Wash Basins.....	112
Ear Syringes.....	112
Splints.....	112
Cast Padding	113
Stockinettes.....	113
Casting Tape.....	113
Compression Stockings	114
Anti-Embolism Stockings.....	115
Arm Slings and Cervical Collars	115
Wrist & Forearm Splints.....	116
Ankle & Foot Supports	116
Post-Op Shoes.....	116
Walkers	117
Podiatry Products	118

MDS090735Z



Alcohol Prep Pads and Swabsticks

- » For removal of oils and residue from skin
- » Soft, absorbent, non-woven pads and swabsticks
- » Ideal for use prior to injections
- » Saturated with 70% isopropyl alcohol, USP
- » This product is intended for use on the body only

Item No.	Description	Pkg
MDS090735Z	Medium, 2-ply Sterile Alcohol Prep Pad	200/bx
MDS090670Z	Large, 2-ply Sterile Alcohol Prep Pad	100/bx
MDS093810Z	Alcohol Swabsticks, Triple Pack	25 pk/bx



MDS098003ZH



Isopropyl Rubbing Alcohol

- » Help prevent infections in minor cuts, scrapes and burns
- » Can be used in small amounts on affected areas one to three times daily
- » Once alcohol dries, wounds may be covered with sterile bandages
- » 70% isopropyl alcohol, USP

Item No.	Description	Pkg
MDS098011H	Bottle: 4-oz./118 ml	1/ea
MDS098003Z	Bottle: 16-oz./473 ml	12/cs
MDS098003ZH	Bottle: 16-oz./473 ml	1/ea
MDS098017	Bottle: 32-oz./946 ml	12/cs
MDS098017H	Bottle: 32-oz./946 ml	1/ea



MDS098001H



Hydrogen Peroxide

- » USP grade 3% hydrogen peroxide is a multi-use first aid antiseptic and oral debriding agent
- » Ideal for prepping skin prior to administering injections and cleaning minor cuts, scrapes and burns
- » This product is intended for use on the body only

Item No.	Description	Pkg
MDS098014H	Bottle: 4-oz./118 ml	1/ea
MDS098014	Bottle: 4-oz./118 ml	24/cs
MDS098015H	Bottle: 8-oz./237 ml	1/ea
MDS098015	Bottle: 8-oz./237 ml	12/cs
MDS098001H	Bottle: 16-oz./473 ml	1/ea
MDS098001Z	Bottle: 16-oz./473 ml	12/cs
MDS098002	1-gal./3.8 L	4/cs



MDS090780H



Nail Polish Remover Pads

- » Effectively remove nail polish when patient monitoring of nail beds is necessary
- » Acetone-Free

Item No.	Description	Pkg
MDS090780H	Nail Polish Remover Pads	100/bx
MDS090780	Nail Polish Remover Pads	1,000/cs



MDS090855H



Adhesive Tape Remover Pads

- » Painlessly and effectively remove barrier films and residue left behind by adhesive tape
- » Pads help reduce patient discomfort while contents help soften the skin

Item No.	Description	Pkg
MDS090855H	Adhesive Tape Remover Pads	100/bx
MDS090855	Adhesive Tape Remover Pads	1,000/cs



MDS094181H



Cleansing Towelettes

- » For quick, general cleaning in both healthcare and personal use

Item No.	Description	Pkg
MDS094181H	Benzalkonium Chloride Towelettes, 0.4% BZK, Contains No Alcohol	100/bx
MDS094184H	Cleansing Towelettes, 0.13% BZK, with Alcohol	100/bx
MDS094186Z	Obstetrical Cleansing Towelettes, 0.4% BZK, Contains No Alcohol	100/bx
MDS094188H	Wash-Up Towelettes	100/bx



MDS093902ZZ



Povidone-Iodine Prep Pads and Swabsticks

- » Broad-range, first aid antimicrobial* that helps prevent infection in minor cuts, scrapes and burns
- » Each pad or swabstick is saturated with a 10% povidone-iodine, USP solution

Item No.	Description	Pkg
MDS093917H	Povidone Iodine Prep Pad, Medium	100/bx
MDS093901ZZ	Povidone Iodine Swabstick, Single	50 pk/bx
MDS093902ZZ	Povidone Iodine Swabstick, Triple	25 pk/bx

*This product is intended for use on the body only.



DYND371603Z



Surgical Scrub Dry Brushes

- » Helps to achieve maximum cleaning with minimal skin irritation
- » Flexible design fits comfortably in most hands
- » Nail cleaner included

Item No.	Description	Pkg
DYND371603Z	Dry Brushes	30 ea/bx



MDS093944H



Povidone-Iodine Solutions

- » Essential for surgical and other procedures
- » Effective for skin and mucous membrane preparation
- » Use as surgical scrubs or first aid antiseptic
- » Prep Solution is 10% weight by volume povidone-iodine, USP solution equivalent to 1% available iodine. Product is intended for use on the body only
- » Antimicrobial Scrub Solution with surfactants is 7.5% weight by volume povidone-iodine, USP solution equivalent to 0.75% available iodine. Product is intended for use on the body only

Item No.	Description	Pkg
Povidone-Iodine Prep Solution		
MDS093940HH	2-fl. oz./60 ml	1/ea
MDS093944H	4-fl. oz./118 ml	1/ea
MDS093943H	8-fl. oz./237 ml	1/ea
MDS093906H	1-pint/473 ml	1/ea
MDS093941H	1-qt./946 ml	1/ea
MDS093904H	1-gal./3.8 L	1/ea

Povidone-Iodine Antimicrobial* Scrub Solution

MDS093945H	4-fl. oz./118 ml	1/ea
MDS093947H	1-pint/473 ml	1/ea
MDS093946H	1-qt./946 ml	1/ea

*This product is intended for use on the body only.

Avoid contamination. Don't forget about Hand Hygiene products on pages 74–77.



MDS032290H



Lubricating Jelly, Sterile

- » Facilitates the entry of diagnostic and therapeutic devices
- » Spreads easily over rubber, metal or synthetic surfaces to form transparent lubricating film
- » Adheres well, yet may be washed away easily
- » Greaseless, water-soluble

Item No.	Description	Pkg
MDS032273H	Foil Pack, 2.7 g	144/bx
MDS032280Z	Foil Pack, 5 g	150/bx
MDS032285H	Flip-top tube, 2-oz.	1/ea
MDS032285Z	Flip-top Tube, 2-oz.	12/bx
MDS032290H	Flip-top tube, 4-oz.	1/ea
MDS032290Z	Flip-top Tube, 4-oz.	12/bx



DAY020536H



HR Pharmaceuticals, Inc.

Surgilube® Lubricating Jelly

- » Sterile, easy-to-use lubrication with bacteriostatic blend of natural water-soluble gums
- » Provides persistent lubrication for easy insertion of catheters, endoscopes, surgical instruments, gloves, finger cots, and more
- » Forms a transparent, crystal clear film that won't cloud optics
- » Spreads readily, adheres well, washes off easily, and is non-irritating
- » Expiration date is stamped on the end of each tube

Item No.	Description	Pkg
DAY020543Z	Foil Packs: 3 g	144/bx
DAY020536H	Tube: 4.25-oz./125.7 ml	1/ea



Ultrasound Gel

- » An excellent conductive agent for use in therapeutic ultrasound, diagnostic imaging and other external procedures
- » Water soluble, non-corrosive, greaseless and non-staining

Item No.	Description	Pkg
MDS092005H	Blue, Squeeze Bottle, Twist-top, 8.5-oz.	1/ea
MDS092010H	Blue, Jug, 5 L	1/ea
CTR001305H	Clear, Squeeze Bottle, Flip-top, 8.5-oz.	1/ea
CTR001201	Clear Foil Packet, 20 g	100/cs



Plastic Zip Closure Bags with White Write-On Block

- » Reclosable bags are ideal for clean and convenient storage
- » Feature a white block for clear labeling
- » Bags are made out of a durable 2 mil plastic
- » Available in a variety of sizes

Item No.	Description	Pkg
NONZIP35	3 x 5"/7.62 x 12.7 cm, 2 mil	1,000/cs
NONZIP35Z	3 x 5"/7.62 x 12.7 cm, 2 mil	100/pk
NONZIP48	4 x 8"/10.16 x 20.32 cm, 2 mil	1,000/cs
NONZIP58	5 x 8"/12.7 x 20.32 cm, 2 mil	1,000/cs
NONZIP58Z	5 x 8"/12.7 x 20.32 cm, 2 mil	100/pk
NONZIP69	6 x 9"/15.24 x 22.86 cm, 2 mil	1,000/cs
NONZIP69Z	6 x 9"/15.24 x 22.86 cm, 2 mil	100/pk
NONZIP810	8 x 10"/20.32 x 25.4 cm, 2 mil	1,000/cs
NONZIP810Z	8 x 10"/20.32 x 25.4 cm, 2 mil	100/pk
NONZIP912	9 x 12"/22.86 x 30.48 cm, 2 mil	1,000/cs
NONZIP912Z	9 x 12"/22.86 x 30.48 cm, 2 mil	100/pk
NONZIP1215	12 x 15"/30.48 x 38.1 cm, 2 mil	1,000/cs
NONZIP1215Z	12 x 15"/30.48 x 38.1 cm, 2 mil	100/pk



MDS202050

Cotton-Tipped Applicators

- » Non-sterile products are available in autoclavable bags for sterile applications

Item No.	Description	Pkg
Non-Sterile Cotton-Tipped Applicator		
MDS202050Z	Wood Shaft, Non-Sterile, 3"/7.6 cm, 100s	100/bg
MDS202055Z	Wood Shaft, Non-Sterile, 6"/15.2 cm, 100s	1,000/bx
MDS202055H	Wood Shaft, Non-Sterile, 6"/15.2 cm, 100s	100/bg

Sterile Cotton-Tipped Applicators

MDS202000Z	Wood Shaft, Sterile, 6"/15.2 cm, 2s	200/bx
MDS202010Z	Wood Shaft, Sterile, 6"/15.2 cm, 10s	200/bx
MDS202010H	Wood Shaft, Sterile, 6"/15.2 cm, 10s	10/bg
MDS202105Z	Wood Shaft, Sterile, 3"/7.6 cm, 2s	200/bx
MDS202095Z	Plastic Shaft, Sterile, 6"/15.2 cm, 2s	200/bx



MDS21460Z

Cotton Balls

- » Soft and highly absorbent 100% cotton
- » Compact and uniform in size and weight
- » Will not come apart when saturated
- » Non-sterile

Item No.	Description	Pkg
MDS21460Z	Medium, 1"/3 cm Diameter	2,000/bg
MDS21462H	Large, 1 1/4"/3 cm Diameter	1,000/bg
MDS21463H	Large, 1 1/4"/3 cm Diameter	50/bg



MDS202090

Specialty Swabs

- » Feature a rayon tip with a plastic shaft
- » Non-sterile

Item No.	Description	Pkg
MDS202090Z	OB/GYN Swab, Rayon-Tip, 8"/20.3 cm	100/bx
MDS202085Z	Procto Swab, Rayon-Tip, 16"/40.6 cm	100/bx



MDS202075



Tongue Depressors

- » Made of high-quality wood, smooth and splinter-free
- » Useful for spreading ointments or stirring medications
- » Sterile and non-sterile available in two different sizes

Item No.	Description	Pkg
----------	-------------	-----

Non-Sterile Tongue Depressors

MDS202065H	Tongue Blade, Non-Sterile, 6"/15.2 cm	500/bx
MDS202077Z	Tongue Blade, Non-Sterile, Wrapped, 6"/15.2 cm	250/bx
MDS202070HH	Tongue Blade, Non-Sterile, 5.5"/14 cm	500/bx
MDS202110Z	Infant Tongue Blade, Non-Sterile, 4.5"/11.4 cm	1,000/bx

Sterile Tongue Depressors

MDS202075Z	Tongue Blade, Sterile, 6"/15.2 cm	100/bx
MDS202073Z	Tongue Blade, Sterile, 5.5"/14 cm	100/bx



MDS202060Z



Wooden Applicator Sticks

- » Ideal for specimen transfer use, mixing rods and other standard laboratory applications

Item No.	Description	Pkg
MDS202060Z	Non-Sterile, 6"/15.2 cm	864/bx



MDS202080Z



Wooden Medical Spoons

- » These spoons are excellent as an aid for dispensing wet medication mixtures

Item No.	Description	Pkg
MDS202080Z	Wood Spoon, Non-Sterile, 3.63"/9.2 cm	1,000/bx



MDS138020



Accu-Therm® Reusable Hot/Cold Gel Pack

- » Therapeutic pack made with non-toxic gel
- » Can be microwaved or frozen
- » Cleans easily

Item No.	Description	Pkg
MDS138020	5 x 10"/12.7 x 25.4 cm	16/cs
MDS138020H	5 x 10"/12.7 x 25.4 cm	1/ea



Instant Hot Pack

- » High-performance pack providing instant and longer lasting therapy
- » Can be applied directly to skin

Item No.	Description	Pkg
MDS139008*	6 x 6"/15.2 x 15.2 cm	24/cs
MDS139009*	6 x 9"/15.2 x 22.9 cm	24/cs
MDS139010*	4 x 9"/10.2 x 22.9 cm	24/cs

*Add "H" to item to order by the each.

MABIS

Sundry/Dressing Jars

- » Ideal for storing for exam room necessities, such as bandages, applicators, tongue depressors, cotton and gauze
- » Clear lid allows easy viewing of contents
- » Labels not included



MBH39803000

Item No.	Description	Pkg
MBH39803000	Clear Unlabeled Sundry Jar	1/ea
MBH39802000	Clear Unlabeled Sundry Jar	5/cs



MDS137000



Standard Instant Cold Packs

- » Activates instantly with fold

Item No.	Description	Pkg
MDS137000	6 x 9"/15.2 x 22.9 cm	24/cs
MDS137000H	6 x 9"/15.2 x 22.9 cm	1/ea



MDS148065



Deluxe Perineal Warm Packs

- » Highly absorbent, soothing pad with adhesive for secure placement

Item No.	Description	Pkg
MDS148065*	4.5 x 14.25"/11.2 x 36.2 cm	24/cs

*Add "H" to item to order by the each.



DYND70400



Pap Smear Tray

- » Complete kit contains a Medline speculum, cervical scraper, two 6" polyester-tipped applicators, pair of vinyl gloves, lubricating jelly, two frosted end slides, slide holder and mailer

Item No.	Description	Pkg
DYND70400	Pap Smear Tray	20/cs



DYND70401S



Light Source Adaptable Disposable Vaginal Specula

- » Cost-effective and efficient
- » Single-use speculums reduce the risk of cross contamination and eliminate sterilization expenses
- » Precision-molded plastic provides patient comfort and ease of use during exams and procedures
- » Corded light source is required
- » Non-sterile

Item No.	Description	Pkg
DYND70401S	Small	25/bx, 4 bx/cs, 100/cs
DYND70401SZ	Small	25/bx
DYND70401M	Medium	25/bx, 4 bx/cs, 100/cs
DYND70401MZ	Medium	25/bx
DYND70401L	Large	25/bx, 4 bx/cs, 100/cs
DYND70401LZ	Large	25/bx



WSD908015ZZ



Endocell® Endometrial Cell Sample

Item No.	Description	Pkg
WSD908015ZZ	Cell, Endometrial	35/bx



DYND70420M



Disposable Vaginal Specula

- » Single use helps to reduce the risk of cross contamination
- » Made of a high-strength, clear non-toxic plastic
- » Easily adjustable with one hand and self-holding in various positions
- » Smooth, rounded plastic helps reduce the risk of pinching or binding tissue during procedures
- » Non-sterile

Item No.	Description	Pkg
DYND70420S	Small	10/bg, 10 bg/cs, 100 ea/cs
DYND70420SZ	Small	10/pk
DYND70420M	Medium	10/bg, 10 bg/cs, 100 ea/cs
DYND70420MZ	Medium	10/pk
DYND70420L	Large	10/bg, 10 bg/cs, 100 ea/cs
DYND70420LZ	Large	10/pk



W-A59001



Disposable Vaginal Speculums and Sheaths

Item No.	Description	Pkg
W-A590XS	Extra Small	24/bx, 4 bx/cs, 96 ea/cs
W-A59000	Small	24/bx, 4 bx/cs, 96 ea/cs
W-A59000Z	Small	24/bx
W-A59001	Medium	24/bx, 4 bx/cs, 96 ea/cs
W-A59001Z	Medium	24/bx
W-A59004CSZ	Large	18/bx
W-A59010	KleenSpec Disposable Sheaths	500/cs
W-A80010	Cordless Illuminator for KleenSpec Vaginal Speculums	1/ea



HCS4511B



Standard O₂ Cannulas

- » Economy nasal cannulas
- » Designed for short-term, post-operative patients with straight, non-flared nasal prongs

Item No.	Description	Pkg
HCS4511B	Adult, Straight Prongs 7' Tubing, Standard	50/cs
HCSU4511B	Adult, Straight Prongs 7' Tubing, Universal	50/cs



HCS4514



Soft Touch O₂ Cannulas

- » Soft, curved and non-flared
- » Available in a wide variety of lengths and sizes with crush resistant tubing

Item No.	Description	Pkg
HCS4514	Adult, 7' Tubing, Standard	50/cs
HCS4516B	Adult, 14' Tubing, Standard	50/cs
HCS4518	Pediatric, 7' Tubing, Standard	50/cs
HCS4517	Infant, 7' Tubing, Standard	50/cs



HCSS4514S



Super Soft O₂ Cannulas

- » Made of soft, flexible and non-tacky material
- » Helps reduce pressure and friction on the ears and cheeks

Item No.	Description	Pkg
HCSU4514S	Adult, 7' Tubing, Universal	50/cs
HCSS4514S	Adult, 7' Tubing, Standard	50/cs
HCSU4518S	Pediatric, 7' Tubing, Universal	50/cs
HCSU4517S	Infant, 7' Tubing, Universal	50/cs



HCS4568



CO₂ Cannulas

- » Provide simultaneous oxygen delivery and CO₂ monitoring on both nasal passages
- » Oxygen is delivered by both nasal prongs, while CO₂ is monitored through lower ports
- » Universal O₂ Connector is designed to fit any flowmeter without the need for adaptors

Item No.	Description	Pkg
HCS4561	Adult, 7', 4" Pigtail Female Luer Connector	25/cs
HCS4568	Adult, 7', 7' Male Luer Connector	25/cs
HCS4568F	Adult, 7', 7' Female Luer Connector	25/cs
HCS4565	Adult, 10', 10' Male Luer Connector	25/cs
HCS4566	Adult, 14', 14' Male Luer Connector	25/cs



HCS4630B



O₂ Delivery/Medication Delivery

- » Lightweight
- » Made of soft, comfortable vinyl
- » Available with or without a swivel connector

Item No.	Description	Pkg
HCS4630B	Adult Elongated Aerosol Mask	50/cs
HCS4630	Adult Elongated Aerosol Mask with Swivel Connector	50/cs
HCS4631B	Pediatric Elongated Aerosol Mask	50/cs
HCS4621B	Pediatric Trach Mask	50/cs



CPRM2216



Manual Resuscitation

- » Textured for a non-slip grip
- » Premium soft-cushion mask allows for effective seal on the patients face

Item No.	Description	Pkg
CPRM1116	Adult Bag with Bag Reservoir	6/cs
CPRM1126	Adult Bag with Tube Reservoir	6/cs
CPRM2216	Pediatric Bag with Bag Reservoir	6/cs
CPRM2226	Pediatric Bag with Tube Reservoir	6/cs
CPRM3312	Infant Bag with Bag Reservoir	12/cs
CPRM3322	Infant Bag with Tube Reservoir	12/cs
HCS64184	Mouth to Mask Resuscitator with Valve	12/cs



Handheld Nebulizer Kits

» Medline offers a wide selection of nebulizer kits for asthma medication delivery with optimal particle size for maximum medicine delivery and minimal residual medication

Item No.	Tee Adaptor	Mouthpiece	Reservoir Tube	Tubing Length	Connector	Mask Type	Disposable	Pkg
HCS4485H				7 foot	Standard	Adult	•	1/ea
HCS4483H	22mm OD/ID	•	•	7 foot	Standard		•	1/ea
HCS4486H				7 foot	Standard	Pediatric	•	1/ea

*Remove H for case quantity 50 ea/cs.



HCS4483

HCS4485




HCS70004



HCS70004H



HCS70004RD



Accessories



Compact Nebulizer Compressor with Kit

» Compact and stable with durable hard plastic to sustain continuous use

» Kits include a disposable nebulizer medicine cup, 7 feet of tubing and a mouthpiece.

Item No.	Description	Pkg
HCS70004	With Disposable Nebulizer Kit	8/cs
HCS70004H	With Disposable Nebulizer Kit	1/ea
HCS70004RD	With Disposable & Reusable Nebulizer Kit	8/cs
HCS70004RDH	With Disposable & Reusable Nebulizer Kit	1/ea
HCSFISHNEB	Character Fish with Disposable Nebulizer Kit	8/cs
HCSFISHNEBH	Character Fish with Disposable Nebulizer Kit	1/ea
HCS70004F	Replacement Filters	5/pk



MDI Spacers

» Simple and economical way to improve effectiveness of inhaler

» Anti-static and compatible with any standard mask with a 22 mm fitting

Item No.	Description	Pkg
HCSA100	MDI Spacer Only	10/bx
HCSA113	MDI Spacer with Pediatric Mask	10/bx
HCSA114	MDI Spacer with Adult Mask	10/bx



HCSPF100

Peak Flow Meters

» Help patients monitor their asthma

» Durable for years of reliable performance and individually calibrated for accuracy

Item No.	Description	Pkg
HCSPF100	Portable Peak Flow Meter, Full Range, 60-720 lpm	12/bx
HCSPF200	Standard Peak Flow Meter, Full Range, 60-850 lpm	12/bx

Finger Pulse Oximeters

Get quick and precise SpO₂ and pulse readings in seconds with our high-quality fingertip pulse oximeters

- » Small and lightweight, these clip sensors automatically take readings when a finger is inserted
- » They fit a wide range of finger sizes, have a long battery life and are solidly constructed for worry-free use
- » All models include AAA batteries, lanyard and warranty



HCSM70



HCSM70P



HCSFSATSPORT



Easy-Grip Finger Pulse Oximeter

- » Hinges adjust to comfortably accommodate larger fingers
- » SpO₂, pulse rate, pulse bar, plethysmograph
- » Backlit LED display with four viewing modes
- » Slip-resistant grip
- » Hinged battery cover
- » 2-year warranty

Item No.	Description	Pkg
HCSM70	Fingerclip, Easy-Grip	1/ea

Pediatric Finger Pulse Oximeter

- » Child-friendly frog design makes kids smile
- » SpO₂, pulse rate, pulse bar, plethysmograph
- » Patented multi-directional display has 6 viewing modes
- » Bright, digital OLED with adjustable readout sizes
- » Soft padding for patient comfort
- » Hinged battery cover
- » 2-year warranty

Item No.	Description	Pkg
HCSM70P	Fingerclip, Pediatric	1/ea

High-Impact Finger Pulse Oximeter

- » Rugged casing withstands accidental drops
- » SpO₂, pulse rate and pulse bar readings
- » Multi-directional display with full-color OLED screen
- » Durable silicone boot construction resists cracking and chipping
- » Hinged battery cover
- » 2-year warranty

Item No.	Description	Pkg
HCSFSATSPORT	Fingertip, Fingersat	1/ea

For a full catalog of Respiratory product options, please visit Medline.com or contact your Medline Representative




MD400 Handheld Oximeters

- » SpO₂, pulse rate, pulse bar, waveform
- » Adjustable audible and visual alarms
- » Two display modes
- » Comfortable finger clip probe
- » Data transfer to computer via USB cable
- » Low battery indicator

Item No.	Description	Pkg
HCSMD400	Handheld Continuous Economy Pulse Oximeter with Adult Finger Probe	1/ea
HCSMD400A	Replacement Probe, Adult, for HCSMD400	1/ea
HCSMD400P	Replacement Probe, Pediatric, for HCSMD400	1/ea



HCSMD400



DYND11216



2-Way Coude-Tip Silicone-Elastomer Coated Latex Foley Catheters

- » Coating helps resist encrustation and aids in patient comfort
- » Large, smooth drainage eyes for maximum comfort
- » 10 ml balloon size
- » For difficult insertions, coude-tip Foley catheters feature a bent tip that helps maneuver around obstructions leading to the bladder

Item No.	Description	Pkg
DYND11214*	14 Fr	12/cs
DYND11216*	16 Fr	12/cs
DYND11218*	18 Fr	12/cs
DYND11220*	20 Fr	12/cs

*Add "H" to item to order by the each.



DYND11756



Silicone-Elastomer Coated Latex Foley Catheter

- » Innovative manufacturing process provides the maximum amount of silicone coating available
- » Silicone-elastomer coating resists encrustation
- » Sterile

Item No.	Description	Pkg
DYND11752*	12 Fr	12/cs
DYND11754*	14 Fr	12/cs
DYND11756*	16 Fr	12/cs
DYND11758*	18 Fr	12/cs
DYND11760*	20 Fr	12/cs
DYND11762*	22 Fr	12/cs

*Add "H" to item to order by the each.



Clear Vinyl Intermittent Catheters

- » Manufactured with a smooth tip for easy insertion and comfort
- » Staggered eyes for maximum flow

Item No.	Description	Pkg
----------	-------------	-----

6"/15.2 cm Female Vinyl Intermittent Self Catheter

DYND10701	10 Fr	30/cs
DYND10702*	12 Fr	30/cs
DYND10703*	14 Fr	30/cs

10"/25.4 cm Pediatric Vinyl Intermittent Self Catheter

DYND10712*	6 Fr	30/cs
DYND10713*	8 Fr	30/cs
DYND10714*	10 Fr	30/cs

16"/40.6 cm Universal Vinyl Intermittent Self Catheter

DYND10721*	10 Fr	30/cs
DYND10722*	12 Fr	30/cs
DYND10723*	14 Fr	30/cs
DYND10700*	14 Fr	100/cs
DYND10724*	16 Fr	30/cs
DYND10725	18 Fr	30/cs

*Add "H" to item to order by the each.



6-inch/15.2 cm Female Catheter



10-inch/25.4 cm Pediatric Catheter



16-inch/40.6 cm Male/Universal Catheter



Hydrophilic Straight Catheters

- » Minimizes friction during catheterization to help prevent discomfort and microtrauma
- » Does not require any additional lubrication
- » 16-inch/40.64 cm Universal

Item No.	Description	Pkg
URO3PVC1016L	10 Fr	30/cs
URO3PVC1216L	12 Fr	30/cs
URO3PVC1416L	14 Fr	30/cs
URO3PVC1616L	16 Fr	30/cs
URO3PVC1816L	18 Fr	30/cs



Vinyl Coude Intermittent Catheters

- » Featuring a curved, coude, tip for help in navigating past strictures or an enlarged prostate
- » 16-inch/40.6 cm length
- » Not made with natural rubber latex

Item No.	Description	Pkg
URO3PVC1016C	10 Fr	30/cs
URO3PVC1216C	12 Fr	30/cs
URO3PVC1416C	14 Fr	30/cs
URO3PVC1616C	16 Fr	30/cs
URO3PVC1816C	18 Fr	30/cs



DYND12574

URO5454

DYND12588

Leg Bags with Drainage Valve

- » These sturdy bags feature a comfort strap, twist flip, or slide-tap valve drainage port, and a sterile fluid pathway
- » Not made with natural rubber latex

Item No.	Description	Pkg
----------	-------------	-----

Leg Bag with Twist Valve

DYND12574*	Medium, 20-oz./600 ml	48/cs
DYND12578*	Large, 30-oz./900 ml	48/cs

Leg Bag with Flip Drain Valve

URO5454*	Medium, 20-oz./600 ml	48/cs
----------	-----------------------	-------

Leg Bag with Slide-Tap Valve

DYND12584*	Medium, 20-oz./600 ml	48/cs
DYND12588*	Large, 30-oz./900 ml	48/cs

Leg Bag Strap with Hook-and-Loop Fasteners

DYND12582H		1/pr
------------	--	------

*Add "H" to item to order by the each.



URO4PVCP12MC

My-Cath Touch-Free Self Catheter System

- » Special closed-system design helps reduce mess, eases disposal and reduces contamination risk by allowing a touch-free catheterization
- » The catheter and introducer are coated with a water-soluble lubricant for ease of insertion
- » Catheter guide allows for simple, protected control of catheter during insertion
- » 1,500 ml collection bag eliminates the need for separate container, reducing the chance of spillage
- » Kit includes: collection bag, catheter, moisture-proof underbuttocks drape, wrapped gloves and PVP swabsticks (3/pk)

Item No.	Description	Pkg
URO4PVCP12MC	12 Fr	60/cs
URO4PVCP12MH	12 Fr	1/ea
DYND10440	14 Fr	60/cs
DYND10440H	14 Fr	1/ea
URO4PVCP16MC	16 Fr	60/cs
URO4PVCP16MH	16 Fr	1/ea



Coude Tip



Straight Tip



DYND13514

Red Rubber All-Purpose Intermittent Catheters

- » Each style has a smooth tip and funnel end. 16-inch/40.64 cm length is suitable for male or female catheterization

Item No.	Description	Pkg
----------	-------------	-----

Coude-Tip Red Rubber Intermittent Catheter

DYND13614	14 Fr	10/cs
-----------	-------	-------

Caution: This product contains natural rubber latex which may cause allergic reactions.

Red Poly-Cath Urethral Catheter

DYND13714	14 Fr*	10/cs
-----------	--------	-------

*Not made with natural rubber latex.

Item No.	Description	Pkg
----------	-------------	-----

Red Rubber Latex Urethral Catheter, Robinson Style

DYND13508	8 Fr	12/cs
DYND13510	10 Fr	12/cs
DYND13512	12 Fr	12/cs
DYND13514**	14 Fr	12/cs
DYND13515	15 Fr	50/cs
DYND13516*	16 Fr	12/cs
DYND13518	18 Fr	12/cs
DYND13520	20 Fr	12/cs
DYND13522	22 Fr	12/cs

Caution: This product contains natural rubber latex which may cause allergic reactions.

**Add "H" to item to order by the each.



DYNC1810

Foley Catheter Insertion Tray

- » Peel-lid tray has the materials required to insert a Foley catheter: a graduated collection basin, fenestrated drape, moisture-proof underbuttocks drape, gloves not made with natural rubber latex, sterile water inflation syringe, triple pack of PVP swabsticks, and lubricating jelly packet
- » Not made with natural rubber latex

Item No.	Description	Pkg
DYNC1810*	Tray with 10 ml Syringe	20/cs

*Add "H" to item to order by the each.

DYND10402



Pre-Connected Intermittent Catheter Trays

- » All the necessary components available in order to complete an intermittent catheterization, with the additional benefit of a pre-attached 30-oz./900 ml drain bag
- » Available with a latex or not made with natural rubber latex catheter
- » All kits include:
 - Aloetouch 3G gloves
 - 3 pack of PVP swabsticks
 - Underbuttocks drape
 - Lubricating jelly
 - Fenestrated drape
 - Specimen container with label

Item No.	Description	Pkg
DYND10402*	Medium, 20-oz./600 ml, 14 Fr Vinyl Catheter	20/cs
DYND10407*	Large, 30-oz./900 ml, 15 Fr Latex Catheter	20/cs

*Add "H" to item to order by the each.



Intermittent Catheterization Trays

- » Provide you with all of the necessary supplies required for catheterization
- » All insertion kits are sterile
- » All kits include:
 - Graduated basin
 - Swabsticks (3/pk)
 - Gloves not made with natural rubber latex
 - Lubricating jelly
 - Underbuttocks drape
 - Specimen container with label
 - Fenestrated drape



DYND10305

Item No.	Peel-Lid Tray	CSR Wrapped Tray	Red Rubber Catheter	Vinyl Catheter	PVP Swabsticks 3/pk	BZK Swabsticks 3/pk	Aloetouch 3G Gloves	Vinyl Gloves	Pkg
DYNC1820*	•			14 Fr	•			•	20/cs
DYNC1816*	•		15 Fr		•			•	20/cs
DYND10300*		•		14 Fr	•		•		20/cs
DYND10305*		•	No Cath	No Cath	•		•		20/cs
DYND10350*		•	15 Fr		•		•		20/cs
DYND10500*		•		14 Fr		•	•		20/cs

*Add "H" to item to order by the each.



DYND15205



Urinary Drain Bag with Anti-Reflux Tower

- » Facilitates urine flow, especially useful for urine with heavy sediment or mucous
- » Designed for easy bag suspension; 50-inch/127 cm premium 1¹/₃₂-inch/0.9 cm STAR™ tubing provides maximum drainage and minimum kinking. Includes bed sheet clamp and slide-tap drainage mechanism
- » Sterile fluid pathway and not made with natural rubber latex

Item No.	Description	Pkg
DYND15205*	Urinary Drain Bag, 2,000 ml	20/cs

*Add "H" to item to order by the each.



DYND12200



DYND16900

Accessories

- » Catheter plug fits securely into funnel end of most Foley catheters
- » Cover protects connector on drainage bag
- » Leg bag extension tubing can be cut to length to meet the patient's needs
- » Elastic catheter strap helps keep indwelling catheters in place

Item No.	Description	Pkg
DYND12550*	Leg Bag Extension Tubing, 18"/45.7 cm, Sterile, Not made with Natural Rubber Latex	50/cs
DYND12200*	Catheter Plug and Drain Tube Protector	100/cs
DYND16900*	Not made with Natural Rubber Latex Catheter Strap, 22 x 2"/56 x 5 cm	10/bx

*Add "H" to item to order by the each.



Speci-Cath Kits

- » Kits provide an easy way to collect a sterile urine sample
- » Minimally invasive catheter is pre-inserted into the specimen tube to ensure sterile technique
- » Sterile
- » Not made with natural rubber latex

Item No.	Catheter Size	Catheter Length	15 ml Specimen Tube with Cap	Specimen Tube Label	Aloetouch 3G Gloves	PVP Swabsticks 3/pk	Lubricating Jelly	Pkg
DYND10805*	8 Fr, Female	6½"/17 cm	•	•	•	•	•	1/ea
DYND10806	8 Fr Mini, Female	6½"/17 cm	•	•	•	•		25/cs
DYND10800	8 Fr, Female	6½"/17 cm	•	•				1/ea
DYND10810*	8 Fr, Male	14"/36 cm	•	•	•	•	•	1/ea
DYND10815*	8 Fr, Pediatric	9"/23 cm	•	•	•	•	•	1/ea
DYND10820*	5 Fr, Neonatal	9"/23 cm	•	•	•	•	•	1/ea

*Add "H" to item to order by the each.



DYND80024H



24-Hour Urine Collection Bottle

- » Easy-grip, extra-wide handle makes this amber-colored, plastic bottle a convenient addition to everyday medical supplies
- » 101.4-oz./3,000 ml graduated bottle
- » Can be used with calculi strainer

Item No.	Description	Pkg
DYND80024H	Collection Bottle	1/ea



DYND4712H



Calculi Strainer

- » Use this strainer with the 24-hour urine collection bottle to provide stability during the collection process
- » Fine mesh filters stones and foreign matter
- » Built-in ridge serves as a handle for simple application and removal

Item No.	Description	Pkg
DYND4712H	Calculi Strainer	1/ea



MDS190510



Infant Urine Collector

- » Versatile infant urine collector can be used with both male and female infants
- » Adhesive is gentle on skin and easy to remove, yet holds securely
- » Folds to form a spillproof specimen transport bag
- » Available in both sterile and non-sterile packaging
- » 5-oz./150 ml

Item No.	Description	Pkg
MDS190505	Pediatric Urine Collector, Non-Sterile	50/bx
MDS190510*	Pediatric Urine Collector, Sterile	50/bx

*Add "H" to item to order by the each.



DYND30330



Specimen Containers, Polypropylene

- » General-use specimen containers with screw-on lids. Easy-to-read graduations in both ounces and milliliters
- » Packaged Individually in Poly-Bag, Sterile Fluid Pathway, Patient ID Label

Item No.	Description	Pkg
DYND30330*	4-oz./120 ml	100/cs
DYND30350*	4.5-oz./130 ml	100/cs

*Add "H" to item to order by the each



DYND36605H



Speci-Pan

- » Designed for efficient collection of urine and stool specimens, these pans feature a measuring device that is graduated to 30.4-oz./900 ml or 4.6-oz./1,200 ml
- » Reinforced ribbing for durability
- » Fits any standard size commode

Item No.	Description	Pkg
DYND36600H	Translucent, 30.4-oz./900 ml	1/ea
DYND36605H	Deluxe, White, 40.6-oz./1,200 ml	1/ea



DYND80235S



Translucent Urinals

- » Helps reduce spills and odor
- » Attached lid helps prevent spilling and reduces odors with odor shield. Contoured handle attaches easily to healthcare beds. Graduated to 32-oz./1,000 ml

Item No.	Description	Pkg
DYND80235S	Male Urinal	48/cs
DYND80235H	Male Urinal	1/ea
DYND80234	Male Urinal, Mini Pack	12/cs



DYND36500



Two-Piece Specimen Collector for Commode

- » Made for security
- » Light blue stabilizing frame and a white bowl
- » Graduated up to 30.4-oz./900 ml

Item No.	Description	Pkg
DYND36500H	Specimen Collector	1/ea



NON21001Z



Sterile Disposable Drapes

- » Paper drapes have a repellent inner layer that helps prevent moisture strike-through

Item No.	Description	Pkg
NON21001Z	Sterile Field, Tissue-Poly-Tissue, 18 x 26"/45.72 x 66 cm	50/bx



MSC281226C



Protection Plus® Disposable Underpads

- » The non-woven facing of these pads is soft against skin
- » Absorbent material aids fluid containment
- » Wide variety of sizes and weights to meet every budget and need
- » Polypropylene backing helps keep underpad in place

Item No.	Description	Pkg
----------	-------------	-----

Standard Weight, Moderate Absorbency

MSC281230	17 x 24"/43 x 61 cm	300/cs
MSC281236P	23 x 36"/58 x 91 cm	10/bg, 15 bg/cs, 150/cs

Economy Weight, Light Absorbency

MSC281224C	17 x 24"/43 x 61 cm	300/cs
MSC281226C	23 x 36"/58 x 91 cm	150/cs
MSCB281245LB	30 x 30"/76 x 76 cm	5/bg, 150/cs



NON21002Z



Sterile Fenestrated Disposable Drapes

- » Combat moisture strike-through with a repellent inner layer

Item No.	Description	Pkg
NON21002Z	Sterile Field, 18 x 26"/45.72 x 66 cm, 2.75"/6.98 cm	50/bx



NON241288

TAMPAX

Tampons

- » Individually wrapped tampons have a cardboard sleeve and flushable applicator
- » Vending machine compatible

Item No.	Description	Pkg
NON241288	For Bulk Vending	200/cs
NON241288H	For Bulk Vending	1/ea



NON241289



Maxi Sanitary Pads

- » Offer maximum absorbency for increased confidence and longer protection
- » Packaged individually or in bulk for easy distribution

Item No.	Description	Pkg
NON241289Z	11"/28 cm, with Wings, Individually Wrapped	12/bg
NON241278	9"/20 cm, Non-Sterile	1/bx, 250/cs
NON241278H	9"/23 cm, Sterile with Adhesive Strip	1/ea
NON241279	9"/20 cm, Non-Sterile	12/bg, 600/cs
NON241279H	9"/23 cm, Sterile	12/bg
NON241289	11"/28 cm, Non-Sterile with Wings, Individually Wrapped	12/bg, 288/cs
NON241286	11"/28 cm, Non-Sterile, Individually Wrapped	12/bg, 288/cs
NON241286Z	11"/28 cm, Non-Sterile, Individually Wrapped	12/bg
MPL1228	7", Panty Liner	28/bx, 336/cs



MSC323950C



Maternity Pads

- » Provides reliable protection for new moms
- » Designed specifically for postpartum flow

Item No.	Description	Pkg
NON241280	Non-Sterile with Tails, Without Adhesive, Bagged: 11"/28 cm	12/bg, 288/cs
NON241281	Sterile with Tails, Without Adhesive, Individually Wrapped: 11"/28 cm	12/bg, 288/cs
NON241283	Sterilizable, with Tails, Without Adhesive, Individually Wrapped: 11"/28 cm	12/bg, 288/cs
MDT214055	Sanitary Belt, for Use with Maternity Pads	6 dz/bx
MSC323950C	Contoured, Peach, Non-Sterile with Adhesive Back: 7 x 14"	12/bg, 240/cs
MSC323000C	Non-Sterile with Adhesive Back: 7 x 17"	12/bg, 240/cs
MSC322505	Non-Sterile with Adhesive Back: 9 x 24"	20/bg, 80/cs



MSC263040H



Fragrance-Free Baby Wipes

- » pH-balanced, hypoallergenic and alcohol-free wet wipes provide exceptional skin care
- » Convenient, high-count packs with standard-weight wipes
- » Plastic, with flip-top lid
- » Soft pack, fragrance-free
- » Do not flush wipes in the toilet

Item No.	Description	Pkg
MSC263040H	Soft Pack, 5.9 x 7.9"/15 x 20 cm	80 ea/pk



SUS750B



Surgical Suture

Item No.	Description	Pkg
Poly Mono		
SUS1020B	4/0, Black, PC31	12/bx
SUS1085B	5/0, Blue, C3	12/bx
Nylon, Black		
SUS1282B	6/0, 18"/45 cm, C3	12/bx
SUS1291B	4/0, 18"/45 cm, C3	12/bx
SUS931B	4/0, 10"/25.4 cm, C17	12/bx
Polyglycolic, Undyed		
SUS493B	5/0, C3	12/bx
SUS494B	4/0, C3	12/bx
Silk Braided, Black		
SUS750B	6/0, 18"/45 cm, C22	12/bx
Nylon Mono Cuticular		
SUS928B	3/0, 18"/45 cm, C6	12/bx
Silk Braided, Black		
SUSSP104BX	4/0, 25-yds./24.7 m, No Needle	12/bx



Covidien-Medtronic® Sutures

Item No.	Description	Pkg
Polysorb™, Undyed		
MSUSL5628GMDL	3/0, 18"/45 cm, P-12	12/bx
MSUSL5627GMDL	4/0, 18"/45 cm, P-12	12/bx
MSUSL5691GMDL	4/0, 18"/45 cm, P-13	12/bx
MSUSL5687GMDL	5/0, 18"/45 cm, P-13	12/bx
Biosyn™, Undyed		
MSUSM5637GMDL	4/0, 30"/76 cm, P-12	12/bx
MSUSM5687GMDL	5/0, 18"/45 cm, P-13	12/bx
Surgipro™, Blue		
MSUSP5698GMDL	5/0, 18"/45 cm, P-13	12/bx
MSUSP5695GMDL	6/0, 18"/45 cm, P-13	12/bx
Monosof™, Black		
MSUSN5667GMDL	4/0, 18"/45 cm, P-12	12/bx
MSUSN5699GMDL	4/0, 18"/45 cm, P-13	12/bx
MSUSN5698GMDL	5/0, 18"/45 cm, P-13	12/bx
MSUSN5696GMDL	6/0, 18"/45 cm, P-13	12/bx



ETH1698G

ETHICON Shaping the future of surgery

J&J Ethicon® Suture

Item No.	Description	Pkg
Ethilon® Suture, Black Monofilament		
ETH1698G	6/0, 18"/45 cm, P-3, 3/8 Circle Needle	12/bx
ETH698G	5/0, 18"/45 cm, P-3, 3/8 Circle Needle	12/bx
ETH699G	4/0, 18"/45 cm, P-3, 3/8 Circle Needle	12/bx

Prolene® Suture, Blue Monofilament

ETH8698G	5/0, 18"/45 cm, P-3, 3/8 Circle Needle	12/bx
ETH8699G	4/0, 18"/45 cm, P-3, 3/8 Circle Needle	12/bx

Vicryl® Suture, Undyed Braided

ETHJ493G	5/0, 18"/45 cm, P-3, 3/8 Circle Needle	12/bx
ETHJ494G	4/0, 18"/45 cm, P-3, 3/8 Circle Needle	12/bx

Dermabond® Topical Skin Adhesive

ETHDHVM12	2-Octyl Cyanoacrylate, 0.25 ml	12/bx
ETHDNX6	Dermabond Advanced, 0.7 ml	6/bx



ETHICON Shaping the future of surgery

Reverse Cutting Suture

» For unprecedented monofilament pliability and smooth passage through tissue. Indicated for use in general soft tissue approximation and/or ligation, but not for use in cardiovascular or neurological tissues

Item No.	Description	Pkg
ETHY494G	4-0, 18"/45 cm, P-3, 3/8 Circle	12/bx

Precision Point, Undyed Monofilament



ETHICON Shaping the future of surgery

Absorbable Suture

» Composed of purified connective tissue (mostly collagen) derived from either serosal layer of beef (bovine) or submucosal fibrous layer of sheep (ovine) intestines

» Sterile

Item No.	Description	Pkg
Chromic, Taper Point		
ETHG122H	3/0, 27"/70 cm, Needle: SH, 1/2 Circle, 26 mm Length, Chromic Gut	36/bx



PAINLOR5LL25



PAINLOR10LL25



Glass Loss of Resistance Syringes

Item No.	Description	Pkg
PAINLOR5LL25	5 ml, Metal Luer Lock Tip	25 ea/cs
PAINLOR5LS25	5 ml, Metal Luer Slip Tip	25 ea/cs
PAINLOR10LL25	10 ml, Metal Luer Lock Tip	25 ea/cs
PAINLOR10LS25	10 ml, Metal Luer Slip Tip	25 ea/cs



Procedural Drapes and Syringes

» Complement Regional Anesthesia and Pain Management trays

Item No.	Description	Pkg
Procedural Drapes		
PAINDRAPE1	Clear Plastic Drape, 24 x 26"/60.9 x 66.0 cm, 4"/10.2 cm, Circle Fenestration, Adhesive Along One Edge	10 ea/cs
PAINDRAPE2	Blue Paper Drape, 18 x 27"/45.7 x 68.6 cm, 3 x 4"/7.6 x 10.2 cm, Rectangle Fenestration, Adhesive Along One Edge	25 ea/cs
PAINDRAPE3	Clear Plastic Drape, 24 x 26"/60.9 x 66.0 cm, 4 x 8"/10.2 x 23.3 cm, Rectangle Fenestration, Adhesive Along One Edge	10 ea/cs
PAINDRAPE4	Clear Plastic Drape, 24 x 26"/60.9 x 66.0 cm, 4"/10.2 cm, Circle Fenestration, Adhesive Patch Around Fenestration and Along One Edge	10 ea/cs

Exactly what you want. Everything you need.

What do you look for in Regional Anesthesia and Pain Management trays? We offer custom trays for maximum efficiency and stock trays for simplicity. Whatever your preference, we have what you need—when you need it. Everything from syringes and needles to surgical drapes and gloves in one, sterile kit.

Create Efficiencies

Eliminate time wasted picking additional products.

Increase Value

Pay for only the products you know you'll use.

Reduce Waste

Avoid throwing away unused components.



Block Trays

Description

Procedural Needles

	DYNJRA9026*	DYNJRA9036	DYNJRA9037	DYNJRA9029*	DYNJRA9031	PAIN9016*	DYNJRA9028*	DYNJRA9030	DYNJRA9032	PAIN0609	PAIN9013S	PAIN9014S
22 g x 3.5"/8.9 cm, Quincke Spinal Needle	•	•										
25 g x 3.5"/8.9 cm, Quincke Spinal Needle			•									
24 g x 4"/10.2 cm, Whitacre Spinal Needle				•	•							
25 g x 3.5"/8.9 cm, Whitacre Spinal Needle						•	•	•	•			
20 g Introducer Needle			•	•	•	•	•	•	•			

Procedural Syringes

5 ml/0.17-fl. oz., Glass Spinal Syringe Luer Slip	•	•	•	•	•	•	•	•	•			
---	---	---	---	---	---	---	---	---	---	--	--	--

Pharmaceuticals

Bupivacaine 0.75% in Dextrose 8.25%, 2 ml/0.07-fl. oz.	•			•	•		•	•	•			
Lidocaine 1% 5 ml/0.17-fl. oz., Glass Ampule	•			•	•		•	•	•			•
Sodium Chloride 0.9% 10 ml/0.34-fl. oz., Glass Ampule												•

Hypodermic Needles

18 g x 1.5"/3.8 cm, Hypodermic Needle	•	•	•	•	•	•	•	•		•	•	•
22 g x 1.5"/3.8 cm, Hypodermic Needle	•	•	•			•	•		•	•	•	•
25 g x 1.0"/2.5 cm, Hypodermic Needle									•			
25 g x 1.5"/3.8 cm, Hypodermic Needle	•	•	•	•	•	•	•	•		•	•	•
19 g x 1.5"/3.8 cm, Filter Needle	•	•	•	•	•	•	•	•		•	•	•
18 g x 1.5"/3.8 cm, Blunt Fill Needle									•			
Foam Needle Pad	•	•	•	•	•	•	•	•	•	•	•	•

Syringes

3 ml/0.1-fl. oz., Syringe Luer Lock	•	•	•	•	•	•	•	•	•		•	•
10 ml/0.34-fl. oz., Syringe Luer Lock	•	•	•			•	•			•	•	•
20 ml/0.7-fl. oz., Syringe Luer Lock	•	•	•			•	•			•	•	•

Miscellaneous Components

Three Sponge Sticks, Angled	•	•	•	•	•	•	•	•	•	•	•	•
Povidone-Iodine Solution	•	•	•	•	•	•	•	•	•	•	•	•
Gauze, 4 x 4"	•	•	•	•	•	•	•	•	•	•	•	•
Blue Paper Drape 18 x 27" with 3 x 4" Fenestration and Adhesive Strip	•	•	•			•	•		•	•	•	•
Clear Plastic Drape 24 x 26" with 4" Circular Fenestration and Adhesive Strip				•	•			•				
Absorbent Towel	•	•	•	•	•	•	•	•	•	•	•	•
Drug Label Sheet Six Pre-Printed and Four Blank Labels	•	•	•	•	•	•	•	•		•	•	•
Adhesive Bandage 1 x 3"									•			
6" Extension Set, Microbore										•		

All kits are packaged in two-deck, pre-formed trays wrapped in CSR. 20 sterile kits per case except for those marked with*, which designates that 10 sterile kits come per case. All kits are not made with natural rubber latex. Sterilized by ethylene oxide. Single use only.



Single Shot Epidural Trays

Description

Item No.	PAIN9019S	PAIN9020S	PAIN9001S	PAIN9002S	PAIN9007S	PAIN9009S	PAIN9010S	PAIN0611	PAIN9003S	PAIN9004S	PAIN9005S	PAIN9006S	PAIN9011S	PAIN9012S
Procedural Needles														
17 g x 3.5"/8.9 cm, Tuohy Epidural Needle with Wings	•	•												
18 g x 3.5"/8.9 cm, Tuohy Epidural Needle with Wings			•	•	•									
20 g x 3.5"/8.9 cm, Tuohy Epidural Needle with Wings						•	•	•	•	•				
Procedural Syringes														
5 ml/0.17-fl. oz., Glass Loss-of-Resistance Syringe Metal Luer Lock	•	•	•	•					•	•	•	•		
7 ml/0.24-fl. oz., Plastic Loss-of-Resistance Syringe Luer Lock					•	•	•	•					•	•
Pharmaceuticals														
Lidocaine 1% 5 ml/0.17-fl. oz., Glass Ampule		•		•			•	•		•		•		•
Sodium Chloride 0.9% 10 ml/0.34-fl. oz., Glass Ampule		•		•			•	•		•		•		•
Hypodermic Needles														
18 g x 1.5"/3.8 cm, Hypodermic Needle	•	•	•	•	•	•	•	•	•	•	•	•	•	•
22 g x 1.5"/3.8 cm, Hypodermic Needle	•	•	•	•	•	•	•	•	•	•	•	•	•	•
25 g x 1.5"/3.8 cm, Hypodermic Needle	•	•	•	•	•	•	•	•	•	•	•	•	•	•
19 g x 1.5"/3.8 cm, Filter Needle	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Foam Needle Pad	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Syringes														
3 ml/0.1-fl. oz., Syringe Luer Lock	•	•	•	•	•	•	•	•	•	•	•	•	•	•
10 ml/0.34-fl. oz., Syringe Luer Lock	•	•	•	•	•	•	•	•	•	•	•	•	•	•
20 ml/0.7-fl. oz., Syringe Luer Lock	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Miscellaneous Components														
Three Sponge Sticks, Angled	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Povidone-Iodine Solution	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Gauze, 4 x 4"	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Blue Paper Drape 18 x 27"	•	•	•	•	•	•	•	•	•	•	•	•	•	•
with 3 x 4" Fenestration and Adhesive Strip	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Absorbent Towel	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Drug Label Sheet, Six Pre-Printed and Four Blank Labels	•	•	•	•	•	•	•	•	•	•	•	•	•	•
6" Extension Set, Microbore								•						

All kits are packaged in two-deck, pre-formed trays wrapped in CSR. 20 sterile kits per case. All kits are not made with natural rubber latex. Sterilized by ethylene oxide. Single use only.



PAIN8005



PAIN8017



PAIN8019

Procedural Needles

Item No.	Description	Pkg
Epidural Needles		
PAIN8003	8 g x 3.5"/8.9 cm Hustead	25 ea/cs
PAIN8004	8 g x 2"/5.08 cm Tuohy	25 ea/cs
PAIN8001	17 g x 3.5"/8.9 cm	25 ea/cs
PAIN8002	17 g x 5"/12.7 cm	25 ea/cs
PAIN8005	18 g x 3.5"/8.9 cm	25 ea/cs
PAIN8006	18 g x 6"/15.2 cm	25 ea/cs
PAIN8007	18 g x 7"/17.8 cm Tuohy	25 ea/cs
PAIN8009	20 g x 3.5"/8.9 cm	25 ea/cs
PAIN8008	20 g x 4.5"/11.43 cm Tuohy	25 ea/cs
PAIN8010	20 g x 6"/15.3 cm	25 ea/cs
PAIN8013	22 g x 3.5"/8.9 cm	25 ea/cs

Item No.	Description	Pkg
Quincke Spinal Needles		
PAIN8021	18 g x 3.5"/8.9 cm	25 ea/cs
PAIN8033	18 g x 6"/15.3 cm	25 ea/cs
PAIN8022	20 g x 3.5"/8.9 cm	25 ea/cs
PAIN8028	20 g x 6"/15.3 cm	25 ea/cs
PAIN8011	22 g x 1.5"/3.8 cm	25 ea/cs
PAIN8012	22 g x 2.5"/6.4 cm	25 ea/cs
PAIN8017	22 g x 3.5"/8.9 cm	25 ea/cs
PAIN8023	22 g x 5"/12.7 cm	25 ea/cs
PAIN8024	22 g x 7"/17.8 cm	25 ea/cs
PAIN8016	25 g x 2.5"/6.4 cm	25 ea/cs
PAIN8034	23 g x 3.5"/8.9 cm	25 ea/cs
PAIN8018	25 g x 3.5"/8.9 cm	25 ea/cs
PAIN8025	25 g x 5"/12.7 cm	25 ea/cs
PAIN8036	23 g x 6"/15.3 cm	25 ea/cs

Item No.	Description	Pkg
Whitacre Spinal Needles		
PAIN8015	22 g x 3.5"/8.9 cm	25 ea/cs
PAIN8014	22 g x 5"/12.7 cm	25 ea/cs
PAIN8030	24 g x 4"/10.0 cm	25 ea/cs
PAIN8019	25 g x 3.5"/8.9 cm	25 ea/cs
PAIN8029	25 g x 5"/12.7 cm	25 ea/cs
PAIN8031	27 g x 3.5"/8.9 cm	25 ea/cs
PAIN8032	27 g x 5"/12.7 cm	25 ea/cs



ORA36M1C



Portable Cryosurgical System

- » For the treatment of warts and benign skin lesions
- » A one-hand procedure takes less than a minute*
- » Minimum treatments per canister varies with product size
- » Average shelf life is two to three years
- » Resolution rate of 93.2% vs. liquid nitrogen at 95%**

Item No.	Description	Pkg
ORA10010292Z	60 Applicators	1/kt
ORA36M1C	36 Applicators	1/kt

*Data on file, available at OraSure Technologies, Inc., Bethlehem, PA.
 **F. Caballero Martínez, MD, et al. Dermatological cryosurgery in primary care with dimethyl ether propane spray in comparison with liquid nitrogen. Translated from: Atención Primaria. 1996; Vol. 18 No. 5 (211, 216).



PMU9034000



CryOmega® Cryosurgical Device

- » Disposable cryosurgical device utilizing liquid nitrous oxide. Internal cartridge provides 240 seconds of spray time for use on multiple patients
- » Approved to treat nine different benign lesions
- » Entire unit is discarded when empty

Item No.	Description	Pkg
PMU9034000	CryOmega Unit	1/ea
PMU9034002	CryOmega Twin Pack	1/ea

DYND80301H



Rectangular Wash Basins

- » The large soaking area equally accommodates limb soaking
- » Rolled rims make carrying easy
- » Slots in rim allow easy clip-on for a soap dish accessory (not included)

Item No.	Description	Pkg
----------	-------------	-----

10 x 13¼ x 3¾"/25 x 34 x 10 cm

DYND80301H	Gold, 8 qt./7.6 L	1/ea
DYND80342H	Graphite, 8 qt./7.6 L	1/ea

9¼ x 11¼ x 4½"/23 x 30 x 11 cm

DYND80347H	Gold, 6-qt./5.7 L	1/ea
------------	-------------------	------

DYND70292



Denture Containers

- » Ideal for safe and convenient bedside storage
- » Clear plastic imprinted lids make identifying contents easy
- » 8-oz./227 g capacity

Item No.	Description	Pkg
----------	-------------	-----

DYND70292Z	Translucent	25/tb
DYND70292H	Translucent	1/ea
DYND70293Z	Turquoise	25/tb
DYND70293H	Turquoise	1/ea



DYND80367H

Emesis Basins

- » Gradations are printed in ounces and cc's
- » Ideal for oral hygiene and ear irrigation

Item No.	Description	Pkg
----------	-------------	-----

16-oz./500 ml, 8½"/22 cm

DYND80321	Gold	1/ea
DYND80322	White	250/cs
DYND80327H	Graphite	1/ea

20-oz./700 ml, 10"/25 cm

DYND80367H	Graphite	1/ea
------------	----------	------



DYND70280

Ear Ulcer Syringes

- » Ideal for irrigating smaller wounds and cleansing nasal passages of newborns
- » Made of PVC
- » Individually packaged

Item No.	Description	Pkg
----------	-------------	-----

DYND70279	1-oz./29.6 ml, Sterile	50/cs
DYND70280	2-oz./59.1 ml, Sterile	50/cs
URO70283	2-oz./59.1 ml, Sterile, Red	24/cs
DYND70277	3-oz./88.7 ml, Sterile	50/cs

*Add H to end of item number to order by the each.



ORT32700

Stackies Kit and Replacements

- » Quick therapy to protect and stabilize the distal fingertip
- » Open nail design with ventilation holes
- » 30 splints in eight sizes
- » Made with polypropylene

Item No.	Description	Pkg
----------	-------------	-----

ORT32700	Assorted Sizes	1/ea
----------	----------------	------



ORT32200S

Frog Finger Splint

- » Folds around finger to stabilize
- » Malleable aluminum lined with blue polyester urethane foam

Item No.	Description	Pkg
----------	-------------	-----

ORT32200S	Small, 2.25 x 2.75"/5.7 x 7 cm	12/cs
ORT32200M	Medium, 2.5 x 3"/6.4 x 7.6 cm	12/cs
ORT32200L	Large, 3 x 3.5"/7.6 x 8.9 cm	12/cs



ORT33100RM

Forearm Splint

- » Pre-molded and malleable for anatomical fit
- » Padded aluminum with blue polyester urethane foam

Item No.	Description	Pkg
----------	-------------	-----

Right, Padded		
ORT33100RM	Medium	1/ea

Left, Padded		
ORT33100LM	Medium	1/ea



ORT33200RM

Metacarpal Splint

- » Pre-molded and malleable for anatomical fit
- » Padded aluminum with blue polyester urethane foam

Item No.	Description	Pkg
----------	-------------	-----

Right, Padded		
ORT33200RS	Small	1/ea
ORT33200RM	Medium	1/ea

Left, Padded		
ORT33200LS	Small	1/ea
ORT33200LM	Medium	1/ea



MDS067003Z



Non-Sterile Synthetic Undercast Padding

- » Synthetic cast padding rolls are ideal for use with all forms of synthetic hard casting material
- » 100% polyester needle-punched construction absorbs excess moisture and dries quickly to retain a comfortable fit

Item No.	Description	Pkg
MDS067003Z	3" x 4-yds./ 7.26 cm x 3.66 m	12/bg
MDS067004Z	4" x 4-yds./ 10.16 cm x 3.66 m	12/bg
MDS067006Z	6" x 4-yds./ 15.24 cm x 3.66 m	6/bg



MDS066004



Non-Sterile Cotton Undercast Padding

- » Constructed of 100% pure cotton fibers to provide a soft environment under a hard cast
- » A low linting product with more loft to provide a better feel for the patient

Item No.	Description	Pkg
MDS066002ZZ	2" x 4-yds./ 5.08 cm x 3.66 m	24/bg
MDS066003Z	3" x 4-yds./ 7.62 cm x 3.66 m	12/bg
MDS066004Z	4" x 4-yds./ 10.16 cm x 3.66 m	12/bg
MDS066006Z	6" x 4-yds./ 15.24 cm x 3.66 m	6/bg



DYNJ066002



Sterile Wytex Undercast Padding

- » Constructed of 100% pure cotton fibers
- » Wytex is low-linting with more loft to provide a better feel for the patient

Item No.	Description	Pkg
DYNJ066002H	2" x 4 yds./ 5.08 cm x 3.66 m	1/ea
DYNJ066003H	3" x 4-yds./ 7.26 cm x 3.66 m	1/ea
DYNJ066004H	4" x 4-yds./ 10.16 cm x 3.66 m	1/ea
DYNJ066006H	6" x 4-yds./ 15.24 cm x 3.66 m	1/ea



MDT221203



Non-Sterile Tubular Stockinettes

- » Heavyweight, premium-grade, 100% unbleached cotton tubular stockinettes are made for use under and over casts, as well as for other dressing applications
- » Not made with natural rubber latex

Item No.	Description	Pkg
MDT221203	3" x 25-yds./8 cm x 22.9 m	1/r
MDT221204	4" x 25-yds./10 cm x 22.9 m	1/r
MDT221206	6" x 25-yds./15 cm x 22.9 m	1/r



CT24BLK



Fiberglass Casting Tape

- » Water-activated resin sets in 20 minutes for patient and staff convenience
- » Dries to a smooth finish; gentle against exposed skin
- » Fresh fragrance and color choice enhance patient experience
- » Durable and highly conformable
- » Allows for penetration by X-rays to monitor bone healing

Item No.	Description	Pkg
CT24BLK	2" x 4-yd., Black	10 rls/bx
CT34BLK	3" x 4-yd., Black	10 rls/bx
CT44BLK	4" x 4-yd., Black	10 rls/bx
CT24BLU	2" x 4-yd., Blue	10 rls/bx
CT34BLU	3" x 4-yd., Blue	10 rls/bx
CT44BLU	4" x 4-yd., Blue	10 rls/bx
CT24WHT	2" x 4-yd., White	10 rls/bx
CT34WHT	3" x 4-yd., White	10 rls/bx
CT44WHT	4" x 4-yd., White	10 rls/bx
CT24PNK	2" x 4-yd., Pink	10 rls/bx
CT34PNK	3" x 4-yd., Pink	10 rls/bx
CT44PNK	4" x 4-yd., Pink	10 rls/bx
CT24PRP	2" x 4-yd., Purple	10 rls/bx
CT34PRP	3" x 4-yd., Purple	10 rls/bx
CT44PRP	4" x 4-yd., Purple	10 rls/bx



Knee High Compression Hosiery

- » Fashionable, hospital-quality compression
- » 20-30 mmHg applies extra controlled pressure to the legs to help to prevent pooling of blood in heavy, aching legs; most commonly chosen treatment for varicose veins



Item No.	Size	Ankle Circumference	Calf Circumference	Length	Pkg
Beige					
MDS1702ATH	A	7-7.875"/18-20 cm	10-13"/25-33 cm	>15"/39 cm	1/pr
MDS1702BTH	B	8-8.875"/20-23 cm	12-15"/31-39 cm	>16"/41 cm	1/pr
MDS1702CTH	C	9-9.875"/23-25 cm	14-17"/36-44 cm	>17"/44 cm	1/pr
MDS1702DTH	D	10-10.875"/25-28 cm	16-19"/41-49 cm	>18"/46 cm	1/pr
MDS1702ETH	E	11-11.875"/28-31 cm	18-21"/46-54 cm	>18"/46 cm	1/pr
MDS1702FTH	F	12-12.875"/31-33 cm	20-23"/51-59 cm	>18"/46 cm	1/pr
MDS1702GTH	G	13-13.875"/33-36 cm	22-26"/56-66 cm	>18"/46 cm	1/pr

Black

MDS1703ABH	A	7-7.875"/18-20 cm	10-13"/25-33 cm	>15"/39 cm	1/pr
MDS1703BBH	B	8-8.875"/20-23 cm	12-15"/31-39 cm	>16"/41 cm	1/pr
MDS1703CBH	C	9-9.875"/23-25 cm	14-17"/36-44 cm	>17"/44 cm	1/pr
MDS1703DBH	D	10-10.875"/25-28 cm	16-19"/41-49 cm	>18"/46 cm	1/pr
MDS1703EBH	E	11-11.875"/28-31 cm	18-21"/46-54 cm	>18"/46 cm	1/pr
MDS1703FBH	F	12-12.875"/31-33 cm	20-23"/51-59 cm	>18"/46 cm	1/pr
MDS1703GBH	G	13-13.875"/33-36 cm	22-26"/56-66 cm	>18"/46 cm	1/pr

8-15 mmHg

Description	Size		
	Small	Medium	Large
Regular Beige	MDS1712ATH	MDS1712BTH	MDS1712CTH
Regular Black	MDS1713ABH	MDS1713BBH	MDS1713CBH
Short Beige	MDS1712ATSH	MDS1712BTSH	MDS1712CTSH
Short Black	MDS1713ABSH	MDS1713BBSH	MDS1713CBSH

20-30 mmHg

Description	Size						
	A	B	C	D	E	F	G
Regular Beige	MDS1702ATH	MDS1702BTH	MDS1702CTH	MDS1702DTH	MDS1702ETH	MDS1702FTH	MDS1702GTH
Regular Black	MDS1703ABH	MDS1703BBH	MDS1703CBH	MDS1703DBH	MDS1703EBH	MDS1703FBH	MDS1703GBH
Short Beige	MDS1702ATSH	MDS1702BTSH	MDS1702CTSH	MDS1702DTSH	MDS1702ETSH	MDS1702FTSH	MDS1702GTSH
Short Black	MDS1703ABSH	MDS1703BBSH	MDS1703CBSH	MDS1703DBSH	MDS1703EBSH	MDS1703FBSH	MDS1703GBSH

Sizing Chart for 8-15 mmHg only

Shoe Size	Height	Size	
6-8.5 Women's	<5'5"/165 cm	Short Regular	S
6-8 Men's	>5'5"/165 cm	Regular	
9-10.5 Women's	<5'6"/167 cm	Short Regular	M
8.5-10 Men's	>5'6"/167 cm	Regular	
11-13 Women's	<5'7"/170 cm	Short Regular	L
10.5-12 Men's	>5'7"/170 cm	Regular	

Sizing Chart for 15-20 mmHg, 20-30 mmHg and 30-40 mmHg

Ankle	Calf	Length	Size	
7"-7.88"/18-20 cm	10"-13"/25-33 cm	<15"/39 cm	Short	A
		>15"/39 cm	Regular	
8"-8.88"/20-23 cm	12"-15"/31-39 cm	<16"/41 cm	Short	B
		>16"/41 cm	Regular	
9"-9.88"/23-25 cm	14"-17"/36-44 cm	<17"/44 cm	Short	C
		>17"/44 cm	Regular	
10"-10.88"/25-28 cm	16"-19"/41-49 cm	<18"/46 cm	Short	D
		>18"/46 cm	Regular	
11"-11.88"/28-31 cm	18"-21"/46-54 cm	<18"/46 cm	Short	E
		>18"/46 cm	Regular	
12"-12.88"/31-33 cm	20"-23"/51-59 cm	<18"/46 cm	Short	F
		>18"/46 cm	Regular	
13"-13.88"/33-36 cm	22"-26"/56-66 cm	<18"/46 cm	Short	G
		>18"/46 cm	Regular	

Ankle Circumference

Measure around the ankle bone.

Calf Circumference

Measure around the widest point of the calf.

Length

Measure from the bend in the knee to the bottom of the heel.

MDS160644



MDS160624



Anti-Embolism Stockings

- » Color-coded toes for easy sizing
- » A tailored, contoured fit with just the right amount of gradient pressure
- » Pressure is greatest at the ankle and gradually decreases, 15-18 mmHg, throughout the length of the leg

Item No.	Size	Calf Circumference	Leg Length*	Pkg
Knee-Length Stockings				
MDS160624	Small, Regular	10 - 13"/25.5-33 cm	14-16"/36-40.5 cm	12 pr/bx
MDS160644	Medium, Regular	13 - 15"/33-38.1 cm	15-17"/38-43 cm	12 pr/bx
MDS160664	Large, Regular	15 - 17"/38.1-43 cm	16-18"/41-45.5 cm	12 pr/bx
MDS160684	X-Large, Regular	17 - 19"/43-48 cm	16-18"/41-45.5 cm	12 pr/bx
MDS160694	XX-Large, Regular	18 - 21"/46-53 cm	16-18"/41-45.5 cm	12 pr/bx
MDS160604	XXX-Large, Regular	21 - 26"/53-66 cm	16-18"/41-45.5 cm	12 pr/bx

Thigh-Length Stockings

MDS160824	Small, Regular	10-13"/25.5-33 cm	29-31"/74-79 cm	6 pr/bx
MDS160844	Medium, Regular	13-15"/33-38.1 cm	29-31"/74-79 cm	6 pr/bx
MDS160864	Large, Regular	15-17"/38.1-43 cm	29-31"/74-79 cm	6 pr/bx
MDS160884	X-Large, Regular	17-19"/43-48 cm	29-31"/74-79 cm	6 pr/bx
MDS160894	XX-Large, Regular	18-21"/46-53 cm	29-31"/41-45.5 cm	6 pr/bx

*Leg length is the measurement from the bend in the knee to heel **Leg length is the measurement from the gluteal fold to the heel.

ORT113005



Deep Pocket Arm Slings

- » Durable, large pocket sling made of poly/cotton easily accommodates bandages and casts
- » Comes with adjustable wide shoulder strap, thumb loop and hook and loop closure for easy use

Item No.	Description	Pkg
ORT11300XS	X-Small, 10 3/4" L x 6 1/2" W/27.3 x 16.5 cm	1/ea
ORT11300S	Small, 13 3/4" L x 7 1/2" W/34.9 x 19.1 cm	1/ea
ORT11300M	Medium, 16 1/2" L x 8 1/2" W/41.9 x 21.6 cm	1/ea
ORT11300L	Large, 20" L x 9 3/4" W/50.8 x 24.8 cm	1/ea
ORT11300XL	X-Large, 23 1/4" L x 9 3/4" W/59.1 x 24.8 cm	1/ea



ORT11100L



Standard Arm Slings

- » Budget friendly, poly/cotton sling constructed with adjustable strap and buckle

Item No.	Description	Pkg
ORT11100XS	X-Small	1/ea
ORT11100S	Small	1/ea
ORT11100M	Medium	1/ea
ORT11100L	Large	1/ea
ORT11100XL	X-Large	1/ea



ORT13300S



Serpentine Cervical Collars

- » Optimum support with firm density support and contour style
- » Graduated sizing and washable synthetic cover

Item No.	Description	Pkg
ORT13300S	Small, 3 1/2" H x 19" L/ 9 x 48 cm	1/ea
ORT13300M	Medium, 4" H x 21" L/ 10 x 53 cm	1/ea
ORT13300L	Large, 5" H x 23 1/2" L/ 13 x 60 cm	1/ea
ORT13300SXL	Small/X-Large, 4" H x 22" L/10 x 56 cm	1/ea
ORT13300MXL	Medium/X-Large, 4" H x 24" L/10 x 61 cm	1/ea



ORT18000L



10"/25.4 cm, Universal Wrist & Forearm Splints

- » Made with plush foam for optimum comfort
- » Easily adjustable and sewn-in aluminum support

Item No.	Description	Pkg
----------	-------------	-----

Right

ORT18000R	4-11"/10.16-27.9 cm, Circumference	1/ea
-----------	------------------------------------	------

Left

ORT18000L	4-11"/10.16-27.9 cm, Circumference	1/ea
-----------	------------------------------------	------



ORT19400RS



8"/20 cm, Wrist Splints

- » Made with soft and breathable nylon lined foam
- » Slip on design with three loop lock closures
- » Adjustable aluminum support

Item No.	Description	Pkg
----------	-------------	-----

Right

ORT19400RXS	X-Small, 4½-5½"/11.4-14 cm	1/ea
-------------	----------------------------	------

ORT19400RS	Small, 5½-6½"/14-16.5 cm	1/ea
------------	--------------------------	------

ORT19400RM	Medium, 6½-7½"/16.5-19.1 cm	1/ea
------------	-----------------------------	------

ORT19400RL	Large, 7½-8½"/19-21.6 cm	1/ea
------------	--------------------------	------

ORT19400RXL	X-Large, 8½-9½"/21.6-24.1 cm	1/ea
-------------	------------------------------	------

Left

ORT19400LXS	X-Small, 4½-5½"/11.4-14 cm	1/ea
-------------	----------------------------	------

ORT19400LS	Small, 5½-6½"/14-16.5 cm	1/ea
------------	--------------------------	------

ORT19400LM	Medium, 6½-7½"/16.5-19.1 cm	1/ea
------------	-----------------------------	------

ORT19400LL	Large, 7½-8½"/19.1-21.6 cm	1/ea
------------	----------------------------	------

ORT19400LXL	X-Large, 8½-9½"/21.6-24.1 cm	1/ea
-------------	------------------------------	------



ORT27200



Air and Foam Ankle Splints

- » Secure stirrup can be inflated or deflated for desired level of compression
- » Fits right or left ankle and adjusts to many widths

Item No.	Description	Pkg
----------	-------------	-----

ORT27200	Universal, 10"H/25.4 cm	1/ea
----------	-------------------------	------

ORT27200T	Trainer, 9½"H/24.1 cm	1/ea
-----------	-----------------------	------

ORT27200S	Small, 8¾"H/22.2 cm	1/ea
-----------	---------------------	------



ORT28350L



Ankle Lacer

- » Lace-up style brace providing stability
- » For mild to moderate ankle sprains
- » Effortless application and removal
- » Desired fit with lacer style

Item No.	Description	Pkg
----------	-------------	-----

ORT28350S	Small	1/ea
-----------	-------	------

ORT28350M	Medium	1/ea
-----------	--------	------

ORT28350L	Large	1/ea
-----------	-------	------

ORT28350XL	X-large	1/ea
------------	---------	------



ORT28400S



Dorsal Night Splint

- » Designed to hold the foot in a neutral position during sleep for added comfort
- » Helps alleviate morning discomfort associated with plantar fasciitis

Item No.	Description	Pkg
----------	-------------	-----

ORT28400S	Small/Medium	1/ea
-----------	--------------	------

ORT28400L	Large/X-Large	1/ea
-----------	---------------	------



ORT27220



Air/Gel Ankle Stirrup

- » Conforms while limiting ankle rotation
- » Can be used for cold therapy
- » Fits right or left ankle
- » Not made with natural rubber latex

Item No.	Description	Pkg
----------	-------------	-----

ORT27220	Universal, 10"H/25.4 cm	1/ea
----------	-------------------------	------

ORT27220T	Trainer, 9½"H/24.1 cm	1/ea
-----------	-----------------------	------

ORT27220S	Small, 8¾"H/22.2 cm	1/ea
-----------	---------------------	------



ORT30300MM



Semi-Rigid Post-Op Shoes

- » Durable, nylon cushioned rocker sole with a heel and toe pitch for natural gait
- » Molded heel collar secures foot placement and reduces heel slipping
- » Loop lock closure for secure fit
- » Removable tongue allows optimum comfort

Item No.	Description	Pkg
----------	-------------	-----

Men's Relative Shoe Size

ORT30300MS	6-8	1/ea
------------	-----	------

ORT30300MM	8.5-10	1/ea
------------	--------	------

ORT30300ML	10.5-12	1/ea
------------	---------	------

ORT30300MXL	12.5-15	1/ea
-------------	---------	------

Women's Relative Shoe Size

ORT30300WXS	4-6	1/ea
-------------	-----	------

ORT30300WS	6-8	1/ea
------------	-----	------

ORT30300WM	8.5-10	1/ea
------------	--------	------

ORT30300WL	10.5-12	1/ea
------------	---------	------



ORT28200M



ORT28100M

Short & Long Leg Walkers

- » Low profile design helps maintain a normal gait
- » Wide distance between the uprights accommodates extra dressings and swellings
- » Nonskid rocker bottom sole
- » Fits left or right leg

Item No.	Description	Pkg
----------	-------------	-----

Long Leg Walkers

ORT28100XS	X-Small, Women's* 3-5	1/ea
ORT28100S	Small, Men's* 5-6, Women's* 5.5-7	1/ea
ORT28100M	Medium, Men's* 6.5-10, Women's* 7.5-11	1/ea
ORT28100L	Large, Men's* 10.5-12, Women's* 11.5-13	1/ea
ORT28100XL	X-Large, Men's* 12.5-13	1/ea

Short Leg Walkers

ORT28200XS	X-Small, Women's* 3-5	1/ea
ORT28200S	Small, Men's* 5-6, Women's* 5.5-7	1/ea
ORT28200M	Medium, Men's* 6.5-10, Women's* 7.5-11	1/ea
ORT28200L	Large, Men's* 10.5-12, Women's* 11.5-13	1/ea
ORT28200XL	X-Large, Men's* 12.5-13	1/ea

*Sizes according to shoe size.



ORT28620M



ORT28610M

Deluxe Pneumatic Short Leg and Ankle Walkers

- » Pump air into walker for desired compression and fit
- » Low profile design helps maintain a normal gait
- » Wide distance between the uprights accommodates extra dressings and swellings
- » Nonskid rocker bottom sole
- » Fits left or right leg

Item No.	Description	Pkg
----------	-------------	-----

Short Leg Walkers

ORT28610XS	X-Small, Women's* 3-5	1/ea
ORT28610S	Small, Men's* 5-6, Women's* 5.5-7	1/ea
ORT28610M	Medium, Men's* 6.5-10, Women's* 7.5-11	1/ea
ORT28610L	Large, Men's* 10.5-12, Women's* 11.5-13	1/ea
ORT28610XL	X-Large, Men's* 12.5-13	1/ea

Ankle Leg Walkers

ORT28620XS	X-Small, Women's* 3-5	1/ea
ORT28620S	Small, Men's* 5-6, Women's* 5.5-7	1/ea
ORT28620M	Medium, Men's* 6.5-10, Women's* 7.5-11	1/ea
ORT28620L	Large, Men's* 10.5-12, Women's* 11.5-13	1/ea
ORT28620XL	X-Large, Men's* 12.5-13	1/ea

*Sizes according to shoe size.



ORT28640L



ORT28630M

Air Shell Walkers

- » Helps increase stability following post-op procedures
- » Air bladders offer custom compression to increase comfort and reduce edema
- » Rocker sole promotes a natural gait
- » Shock absorbing insole helps reduce impact of the heel strike while walking
- » Toe guard for added protection

Item No.	Description	Pkg
----------	-------------	-----


Short Leg Walker

ORT28630XS	X-Small	1/ea
ORT28630S	Small	1/ea
ORT28630M	Medium	1/ea
ORT28630L	Large	1/ea
ORT28630XL	X-Large	1/ea

Ankle Walker

ORT28640XS	X-Small	1/ea
ORT28640S	Small	1/ea
ORT28640M	Medium	1/ea
ORT28640L	Large	1/ea
ORT28640XL	X-Large	1/ea

**More podiatric supplies are available.
Ask your Medline Sales Representative for our Foot & Ankle catalog.**





POD14201

Gel Hallux Bunion Guard

» Made for the big toes to cushion, protect and help reduce pressure on bunions

Item No.	Description	Pkg
POD14201	One Size Fits Most	1/ea



POD14213L

Gel Toe Separators

» Help align and straighten toes with gentle, constant pressure

Item No.	Description	Pkg
POD14213S	Small	15/pkg
POD14213L	Large	15/pkg



POD14207SM

Ribbed Gel

» Comfortable and stretchable fabric caps fully coated with gel

Item No.	Description	Pkg
POD14203	Small/Medium	1/ea
POD14204	Large/X-Large	1/ea
POD14207SM	Small/Medium	6/pkg
POD14207LXL	Large/X-Large	6/pkg



Knit Gel Toe Cap

» Comfortable and stretchable fabric caps fully coated with gel

Item No.	Description	Pkg
POD14208SM	Small/Medium	6/pkg



Gel Toe Straightener

» A soft gel spacer with two loops for better spacing and positioning

Item No.	Description	Pkg
POD14206	One Size Fits Most	4/pkg



POD14211N



Ribbed Gel Digital Tube

» A comfortable fabric sleeve fully lined with medical grade mineral oil gel to provide all around digital protection, soften to help absorb shock and pressure

Item No.	Description	Pkg
POD14211N	Narrow	3/pkg
POD14211W	Wide	3/pkg



POD14209LXL



Gel Toe Cap

» Helps condition and soften the skin while protecting the toes

Item No.	Description	Pkg
POD14209SM	Small/Medium	10/pkg
POD14209LXL	Large/X-Large	10/pkg



POD142125M



Gel Toe Pads

» A comfortable fabric sleeve with a medical grade mineral oil pad that releases conditioning oils onto the skin to help soften and moisturize

Item No.	Description	Pkg
POD142125M	Small/Medium	12/pkg
POD14212LXL	Large/X-Large	12/pkg



POD14210S



Gel Toe Spacers

» Help align and straighten toes with gentle, constant pressure

Item No.	Description	Pkg
POD14210S	Small	4/pkg
POD14210M	Medium	4/pkg
POD14210L	Large	4/pkg



POD14202



Gel Hammer Toe Crest

» Loops around the second toe; made to cushion, helps protect and reduce foot pain and pressure

Item No.	Description	Pkg
POD14202	One Size Fits Most	1/ea

PHARMACEUTICALS

Contrast	120-122
Inhalation	122-123
Injectables	123-130
Miscellaneous	131
Ophthalmics.....	131-132
Oral Tablets, Capsules & Liquids.....	132-134
Otics.....	135
RX/OTC Topicals & Ointments.....	135-136
Vaccines	136-138



Contrast

Item No.	Description	Qty	Manufacturer
BCC900004	Barium Sulfate, Honey	12/cs	Bracco Diagnostics
BCC900002	Barium Sulfate, Nectar	24/cs	Bracco Diagnostics
BCC900006	Barium Sulfate, Pudding	12/cs	Bracco Diagnostics
BCC450701	Barium Sulfate, Suspension	24/cs	Bracco Diagnostics
BCC900001	Barium Sulfate, Thin Suspension	24/cs	Bracco Diagnostics
BCC900001H	Barium Sulfate, Thin Suspension	1/ea	Bracco Diagnostics
BCC900005	Barium Sulfate, Varibar, Honey	12/cs	Bracco Diagnostics
BCC901702	Barium Thick E-Z-Hd 12-oz.	24/cs	Bracco Diagnostics
BCC901901	Barium Thin E-Z-Paque	24/cs	Bracco Diagnostics
BCC450104	Barium, Read-Cat 2, 450 ml, 2%	24/cs	Bracco Diagnostics
BCC902001	Baro Effervescent Granules	50/cs	Bracco Diagnostics
BCC450301	Baruim Read-Cat 2% 450 ml Berry	24/cs	Bracco Diagnostics
BCC450304	Baruim Readicat 450 ml, Banana	24/cs	Bracco Diagnostics
BCC450305	Baruim Readicat 450 ml, Creamy Vanilla	24/cs	Bracco Diagnostics
BCC450307	Baruim, Readicat 2% Sol, 450 ml, Mochaccino	24/cs	Bracco Diagnostics
LBI095350	Conray, 12 x 150 ml Btl	12/cs	Guerbet
BCC270131652	Contrast Isovue 370 75 ml	10/cs	Bracco Diagnostics
BCC014957	Cysto-Graphin, 300 ml	10/cs	Bracco Diagnostics
BCC141030	Cystograffin Diluent	10/cs	Bracco Diagnostics
BCC902102HH	E-Z-Disk™, 700 mg, Tblt	1/bt	Bracco Diagnostics
BCC902801	E-Z-Paque®, 12-oz. Liquid	24/cs	Bracco Diagnostics
BCC902801H	E-Z-Paque®, 12-oz. Liquid	1/ea	Bracco Diagnostics
BCC901407	Entero, Vu, 600 ml	12/cs	Bracco Diagnostics
50419-320-05	Eovist 0.25 Mol VI 5 x 10 ml	5/bx	Bayer Healthcare
BCC902101	Esophotrast, Paste, E-Z	12/cs	Bracco Diagnostics
50419-325-28	Gadavist Pfs 5 x 10 ml	5/bx	Bayer Healthcare
50419-325-29	Gadavist Pfs 5 x 15 ml	5/bx	Bayer Healthcare
50419-325-27	Gadavist Pfs 5 x 7.5 ml	5/bx	Bayer Healthcare
5041932512CS	Gadavist Sdv 20 x 10 ml	20/cs	Bayer Healthcare
5041932513CS	Gadavist Sdv 20 x 15 ml	20/cs	Bayer Healthcare
5041932511CS	Gadavist Sdv 20 x 7.5 ml	20/cs	Bayer Healthcare
BCC44535	Gastrograffin, 30 ml	24/cs	Bracco Diagnostics
BCC044540	Gastrografen, 120 ml	12/cs	Bracco Diagnostics
BCC044540H	Gastrografen, 120 ml	1/ea	Bracco Diagnostics
NYOY101	Iohexol, Omnipaque, 180 mg, 10 ml	10/bx	GE Healthcare Technologies
NYOY102	Iohexol, Omnipaque, 180 mg, 20 ml	10/bx	GE Healthcare Technologies
NYOY203	Iohexol, Omnipaque, 240 mg, 10 ml	10/bx	GE Healthcare Technologies
NYOY524	Iohexol, Omnipaque, 240 mg, 150 M	10/bx	GE Healthcare Technologies
NYOY306	Iohexol, Omnipaque, 300 mg, 10 ml	10/bx	GE Healthcare Technologies
NYOY220	Iohexol, Omnipaque, 240 mg, 200 ml	10/bx	GE Healthcare Technologies
NYOY542	Iohexol, Omnipaque, 350 mg, 100 ml	10/bx	GE Healthcare Technologies
NYOY420	Iohexol, Omnipaque, 350 mg, 125 ml	10/bx	GE Healthcare Technologies
NYOY548B	Iohexol, Omnipaque, 350 mg, 500 ml	10/bx	GE Healthcare Technologies
NYO548Y	Iohexol, Omnipaque, 350 mg, 500 ml	10/bx	GE Healthcare Technologies
NYOY540	Iohexol, Omnipaque, 350 mg, Plastic	10/bx	GE Healthcare Technologies
NYOY530	Iohxl, Omnipaque, 300 mg, 50 ml, Plastic	10/bx	GE Healthcare Technologies
NYO532Y	Iohxl, Omnpqu, 300 mg, 100 ml, Nova	10/bx	GE Healthcare Technologies
NYO544Y	Iohxl, Omnpqu, Nova, Plastic	10/bx	GE Healthcare Technologies
BCC141125	Iopamidol, Isovue-M, 200, 20 ml	10 bt/cs	Bracco Diagnostics
BCC141215H	Iopamidol, Isovue-M, 300, 15 ml	1 bt/ea	Bracco Diagnostics
BCC141215	Iopamidol, Isovue-M, 300, 15 ml	10 bt/cs	Bracco Diagnostics

Item No.	Description	Qty	Manufacturer
BCC131709	Iopamidol, Isovue, 150 ml	10/bx	Bracco Diagnostics
BCC131535	Iopamidol, Isovue, 300, 100 ml	10/cs	Bracco Diagnostics
BCC131550	Iopamidol, Isovue, 300, 150 ml	10/cs	Bracco Diagnostics
BCC131530	Iopamidol, Isovue, 300, 50 ml	10/cs	Bracco Diagnostics
BCC131598	Iopamidol, Isovue, 300, 500 ml	6/cs	Bracco Diagnostics
BCC131635	Iopamidol, Isovue, 370, 100 ml	10/cs	Bracco Diagnostics
BCC131637	Iopamidol, Isovue, 370, 150 ml	10/cs	Bracco Diagnostics
BCC0270131630	Iopamidol, Isovue, 370, 50 ml	10/cs	Bracco Diagnostics
BCC131698	Iopamidol, Isovue, 500 ml	6/bx	Bracco Diagnostics
BCC131641	Iopamidol, Isovue, Multi Pack	10/cs	Bracco Diagnostics
LBI318350	Iothalamate Meglumine, Usp 43%	12/cs	Guerbet
LBI318305	Iothalamate Meglumine, Usp 43%	25/cs	Guerbet
LBI095305	Iothalamate Meglumine, Usp 60%, 25 x 50 ml	25/cs	Guerbet
LBI132406	Ioversol, 51%, 240 bt, 25 x 50 ml	25/cs	Guerbet
LBI132306	Ioversol, Optiray, 68%, 320 bt, 25 x 50 ml	25/cs	Guerbet
LBI133306	Ioversol, Optiray, 74%, 350 mg/ml	25/cs	Guerbet
LBI133311	Ioversol, Optiray, 74%, 100 ml	12/cs	Guerbet
BCC131525	Isovue 300, 30 ml	10/cs	Bracco Diagnostics
BCC131547	Isovue 300, 75 ml	10/cs	Bracco Diagnostics
BCC131604	Isovue 370, 125 ml	10/cs	Bracco Diagnostics
BCC131541	Isovue Multipack -300 200 ml	10/cs	Bracco Diagnostics
BCC141111	Isovue-M 200, VI 41% 10 ml	10/bx	Bracco Diagnostics
BCC131702	Isovue, 250 100 ml Vials	10/cs	Bracco Diagnostics
BCC131705	Isovue, 250 50 ml Vials	10/cs	Bracco Diagnostics
50419-188-36	Magnevist Inj Syr 5 x 10 ml	5/bx	Bayer Healthcare
50419-188-37	Magnevist Inj Syr 5 x 15 ml	5/bx	Bayer Healthcare
50419-188-38	Magnevist Inj Syr 5 x 20 ml	5/bx	Bayer Healthcare
50419-188-58	Magnevist Pbp Vial 10 x 50 ml	10/bx	Bayer Healthcare
50419-188-11	Magnevist Pbp VI 10 x 100 ml	10/bx	Baxter Healthcare
50419-188-01	Magnevist VI 10 ml	1/ea	Bayer Healthcare
50419-188-15	Magnevist VI 15 ml	1/ea	Bayer Healthcare
50419-188-05	Magnevist VI 5 ml	1/ea	Baxter Healthcare
LBI481604	Meglumine Diatrizoate, 30 ml	25/cs	Guerbet
BCC516412BX	Multihance 5 ml	5/bx	Bracco Diagnostics
BCC516413BX	Multihance 10 ml	5/bx	Bracco Diagnostics
BCC516414BX	Multihance 15 ml	5/bx	Bracco Diagnostics
BCC516415BX	Multihance 20 ml	5/bx	Bracco Diagnostics
BCC526417	Multihance, Multipack, 100 ml	5/cs	Bracco Diagnostics
BCC526416	Multihance, Multipack, 5 x 50 ml	5/bx	Bracco Diagnostics
NYOY532	Omnipaque 100 ml 300 mg, Plastic	10/bx	GE Healthcare Technologies
NYOY522	Omnipaque 240 mg 100 ml-Plastic Bottle	10/bx	GE Healthcare Technologies
NYOY534	Omnipaque 300 mg 150 ml Pls Pk, Plastic	10/bx	GE Healthcare Technologies
NYOY538B	Omnipaque 300 mg 500 ml Plus Pak	10/bx	GE Healthcare Technologies
NYO420Y	Omnipaque 350 mg 125 ml (Nova)	10/bx	GE Healthcare Technologies
NYOY546	Omnipaque 350 mg 200 ml Plus Pak	10/bx	GE Healthcare Technologies
NYOY526	Omnipaque, 240 mg, 200 ml, Bottle	10/cs	GE Healthcare Technologies
NYOY503	Omnipaque, 300 mg, 30 ml In 50 ml +Pluspak	10/bx	GE Healthcare Technologies
NYOY510	Omnipaque, 50 ml, D, Rx	10/bx	GE Healthcare Technologies
NYO520Y	Omnipaque, Novaplus, Pluspak 240 mg, 50 ml	10/bx	GE Healthcare Technologies
NYO540Y	Omnipaque, Novatn, Plastic	10/bx	GE Healthcare Technologies
NYOJ100	Omniscan 287 mg 10 ml	10/bx	GE Healthcare Technologies

Contrast Continued

Item No.	Description	Qty	Manufacturer
NYOJ120	Omniscan 287 mg 15 ml	10/bx	GE Healthcare Technologies
NYOJ140	Omniscan 287 mg 20 ml	10/bx	GE Healthcare Technologies
NYOJ640	Omniscan Prefilled Plus 15 ml	10/bx	GE Healthcare Technologies
NYO140J	Omniscan, 287 mg, 20 ml, (Nova)	10/bx	GE Healthcare Technologies
LBI132411	Optiray 240, 12 x 100 ml	12/cs	Guerbet
LBI132487	Optiray 240, 20 x 125 ml Pi Syr	20/cs	Guerbet
LBI133200	Optiray 300, 100 ml Pi Syrg Rfid	20/cs	Guerbet
LBI133221	Optiray 300, 200 ml, Bottle	12/cs	Guerbet
LBI133206	Optiray 300, 50 ml Bottle	25/cs	Guerbet
LBI132321	Optiray 320 12 x 200 ml Btl	12/cs	Guerbet
LBI132327	Optiray 320 125 ml Pi Syrg Rfid	20/cs	Guerbet
LBI132316	Optiray 320 150 ml Bottle	12/cs	Guerbet
LBI132390	Optiray 320 20 x 100 ml Pi Syr	20/cs	Guerbet
LBI132387	Optiray 320 20 x 125 ml Pi Syr	20/cs	Guerbet
LBI132352	Optiray 320 20 x 50 ml Powr Inject	20/cs	Guerbet
LBI132395	Optiray 320 20 x 75 ml Pi Syr	200/cs	Guerbet
LBI132361	Optiray 320, 6 x 500 ml Pbp Multi-Dose	6/cs	Guerbet
LBI133316	Optiray 350 12 x 150 ml Btl	12/cs	Guerbet
LBI133321	Optiray 350 12 x 200 ml Btl	12/cs	Guerbet
LBI133387	Optiray 350 20 x 125 ml Pi	20/cs	Guerbet
LBI133395	Optiray 350 20 x 75 ml Pi Syr	20/cs	Guerbet
LBI133390	Optiray 350, 20 x 100 ml, Pi	20/cs	Guerbet
LBI133361	Optiray 350, 6 x 500 ml Pbp Multi-Dose	6/cs	Guerbet
LBI133216	Optiray-300 12 x 150 ml Bottle	12/cs	Guerbet
BCC111104BX	Prohance Gadolinium, Nonionic, 5 ml	5/bx	Bracco Diagnostics
BCC111101BX	Prohance, Gadoteridol, 10 ml, Vial	5/bx	Bracco Diagnostics
BCC111102BX	Prohance, Gadoteridol, 15 ml, Vial	5/bx	Bracco Diagnostics
BCC111103BX	Prohance, Gadoteridol, 20 ml, Vial	5/bx	Bracco Diagnostics
50419-342-10	Ultravist 240 mg/ml Vial 10 x 100 ml	10/bx	Bayer Healthcare
50419-344-58	Ultravist 300 mg/ml Pbp Vial 8 x 500 ml	8/bx	Bayer Healthcare
50419-344-10	Ultravist 300 mg/ml Vial 10 x 100 ml	10/bx	Bayer Healthcare
50419-344-15	Ultravist 300 mg/ml Vial 10 x 150 ml	10/bx	Bayer Healthcare
50419-344-05	Ultravist 300 mg/ml Vial 10 x 50 ml	10/bx	Bayer Healthcare
50419-346-58	Ultravist 370 mg/ml Pbp Vial 8 x 500 ml	8/bx	Bayer Healthcare
50419-346-10	Ultravist 370 mg/ml V1 10 x 100 ml	10/bx	Bayer Healthcare
NYOV552	Visipaque 270 mg, Plastic, 10 x 100 ml	10/bx	GE Healthcare Technologies
NYOV554	Visipaque 270 mg, Plastic, 10 x 150 ml	10/bx	GE Healthcare Technologies
NYOV562	Visipaque 320 mg 100 ml	10/bx	GE Healthcare Technologies
NYOV568B	Visipaque 320 mg 500 ml	10/bx	GE Healthcare Technologies
NYOV564	Visipaque 320 mg, 150 ml, Plastic	10/bx	GE Healthcare Technologies
NYOV566	Visipaque 320 mg, Plastic, 10 x 200 ml	10/bx	GE Healthcare Technologies
NYOV560	Visipaque 320 mg/ml, 50 ml	10/bx	GE Healthcare Technologies

Inhalation

Item No.	Description	Qty	Manufacturer
0487-9501-60	Albuterol Inh 0.83 mg/ml Sol 60 x 3 ml	60/bx	Nephron Pharmaceuticals
0487-9501-25	Albuterol Inh, 0.83 mg/ml Sol, 3 ml	25/bx	Nephron Pharmaceuticals
50383-741-20	Albuterol Inh, 5 mg/ml Sol, 20 ml	1/ea	Hi-Tech Pharmaco Co. Inc

Item No.	Description	Qty	Manufacturer
0487-9904-25	Albuterol Sulf 1.25 mg/3 ml Inh Sol 3 ml	25/bx	Nephron Pharmaceuticals
76204-600-30	Albuterol/Ipratropium Inh Soln 3 ml	30/bx	Ritedose Pharmaceuticals Llc
10019-360-60	Forane Inh 6 x 250 ml	6/bx	Baxter Healthcare
10019-360-40	Forane Liq 6 x 100 ml	6/bx	Baxter Healthcare
0487-9801-30	Ipratrop Brom, 0.02% Sol, 30 x 2.5 ml	30/bx	Nephron Pharmaceuticals
66794-017-10	Isoflurane Inh 100 ml	1/ea	Piramal Critical Care Inc
66794-017-25	Isoflurane Inh 250 ml	1/ea	Piramal Critical Care Inc
66993-022-27	Levalbuterol 0.63 mg/3 ml Inh Soln 24 x 3 ml	24/bx	Prasco Laboratories
66993-023-27	Levalbuterol 1.25 mg/3 ml Sol 24 x 3 ml	24/bx	Prasco Laboratories
63323-877-15	Nebupent 300 mg Sdv (For Inh Soln)	1/ea	APP Pharm
0186-1990-04	Pulmicort Res 1 mg/2ml Res 30 x 2 ml	30/bx	Astra Zeneca LP
0487-5901-99	S-2 Racemie E 2.25% Brc Sol 30 x 0.5 ml U	30/bx	Nephron Pharmaceuticals
10019-651-64	Sevoflurane Inh 6 x 250 ml	6/bx	Baxter Healthcare
0378-6987-89	Sod Chl 0.9 % Sol 50 x 15 ml	50/bx	Mylan
0487-9301-03	Sodium Chloride, 0.9%, 3 ml	100/bx	Nephron Pharmaceuticals
10019-641-34	Suprane Inh Liq 6 x 240 ml	6/bx	Baxter Healthcare
0173-0682-20	Ventolin HFA INH 18 gm (200 Actuations)	1/ea	Glaxo Smithkline Consumer



0009-0280-03



0409-0805-01



0409-1746-30H



0409-1761-02



0409-3182-01



61703-324-18



0409-3182-02



0409-4278-01

Injectables

Item No.	Description	Qty	Manufacturer
39822-0190-1	Acetazolamide 500 mg Vial	1/ea	X-Gen Pharmaceuticals
0409-3308-03	Acetylcysteine 20% Sol TTV 3 x 30 ml	3/bx	Hospira Worldwide Inc
50242-044-13	Activase 50 mg Vial	1/ea	Genentech, Inc
17478-542-02	Adenosine 3 mg/ml Vial 10 x 2 ml	10/bx	Akorn Inc
0409-1932-02	Adenosine 3 mg/ ml Sdv 30 ml	1/ea	Pfizer Inc/Consumr Hlthcr
42023-168-01	Adrenalin 1 mg/ml Mdv 30 ml	1/ea	Par Pharmaceutical Companies, Inc
17478-503-05	Alcohol Dehydrated 98% Vial 5 ml	10/bx	Akorn Inc
0409-8062-01	Amidate 2 mg/ml Amp 5 x 10 ml	5/bx	Hospira Worldwide Inc
0409-6695-01	Amidate 2 mg/ml SDV 10 x 10 ml	10/bx	Hospira Worldwide Inc
0409-6695-02	Amidate 2 mg/ml SDV 10 x 20 ml	10/bx	Hospira Worldwide Inc
52584-921-01	Aminophylline 25 mg/ml SDV 10 ml	1/ea	Hospira Worldwide Inc
0409-5921-01	Aminophylline 25 mg/ml SDV 25 x 10 ml	25/bx	Hospira Worldwide Inc
0409-5922-01	Aminophylline 25 mg/ml SDV 25 x 20 ml	25/bx	Hospira Worldwide Inc


Injectables Continued



61703-341-06



63323-545-01



76329-1469-5



0517-0031-25



0641-0928-25



0409-4688-23



0409-6482-01



61703-344-18

Item No.	Description	Qty	Manufacturer
63323-616-09	Amiodarone 50 mg/ml Vial 10 x 9 ml	10/bx	APP Pharm
55150-114-20	Ampicillin 2 gm Vial	10/bx	Auromedics
63323-708-00	Ampicillin 500 mg Vial	10/bx	Fresenius Kabi
67457-118-50	Ascorbic Acid 500 mg/ml Vial 50 ml	1/ea	Bioniche Life Sciences Inc
0409-4911-34	Atropine Sul 0.1 mg/ml PFs 10 x 10 ml	10/bx	Hospira Worldwide Inc
0409-4910-34	Atropine Sulf 0.1 mg/ml PFs 10 x 5 ml	10/bx	Hospira Worldwide Inc
0517-1010-25	Atropine Sulfate, 1 mg/ml, 1 ml Vial	25/bx	American Regent Inc
52584-044-25	Atropine Sulfate, 1 mg/ml, 1 ml Vial	1/ea	American Regent Inc
63323-329-31	Bacitracin 50,000U/Vial	10/bx	Fresenius Kabi
58914-080-52	Bentyl 10 mg Amp 5 x 2 ml	5/bx	Aptalis Pharma Us Inc
0517-0720-01	Betamethasone, 6 mg/ml, MDV, 5 ml	1/ea	American Regent Inc
60793-600-10	Bicillin Cr 1200 Mu A PFs, 10 x 2 ml	10/bx	King Health Care
60793-602-10	Bicillin Cr 900/ 300 P Tbx 10 x 2 ml	10/bx	Pfizer Laboratories Div Pfizer Inc.
60793-701-10	Bicillin La 1200 Mu PFs 10 x 2 ml	10/bx	Pfizer Inc/Consumr Hlthcr
60793-702-10	Bicillin La 2400Mu PFs 10 x 4 ml	10/bx	Pfizer Inc/Consumr Hlthcr
60793-700-10F	Bicillin La 600 Mu PFs 10 x 1 ml	10/cs	Pfizer Inc/Consumr Hlthcr
0004-0191-09	Boniva 3 mg PFs 3 ml	1/ea	Genentech, Inc
0409-1412-10	Bumetanide 0.25 mg/ml MDV 10 x 10 ml	10/bx	Hospira Worldwide Inc
0409-1412-04	Bumetanide 0.25 mg/ml SDV 10 x 4 ml	10/bx	Hospira Worldwide Inc
0409-1160-01	Bupivacaine 0.25% MDV 25 x 50 ml	25/bx	Hospira Worldwide Inc
0409-1159-01	Bupivacaine 0.25% PF SDV 25 x 10 ml	25/bx	Hospira Worldwide Inc
0409-1159-02	Bupivacaine 0.25% PF SDV 25 x 30 ml	25/bx	Hospira Worldwide Inc
0409-9043-01	Bupivacaine 0.25%+Epi 1:200,000 FTV 25 x 50	25/cs	Hospira Worldwide Inc
0409-9042-01	Bupivacaine 0.25%+Epi 1:200000 SDV 10 x 10 ml	10/bx	Hospira Worldwide Inc
0409-1163-01	Bupivacaine 0.5% MDV 25 x 50 ml	25/bx	Hospira Worldwide Inc
0409-1162-01	Bupivacaine 0.5% PF SDV 25 x 10 ml	25/bx	Hospira Worldwide Inc
0409-1162-02	Bupivacaine 0.5% PF SDV 25 x 30 ml	25/bx	Hospira Worldwide Inc
0409-9046-01	Bupivacaine 0.5%+Epi 1:200,000 FTV 25 x 50 ml	25/cs	Hospira Worldwide Inc
17478-931-01	Calcitriol 1 mcg/ml Amp 1 ml	10/bx	Akorn Inc
0409-1631-10	Calcium Chloride 100 mg/ml PFs 10 x 10	10/bx	Hospira Worldwide Inc

Item No.	Description	Qty	Manufacturer
0409-4928-34	Calcium Chloride 100 mg/ml PFs 10 x 10	10/bx	Hospira Worldwide Inc
59584-138-01	Candin Mdv 1 ml	1/ea	Nielsen Biosciences, Inc
0409-1038-50	Carbocaine 1% FTV 1 x 50 ml	1/ea	Hospira Worldwide Inc
0409-1036-30	Carbocaine 1% PF SDV 1 x 30 ml	1/ea	Hospira Worldwide Inc
0409-1041-30	Carbocaine 1.5% PF SDV 1 x 30 ml	1/ea	Hospira Worldwide Inc
50242-041-64	Cathflo Activase 2 mg Vial 2 ml	1/ea	Genentech, Inc
0143-9924-90	Cefazolin (Ancef), 1 gm 10 ml Vial	25/bx	West-Ward Pharm Corp
52584-023-10	Cefazolin (Ancef), 1 gm 15 ml Vial	1/ea	West-Ward Pharm Corp
0409-0805-01	Cefazolin 1G SDV 25 x 1G	25/bx	Hospira Worldwide Inc
0143-9923-90	Cefazolin, 500 mg Vial	25/bx	West-Ward Pharm Corp
25021-121-20	Cefepime 1 gm Vial	10/bx	Sagent Pharmaceuticals
25021-109-10	Cefoxitin 1 gm 10 Vial	10/bx	Sagent Pharmaceuticals
25021-110-20	Cefoxitin 2 gm 10 Vial	10/bx	Sagent Pharmaceuticals
25021-106-10	Ceftriaxone 1 gm Vial 25	25/bx	Sagent Pharmaceuticals
0409-7334-10	Ceftriaxone 10G FTV 1 x 10G	1/ea	Hospira Worldwide Inc
0409-7332-01	Ceftriaxone 1G SDV 10 x 1G	10/bx	Hospira Worldwide Inc
0781-3206-85	Ceftriaxone 250 mg SDV	1/ea	Sandoz
0409-7337-01	Ceftriaxone 250 mg SDV 10 x 250	10/bx	Hospira Worldwide Inc
0409-7337-01H	Ceftriaxone 250 mg SDV 10 x 250	1/ea	Hospira Worldwide Inc
0409-7335-03	Ceftriaxone 2G SDV 10 x 2G	10/bx	Hospira Worldwide Inc
0409-7338-01	Ceftriaxone 500 mg SDV 10 x 500	10/bx	Hospira Worldwide Inc
0143-9857-25	Ceftriaxone, 1 gm Vial	25/bx	West-Ward Pharm Corp
0641-1397-35	Chlorpromazine 25 mg/ml Amp 25 x 1 ml	25/bx	Baxter Healthcare
0409-4777-02	Ciprofloxacin 2 mg/ml 5% Dex 24 x 200 ml	24/cs	Hospira Worldwide Inc
0009-0775-26	Cleocin Phosphate 600 Mg SDV 4 ml	25/bx	Pfizer Inc/Consumr Hlthcr
63323-282-60	Clindamycin 150 mg/ml Vial 60 ml	1/ea	APP Pharm
25021-115-02	Clindamycin 150 mg/ml SDV 2 ml	25/bx	Sagent Pharmaceuticals
25021-115-06	Clindamycin 150 mg/ml SDV 6 ml	25/bx	Sagent Pharmaceuticals
63323-282-04	Clindamycin 600 mg SDV 25 x 4 ml	25/bx	APP Pharm
52584-399-10	Cyanocobalamin 1,000 mcg/ml MDV 10 ml	1/ea	West-Ward Pharmaceuticals
0517-0031-25	Cyanocobalamin 1,000 mcg/ml 25 x 1 ml	25/bx	American Regent Inc
0517-0032-25	Cyanocobalamin 1,000 mcg/ml 25 x 10 ml	25/bx	American Regent Inc
52584-080-01	Cyanocobalamin 1,000 mcg/ml 30 ml	1/ea	American Regent Inc
0517-0130-054	Cyanocobalamin 1,000 mcg/ml 5 x 30 ml	5/bx	American Regent Inc
63323-044-01	Cyanocobalamin 1,000 mcg/ml Vial 25 x 1 M	25/bx	APP Pharm
10019-955-01	Cyclophosphamide 500 mg SDV	1/ea	Baxter Healthcare
42023-123-06	Dantrium lv 20 mg/Vial	6/bx	JHP Pharmaceuticals
42023-111-01	Delestrogen, 20 mg/ml Vial, 5 ml	1/ea	JHP Pharmaceuticals
0009-0271-01	Depo Estradio, 5 mg/ml,Vial, 5 ml	1/ea	Pharmacia Deltec Inc
0009-3073-03	Depo Medrol 40 mg/ml Vial 25 x 1 ml	25/bx	Pfizer Inc/Consumr Hlthcr
0009-0274-01	Depo Medrol, 20 mg/ml, Vial, 5 ml	1/ea	Pfizer Inc/Consumr Hlthcr
0009-3073-01	Depo Medrol, 40 mg, 1 ml, SDV (5-6 Mos Exp)	1/ea	Pfizer Inc/Consumr Hlthcr
0009-0280-02	Depo Medrol, 40 mg/ml, 5 ml, Vial	1/ea	Pfizer Inc/Consumr Hlthcr
0009-0626-01	Depo Provera 400 mg Vial 2.5 ml	1/ea	Pfizer Inc/Consumr Hlthcr
0009-0280-03	Depo-Medrol 40 mg/ml MDV 10 ml	1/ea	Pfizer Inc/Consumr Hlthcr
0009-0280-51	Depo-Medrol 40 mg/ml Vial 25 x 5 ml	25/bx	Pfizer Inc/Consumr Hlthcr
0009-3475-03	Depo-Medrol 80 mg/ml 1 ml 25 x 1 ml Vial	25/bx	Pfizer Inc/Consumr Hlthcr
0009-3475-01	Depo-Medrol 80 mg/ml Sdv 1 ml (5-6 Mos Exp)	1/ea	Pfizer Inc/Consumr Hlthcr
0009-0306-02	Depo-Medrol, 80 mg/ml, 5 ml, MDV	1/ea	Pfizer Inc/Consumr Hlthcr
0009-0746-30	Depo-Provera 150 mg/ml 1 ml SDV	1/ea	Pfizer Inc/Consumr Hlthcr
0009-4709-13	Depo-Subq Pro 104 mg PFs 0.65 ml	1/ea	Pfizer Inc/Consumr Hlthcr

Injectables Continued

Item No.	Description	Qty	Manufacturer
0641-0367-25	Dexamethasone 10 mg/ml 25 x 1 ml	25/bx	Baxter Healthcare
63323-506-01	Dexamethasone 10 mg/ml PF SDV 25 x 1 ml	25/bx	APP Pharm
63323-516-10	Dexamethasone 10 mg/ml Vial 10 x 10 ml	10/bx	APP Pharm
63323-165-05	Dexamethasone 4 mg/ml 25 x 5 ml Vial	25/bx	APP Pharm
63323-165-30	Dexamethasone 4 mg/ml 30 ml Vial	1/bx	APP Pharm
55150-237-01	Dexamethasone 4 mg/ml Sdv 1 ml	25/bx	Auromedics
0409-7517-16	Dextrose 50% Ansyrr Syringe 10 x 50 ml	10/bx	Hospira Worldwide Inc
0409-4902-34	Dextrose 50% PFs 10 x 50 ml	10/bx	Hospira Worldwide Inc
52584-902-34	Dextrose 50% PFs 50 ml	1/ea	Hospira Worldwide Inc
0574-0850-05	Dihydroergot 1 mg/ml Amp 5 x 1 ml	5/bx	Paddock Laboratories Inc
0574-0850-10	Dihydroergotamine 1 mg/ml Amp 10 x 1 ml	10/bx	Paddock Laboratories Inc
0641-6013-10	Diltiazem 5 mg/ml SDV 10 x 5 ml	10/bx	Baxter Healthcare
0409-1171-01	Diltiazem HCL 5 mg/ml SDV 10 x 5 ml	10/bx	Hospira Worldwide Inc
0409-2290-31	Diphenhydramine HCL 50 mg/ml PFs 10 x 1 ml	10/bx	Hospira Worldwide Inc
0641-0376-25	Diphenhydramine, 50 mg/ml, 1ml	25/bx	Baxter Healthcare
52584-664-01	Diphenhydramine, 50 mg/ml, 1ml	1/ea	Baxter Healthcare
0641-2569-44	Dipyridamole, 5 mg/ml, 10 ml, SDV	5/bx	Baxter Healthcare
0409-2347-32	Dobutamine +5% Dext 500 mg FC 12 x 250 ml	12/cs	Hospira Worldwide Inc
0409-2344-02	Dobutamine 12.5 mg/ml Sdv 20 ml	10/bx	Pfizer Inc/Consumr Hlthcr
0338-1073-02	Dobutamine 250 mg Sol 18 x 250 ml	18/bx	Baxter Healthcare
0409-9104-20	Dopamine HCL 40 mg/ml SDV 25 x 10 ml	25/bx	Hospira Worldwide Inc
0409-5820-01	Dopamine HCL 40 mg/ml SDV 25 x 5 ml	25/bx	Hospira Worldwide Inc
0409-2122-01	Enalaprilat 1.25 mg/ml FTV 1 x 1 ml	1/ea	Hospira Worldwide Inc
0409-2122-02	Enalaprilat 1.25 mg/ml FTV 1 x 2 ml	1/ea	Hospira Worldwide Inc
0781-3224-64	Enoxaparin 40 mg Syg 10	10/bx	Sandoz
0409-4921-34	Epinephrine 0.1 mg/ml PFs 10 x 10 ml	10/bx	Hospira Worldwide Inc
52584-004-34	Epinephrine 0.1 mg/ml PFs 10 ml	1/ea	Pfizer Inc/Consumr Hlthcr
54288-103-10	Epinephrine 1 mg/ml Amp 10 x 1 ml	10/bx	BPI Labs
49502-500-02	Epipen 0.3 mg 2-Pak (<12 Mo Exp Date)	1/ea	Dey Pharma
49502-501-02	Epipen Jr 0.15 mg 2-Pak (<12 Mo Exp Date)	1/ea	Dey Pharma
55566-4100-1	Euflexxa Syg 3 x 2 ml	3/bx	Ferring Pharmaceuticals
0641-6022-25	Famotidine 10 mg/ml 2 ml SDV 25 x 2 ml	25/bx	West-Ward Pharm Corp
0024-2792-10	Ferrlecit 62.5 mg Vial 10 x 5 ml	10/bx	Sanofi-Adventis
0143-9899-91	Fluconazole 200 mg Vial 6 x 100 ml	6/bx	West-Ward Pharm Corp
0409-4688-02	Fluconazole+Sod Cl 2 mg/ml Fc 200 ml	6/bx	Pfizer Inc/Consumr Hlthcr
63323-424-10	Flumazenil 0.1 mg/ml MDV 10 ml	1/ea	APP Pharm
0143-9784-10	Flumazenil 0.1 mg/ml Vial 10 x 5 ml	10/bx	West-Ward Pharm Corp
63323-117-10	Fluorouracil 500 mg PF Vial 10 x 10 ml	10/bx	APP Pharm
0409-1639-10	Furosemide 10 mg/ml PFs 10 x 10 ml	10/bx	Hospira Worldwide Inc
0409-9631-04	Furosemide 10 mg/ml PFs 10 x 4 ml	10/bx	Hospira Worldwide Inc
0409-6102-10	Furosemide 10 mg/ml SDV 25 x 10 ml	25/bx	Hospira Worldwide Inc
0409-6102-02	Furosemide 10 mg/ml SDV 25 x 2 ml	25/bx	Hospira Worldwide Inc
0409-6102-04	Furosemide 10 mg/ml SDV 25 x 4 ml	25/bx	Hospira Worldwide Inc
52584-102-02	Furosemide 10 mg/ml Sdv 2 ml	1/ea	Pfizer Inc/Consumr Hlthcr
63323-010-02	Gentamicin 40 mg/ml MDV 25 x 2 ml	25/bx	APP Pharm
0409-1207-03	Gentamicin 40 mg/ml SDV 25 x 2 ml	25/bx	Hospira Worldwide Inc
0002-8031-01	Glucagon Emerg 1 mg Kit 1 ml	1/ea	Eli Lilly and Company
0069-0058-01	Heparin 1,000U/ml MDV 25 x 10 ml	25/bx	Pfizer Inc/Consumr Hlthcr
0409-2720-03	Heparin 1,000U/ml FTV 25 x 30 ml	25/bx	Hospira Worldwide Inc
0409-2720-01	Heparin 1,000U/ml SDV 25 x 1 ml	25/bx	Hospira Worldwide Inc
0409-2721-01	Heparin 10,000U/ml SDV 25 x 1 ml	25/bx	Hospira Worldwide Inc

Item No.	Description	Qty	Manufacturer
0409-7650-62	Heparin 100U/ml+0.45% Sod Chl 24 x 250 ml	24/cs	Hospira Worldwide Inc
0409-7793-62	Heparin 100U/ml+5% Dex FC 24 x 250 ml	24/cs	Hospira Worldwide Inc
0409-7620-03	Heparin 2U/ml+0.9 Sod FC 18 x 500 ml	18/cs	Hospira Worldwide Inc
0409-2723-01	Heparin 5,000U/ml SDV 25 x 1 ml	25/bx	Hospira Worldwide Inc
63323-545-01	Heparin Lf 100 Un/ml Vial 25 x 1 ml	25/bx	APP Pharm
63323-545-05	Heparin Lock Flush 100Un/ml VI 5 ml	25/bx	Fresenius Kabi
63323-0542-07	Heparin Sod 10 M Un/ml MDV 25 x 5 ml	25/bx	APP Pharm
63323-0915-01	Heparin Sod 20 M Un MDV 25 x 1 ml	25/bx	APP Pharm
64253-333-35	Heparin, 100U /ml, 12 Syg, 60 x 5 ml	60/bx	Medefil
0002-8715-01	Humulin 70/ 30 Vial 10 ml	1/ea	Eli Lilly And Company
0002-8315-01	Humulin N 100U /ml Vial 10 ml	1/ea	Eli Lilly And Company
0002-8215-01	Humulin R 100U /ml Vial 10 ml	1/ea	Eli Lilly And Company
89122-724-20	Hyalgan 20 mg Syg 2 ml	1/ea	Fidia Pharma USA Inc
0591-2888-30	Hydroxocobalamin 1,000 mcg/ml MDV 30 ml	1/ea	Watson Pharmaceuticals
0517-5601-25	Hydroxyzine 50 mg/ml Vial 25 x 1 ml	25/bx	American Regent Inc
18657-117-04	Hylenex 150U/ml VI 1 ml	4/bx	Halozyme Therapeutics, Inc
17478-701-02	Ic-Green 25 mg Kit (6 Vials with Solvent)	1/bx	Akorn Inc
0006-3843-71	Invanz 1 Gm Vial	10/bx	Merck & Co Inc
0003-0494-20	Kenalog-10, 10 mg/ml, 5 ml Vial	1/ea	Bristol-Myers Squibb
0003-0293-28	Kenalog-40 40 mg/ml 10 ml Vial	1/ea	Bristol-Myers Squibb
0003-0293-20	Kenalog-40 40 mg/ml 5 ml Vial	1/ea	Bristol-Myers Squibb
0003-0293-05	Kenalog-40, 40 mg/ml, 1 ml Vial	1/ea	Bristol-Myers Squibb
52584-584-93	Ketorolac Tro 30 mg/ml SDV 1 ml	1/ea	Pfizer Inc/Consumr Hlthcr
0409-3793-01	Ketorolac Trometh 15 mg/ml SDV 25 x 1 ml	25/bx	Hospira Worldwide Inc
0409-2287-31	Ketorolac Trometh 30 mg/ml PFs 10 x 1 ml	10/bx	Hospira Worldwide Inc
0409-3795-01H	Ketorolac Trometh 30 mg/ml SDV 1 ml	1/ea	Hospira Worldwide Inc
0409-3795-01	Ketorolac Trometh 30 mg/ml SDV 25 x 1 ml	25/bx	Hospira Worldwide Inc
0409-3796-01H	Ketorolac Trometh 60 mg/2ml SDV 2 ml	1/ea	Hospira Worldwide Inc
0409-3796-01	Ketorolac Trometh 60 mg/2ml SDV 25 x 2 ml	25/bx	Hospira Worldwide Inc
0409-2267-20	Labetalol HCL 5 mg/ml FTV 1 x 20 ml	1/ea	Hospira Worldwide Inc
0409-2267-54	Labetalol HCL 5 mg/ml FTV 1 x 40 ml	1/ea	Hospira Worldwide Inc
0409-2339-34	Labetalol HCL 5 mg/ml PFs 10 x 4 ml	10/bx	Hospira Worldwide Inc
17478-107-20	Levofloxacin 25 mg/ml Vial 20 ml	1/ea	Akorn Inc
17478-107-30	Levofloxacin 25 mg/ml Vial 30 ml	1/ea	Akorn Inc
0409-1443-04	Levophed 1 mg/ml Amp 10 x 4 ml	10/bx	Hospira Worldwide Inc
0469-6501-89	Lexiscan 0.4 mg/5 ml Syr	1/ea	Astellas Pharma Us, Inc
0409-4278-01	Lidocaine 0.5% PF SDV 25 x 50 ml	25/bx	Hospira Worldwide Inc
0409-3177-01	Lidocaine 0.5%+Epi 1:200,000 MDV 25 x 50 ml	25/bx	Hospira Worldwide Inc
52584-276-01	Lidocaine 1% MDV 20 ml	25/bx	Hospira Worldwide Inc
0409-4276-01	Lidocaine 1% MDV 25 x 20 ml	25/bx	Hospira Worldwide Inc
0409-4276-02	Lidocaine 1% MDV 25 x 50 ml	25/bx	Hospira Worldwide Inc
0409-4713-02	Lidocaine 1% PF Amp 25 x 5 ml	25/bx	Hospira Worldwide Inc
0409-4713-32	Lidocaine 1% PF Amp 2 ml	50/bx	Pfizer Inc/Consumr Hlthcr
0409-4713-32	Lidocaine 1% PF Amp 50 x 2 ml	50/bx	Hospira Worldwide Inc
0409-9137-05	Lidocaine 1% PF PFs 10 x 5 ml	10/bx	Hospira Worldwide Inc
0409-4279-02	Lidocaine 1% PF SDV 25 x 30 ml	25/bx	Hospira Worldwide Inc
0409-3178-01	Lidocaine 1%+Epi 1:100,000 MDV 25 x 20 ml	25/bx	Hospira Worldwide Inc
0409-3178-03	Lidocaine 1%+Epi 1:100,000 MDV 25 x 50 ml	25/bx	Hospira Worldwide Inc
52584-178-02	Lidocaine 1%+Epi 1:100,000 MDV 30 ml	1/ea	Hospira Worldwide Inc
52584-178-03	Lidocaine 1%+Epi 1:100,000 MDV 50 ml	1/ea	Hospira Worldwide Inc
0409-3178-02	Lidocaine 1%+Epi 100,000 MDV 25 x 30 ml	25/bx	Hospira Worldwide Inc

Injectables Continued

Item No.	Description	Qty	Manufacturer
0409-4277-02	Lidocaine 2% FTV 25 x 50 ml	25/bx	Hospira Worldwide Inc
52584-277-02	Lidocaine 2% FTV 50 ml	1/ea	Hospira Worldwide Inc
52584-277-01	Lidocaine 2% MDV 20 ml	1/ea	Hospira Worldwide Inc
0409-4277-01	Lidocaine 2% MDV 25 x 20 ml	25/bx	Hospira Worldwide Inc
0409-4282-02	Lidocaine 2% PF Amp 25 x 10 ml	25/bx	Hospira Worldwide Inc
0409-4282-01	Lidocaine 2% PF Amp 25 x 2 ml	25/bx	Hospira Worldwide Inc
0409-4903-34	Lidocaine 2% PF PFs 10 x 5 ml	10/bx	Hospira Worldwide Inc
0409-1323-05	Lidocaine 2% PF PFs 10 x 5 ml	10/bx	Hospira Worldwide Inc
0409-2066-05	Lidocaine 2% PF SDV 10 x 5 ml	10/bx	Hospira Worldwide Inc
0409-3182-01	Lidocaine 2%+Epi 1:100,000 FTV 25 x 20 ml	25/bx	Hospira Worldwide Inc
0409-3182-02	Lidocaine 2%+Epi 1:100,000 FTV 25 x 30 ml	25/bx	Hospira Worldwide Inc
0409-3182-03	Lidocaine 2%+Epi 1:100,000 FTV 25 x 50 ml	25/bx	Hospira Worldwide Inc
52584-182-03	Lidocaine 2%+Epi 1:100,000 FTV 50 ml	1/ea	Hospira Worldwide Inc
0409-3183-01	Lidocaine 2%+Epi 1:200,000 SDV 5 x 20	5/bx	Hospira Worldwide Inc
0409-4283-01	Lidocaine 4% PF Amp 25 x 5 ml	25/bx	Hospira Worldwide Inc
0009-0555-02	Lincocin, 300 mg/ml, Vial, 10 ml	1/ea	Pfizer Inc/Consumr Hlthcr
67457-134-50	Magnesium Chl 200 mg/ml Vial 50 ml	1/ea	Mylan
0409-6729-09	Magnesium Sul 40 mg/ml FC 12 x 1000 ml	12/cs	Hospira Worldwide Inc
63323-064-03	Magnesium Sul 50% SDV 25 x 2 ml	25/bx	Fresenius Kabi
0409-1587-50	Marcaine 0.25% FTV 1 x 50 ml	1/ea	Hospira Worldwide Inc
0409-1559-10	Marcaine 0.25% SDV 10 x 10 ml	10/bx	Hospira Worldwide Inc
0409-1559-30	Marcaine 0.25% SDV 10 x 30 ml	10/bx	Hospira Worldwide Inc
0409-1752-50	Marcaine 0.25%+Epi 1:200,000 MDV 1 x 50 ml	1/ea	Hospira Worldwide Inc
0409-1746-10	Marcaine 0.25%+Epi 1:200,000 SDV 10 x 10 ml	10/bx	Hospira Worldwide Inc
0409-1746-30	Marcaine 0.25%+Epi 1:200,000 SDV 10 x 30 ml	10/bx	Hospira Worldwide Inc
0409-1610-50	Marcaine 0.5% FTV 1 x 50 ml	1/ea	Hospira Worldwide Inc
0409-1560-10	Marcaine 0.5% SDV 10 x 10 ml	10/bx	Hospira Worldwide Inc
0409-1560-29	Marcaine 0.5% SDV 10 x 30 ml	10/bx	Hospira Worldwide Inc
0409-1755-50	Marcaine 0.5%+Epi 1:200,000 FTV 1 x 50 ml	1/ea	Hospira Worldwide Inc
0409-1749-10	Marcaine 0.5%+Epi 1:200,000 SDV 10 x 10 ml	10/bx	Hospira Worldwide Inc
0409-1749-29	Marcaine 0.5%+Epi1:200000 Sdv 30 ml	10/bx	Pfizer Inc/Consumr Hlthcr
0409-1582-29	Marcaine 0.75% SDV 10 x 30 ml	10/bx	Hospira Worldwide Inc
0409-1582-10	Marcaine 0.75% SDV, 10 x 10 ml	10/bx	Hospira Worldwide Inc
0409-1761-02	Marcaine Spinal 0.75% Amp 10 x 2 ml	10/bx	Hospira Worldwide Inc
0548-5701-00	Medroxyprogesterone 150 mg/ml PFS 1 ml	1/ea	Amphastar Pharmaceuticals
59762-4537-1	Medroxyprogesterone 150 mg/ml Vial 1 ml	1/ea	Greenstone LLC
0548-5400-00	Medroxyprogesterone 150 mg/ml SDV 1 ml	1/ea	Amphastar Pharmaceuticals
67457-273-10	Methocarbamol 100 mg/ml SDV 25 x 10 ml	25/bx	Mylan
17478-504-10	Methylene Blue 1% Sdv 10 ml	10/bx	Akorn Inc
0409-3414-01	Metoclopramide 5 mg/ml SDV 25 x 2 ml	25/bx	Hospira Worldwide Inc
63323-614-01	MBO-Hydralazine 20 mg Vial 25 x 1 ml	25/bx	Fresenius Kabi
63323-064-11	MBO-Magnesium Sul 50% SDV 25 x 10 ml	25/bx	Fresenius Kabi
0641-1495-35	Mbo-Promethazine 25 mg/ml Amp 25 x 1 ml	25/bx	Baxter Healthcare
0409-1778-35	Metoprolol Tartrate 1 mg/ml PFs 3 x 5 ml	3/bx	Hospira Worldwide Inc
0409-1778-05	Metoprolol Tartrate 1 mg/ml SDV 10 x 5 ml	10/bx	Hospira Worldwide Inc
0409-2285-05	Metoprolol Tartrate 1 mg/ml Amp 12 x 5 ml	12/bx	Hospira Worldwide Inc
16729-108-11	Mitomycin 20 mg Vial	1/ea	Accord Healthcare Inc
16729-116-38	Mitomycin 40 mg Vial	1/ea	Accord Healthcare Inc
63323-330-60	NaFCillin 10 gm Vial	1/ea	APP Pharm
0409-1464-01	Nalbuphine HCL 10 mg/ml FTV 25 x 10 ml	25/cs	Hospira Worldwide Inc
0409-1467-01	Nalbuphine HCL 20 mg/ml FTV 25 x 10 ml	25/cs	Hospira Worldwide Inc

Item No.	Description	Qty	Manufacturer
76329-3369-1	Naloxone 2 mg/2ml LI Syg 10 x 2 ml	10/bx	International Medication Systems
0409-1219-01	Naloxone HCL 0.4 mg/ml MDV 25 x 10 ml	25/cs	Hospira Worldwide Inc
0409-1782-69	Naloxone HCL 0.4 mg/ml PFs 10 x 1 ml	10/bx	Hospira Worldwide Inc
0409-1215-01	Naloxone HCL 0.4 mg/ml SDV 10 x 1 ml	10/bx	Hospira Worldwide Inc
63323-477-27	Nesacaine-MPF 2% SDV 25 x 20 ml	25/bx	APP Pharm
0517-4810-25	Nitroglycerin 5 mg/ml Vial 25 x 10 ml	25/bx	American Regent Inc
42998-302-02	Nitropress 25 mg/ml SDV 2 ml	1/ea	Marathon Pharmaceuticals LLC
55566-1501-1	Novarel 10 M Un + Dil Mdv 10 ml	1/kt	Ferring Pharmaceuticals Inc
0169-1837-11	Novolin 70/30 100U/ml Vial 10 ml	1/ea	Novo Nordisk Inc
0169-1834-11	Novolin N 100U/ml Vial 10 ml	1/ea	Novo Nordisk Inc
0169-1833-11	Novolin R 100U/ml Vial 10 ml	1/ea	Novo Nordisk Inc
43825-102-01	Ofirmev 10 mg/ml Inj 24 x 100 ml	24/bx	Cadence Pharmaceuticals
0143-9890-01	Ondansetron 2 mg/ml 20ml MDV	1/ea	West-Ward Pharm Corp
0409-4759-01	Ondansetron 2 mg/ml FTV 1 x 20 ml	1/ea	Hospira Worldwide Inc
0409-4755-03	Ondansetron 2 mg/ml SDV 25 x 2 ml	25/bx	Hospira Worldwide Inc
63323-373-02	Ondansetron 2 mg/ml PF SDV 25 x 2 ml	25/bx	Fresenius Kabi
0591-3222-47	Orphenad Cit 30 mg/ml Amp 10 x 2 ml	10/bx	Watson Pharmaceuticals
17478-538-02	Orphenad Cit, 30 mg/ml Vial, 10 x 2 ml	10/bx	Akorn Inc
59676-360-01	Orthovisc 15 mg/ml Syg 2 ml	1/ea	Anika Therapeutics, Inc.
63323-012-10	Oxytocin 10 Un MDV 25 x 10 ml	25/bx	APP Pharm
63323-012-01	Oxytocin 10 Un Vial 25 x 1 ml	25/bx	APP Pharm
0409-4646-01	Pancuronium Bromide Inj 1 mg/ml 10 ml	25/bx	Hospira Worldwide Inc
0517-4010-01	Papaverine 30 mg/ml Vial 10 ml	1/ea	American Regent Inc
0781-6135-95	Penicillin G Potassium 5 Mmu Vial	10/bx	Sandoz
60793-131-10	Penicillin G 1.2M Un PFs 10 x 2 ml	10/bx	King Pharmaceuticals
60793-130-10	Penicillin G 600M Un Tbx Syg 10 x 1 ml	10/bx	Pfizer Inc/Consumr Hlthcr
0641-6082-25	Phenergan 25 mg Amp 25 x 1 ml	25/bx	Baxter Healthcare
0641-6142-01	Phenylephrine HCL 10 mg/ml Sdv 1 ml	1/ea	West-Ward Pharmaceuticals
67457-522-37	Piperacillin/Tazobactam 3.375G Sdv	10/bx	Mylan
0409-6651-06	Potassium Cl 20 Meq SDV 25 x 10 ml	25/bx	Hospira Worldwide Inc
0409-6653-05	Potassium Cl 40 Meq SDV 25 x 20 ml	25/bx	Hospira Worldwide Inc
0409-1902-01	Procainamide HCL 100 mg/ml FTV 25 x 10 ml	25/cs	Hospira Worldwide Inc
0409-1903-01	Procainamide HCL 500 mg/ml FTV 25 x 2 ml	25/bx	Hospira Worldwide Inc
23155-294-42	Prochlorperazine 5 mg/ml Vial 2 ml	10/bx	Heritage Labs
0591-3128-79	Progesterone, 50 mg Oil Vial, 10 ml	1/ea	Watson Pharmaceuticals
0641-0928-25	Promethazine 25 mg/ml Vial 25 x 1 ml	25/bx	West-Ward Pharm Corp
0641-1496-35	Promethazine HCL, 50 mg/ml, 1 ml Amp	25/bx	Baxter Healthcare
63323-604-01	Propranolol HCL 1 mg/ml SDV 10 x 1 ml	10/bx	Fresenius Kabi
63323-229-05	Protamine Sul 50 mg Vial 25 x 5 ml	25/bx	APP Pharm
63323-180-01	Pyridoxine 100 mg/ml Vial 25 x 1 ml	25/bx	APP Pharm
0409-6629-02	Quelicin/Succinyl 20 mg/ml Mdv 25 x 10 ml	25/bx	Pfizer Inc/Consumr Hlthcr
68382-423-06	Ranitidine 25 mg/ml Mdv 6 ml	1/ea	Zydus Pharmaceuticals USA Inc
0078-0435-61	Reclast 5 mg Sol 100 ml	1/ea	Novartis Consumer Health
57894-030-01	Remicade 100 mg Vial	1/ea	Janssen Pharmaceuticals,
44206-300-01	Rhophylac 300 mcg	1/ea	Castle
44206-300-10	Rhophylac 300 mcg	10/bx	CSL Behring
0078-0180-01	Sandostatin 50 mcg/ml Amp 10 x 1 ml	10/bx	Novartis Consumer Health
63323-467-57	Sensorcaine 0.5 % Vial 25 x 50 ml	25/bx	APP Pharm
63323-466-17	Sensorcaine-MPF 0.5% SDV 25 x 10 ml	25/bx	APP Pharm
63323-187-30	Sod Chl Conc 23.4 % Vial 25 x 30 ml	25/bx	APP Pharm
0409-4916-34	Sodium Bicarbonate 7.5% PFs 10 x 50 ml	10/bx	Hospira Worldwide Inc

Injectables Continued

Item No.	Description	Qty	Manufacturer
0409-6637-34	Sodium Bicarbonate 8.4% PFs 10 x 50 ml	10/bx	Hospira Worldwide Inc
52584-637-34	Sodium Bicarbonate 8.4% PFs 50 ml	1/ea	Hospira Worldwide Inc
0409-6625-02	Sodium Bicarbonate 8.4% SDV 25 x 50 ml	25/bx	Hospira Worldwide Inc
0409-5534-34	Sodium Bicarbonate Inft 4.2% PFs 10 x 10 ml	10/bx	Hospira Worldwide Inc
0409-4900-34	Sodium Bicarbonate Ped 8.4% PFs 10 x 10 ml	10/bx	Hospira Worldwide Inc
63323-088-63	Sodium Chloride 23.4% Maxivial 200 ml	20/cs	APP Pharm
0009-0825-01	Solu Cortef 100 mg Vial	1/ea	Pfizer Inc/Consumr Hlthcr
0009-0013-05	Solu Cortef 250 mg Aov 2 ml	1/ea	Pfizer Inc/Consumr Hlthcr
0009-0013-06	Solu Cortef 250 mg Aov 25 x 2 ml	25/bx	Pfizer Inc/Consumr Hlthcr
0009-0018-20	Solu-Medrol 1 gm Aov 8 ml	1/ea	Pfizer Inc/Consumr Hlthcr
52584-047-25	Solu-Medrol 125 mg Aov 2 ml	1/ea	Pfizer Inc/Consumr Hlthcr
0009-0047-22	Solu-Medrol 125 mg Aov 2 ml 25 x 2 ml	25/bx	Pfizer Inc/Consumr Hlthcr
0009-0039-28	Solu-Medrol 40 mg/ml Aov	25/bx	Pfizer Inc/Consumr Hlthcr
52584-039-30	Solu-Medrol 40 mg/ml Aov	1/ea	Pfizer Inc/Consumr Hlthcr
0009-0003-02	Solu-Medrol 500 mg PF Aov SDV 4 ml	1/ea	Pfizer Inc/Consumr Hlthcr
63323-185-10	Sterile Water For Inj SDV 10 ml	25/bx	Fresenius Kabi
58468-090-01	Synvisc Syg 3 x 2 ml	3/bx	Teleflex Medical
58468-090-03	Synvisc-One 10 ml Syg 6 ml	1/ea	Teleflex Medical
63323-013-02	Thiamine 100 mg/ml Vial 25 x 2 ml	25/bx	APP Pharm
60793-215-05	Thrombin Jmi 5Mun /Dil Vial	1/ea	Pfizer Inc/Consumr Hlthcr
42023-118-01	Tigan, 100 mg/ml Vial, 20 ml	1/ea	JHP Pharmaceuticals
0409-3578-01	Tobramycin 40 mg/ml FTV 25 x 2 ml	25/cs	Hospira Worldwide Inc
0338-0644-04	Travasol 10 % Sol 12 x 1000 ml	12/bx	Baxter Healthcare
52544-092-76	Trelstar 22.5 mg Inj	1/ea	Watson Pharmaceuticals
52544-188-76	Trelstar La 11.2 5 mg Inj	1/ea	Watson Pharmaceuticals
0008-4990-20	Tygacil 50 mg SDV 10 x 10 ml	10/bx	Wyeth Consumer
0143-9785-10	Valproate 500 mg/5 ml Vial 10 x 5 ml	10/bx	West-Ward Pharm Corp
63323-494-05	Valproate Sod 500 mg/5 ml Vial 10 x 5 ml	10/bx	APP Pharm
0409-6533-01	Vancomycin HCL 1G SDV 10 x 1 g	10/bx	Hospira Worldwide Inc
0409-4332-01	Vancomycin HCL 500 mg SDV 10 x 500 mg	10/bx	Hospira Worldwide Inc
42023-164-25	Vasostrict 20U/ml Mdv 1 ml	25/bx	Par Pharmaceutical Companies, Inc
0517-2340-10	Venofer 20 mg/ml SDV 10 x 5 ml	10/bx	American Regent Inc
0409-4011-01	Verapamil HCL 2.5 mg/ml Amp 5 x 2 ml	5/bx	Hospira Worldwide Inc
0409-1144-05	Verapamil HCL 2.5 mg/ml SDV 25 x 2 ml	25/bx	Hospira Worldwide Inc
0409-1144-02	Verapamil HCL 2.5 mg/ml SDV 5 x 4 ml	5/bx	Hospira Worldwide Inc
MEN6065	Vial, Astringyn 8 g	12/bx	Cooper Surgical
61703-309-06	Vincristine Sulf 1 mg/ml PF SDV 240 x 1 ml	1/ea	Hospira Worldwide Inc
67457-146-30	Vit B-Complex, 100 Inj, 30 ml	1/ea	Bioniche Life Sciences Inc
0409-9157-01	Vitamin K1 1 mg/0.5 ml Amp 25 x 0.5 ml	25/bx	Hospira Worldwide Inc
0409-9158-01	Vitamin K1 10 mg/ml Amp 25 x 1 ml	25/bx	Hospira Worldwide Inc
24208-002-02	Vitrase Ov 200U/ml 1.2 ml(<12 Mo exp)	2/bx	Bausch & Lomb
63323-484-57	Xylocaine 0.5 % MDV 25 x 50 ml	25/bx	APP Pharm
52584-485-57	Xylocaine 1 % MDV 50 ml	1/ea	APP Pharm
63323-492-27	Xylocaine 1 % Vial 25 x 2 ml	25/bx	APP Pharm
63323-485-57	Xylocaine 1 % Vial 25 x 50 ml	25/bx	APP Pharm
63323-486-27	Xylocaine 2 % Vial 25 x 20 ml	25/bx	APP Pharm
63323-486-57	Xylocaine 2 % Vial 25 x 50 ml	25/bx	APP Pharm
52584-489-27	Xylocaine-MPF 2% & Epi 1:200000 SDV 20 ml	1/ea	APP Pharm

Miscellaneous

Item No.	Description	Qty	Manufacturer
OTC073230	Acetaminophen 120 mg Suppository 12/bx	1/ea	Perrigo
OTC211402	Acetaminophen 80 mg Infant Supp	6/bx	Taro Pharmaceutical Industries Ltd
0132-0202-20	Fleet Enema Child Enm 2.25-oz.	1/ea	CB Fleet
OTC211602	MBO Acetaminophen 325 mg Suppository 6/bx	1/ea	Taro Pharmaceutical Industries Ltd
OTC010015	Neo-Synephrine 1% X/S Nasal Spray 15 ml	1/ea	Foundation Consumer Healthcare LLC
0003-0526-20	Nitratetest Dis 1-Roll	1/ea	Bristol-Myers Squibb
0884-6297-30	Phenol Ez Swb 30	1/ea	Pedinol Pharmacal Inc
59365-5030-0	Trimo-San W/Apppl Jel 4-oz.	1/ea	Cooper Surgical

Ophthalmics

Item No.	Description	Qty	Manufacturer
17478-716-10	Apraclonidine 0.5% Oph Soln 5 ml	1/ea	Akorn Inc
17478-215-05	Atropine Sulfate 1% Oph Soln 5 ml	1/ea	Akorn Inc
0574-4022-35	Bacitracin 500U/gm Oph Oint 3.5 gm	1/ea	Perrigo
BLP39182	Bio Glo Fluorescein Sodium Ophthalmic Strips	100/bx	Bausch & Lomb
0065-0795-50	Bss Opht Glass Sol 6 x 500 ml	6/bx	Alcon Labs
0065-0795-15	Bss Opth 36 x 15 ml	36/bx	Alcon Labs
61314-656-05	Ciprofloxacin 0.3% Opth Drp 5 ml	1/ea	Sandoz
OTC000252	Eye Wash Soln 4-oz.	1/ea	Bausch & Lomb
0065-0359-05	Cyclomydril Drp 5 ml	1/ea	Alcon Labs
17478-100-12	Cyclopentolate 1% Oph Soln 15 ml	1/ea	Akorn Inc
17478-100-02	Cyclopentolate 1% Oph Soln 2 ml	1/ea	Akorn Inc
17478-070-35	Erythromycin 0.5% Oph Oint 3.5 gm	1/ea	Akorn Inc
60758-880-15	Fluorometh 0.1% Drp 15 ml	1/ea	Pacific Pharma, Inc.
17478-284-35	Gentak 0.3% O/O 3.5 gm	1/ea	Akorn Inc
60758-188-05	Gentamicin 0.3% Drp 5 ml	1/ea	Pacific Pharma, Inc.
24208-580-60	Gentamicin 0.3% Ophthalmic Soln 5 ml	1/ea	Bausch & Lomb
17478-064-12	Gonak 2.5% Sol 15 ml	1/ea	Akorn Inc
OTC071710	Ketotifen 0.025% Oph Soln 5 ml	1/ea	Akorn Inc
24208-780-55	Neo/Poly/Bac O/O, 3.5 gm	1/ea	Bausch & Lomb
61314-630-06	Neo/Poly/Dex Drp 5 ml	1/ea	Sandoz
61314-631-36	Neo/Poly/Dex O/O 3.5 gm	1/ea	Sandoz
24208-631-10	Neomycin/Polymyxin B Hc Oph Sol 10 ml	1/ea	Bausch & Lomb
17478-713-10	Ofloxacin 0.3% O/S 5 ml	1/ea	Akorn Inc
42702-102-15	Phenylephrine Hcl 2.5% Oph Soln 15 ml	1/ea	Paragon Biotech, Inc.
17478-223-1	Pilocarpine 1% Oph Soln 15 ml	1/ea	Akorn Inc
61314-204-15	Pilocarpine 2% Drp 15 ml	1/ea	Sandoz
17478-703-11	Polymyxin B & Trimethoprim Oph Sol 10 ml	1/ea	Akorn Inc
0023-0106-05	Pred-G 0.3% /1% Oph Susp 5 ml	1/ea	Allergan, Inc.
61314-637-05	Prednisolone Acetate 1% Oph Susp 5 ml	1/ea	Falcon Safety
60758-119-05	Prednisolone Acetate 1% Oph Susp 5 ml	1/ea	Pacific Pharma, Inc.
17478-263-12	Proparacaine 0.5% Drp 15 ml	1/ea	Akorn Inc
0023-4554-30	Refresh Cellu Visc Sc O/S 30	30/bx	Allergan, Inc.
0023-3240-15	Refresh Optiv Ophthalmic Drop 15 ml	1/ea	Allergan, Inc.
0023-9163-30	Restasis 0.05% VI 30 x 0.4 ml	30/bx	Allergan, Inc.
61314-701-01	Sod Sulfacet 10% Drp 15 ml	1/ea	Sandoz
17478-404-01	Strip, Ful-Glo, 1 mg Stp, 100	1/ea	Akorn Inc

Ophthalmics Continued

Item No.	Description	Qty	Manufacturer
17478-045-32	Tetracaine 1% Amp 25 x 2 ml	25/bx	Akorn Inc
0065-0741-14	Tetracaine HCL 0.5% Oph Soln 12 x 4 ml UD	1/ea	Alcon Labs
61314-647-25	Tobramyc/Dexa 0.3/ 0.1% Sus 2.5 ml	1/ea	Sandoz
61314-643-05	Tobramycin 0.3% Drp 5 ml	1/ea	Sandoz
61314-647-05	Tobramycin/Dexamethasone Oph Susp 5 ml	1/bt	Falcon Safety
17478-101-12	Tropicacyl, 0.5% O/S, 15 ml	1/ea	Akorn Inc
17478-102-12	Tropicamide 1% Oph Soln 15 ml	1/ea	Akorn Inc

Oral Tablets, Capsules & Liquids

Item No.	Description	Qty	Manufacturer
OTC664524CS	Acetamin Jr 160 mg Chew Tab 24/bx	48/cs	Major Pharmaceutical
0527-1050-01	Acetazolamide 250 mg Tab	100/bt	Lannett Company
68382-261-01	Acetazolamide 500 mg Er Cap	100/bt	Zyodus Pharmaceuticals
0536-7180-01	Acidophilus Pectin Cpl	100/bt	Watson Pharmaceuticals
00591-5543-01	Allopurinol 100 mg Tab	100/bt	Watson Pharmaceuticals
0603-2116-21	Allopurinol 300 mg Tab	100/bt	Qualitest Pharmaceuticals
0591-5544-01	Allopurinol 300 mg Tab	100/bt	Watson Pharmaceuticals
67877-198-90	Amlodipine Besylate 5 mg Tab 90/bt	1/ea	Ascend Laboratories
0781-2020-01	Amoxicillin 250 mg Cap 100	100/bt	Sandoz
0093-2268-01	Amoxicillin 250 mg Chewable Tab	100/bt	Teva Pharmaceuticals USA
0781-6041-46	Amoxicillin 250 mg/5ml Sus 100 ml	1/ea	Sandoz
0143-9887-01	Amoxicillin 400 mg/5ml Oral Susp 100 ml	1/ea	West-Ward Pharm Corp
0781-2613-01	Amoxicillin 500 mg Cap	100/bt	Sandoz
0781-5061-01	Amoxicillin 875 mg Tab	100/bt	Sandoz
0781-2613-05	Amoxicillin, 500 mg Cap	500/bt	Sandoz
0093-0752-01	Atenolol 50 mg Tab	100/bt	Teva Pharmaceuticals USA
68180-160-06	Azithromycin 250 mg Tab	30/bt	Lupin Pharmaceuticals Inc
68180-160-13	Azithromycin 250 mg Tab Uou 3 x 6 Tabs	18/bt	Lupin Pharmaceuticals Inc
59762-3051-2	Azithromycin, 1 gm Sdp Pkt	3/bx	Greenstone Llc
65597-103-90	Benicar 20 mg Tab	90/bt	Daiichi Sankyo
0054-3117-63	Calcium Carb 1250 mg/5 ml Sus 500 ml	1/ea	Roxane Laboratories Inc
0869-2452-10	Castor Oil 2-oz.	1/ea	Cumberland-Swan Inc
0025-1520-31	Celebrex 100 mg Cap	100/bt	Pfizer Inc/Consumr Hlthcr
0093-3145-01	Cephalexin 250 mg Cap	100/bt	Teva Pharmaceuticals USA
68180-122-01	Cephalexin 500 mg Cap	100/bt	Lupin Pharmaceuticals Inc
68180-122-02	Cephalexin,500 mg Cap	500/bt	Lupin Pharmaceuticals Inc
OTC27716	Cetirizine, 1 mg/ml, 4-oz. Syrup, Perrigo	1/bt	Perrigo
65862-077-01	Ciprofloxacin 500 mg Tab	100/bt	Aurobindo Pharma
0591-5708-01	Clindamycin 150 mg Cap	100/bt	Watson Pharmaceuticals
0603-2959-21	Clonidine 0.3 mg Tab	100/bt	Qualitest Pharmaceuticals
52584-656-61	Clonidine Hcl 0.1 mg UD Tab (1 Tablet)	1/ea	Major Pharmaceutical
68001-238-00	Clonidine Hcl 0.2 mg Tab	100/bt	Bluepoint Laboratories
68001-237-00	Clonidine Hcl 0.1 mg Tab 100/bt	10/bx	Bluepoint Laboratories
64764-119-01	Colcrys 0.6 mg	100/bt	Takeda Pharmaceuticals
10702-007-01	Cyclobenzaprine Hcl 10 mg Tab	100/bt	Kvk-Tech, Inc.
0054-3177-63	Dexameth Oral 0.5 mg/5 ml Sol 500 ml	1/ea	Roxane Laboratories Inc

Item No.	Description	Qty	Manufacturer
0054-4184-25	Dexamethasone 4 mg Tab	100/bt	West-Ward Pharmaceuticals
0093-3125-01	Dicloxacillin 500 mg Cap	100/bt	Teva Pharmaceuticals USA
0378-1610-01	Dicyclomine 10 mg Cap	100/bt	Mylan
0078-0359-34	Diovan 160 mg Tab	90/bt	Novartis Consumer Health
68382-187-01	Dipyridamole 25 mg Tab	100/bt	Zydus Pharmaceuticals
64980-171-01	Disulfiram 250 mg Tab	100/bt	Rising Pharmaceuticals,
0409-2344-02	Dobutamine 12.5 mg/ml Sdv 20 ml	10/bx	Pfizer Inc/Consumr Hlthcr
49884-093-03	Doxycycl Mono 100 mg Tab	50/bt	Par Pharmaceutical
0143-3142-50	Doxycycline 100 mg Cap	50/bt	West-Ward Pharm Corp
0143-3142-05	Doxycycline Hyclate 100 mg Cap	500/bt	West-Ward Pharmaceuticals
0143-2112-50	Doxycycline, 100 mg Tab	50/bt	West-Ward Pharm Corp
24338-126-13	Ery-Tab Dr, 500 mg, 100 Tab	100/bt	Arbor Pharma
0378-1454-01	Estradiol 1 mg Tab	100/bt	Mylan
51079-966-20	Famotidine 20 mg Ud Tab (Rx)	100/bx	Mylan
45802-425-78	Fexofenadine 60 mg Tab	100/bt	Perrigo
62584-897-01	Folic Acid 1 mg Tab	100/bx	American Health Prod
51079-105-20	Folic Acid 1 mg Tab	100/bx	Udl Laboratories Inc
65162-361-10	Folic Acid 1 mg Tab 100/bt	1/ea	Amneal Pharmaceuticals
0054-4299-25	Furosemide 40 mg Tab	100/bt	Roxane Laboratories Inc
16729-183-01	Hydrochlorothiazide 25 mg Tab	100/bt	Accord Healthcare Inc
0185-0674-01	Hydroxyzine Pamoate 25 mg Capsule	100/bt	Sandoz
0006-0717-54	Hyzaar 50 mg/12.5 mg Tab	90/bt	Merck & Co Inc
0904-7915-40	Ibuprofen 200 mgbrn Tab	500/bt	Major Pharmaceutical
67877-319-05	Ibuprofen 400 mg Tab 500/bt	1/bt	Ascend Laboratories
0904-5854-40	Ibuprofen 600 mg Tab	500/bt	Major Pharmaceutical
67877-321-01	Ibuprofen 800 mg Tab 100/bt	1/bt	Ascend Laboratories
0904-5309-09	Ibuprofen Chld Berry Sus 8-oz.	1/ea	Major Pharmaceutical
50383-584-04	Ibuprofen 100 mg/5 ml Oral Susp 4-oz. (Rx)	1/ea	Akorn Inc
0093-4029-01	Indomethacin 25 mg Cap	100/bt	Teva Pharmaceuticals USA
0187-0746-33	Insta Glucose Tri-Pack Gel 3 x 31 gm	3/bx	Valeant Pharmaceuticals
61748-013-10	Isoniazid 300 mg Tab	1000/bt	Versapharm, Inc
52544-892-28	Jolivette-28 0.35 mg Tab 6 x 28	6/bx	Watson Pharmaceuticals
0555-9027-42	Junel 21 1.5 M /30 mcg Tab 3 x 21	3/bx	Teva Pharmaceuticals USA
51672-4026-6	Ketoconazole, 200 mg Tab	30/bt	Taro Pharmaceutical
66220-719-30	Kristalose Pkt 30 x 10 Gm	1/bx	Cumberland
0527-1343-01	Levothyroxine, 75 mcg Tab	100/bt	Lannett Company
0185-0074-01	Lovastatin 40 mg Tab	100/bt	Sandoz
68382-050-01	Meloxicam 7.5 mg Tab	100/bt	Zydus Pharmaceuticals
31722-534-01	Methocarbamol 750 mg Tab	100/bt	Camber Pharmaceuticals
0555-0572-02	Methotrexate 2.5 mg Tab	100/bt	Teva Pharmaceuticals USA
57664-477-52	Metoprolol Tartrate 50 mg Tab	100/bt	Sun Pharmaceutical Industries
50111-333-01	Metronidazole 250 mg Tab	100/bt	Teva Pharmaceuticals USA
50111-334-01	Metronidazole 500 mg Tab	100/bt	Teva Pharmaceuticals USA
12843-158-59	Midol Complete Caplet	24/bt	Bayer Healthcare
16477-505-01	Millipred 5 mg Tab	100/bt	Laser Pharmaceuticals
59762-5008-1	Misoprostol 200 mcg Tab	60/bt	Greenstone LLC
0004-6316-01	Naprosyn 500 mg Tab	100/bt	Genentech, Inc
53746-190-10	Naproxen 500 mg Tab	1000/bt	Amneal Pharmaceuticals
68462-190-01	Naproxen, 500 mg Tab	100/bt	Glenmark Generics Limited

Oral Tablets, Capsules & Liquids Continued

Item No.	Description	Qty	Manufacturer
0536-4078-10	Niacin 500 mg Tab	1000/bt	Watson Pharmaceuticals
0228-2530-10	Nifedipine 20 mg Cap	100/bt	Amerisource
0378-0353-01	Nifedipine Er 30 mg Cap	100/bt	Amerisource
51079-348-20	Nitrofurantoin Mono/Hyd, 100 mg, Tablet	100/bx	Udl Laboratories Inc
0185-0122-01	Nitrofurantoin Monohydr 100 mg Cap	100/bt	Sandoz
70003-300-20	Nitroglycerin Lingual 0.4 mg Spray 12 Gm	1/ea	Jacksonville Pharmaceuticals, Inc.
70003-300-65	Nitroglycerin Lingual 0.4 mg Spray 4.9 Gm	1/ea	Jacksonville Pharmaceuticals, Inc.
52584-418-13	Nitrostat Tabs, 0.4 mg	25/bt	Pfizer Inc/Consumr Hlthcr
0071-0418-13	Nitrostat Tabs, 0.4 mg, 4 x 25/Bottle	4/bx	Pfizer Inc/Consumr Hlthcr
0555-9131-67	Ocella 3/ 0.03 mg Tab 3 x 28	3/bx	Teva Pharmaceuticals USA
68462-396-10	Omeprazole 20 mg Dr Cap 1,000/bt (Rx)	1/ea	Glenmark Generics Limited
62175-118-32	Omeprazole 20 mg Dr Cap (Rx)	30/bt	Kremers Urban
0378-0315-93	Ondansetron 4 mg Tab	30/bt	Mylan
68001-246-17	Ondansetron 4 mg Odt Tab	30/bx	Bluepoint Laboratories
45963-539-30	Ondansetron 8 mg Tab	30/bt	Actavis Us
50458-197-15	Ortho-Cyclen 28 Day Tab 6 x 28	6/bx	Janssen Pharmaceuticals, Inc.
31722-713-90	Pantoprazole 40 mg Dr Tab	90/bt	Camber Pharmaceuticals
65862-175-01	Penicillin Vk 250 mg Tab	100/bt	Aurobindo Pharma
0781-5710-01	Potassium Chl 10 Meq Tab	100/bt	Sandoz
0121-0759-08	Prednisolone Sod Phos 15/5 ml Soln 237 ml	1/ea	Pharmaceutical Associates Inc
50383-042-24	Prednisolone Sol, 8-oz.	1/ea	Hi-Tech Pharmcal Co. Inc
0054-4741-25	Prednisone 1 mg Tab	100/bt	Roxane Laboratories Inc
0603-5338-21	Prednisone 10 mg Tab	100/bt	Qualitest Pharmaceuticals
0054-0017-25	Prednisone 10 mg Tab	100/bt	West-Ward Pharmaceuticals
0054-4742-25	Prednisone 2.5 mg Tab	100/bt	Roxane Laboratories Inc
0143-9738-01	Prednisone 20 mg Tab	100/bt	West-Ward Pharmaceuticals
0603-5337-21	Prednisone 5 mg Tab	100/bt	Par Pharmaceutical Companies, Inc
59746-115-06	Prochlorperazine 10 mg Tablet	100/bt	Jubilant Cadista Pharmaceuticals
68382-041-01	Promethazine 25 mg Tab	100/bt	Zydus Pharmaceuticals USA Inc
68094-204-62	Ranitidine 150 mg Syr 30 x 10 ml Ud	30/bx	American Health Packaging
64380-803-07	Ranitidine 150 mg Tab	500/bt	Strides Pharma, Inc
61748-018-01	Rifampin 300 mg Cap	100/ea	Versapharm, Inc
OTC589430	Simethicone 20 mg/0.3 ml Oral Susp 30 ml	1/ea	Major Pharmaceutical
0223-1721-01	Sodium Bicarbonate 10 Gr Tab	100/bt	Consolidated Midland Corporation
53746-272-01	Sulfameth/Trimeth 800/160 mg Tab	100/bt	Amneal Pharmaceuticals
63304-098-19	Sumatriptan 50 mg Tab 9	1/ea	Ranbaxy Laboratories Limited
0004-0801-85	Tamiflu 45 mg Blisterpak	10/bx	Genentech, Inc
0004-0800-85	Tamiflu 75 mg 10 Uou Cap	10/bx	Genentech, Inc
24208-920-64	Tetracaine 0.5% Oph Solution 15 ml	1/ea	Bausch & Lomb
0781-1123-01	Triamter/Hctz 37.5 /25 mg Tab	100/bt	Sandoz
0781-1123-05	Triamter/Hctz 37.5 /25 mg Tab	500/bt	Sandoz
61958-701-01	Truvada Tab	30/bt	Gilead Sciences, Inc.
68084-485-01	Venlafaxine 75 mg Er Cap Ud	100/bt	American Health Packaging
0536-4406-01	Vit B-6 25 mg Tab	100/bt	Watson Pharmaceuticals
0904-5986-60	Vitamin D-3 5000 Iu Cap	100/bt	Major Pharmaceutical
50383-917-50	Vitamin D-3 Drops 400 Iu/ml 50 ml	100/bt	Hi-Tech Pharmcal Co. Inc
66582-313-54	Vytorin 10/40 Tab	90/bt	Merck & Co Inc

Otics

Item No.	Description	Qty	Manufacturer
50383-901-10	Acetic Acid Hc Drp 10 ml	1/ea	Hi-Tech Pharmacial Co. Inc
0065-8533-02	Ciprodex Otic Susp 7.5 ml	1/ea	Alcon Labs
65162-702-94	Fluocinolone Acet 0.01% Oil Ear Drp 20 ml	1/bt	Amneal Pharmaceuticals
61314-645-11	Neo/Poly/Hc Otic Af Sus 10 ml	1/ea	Sandoz
6131-4646-10	Neo/Poly/Hc/Otic Soln 10 ml	1/ea	Sandoz

RX/OTC Topicals & Ointments

Item No.	Description	Qty	Manufacturer
0235-0378-82	Americaine 20% Aerosol 2-oz.	1/ea	Emerson Healthcare
OTC052547	Ayr Saline Nasal Gel with Aloe 0.5-oz.	1/ea	B F Ascher & Co Inc
0713-0280-31	Bacitracin 500 Un/gm Ont 1-oz.	1/ea	G & W Laboratories Inc
CUR110556H	Bacitracin+Zinc Oin, 500Un/G, 1-oz.Tube	1/ea	Medline
0501-3203-02	Benadryl M/S Spy 2-oz.	1/ea	Johnson & Johnson
OTC010002	Benzoin Compound Tincture 2-oz.	1/ea	Geritrex Corporation
OTC105656	Benzoyl Peroxide 10% Gel 1.5-oz.	1/ea	Rugby
45802-021-46	Betameth Dip 0.05% Lot 60 ml	1/ea	Perrigo
0168-0040-15	Betameth Val 0.1% Cream 15 gm	1/ea	Fougera & Company
8307-8133-14	Carmex Bp Sqz Tube 12 x 0.35-oz.	12/bx	Carma Laboratories Inc
10223-201-03	Cetacaine Spray 20 Gm	1/ea	Cetylite Industries
0168-0258-15	Clotrim/Betam Dibr 1/.05% Cream 15 gm	1/ea	Fougera & Company
0168-0080-31	Cream, Hydrocortisone, 2.5%, 30 gm	1/bx	Fougera & Company
54162-011-01	Dibucaine 1% Ointment 1 Oz	1/ea	Geritrex Corporation
IIS2282906	Drysol 20% 37.5 ml	1/ea	Person & Covey, Inc.
0096-0707-60	Drysol 20% Sol 60 ml	1/ea	Person And Covey Inc
0096-0707-35	Drysol Daboma Sol, 35 ml	1/ea	Person And Covey Inc
0096-0707-37	Drysol Sol, 37.5 ml	1/ea	Person And Covey Inc
IMS3015SP	Gel 2% Lidocaine 20 ml	25/bx	International Medical
45802-056-35	Gentamic Top, 0.1% Cream, 15 gm	1/ea	Perrigo
45802-046-35	Gentamicin Top 0.1% Ont 15 gm	1/ea	Perrigo
0386-0004-04	Gerbauers Spray & Stretch Spry 3.5-oz.	1/ea	Gerbauer Company
24385-003-46	Gnp Gentian 1% Violet Liq 2-oz.	1/ea	Humco
0283-0871-31	Hurricaine 20% Gel 1-oz.	1/ea	Beutlich Lp
0283-1185-20	Hurricaine Ext Tube 200	200/bx	Beutlich Lp
45802-980-64	Hydroquinone 4% Cream 30 gm	1/ea	Perrigo
76329-3012-5	Lidocaine 2% Jel 25 x 5 ml	25/bx	International Medication
52584-035-05	Lidocaine 2% Jelly Urojet Pfs 20 ml	1/ea	International Medication
76329-3015-5	Lidocaine 2% Jelly Urojet Pfs 25 x 20 ml	25/bx	International Medication
76329-3013-05	Lidocaine 2% Urojet Jel 25 x 10 ml	25/bx	International Medication
0054-3505-47	Lidocaine 4% Topical Soln 50 ml	1/ea	West-Ward Pharmaceuticals
50383-933-55	Lidocaine 5% Ointment 50Gm Jar	1/ea	Hi-Tech Pharmacial Co. Inc
17478-711-10	Lidocaine Hcl 2% Jel 10 x 5 ml	10/bx	Akorn Inc
17478-711-30	Lidocaine Hcl, 2% Jel 30 ml	1/ea	Akorn Inc
0168-0204-37	Lidocaine Top 5% Ont 35 gm	1/ea	Fougera & Company
50383-775-04	Lidocaine Vis, 2% Sol, 100 ml	1/ea	Hi-Tech Pharmacial Co. Inc
50383-667-30	Lidocaine/Priloc, 2.5% Cream, 30 gm	1/ea	Hi-Tech Pharmacial Co. Inc
0496-0882-15	Lmx4 4% Cream 15 gm	1/ea	Eloquest Healthcare
0496-0882-30	Lmx4 4% Cream 30 gm	1/ea	Eloquest Healthcare

RX/OTC Topicals & Ointments Continued

Item No.	Description	Qty	Manufacturer
0496-0883-30	Lmx5 5% Cream 30 gm	1/ea	Eloquest Healthcare
51672-4116-6	Metronidazole 0.75% Top Gel 45 gm	1/ea	Taro Pharmaceutical Industries Ltd
68462-564-17	Mupirocin 2% Crm 15 gm	1/ea	Glenmark Generics Limited
68462-180-22	Mupirocin 2% Oint 22 gm	1/ea	Bionpharma Inc
0501-3712-05	Neosporin Plus Cream 15 gm	1/ea	Johnson & Johnson
1254-7237-79	Neosporin Plus Cream 1-oz.	1/ea	Johnson & Johnson
0281-0326-30	Nitro-Bid 2% Ont 30 gm	1/ea	Savage Labs
0281-0326-60	Nitro-Bid 2% Ont 60 gm	1/ea	Savage Labs
51672-1263-1	Nys/Triam/Ace Cream 15 gm	1/ea	Taro Pharmaceutical
0472-0166-30	Nystatin 100M Un/gm Ont 30 gm	1/ea	Actavis Us
68308-152-15	Nystatin 100M Un/gm Pwd 15 gm	1/ea	Hi-Tek Laboratories Inc
OTC571135	Oxymetazone 0.05% Nasal Spry 15 ml	1/ea	Major Pharmaceutical
0386-0008-03	Pain Ease Med Spry Non-Flam Liq 3.5-oz.	1/ea	Gebauer Company
0574-0601-15	Podocon-25 Liq 15 ml	1/ea	Paddock Laboratories Inc
OTC006009	Polysporin Oint 144 x 0.9 gm. Ud	144/bx	Johnson & Johnson
0501-3798-01	Polysporin Oint, 1-oz. Tube	1/ea	Johnson & Johnson
OTC110231	Povidone Iodine 10 % Oint 1-oz.	1/ea	Major Pharmaceutical
0743-0042-11	Rid For Lice/Scabies Control Spray 5-oz.	1/ea	Bayer Healthcare
50484-010-30	Santyl 250U/gm Oint 30 gm	1/ea	Smith & Nephew
67877-124-40	Silver Sulfadiazine 1% Cream 400 gm	1/ea	Ascend Laboratories
IMS6300	Solution, Lidocaine HCL 4%, Topical	25/cs	International Medical
43598-0210-40	Ssd 1% Jar Cream 400 gm	1/ea	Dr. Reddys Labs
43598-210-55	Ssd 1% Jar Cream 50 gm	1/ea	Dr. Reddys Labs
43598-210-25	Ssd 1% Tube Cream 25 gm	1/ea	Dr. Reddys Labs
IMS3013SP	Syringe, Prefilled Lidocaine HCl Jelly 2	25/bx	International Medical
0536-5150-26	Tolnaftate 1% Pwd 45 gm	1/ea	Watson Pharmaceuticals
0168-0004-80	Triamcino Ace 0.1% Cream 80 gm	1/ea	Fougera & Company
0168-0006-15	Triamcino Ace 0.1% Ont 15 gm	1/ea	Fougera & Company
10481-3008-1	Trichloroacetic Acid Liq, 15 ml	1/ea	Gordon Laboratories
0884-0449-04	Ureacin 20 20% Cream 4-oz.	1/ea	Pedinol Pharmcal Inc
0521-2326-00	Vaseline Jar Jel 3.75-oz.	1/ea	Unilever

Vaccines

Item No.	Product	Description	Size
49281-545-03	ActHIB Vacc	Haemophilus B Conjugate (Hib)	5/bx
49281-400-10	Adacel SDV	Tetanus, Diphtheria, Acellular Pertussis (TDaP) Booster for Adolescents/Adults	10 x 0.5 ml
49281-400-20	Adacel® PFS	Tetanus, Diphtheria, Acellular Pertussis (TDaP) Booster for Adolescents/Adults	5 x 0.5 ml
42023-104-01	Aplisol 5TU	Tuberculin Purified Protein Derivative (PPD), Mantoux, Tuberculin Skin Test (TB)	1 ml 10-Test Vial
42023-104-05	Aplisol® 5TU	Tuberculin Purified Protein Derivative (PPD), Mantoux, Tuberculin Skin Test (TB)	5 ml 50-Test Vial
58160-976-06	Bexsero PFS	Meningococcal Group B	1 x 0.5 ml
58160-976-02	Bexsero PFS	Meningococcal Group B	10 x 0.5 ml
58160-842-11	Boostrix SDV	Tetanus, Diphtheria, Acellular Pertussis (TDaP) Booster for Adolescents/Adults	10 x 0.5 ml
58160-842-52	Boostrix® PFS	Tetanus, Diphtheria, Acellular Pertussis (TDaP) Booster for Adolescents/Adults	10 x 0.5 ml
49281-286-10	Daptacel SDV	Diphtheria, Tetanus Acellular Pertussis (TDaP)-Infants & Children 6 wks-6 yrs	10 x 0.5 ml

Item No.	Product	Description	Size
49281-225-10	Diphtheria & Tetanus Toxoid	Adsorbed (DT)-(Infants & Children 6 wks-6 yrs)	10 x 0.5 ml
58160-821-52	Engerix-B 20 mcg PFS	Hepatitis B (Adult)	10 x 1 ml
58160-821-11	Engerix-B 20 mcg SDV	Hepatitis B (Adult)	10 x 1 ml
58160-820-52	Engerix-B® 10 mcg PFS	Hepatitis B (Pediatric/Adolescent)	10 x 0.5 ml
0006-4119-03	Gardasil SDV	Human Papillomavirus (HPV)-9 valent types (Types 6,11,16,18,31,33,45,52,58)	10/bx
0006-4121-02	Gardasil PFS	Human Papillomavirus (HPV)-9-Valent (Types 6,11,16,18,31,33,45,52,58)	10/bx
58160-825-52	Havrix 720 ELU/0.5 ml PFS	Hepatitis A (Pediatric/Adolescent)	10 x 0.5 ml
58160-825-11	Havrix 720 ELU/0.5 ml SDV	Hepatitis A (Pediatric/Adolescent)	10 x 0.5 ml
58160-826-52	Havrix® 1440 ELU/ml PFS	Hepatitis A (Adult)	10 x 1 ml
43528-002-05	Heplisav B- 20mcg/0.5 ml/ SDV	Hepatitis B	5 x 0.5 ml
58160-818-11	Hiberix Vaccine W/Diluent Sdv	Haemophilus B Conjugate	10/bx
49281-250-51	Imovax Rabies 1 ml	Single Dose Kit; Rabies	1/kt
58160-810-52	Infanrix PFS	Diphtheria, Tetanus, Acellular Pertussis (DTaP)-(Infants & Children 6 wks-7 yrs)	10 x 0.5 ml
58160-810-11	Infanrix® SDV	Diphtheria, Tetanus, Acellular Pertussis (DTaP)-(Infants & Children 6 wks-7 yrs)	10 x 0.5 ml
49281-860-10	IPOL® 10-Dose MDV	Polio Virus (IPV)	5 ml
42515-002-01	Ixiaro PFS 0.5 ml	Japanese Encephalitis	1/bx
58160-812-52	Kinrix PFS	Diphtheria Tetanus, Acellular Pertussis, Polio (DTaP/IPV)	10/bx
58160-812-11	Kinrix® SDV	Diphtheria, Tetanus, Acellular Pertussis, Polio (DTaP/IPV)	10 x 0.5 ml
49281-589-05	Menactra® SDV	Meningococcal	5 x 0.5 ml
46028-208-01	Menveo SDV	Meningococcal	5/kt
0006-4681-00	MMR SDV	Measles, Mumps, Rubella (MMR)	10 x 0.5 ml
58160-811-52	Pediarix® PFS	Diphtheria, Tetanus, Acellular Pertussis, Hepatitis B, Polio (DTaP/Hep B/IPV)	10 x 0.5 ml
0006-4897-00	PedvaxHIB® SDV	Haemophilus B Conjugate (Hib)	10 x 0.5 ml
49281-510-05	Pentacel® SDV	Diphtheria, Tetanus, Acellular Pertussis, Polio, Haemophilus Influenzae B Conjugate (DTaP/IPV/Hib)	5/bx
0006-4837-03	Pneumovax 23 PFs 0.5 ml	Pneumococcal	10/bx
0006-4943-00	Pneumovax SDV	Pneumococcal	10 x 0.5 ml
0005-1971-02	Prevnar 13® PFS	Pneumococcal Conjugate	10 x 0.5 ml
0005-1971-05	Prevnar 13 PFS	Pneumococcal Conjugate	1/ea
0006-4171-00	Proquad SDV	Measles, Mumps, Rubella, Varicella (MMRV)	10/bx
49281-562-10	Quadracel SDV	Diphtheria, Tetanus, Acellular Pertussis, Polio (DTaP/IPV)	10/bx
58160-964-12	Rabavert 1 ml Vial	Rabies	1/kt
0006-4093-02	Recombivax 5 mcg/0.5 ml	Hepatitis B (Pediatric/Adolescent);	10/bx
0006-4992-00	Recombivax HB 40 mcg/ml SDV	Hepatitis B (Dialysis Formulation)	1 ml
0006-4047-41	RotaTeq® (For Oral Use Only)	Rotavirus	10 x 2 ml Single Dose Tubes
0006-4047-20	RotaTeq (For Oral Use Only)	Rotavirus	25 x 2 ml
58160-819-12	Shingrix SDV	Shingles	1/kt
58160-823-11	Shingrix SDV	Shingles	10/kt
49281-215-10	Tenivac Vial	Tetanus & Diphtheria Toxoids Adsorbed (Td)	10 x 0.5ml
49281-215-15	Tenivac® PFS	Tetanus & Diphtheria Toxoids Adsorbed (Td)	10 x 0.5 ml
13533-131-01	Tetanus & Diphtheria Tox Ad SDV 0.5 ml	Tetanus & Diphtheria Toxoids Adsorbed (Td)	10/bx
0005-0100-05	Trumemba PFS 5 x 0.5 ml	Meningococcal Group B	5 x 0.5 ml
0005-0100-10	Trumemba PFS 10 x 0.5 ml	Meningococcal Group B	10 x 0.5 ml
49281-752-21	Tubersol 5TU	Tuberculin Purified Protein Derivative (PPD), Mantoux, Tuberculin Skin Test (TB)	1 ml 10-Test Vial
49281-7522-2	Tubersol® 5TU	Tuberculin Purified Protein Derivative (PPD), Mantoux, Tuberculin Skin Test (TB)	5 ml 50-Test Vial


Vaccines Continued

Item No.	Product	Description	Size
58160-815-52	Twinrix® PFS	Hepatitis A, Hepatitis B	10 x 1 ml
49281-790-20	Typhim Vi MDV	Typhoid	1/ea
49281-790-51	Typhim VI® PFS	Typhoid Vaccine	0.5 ml
0006-4827-00	Varivax SDV	Varicella	10 x 0.5 ml
0006-4096-02	Vaqta 50 u/PFS	Hepatitis A (Adult)	10 x 0.5 ml
0006-4095-02	Vaqta 25 u/PFS	Hepatitis A (Pediatric/Adolescent)	10 x 0.5 ml
0006-4841-00	Vaqta 50 u/SDV	Hepatitis A (Adult)	1/ea
0006-4841-41	Vaqta 50 u/SDV	Hepatitis A (Adult)	10 x 1 ml
0006-4831-41	Vaqta 25 u/0.5 ml SDV	Hepatitis A (Pediatric/Adolescent)	10 x 0.5 ml
69401-000-02	Vivotif® Cap 4	Typhoid Vaccine	4 Oral Capsules
0006-4963-00	Zostavax SDV 0.65 ml	Zoster (Shingles)	1/ea
0006-4963-41	Zostavax SDV 10 x 0.65 ml	Zoster (Shingles)	10/bx



SYRINGES & BLOOD DRAW

Needles/Syringes.....	140-143
Tourniquets	143
Blood Draw	143-146
Infusion Sets	146
IV Catheters.....	147
Blood Collection Tubes	148
Vacutainers	148
IV Tubing	149
Sterile Solutions.....	150
Heel Incision Devices.....	151
Safety Lancets.....	151
Sharps Containers	152



Hypodermic Needles & Syringes

- » High-quality hypodermic needles and syringes
- » Sterile and individually packaged
- » Ready for use for dependable drug delivery

Item No.	Description	Pkg
Luer Lock Syringe and Cap		
SYR101010*	1 ml	800/cs
SYR103010*	3 ml	800/cs
SYR105010*	5 ml	400/cs
SYR110010*	10 ml	400/cs
SYR120010*	20 ml	160/cs
SYR130010*	30 ml	160/cs
SYR160010*	60 ml	160/cs
SYRTCAPZ	Luer Lock Tip Caps	200/bx

Luer Slip Syringe Only

SYR101020*	1 ml	800/cs
SYR103020*	3 ml	800/cs
SYR105020*	5 ml	400/cs
SYR110020*	10 ml	400/cs
SYR120020*	20 ml	160/cs
SYR130020*	30 ml	160/cs
SYR160020*	60 ml	160/cs

Standard Syringe with Needle

SYR101272*	1 ml, TB, 27 Gauge x 1/2"	800/cs
SYR101282*	1 ml, TB, 28 Gauge x 1/2"	800/cs
SYR103205*	3 ml, LL, 20 Gauge x 1"	800/cs
SYR103215*	3 ml, LL, 21 Gauge x 1"	800/cs
SYR103217*	3 ml, LL, 21 Gauge x 1 1/2"	800/cs
SYR103225*	3 ml, LL, 22 Gauge x 1"	800/cs
SYR103227*	3 ml, LL, 22 Gauge x 1 1/2"	800/cs
SYR110235*	3 ml, LL, 23 Gauge x 1"	800/cs
SYR103253*	3 ml, LL, 25 Gauge x 5/8"	800/cs
SYR105205*	5 ml, LL, 20 Gauge x 1"	400/cs

Item No.	Description	Pkg
Standard Syringe with Needle (Cont.)		
SYR105207*	5 ml, LL, 20 Gauge x 1 1/2"	400/cs
SYR105217*	5 ml, LL, 21 Gauge x 1 1/2"	400/cs
SYR110205*	10 ml, LL, 20 Gauge x 1"	400/cs
SYR110207*	10 ml, LL, 20 Gauge x 1 1/2"	400/cs
SYR110215*	10 ml, LL, 21 Gauge x 1"	400/cs
SYR110217*	10 ml, LL, 21 Gauge x 1 1/2"	400/cs
SYR110225*	10 ml, LL, 22 Gauge x 1"	400/cs
SYR110227*	10 ml, LL, 22 Gauge x 1 1/2"	400/cs

Standard Hypodermic Needles

SYR100165*	16 Gauge x 1"	1,000/cs
SYR100167*	16 Gauge x 1 1/2"	1,000/cs
SYR100185*	18 Gauge x 1"	1,000/cs
SYR100187*	18 Gauge x 1 1/2"	1,000/cs
SYR100195*	19 Gauge x 1"	1,000/cs
SYR100197*	19 Gauge x 1 1/2"	1,000/cs
SYR100205*	20 Gauge x 1"	1,000/cs
SYR100207*	20 Gauge x 1 1/2"	1,000/cs
SYR100215*	21 Gauge x 1"	1,000/cs
SYR100225*	22 Gauge x 1"	1,000/cs
SYR100227*	22 Gauge x 1 1/2"	1,000/cs
SYR100235*	23 Gauge x 1"	1,000/cs
SYR100252*	25 Gauge x 1/2"	1,000/cs
SYR100253*	25 Gauge x 5/8"	1,000/cs
SYR100255*	25 Gauge x 1"	1,000/cs
SYR100257*	25 Gauge x 1 1/2"	1,000/cs
SYR100302*	30 Gauge x 1/2"	1,000/cs

Safety Needles

SYRS100185*	18 Gauge x 1"	800/cs
SYRS100187*	18 Gauge x 1 1/2"	800/cs
SYRS100197*	19 Gauge x 1 1/2"	800/cs
SYRS100205*	20 Gauge x 1"	800/cs
SYRS100207*	20 Gauge x 1 1/2"	800/cs
SYRS100215*	21 Gauge x 1"	800/cs
SYRS100217*	21 Gauge x 1 1/2"	800/cs
SYRS100225*	22 Gauge x 1"	800/cs

Item No.	Description	Pkg
Safety Needles (Cont.)		
SYRS100227*	22 Gauge x 1 1/2"	800/cs
SYRS100235*	23 Gauge x 1"	800/cs
SYRS100237*	23 Gauge x 1 1/2"	800/cs
SYRS100255*	25 Gauge x 1"	800/cs
SYRS100257*	25 Gauge x 1 1/2"	800/cs
SYRS100253*	25 Gauge x 5/8"	800/cs
SYRS100273*	27 Gauge x 5/8"	800/cs
SYRS100272*	27 Gauge x 1/2"	800/cs
SYRS100276*	27 Gauge x 1 1/4"	800/cs
SYRS100302*	30 Gauge x 1/2"	800/cs

Safety Syringes with Needle

SYRSI193292*	Insulin, 0.5 ml, 29 Gauge x 1/2"	500/cs
SYRS101235F*	1 ml, 23 Gauge x 1"	1,200/cs
SYRS101255F*	1 ml, 25 Gauge x 1"	1,200/cs
ULT101253*	TB, 1 ml, 25 Gauge x 5/8"	500/cs
SYRS101272F*	1 ml, 27 Gauge x 1/2"	1,200/cs
SYRSI101292*	Insulin, 1 ml, 29 Gauge x 1/2"	500/cs
SYRS101392F*	1 ml, 30 Gauge x 1/2"	1,200/cs
SYRS103205*	3 ml, 20 Gauge x 1"	1,200/cs
SYRS103207*	3 ml, 20 Gauge x 1 1/2"	1,200/cs
SYRS103215*	3 ml, 21 Gauge x 1"	1,200/cs
SYRS103217*	3 ml, 21 Gauge x 1 1/2"	1,200/cs
SYRS103225*	3 ml, 22 Gauge x 1"	1,200/cs
SYRS103227*	3 ml, 22 Gauge x 1 1/2"	1,200/cs
SYRS103235*	3 ml, 23 Gauge x 1"	1,200/cs
SYRS103253*	3 ml, 25 Gauge x 5/8"	1,200/cs
SYRS103255*	3 ml, 25 Gauge x 1"	1,200/cs
SYRS105205*	5 ml, 20 Gauge x 1"	1,200/cs
SYRS105207*	5 ml, 20 Gauge x 1 1/2"	1,200/cs
SYRS105215*	5 ml, 21 Gauge x 1"	1,200/cs
SYRS105227*	5 ml, 22 Gauge x 1 1/2"	1,200/cs
SYRS110215*	10 ml, 21 Gauge x 1"	900/cs
SYRS110227*	10 ml, 22 Gauge x 1 1/2"	900/cs

*Add "Z" to item to order by the box.



SYRA101272



Allergy Syringe Trays

- » Designed to deliver small doses of medication easily and accurately for efficient patient care
- » Smooth-gliding rubber stopper promotes accurate medication measurement
- » Low dead-space minimizes allergen waste
- » Easy-to-read number scale helps you prepare precise dose
- » Sterile packaged 25 per tray for easy accessibility; trays stack for ease and convenience

Item No.	Description	Pkg
SYRA101272	27 Gauge x 0.5", 1 ml	40/cs
SYRA101261	26 Gauge x 3/8", 1 ml	40/cs



EMZ11240



EMZE060055



ZR® Saline Flush Syringes

- » The ZR difference is clearly, part of the solution
- » Utilizes anti-reflux technology
- » Enhanced label helps lower the risk of medication error
- » Syringe clarity improves visual inspection of solution
- » All ZR syringes are terminally sterilized

Item No.	Description	Pkg
EMZ111240	10 ml Saline in 10 ml Syringe	240/cs, 30/bx
EMZ151240	5 ml Saline in 10 ml Syringe	240/cs, 30/bx
EMZ131240	3 ml Saline in 10 ml Syringe	240/cs, 30/bx
EMZE060055	100 u/ml Heparin 5 ml in 10 ml	600/cs
EMZE050055	10 u/ml Heparin 5 ml in 10 ml	600/cs



B-D309657



Standard Syringes, Graduated

- » Features clear barrels with bold scale markings, tapered plunger rods for ease of aspiration and positive plunger rod stops
- » Luer-Lok thread for increased secure connection
- » Sterile
- » Not made with natural rubber latex

Item No.	Description	Pkg
B-D309657Z	3 ml, LL, 1/10 ml	200/bx
B-D309646Z	5 ml, LL, 1/5 ml	125/bx
B-D302995Z	10 ml, LL, 1/5 ml	200/bx
B-D302832Z	30 ml, LL, 1 ml	56/bx



B-D309626Z



Tuberculin Syringe with Detachable Needle, Slip Tip

- » Helps improve dosage accuracy and virtually eliminates needle pop-off
- » Bold, clear scale markings make easy and accurate dosage measurements

Item No.	Description	Pkg
B-D309626Z	1 ml, 25 Gauge x 5/8"	100/bx



BHL2N3379Z



Baxter® INTERLINK® Injection Site, Male Luer Lock Adapter

- » Sterile, nonpyrogenic
- » Not made with natural rubber latex

Item No.	Description	Pkg
BHL2N3379Z	Injection Site Adapter	200/bx



B-D305500Z



1 ml Allergy Syringe with PrecisionGlide® Permanently Attached Needle

Item No.	Description	Pkg
B-D305500Z	28 Gauge x 1/2" Regular Bevel	100/bx



B-D309570Z



Standard Syringes with Needle

Item No.	Description	Pkg
B-D309570Z	3 ml, 25 Gauge x 5/8", LL	100/bx
B-D309581Z	3 ml, 25 Gauge x 1"	100/bx
B-D309571Z	3 ml, 23 Gauge x 1", LL	100/bx
B-D309589Z	3 ml, 23 Gauge x 1 1/2", Intramuscular, Thin Wall	100/bx
B-D309572Z	3 ml, 22 Gauge x 1", LL	100/bx
B-D309643Z	10 ml, 21 Gauge x 1 1/2"	100/bx



B-D305128Z



Regular Beveled Needles

Item No.	Description	Pkg
B-D305128Z	30 Gauge x 1"	100/bx
B-D305136Z	27 Gauge x 1 1/4"	100/bx
B-D305122Z	25 Gauge x 5/8"	100/bx
B-D305127Z	25 Gauge x 1 1/2"	100/bx
B-D305155Z	22 Gauge x 1"	100/bx
B-D305156Z	22 Gauge x 1 1/2"	100/bx
B-D305195Z	18 Gauge x 1"	100/bx
B-D305196Z	18 Gauge x 1 1/2"	100/bx



B-D305124b



PrecisionGlide™ Needles

- » Needle cannula is made of stainless steel
- » Polypropylene hub is designed to fit any slip tip, reusable or single-use
- » Strips of five sterile, fiber-free, peel-apart packages
- » Shield protecting point and cannula
- » Package is color-coded according to needle gauge
- » Lot control number on package

Item No.	Description	Pkg
B-D305106Z	30 Gauge x ½"	100/bx
B-D305125	25 Gauge x 1"	1,000/cs
B-D305125Z	25 Gauge x 1"	100/bx
B-D305124Z	25 Gauge x ¾"	100/bx



B-D305758Z



Luer-Lok® Syringe with Eclipse Needles

Item No.	Description	Pkg
1 ml with Detachable Needle		
B-D305778Z	30 Gauge x ½"	50/bx
B-D305789Z	27 Gauge x ½"	50/bx
B-D305780Z	25 Gauge x ⅝"	50/bx
3 ml with Detachable Needle		
B-D305781Z	25 Gauge x ⅝"	50/bx
B-D305782Z	23 Gauge x 1"	50/bx
B-D305783Z	22 Gauge x 1 ½"	50/bx



B-D305758Z



Eclipse™ Safety Needles

- » Pivoting-shield needle allows single finger stroke activation
- » No hard surface required

Item No.	Description	Pkg
B-D305761Z	25 Gauge x 1", Use LL	100/bx
B-D305767Z	25 Gauge x 1 ½", Use LL	100/bx
B-D305762Z	23 Gauge x 1"	100/bx
B-D305763Z	22 Gauge x 1 ½"	100/bx
B-D305766Z	18 Gauge x 1 ½", Use LL	100/bx



B-D305950



SafetyGlide® 1 ml Allergist Tray with Permanently Attached Needle

Item No.	Description	Pkg
B-D305950	27 Gauge x ½" Regular Bevel	1,000/cs
B-D305950Z	27 Gauge x ½" Regular Bevel	25/tr



B-D305900



SafetyGlide Needles

- » Activation-Assist™ technology allows for fast, easy, single-handed shielding


Item No.	Description	Pkg
B-D305900Z	22 Gauge x 1 ½"	50/bx
B-D305901Z	25 Gauge x ⅝"	50/bx
B-D305902Z	23 Gauge x 1"	50/bx
B-D305916Z	25 Gauge x 1"	50/bx

1 ml Syringe with Needle

Item No.	Description	Pkg
B-D305903Z	25 Gauge x ⅝" Subcutaneous Needle	50/bx

3 ml Syringe with Needle

Item No.	Description	Pkg
B-D305904ZZ	25 Gauge x ⅝" Subcutaneous Needle	50/bx
B-D305906Z	22 Gauge x 1 ½" Intramuscular Needle	50/bx



B-D305270



Integra® 3 ml Retracting Syringes with Detachable Needle

- » Syringes with Tru-Lok technology are designed to have significantly lower waste space volume as compared to traditional luer connections

Item No.	Description	Pkg
B-D305269Z	25 Gauge x ⅝"	100/bx
B-D305270Z	25 Gauge x 1"	100/bx
B-D305271Z	23 Gauge x 1"	100/bx



DYNDE75020

Disposable Tourniquets

- » Offers an unbeatable combination of value, quality and ease of use
- » Features a comfortable texture that allows the tourniquet to remain secure while in use
- » Promotes a single-use policy
- » DYND75020, DYND75020H, DYND75021 and JPS16332 are not made with natural rubber latex

Item No.	Description	Pkg
----------	-------------	-----

Rolled Tourniquets

DYND75020	1 x 18"/2.5 x 45.7 cm, Blue	100/cs
DYND75020H	1 x 18"/2.5 x 45.7 cm, Blue	1/ea
JPS16332	1 x 18"/2.5 cm x 45.7 cm, White, Bulk	250/cs
DYND75026	1 x 18"/2.5 cm x 45.7 cm, Orange, Bulk	50/bx, 20 bx/cs, 1,000/cs
DYND75021	¾ x 18"/1.9 x 45.7 cm, Blue	100/cs
DYNDE75020	1 x 18"/2.5 cm x 45.7 cm, Blue, Dispenser Box	50/bx, 20 bx/cs, 1,000/cs

Flat Tourniquets

JPS16332	1 x 18"/2.5 x 45.7 cm, White, Bulk	250/cs
----------	------------------------------------	--------

Disposable Latex Tourniquets

NON494109	1 x 18"/2.5 x 45.7 cm, Bulk	100/cs
NON494109H	1 x 18"/2.5 x 45.7 cm, Bulk	1/ea



B-D305946



SafetyGlide Tuberculin Syringes

- » Offers single-handed activation with minimal change in technique

Item No.	Description	Pkg
B-D305945Z	1 ml, 26 Gauge x ½"	100/bx
B-B-D305946Z	26 Gauge x ¾"	100/bx



NPMNM20G25Z



Multi-Sample Blood Collection Needles

- » Multi-sample needles for essential venipuncture
- » Fits most phlebotomy tube holders

Item No.	Description	Pkg
NPMNM22G38Z	22 Gauge x 1 ½"	100/bx
NPMNM22G25Z	22 Gauge x 1"	100/bx
NPMNM21G38Z	21 Gauge x 1 ½"	100/bx
NPMNM21G25Z	21 Gauge x 1"	100/bx
NPMNM20G38Z	20 Gauge x 1 ½"	100/bx
NPMNM20G25Z	20 Gauge x 1"	100/bx



DYNDTN0657



Blunt Tip Cannulas

Item No.	Description	Pkg
DYNDTN0657	Soft Pack	100/cs



SWD216066H



Monoject Blood Collection Needles, Thin Wall

Item No.	Description	Pkg
SWD216066	22 Gauge x 1"	1,000/cs
SWD216033	21 Gauge x 1"	1,000/cs



B-D364815



Blood Collection Needles/Tube Holders

Item No.	Description	Pkg
B-D364815Z	BD Vacutainer® One-Use, Non-Stackable Holder	250/bx
SWD610102	Monoject, Regular 13 mm Diameter	100/bx



OnGuard Blood Collection Holder

- » Compatible with most blood collection sets, needles and tubes

Item No.	Description	Pkg
DYNDBCHLDR	Blood Collection Tube Holder	1,000/cs



B-D368607Z



Vacutainer® Eclipse Safety Blood Collection Needles

- » Single-handed activation and immediate protection
- » Needles with luer adapters and pre-attached holders are designed to help your institution comply with OSHA standards
- » BD PrecisionGlide needle technology ensures highest quality and sharpness
- » Each holder is engraved with the words "DO NOT REUSE" and "SINGLE USE ONLY" to reinforce single-use compliance


Item No.	Description	Pkg
B-D368608CSZ	22 Gauge x 1 ¼"	48/bx
B-D368607Z	21 Gauge x 1 ¼"	48/bx
B-D368651CS	22 Gauge x 1 ¼" Pre-Attached Holder	100/cs
B-D368650CS	21 Gauge x 1 ¼" Pre-Attached Holder	100/cs



Multiple Sample Luer Adapters

- » Provides the security of a threaded, locking luer connection, replacing a luer slip device
- » Compatible with a female luer connection
- » Sterile, secure and safe specimen sampling

Item No.	Description	Pkg
B-D367290	Vacutainer Multiple Sample Luer Adapters	1,000/cs
SWD225257Z	Monoject Luer Adapter	100/bx



B-D367297Z



Safety-Lok Blood Collection and Infusion Sets

- » Designed to shield healthcare workers from exposed needles
- » Manually activated safety shield fully covers needle
- » Not made with natural rubber latex

Item No.	Description	Pkg
B-D367285Z	25 Gauge x ¾", with 12" Tubing	50/bx
B-D367297Z	23 Gauge x ¾", with 12" Tubing, Without Luer Adapter	50/bx
B-D367292Z	23 Gauge x ¾", with 7" Tubing, with Luer Adapter	50/bx




SWD225281



Blood Collection Sets

» Blood collection sets with female luer adapter attached are used for routine venipuncture applications

Item No.	Description	Pkg
B-D367281Z	BD Safety-Lok 21 Gauge x 3/4", 12"/31 cm Tubing, Green Wing	50/bx
B-D367283Z	BD Safety-Lok 23 Gauge x 3/4", 12"/31 cm Tubing, Lt. Blue Wing	50/bx
SWD225299	Monoject 21 Gauge x 3/4", 12"/31 cm Tubing, Green Wing with Multi-Sample Luer	50/cs
SWD225190	ANGEL WING® 23 Gauge x 3/4", 12"/31 cm Tubing, Blue Wing with Female Luer	50/cs
SWD225281	Monoject 19 Gauge x 3/4", 12"/31 cm Tubing, Cream Wing	50/cs




KAWDBMA21G



K-Shield Advantage Safety Winged Blood Collection Sets

- » Fast and simple one-step safety activation with an "audible" click
- » Needle is covered during removal, ensuring confidence in safety
- » Kink resistant and memory free tubing
- » Multi-sample adapter luer locks securely onto female connector for ease of removal with a syringe draw
- » Sharper, double beveled venipuncture needle with a silicone coating for added patient comfort and less pain

Item No.	Description	Pkg
KAWDBMA21G	21 Gauge Needle, Green Wing Color, 3/4" Needle Length, 12" Tubing	200/cs
KAWDBMA21GZ	21 Gauge Needle, Green Wing Color, 3/4" Needle Length, 12" Tubing	50/bx
KAWDBMA23G	23 Gauge Needle, Light Blue Wing Color, 3/4" Needle Length, 12" Tubing	200/cs
KAWDBMA23GZ	23 Gauge Needle, Light Blue Wing Color, 3/4" Needle Length, 12" Tubing	50/cs
KAWDBMA25G	25 Gauge Needle, Dark Blue Wing Color, 3/4" Needle Length, 12" Tubing	200/cs
KAWDBMA25GZ	25 Gauge Needle, Dark Blue Wing Color, 3/4" Needle Length, 12" Tubing	50/cs



KAWDBMH21G



K-Shield Advantage Winged Blood Collection Sets with Pre-Attached Tube Holder

- » Includes pre-attached tube holder for added convenience
- » Needle is covered during removal, ensuring confidence in safety
- » Fast and simple one-step safety activation with an "audible" click
- » Kink resistant and memory free tubing

Item No.	Description	Pkg
KAWDBMH21G	21 Gauge Needle, Green Wing Color, 3/4" Needle Length, 12" Tubing	100/cs
KAWDBMH21GZ	21 Gauge Needle, Green Wing Color, 3/4" Needle Length, 12" Tubing	50/bx
KAWDBMH23G	23 Gauge Needle, Light Blue Wing Color, 3/4" Needle Length, 12" Tubing	100/cs
KAWDBMH23GZ	23 Gauge Needle, Light Blue Wing Color, 3/4" Needle Length, 12" Tubing	50/cs
KAWDBMH25G	25 Gauge Needle, Dark Blue Wing Color, 3/4" Needle Length, 12" Tubing	100/cs
KAWDBMH25GZ	25 Gauge Needle, Dark Blue Wing Color, 3/4" Needle Length, 12" Tubing	50/cs



B-D367324Z



Push-Button Safety

- » Safety at the push of a button provides split-second protection against needle-stick injuries
- » In-vein activation irreversibly retracts the needle, locking it safely out of harms way

Item No.	Description	Pkg
B-D367324Z	Vacutainer 23 Gauge x 3/4", with 12"/31 cm Tubing, without Luer	50/bx
B-D367342Z	Vacutainer 23 Gauge x 3/4", with 12"/31 cm Tubing, Lt. Blue Wing with Luer Adapter	50/bx
B-D367344Z	Vacutainer 21 Gauge x 3/4", with 12"/31 cm Tubing, Lt. Green Wing with Luer Adapter	50/bx



BMG471690



Huber Safety Needles

- » Female luer connector, on/off clamp, 90° WHIN Safe Huber needle
- » Lipid-compatible, pressure rated
- » Priming vol: 0.23 ml, length: 7-inch/17.8 cm
- » DEHP-free, not made with natural rubber latex

Item No.	Description	Pkg
BMG573114	20 Gauge x 1"	100/cs



KAWD3K223GZ



Small Vein Infusion Sets

- » Thin wall needles for adequate flow rates
- » Ideal for use on children or adults with smaller veins
- » Flexible wings with textured surface easily pinch together to assure a secure grip
- » Kink-resistant tubing with female luer lock connector

Item No.	Description	Pkg
KAWD3K223GZ	23 Gauge x 3/4"/2 cm Needle; 12"/31 cm Tubing	50/bx



Safety Small Vein Butterfly Infusion Sets

- » Small vein infusion sets are designed for effective use on children or adults with smaller veins
- » Flexible wings with textured surface easily pinch together to assure a secure grip
- » Kink-resistant tubing

Item No.	Description	Pkg
NPMSPR21G19	21 Gauge x 3/4" Safety	500/cs, 50/bx
NPMSPR23G19	23 Gauge x 3/4" Safety	500/cs, 50/bx
NPMSPR25G19	25 Gauge x 3/4" Safety	500/cs, 50/bx



Standard IV Catheters

- » Features a sharp bevel, thin-walled inner needle that assures an easy puncture
- » Color-coded translucent hub allows simple identification of catheter size and detection of blood flashback

Item No.	Description	Pkg
NPMCI1450*	14 Gauge x 2"	200/cs, 50/bx
NPMCI1632*	16 Gauge x 1 1/4"	200/cs, 50/bx
NPMCI1650*	16 Gauge x 2"	200/cs, 50/bx
NPMCI1832*	18 Gauge x 1 1/4"	200/cs, 50/bx
NPMCI1850*	18 Gauge x 2"	200/cs, 50/bx
NPMCI2025*	20 Gauge x 1"	200/cs, 50/bx
NPMCI2032*	20 Gauge x 1 1/4"	200/cs, 50/bx
NPMCI2050*	20 Gauge x 2"	200/cs, 50/bx
NPMCI2225*	22 Gauge x 1"	200/cs, 50/bx
NPMCI2232*	22 Gauge x 1 1/4"	200/cs, 50/bx
NPMCI2420*	24 Gauge x 3/4"	200/cs, 50/bx
NPMCI2425*	24 Gauge x 1"	200/cs, 50/bx

*Add "Z" to item to order by the box.



DYNCSA20100



DYNCSA22100



SureSite IV™ Auto Safety IV Catheter

- » Push-button safety activation
- » Fully encapsulating needle shield
- » Notched needle for initial flash observation in sizes 20 - 24 gauge
- » Flashback chamber for clear visual confirmation
- » Manufactured with PUR (Polyurethane)

Item No.	Description	Pkg
DYNCSA14175	14 Gauge x 1.75", Orange	200/cs, 50/bx
DYNCSA16177	16 Gauge x 1.77", Gray	200/cs, 50/bx
DYNCSA18116	18 Gauge x 1.16", Green	200/cs, 50/bx
DYNCSA18188	18 Gauge x 1.88", Green	200/cs, 50/bx
DYNCSA20100	20 Gauge x 1", Pink	200/cs, 50/bx
DYNCSA20116	20 Gauge x 1.16", Pink	200/cs, 50/bx
DYNCSA20188	20 Gauge x 1.88", Pink	200/cs, 50/bx
DYNCSA22100	22 Gauge x 1", Blue	200/cs, 50/bx
DYNCSA24075	24 Gauge x 0.75", Yellow	200/cs, 50/bx



DYNCS24340



DYNCS22100



SureSite IV Slide Safety IV Catheter

- » Clicking lock safety mechanism
- » Self-sheathing needle shield
- » Flashback chamber for clear visual confirmation
- » Manufactured with PUR (Polyurethane)

Item No.	Description	Pkg
DYNCS214114	14G x 1 1/4", Orange	200/cs, 50/bx
DYNCS216114	16G x 1 1/4", Grey	200/cs, 50/bx
DYNCS218114	18G x 1 1/4", Green	200/cs, 50/bx
DYNCS220114	20G x 1 1/4", Pink	200/cs, 50/bx
DYNCS220100	20G x 1", Pink	200/cs, 50/bx
DYNCS222100	22G x 1", Blue	200/cs, 50/bx
DYNCS224340	24G x 3/4", Yellow	200/cs, 50/bx



Insyte™ Autoguard™ Shielded IV Catheters

- » With the push of a button, the needle instantly retracts, reducing the risk of accidental needle-stick injuries
- » Made of Vialon™, a unique proprietary biomaterial proven to enhance clinical and economic outcomes
- » Gives the clinician control over when to activate safety mechanism
- » Non-winged, sterile

Item No.	Description	Pkg
B-D381444Z	18 Gauge x 1.16"/1.3 x 30 mm, Retractable	50/bx
B-D381423Z	22 Gauge x 1"/0.9 x 25 mm, Notched Needle	50/bx
B-D381412Z	24 Gauge x 3/4"/0.7 x 19 mm, Notched Needle	50/bx



Vacutainer® Blood Collection Tubes

Item No.	Tube Type	Tube Size (mm)	Draw Volume	Color	Additive/Concentration	Pkg
Serum						
B-D367986Z	Plastic	13 x 100	5 ml	Gold	Clot Activator/Polymer Gel	100/bx
B-D367988Z	Plastic	16 x 100	8.5 ml	Red/Gray	Clot Activator/Polymer Gel	100/bx
B-D367815Z	Plastic	13 x 100	6 ml	Red/Gray	Clot Activator, Silicone-Coated	100/bx
B-D366430Z	Glass	16 x 100	10 ml	Red	Silicone-Coated	100/bx
Plasma						
B-D367960Z	Plastic	13 x 75	3 ml	Lt Green	Lithium Heparin 56 USP, Polymer Gel	100/bx
B-D367961Z	Plastic	13 x 100	3.5 ml	Lt Green	Lithium Heparin 65 USP, Polymer Gel	100/bx
B-D367962Z	Plastic	13 x 100	4.5 ml	Lt Green	Lithium Heparin 84 USP, Polymer Gel	100/bx
Fluoride						
B-D367587Z	Plastic	13 x 75	2 ml	Gray	Sodium Fluoride 3 mg, Na2EDTA 6 mg	100/bx
B-D367925Z	Plastic	13 x 100	6 ml	Gray	Sodium Fluoride 15 mg, Potassium Oxalate 12 mg	100/bx
Heparin						
B-D367871Z	Plastic	13 x 75	4 ml	Green	Sodium Heparin 75 USP	100/bx
Whole Blood						
B-D367856Z	Plastic	13 x 75	3 ml	Lavender	K2EDTA 5.4 mg	100/bx
B-D367844Z	Plastic	13 x 75	4 ml	Lavender	K2EDTA 7.2 mg	100/bx
B-D367861Z	Plastic	13 x 75	4 ml	Lavender	K2EDTA 7.2 mg	100/bx
B-D366450Z	Plastic	13 x 100	7 ml	Lavender	K3EDTA 12.15 mg 15% Sol, 0.081 ml	100/bx
Citrate						
B-D369714Z	Glass	13 x 75	4.5 ml	Lt Blue	Buff Sodium Citrate 105 M, 3.2%	100/bx
B-D363083Z	Plastic	13 x 75	2.7 ml	Lt Blue	Citrate	100/bx



greiner bio-one

VACUETTE® Low-Volume Blood Collection Tubes

- » Standard PET plastic tubes with color-coded safety caps
- » Availability of additives required for laboratory testing
- » Optimal blood to additive ratios for reduced blood volumes
- » Physician office friendly sizes available

Item No.	Description	Pkg
GRI454236	2 ml Z Serum Clot Activator, 13 x 75 Red Cap-White Ring	1,200/cs
GRI454081	1 ml LH Lithium Heparin, 13 x 75 Green Cap-White Ring, PREMIUM	1,200/cs
GRI454237	2 ml LH Lithium Heparin, 13 x 75 Green Cap-White Ring	1,200/cs
GRI454302	2 ml NH Sodium Heparin Separator, 13 x 75 Green Cap-Green Ring	1,200/cs
GRI454052	1 ml K2E K2EDTA, 13 x 75 Lavender Cap-White Ring	1,200/cs
GRI454428CS	2 ml K2E K2EDTA, 13 x 75 Lavender Cap-White Ring	1,200/cs
GRI454222	2 ml K3E K3EDTA, 13 x 75 Lavender Cap-White Ring	1,200/cs
GRI454320	1 ml 9NC Coagulation Sodium Citrate 3.2%, 13 x 75 Blue Cap-White Ring, PREMIUM	1,200/cs
GRI454322	2 ml 9NC Coagulation Sodium Citrate 3.2%, 13 x 75 Blue Cap-White Ring	1,200/cs
GRI454238	2 ml FX Sodium Fluoride/Potassium Oxalate, 13 x 75 Grey Cap-White Ring	1,200/cs



DYNDTN1512



DYNDTN1812



Gravity IV Administration Sets

- » Crystal clear tubing and drip chambers
- » Flat, clear needle-free Y-sites
- » Not compatible with infusion pumps

Item No.	Description	Pkg
DYNDTN1529	92", 1 Needle-Free Y-Site, 15 drp/ml	50/cs
DYNDTN1512	96", 2 Needle-Free Y-Sites, 15 drp/ml	50/cs
DYNDTN1545	109", 3 Needle-Free Y-Sites, 15 drp/ml	50/cs
DYNDTN0529	92", 1 Needle-Free Y-Site, 60 drp/ml	50/cs
DYNDTN0512	96", 2 Needle-Free Y-Sites, 60 drp/ml	50/cs
DYNDTN0555	111", 3 Needle-Free Y-Sites, 60 drp/ml	50/cs
DYNDTN1812	96", 2 Split Septum Y-Sites, 15 drp/ml	50/cs
DYNDTN1530	92", 1 Injection Y-Site, 15 drp/ml	50/cs



Secondary IV Administration Sets

Item No.	Description	Pkg
DYNDTN1540A	40" IV Secondary Set with Hanger	50/cs



DYNDTC2000A



IV Accessories

Item No.	Description	Pkg
DYNDTC2000A	SmartSite® Needle-Free Valve	100/cs
DYNDTN0606	Mini Vented Transfer Pin	50/cs
DYND75030	Intermittent Injection Site	200/cs, 50/bx
DYNDTN0657	Blunt Tip Cannula	100/cs



DYNDTC5081



DYNDTN0003



IV Extension Sets

- » Clear, flexible tubing
- » Male and female luer locks
- » Not made with DEHP

Item No.	Description	Pkg
DYNDTC5077	6.5" Set with Non-Bonded Needle-Free Valve, Non-Removable Slide Clamp, Non-Rotating Male Luer Lock, 0.8 ml Priming Volume	100/cs
DYNDTC5081	8" Set with Non-Bonded Needle-Free Valve, Non-Removable Slide Clamp, Rotating Male Luer Lock, 0.3 ml Priming Volume	100/cs
DYNDTC5002	7" Set with Bonded Needle-Free Valve, Non-Removable Slide Clamp, Rotating Male Luer Lock, 0.1 ml Priming Volume	100/cs
DYNDTC5082	8" Set with Non-Bonded Needle-Free Valve, Needle-Free T-Port Connector, Non-Removable Slide Clamp, 0.4 ml Priming Volume	100/cs
DYNDTN1002	6" Set with Non-Removable Slide Clamp, Non-Rotating Male Luer Lock, 0.2 ml Priming Volume	100/cs
DYNDTN0001	7" Set with Removable Slide Clamp, Non-Rotating Male Luer Lock, 0.2 ml Priming Volume	100/cs
DYNDTN0003	12" Set with Non-Removable Slide Clamp, Non-Rotating Male Luer Lock, 0.3 ml Priming Volume	100/cs
DYNDTN0007	60" Set with Non-Removable Slide Clamp, Non-Rotating Male Luer Lock, 0.6 ml Priming Volume	100/cs



BMGL850101

BHL2B0061H



Sterile IV Solutions

Item No.	Description	Pkg
ABB079830262H	Hospira®, Sodium Chloride 0.9%, 8.45-oz./250 ml	1/ea
ABB79830939H	Hospira, Sodium Chloride 0.9%, 33.81-oz./1,000 ml	1/ea
BHL2B1324XH	Sodium Chloride Solution 0.9%, 33.81-oz./1,000 ml	1/ea
BHL2B0062QH	5% Dextrose Injection USP, 8.45-oz./250 ml	1/ea
BHL2B0088	5% Dextrose Injection USP, 1.69-oz./50 ml	96 bg/cs
BHL2B0089	5% Dextrose Injection USP, 3.38-oz./100 ml, Multi-Pack	96 bg/cs
BHL2B1308Z	0.9% NaCl Injection USP, 1.69-oz./50 ml, Multi-Pack	16/bx
BHL2B1309Z	0.9% Sodium Chloride, 3.38-oz./100 ml, Multi-Pack VI	16/bx
BHL2B2322Q	Lactated Ringer's Injection, USP, 8.45-oz./250 ml	36/cs
BHL2B1322QH	0.9% Sodium Chloride Injection, USP, 8.45-oz./250 ml	1/ea
BMGL8001H	B Braun®, Injection, Sodium Chloride 0.9%, 16.91-oz./500 ml	1/ea
BMGL850101	B Braun, IV Solution, H ₂ O, Sterile, 16.91-oz./500 ml	24/cs
BMGR520101H	B Braun, Sterile Saline (0.9% NaCl), 16.91-oz./500 ml	1/ea
BMGR520001H	B Braun, Sterile Saline (0.9% NaCl), 33.81-oz./1,000 ml	1/ea
BMGL8502H	B Braun, Injection, Sterile Water, 16.91-oz./500 ml	1/ea
BMGL8500H	B Braun, Injection, Sterile Water, 33.81-oz./1,000 ml	1/ea



BHL2F7113HH

Baxter

Sterile Irrigation Solutions

Item No.	Description	Pkg
Sterile Water for Irrigation		
BHL2F7113CS	USP, 16.91-oz./500 ml	18/cs
BHL2F7113HH	USP, 16.91-oz./500 ml	1/ea
BHL2B7114X	33.81-oz./1,000 ml	14/cs
BHL2B7114XH	33.81-oz./1,000 ml	1/ea
Sodium Chloride 0.9%		
BHL2F7122HH	8.45-oz./250 ml	1/ea
BHL2F7123HH	16.91-oz./500 ml	1/ea



Sterile Water and Saline Solutions

- » Prepackaged, sterile and ready-to-use containers provide just the right amount of solution necessary for procedures
- » For moistening of wound dressings, wound debridement, and device irrigation
- » Not for irrigation of wound or body cavities
- » Not for injection

Item No.	Description	Pkg
Sterile Water		
RDI30295	Screw-Top Bottle: 3.38-oz./100 ml	48/cs
DYND40570	Foil-Lid Cup: 3.72-oz./110 ml	48/cs
PCS1550	Screw-Top Bottle: 8.45-oz./250 ml	24/cs
DYND500MLW	Bottle: 16.91-oz./500 ml	18/cs

Sterile Saline 0.9%		
RDI30296	Screw-Top Bottle: 3.38-oz./100 ml	48/cs
DYND40540	Foil-Lid Cup: 3.72-oz./110 ml	48/cs
PCS1650	Screw-Top Bottle: 8.45-oz./250 ml	24/cs
DYND500MLS	Bottle: 16.91-oz./500 ml	18/cs

Glacier Distilled Water		
SAL5415047	1-gal./3.8 L	3/cs



Heel Incision Devices

- » High-quality stainless steel blades make small, precise incisions for reduced bruising
- » Disposable, one-time use devices
- » Easy trigger activates blade quickly, safely and precisely
- » Color-coded for easy size identification of devices

Item No.	Description	Pkg
MPHMPREEMIE2	Micro-Preemie Heel Lancet, Light Blue	50/bx, 4 bx/ cs
MPHMPREEMIE2Z	Micro-Preemie Heel Lancet, Light Blue	50/bx
MPHPREEMIE2	Preemie Heel Lancet, White	50/bx, 4 bx/ cs
MPHPREEMIE2Z	Preemie Heel Lancet, White	50/bx
MPHNEWBORN2	Newborn Heel Lancet, Pink	50/bx, 4 bx/ cs
MPHNEWBORN2Z	Newborn Heel Lancet, Pink	50/bx
MPHTODDLER2	Toddler Heel Lancet, Purple	50/bx, 4 bx/ cs
MPHTODDLER2Z	Toddler Heel Lancet, Purple	50/bx



MPHADJDEPTH



Adjustable-Depth Safety Lancets

- » Innovative rotating dial allows for three easy depth adjustments: 1.3 mm, 1.8 mm, 2.3 mm
- » Pressure activated
- » Compares to the Accu-Chek Safe-T-Pro Plus

Item No.	Description	Pkg
MPHADJDEPTH	Adjustable-Depth Safety Lancet	1,000/cs
MPHADJDEPTHZ	Adjustable-Depth Safety Lancet	200/bx



Safety Lancets

- » Simple and safe for healthcare professionals and patients
- » Needle retracts automatically and cannot be re-used
- » Different sizes allow healthcare professionals to obtain the proper amount of blood for capillary testing
- » Double spring action and controlled needle depth provides greater patient comfort
- » Quick and easy-to-use

Item No.	Description	Pkg
Push Button		
MPHSAFETY28	28 Gauge, 1.6 mm	200/bx, 18 bx/cs
MPHSAFETY28Z	28 Gauge, 1.6 mm	200/bx
MPHSAFETY285	28 Gauge, 1.6 mm	50/bx, 24 bx/cs
MPHSAFETY285Z	28 Gauge, 1.6 mm	50/bx
MPHSAFETY281	28 Gauge, 1.6 mm	100/bx, 20 bx/cs
MPHSAFETY281Z	28 Gauge, 1.6 mm	100/bx
MPHSAFETY23	23 Gauge, 1.8 mm	200/bx, 10 bx/cs
MPHSAFETY23Z	23 Gauge, 1.8 mm	200/bx
MPHSAFETY21	21 Gauge, 1.8 mm	200/bx, 18 bx/cs
MPHSAFETY21Z	21 Gauge, 1.8 mm	200/bx
MPHSAFETYBL	1.6 mm Blade; 1.5 mm Depth	200/bx, 18 bx/cs
MPHSAFETYBLZ	1.6 mm Blade; 1.5 mm Depth	200/bx

Pressure Activated

MPHST28	28 Gauge, 1.8 mm	200/bx, 10 bx/cs
MPHST28Z	28 Gauge, 1.8 mm	200/bx
MPHST21	21 Gauge, 2.2 mm	200/bx, 10 bx/cs
MPHST21Z	21 Gauge, 2.2 mm	200/bx

Side Button

MPHSTSIDE28	28 Gauge, 1.8 mm	200/bx, 10 bx/cs
MPHSTSIDE28Z	28 Gauge, 1.8 mm	200/bx
MPHSTSIDE21	21 Gauge, 2.4 mm	200/bx, 10 bx/cs
MPHSTSIDE21Z	21 Gauge, 2.4 mm	200/bx

Contact Activated

MHPRESS30	30 Gauge, 1.5 mm	100/bx, 20 bx/cs
MHPRESS30Z	30 Gauge, 1.5 mm	100/bx
MHPRESS28	28 Gauge, 1.8 mm	100/bx, 20 bx/cs
MHPRESS28Z	28 Gauge, 1.8 mm	100/bx
MHPRESS26	26 Gauge, 1.8 mm	100/bx, 20 bx/cs
MHPRESS26Z	26 Gauge, 1.8 mm	100/bx
MHPRESS23	23 Gauge, 2.2 mm	100/bx, 20 bx/cs
MHPRESS23Z	23 Gauge, 2.2 mm	100/bx
MHPRESS21	21 Gauge, 2.2 mm	100/bx, 20 bx/cs
MHPRESS21Z	21 Gauge, 2.2 mm	100/bx



MDS705110H



Phlebotomy Sharps Containers

- » Designed to fit into blood drawing trays for point-of-use disposal
- » Square base for added stability
- » Translucent top allows easy identification of contents
- » Accommodates a variety of sharps
- » Containers lock for final disposal and safe transport
- » Brackets available

Item No.	Description	Pkg
----------	-------------	-----

Phlebotomy Container

MDS705110H	1-qt./0.9 L, 6.25 H x 4.25 L x 4.5"W/159 x 108 x 114 mm	1/ea
MDS707011	Bracket for MDS705110H	5 ea/cs

Transportable Container

MDS705115H	1.5-qt./1.4 L, Screw Top Lid	1/ea
------------	------------------------------	------



MDS705153H Shown with
Cabinet MDS707953H



Patient Room Sharps Containers

- » Designed for use in patient rooms, on med carts and in high-traffic areas
- » Lid design facilitates one-handed operation
- » Horizontal drop design for maximum capacity to fill line
- » Easy-to-see fill levels

Item No.	Description	Pkg
----------	-------------	-----

MDS705153H	5-qt./4.7 L, Red, 11 H x 4.75 L x 10.75"W/279 x 121 x 273 mm	1/ea
MDS705154H	5-qt./4.7 L, Clear, 11 H x 4.75 L x 10.75"W/279 x 121 x 273 mm	1/ea
MDS707953H	5-qt./4.7 L with Locking Cabinet	1/ea
MDS705203H	3-gal./11.4 L, Red, 20.5 H x 6 L x 14"W/521 x 152 x 356 mm	1/ea
MDS707930H	3-gal./11.4 L, Locking Cabinet	1/ea



MDS705308



Large PG-II Containers

- » Provides safe disposal for regulated medical waste and they comply with D.O.T. requirements to eliminate the need for secondary packaging
- » Gasketed lid resists leakage
- » Integrated side handles make the container easy to carry

Item No.	Description	Pkg
MDS705208	8-gal./30.3 L, Flat Lid, Red 13.5 H x 13 L x 17.3"W/343 x 330 x 439 mm	9/cs
MDS705308	8-gal./30.3 L, Port Lid, Red 13.5 H x 13 L x 17.3"W/343 x 330 x 439 mm	9/cs



MDS705202H



Multipurpose Sharps/Biohazard Containers

- » Designed for use in restricted access areas
- » Containers lock for final disposal
- » Nestable containers save space and reduce shipping costs
- » Sliding lid design accommodates a variety of sharps sizes
- » Horizontal drop maximizes use of container volume
- » Brackets available

Item No.	Description	Pkg
MDS705201H	1-gal./3.8 L, Sliding Lid	1/ea
MDS705202H	2-gal./7.6 L, Sliding Lid	1/ea
MDS705208SLH	8-gal./30.3 L, Sliding Lid	1/ea



ATT1G1V4H



Mail Back Sharps Containers

Item No.	Description	Pkg
ATT1Q1V4H	1-qt. /0.9 L	1/ea
ATT1G1V4H	1-gal./3.8 L	1/ea
ATT2G1V4H	2-gal./7.6 L	1/ea



TRADITIONAL WOUND CARE

First Aid Ointments & Creams.....	154
Sponges	155-156
Abdominal/Combine Pads.....	157
Gauze Dressings.....	158-159
Elastic Bandages.....	160
Self-Adherent Wrap.....	161
Adhesive Bandages	162-163
Tape	163-164
Wound Closure Strips	164
Casting Tape.....	164
Additional Supplies	165

ADVANCED WOUND CARE

Wound Bed Prep.....	166
Skin Substitute	166
Moisture Balance	166
Infection/Inflammation	167
Post-Op.....	168
Superabsorbent Dressings.....	168
Fiber Dressings	169
Alginate Dressings	169
Compression Systems.....	169-170
Antimicrobial Gels/Powders	170
Hydrocolloids/Collagen Dressings.....	170
Honey Dressings.....	171
Skin Protectants.....	171
Additional Supplies	171-172



CUR5305



Germ Shield™ Antimicrobial Wound Gel

- » Antimicrobial gel with silver as the active ingredient
- » Formulated to offer broad spectrum protection against bacteria
- » Applies easily and absorbs thoroughly
- » Non-stinging, non-staining barrier that helps hydrate the wound area

Item No.	Description	Pkg
CUR5305	0.5-oz./14.2 g Antimicrobial Gel	1/bx, 12 bx/cs



CUR01496



Mediplast® Wart Pads

- » The active ingredient salicylic acid (40%), is widely recommended for the removal of corns, calluses and warts

Item No.	Description	Pkg
CUR01496	2 x 3"/5.08 x 7.62 cm	150/cs
CUR01496H	2 x 3"/5.08 x 7.62 cm	1/ea



CUR00354Z



Vitamin A & D Ointment

- » Active ingredients: white petrolatum 93.5%. Inactives: corn oil, light mineral oil, vitamin A palmitate and vitamin D
- » Not made with natural rubber latex

Item No.	Description	Pkg
CUR003501	2-oz./56.6 g Tube	12/cs
CUR00354Z	0.18-oz./5 g Packets	144/bx



CUR001209Z



Triple Antibiotic Ointment

- » Helps prevent infection in minor cuts, scrapes and burns
- » This product is intended for use on the body only
- » Active: Bacitracin Zinc 400 units, Neomycin Sulfate 5 mg, Polymyxin B Sulfate 5,000 units
- » Not made with natural rubber latex

Item No.	Description	Pkg
CUR001231	1-oz./28.4 g Tube	12/cs
CUR001209Z	0.03-oz./0.9 g Packets	144/bx



CUR005331



Petroleum Jelly

- » Temporarily protects minor cuts, scrapes and burns
- » Active: White Petrolatum 100%
- » Not made with natural rubber latex

Item No.	Description	Pkg
CUR005331	1-oz./28.4 g Tube	12/cs
CUR005345Z	0.18-oz./5 g Packets	144/bx



CUR110556



Bacitracin with Zinc

- » Helps prevent infection in minor cuts, scrapes and burns
- » Active ingredients: bacitracin zinc 500 units
- » Not made with natural rubber latex

Item No.	Description	Pkg
CUR110556	1-oz./28.4 g Tube	12/cs
CUR001109Z	0.03-oz./0.9 g Packets	144/cs



Woven Gauze Sponges

- » 100% cotton
- » Ideal for wound dressings, wound packing and general wound care
- » Strictly controlled manufacturing assures clean, debris-free packing with folded edges to prevent unraveling
- » Sterile 2's packed in coated paper envelopes to reduce fiber debris
- » Easy-open envelopes carefully sealed to prevent the intrusion of dust and contaminants
- » Not made with natural rubber latex

Premium

- » "H fold" design offers a premium approach to minimizing loose threads and lint

Item No.	Description	Pkg
Sterile 2's USP Type VII		
NON21420H	2 x 2"/5.1 x 5.1 cm, 8-ply	50 pk/bx
NON21422H	3" x 3"/7.6 x 7.6 cm, 12-ply	40 pk/bx
NON21448H	4 x 4"/10.2 x 10.2 cm, 8-ply	50 pk/bx
NON21424H	4 x 4"/10.2 x 10.2 cm, 12-ply	25 pk/bx
NON21842Z	8 x 4"/20.3 x 10.2 cm, 12-ply	25 pk/bx

Sterile 10's (Trays) USP Type VII

NON21425Z	4 x 4"/10.2 x 10.2 cm, 8-ply	10 ea/tray
NON21426H	4 x 4"/10.2 x 10.2 cm, 12-ply	10 ea/tray
NON21428H	4 x 4"/10.2 x 10.2 cm, 16-ply	10 ea/tray

Non-Sterile (Bulk)

NON25208H	2 x 2"/5.1 x 5.1 cm, 8-ply	200/sleeve
NON25212H	2 x 2"/5.1 x 5.1 cm, 12-ply	200/sleeve
NON25312H	3 x 3"/7.6 x 7.6 cm, 12-ply	200/sleeve
NON25408H	4 x 4"/10.2 x 10.2 cm, 8-ply	200/sleeve
NON25412H	4 x 4"/10.2 x 10.2 cm, 12-ply	200/sleeve
NON25416Z	4 x 4"/10.2 x 10.2 cm, 16-ply	200/sleeve
NON25812ZZ	8 x 4"/20.3 x 10.2 cm, 12-ply	200/sleeve



Medline Essentials

- » "C fold" design minimizes loose threads and lint

Item No.	Description	Pkg
Sterile 2's		
PRM2208Z	2 x 2"/5.1 x 5.1 cm, 8-ply	2/pk, 50 pk/bx
PRM3312Z	3 x 3"/7.6 x 7.6 cm, 12-ply	2/pk, 40 pk/bx
PRM4408Z	4 x 4"/10.2 x 10.2 cm, 8-ply	2/pk, 50 pk/bx
PRM4412H	4 x 4"/10.2 x 10.2 cm, 12-ply	2/pk

Sterile 1's

PRM21419	2 x 2"/5.1 x 5.1 cm, 12-ply	100 pk/bx, 24 bx/cs, 2,400 pk/cs
PRM21423	4" x 4"/10.2 x 10.2 cm, 12-ply	100 pk/bx, 12 bx/cs, 1,200 pk/cs

Non-Sterile

PRM21208CZ	2 x 2"/5.1 x 5.1 cm, 8-ply	200/sleeve
PRM21212CZ	2 x 2"/5.1 x 5.1 cm, 12-ply	200/sleeve
PRM21312CZ	3 x 3"/7.6 x 7.6 cm, 12-ply	200/sleeve
PRM21408CZ	4 x 4"/10.2 x 10.2 cm, 8-ply	200/sleeve
PRM21412CZ	4 x 4"/10.2 x 10.2 cm, 12-ply	200/sleeve
PRM21416CZ	4 x 4"/10.2 x 10.2 cm, 16-ply	200/sleeve

Cotton Filled Dental Sponges

PRM22DS	2 x 2"/5.1 x 5.1 cm	200/sleeve, 25 pk/cs, 5,000/cs
---------	---------------------	-----------------------------------



Cotton Gauze Sponges

- » Sterile, 100% cotton, woven gauze sponges are ideal for wound dressings, wound packing and general wound care
- » "C fold" design helps minimize loose threads and lint
- » Easy-open envelopes are carefully sealed to prevent the intrusion of dust and contaminants
- » Not made with natural rubber latex

PRM21426



Item No.	Description	Pkg
Sterile 10's (Soft Pouch) USP Type VII		
PRM21426	4 x 4"/10.16 x 10.16 cm, 12-Ply	10 sponges/pk, 128 pk/cs, 1,280/cs
PRM21428	4 x 4"/10.16 x 10.16 cm, 16-Ply	10 sponges/pk, 128 pk/cs, 1,280/cs



Non-Woven Gauze Sponges

- » Significantly faster wicking ability promotes a healthy healing environment
- » Discover the perfect balance between strength and absorbency
- » Rayon and polyester blend
- » Produces less lint and is less traumatic to fragile granulating wound tissue than traditional woven cotton gauze
- » Virtually non-adherent
- » Not made with natural rubber latex



Premium

Item No.	Description	Pkg
Sterile 2's		
NON21224H	2 x 2"/5.1 x 5.1 cm, 4-ply	50 pk/bx, 100 ea/bx
NON21334Z	3 x 3"/7.6 x 7.6 cm, 4-ply	40 pk/bx, 80 ea/bx
NON21444H	4 x 4"/10.2 x 10.2 cm, 4-ply	50 pk/bx, 100 ea/bx
NON21446Z	4 x 4"/10.2 x 10.2 cm, 6-ply	25 pk/bx, 50 ea/bx
Sterile 10's (Trays)		
NON21449H	4 x 4"/10.2 x 10.2 cm, 4-ply	10/tray
Non-Sterile (Bulk)		
NON25223Z	2 x 2"/5.1 x 5.1 cm, 3-ply	200/sleeve
NON25224Z	2 x 2"/5.1 x 5.1 cm, 4-ply	200/sleeve
NON25334Z	3 x 3"/7.6 x 7.6 cm, 4-ply	200/sleeve
NON25443Z	4 x 4"/10.2 x 10.2 cm, 3-ply	200/sleeve
NON25444H	4 x 4"/10.2 x 10.2 cm, 4-ply	200/sleeve
NON25449	4 x 4"/10.2 x 10.2 cm, 4-ply	25/sleeve, 100 sleeve/cs

Medline Essentials

Item No.	Description	Pkg
Sterile 2's		
PRM21224	2 x 2"/5.1 x 5.1 cm, 4-ply	50 pk/bx, 30 bx/cs, 1500 pk/cs
PRM21444	4 x 4"/10.2 x 10.2 cm, 4-ply	50 pk/bx, 12 bx/cs, 600 pk/cs
PRM21444H	4 x 4"/10.16 x 10.16 cm, 4-Ply	50 pk/bx
Non-Sterile (Bulk)		
PRM25224	2 x 2"/5.1 x 5.1 cm, 4-ply	200/sleeve, 40 sleeves/cs
PRM25224H	2 x 2"/5.08 x 5.08 cm, 4-Ply	200/sleeve, 40 sleeve/cs
PRM25444	4 x 4"/10.2 x 10.2 cm, 4-ply	200/sleeve, 10 sleeves/cs
PRM25444H	4 x 4"/10.16 x 10.16 cm, 4-Ply	200/sleeve, 10 sleeves/cs



Drain Sponges

- » Rayon/polyester blend, 6-ply, highly absorbent, low-linting and fast wicking
- » Precisely fenestrated for easy and consistent application
- » Unique fenestration allows sponge to lie flat when wrapped around a drainage tube
- » Pre-cut sponge eliminates the need for sterile scissors and conserves nursing time
- » All the benefits of non-woven
- » Not made with natural rubber latex

Premium

Item No.	Description	Pkg
NON256001Z	2" x 2"/5.1 x 5.1 cm, Sterile 2's	2/pk, 35/pk/bx
NON256000H	4" x 4"/10.2 x 10.2 cm, Sterile 2's	2/pk, 35/pk/bx

Medline Essentials

Item No.	Description	Pkg
PRM256001H	2" x 2"/5.1 x 5.1 cm	2/pk, 35/pk/bx
PRM256000H	4" x 4"/10.2 x 10.2 cm	2/pk, 35/pk/bx



Abdominal (ABD) Pads

- » Features extra absorbency and wicking for heavily draining wounds
- » Soft non-woven outer layer quickly wicks fluid to a cellulose center
- » Cellulose quickly absorbs and disperses fluids laterally to prevent pooling
- » Hydrophobic back impedes strike-through
- » Sealed edges prevent linting
- » Not made with natural rubber latex

Premium

Item No.	Description	Pkg
Sterile 1's		
NON21450H	5 x 9"/12.7 x 22.9 cm	25/bx
NON21453H	8 x 7.5"/20.3 x 19.1 cm	20/bx
NON21454H	8 x 10"/20.3 x 25.4 cm	18/bx
Non-Sterile (Bulk)		
NON21451	5 x 9"/12.7 x 22.9 cm	576/cs
NON21452	8 x 7.5"/20.3 x 19.1 cm	576/cs
NON21456	8 x 10"/20.3 x 25.4 cm	432/cs
NON21457	12 x 16"/30.5 x 40.6 cm	144/cs



NON21450



PRM21450

Medline Essentials

Item No.	Description	Pkg
Sterile 1's		
PRM21450Z	5 x 9"/12.7 x 22.9 cm	25/bx
PRM21450H	5 x 9"/12.7 x 22.9 cm	1/ea
PRM21453Z	8 x 7.5"/20.3 x 19.1 cm	20/bx
PRM21453H	8 x 7.5"/20.3 x 19.1 cm	1/ea
PRM21454Z	8 x 10"/20.3 x 25.4 cm	16/bx
Non-Sterile (Bulk)		
PRM21451	5 x 9"/12.7 x 22.9 cm	576/cs
PRM21456	8 x 10"/20.3 x 25.4 cm	432/cs



Sterile Non-Adherent Pads

Compare to Covidien Telfa®*

- » This sterile, absorbent cotton filler and polyester outer layer will not interfere with granulating tissue by sticking to the wound
- » Soft perforated film allows air to circulate and the absorbent pad allows the fluid to collect
- » Can easily be cut to fit any shape without separating
- » Not made with natural rubber latex

Item No.	Description	Pkg
Non-Adherent Pads		
NON25700H	2 x 3"/5.08 x 7.62 cm	1/pk, 100/bx
NON25710Z	3 x 4"/7.62 x 10.16 cm	1/pk, 100/bx
NON25720Z	3 x 8"/7.62 x 20.32 cm	1/pk, 50/bx



NON25700H

Item No.	Description	Pkg
Non-Adherent Pads with Adhesive		
NON25700A	2 x 3"/5.08 x 7.62 cm	1/pk, 100 pk/bx, 12 bx/cs, 1,200, ea/cs
NON25710A	3 x 4"/7.62 x 10.16 cm	1/pk, 100 pk/bx, 12 bx/cs, 1,200 ea/cs

*This product is not manufactured or distributed by Covidien LP, makers of Telfa® Non-Adherent Pads.

For even more selection, including Post-Op Gauze Sponges or x-ray detectable gauze, visit www.medline.com/physician for more information.



Conforming Bandages

- » Rayon and polyester blend
- » Gently secures to body contours without restricting movement
- » Accommodates an increase or decrease in swelling
- » Uniquely constructed weave and finished edges help eliminate unraveling and lint
- » 1-ply bandages are ideal for secondary dressing applications and provides stretch memory when applied
- » Not made with natural rubber latex

Premium

Compare to Kendall Conform®*

Item No.	Description	Pkg
Sterile 1-Ply (Singles), Relaxed		
NON254955H	1 x 75"/2.5 x 190.5 cm	12/bx
NON25496H	2 x 75"/5.1 x 190.5 cm	12/bx
NON25497H	3 x 75"/7.6 x 190.5 cm	12/bx
NON25498H	4 x 75"/10.2 x 190.5 cm	12/bx
NON25499H	6 x 80"/15.2 x 203.2 cm	6/bx

Non-Sterile 1-Ply (Bulk), Relaxed

NON25491H	1 x 75"/2.5 x 190.5 cm	12/bg
NON25492H	2 x 75"/5.1 x 190.5 cm	12/bg
NON25493H	3 x 75"/7.6 x 190.5 cm	12/bg
NON25494H	4 x 75"/10.2 x 190.5 cm	12/bg
NON25495H	6 x 80"/15.2 x 203.2 cm	6/bg

*This product is not manufactured or distributed by Medtronic, makers of Conform® gauze bandages.

Medline Essentials

Item No.	Description	Pkg
Non-Sterile		
PRM254192	2 x 4.1-yd./5.1 x 3.8 m	8 bx/cs
PRM254193	3 x 4.1-yd./7.6 x 3.8 m	8 bx/cs
PRM354194	4 x 4.1-yd./10.2 x 3.8 m	8 bx/cs
PRM354195	6 x 4.1-yd./15.2 x 3.8 m	8 bx/cs
Sterile		
PRM354196	2 x 4.1-yd./5.1 x 3.8 m	8 bx/cs
PRM354197	3 x 4.1-yd./7.6 x 3.8 m	8 bx/cs
PRM354198	4 x 4.1-yd./10.2 x 3.8 m	8 bx/cs
PRM354199	6 x 4.1-yd./15.2 x 3.8 m	8 bx/cs



100% Cotton Gauze Bandages

- » Ideal as a primary or secondary dressing
- » Open weave design provides fast wicking, absorbency and aeration
- » Tightly finished edge minimizes lint
- » Increased crimps for improved cushioning and patient comfort
- » 6-ply bandage stretches and conforms for a variety of patient needs, especially burn care
- » Not made with natural rubber latex

Premium Bulkee II

Compare to Kerlix®*

Item No.	Description	Pkg
Sterile 6-Ply (Singles), Stretched		
NON25865H	4.5" x 4.1-yd./11.4 cm x 3.8 m	1/ea
NON25861H	3.4" x 3.6-yd./8.6 cm x 3.3 m	1/ea
NON25850H	2.25" x 3-yd./5.7 cm x 2.7 m	1/ea

Non-Sterile 6-Ply, Stretched

NON25855H	4.5" x 4.1-yd./11.4 cm x 3.8 m	1/ea
NON25855	4.5" x 4.1-yd./11.4 cm x 3.8 m	100 rolls/cs

*This product is not manufactured by Medtronic, makers of Kerlix® dressings, bandages and sponges.

Medline Essentials

Item No.	Description	Pkg
Sterile 6-Ply		
PRM25865	4.5" x 4.1-yd./11.4 cm x 3.8 m	100 rolls/cs
PRM25865P	4.5" x 4.1-yd./11.4 cm x 3.8 m	5 rolls/pk
PRM25865H	4.5" x 4.1-yd./11.4 cm x 3.8 m	1/ea
Non-Sterile 6-Ply		
PRM25855	4.5" x 4.1-yd./11.4 cm x 3.8 m	100 rolls/cs



Bulkee II® Super Fluff Sponges, Sterile

Compare to Kerlix® Fluff Sponges*

- » All the great features of Bulkee II in a sponge
- » An excellent primary dressing with ideal loft, bulk and absorbency
- » Effectively manages exudate in and around the wound
- » Ideal for pre-op prepping and cleansing, debriding and packing
- » Not made with natural rubber latex

Item No.	Description	Pkg
NON25852Z	6 x 6.75"/15.2 x 17.1 cm, Diagonally	40/bx
NON25853H	6 x 6.75"/15.2 x 17.1 cm, Diagonally	5/tray
NON25854H	6 x 6.75"/15.2 x 17.1 cm, Diagonally	10/tray

*This product is not manufactured by Medtronic, makers of Kerlix® dressings, bandages and sponges.

More gauze dressings are available on www.medline.com/physician



Petrolatum Gauze, Sterile

- » An effective, non-adherent primary dressing that maintains a moist wound bed for autolytic healing
- » USP cotton gauze is evenly impregnated with USP white petrolatum to create a soothing, conforming dressing that won't stick to the wound site
- » Ideal for use on non-draining wounds
- » Not made with natural rubber latex

Item No.	Description	Pkg
CUR251272Z	0.5 x 72"/1.25 x 1.83 cm	12/bx
CUR251180Z	1 x 8"/2.54 x 20.32 cm	50/bx
CUR251136Z	1 x 36"/2.54 x 91.44 cm	12/bx
CUR251390Z	3 x 9"/7.62 x 22.86 cm	50/bx
CUR251318Z	3 x 18"/7.62 x 45.72 cm	12/bx
CUR251336Z	3 x 36"/7.62 x 91.44 cm	12/bx
CUR251636Z	6 x 36"/15.24 x 91.44 cm	12/bx



CUR253180Z



Xeroform Gauze Dressings

- » Ideal for use on non-draining to lightly draining wounds
- » Petrolatum gauze with the addition of 3% Xeroform (Bismuth Tribromophenate), a deodorizing agent
- » Contraindicated for heavily draining wounds
- » Not made with natural rubber latex

Item No.	Description	Pkg
CUR253180WZ	1 x 8"/2.54 x 20.32 cm	50/bx
CUR253590WZ	5 x 9"/12.7 x 22.86 cm	50/bx

Xeroform Petrolatum Overwrap Gauze Dressing, Sterile

Item No.	Description	Pkg
CUR253180Z	1 x 8"/2.54 x 20.32 cm	50/bx
CUR253220Z	2 x 2"/5.08 x 5.08 cm	25/bx
CUR253440Z	4 x 4"/10.16 x 10.16 cm	25/bx
CUR253590Z	5 x 9"/12.7 x 22.86 cm	50/bx
CUR253200H	4" x 3-yd./10.16 cm x 2.74 m	1/ea



Packing Strips

- » Tightly woven, 100% cotton gauze packing strips are ravel resistant and won't fray or lint
- » Strips are available plain or impregnated with iodoform, an antiseptic agent
- » Packaged as 5-yard strip per bottle; Cut to desired length
- » Not made with natural rubber latex

Item No.	Description	Pkg
Plain, Sterile		
NON255145H	¼" x 5-yd./0.625 cm x 4.57 m	1/ea
NON255125H	½" x 5-yd./1.25 cm x 4.57 m	1/ea
NON255015H*	1" x 5-yd./2.5 cm x 4.57 m	1/ea
NON255025H	2" x 5-yd.	1/ea

Iodoform, Sterile

NON256145H	¼" x 5-yd./0.625 cm x 4.57 m	1/ea
NON256125H	½" x 5-yd./1.25 cm x 4.57 m	1/ea
NON256015H	1" x 5-yd./2.5 cm x 4.57 m	1/ea
NON256025H	2" x 5-yd.	1/ea



CUR250383Z



Oil Emulsion Gauze Dressing, Steriles

- » Ideal for most draining wounds
- » Knitted, high porosity cellulose acetate is evenly impregnated with USP white petrolatum in an emulsion blend
- » Permits the flow of exudates without adhering to granulating tissue and minimizing skin maceration
- » Highly conformable dressing
- » Use directly on the wound, beneath an absorbent secondary dressing
- » Not made with natural rubber latex

Item No.	Description	Pkg
CUR250383Z	3 x 8"/7.62 x 20.3 cm	3/pk, 108/bx
CUR250330Z	3 x 3"/7.62 x 7.62 cm	50/bx
CUR250381Z	3 x 8"/7.62 x 20.3 cm	1/pk, 24/bx
CUR250317Z	3 x 16"/7.62 x 40.64 cm	36/bx



Matrix® Elastic Bandages

- » Ideal for medium compression situations, the polyester/cotton weave provides excellent elasticity and recoverability
- » Features an easy-to-use hook-and-loop closure at both ends
- » Color-coded for easy size identification on both ends of roll
- » Not made with natural rubber latex

Item No.	Description	Pkg
----------	-------------	-----

Sterile (Singles), Stretched

DYNJ05152LFH	2" x 5-yd./ 5.08 cm x 4.57 m	1/ea
DYNJ05153LFH	3" x 5-yd./ 7.62 cm x 4.57 m	1/ea
DYNJ05154LFH	4" x 5-yd./ 10.16 cm x 4.57 m	1/ea
DYNJ05155LFH	4" x 10-yd./ 10.16 cm x 9.14 m	1/ea
DYNJ05156LFH	6" x 5-yd./ 15.2 cm x 4.57 m	1/ea
DYNJ05157LFH	6" x 10-yd./ 15.2 cm x 9.14 m	1/ea
DYNJ05158LFH	6" x 15-yd./ 15.2 cm x 13.72 m	1/ea

Non-Sterile, Stretched

MDS087002LFZ	2" x 5-yd./ 5.08 cm x 4.57 m	10/bx
MDS087002LFH	2" x 5-yd./ 5.08 cm x 4.57 m	1/ea
MDS087003LFZ	3" x 5-yd./ 7.62 cm x 4.57 m	10/bx
MDS087003LFH	3" x 5-yd./ 7.62 cm x 4.57 m	1/ea
MDS087004LFZ	4" x 5-yd./ 10.16 cm x 4.57 m	10/bx
MDS087004LFH	4" x 5-yd./ 10.16 cm x 4.57 m	1/ea
MDS087006LFZ	6" x 5-yd./ 15.2 cm x 4.57 m	10/bx
MDS087006LFH	6" x 5-yd./ 15.2 cm x 4.57 m	1/ea
MDS087104LFZ	4" x 10-yd./ 10.16 cm x 9.14 m	10/bx
MDS087106LFZ	6" x 10-yd./ 15.2 cm x 9.14 m	10/bx
MDS087106LFH	6" x 10-yd./ 15.2 cm x 9.14 m	1/ea
MDS087156LFH	6" x 15-yd./ 15.2 cm x 13.72 m	1/ea



Swift-Wrap® Elastic Bandages

- » Convenient hook-and-loop closure, eliminating the need for clips or tape
- » Cotton/polyester knit construction, ideal for high-compression application on a variety of conditions
- » Not made with natural rubber latex

Item No.	Description	Pkg
----------	-------------	-----

Sterile (Singles), Stretched

DYNJ05145H	3" x 5-yd./ 7.62 cm x 4.57 m	1/ea
DYNJ05146H	4" x 5-yd./ 10.16 cm x 4.57 m	1/ea
DYNJ05147H	6" x 5-yd./ 15.2 cm x 4.57 m	1/ea

Non-Sterile, Stretched

MDS077002Z	2" x 5-yd./ 5.08 cm x 4.57 m	20/cs
MDS077002H	2" x 5-yd./ 5.08 cm x 4.57 m	1/ea
MDS077003Z	3" x 5-yd./ 7.62 cm x 4.57 m	20/cs
MDS077003H	3" x 5-yd./ 7.62 cm x 4.57 m	1/ea
MDS077004Z	4" x 5-yd./ 10.16 cm x 4.57 m	20/cs
MDS077004H	4" x 5-yd./ 10.16 cm x 4.57 m	1/ea
MDS077006Z	6" x 5-yd./ 15.2 cm x 4.57 m	20/cs
MDS077006H	6" x 5-yd./ 15.2 cm x 4.57 m	1/ea

To view more elastic bandages visit www.medline.com/physician for more information.



Sure-Wrap® Elastic Bandages

- » Cost-effective elastic bandage
- » Delivers consistent pressure and is appropriate for applications requiring higher compressions.
- » Clips included
- » 10% cotton/90% polyester
- » Not made with natural rubber latex

Item No.	Description	Pkg
----------	-------------	-----

MDS055002Z	2" x 5-yd./ 5.08 cm x 4.57 m	20/cs
MDS055002H	2" x 5-yd./ 5.08 cm x 4.57 m	1/ea
MDS055003Z	3" x 5-yd./ 7.62 cm x 4.57 m	20/cs
MDS055003H	3" x 5-yd./ 7.62 cm x 4.57 m	1/ea
MDS055004Z	4" x 5-yd./ 10.16 cm x 4.57 m	20/cs
MDS055004H	4" x 5-yd./ 10.16 cm x 4.57 m	1/ea
MDS055006Z	6" x 5-yd./ 15.20 cm x 4.57 m	20/cs
MDS055006H	6" x 5-yd./ 15.20 cm x 4.57 m	1/ea



Soft-Wrap® Elastic Bandages

- » Knitted, breathable, suited for general use support and compression
- » Serves as an ideal secondary bandage
- » Not made with natural rubber latex
- » 90% cotton/10% polyester
- » Clips included
- » Ideal for low compression

Item No.	Description	Pkg
----------	-------------	-----

MDS046002ZZ	2" x 5-yd./ 5.08 cm x 4.57 m	10/bx
MDS046003ZZ	3" x 5-yd./ 7.62 cm x 4.57 m	10/bx
MDS046004ZZ	4" x 5-yd./ 10.16 cm x 4.57 m	10/bx
MDS046006ZZ	6" x 5-yd./ 15.20 cm x 4.57 m	10/bx



CoFlex®

Compare to 3M Coban®*

- » Soft foam layer
- » Performs better than latex
- » Features patented hand tear technology
- » Will not constrict
- » Not made with natural rubber latex



Item No.	Description	Pkg
----------	-------------	-----

CoFlex LF², Tan, Stretched, Non-Latex, Sterile

DYNJ089002	2" x 5-yd./5.08 cm x 4.57 m	20 rls/cs
DYNJ089003	3" x 5-yd./7.62 cm x 4.57 m	20 rls/cs
DYNJ089003H	3" x 5-yd./7.62 cm x 4.57 m	1/ea
DYNJ089004	4" x 5-yd./10.16 cm x 4.57 m	20 rls/cs
DYNJ089004H	4" x 5-yd./10.16 cm x 4.57 m	1/ea
DYNJ089006	6" x 5-yd./15.2 cm x 4.57 m	20 rls/cs
DYNJ089006H	6" x 5-yd./15.2 cm x 4.57 m	1/ea

CoFlex LF², Tan, Stretched, Non-Latex, Non-Sterile

MDS089001	1" x 5-yd./2.54 cm x 4.57 m	30 rls/cs
MDS089001H	1" x 5-yd./2.54 cm x 4.57 m	1/ea
MDS089002	2" x 5-yd./5.08 cm x 4.57 m	36 rls/cs
MDS089002H	2" x 5-yd./5.08 cm x 4.57 m	1/ea
MDS089003	3" x 5-yd./7.62 cm x 4.57 m	24 rls/cs
MDS089003H	3" x 5-yd./7.62 cm x 4.57 m	1/ea
MDS089004	4" x 5-yd./10.16 cm x 4.57 m	18 rls/cs
MDS089004H	4" x 5-yd./10.16 cm x 4.57 m	1/ea
MDS089006	6" x 5-yd./15.2 cm x 4.57 m	12 rls/cs

Item No.	Description	Pkg
----------	-------------	-----

CoFlex LF², Color Pack, Stretched, Non-Latex, Non-Sterile

MDS089003CP	3" x 5-yd./7.62 cm x 4.57 m	24 rls/cs
MDS089006CP	6" x 5-yd./15.2 cm x 4.57 m	12 rls/cs

Co-Flex Med, Beige, Latex, Sterile

DYNJ086002	Med, 2" x 5-yd./5.08 cm x 4.57 m	20/cs
DYNJ086002H	Med, 2" x 5-yd./5.08 cm x 4.57 m	1/ea
DYNJ086003	Med, 3" x 5-yd./7.62 cm x 4.57 m	20/cs
DYNJ086003H	Med, 3" x 5-yd./7.62 cm x 4.57 m	1/ea
DYNJ086004	Med, 4" x 5-yd./10.16 cm x 4.57 m	20/cs
DYNJ086004H	Med, 4" x 5-yd./10.16 cm x 4.57 m	1/ea
DYNJ086006	Med, 6" x 5-yd./15.2 cm x 4.57 m	20/cs
DYNJ086006H	Med, 6" x 5-yd./15.2 cm x 4.57 m	1/ea

Co-Flex Med, Latex, Non-Sterile

MDS086001B	Med, 1" x 5-yd./2.54 cm x 4.57 m, Blue	30/ea
MDS086002	Med, 2" x 5-yd./5.08 cm x 4.57 m, Tan	36/ea
MDS086002H	Med, 2" x 5-yd./5.08 cm x 4.57 m, Tan	1/ea
MDS086003	Med, 3" x 5-yd./7.62 cm x 4.57 m, Tan	24/ea
MDS086003H	Med, 3" x 5-yd./7.62 cm x 4.57 m, Tan	1/ea
MDS086004	Med, 4" x 5-yd./10.16 cm x 4.57 m, Tan	18/ea
MDS086004H	Med, 4" x 5-yd./10.16 cm x 4.57 m, Tan	1/ea

CoFlex AFD, Tan, Sterile

DYNJ089202	2" x 2.5 yds./5.08 cm x 2.27 m	16 rls/cs
DYNJ089204	4" x 2.5 yds./10.16 cm x 2.27 m	8 rls/cs
DYNJ089207	4" x 7.5 yds./10.16 cm x 6.86 m	6 rls/cs

CoFlex AFD, Tan, Non-Sterile

MDS089202	2" x 2.5 yds./5.08 cm x 2.27 m	16 rls/cs
MDS089204	4" x 2.5 yds./10.16 cm x 2.27 m	8 rls/cs
MDS089207	4" x 7.5 yds./10.16 cm x 6.86 m	6 rls/cs

*This product is not manufactured or distributed by 3M company, the owner of the registered trademark Coban.



Self-Adherent Wrap

- » Self-adherent (cohesive) wrap sticks to itself but not to skin, hair or dressings
- » Secondary dressing that functions like a tape to retain dressings without tension around digits or limbs
- » Can be applied to provide light compression if desired



Item No.	Description	Pkg
----------	-------------	-----

Tan, Stretched, Latex

PRM086001	1" x 5-yd./2.54 cm x 4.57 m	30 rls/cs
PRM086002	2" x 5-yd./5.08 cm x 4.57 m	36 rls/cs
PRM086003	3" x 5-yd./7.62 cm x 4.57 m	24 rls/cs
PRM086004	4" x 5-yd./10.16 cm x 4.57 m	18 rls/cs
PRM086006	6" x 5-yd./15.24 cm x 4.57 m	12 rls/cs

Item No.	Description	Pkg
----------	-------------	-----

Tan, Stretched, Not Made with Natural Rubber Latex

PRM088001	1" x 5-yd./2.54 cm x 4.57 m	30 rls/cs
PRM088002	2" x 5-yd./5.08 cm x 4.57 m	36 rls/cs
PRM088003	3" x 5-yd./7.62 cm x 4.57 m	24 rls/cs
PRM088004	4" x 5-yd./10.16 cm x 4.57 m	18 rls/cs
PRM088006	6" x 5-yd./15.24 cm x 4.57 m	12 rls/cs

Color Pack, Not Made with Natural Rubber Latex

PRM088001CP	1" x 5-yd./2.54 cm x 4.57 m	30 rls/cs
PRM08802CP	2" x 5-yd./5.08 cm x 4.57 m	36 rls/cs
PRM08803CP	3" x 5-yd./7.62 cm x 4.57 m	24 rls/cs



Kid's Neon Colors Bandages

- » Latex-free adhesive bandages are designed to fit kid-sized wounds

Item No.	Description	Pkg
NON256131Z	0.75 x 3"/1.9 x 7.62 cm	50/bx



Medtoons® Cartoon Bandages

- » Various designs
- » Durable acrylic adhesive is skin friendly and absorbent
- » Non-adherent island pad
- » Convenient dispenser box
- » Not made with natural rubber latex

Item No.	Description	Pkg
NON256130Z	0.75 x 3"/1.9 x 7.62 cm	50/bx



Cupcake Covers Bandages

- » Colorful and kid-friendly
- » Not made with natural rubber latex


Item No.	Description	Pkg
CUR00002RB	20 Count/Assorted Sizes	24/cs



Busy Bugs Adhesive Bandages

- » Colorful and kid-friendly
- » Not made with natural rubber latex

Item No.	Description	Pkg
CUR00001RB	20 Count/Assorted Sizes	24/cs



CUR45702RB



CAMO Fabric Adhesive Bandages

- » Designed to fit kid-sized wounds
- » Not made with natural rubber latex

Item No.	Description	Pkg
CUR45701RB	Green Camo: 0.75 x 3"/ 1.9 x 7.62 cm	24/bx, 25 bx/cs, 600/cs
CUR45702RB	Pink and Blue Camo: 0.75 x 3"/ 1.9 x 7.62 cm	24/bx, 25 bx/cs, 600/cs



CUR45701RB



Plastic Sterile Adhesive Bandages

- » Constructed with PVC-free durable polyethylene and are perforated for breathability
- » Features an island wound pad allowing for a four-sided seal

Item No.	Description	Pkg
NON25509Z	0.375 x 1.5"/1 x 3.8 cm, Junior	100/bx
NON25500Z	0.75 x 3"/1.9 x 7.62 cm	100/bx
NON25501H	0.875"/2.2 cm, Spot	100/bx
NON25600Z	1 x 3"/2.54 x 7.62 cm	100/bx
NON25504Z	2 x 4"/5.08 x 10.16 cm, X-Large	50/bx



NON25650Z

Flex-Fabric™ Sterile Adhesive Bandages

- » Woven to allow them to stretch and conform to the wound, material naturally breathes
- » Features an island wound pad allowing for a four-sided seal

Item No.	Description	Pkg
NON25650Z	0.75 x 3"/1.9 x 7.62 cm	100/bx
NON25502Z	0.875"/2.2 cm, Spot	100/bx
NON25660Z	1 x 3"/2.54 x 7.62 cm	100/bx
NON25510Z	1.5 x 3"/3.8 x 7.62 cm, Knuckle	100/bx
NON25513Z	1.75 x 2"/4.44 x 5.08 cm, Fingertip	100/bx
NON25524Z	2 x 4"/5.08 x 10.16 cm	50/bx



NON25750

Comfort Fabric Sterile Adhesive Bandages

- » Constructed of a soft non-woven gentle material that is breathable, comfortable, and durable
- » Offers even adhesion that stretches with your skin
- » Features an island wound pad allowing for a four-sided seal
- » Not made with natural rubber latex

Item No.	Description	Pkg
NON25750	0.75 x 3"/1.9 x 7.62 cm	100/bx, 12 bx/cs, 1200/cs
NON25760	1 x 3"/2.54 x 7.62 cm	100/bx, 12 bx/cs, 1200/cs
NON25744	2 x 4"/5.08 x 10.16 cm, X-Large	50/bx, 12 bx/cs, 600/cs



PRM25600H



Adhesive Bandages, Sterile

- » Ideal for bandaging minor wounds and scrapes
- » Tan plastic backing is water-resistant
- » Long-lasting adhesive with minimal residue
- » Not made with natural rubber latex

Item No.	Description	Pkg
PRM25500H	0.75 x 3"/1.91 x 7.62 cm	100/bx
PRM25600H	1 x 3"/2.54 x 7.62 cm	100/bx



NON270001Z



Paper Tape Rolls

Compare to Micropore®*

- » A breathable tape with a gentle adhesion for sensitive skin, recommended especially for pediatric and geriatric patients
- » Single-use rolls are great for reducing waste and reducing cross contamination
- » Adheres well to damp skin
- » Not made with natural rubber latex

Item No.	Description	Pkg
NON270012Z	0.5" x 10-yd./ 1.25 cm x 9.14 m	24 rls/bx
NON270001Z	1" x 10-yd./ 2.54 cm x 9.14 m	12 rls/bx
NON270002Z	2" x 10-yd./ 5.08 cm x 9.14 m	6 rls/bx
NON270003Z	3" x 10-yd./ 7.62 cm x 9.14 m	4 rls/bx

With Dispenser

NON260001DZ	1" x 10-yd./ 2.54 cm x 9.14 m	12/bx
NON260002DZ	2" x 10-yd./ 5.08 cm x 9.14 m	6/bx

*This product is not manufactured or distributed by 3M Company, makers of Micropore®.

NON270201Z



Transparent Tape Rolls

Compare to Transpore®*

- » A breathable tape that permits skin examination without tape removal
- » Provides a clean and easy bilateral tear for custom widths
- » Single-use rolls are great for reducing waste and reducing cross contamination
- » Holds facial dressings
- » Not made with natural rubber latex

Item No.	Description	Pkg
NON270212Z	0.5" x 10-yd./ 1.25 cm x 9.14 m	24/bx
NON270201SZ	1" x 1½-yd./ 2.54 cm x 1.37 m	100/bx
NON270201Z	1" x 10-yd./ 2.54 cm x 9.14 m	12/bx
NON270202SZ	2" x 1.5-yd./ 5.08 cm x 4.57 m	50/bx
NON270202Z	2" x 10-yd./ 5.08 cm x 9.14 m	6/bx
NON270202SZ	2" x 1½-yd./ 5.08 cm x 1.37 m	50/bx
NON270203Z	3" x 10-yd./ 7.62 cm x 9.14 m	4/bx

*This product is not manufactured or distributed by 3M Company, makers of Transpore®.



Cloth "Silk" Tape Rolls

Compare to Durapore®*

- » Strong woven construction
- » Water resistant and easy to tear
- » Not made with natural rubber latex

Item No.	Description	Pkg
NON270112Z	0.5" x 10-yd./ 1.25 cm x 9.14 m	24/bx
NON270101Z	1" x 10-yd./ 2.54 cm x 9.14 m	12/bx
NON270102Z	2" x 10-yd./ 5.08 cm x 9.14 m	6/bx
NON270103Z	3" x 10-yd./ 7.62 cm x 9.14 m	4/bx

*This product is not manufactured or distributed by 3M Company, makers of Durapore®.



PRM260001



Paper Tape Rolls

- » Breathable, with a gentle adhesion for sensitive skin
- » Gentle, air-permeable, non-woven material
- » Especially recommended for geriatric patients or others with fragile skin
- » Conforms and adheres even to damp skin
- » Secures lightweight dressings and tubing, ostomy pouches and facial dressings
- » Not made with natural rubber latex

Item No.	Description	Pkg
PRM260001H	1" x 10-yd./ 2.54 cm x 9.14 m	12 rls/bx
PRM260001HH	1" x 10-yd./ 2.54 cm x 9.14 m	1/ea
PRM260002H	2" x 10-yd./ 5.08 cm x 9.14 m	6 rls/bx
PRM260003H	3" x 10-yd./ 7.62 cm x 9.14 m	4 rls/bx



PRM260201H



Transparent Tape Rolls

- » Permits skin examination without tape removal
- » Clear, breathable plastic material
- » Micro-perforation allows for bidirectional tearing
- » Excellent for IV sites and tubing retention
- » Not made with natural rubber latex

Item No.	Description	Pkg
PRM260201H	1" x 10-yd./ 2.54 cm x 9.14 m	12 rls/bx
PRM260202H	2" x 10-yd./ 5.08 cm x 9.14 m	6 rls/bx



Cloth "Silk" Tape Rolls

- » Sturdy rayon/acetate fabric helps secure dressings in place
- » Serrated edges facilitate tearing
- » Durable adhesion with minimal residue on skin
- » Water repellent
- » Not made with natural rubber latex

Item No.	Description	Pkg
PRM260101H	1" x 10-yd./ 2.54 cm x 9.14 m	12 rls/bx
PRM260102H	2" x 10-yd./ 5.08 cm x 9.14 m	6 rls/bx
PRM260103	3" x 10-yd./ 7.62 cm x 9.14 m	4 rls/bx, 12 bx/cs, 48 rls/cs



NON260302Z

Ortho-Porous Sports Adhesive Tape Rolls

- » Strong woven cotton backing provides high tensile strength and support to highly active patients
- » Powerful adhesive is designed for long-lasting adhesion with minimal residue
- » Highly porous construction allows the skin to breathe and permits moisture vapor transmission
- » Not made with natural rubber latex

Item No.	Description	Pkg
NON260301H	1" x 10-yd./ 2.5 cm x 9.1 m	1/ea
NON260315H	1.5" x 10-yd./ 3.8 cm x 9.1 m	1/ea
NON260302H	2" x 10-yd./ 5.08 cm x 9.1 m	1/ea
NON260302Z	2" x 10-yd./ 5.08 cm x 9.14 m	6/bx
NON260303HH	3" x 10-yd./ 7.5 cm x 9.1 m	1/ea



Medi-Strips® Reinforced Wound Closure Strips

Compare to 3M® Steri-Strip®*

- » Adhesive wound closures with fiber backing use reinforced filaments for high tensile strength
- » Medical grade adhesive affords excellent adhesion while being gentle to the skin
- » Use in addition to, or in place of, sutures
- » Each size is in a color-coded package for quick and easy size identification
- » Not made with natural rubber latex



Item No.	Description	Pkg
NON250318Z	0.125 x 3"/0.32 x 7.62 cm	5/pk, 50 pk/bx, 250/bx
NON250114Z	0.25 x 1.5"/0.6 x 3.8 cm	6/pk, 50 pk/bx, 300/bx
NON250314Z	0.25 x 3"/0.6 x 7.62 cm	3/pk, 50 pk/bx, 50/bx
NON250414Z	0.25 x 4"/0.6 x 10.16 cm	10/pk, 50 pk/bx, 500/bx
NON250412Z	0.5 x 4"/1.3 x 10.16 cm	6/pk, 50 pk/bx, 300/bx
NON250501Z	1 x 5"/2.54 x 12.7 cm	4/pk, 25 pk/bx, 300/bx

*This product is not manufactured or distributed by 3M Company, makers of Steri-Strip®.




Fiberglass Casting Tape

- » Easy application for patient and staff convenience
- » Dries to a smooth finish; gentle against exposed skin
- » Fresh fragrance and color choice enhance patient experience
- » Conforms well for a comfortable fit
- » Allows for penetration by X-rays to monitor bone healing

Item No.	Color	Size	Pkg
CT24BLK	Black	2" x 4-yd.	10 rolls/bx
CT34BLK		3" x 4-yd.	10 rolls/bx
CT44BLK		4" x 4-yd.	10 rolls/bx
CT24BLU	Blue	2" x 4-yd.	10 rolls/bx
CT34BLU		3" x 4-yd.	10 rolls/bx
CT44BLU		4" x 4-yd.	10 rolls/bx
CT24WHT	White	2" x 4-yd.	10 rolls/bx
CT34WHT		3" x 4-yd.	10 rolls/bx
CT44WHT		4" x 4-yd.	10 rolls/bx
CT24PNK	Pink	2" x 4-yd.	10 rolls/bx
CT34PNK		3" x 4-yd.	10 rolls/bx
CT44PNK		4" x 4-yd.	10 rolls/bx
CT24PRP	Purple	2" x 4-yd.	10 rolls/bx
CT34PRP		3" x 4-yd.	10 rolls/bx
CT44PRP		4" x 4-yd.	10 rolls/bx





Durable and MRI-Safe Arm Boards

- » Available as reusable or disposable
- » Soft, spongy core for patient comfort and a firm covering for stability
- » Made of durable vinyl

Item No.	Description	Pkg
Reusable		
NON24285	Pediatric, with Hook-and-Loop, 3 x 9"/7.6 x 22.9 cm	50/cs
NON24286	Infant, with Hook-and-Loop, 2 x 6"/5.1 x 15.2 cm	50/cs
NON24292	Adult, 3.5 x 9"/8.9 x 22.9 cm	24/cs
NON24290	Adult, 3.5 x 18"/8.9 x 45.7 cm	24/cs

Disposable

NON24288	IV, Infant, 1 x 4"/2.5 x 10.2 cm	50/cs
NON24288H	IV, Infant, 1 x 4"/2.5 x 10.2 cm	1/ea
NON24289	IV, Infant, 2 x 4"/5.1 x 10.2 cm	50/cs
NON24289H	IV, Infant, 2 x 4"/5.1 x 10.2 cm	1/ea
NON24297	IV, Infant, 2 x 6"/5.1 x 15.2 cm	50/cs
NON24297H	IV, Infant, 2 x 6"/5.1 x 15.2 cm	1/ea
NON24294	IV, Infant, 2 x 9"/5.1 x 22.9 cm	50/cs
NON24293	IV, Infant, 3 x 9"/7.6 x 22.9 cm	50/cs
NON24293H	IV, Infant, 3 x 9"/7.6 x 22.9 cm	1/ea
NON24291	IV, Adult, 3 x 17.5"/7.6 x 44.5 cm	100/cs
NON24291H	IV, Adult, 3 x 17.5"/7.6 x 44.5 cm	1/ea

View Compression Stockings in the Patient Care section on page 114, or visit www.medline.com/physician for more information.



CURSALINE7



Saline Wound Wash

- » Gently cleanses minor scrapes and wounds
- » Pressurized flow effectively removes dirt and debris with better control and less waste
- » Dispenses a controlled spray at any angle, even upside down

Item No.	Description	Pkg
CURSALINE7	CURAD Sterile Saline Wound Wash, 7.1 fl. oz.	12/cs



Emesis Bags

- » Tie bag off at the ring to reduce contamination and smell
- » More private, more discreet for patient dignity
- » Contaminants and pathogens are isolated



NON80329



NON80328

Item No.	Description	Pkg
Bag, Emesis, Sickness Clean-Up, Clean Sack		
NON80328	Fits in Dispenser, NON80328D	144/cs
NON80328Z	Fits in Dispenser, NON80328D	24/pk
NON80329	Fits in Dispenser, NON80328D	144/cs
NON80329Z	Fits in Dispenser, NON80328D	24/pk

Bag, Vomit, Emesis, Sickness, with Gradations		
NON70600	Individually Packaged	500/cs
NON70600H	Individually Packaged	1/ea

Dispenser, Sickness Bag		
NONEMBGDISP	Fits Items NON80328 and NON80329	6/cs
NONEMBGDISPH	Fits Items NON80328 and NON80329	1 ea



PGL050 PGL400 PGL050TUBE

PluroGel Burn and Wound Dressing

- » 100% water-soluble, biocompatible, cell-friendly, and non-ionic
- » Softer wound debris rinses off easily at dressing change
- » Helps to create and maintain an optimal moist wound environment
- » Thickens in warm environment and remains on the wound bed for better protection

Item No.	Description	Pkg
PGL050	PluroGel, 50 g jar	24/cs
PGL050H	PluroGel, 50 g jar	1 ea
PGL050TUBE	PluroGel, 50 g tube	12/cs
PGL050TUBEH	PluroGel, 50 g tube	1 ea
PGL400	PluroGel, 400 g jar	6/cs
PGL400H	PluroGel, 400 g jar	1 ea



MSS4048 MSS4044 MSS4022

Hyalomatrix Hyaluronic Acid Wound Device

- » 100% biocompatible HA is not derived from human or animal tissue (no immunorejection)
- » Primary layer, composed entirely of HYAFF 11, acts as a scaffold for cellular colonization and capillary growth
- » Facilitates ordered reconstruction of dermal tissue
- » Semipermeable silicone membrane controls MVTR (moisture vapor transmission rate)

Item No.	Description	Pkg
MSS4011	Hyalomatrix, 2.5 x 2.5 cm	10/bx
MSS4011H	Hyalomatrix, 2.5 x 2.5 cm	1 ea
MSS4022	Hyalomatrix, 5 x 5 cm	10/bx
MSS4022H	Hyalomatrix, 5 x 5 cm	1 ea
MSS4044	Hyalomatrix, 10 x 10 cm	1 ea
MSS4048	Hyalomatrix, 10 x 20 cm	1 ea



Optifoam Gentle EX

- » Five layer silicone bordered foam dressing with superabsorbent core


Item No.	Description	Pkg
MSCEX33EP	3 x 3"/7.6 x 7.6 cm, Pad: 1.75 x 1.75"/3.8 x 3.8 cm	10/bx, 10 bx/cs
MSCEX44EP	4 x 4"/10.2 x 10.2 cm, Pad: 2.5 x 2.5"/6.4 x 6.4 cm	10/bx, 10 bx/cs
MSCEX66EP	6 x 6"/15.2 x 15.2 cm, Pad: 4.5 x 4.5"/11.4 x 11.4 cm	10/bx, 10 bx/cs
MSCEX77EP	Sacral: 7 x 7"/17.8 x 17.8 cm	10/bx, 10 bx/cs
MSCEX99EP	Sacral: 10 x 9"/25.4 x 22.9 cm	10/bx, 10 bx/cs



Optifoam Gentle

- » Silicone bordered foam dressing

Item No.	Description	Pkg
MSC20162EP	1.6 x 2"/4 x 5 cm, Pad: 0.8 x 1.2"/2.0 x 3.0 cm	10/bx, 7 bx/cs
MSC2033EP	3 x 3"/7.6 x 7.6 cm, Pad: 1.5 x 1.5"/3.8 x 3.8 cm	10/bx, 10 bx/cs
MSC2044EP	4 x 4"/10.2 x 10.2 cm, Pad: 2.5 x 2.5"/6.4 x 6.4 cm	10/bx, 10 bx/cs
MSC2066EP	6 x 6"/15.2 x 15.2 cm, Pad: 4.5 x 4.5"/11.4 x 11.4 cm	10/bx, 10 bx/cs
MSC2077EP	Sacral: 7 x 7"/17.8 x 17.8 cm	5/bx, 10 bx/cs
MSC20109EP	Sacral: 10 x 9"/25.4 x 22.9 cm	5/bx, 10 bx/cs



Optifoam Gentle Non-Bordered

- » Non-bordered silicone faced foam dressing

Item No.	Description	Pkg
MSC2244EP	4 x 4"/10.2 x 10.2 cm	10/bx, 10 bx/cs
MSC2266EP	6 x 6"/15.2 x 15.2 cm	10/bx, 10 bx/cs
MSC2288EP	8 x 8"/20.4 x 20.4 cm	5/bx, 10 bx/cs

Ordering Information

To order by the box, add "Z" to the end of the item number.
To order by the dressing, add "H" to the end of the item number.



Optifoam Adhesive

» Adhesive bordered foam dressing

Item No.	Description	Pkg
MSC1044EP	4 x 4"/10.2 x 10.2 cm, Pad: 2.5 x 2.5"/6.4 x 6.4 cm	10/bx, 10 bx/cs
MSC1066EP	6 x 6"/15.2 x 15.2 cm, Pad: 4.5 x 4.5"/11.4 x 11.4 cm	10/bx, 10 bx/cs
MSC1065EP	Sacrum: 6 x 5.5"/15.5 x 14.2 cm	10/bx, 10 bx/cs



Optifoam Non-Adhesive

» Non-bordered foam dressing

Item No.	Description	Pkg
MSC1244EP	4 x 4"/10.2 x 10.2 cm	10/bx, 10 bx/cs
MSC1266EP	6 x 6"/15.2 x 15.2 cm	10/bx, 10 bx/cs



Optifoam Heel

» Heel-shaped non-adhesive foam dressing

Item No.	Description	Pkg
MSC1200EP	Heel Shaped	5/bx, 8 bx/cs



Optifoam Basic

» Foam dressing without polyurethane layer

Item No.	Description	Pkg
MSC1133	3 x 3"/7.6 x 7.6 cm	10/bx, 10 bx/cs
MSC1133F	3 x 3"/7.6 x 7.6 cm with Fenestration	10/bx, 10 bx/cs
MSC1145	4 x 5"/10.2 x 12.7 cm	10/bx, 10 bx/cs



Optifoam Site

» Circular site dressing with radial slit and starburst opening


Item No.	Description	Pkg
MSC1104	Adhesive Dressing: 4"/10.2 cm, Pad: 2"	30/bag, 4 bags/cs
MSC1102	Non-Adhesive Pad: 2"/5.1 cm	10/bx, 10 bx/cs



Optifoam Thin

» Thin adhesive faced dressing

Item No.	Description	Pkg
MSC1523EP	2 x 3"/5.1 x 7.6 cm	10/bx, 10 bx/cs
MSC1544EP	4 x 4"/10.2 x 10.2 cm	10/bx, 10 bx/cs



Optifoam Gentle Ag+ Silicone Faced and Bordered

» Antimicrobial foam dressing with superabsorbent core

Item No.	Description	Pkg
MSC9644EP	4 x 4"/10.2 x 10.2 cm Pad: 2.5 x 2.5"/6.4 x 6.4 cm	10/bx, 10 bx/cs
MSC9666EP	6 x 6"/15.2 x 15.2 cm Pad: 4.5 x 4.5"/11.4 x 11.4 cm	10/bx, 10 bx/cs



Optifoam Gentle Ag+ Non-Bordered Silicone Faced

» Antimicrobial non-bordered silicone faced foam dressing

Item No.	Description	Pkg
MSC9544EP	4 x 4"/10.2 x 10.2 cm	10/bx, 10 bx/cs
MSC9566EP	6 x 6"/15.2 x 15.2 cm	10/bx, 10 bx/cs
MSC9588EP	8 x 8"/20.3 x 20.3 cm	10/bx, 10 bx/cs



Optifoam Gentle Ag+ Silicone Bordered

» Antimicrobial silicone bordered foam dressing



Item No.	Description	Pkg
MSC9744EP	4 x 4"/10.2 x 10.2 cm, Pad: 2.5 x 2.5"/6.4 x 6.4 cm	10/bx, 10 bx/cs
MSC9766EP	6 x 6"/15.2 x 15.2 cm, Pad: 4.5 x 4.5"/11.4 x 11.4 cm	10/bx, 10 bx/cs



Optifoam Ag+ Adhesive

» Antimicrobial silver bordered foam



Item No.	Description	Pkg
MSC9604EP	4 x 4"/10.2 x 10.2 cm, Pad: 2.5 x 2.5"/6.4 x 6.4 cm	10/bx, 10 bx/cs
MSC9606EP	Sacrum: 6 x 5.5"/15.2 x 14.2 cm	10/bx, 10 bx/cs



Optifoam Ag+ Non-Adhesive

» Antimicrobial non-bordered foam



Item No.	Description	Pkg
MSC9614EP	4 x 4"/10.2 x 10.2 cm	10/bx, 10 bx/cs



Optifoam Gentle Ag+ SuperAbsorbent Post-Op

» Antimicrobial silicone bordered and faced post operative foam dressing with superabsorbent core and flexible design



Item No.	Description	Pkg.
MSC9746	4 x 6"/10.2 x 15.2 cm, Pad: 2 x 4"/5.08 x 10.16 cm	10/bx, 10 bx/cs
MSC9748	4 x 8"/10.2 x 20.3 cm, Pad: 2 x 6"/5.1 x 15.2 cm	10/bx, 10 bx/cs
MSC97410	4 x 10"/10.2 x 25.4 cm, Pad: 2 x 8"/5.1 x 20.3 cm	10/bx, 10 bx/cs
MSC97412	4 x 12"/10.2 x 30.5 cm, Pad: 2 x 10"/5.1 x 25.4 cm	10/bx, 10 bx/cs
MSC97414	4 x 14"/10.2 x 35.6 cm, Pad: 2 x 12"/5.1 x 30.5 cm	10/bx, 10 bx/cs

Ordering Information

To order by the box, add "Z" to the end of the item number.

To order by the dressing, add "H" to the end of the item number.



Optifoam Gentle Ag+ Post-Op

» Antimicrobial silicone bordered post-operative foam dressing



Item No.	Description	Pkg.
MSC9736	3.5 x 6"/8.9 x 15.2 cm, Pad: 1.5 x 4"/3.8 x 10.2 cm	10/bx, 10 bx/cs
MSC97310	3.5 x 10"/8.9 x 25.4 cm, Pad: 1.5 x 8"/3.8 x 20.3 cm	10/bx, 10 bx/cs
MSC97314	3.5 x 14"/8.9 x 35.6 cm, Pad: 1.5 x 12"/3.8 x 30.5 cm	10/bx, 10 bx/cs



Optifoam Ag+ Post-Op

» Antimicrobial silver foam post-operative dressing




Item No.	Description	Pkg.
MSC9636	3.5 x 6"/8.9 x 15.2 cm, Pad: 1.5 x 4"/3.8 x 10.2 cm	10/bx, 10 bx/cs
MSC96310	3.5 x 10"/8.9 x 25.4 cm, Pad: 1.5 x 8"/3.8 x 20.3 cm	10/bx, 10 bx/cs
MSC96314	3.5 x 14"/8.9 x 35.6 cm, Pad: 1.5 x 12"/3.8 x 30.5 cm	10/bx, 10 bx/cs



Quick Superabsorbent Aquaconductive™ Dressing

» Superabsorbent conformable dressing



Item No.	Description	Pkg.
MSC5822	2 x 2"/5.1 x 5.1 cm	10/bx, 5 bx/cs
MSC5844	4.25 x 4"/10.8 x 10.2 cm	10/bx, 5 bx/cs
MSC5868	6.125 x 8"/15.56 x 20.321 cm	10/bx, 5 bx/cs



Optilock Superabsorbent Dressing

» Non-adherent and superabsorbent dressing



Item No.	Description	Pkg.
MSC6433EP	3 x 3"/7.5 x 7.5 cm	10/bx, 10 bx/cs
MSC6444EP	4 x 4"/10.2 x 10.2 cm	10/bx, 10 bx/cs
MSC6455EP	5 x 5.5"/12.7 x 14 cm	10/bx, 10 bx/cs
MSC64610EP	6.5 x 10"/16.5 x 25.4 cm	10/bx, 5 bx/cs
MSC64812EP	8 x 12"/20.3 x 30.5 cm	10/bx, 8 bx/cs



Opticell® AG+ Gelling Fiber with Silver

- » Unique gelling fiber ideal for dry to heavily draining wounds
- » Maintains optimal moist wound healing environment
- » Broad-spectrum antibacterial silver

Item No.	Description	Pkg
MSC9822EPZ	2 x 2"/5.1 x 5.1 cm	10/bx
MSC9845EPZ	4 x 5"/10.2 x 12.7 cm	10/bx
MSC9866EPZ	6 x 6"/15.2 x 15.2 cm	5/bx
MSC98812EPZ	8 x 12"/20.3 x 15.2 cm	5/bx
MSC9818RZ	1 x 18"/2.5 x 46 cm	5/bx



Maxorb® II AG+ Dressings

- » Ideal for moderate to heavily draining wounds
- » Maintains optimal moist wound healing environment
- » Broad-spectrum antibacterial silver
- » 21-day extended wear time

Item No.	Description	Pkg
MSC9922EPZ	2 x 2"/5.1 x 5.1 cm	10/bx
MSC9945EPZ	4 x 4 3/4"/10.2 x 12.1 cm	10/bx
MSC9948EPZ	4 x 8"/10.2 x 20.3 cm	5/bx
MSC9966EPZ	6 x 6"/15.2 x 15.2 cm	5/bx
MSC99812EPZ	8 x 12"/20.3 x 30.5 cm	5/bx
MSC9912EPZ	1 x 12"/2.5 x 46 cm	5/bx



Maxorb II Alginate Dressings

- » Ideal for moderate to heavily draining wounds
- » 100% calcium alginate designed for superior absorbency and wet strength
- » Maintains optimal moist wound healing environment

Item No.	Description	Pkg
MSC7322EPZ	2 x 2"/5 x 5 cm	10/bx
MSC7344EPZ	4 x 4"/10 x 10 cm	10/bx
MSC7348EPZ	4 x 8"/10 x 20 cm	5/bx
MSC7366EPZ	6 x 6"/15 x 15 cm	5/bx
MSC7312EPZ	1 x 12"/2.5 x 30 cm	10/bx
MSC7318EPZ	1 x 18"/2.5 x 46 cm	10/bx



NONUNNAS130

Unna-Z® Stretch Unna Boots

- » Maintains moist healing environment and provides light compression for venous leg ulcers
- » Impregnated with zinc oxide or zinc oxide and calamine
- » Inner plastic core for easy application

Item No.	Description	Pkg
Unna Boot with Zinc and Calamine		
NONUNNAS30	3" x 10-yd./ 7.62 cm x 9.14 m	12/cs
NONUNNAS40	4" x 10-yd./ 10.16 cm x 9.14 m	12/cs
Unna Boot with Zinc		
NONUNNAS130	3" x 10-yd./ 7.62 cm x 9.14 m	12/cs
NONUNNAS140	4" x 10-yd./ 10.16 cm x 9.14 m	12/cs



MSC4300



Threeflex® Compression Bandage System

- » Lighter three-layer compression for chronic venous insufficiency and associated edema
- » Effective appropriate compression
- » Extended wear time

Item No.	Description	Pkg
MSC4300	Threeflex Kit	8 kits/cs



MSC4400

Fourflex® Compression Bandage System

- » For the treatment of chronic venous insufficiency and associated edema
- » Effective appropriate compression
- » Extended wear time

Item No.	Description	Pkg
MSC4400	Fourflex Kit	8 kits/cs
MSC4400XL	Fourflex XL Kit	8 kits/cs



MSC9500



Medigrip™ Tubular Bandages

- » Provides excellent support for joints
- » Can provide mild compression when doubled
- » Wide range of applications

Item No.	Description	Pkg
MSC9500	Size A, 1 3/4" x 11-yd./ 4.5 cm x 10 m	1 roll/bx
MSC9501	Size B, 2 1/2" x 11-yd./ 6.3 cm x 10 m	1 roll/bx
MSC9502	Size C, 2 5/8" x 11-yd./ 6.8 cm x 10 m	1 roll/bx
MSC9503	Size D, 3" x 11-yd./ 7.5 cm x 10 m	1 roll/bx
MSC9504	Size E, 3 1/2" x 11-yd./ 8.75 cm x 10 m	1 roll/bx
MSC9505	Size F, 4" x 11-yd./ 10 cm x 10 m	1 roll/bx
MSC9506	Size G, 4 3/4" x 11-yd./ 12 cm x 10 m	1 roll/bx
MSC9507	Size J, 6 3/4" x 11-yd./ 17.15 cm x 10 m	1 roll/bx
MSC9508	Size K, 8 1/4" x 11-yd./ 20.96 cm x 10 m	1 roll/bx
MSC9504YDH	3 1/2" x 1-yd./ 8.75 cm x 0.91 m	1/bx, 30 bx/cs
MSC9505YDH	4" x 1-yd./ 10 cm x 0.91 m	1/bx, 30 bx/cs
MSC9504YDH	4 3/4" x 1-yd./ 12 cm x 0.91 m	1/bx, 30 bx/cs



Ag+ Powder

- » Silver powder for hard to dress wounds
- » Provides continuous antimicrobial protection
- » Can be used for up to five days



Item No.	Description	Pkg
MSC9405Z	Bottle: 5 g	5/bx
MSC9410Z	Bottle: 10 g	5/bx



MSC903025EP



SilvaSorb® Antimicrobial Wound Gel

- » For lightly draining wounds in need of an antimicrobial barrier
- » Helps manage bacterial burden and provides continuous antimicrobial protection
- » Non-staining, gentle for the patient and extended wear time

Item No.	Description	Pkg
MSC93025EP	Tube: 0.25-oz./7.09 g	25/bx
MSC9301EP	Tube: 1.5-oz./42.52 g	12/cs
MSC9303EP	Tube: 3-oz./85.05 g	12/cs
MSC9308EP	Tube: 8-oz./226.8 g	6/cs



MSC5444



Exuderm® Satin Hydrocolloid Dressing

- » Satin finish backing designed to reduce friction and roll-up
- » Ideal for low to moderately draining wounds

Item No.	Description	Pkg
Square		
MSC5422	2 x 2"/5.08 x 5.08 cm	20/bx
MSC5444	4 x 4"/10.16 x 10.16 cm	10/bx
MSC5466	6 x 6"/15.2 x 15.2 cm	5/bx
MSC5488	8 x 8"/20.3 x 20.3 cm	5/bx

Sacral

MSC5470	4 x 3 3/5"/10.2 x 9.1 cm	10/bx
MSC5475	6 2/5 x 6 1/2"/16.3 x 16.5 cm	5/bx



MSC8622EPZ



Puracol® Plus Collagen Dressing

- » Ideal for wounds that are chronic or stalled
- » 100% native collagen designed to help manage stalled chronic wounds
- » Features unique MicroScaffold™ structure for longer wear time (up to 7 days)

Item No.	Description	Pkg
MSC8622EPZ	2 x 2 1/2"/5.1 x 5.7 cm	10/bx
MSC8644EPZ	4 1/5 x 4 1/2"/ 10.8 x 11.4 cm	10/bx
MSC861X8EPZ	Rope: 1 x 8'/ 2.5 x 20.3 cm	10/bx



MSC8722EPZ



Puracol® Plus AG+

- » Ideal for stalled wounds when antimicrobial properties of silver are desired
- » 100% native collagen impregnated with ionic silver designed to help jump-start stalled chronic wounds and manage bioburden
- » Helps promote natural healing

Item No.	Description	Pkg
MSC8722EPZ	2 x 2"/5.1 x 5.1 cm	10/bx
MSC8744EPZ	4 x 4"/10.2 x 10.2 cm	10/bx
MSC871X8EPZ	Rope: 1 x 8'/ 2.5 x 20.3 cm	10/bx
MSC8488	8 x 8'/ 20.3 x 20.3 cm	50/cs



TheraHoney® Honey Dressings

- » Made of 100% medical grade honey (*Leptospermum*) with no additives or fillers
- » Use when moist wound healing is indicated
- » Promotes autolytic debridement, pair with secondary absorbent dressing

Item No.	Description	Pkg
MNK0005	Gel Tube: 0.5-oz./14.2 g	10/bx
MNK0015	Gel Tube: 1.5-oz./42.5 g	12/cs
MNK0077	Sheet: 4 x 5"/10.2 x 12.7 cm	10/bx
MNK0082	HD Sheet: 2 x 2"/5.1 x 5.1 cm	10/bx
MNK0087	HD Sheet: 4 x 5"/10.2 x 12.7 cm	10/bx
MNK0089	HD Ribbon: 1 x 12"/1.9 x 30.5 cm	10/bx
MNK1344	Foam Flex: 4 x 4"/10.2 x 10.2 cm	10/bx



MSC1500Z

Sureprep®

- » Fast-drying and economical skin protectant

Item No.	Description	Pkg
MSC1500Z	Sureprep Skin Protectant	50/bx



MSC1510Z

MSC1513Z

MSC1505Z

Sureprep® No-Sting

- » Alcohol-free, no-sting skin protectant
- » Ideal for damaged or fragile skin to protect from adhesive stripping or moisture-related damage

Item No.	Description	Pkg
MSC1505Z	Wipe	50/bx
MSC1510Z	Wand: 1 ml	25/bx
MSC1513Z	Wand: 3 ml	25/bx



MSC093005

Marathon® No-Sting Skin Protectant

- » Robust, flexible and long-lasting skin protection against friction and moisture
- » Clinically proven to protect against moisture-associated skin damage (MASD)*
- » Can be used to treat skin tears and minor abrasions

Item No.	Description	Pkg
MSC093005	0.02-oz./0.5 g Applicator	10/bx
MSC093001	0.02-oz./0.5 g Applicator	5/bx

*Independent lab test results on file.



MSC1823EPZ

Versatel® One Contact Layer Dressing

- » One-sided contact layer coated with silicone
- » Allows easy application by preventing gloves and scissors from sticking to the dressing
- » Adheres gently to skin but won't stick to moist wound beds, helping minimize pain during removal

Item No.	Description	Pkg
MSC1823EPZ	2 x 3"/5.1 x 7.6 cm	10/bx
MSC1834EPZ	3 x 4"/7.6 x 10.2 cm	10/bx
MSC1845EPZ	4 x 5"/10.2 x 12.7 cm	10/bx
MSC1847EPZ	4 x 7"/10.2 x 17.8 cm	10/bx
MSC18812EPZ	8 x 12"/20.3 x 30.5 cm	5/bx



MSC4004

Medfix™ Dressing Retention Tape

- » Low sensitivity adhesive, gentle for the patient
- » Flexible design with customized sizing
- » Water resistant

Item No.	Description	Pkg
MSC4002	2" x 11-yd./5.1 x 10 m	1/bx
MSC4004	3" x 11-yd./7.6 x 10 m	1/bx
MSC4006	6" x 11-yd./15.2 x 10 m	1/bx



Medline has a long history of rapid innovation in wound care and an extensive product line.

Don't see what you need?
Contact us to learn more.



MSC2302

Suresite® Window

» Transparent film with traditional window frame application method

Item No.	Description	Pkg
MSC2302	2.4 x 2.8"/6 x 7 cm	100/bx
MSC2304	4 x 4.5"/10.2 x 11.4 cm	50/bx



MSC2701Z

Suresite 123 Transparent Film Dressings

» Transparent film dressing with easy 123 application method

Item No.	Description	Pkg
MSC2701Z	1.5 x 1.5"/3.9 x 3.9 cm	100/bx
MSC2703Z	2.4 x 2.8"/6.1 x 7.1 cm	100/bx
MSC2705Z	4 x 4.8"/10.2 x 12.2 cm	50/bx
MSC2706ZZ	6 x 8"/15.2 x 20.3 cm	25/bx
MSC2710ZZ	4 x 10"/10.2 x 25.4 cm	25/bx
MSC2712ZZ	8 x 12"/20.3 x 30.5 cm	25/bx



MSC6008H

Skintegrity® Wound Cleanser

» Gently cleanses wounds before and between dressing changes
 » Specially formulated to allow fast, thorough and gentle cleansing of wounds in all stages
 » Mist or stream delivery

Item No.	Description	Pkg
MSC6008H	Spray Bottle: 8-oz.	1/ea
MSC6016H	Spray Bottle: 16-oz.	1/ea



Bordered Gauze

» Non-woven adhesive border with absorbent pad
 » Water resistant backing

Item No.	Description	Pkg
MSC3222Z	Dressing: 2 x 2"/5.1 x 5.1 cm, Pad: 1 x 1" /2.6 x 2.6 cm	15/bx
MSC3236Z	Dressing: 3 x 6"/7.6 x 15.2 cm Pad: 1.5 x 4"/3.8 x 10.1 cm	15/bx
MSC3244Z	Dressing: 4 x 4"/10.2 x 10.2 cm, Pad: 2 ½ x 2 ½"/6.4 x 6.4 cm	15/bx
MSC3245Z	Dressing: 4 x 5"/10.2 x 12.7 cm, Pad: 2 x 2 ½"/5.1 x 6.4 cm	15/bx
MSC3248Z	Dressing: 4 x 8"/10.2 x 20.3 cm, Pad: 2 x 6"/5.1 x 15.2 cm	15/bx
MSC32410Z	Dressing: 4 x 10"/10.2 x 25.4 cm, Pad: 2 x 8"/5.1 x 20.3 cm	15/bx
MSC32414Z	Dressing: 4 x 14"/10.2 x 35.6 cm, Pad: 2 x 12"/5.1 x 30.5 cm	15/bx
MSC3266Z	Dressing: 6 x 6"/15.2 x 15.2 cm, Pad: 4 x 4"/10.2 x 10.2 cm	15/bx



Stratasorb® Composite Dressings

» Waterproof composite dressing with absorbent pad
 » Four distinct layers designed to cover all stages of wound care

Item No.	Description	Pkg
MSC30410Z	Dressing: 4 x 10"/10.2 x 25.4 cm, Pad: 2 x 8"/5.1 x 20.3 cm	10/bx
MSC30414Z	Dressing: 4 x 14"/10.2 x 35.6 cm, Pad: 2 x 12"/5.1 x 30.5 cm	10/bx
MSC3044Z	Dressing: 4 x 4"/10.2 x 10.2 cm, Pad: 2 ½ x 2"/6.4 x 5.1 cm	10/bx
MSC3066Z	Dressing: 6 x 6"/15.2 x 15.2 cm, Pad: 4 x 4"/10.2 x 10.2 cm	10/bx
MSC3068Z	Dressing: 6 x 7.5"/15.2 x 19.1 cm, Pad: 4 x 6"/10.2 x 15.2 cm	10/bx

CAPITAL EQUIPMENT

Storage	174
Cabinetry	174
Medical Carts	174-175
Patient Mobility	176-177
Patient Scales	179-181
Measuring Instruments	181-182
Sterilization Equipment	182-183
Electrosurgery Generators	184
Stools	185-186
Hamper Stands	186-187
Mayo Stands	187
Lighting	188-189
IV Poles	190
Exam Tables	190-193
Exam Tables Accessories	194
Treatment Tables	194
Chairs	195
Disposable Cubicle Curtains	196
Reusable Cubicle Curtains	196-197



Ultra Stackable and Hanging Bins

- » Quantum® is one of the industry-leading hanging and stacking plastic storage bin
- » Heavy-duty, industrial-grade polymer bins stack together or hang from Quantum louvered racks and panels, rail racks, storage cabinets and louvered carts for high-density storage
- » Quantum are corrosion-proof and unaffected by weak acids and alkalis
- » Dividers, lids and a wide range of labeling options are available
- » For clear options, contact your Medline sales representative



Item No.	Description	Pkg
MDRQUS210C	5.375 W x 4.125 D x 3"H/13.7 x 10.5 x 7.6 cm	24
MDRQUS220C	7.375 W x 4.125 D x 3"H/18.7 x 10.5 x 7.6 cm	24
MDRQUS224C	10.875 W x 4.125 D x 4"H/27.6 x 10.5 x 10.2 cm	12
MDRQUS230C	10.875 W x 5.5 D x 5"H/27.6 x 14 x 12.7 cm	12
MDRQUS234C	14.75 W x 5.5 D x 5"H/37.5 x 14 x 12.7 cm	12
MDRQUS235C	10.875 W x 11 D x 5"H/27.6 x 27.9 x 12.7 cm	6
MDRQUS239C	10.75 W x 8.25 D x 7"H/27.3 x 21 x 17.8 cm	6
MDRQUS240C	14.75 W x 8.25 D x 7"H/37.5 x 21 x 17.8 cm	12
MDRQUS250C	14.75 W x 16.5 D x 7"H/37.5 x 41.9 x 17.8 cm	6
MDRQUS245C	10.875 W x 16.5 D x 5"H/27.6 x 41.9 x 12.7 cm	6
MDRQUS260C	18 W x 11 D x 10"H/45.7 x 27.9 x 25.4 cm	4
MDRQUS265C	18 W x 8.25 D x 9"H/45.7 x 21 x 22.9 cm	6
MDRQUS270C	18 W x 16.5 D x 11"H/45.7 x 41.9 x 27.9 cm	3

ESYDC150



WarmRight® Blanker Warmers

- » Multiple zone warming technology, provides heat where it is needed
- » Glass door for easy viewing
- » Warm, safe zone heating technology
- » Hands-free push button on door

Item No.	Description	Pkg
ESYDC150	DC150, 1.5 cu. ft. (3-4 Blankets Recommended)	1/ea
ESYDC250	DC250, 2.5 cu. ft. (4-6 Blankets Recommended)	1/ea
ESYDC350	DC350, 3.5 cu. ft. (10-12 Blankets Recommended)	1/ea
ESYDC750	DC750, 7.5 cu. ft. (20-25 Blankets Recommended)	1/ea



ivNOW® Fluid Warmers

- » Quickly heats and maintains safe temperatures of intravenous and irrigation fluids
- » Modular configuration
- » Anodized aluminum heat plate
- » Lightweight, medical-grade plastic housing
- » Intelligent bag sensors
- » Dual temperature sensors
- » Heats to 40°C/104°F in 30 minutes or less
- » Accurate to +0°/-2°C (+0°/-3°F)



ESYIVNOW4

Item No.	Description	Pkg	Item No.	Description	Pkg
ESYIVNOW1	One Pod	1/ea	ESYIVNOW4	Four Pods	1/ea
ESYIVNOW2	Two Pods	1/ea	ESYIVNOW5	Five Pods	1/ea
ESYIVNOW3	Three Pods	1/ea	ESYIVNOW6	Six Pods	1/ea

MPH01WNCS2



Narcotics Cabinet

- » Double door locks helps provide tight security for medications
- » Durable zinc-plated steel construction



Item No.	Description	Pkg
MPH01WNCS2	Small, Double Doors, 12 W x 8 D x 15"H/30.5 x 20.3 x 38.1 cm with Two Locks	1/ea



Pediatric Cart

- » Available in lightweight aluminum and traditional steel
- » Aluminum Overall: 29 W x 24.5 D x 42"H/73.7 x 62.2 x 106.7 cm
- » Steel Overall: 29 W x 24.5 D x 43.5"H/73.7 x 62.2 x 110.5 cm
- » Drawer: 22 W x 16.5 D x 4"H/55.9 x 41.9 x 10.2 cm



MPH01WMLA9P

Item No.	Description	Pkg
MPH01WMLA9P	Aluminum Construction, 8-3", 1-6" Drawers	1/ea
MPH01WML9PED	Steel Construction, 9-3" Drawers	1/ea

California customers—For all items on this page, **WARNING:** Cancer and Reproductive Harm—www.P65Warnings.ca.gov

MPH01WML3



Three-Drawer Medical Carts

- » Aluminum Overall: 29 W x 24.5 D x 36"H/ 73.7 x 62.2 x 91.4 cm
- » Steel Overall: 29 W x 24.5 D x 37.5"H/ 73.7 x 62.2 x 95.3 cm
- » Drawer: 22 W x 16.5"D/ 55.9 x 41.9 cm

Item No.	Description	Pkg
MPH01WMLA1	Aluminum Construction, 1-3"/2.5-7.6 cm, 2-9"/5.1-22.9 cm Drawers	1/ea
MPH01WML3	Steel Construction, 1-3"/2.5-7.6 cm, 2-9"/5.1-22.9 cm Drawers	1/ea



MPH01WML4



Four-Drawer Medical Carts

- » Aluminum Overall: 29 W x 24.5 D x 36"H/ 73.7 x 62.2 x 91.4 cm
- » Steel Overall: 29 W x 24.5 D x 37.5"H/ 73.7 x 62.2 x 95.3 cm
- » Drawer: 22 W x 16.5"D/ 55.9 x 41.9 cm

Item No.	Description	Pkg
MPH01WMLA2	Aluminum Construction, 2-3"/5.1-7.6 cm, 1-6"/2.5-15.2 cm, 1-9"/2.5-22.9 cm Drawers	1/ea
MPH01WML4	Steel Construction, 3"/2.5-7.6 cm, 1-6"/2.5-15.2 cm, 2-3" Drawers	1/ea



MPH01WMLA5



Five-Drawer Medical Carts

- » Aluminum Overall: 29 W x 24.5 D x 42"H/ 73.7 x 62.2 x 106.7 cm
- » Steel Overall: 29 W x 24.5 D x 43.5"H/ 73.7 x 62.2 x 110.5 cm
- » Drawer: 22 W x 16.5"D/ 55.9 x 41.9 cm

Item No.	Description	Pkg
MPH01WMLA5	Aluminum Construction, 2-3"/5.1-7.6 cm, 2-6"/5.1-15.2 cm, 1-9"/2.5-22.9 cm Drawers	1/ea
MPH01WML5	Steel Construction, 3"/7.6-7.6 cm, 4-3" Drawers	1/ea



MPH01WML6



Six-Drawer Medical Carts

- » Aluminum Overall: 29 W x 24.5 D x 42"H/ 73.7 x 62.2 x 106.7 cm
- » Steel Overall: 29 W x 24.5 D x 43.5"H/ 73.7 x 62.2 x 110.5 cm
- » Drawer: 22 W x 16"D/55.9 x 40.6 cm

Item No.	Description	Pkg
MPH01WMLA9	Aluminum Construction, 4-3"/10.2-7.6 cm, 1-6"/2.5-15.2 cm, 1-9"/2.5-22.9 cm Drawers	1/ea
MPH01WML6	Steel Construction, 6-3" Drawers	1/ea

MPH01WML6



Traditional Steel

- » Antimicrobial* properties help inhibit the growth of microorganisms
- » Removable ABS plastic-molded top for protection and easy cleanup
- » Full-extension ball-bearing slides for smooth, quiet drawers
- » New dual push handles for ease while pushing
- » Available with three locking options; breakaway (lever) lock and key lock are standard or mechanical push-button lock for an additional charge
- » Premium 5-inch/12.7 cm easy-roll casters mount to our non-marking high-impact plastic stabilizer base
- » Pull-out side shelf (single)

Item No.	Description	Pkg
MPH01WML6	Steel Cart, 6 Drawer, 4-3"/10.2-7.6 cm, 1-6"/2.5-15.2 cm, 1-9"/2.5-22.9 cm Drawers	1/ea

*This product is treated with an antimicrobial agent to protect the product. This product does not protect users or others against disease-causing bacteria.



Medical Carts Catalog

See our complete line of crash, isolation, anesthesia, bedside/slim, treatment/procedure and medication carts. Medline is your source for professional medical carts that store and transport vital equipment, materials and medications.

California customers—For all items on this page, **WARNING:** Cancer and Reproductive Harm—www.P65Warnings.ca.gov



MDS806700



Extra-Wide Wheelchairs

- » Ideal multi-person use chair built to last long and carry heavier weight capacity
- » Threaded seat-screw inserts minimize stripping from upholstery changes
- » Padded and upholstered armrests are standard
- » Carbon steel frame with rust- and chip-resistant chrome plating
- » Easy-to-clean navy vinyl upholstery
- » Chart pocket on back
- » 500-lb./226.8 kg weight capacity

Item No.	Description	Pkg
MDS806700	Desk-Length Removable Arm, Detachable Footrest, 20 x 18"/50.8 x 45.7 cm	1/ea
MDS806800	Detachable Elevated Legrest, Seat 22 x 18"/45.7 x 40.6 cm	1/ea
MDS806900	Detachable Elevated Legrest, Seat 24 x 18"/45.7 x 40.6 cm	1/ea

Accessories

MDS85196	Anti-Fold Device	1/ea
MDS85197A	Overhead Anti-Theft Device for 18–22"/45.7–55.9 cm Models	1/ea
MDS85183	IV Pole	1/ea
WCA6100	Oxygen Holder	1/ea
MDS85190	Oxygen/IV Combo	1/ea
MDS85189	Anti-Tip Device	1/ea
MDS85192	22"/55.9 cm Pushbar	1/ea
WCA806991	Tinnerman Legrest Locks	1/ea
MDS851965SH	Anti-Fold, 24"/61 cm Models	1/ea

Full-length arms available, just add "FLA" to the item number.



MDS806300D




Excel® 2000 Wheelchairs

- » Ideal multi-person use chair built to last long and adjust to different heights
- » Padded and upholstered armrests and calf pads
- » Threaded seat-screw inserts reduce stripping from upholstery changes
- » Hemi height-adjustable axles raise and lower seat 2-inches/5.08 cm
- » Desk-length removable arm
- » Chart pocket on back
- » Easy-to-clean vinyl upholstery
- » Weighs: 38-lbs./17 kg; weight capacity: 300-lbs./135 kg
- » Warranty: Lifetime on the frame; one year on all parts

Item No.	Description	Pkg	Item No.	Description	Pkg
MDS806400	Detachable Footrest, Seat 20 x 16"/50.8 x 40.6 cm	1/ea	Accessories		
MDS806250D	Detachable Footrest, Seat 18 x 16"/45.7 x 40.6 cm	1/ea	WCA6100	Oxygen Holder	1/ea
MDS806300D	Detachable Elevated Legrest, Seat 18 x 16"/45.7 x 40.6 cm	1/ea	MDS85196	Anti-Fold Device	1/ea
			MDS85183	IV Pole	1/ea
			MDS85197A	Overhead Anti-Theft Device	1/ea
			MDS85190	Oxygen/IV Combo	1/ea
			MDS85189	Anti-Tip Device	1/ea
			WCA806991	Tinnerman Legrest Locks	1/ea

Full-length arms available by adding "FLA" to the end of the item number.



MDS51514-10




Basic Crutches

- » Designed to be budget friendly
- » 250-lb./113.4 kg weight capacity

Item No.	Description	Pkg
MDS50514-10	Tall, Patient Height 5'10"–6'6", Height Adjustment 52–60.5"/132.1–153.7 cm	10 pr/cs
MDS51514-10	Adult, Patient Height 5'2"–5'10", Height Adjustment 44.5–52.5"/113–133.4 cm	10 pr/cs
MDS52514-10	Youth, Patient Height 4'6"–5'2", Height Adjustment 36.5"–44.5"/92.7–113 cm	10 pr/cs

California customers—For all items on this page, **⚠️ WARNING: Cancer and Reproductive Harm—www.P65Warnings.ca.gov**



I-beam post

MDSV80534



Standard Crutches

- » The I-beam post helps reduce the wobbliness one sometimes feels when using crutches
- » 300-lb./136.1 kg weight capacity

Item No.	Description	Pkg
MDSV80534	Adult, Tall, Patient Height 5'10"-6'6", Height Adjustment 52.5-60.5"/133.4153.7 cm	8 pr/cs
MDSV80535	Adult, Medium, Patient Height 5'1"-5'10", Height Adjustment 44.5-52.5"/113-133.4 cm	8 pr/cs
MDSV80536	Youth, Patient Height 4'6"-5'2", Height Adjustment 36.5-44.5"/92.7-113 cm	8 pr/cs
MDS80337Z	Economy, 250-lb./113 kg Weight Capacity, Patient Height 4'0"-4'6", Height Adjustment 31-37"/78-94 cm	2 pr/cs
G53314-8	Child, 275-lb./125 kg Weight Capacity, Patient Height 3'9"-4'7", Height Adjustment 31.5-37.5"/80-95.3 cm	8 pr/cs



Crutch Pillows

- » Padded fleece cushions for more comfort during crutch use



Item No.	Description	Pkg
MDSPPC107	Charcoal	6 sets/cs

3-in-1 Quick Fit
One crutch for
children through
adult sizes!



3-in-1 Quick-Fit Aluminum Crutches

- » One crutch fits three sizes Youth, Adult and Tall Adult
- » Incorporates the patented Red Dot® easy handgrip adjustment and the I-Beam height adjustment
- » 300-lb./136 kg weight capacity

Item No.	Description	Pkg
MDS80540	Quick-Fit, Patient Height 4'7"-6'7", Height Adjustment 38.5-59.5"/97.8-151.1 cm	2 pr/cs
MDS80540-8	Quick-Fit, Patient Height 4'7"-6'7", Height Adjustment 38.5-59.5"/97.8-151.1 cm	8 pr/cs

G30754B



Guardian® Bariatric Folding Walkers

- » Higher weight capacities help keep larger patients active and moving
- » Two-button folding capabilities let patients receive support through narrow spaces
- » Folds down to only 4.5-inches/11.4 cm

Item No.	Description	Pkg
G30754B	Extra-Wide Bariatric Walker, 650-lb./294.8 kg Weight Capacity, Inside Handgrip 22"W/50.8 cm, Height Adjustment 32-39"/81.3-99.1 cm	2/cs
MDS86410XW	Adult Bariatric Walker, 500-lb./226.8 kg Weight Capacity, Inside Handgrip 20"W/50.8 cm, Height Adjustment 31-41"/78.7-104.1 cm	1/ea
MDS86410XWW	Adult Bariatric Walker, 500-lb./226.8 kg Weight Capacity, Inside Handgrip 20"W/50.8 cm, Height Adjustment 31-41"/78.7-104.1 cm, with 5"/12.7 cm Dual Wheels	1/ea



MDS864104

Guardian® Two-Button Folding Walkers

- » Lightweight aluminum-frame and cross braces for extra durability
- » Easy-to-use, push-button mechanisms lets walker fold easily
- » Each side folds independently, allowing easy movement through narrow spaces
- » Weight capacity: 300-lbs./136 kg
- » Folded width: 4.5-inches/11 cm
- » Push-button, height-adjustable legs
- » Fits users 5'5"-6'4"/165-193 cm

Item No.	Description	Pkg
MDS864104	Adult	4/cs
MDS864104H	Adult	1/ea
MDS86410W4	Adult, 3"/7.6 cm Wheels	4/cs
MDS86410W54	Adult, 5"/12.7 cm Wheels	4/cs
MDS86410J4	Junior	4/cs
MDS86410JW4	Junior, 3"/7.6 cm Wheels	4/cs
MDS86410JW54	Junior, 5"/12.7 cm Wheels	4/cs

California customers—For all items on this page, **⚠️ WARNING: Cancer and Reproductive Harm—www.P65Warnings.ca.gov**



MEDLINE



Assist in Strengthening your facility.

Help improve patient outcomes with innovative, cost-effective therapy and rehab solutions.

Discover everything you need to elevate care and achieve results. We offer a wide assortment of therapy and rehab essentials, including innovative and exclusive products, well-known national brands, and quality Medline brand items.



Equipment

Advanced clinical equipment
to help set your facility apart



Modalities

To help manage pain
for **effective outcomes**



Supplies

Over **20,000** quality and value-driven
products to choose from



» Visit [medline.com/go/therapy](https://www.medline.com/go/therapy)
to view our full selection.



MDR141447



Digital Wheelchair Scales

- » Handrail-mounted or remote-mounted indicator designs to provide quality weighing solutions for a variety of patient weighing applications
- » Weight capacity: 1,000-lbs./453.6 kg
- » Graduation: 0.2-lb./0.1 kg
- » Power supply: 9v DC, provided by six AA alkaline batteries (included) or AC adapter (optional)
- » Unit of measure: lb. only, kg only, lb./kg
- » Functions: BMI, tare, hold, reweigh
- » EMR compatible, USB
- » Two-year warranty

Item No.	Description	Pkg
MDR141447	Digital Wheelchair Scale, Handrail (Dual Ramp), Platform Dimensions: 29.5 W x 32"D/ 74.9 x 81.3 cm	1/ea
MDR133077	AC Adapter	1/ea



CNS2101KL



Digital Platform Scales

- » Weight capacity: 1,000-lbs./453.6 kg
- » Resolution: 0.2-lb./0.09 kg
- » Connectivity via USB
- » 0.75-inch/1.9 cm LCD display, with 180° swivel for ease of viewing
- » Functions: lb/kg conversion, lb./kg lockout, body mass index (bmi), zero out/tare, hold/release, reweigh, auto zero, auto off
- » Handrails may be used during weighing for patient stability
- » Two wheels for ease of movement
- » Uses 6 AA batteries or 120V adapter (Includes adapter. Batteries sold separately.)
- » Two-year limited warranty

Item No.	Description	Pkg
CNS1101KL	22 W x 15.75"D/ 55.9 x 40 cm	1/ea
CNS2101KL	26 W x 22"D/66 x 55.9 cm	1/ea



Multifunctional—height and weight are indicated in a one-task operation. Weight information can be read from the front and back.



Heavy-duty wheels allow easy portability

MDR102613



Mechanical Physician Scale

- » Durable, accurate and easy to operate, move and store
- » Weight can be read from either the front or back
- » Slip-resistant base cover is removable for easy cleaning
- » Height rod: 24–84-inches/ 61–213.4 cm
- » Graduation: 0.25-lb./11.3 kg
- » Unit of measure: lbs./kg
- » Rear wheels for ease of portability
- » Two-year warranty

Item No.	Description	Pkg
MDR102613	Platform Dimensions 10.75 W x 14.75"D/ 27.3 x 37.5 cm, 440-lb./ 200 kg Weight Capacity	1/ea



MDR110611



Eye-Level Digital Physician Scale

- » Lightweight and easy to move, with two heavy-duty rear wheels
- » Height measuring rod: 27.5–82 x 0.125-inches/70–208 cm x 1 mm
- » Weight capacity: 550-lbs./250 kg
- » Graduation: 0.2-lb./0.1 kg
- » Power supply: 9V DC provided by six AA batteries (included) or AC adapter (optional)

- » Unit of measure: lbs. only, kg only, lbs./kg
- » Functions: BMI, hold, low-battery indication, auto shutoff
- » Two-year warranty
- » EMR capable, USB

Item No.	Description	Pkg
MDR110611	Digital Eye-Level Scale, Platform 14.5 W x 14"D/ 36.8 x 35.6 cm	1/ea
MDR133077	9V AC Adapter, 120v/ 230v Auto Switching	1/ea

California customers—For all items on this page, **⚠ WARNING: Cancer and Reproductive Harm—www.P65Warnings.ca.gov**



MPH07HP4



Eye- and Waist-Level Digital Physician Scales

- » Weight capacity: 500-lbs./226.8 kg
- » Resolution: 0.2-lb./0.1 kg
- » EMR connectivity via USB
- » 1-inch/2.5 cm LCD display
- » Functions: lb/kg conversion, lb./kg lock out, body mass index (bmi), zero out/tare, hold/release, auto zero, auto off
- » Height Rod: 30–84 ¼-inches/76–214 cm on MPH07HP4 only
- » Two wheels for ease of movement
- » Six AA batteries included, optional 120V adapter (CNSADPT31 AC adapter)
- » Two-year limited warranty

Item No.	Description	Pkg
MPH07HP4	Eye-Level Digital Scale, Platform 13.75 W x 16.5"D/ 34.9 x 41.9 cm	1/ea
MPH07H499KL	Waist-Level Digital Scale, Platform 13.75 W x 16.5"D/ 34.9 x 41.9 cm	1/ea
CNSADPT31	AC Adapter	1/ea



MDR440FD

Digital Step-On Scale

- » Provides precision and ease of use with its auto-on function
- » Backlit LCD display: 3 5/8 x 2.125-inches/9.2 x 5.4 cm
- » Display shows both lbs. and kg
- » Extra-wide platform surface: 14 x 11 7/8-inches/35.5 x 30.2 cm
- » Only one inch from the floor
- » Non-skid feet help prevent sliding
- » Weight capacity: 440-lbs./200 kg
- » One-year warranty

Item No.	Description	Pkg
MDR440FD	Digital Step-On Scale	1/ea

Power source requires four AA batteries (not included).



MDR121304



Digital Floor Level Physician Scale

- » Weight capacity: 550-lbs./249.5 kg
- » Graduation: 0.2-lb./0.1 kg
- » Power supply: 9V DC provided by six AA batteries (included) or AC adapter (optional)
- » Unit of measure: lbs only, kg only, lbs./kg
- » Functions: BMI (body mass index), hold, low-battery auto shutoff
- » EMR capable, USB
- » Two-year warranty

Item No.	Description	Pkg
MDR121304	Platform Dimensions: 14 W x 14.5 D x 3" H/ 35.6 x 36.8 x 7.6 cm	1/ea

6.5-foot/1.98 m cable between the scale base and indicator.



MDR75470



iPC Digital Dietary Scale

- » Compact and versatile with a sleek design
- » Simple operation and pinpoint accuracy
- » Convenient four-button operation
- » Easy to move and set up at remote locations
- » Six-digit LCD
- » Powered by two "D" cell batteries with up to 500 hours operation (batteries not included)
- » Functions: Auto power-off, low-battery indicator, overload stops, auto-ranging
- » Two-year limited warranty

Item No.	MDR75454	MDR75457
Capacity	6-lb.	30-lb.
Graduations	0.002-lb. from 0–3-lb./ 0–1.4 kg	0.01-lb. from 0–15-lb./ 0–6.8 kg
	0.005-lb. from 3–6-lb./ 1.4– 2.7 kg	0.02-lb. from 15–30-lb./ 6.7– 13.6 kg

Item No.	Description	Pkg
----------	-------------	-----

Accessories

MDR75470	Stainless Steel Platter Cover For 100-oz./ 2.8 kg, 6-lb./2.7 kg, and 15-lb./6.8 kg Capacity Scales	1/ea
MDR75473	AC Adapter (Batteries and AC Adapter Not Included)	1/ea



MPH07SD1

seca

Digital Diaper Scale

- » This scale is ideal for measuring a baby's urinary output
- » Scale and bowl are perfectly matched
- » Weight capacity: 6.6-lb./3 kg
- » Graduation: 0.05-oz./1 g
- » Dimensions: 8.3 W x 5.5 H x 8.3"D/ 210 x 140 x 210 mm
- » Functions: Tare, lb/kg conversion, automatic switch-off

Item No.	Description	Pkg
MPH07SD1	Digital	1/ea



MDR107423



Digital Baby Scale

- » Weight capacity: 44-lbs./20 kg
- » Graduation: 0.5-oz./14 g
- » Sleek, extremely durable, easy to clean and maintain
- » Platform dimensions: 11 x 23½-inches/28 x 60 cm
- » Functions: lb/kg conversion, recall key, auto shutoff
- » Large, easy-to-read 1-inch/2.5 cm LCD display
- » Built-in measuring tape 0–22½-inches/0–57.2 cm
- » Power: four AA batteries (not included) or AC adapter (included)
- » Two-year limited warranty

Item No.	Description	Pkg
MDR107423	Digital Baby Scale	1/ea
MDR112570	Carry Case	1/ea



SCA206

seca

Mechanical Measuring Tape with Wall Attachment

- » Attaches to the wall with just one screw
- » Durable metal tape
- » Display window in the headpiece
- » Measuring range: 0–87-inches/0–220 cm
- » Two-year warranty

Item No.	Description	Pkg
SCA206	Range: 0–87"/0–220 cm, Graduation: 0.125"/3.2 cm	1/ea



SCA201

seca

Ergonomic Circumference Measuring Tape

- » Measures girth with millimeter precision
- » Easy display unwinds and locks into place
- » Housing fits comfortably into user's hand and is tough enough to withstand accidental dropping

Item No.	Description	Pkg
SCA201	Range: 0–80"/0–205 cm, Graduation: 0.125"/3.2 cm	1/ea



MPH07SM1

seca

Measuring Mat for Babies

- » Allows precise measurement of babies and toddlers while they're lying down
- » Measuring range: 4–39-inches/10–99 cm
- » Washable and foldable
- » Fixed headpiece and sliding foot positioner
- » Easy-carry handle

Item No.	Description	Pkg
MPH07SM1	Measures in cm, Range: 10–99 cm, Graduation: 5 cm	1/ea



Growth Charts for Boys and Girls

- » Track height, weight, head circumference and BMI by age
- » Contains 100 sheets

Item No.	Description	Pkg
MDR158463	Birth to 36 Months, Boy	1/ea
MDR158464	Birth to 36 Months, Girl	1/ea
MDR158465	2–20 Years, Boy	1/ea
MDR158466	2–20 Years, Girl	1/ea

California customers—For all items on this page, **WARNING:** Cancer and Reproductive Harm—www.P65Warnings.ca.gov



SCA213



Freestanding Portable Stadiometer

- » A large measuring range of 8"-6'9"/20-205 cm
- » Graduation: 0.125-inch/3.175 mm
- » Measuring rod can be dismantled and set up easily and quickly
- » The spacer keeps the rod straight and stable for precise measuring
- » Quick and easy setup
- » Integrated handle for easy transport
- » Two-year warranty

Item No.	Description	Pkg
SCA213	Range: 8-81"/20-205 cm, Graduation: 0.125"/3.175 cm	1/ea



SCA217



Freestanding Stadiometer

- » Features an exceptionally high degree of stability
- » Measuring range: 8"-6'9"/20-205 cm
- » Graduation: 0.125-inch/3.175 mm
- » Robust base plate guarantees a sturdy stance
- » Integrated handle for easy transport
- » Two-year warranty

Item No.	Description	Pkg
SCA217	Range 8-81"/20-205 cm, Graduation: 0.125"/3.175 cm	1/ea

CNS201HR



Mechanical Wall-Mounting Measuring Rod

- » Features a variable scale suitable for measuring adults and children
- » Measuring range: 1.4-90-inches/3.5-230 cm
- » Graduation: 1/16"/1 mm
- » Frame can be easily mounted on the wall at desired height
- » Results displayed in side window
- » Locking screw to keep the height display in place
- » Two-year warranty

Item No.	Description	Pkg
CNS201HR	Range: Selectable from 1.4-90"/3.5-230 cm, Graduation: 0.0625"/1 mm	1/ea



EMS Sterilizer Monitoring Service

- » Includes one spore strip for mail-in laboratory certification
- » Test results are available 24 hours after packet is received at lab
- » Test results are available 24/7 online at www.checkyourtest.com

Item No.	Description	Pkg
SDCEMS012BX	EMS Sterilizer Monitoring Service	12/bx
CCYEMS052	EMS Sterilizer Monitoring Service	52/bx



SDCPP012

PassPort and PassPort®Plus Sterilizer Monitoring Services

- » PassPort includes two spore strips for mail-in laboratory certification - for separate test results for front and back of sterilizer
- » The PassPort Plus is the only mail-in spore test system that includes an immediate readout integrator with the spore test strip and one spore strip for mail-in laboratory certification
- » Certified for use in dry heat, chemical vapor, and EO gas processes
- » Test results are available 24 hours after packet is received at lab and 24/7 online at www.checkyourtest.com

Item No.	Description	Pkg
SDCPS012BX	PassPort Sterilizer Monitoring Service	12/bx
SDCPS052BX	PassPort Sterilizer Monitoring Service	52/bx
SDCPP012	PassPort Plus Sterilizer Monitoring Service	12/bx
SDCPP052	PassPort Plus Sterilizer Monitoring Service	52/bx

Don't forget about sterilization pouches and indicators. Check out our Instruments section on pages 79-92.



Ritter™ UltraClave Automatic Sterilizers

- » The Ritter M9 couldn't be easier and is time tested for reliability. You load it... set it... and forget it. Select your cycle and the time is set for you. Press "Start" and the Ritter M9 fills with the correct amount of water. Once the sterilization process is complete, the door opens automatically and quietly to dissipate the steam to give you fast drying time for your instruments. It's real convenience
- » The Ritter M9 sterilizer also has easy front access for filling and draining and controls that are easy to understand
- » Programmable cycle features allow custom creation of different cycle parameters for special sterilization needs. Programmable buttons allow you to change the following: Time (3 to 90 minutes), Temperature (230 to 275°F), Dry Time (0 to 60 minutes), Vent (slow or fast) and two program buttons allow you to store the custom cycles
- » Preprogrammed sterilization cycle controls: Unwrapped (132°C or 270°F for 3 minutes), Pouched (132 °C or 270°F for 5 minutes), and Packs (121 °C or 250 °F for 30 minutes)
- » The display indicates cycle selected, cycle temperature and exposure time for selected cycle. During the cycle, it shows messages describing status of cycle. When the cycle enters sterilization mode and dry mode, remaining cycle time is displayed as well as temperature and pressure



MIMM9022



Ritter M9D UltraClave Automatic Sterilizer

- » Steam-Flush Pressure Pulse Air Removal System helps ensure proper steam penetration
- » Simply load it, set it and forget it
- » Four pre-programmed and two programmable cycles
- » Once the sterilization process is complete, the door opens automatically and quietly to dissipate the steam
- » Four trays: two small, two large sizes

Item No.	Description	Pkg
MIMM9022	3.5-gal./1.4 Liters Automatic Door, 73-lbs.	1/ea
MIMM9D022	3.5-gal./1.4 Liters Manual Door, 73-lbs.	1/ea

Specifications:

Overall length with plug: 20.13-inches/51.13 cm
 Overall width: 15.3-inches/38.86 cm
 Overall height with printer: 15.8-inches/40.13 cm
 Countertop area: 15.3 x 20.38-inches/ 38.86 x 51.77 cm
 Chamber: 9-inch/22.9 cm diameter, 15-inch/38.1 cm depth, 3.5 gallon/13.2 L usable volume
 Trays: 2 Large
 7.31 W x 12 L x 0.875-inch D/ 18.57 x 30.48 x 2.22 cm;
 2 Small
 5.625 W x 12 L x 0.875-inch D/ 14.29 x 30.48 x 2.22 cm
 Weight: 73-lbs./33.1 kg
 Water reservoir capacity: 1.1 gallon/4.16 L
 Electrical: 115 VAC, 50/60 Hz, 15 AMPS
 Single Phase; 115 VAC has 1425-watt tubular immersion heater
 UL and CUL Listed; CSA Certified



Ritter M11 UltraClave Automatic Sterilizer

- » Features a door that opens automatically and quietly, making sterilization a set-it-and-forget-it procedure
- » 11-inch/27.9 cm by 18-inch/45.7 cm chamber
- » Perfect for larger packs and instruments
- » Passes the air removal tests required on more complex vacuum steam sterilizer systems
- » Four pre-programmed sterilization cycle controls: Unwrapped (132°C or 270°F for 3 minutes), Pouched (132°C or 270°F for 5 minutes), Packs (121°C or 250°F for 30 minutes), Handpiece (Dental Handpiece use only); two programmable controls to store time, temp, dry and vent times



MIMM11022

Item No.	Description	Pkg
MIMM11022	Ritter M11 UltraClave Automatic Sterilizer	1/ea



Valueklave® 1730

- » Provides basic sterilization needs at an economical price
- » Fast cycle time (11 min.) and heated drying deliver quick turnaround for wrapped or unwrapped instruments
- » Auto shutoff at the end of sterilization and dry cycles
- » Double safety locks prevent door from opening while chamber is pressurized
- » Compact 7-inch/18 centimeter chamber
- » Reliable with simple controls and high speed

Item No.	Description	Pkg
TUN1730	Valueklave 1730	1/ea



TUN1730

California customers—For all items on this page, **⚠️ WARNING: Cancer and Reproductive Harm—www.P65Warnings.ca.gov**



ARNA952

Electrosurgery Generators/Desiccators

» Bovie has you covered for all your electrosurgery generators and accessory needs in the office

Item No.	Description	Pkg
ARNA952	Aaron 940™ High-Frequency Desiccator, 40-Watt High-Frequency Dessicator with Power Control Handpiece	1/ea
ARNA950	Aaron 950™ High-Frequency Electrosurgical Generator/Desiccator, 60 Watts High-Frequency Dessication with Cut	1/ea
ARNA1250S	Aaron 1250 High-Frequency Electrosurgical Generator, Bovie Specialist PRO 120 Watt Multipurpose Electrosurgical Generator	1/ea
ARNA1250SG	Aaron 1250 OB/GYN Total System Solution, Bovie Specialist PRO-G Electrosurgery System with Smoke Evacuation	1/ea



ARNA805

Electrosurgery Accessories

Item No.	Description	Pkg
ARNA804	Sharp Tip, Non-Sterile	100/bx
ARNA805	Sharp Tip, Sterile	50/bx
ARNA806	Blunt Tip, Non-Sterile	100/bx
ARNA807	Blunt Tip, Sterile	50/bx
ARNA910	Sheath for ES Pencil	100/bx
ARNA910ST	Sterile Sheath for ES Pencil and Cord	25/bx
ARNES09	10 x 10 mm Loop Electrode	5/bx
ARNES53	15 x 12 mm Loop Electrode	5/bx
ARNES42	20 x 12 mm Loop Electrode	5/bx
ARNES13	20 x 15 mm Loop Electrode	5/bx
ARNES14	5 x 5 mm Square Electrode	5/bx
ARNES16	10 x 8 mm Square Electrode	5/bx



BOVCRY8

Freezpoint™ Cryosurgery Device

- » Safe, simple, and exceptionally cold device, which delivers an improved, low-cost-per-procedure, cryosurgical option
- » Freezpoint device application is able to reach -89°C and is up to 62% colder than common multiuse disposable systems
- » The colder temperature may provide a more effective application compared to the higher temperature produced by applicator-based systems

Item No.	Description	Pkg
BOVCRY8	Freezpoint Cryosurgical Device, 8 g	1/ea
BOVCRY16	Freezpoint Cryosurgical Device, 16 g	1/ea



BOVSE02

Smoke Evacuator

- » Designed to be easy-to-use and extremely effective at smoke and particle removal at the operative site
- » For use with all electrosurgical generators, the unit is also excellent for use in controlling laser and Harmonic® Scalpel smoke plume

Item No.	Description	Pkg
BOVSE02	Smoke Shark™ II Smoke Evacuator Unit with 35-Hr. Filter	1/ea



Pneumatic Exam Stools

- » Offer smooth height adjustment via a pneumatic lift
- » Pneumatic lift mechanism operates for optimum height range
- » Washable vinyl seat
- » Dense foam padding
- » Five brushed aluminum, dual-wheeled nylon casters on base for easy maneuverability



Item No.	Description	Height Adjustment	Seat Dimensions	Pkg
MDR505	Stool Without Back	17.75–24.5"/45.1–62.2 cm	16 x 4"/40.6 x 10.2 cm	1/ea
MDR505B	Stool with Back	17.75–24.5"/45.1–62.2 cm	16 x 4"/40.6 x 10.2 cm	1/ea
MDR505LC	Aluminum Base, Locking Caster	18.25–23.75"/46.4–60.3 cm	16 x 4"/40.6 x 10.2 cm	1/ea



MDR712155S	Stool Without Back	18.25–23.75"/46.4–60.3 cm	16 x 4.5"/40.6 x 11.4 cm	1/ea
MDR712155BS	Stool with Back	18.25–23.75"/46.4–60.3 cm	16 x 4.5"/40.6 x 11.4 cm	1/ea
MDR712135SH	Composite Base	19–24.5"/48.3–62.2 cm	16 x 4.5"/40.6 x 11.4 cm	1/ea
MDR712135B	Composite Base with Back	19–24.5"/48.3–62.2 cm	16 x 4.5"/40.6 x 11.4 cm	1/ea



Screw Lift Height-Adjustable Stools

- » Screw lift manually adjusts the seat height on these five-leg stools
- » Black nylon-composite base



Traditional Exam Stools

- » Manually adjust in height from 18.5–25.5-inches/ 47–65 cm using a screw lift
- » Screw lift height adjustment from 18.5–25.5-inches/ 47–65 cm
- » Seat is 14-inches/35.6 cm diameter with 4.5-inch/ 11.4 cm thick high-density foam
- » Four-leg base with rubber casters

Item No.	Description	Screw Lift Height Adjustment	Pkg
MDR712130	Stool without Back	18.5–22.75"/47–57.8 cm	1/ea
MDR712130B	Stool with Back	18.5–22.75"/47–57.8 cm	1/ea

Item No.	Description	Screw Lift Height Adjustment	Pkg
PEDP36	Square Foot Ring	17 7/8–26 1/4"/45.4–66.7 cm	1/ea
MDR712102	Round Foot Ring	18.5–25.5"/47–65 cm	1/ea


California customers—For all items on this page, **WARNING:** Cancer and Reproductive Harm—www.P65Warnings.ca.gov



Foot Stools

- » These stools have a handrail for safety
- » Rubber matting for comfort and traction
- » Non-slip rubber tips for stability
- » Durable, with 17 distinct weld points for enhanced strength

Item No.	Description	Pkg
Chrome		
MDR776250	600-lb. Cap w/Handrail 12.25 x 17	1/ea
MDR717150	350-lb. Cap w/Handrail 7/8" Diameter, 36" High	1/ea
Stainless Steel		
MPHSS150	With Handrail 11.25 X 14.25 x 9"H	1/ea
BLK7757SSHR	With 41"/104.1 cm Handrail, 18 W x 12 D x 8"H/ 45.7 x 30.5 x 20.3 cm	1/ea



Stacking Stools

- » Heavy-duty construction
- » Stackable for easy storage
- » Multiple stools lock side-to-side and end-to-end to create a wide variety of configurations
- » Optional wheeled dolly for easy transport and storage

Item No.	Description	Pkg
MPH08SS01C	18 W x 12 D x 6"H/45.7 x 30.5 x 15.2 cm	1/ea
MPH08SD01C	Dolly, Fits all 18 W x 12 D x 6"H/ 45.7 x 30.5 x 15.2 cm Stacking Stools	1/ea



Step Stool

- » Rubber matting for comfort and traction
- » Non-slip rubber tips provide additional safety
- » Features 17 weld points for enhanced strength

Item No.	Description	Pkg
MDS80430I	Chrome, 13 W x 17 D x 9.8"H/33 x 43.2 x 24.9 cm	1/ea



Chrome Hamper Stands

- » Chrome-plated, stainless steel and PVC stands
- » Foot and manual operated

Item No.	Description	Pkg
MDS80529	Chrome-Plated Steel with Painted Lid, 36-42-gal./136.3-159 L Capacity, 18.33 W x 19.5 D x 37.75"H/46.6 x 49.5 x 95.9 cm	1/ea
MDS80529P	Chrome-Plated Steel with Plastic Lid, 36-42-gal./136.3-159 L Capacity, 18.33 W x 19.5 D x 37.75"H/46.6 x 49.5 x 95.9 cm	1/ea
MDS80525	X-Frame, No Lid, 36-42-gal./136.3-159 L Capacity, 19 W x 19 D x 42"H/ 48.3 x 48.3 x 106.7 cm	1/ea



Basic X-Frame Hamper Stand

- » Open steel base frame
- » No lid
- » 2-inch/5 cm hard rubber casters
- » Bag diameter: 18, 22 or 25-inches/45.7, 55.9 or 63.5 cm
- » 42-inch/107 cm height

Item No.	Description	Pkg
MDS80525	Basic X-Frame Hamper	1/ea



Stainless Steel Hamper Stands

- » Constructed from 1-inch/2.5 cm stainless steel tubing
- » Solid steel bag-rest supports hamper bags
- » Non-slip pedal allows for hands-free operation
- » Hampers easily move with 2-inch/ 5.1 cm hard rubber casters

Item No.	Description	Pkg
MDS80530KD	Triangle Stainless Steel Lid with Foot Pedal, 18 D x 28–34”H/ 45.7 x 71.1–86.4 cm	1/ea
MPH08301	Triangular Shape with Lid, Non-Slip Pedal, 18 D x 34”H/ 68.1 x 86.4 cm	1/ea



Foot-Operated Stainless Steel Mayo Stands

- » Foot-pedal lock and “safe descend” help reduce the chance of accidents to ensure patient and staff security
- » Adjust in height from 36–62-inch/ 91–157 cm with foot pedal
- » Low-profile base slides easily under equipment

Item No.	Description	Pkg
MPH08MS02	21 W x 16”D/53.3 x 40.6 cm	1/ea
MPH08MS03	25 W x 20”D/63.5 x 50.8 cm	1/ea



Knob-Operated Stainless Steel Mayo Stands

- » Feature a Tru-Loc™ friction knob for manual locking at the desired height
- » Constructed on 1-inch/2.5 cm stainless steel tubing
- » Tray removes for easy cleaning

Item No.	Description	Pkg
MPH08741	19 W x 13”D/48.3 x 33 cm, Adjusts 39.5–62”H/ 100.3– 157.5 cm	1/ea
MPH08752	12.63 W x 19.13”D/ 32.1 x 48.6 cm, Adjusts 37.25–64”H/ 94.6–162.6 cm	1/ea



Stainless Steel Solution Stands

- » Feature rubber bumpers inside the ring to help reduce noise from metal-to-metal contact
- » Single or double basins with supporting frame
- » Connected by your choice of a bottom shelf or H-brace
- » Large 8.5-quart/8 Liter capacity

Item No.	Description	Pkg
MPH08017	Single Basin with Shelf, 15 W x 15 D x 33”H/ 38.1 x 38.1 x 83.8 cm	1/ea
MPH08015	Double Basin with Shelf, 29 W x 15 D x 33”H/ 73.7 x 38.1 x 83.8 cm	1/ea
MPH0807SSHB	Single Basin with H-Brace, 15 W x 15.25 D x 33”H/ 38.1 x 38.7 x 83.8 cm	1/ea
MPH0808SSHB	Double Basin with H-Brace, 14 W x 32 D x 33”H/ 35.6 x 81.3 x 83.8 cm	1/ea

Check out our Exam Room Supplies section beginning on page 45 for hamper liners and drapes.



MDR7211216FL

burton®

Outpatient™ II Diagnostic Procedure Lights

- » Provides maximum intensity and shadow-free illumination
- » Focusable flood-type light
- » Minimal heat output and high color temperature
- » Low-maintenance

Item No.	Description	Pkg
MDR7211216FL	Outpatient II Procedure Light, Floor Stand	1/ea
MDR7211122	Outpatient II Procedure Light, Single Ceiling Mount	1/ea
MDR0224132	Outpatient II Procedure Light, Double Ceiling Mount	1/ea
MDR7211142	Outpatient II Procedure Light, Wall Mount	1/ea



MPHM191M

Brewer

LED Gooseneck Exam Light

- » Bright white light supports clinical efficiency
- » Eco-friendly, long-life LED bulb for lower cost of ownership
- » Easy to maneuver with smooth, easy-to-clean head and flexible arm
- » Cool operation for extended use
- » Beige color complements your decor
- » Available with a mobile five-caster base
- » One-year standard warranty

Item No.	Description	Pkg
MPHM191M	LED Gooseneck Exam Light	1/ea
MPHM191M	LED Gooseneck Exam Light with Mobile Base	1/ea



MDR721607

DAZOR
LIGHTING TECHNOLOGY

Swivel Reflector Incandescent Exam Lamp

- » Swivel joint for easy rotation and is air cooled for comfortable handling
- » Anodized double-chambered reflector protects the reflective surface
- » 60-watt bulb (not included) is recommended, but UL-listed for 100-watt bulb for longer life and cooler operation
- » Five-year warranty

Item No.	Description	Pkg
MDR721607	Mobile Incandescent Exam Lamp, Maximum Extended Height is 69.5"/176.5 cm, Base Diameter is 24"/61 cm	1/ea
MDR721606	Stationary Incandescent Exam Lamp, Maximum Extended Height is 68.5"/174 cm, Base Diameter is 15"/38.1 cm	1/ea



MDR721010

CLINTON
LIGHTING

Incandescent Exam Lamp

- » Features a 12-inch/30 cm gooseneck and 50-74-inch/137-188 cm vertical height-adjustment
- » Shade concentrates light with minimal heat buildup
- » Rests firmly on a round 10-inch/25 cm diameter weighted base
- » 12-foot/3.6 m, hospital-grade, UL-approved cord included
- » Accepts up to a 60-watt bulb (not included)

Item No.	Description	Pkg
MDR721010	Incandescent Exam Lamp	1/ea



MDRSE50FL

burton®

Super Exam 50 LED

- » Long life
- » Low energy consumption
- » Versatile

Item No.	Description	Pkg
MDRSE50LSC	Ceiling Mount	1/ea
MDRSE50FL	Floor Stand	1/ea

California customers—For all items on this page, **WARNING:** Cancer and Reproductive Harm—www.P65Warnings.ca.gov



ARNCS103TLED



Colposcopes

- » Colpo-Master™ I is perfect for the OB/GYN performing both exams and procedures
- » The swing-arm design allows easy access to the vaginal canal aiding procedure efficiency

Item No.	Description	Pkg
ARNCS103TLED	Colpo-Master I Swing-Arm Colposcope 110V, 45° Camera-Ready Trinocular Zoom Head, 3 Leg-Base	1/ea
ARNCS103TDV	Swing-Arm Digital Video Package, Includes CC-103T-LED and CS-8185	1/ea

burton®

Handheld Ultraviolet Woods Exam Light

- » Portable exam lamp uses an optically ground and polished magnifier lens to provide a clear viewing area
- » Handheld flexibility
- » Three-diopter lens
- » Focal length of 8-inches/0.2 cm
- » Five-year limited warranty



MDR721109

Item No.	Description	Pkg
MDR721109	UV Light with Magnifier, Two White and Two UV Bulbs	1/ea
MDR7211602	UV Light with Magnifier, Two UV Bulbs	1/ea

burton®

Super Exam 50® Ceiling Light

- » Delivers higher illumination in a smaller area
- » Ideal for obstetrics, gynecology and emergency rooms
- » Uses halogen technology for maximum light output and color rendering
- » A dual lens optical system collects and concentrates the light in a precise beam



MDRSE50FL

Item No.	Description	Pkg
MDRSE50SC	Ceiling Mount	1/ea
MDRSE50FL	Floor Stand	1/ea



Medical Lighting

- » LED technology that improves the life of the light and the comfort of the patient, practitioner, and medical staff
- » Versatile—comes with a variety of mounting options

Item No.	Description	Pkg
ARNXLDPSC	MI 500 Single Ceiling 100V-240V Dual	1/ea
ARNXLDPSC	MI 750 Single Ceiling 100V-240V	1/ea
ARNXLDSC	MI 1000 Single Ceiling 100V-240V	1/ea



ARNXLDSC

California customers—For all items on this page, ⚠️ **WARNING:** Cancer and Reproductive Harm—www.P65Warnings.ca.gov



Heavy-Duty IV Stand

- » 21.5-inch/54.6 cm diameter with a low center of gravity
- » 5-leg cast iron base
- » Height adjusts 48–80-inches/121.9–203.2 cm
- » 50-lb. maximum weight capacity
- » Overall weight 18-lbs./8.2 kg

Item No.	Description	Pkg
MDS80500	Four-Hook IV Stand	2/cs



Deluxe IV Stand

- » 24-inch/61 cm diameter stand and a five-leg aluminum base
- » 2-inch/5 cm heavy-duty casters
- » Height adjusts from 56.75–96-inches/142–244 cm
- » Overall weight 8.3-lb./3.8 kg

Item No.	Description	Pkg
MDS80494	Four-Hook IV Stand	2/cs



Aluminum IV Pole

- » Tool-free push-button design makes adjusting easy
- » 5 wheel plastic base ensures stability
- » Accommodates two pumps up to 14 lbs. each and two 3L bags of fluid
- » Overall weight 5.8-lbs./2.6 kg

Item No.	Description	Pkg
MDS129CS	Two-Hook, Five-Leg	3/cs



UMF ProGlide 300+ Power Table

- » ADA height
- » 650-lb./294.8 kg patient weight capacity
- » Programmable hand control with two memory functions and a return-to-home
- » OneTouch WheelBase® raises and rolls for easy cleaning underneath
- » Seamless upholstery
- » Two storage drawers
- » Foot control
- » Drawer warmer
- » Receptacle

Item No.	Description	Pkg
UME4040650100	ProGlide Power 100, Power Back	1/ea
UME4040650200	ProGlide Power 200, Power Back	1/ea
UME4040650300	ProGlide Power 300, Power Back	1/ea
UME4070650100	ProGlide Manual 100, Manual Back	1/ea
UME4070650200	ProGlide Manual 200, Manual Back	1/ea
UME4070650300	ProGlide Manual 300, Manual Back	1/ea

*Choose from the upholstery color options shown on page 191. For more color options contact your Medline representative.

California customers—For all items on this page, **⚠️ WARNING: Cancer and Reproductive Harm—www.P65Warnings.ca.gov**

Brewer

Brewer Low-Access Power Table

- » Innovative chair design allowing for back support for disabled patients
- » Over 5.6 cubic feet of additional storage
- » Convenient pass-through work surface
- » Patented pass-through side drawers
- » Standard grab bars for added patient security
- » Chair-to-Table design with one-touch return to chair on power back models
- » Height: 18-inches/45.7 cm lowest, 28-inches/71.1 cm when flat
- » Length: 57-inches/144.8 cm extends to 78-inches/198.1 cm
- » 450-lb./204.1 kg weight capacity
- » 700-lb./317.5 kg weight capacity also available on power back models
- » Three-year standard warranty


Item No.	Description	Pkg
BRW5000	450-lb./204.1 kg Capacity, Manual Back	1/ea
BRW5001	450-lb./204.1 kg Capacity, Manual Back with Pelvic Tilt and Drawer Warmer	1/ea
BRW6501	450-lb./204.1 kg Capacity, Power Back with Pelvic Tilt and Drawer Warmer	1/ea
BRW6800	700-lb./317.5 kg Capacity, Power Back	1/ea
BRW6801	700-lb./317.5 kg Capacity, Power Back with Pelvic Tilt and Drawer Warmer	1/ea

*Choose from the upholstery color options shown on this page. For more color options contact your Medline representative.



BRW5001

Upholstery Color Options



Fabric colors are represented as closely as possible given printing limitations. For accurate color verification, please request an actual fabric swatch BEFORE ordering. All upholstered items are custom-made and non-returnable.



M5240



Ultra Comfort Exam Table

- » Features a vacuum-formed and seamless top
- » Easy-glide comfort stirrups: 4 lateral positions and friction lock length adjustment
- » Upholstered, removable leg rest
- » Pull-out, all-steel, non-slip step: 7-inches/17.8 cm high
- » Paper roll holder accepts 18- or 21-inch/45.7 or 53.3 cm paper
- » Plastic 9 x 9-inch/22.9 x 22.9 cm drip pan
- » Paper roll storage drawers
- » Two front and two side storage drawers (field reversible)
- » 500-lb./226.8 kg evenly distributed weight capacity
- » Five-year limited warranty

Specifications

Overall Dimensions: 27 W x 58.5–77.5 L x 31.75-inch H/
68.6 x 148.6–196.9 x 80.6 cm

infinitely adjustable over length of the bar

Front Drawers: (2) 18 W x 15.5 D x 4-inch H/45.7 x 39.4 x 10.2 cm

Side Drawers: (2) 8 W x 15.5 D x 4-inch H/45.7 x 39.4 x 10.2 cm

Item No.	Description	Pkg
M5240*	Ultra Comfort Exam Table	1/ea

*Choose from the upholstery color options shown on this page.

Upholstery Color Options



California customers—For all items on this page, **WARNING: Cancer and Reproductive Harm**—www.P65Warnings.ca.gov



Exam Table

- » One-piece upholstered top promotes quick, efficient cleaning after each patient
- » Steel frame for added strength
- » Multi-position stirrups adjust to standard and extra-wide angles for exceptional patient comfort, with a extra-long pull-out step
- » Independent operation of extra-wide leg extension and irrigation basin allows for quick transition during procedures
- » Five large storage drawers keep supplies conveniently located (drawers will be on the patient's right side unless otherwise specified)
- » UL-listed, hospital-grade, dual electrical outlets located on back side of table; circuit breaker system activates if power is interrupted
- » Paper roll holder accommodates 18–21-inch/ 45.7–53.3 cm paper
- » Back adjustment: horizontal to 80 degrees; 32-inch/ 81.3 cm height; 54-inch/137.2 cm length extends to 72-inches/182.9 cm; 27-inches/68.6 cm wide
- » 600-lb./272.2 kg weight capacity
- » Back adjusts by a counterbalance system, no gas cylinder
- » Four-year warranty

Item No.	Description	Pkg
MDR741700*	General Purpose Exam Table	1/ea
MDR741701*	General Purpose Exam Table with Pelvic Tilt and Warmer Drawer	1/ea

*Choose from cream white or pearl gray for base and beige, gray or natural oak for drawer fronts. Vinyl color options include: windsor wine, cambridge blue, emerald, ocean jade, ninja, smoke, fawn, and earthen brown.



Economy Exam Table

- » Quality exam table offering all the necessary components to make patients feel welcome and staff supplied
- » Large storage areas with doors on either side and two drawers in front
- » Concealed paper dispenser, pull-out leg extension, and nylon-coated, two-position hide-away stirrups
- » Ships as two components for easy single-person assembly
- » 400-lb./181.4 kg weight capacity

Specifications

Length: 56-inches/142.2 cm
 Width: 27-inches/68.6 cm
 Leg Rest: 72-inches/182.9 cm at full extension

Item No.	Description	Pkg
MDR748870*	Family Practice Table	1/ea

*For more color options, contact your Medline Representative..



Pediatric Table

- » Fun, kid friendly table featuring playful graphics
- » Easy-to-clean surfaces
- » Storage with convenient adjustable shelves
- » 350-lb./158.8 kg load capacity

Item No.	Description	Pkg
MDR741202	Pediatric Table	1/ea

California customers—For all items on this page, **⚠️ WARNING: Cancer and Reproductive Harm—www.P65Warnings.ca.gov**



Ritter 204 Manual Examination Tables

- » Enlarged footstep (over 40% larger than industry standard) and ergonomically designed top provides patients' easy access so they feel stable and secure during an exam

Additional Features:

- » Three side drawers for optimal supply management
- » Pass-through drawer design
- » Duplex, hospital-grade electrical receptacle on left side
- » Comfort-design stirrups help make lithotomy position easier for patients
- » Modern styling to suit any exam room
- » Removable top allows for thorough cleaning or quick upholstery change
- » Paper roll holder
- » 500-lb./226.8 kg capacity
- » Soft-touch stitched upholstery

Item No.	Description	Pkg
MIM204011	Exam Table, Base Only	1/ea
MIM204012	Exam Table with Pelvic Tilt/Drawer Heater, Base Only	1/ea



Upholstery Options for Ritter Examination Tables


- » Seamless upholstery helps eliminate places where dirt and bacteria can accumulate so cleanup is easier, faster, and more effective

Item No.	Description	Pkg
Styled Seamless		
MIM0210138846	Curative Copper	1/ea
MIM0210138843	Oasis	1/ea
MIM0210138853	Shaded Garden	1/ea
MIM0210138814	Stone	1/ea
MIM0210138847	Dream	1/ea
MIM0210138860	Latte	1/ea
MIM0210138855	Soothing Blue	1/ea
MIM0210138857	Obsidian	1/ea

- » An economical option that provides a soft and functional alternative for your practice

Item No.	Description	Pkg
Soft-Touch Stitched		
MIM0210137846	Curative Copper	1/ea
MIM0210137843	Oasis	1/ea
MIM0210137853	Shaded Garden	1/ea
MIM0210137814	Stone	1/ea
MIM0210137847	Dream	1/ea
MIM0210137860	Latte	1/ea
MIM0210137855	Soothing Blue	1/ea
MIM0210137857	Obsidian	1/ea

Upholstery Color Options



Color may vary based on monitor settings.

California customers—For all items on this page, **⚠️ WARNING: Cancer and Reproductive Harm—www.P65Warnings.ca.gov**



Brewer

Exam Patient Lift

- » Large non-slip platform with grab bars
- » Two emergency stop buttons
- » Four auto-lock rubber feet for stability
- » Four retractable swivel casters prevent rolling during use
- » 500-lb./226.8 kg weight capacity

Item No.	Description	Pkg
BRW3000	Exam Patient Lift	1/ea



MDREX23



MDREX21



MDREX27



Room Signs

- » Show little ones the way with these fun designs

Item No.	Description	Pkg
MDREX23	Doctor's Office Sign	1/ea
MDREX21	Exit Sign	1/ea
MDREX27	Lavatory Sign	1/ea



MDRT77



Sundry Jars

- » Keep everyday essentials handy with cartoon-labeled sundry jars

Item No.	Description	Pkg
MDRT77	Clear, Labeled Sundry Jar Set	1/pk



Treatment Tables

- » Feature solid hardwood legs and steel connecting support
- » 400-lb./181.4 kg weight capacity



MDR731024

Item No.	Description	Pkg
With Laminate Shelf		
MDR732024	24 W x 72 L x 31"H/61 x 182.9 x 78.7 cm	1/ea
MDR732027	27 W x 72 L x 31"H/61 x 182.9 x 78.7 cm	1/ea
MDR73102030	30 W x 72 L x 31"H/76.2 x 182.9 x 78.7 cm	1/ea

Without Shelf

MDR731024	24 W x 72 L x 31"H/61 x 182.9 x 78.7 cm	1/ea
MDR731027	27 W x 72 L x 31"H/68.6 x 182.9 x 78.7 cm	1/ea
MDR731030	30 W x 72 L x 31" H/76.2 x 182.9 x 78.7 cm	1/ea



MDR738505



Extra-Strength Steel-Leg Treatment Tables

- » 450-lb./181.4 kg weight capacity

Item No.	Description	Pkg
MDR7330102	Steel Frame, H-Brace, 27"W	1/ea
MDR301327	Steel Frame, Drawers, 27"W	1/ea
MDR3020270	Steel Frame, Shelf, 27"W	1/ea



MDR73856S



Treatment Tables with Adjustable Back

- » Back adjusts from nearly horizontal to almost upright
- » High-strength steel construction
- » High-impact leveling feet
- » 2-inch/5.1 cm thick, medium-density foam pad
- » Includes paper holder for 21 x 3-inch/53.3 x 7.6 cm roll
- » 400-lb./181.4 kg weight capacity

Item No.	Description	Pkg
MDR73856S	With Shelf, 28 W x 72 L x 30"H/71.1 x 182.9 x 76.2 cm	1/ea
MDR73857	With H-Brace, 28 W x 72 L x 30"H/71.1 x 182.9 x 76.2 cm	1/ea

California customers—For all items on this page, **⚠️ WARNING: Cancer and Reproductive Harm—www.P65Warnings.ca.gov**



Bariatric Phlebotomy Chairs

- » Extra-wide blood draw chair accommodates larger patients or a child with parent
- » Padded and upholstered seat, backrest and armrest for patient comfort
- » All-welded, heavy-duty frame for strength and stability
- » Both arms flip and adjust in height 29–34-inches/ 73.7–86.4 cm
- » Convenient “clean-out” between the seat and back for easy cleaning
- » 400-lb./181.4 kg weight capacity

Item No.	Description	Pkg
MDR66000	Seat: 34 L x 18 W x 20”H/86.4 x 45.7 x 50.8 cm Overall: 46 L x 26 W x 20”H/116.8 x 66 x 50.8 cm	1/ea



Padded Phlebotomy Chair

- » Heavy-gauge steel construction and welded frame with powder coating for extra durability
- » Padded and upholstered seat, backrest and armrests for patient comfort
- » Arms adjust in height 27–33-inches/68.6–83.8 cm
- » Optional one-piece plastic storage drawer
- » 300-lb./136.1 kg weight capacity

Item No.	Description	Pkg
MDR66010	Padded with Flip Arm	1/ea
MDR66020	Padded with Flip Arm, Storage Drawer	1/ea



Newport Collection

- » Durable and affordable, geared for high-use reception areas
- » Upward-sloping arms for easy egress
- » Matrex® seat support system for long-life comfort
- » Guest rated at 275-lb./124.7 kg, oversize at 400-lb./181.4 kg, and bariatric at 750-lb./340.2 kg capacity
- » Wall-saver design
- » Optional urethane arm caps available
- » Cleanout in back of seat
- » Chairs and loveseats are available for ganging
- » Ships UPS prepaid freight directly to your customer
- » Product requires minimal assembly
- » Limited lifetime warranty on frame

Item No.	Description	Pkg.
Guest Chair, Wall-Saver		
MDRNEWP1S	23 W x 25 D x 35”H/58.4 x 63.5 x 88.9 cm	1/ea
MDRNEWP1W	Wide, 27 W x 25 D x 35”H/68.6 x 63.5 x 88.9 cm	1/ea
MDRNEWPB33	Bariatric, 33 W x 25 D x 35”H/83.8 x 63.5 x 88.9 cm	1/ea
Loveseat		
MDRNEWPL	43 W x 25 D x 35”H/109.2 x 63.5 x 88.9 cm	1/ea

California customers—For all items on this page, **WARNING:** Cancer and Reproductive Harm—www.P65Warnings.ca.gov



Shown in Encinitas
CCD144X99ENC




Disposable Traditional Curtains

- » The look of reusable curtains with the convenience of a disposable curtain
- » Choose from a variety of colors: Tan, Coronado, Encinitas and Pet Parade
- » Traditional track attachment
- » No laundering required

Item No.	Dimensions	Pattern/Color	Pkg
CCD144X99TAN	144"W x 99"H/ 366cm x 251cm	Tan	4 ea/cs
CCD144X99COR	144"W x 99"H/ 366cm x 251cm	Coronado	4 ea/cs
CCD144X99ENC	144"W x 99"H/ 366cm x 251cm	Encinitas	4 ea/cs
CCD144X99PP	144"W x 99"H/ 366cm x 251cm	Pet Parade	4 ea/cs
CCD240X99TAN	240"W x 99"H/ 610cm x 251cm	Tan	4 ea/cs
CCD240X99COR	240"W x 99"H/ 610cm x 251cm	Coronado	4 ea/cs
CCD240X99ENC	240"W x 99"H/ 610cm x 251cm	Encinitas	4 ea/cs
CCD240X99PP	240"W x 99"H/6 10cm x 251cm	Pet Parade	4 ea/cs

Color Options



Tan
TAN

Coronado
COR

Encinitas
ENC

Pet Parade
PP



Shown in Spring Frost Cookie
MITQCSF109COK



Spring Frost® Cubicle Curtains


- » Basket weave pattern
- » Each includes 20-inches/50.8 cm of white mesh
- » Treated with an antimicrobial agent to inhibit the growth of bacteria and fungi on the product**

Item No.	Description	Pkg
MITQCSF101*	72 W x 82"H/182.9 x 208.3 cm	1/ea
MITQCSF102*	72 W x 88"H/182.9 x 223.5 cm	1/ea
MITQCSF103*	72 W x 94"H/182.9 x 238.8 cm	1/ea
MITQCSF104*	108 W x 82"H/274.3 x 208.3 cm	1/ea
MITQCSF105*	108 W x 88"H/274.3 x 223.5 cm	1/ea
MITQCSF106*	108 W x 94"H/274.3 x 238.8 cm	1/ea
MITQCSF107*	144 W x 82"H/365.8 x 208.3 cm	1/ea
MITQCSF108*	144 W x 88"H/365.8 x 223.5 cm	1/ea
MITQCSF109*	144 W x 94"H/365.8 x 238.8 cm	1/ea
MITQCSF110*	180 W x 82"H/457.2 x 208.3 cm	1/ea
MITQCSF111*	180 W x 88"H/457.2 x 223.5 cm	1/ea
MITQCSF112*	180 W x 94"H/457.2 x 238.8 cm	1/ea
MITQCSF113*	216 W x 82"H/548.6 x 208.3 cm	1/ea
MITQCSF114*	216 W x 88"H/548.6 x 223.5 cm	1/ea
MITQCSF115*	216 W x 94"H/548.6 x 238.8 cm	1/ea

*When ordering, add the color code to the end of the item number. All sizes listed above are cut sizes.

**This product doesn't protect users or others against bacteria, viruses, germs or other disease organisms.

Color Options



Blue Moon
BLM

Cookie
COK

Oatmeal
OA

Spa
SPA

California customers—For all items on this page, **WARNING:** Cancer and Reproductive Harm—www.P65Warnings.ca.gov



Shown in Uptown Eden
SNAP66PUTEDEN



Mesh Panel



Snap Panel System

- » For easy installation, storage and changes. Order mesh and fabric panels separately
- » Antimicrobial mesh: white only
- » Comes in six different sizes
- » Panels available in variety of colors
- » Number of panels determined by mesh width and track length, see chart below

Snap Panel Fabric

Item No.	Dimensions	Pattern/Color	Pkg
SNAP66PSFBLM	66"W x 66"H/168cm x 168cm	Spring Frost Blue Moon	6 ea/cs
SNAP66PSFOAT	66"W x 66"H/168cm x 168cm	Spring Frost Oatmeal	6 ea/cs
SNAP66PSFSPA	66"W x 66"H/168cm x 168cm	Spring Frost Spa	6 ea/cs
SNAP66PUTEDEN	66"W x 66"H/168cm x 168cm	Uptown Eden	6 ea/cs

Mesh Panels**

» All mesh panels use white mesh with white snaps

Item No.	Mesh Width	Mesh Height	Pkg	Item No.	Mesh Width	Mesh Height	Pkg
MIT6M066X20V	66"	20"	6 ea/cs	MIT6M066X28V	66"	28"	6 ea/cs
MIT6M124X20V	124"	20"	6 ea/cs	MIT6M124X28V	124"	20"	6 ea/cs
MIT6M184X20V	184"	20"	6 ea/cs	MIT6M184X28V	184"	20"	6 ea/cs
MIT6M242X20V	242"	20"	6 ea/cs	MIT6M242X28V	242"	20"	6 ea/cs
MIT6M300X20V	300"	20"	6 ea/cs	MIT6M300X28V	300"	20"	6 ea/cs


» Determine the number of panels needed for your track length

Track Length	50"-60"	61"-112"	113"-167"	168"-220"	221"-272"
Mesh Width	66"	124"	184"	242"	300"
Number of Panels	1 Panel	2 Panels	3 Panels	4 Panels	5 Panels

*Data on file. This product is treated with an antimicrobial additive to protect the product. This product does not protect users or others against bacteria, germs or other disease organisms.

**A reusable top mesh is required for Snap Panel Systems.


Color Options



Spring Frost Blue Moon
SFBLM



Spring Frost Oatmeal
SFOA



Spring Frost Spa
SFSPA



Uptown Eden
UTEDEN



Supreme Cubicle Track Items



- » Easy-to-install track features a durable reinforced spine and self-lubricating wheels for smooth operation
- » Track includes three carriers per foot (for additional carriers order Supreme Carrier, item number MCD114130)

Item No.	Description	Pkg
Straight		
SUP2FTSTR	24"/60.96 cm	1/ea
SUP3FTSTR	36"/91.44 cm	1/ea
SUP4FTSTR	48"/121.92 cm	1/ea
SUP5FTSTR	60"/152.4 cm	1/ea
SUP6FTSTR	72"/182.88 cm	1/ea
SUP7FTSTR	84"/213.36 cm	1/ea
SUP10FTSTR	94"/238.76 cm	1/ea

With 45° Bend		
SUP2X2FTW45	24 x 24"/60.96 x 60.96 cm	1/ea

With 90° Bend		
SUP2X2FTW90	24 x 24"/60.96 x 60.96 cm	1/ea

Tracks of 8 feet/2.44 m or shorter lengths come unspliced in one piece. All tracks of 8 feet/2.44 m or greater lengths are shipped in two pieces. For uncut pieces, add "U" to the end of the item number (unspliced). Special freight and handling charges will apply. Call for shipping quote.



Medline offers a full line of attractive and elegant cubicle curtain fabrics. Additional fabric options can be found on this catalog, as well as options for curtain tracks and disposable items in our Stocked Interiors catalog.



PEACE OF MIND, DELIVERED.

You found the capital equipment you need. Now, how do you get it delivered to your facility, quickly and efficiently? We can help with world-class services to meet your unique needs.

For more information, or to get a quote from our team of logistics experts, email capitalfreight@medline.com or contact your Medline Representative.

Prior to delivery

- » Warehousing and consolidation
- » Preassembly

During receiving

- » Dock-to-dock delivery*
- » Lift-gate delivery**
- » Inside delivery to specified location

Post-delivery

- » Unpacking and assembly
- » Wall-mounting
- » Old equipment removal
- » Garbage removal with dumpster

*Let us know if our 53-ft. carrier won't fit in your dock. Additional charges apply if we have to return to our hub, put your equipment on a smaller truck and attempt a second delivery.
**Lift gates must be requested when ordering as not all our trucks have them.



OFFICE SUPPLIES

Disposable Serveware and Supplies	200-205
Breakroom	206
Office Supplies	207-208



Plastic Cups

- » Perfect for serving a variety of cold beverages
- » Made of durable plastic
- » Recyclable in many recycling programs

Item No.	Description	Pkg	Compatible Lids
NON030035	3.5-oz./103.5 ml	100/bg, 25 bg/cs, 2,500 ea/cs	—
NON03005	5-oz./147.9 ml	100/bg, 25 bg/cs, 2,500 ea/cs	NON0357L
NON03007	7-oz./207 ml	100/bg, 25 bg/cs, 2,500 ea/cs	NON0357L
NON03009	9-oz./266.2 ml	100/bg, 25 bg/cs, 2,500 ea/cs	NON0309L
NON03010BAR	10-oz./295.7 ml, Graduated	50/bg, 20 bg/cs, 1,000 ea/cs	NON03010L
NON03012	12-oz./354.9 ml	50/bg, 20 bg/cs, 1,000 ea/cs	NON03012L
NON03014	14-oz./414 ml	50/bg, 20 bg/cs, 1,000 ea/cs	NON03012L
NON03016	16-oz./473.2 ml	50/bg, 20 bg/cs, 1,000 ea/cs	NON03016L



Paper Cups

- » Perfect for serving a variety of hot or cold beverages
- » Barium graduated cups have easy-to-read markings in both imperial and metric units

Item No.	Description	Pkg	Compatible Lids
Cold			
NON05003	3-oz./88.7 ml	125/pk, 40 pk/cs, 5,000 ea/cs	—
NON05003BAR	3-oz./88.7 ml, Graduated	125/pk, 40 pk/cs, 5,000 ea/cs	—
NON05005	5-oz./147.9 ml	100/pk, 30 pk/cs, 3,000 ea/cs	—
NON05007	7-oz./207 ml	100/pk, 20 pk/cs, 2,000 ea/cs	—
NON05009	9-oz./266.2 ml	100/pk, 20 pk/cs, 2,000 ea/cs	NON0309L
NON05012	12-oz./354.9 ml	50/pk, 40 pk/cs, 2,000 ea/cs	NONL16BL
NON05016	16-oz./473.2 ml	50/pk, 40 pk/cs, 2,000 ea/cs	NONL16BL
NON05016BAR	16-oz./473.2 ml, Graduated	50/bg, 20 bg/cs, 1,000 ea/cs	NONL16BL
NON05020	20-oz./591.5 ml	50/bg, 20 bg/cs, 1,000 ea/cs	NONL16BL
Hot			
NON7008HDL	8-oz./236.6 ml, With Handle	50/bg, 20 bg/cs, 1,000 ea/cs	DTM8FTL, NONPLCLID8, BEZBEHL8PLA
NONDWHC8PE	8-oz./236.6 ml, Double-Wall	20/pk, 25 pk/cs, 500 ea/cs	DTM8FTL, NONPLCLID8, BEZBEHL8PLA



Foam Cups

- » Polystyrene cups are made from insulating thin-wall foam that keeps hot beverages hot and cold beverages cold
- » Coordinating lids prevent spills and heat from escaping



WNCPC6



WNCPC1618CC

Item No.	Description	Pkg	Compatible Lids
WNCPC6	6-oz./177.4 ml	50/bg, 20 bg/cs, 1,000 ea/cs	—
DTM8J8	8-oz./236.6 ml	25/bg, 40 bg/cs, 1,000 ea/cs	DTM8SL, DTM8JL, DTM8FTL
DTM10J10	10-oz./295.7 ml	25/bg, 40 bg/cs, 1,000 ea/cs	DTM10SL
DTM12J12	12-oz./354.9 ml	25/bg, 40 bg/cs, 1,000 ea/cs	DTM16SL, DTM12JL, DTM16FTLS, DTM16RCL
DTM14J16	14-oz./414 ml	25/bg, 40 bg/cs, 1,000 ea/cs	DTM16SL, DTM12JL, DTM16FTLS, DTM16RCL
WNCPC1618	16-oz./473.2 ml	25/bg, 20 bg/cs, 500 ea/cs	DTM16SL, DTM12JL, DTM16FTLS, DTM16RCL
WNCPC1618CC	16-oz./473.2 ml, Graduated	25/bg, 20 bg/cs, 500 ea/cs	DTM16SL, DTM12JL, DTM16FTLS, DTM16RCL
WNCPC20J16	20-oz./591.5 ml	20/bg, 25 bg/cs, 500 ea/cs	DTM16SL, DTM12JL, DTM16FTLS, DTM16RCL
WNCPC24C18	24-oz./709.8 ml	15/bg, 20 bg/cs, 300 ea/cs	DTM16SL, DTM12JL, DTM16FTLS, DTM16RCL
WNCPC32TJ32	32-oz./946.4 ml	20/bg, 25 bg/cs, 500 ea/cs	DTM32SL



Compostable Cups

- » Compostable cups and lids provide an environmentally-friendly option
- » Classic hot cups are made from all natural plant sugars that can withstand boiling water and contain double sealed seams for exceptional strength
- » Crystal clear cold cups are made of a “green” and environmentally friendly PLA resin



BEZBEHC8PLA



NONPLC8



NONPLC8DW

Item No.	Description	Pkg	Compatible Lids
Hot, Single Wall			
BEZBEHC4PLA	4-oz./118.3 ml	50/bg, 20 bg/cs, 1,000 ea/cs	—
BEZBEHC8PLA	8-oz./236.6 ml	50/bg, 20 bg/cs, 1,000 ea/cs	BEZBEHL8PLA, NONPLCLID8
BEZBEHC10PLA	10-oz./295.7 ml	50/bg, 20 bg/cs, 1,000 ea/cs	BEZBEHL1024PL, NONPLCLID
BEZBEHC12PLA	12-oz./354.9 ml	50/bg, 20 bg/cs, 1,000 ea/cs	BEZBEHL1024PL, NONPLCLID
BEZBEHC16PLA	16-oz./473.2 ml	50/bg, 20 bg/cs, 1,000 ea/cs	BEZBEHL1024PL, NONPLCLID
BEZBEHC20PLA	20-oz./591.5 ml	50/bg, 20 bg/cs, 1,000 ea/cs	BEZBEHL1024PL, NONPLCLID
NONPLC8	8-oz./236.6 ml	50/bg, 20 bg/cs, 1,000 ea/cs	BEZBEHL8PLA, NONPLCLID8
NONPLC12	12-oz./354.9 ml	50/bg, 20 bg/cs, 1,000 ea/cs	BEZBEHL1024PL, NONPLCLID
NONPLC16	16-oz./473.2 ml	50/bg, 20 bg/cs, 1,000 ea/cs	BEZBEHL1024PL, NONPLCLID
NONPLC20	20-oz./591.5 ml	25/bg, 20 bg/cs, 500 ea/cs	BEZBEHL1024PL, NONPLCLID

Hot, Double Wall

BEZBEDWHC8PLA	8-oz./236.6 ml	25/bg, 40 bg/cs, 500 ea/cs	BEZBEHL8PLA, NONPLCLID8
BEZBEDWH12PLA	12-oz./354.9 ml	25/bg, 40 bg/cs, 500 ea/cs	BEZBEHL1024PL, NONPLCLID
BEZBEDWH16PLA	16-oz./473.2 ml	25/bg, 40 bg/cs, 500 ea/cs	BEZBEHL1024PL, NONPLCLID
BEZBEDWH20PLA	20-oz./591.5 ml	25/bg, 40 bg/cs, 500 ea/cs	BEZBEHL1024PL, NONPLCLID
NONPLC8DW	8-oz./236.6 ml	50/bg, 20 bg/cs, 1,000 ea/cs	BEZBEHL8PLA, NONPLCLID8
NONPLC12DW	12-oz./354.9 ml	40/bg, 25 bg/cs, 1,000 ea/cs	BEZBEHL1024PL, NONPLCLID
NONPLC16DW	16-oz./473.2 ml	40/bg, 15 bg/cs, 600 ea/cs	BEZBEHL1024PL, NONPLCLID
NONPLC20DW	20-oz./591.5 ml	40/bg, 15 bg/cs, 600 ea/cs	NONPLCLID20DW

Cold, Clear Plastic

NONPLA2	2-oz./59.1 ml	50/bg, 40 bg/cs, 2,000 ea/cs	NONPLASL
NONPLA4	4-oz./118.3 ml	50/bg, 20 bg/cs, 1,000 ea/cs	NONPLASL
NONPLA9	9-oz./266.2 ml	50/bg, 20 bg/cs, 1,000 ea/cs	NONPLAFL, NONPLADL, NONPLAYL
NONPLA12	12-oz./354.9 ml	50/bg, 20 bg/cs, 1,000 ea/cs	NONPLAFL, NONPLADL, NONPLAYL
NONPLA16	16-oz./473.2 ml	50/bg, 20 bg/cs, 1,000 ea/cs	NONPLAFL, NONPLADL, NONPLAYL
NONPLA20	20-oz./591.5 ml	50/bg, 20 bg/cs, 1,000 ea/cs	NONPLAFL, NONPLADL, NONPLAYL
NONPLA24	24-oz./709.8 ml	50/bg, 20 bg/cs, 1,000 ea/cs	NONPLAFL, NONPLADL, NONPLAYL



NON024214

NON034215

Portion/Souffle Cups

- » Ideal for use as medicine cups

Item No.	Description	Pkg
Souffle		
NON024214	0.50-oz./14.8 ml, Paper	250/bg, 20 bg/cs, 5,000 ea/cs
NON024215	0.75-oz./22.2 ml, Paper	250/bg, 20 bg/cs, 5,000 ea/cs
NON024220	1-oz./29.6 ml, Paper	250/bg, 20 bg/cs, 5,000 ea/cs
NON024230	3.5-oz./103.5 ml, Paper	100/bg, 25 bg/cs, 2,500 ea/cs

Item No.	Description	Pkg
Portion		
NON034215	0.75-oz./22.2 ml, Plastic	250/bg, 20 bg/cs, 5,000 ea/cs
DTM100PC	1-oz./29.6 ml, Plastic	125/bg, 20 bg/cs, 2,500 ea/cs
DTM200PC	2-oz./59.1 ml, Plastic	125/bg, 20 bg/cs, 2,500 ea/cs
DTM400PC	4-oz./118.3 ml, Plastic	125/bg, 20 bg/cs, 2,500 ea/cs

Clear Portion Cup Lids

DTMPL1N	1-oz./29.6 ml, Plastic	125/bg, 20 bg/cs, 2,500 ea/cs
DTMPL2N	2-oz./59.1 ml, Plastic	125/bg, 20 bg/cs, 2,500 ea/cs
DTMPL4N	4-oz./118.3 ml, Plastic	125/bg, 20 bg/cs, 2,500 ea/cs



DTM8JL



DTM8FTL



DTM16FTLS



DTM8SL/DTM10SL



DTM12JL



DTM16RCL



DTM16SL/DTM32SL



NONL16BL



NONPLADL



NONPLCLID

Lids

- » Helps keep hot beverages hot and cold beverages cold, in addition to preventing spills
- » Available in plastic or compostable material with many styles

Item No.	Description	Pkg
----------	-------------	-----



Compostable Cup Lids

BEZBEHL8PLA	8-oz./236.6 ml, Compostable, Translucent, Hot, Sip-Thru	50/bg, 1,000 ea/cs
BEZBEL1024PLA	10-24 oz./295.7-709.8 ml, Compostable, Translucent, Hot, Sip-Thru	50/bg, 1,000 ea/cs
NONPLCLID8	8-oz./236.6 ml, Compostable, Translucent, Hot, Sip-Thru	50/bg, 1,000 ea/cs
NONPLCLID	10-20 oz./295.7-591.5 ml, Compostable, Translucent, Hot, Sip-Thru	50/bg, 1,000 ea/cs
NONPLCLID20DW	20-oz./591.5 ml, Compostable, Translucent, Hot DW, Sip-Thru	50/bg, 600 ea/cs



Foam Cup Lids

DTM8FTL	8-oz./236.6 ml, Plastic, White, Hot, Pull-Back	100/bg, 1,000 ea/cs
DTM8JL	8-oz./236.6 ml, Plastic, White, Hot, Vented	100/bg, 1,000 ea/cs
DTM8SL	8-oz./236.6 ml, Plastic, Translucent, Cold, Straw Slot	100/bg, 1,000 ea/cs
DTM10SL	10-oz./295.7 ml, Plastic, Translucent, Cold, Straw Slot	100/bg, 1,000 ea/cs
DTM12JL	12-24-oz./354.9-709.8 ml, Plastic, Translucent, Cold, Vented	100/bg, 1,000 ea/cs
DTM16FTLS	12-24-oz./354.9-709.8 ml, Plastic, Translucent, Pull-Back/Straw Slot	100/bg, 1,000 ea/cs
DTM16RCL	12-24-oz./354.9-709.8 ml, Plastic, White, Hot, Reclosable	100/bg, 1,000 ea/cs
DTM16SL	12-24-oz./354.9-709.8 ml, Plastic, Translucent, Cold, Straw Slot	100/bg, 1,000 ea/cs
DTM32SL	32-oz./946.4 ml, Plastic, Translucent, Cold, Straw Slot	100/bg, 500 ea/cs



Plastic/Paper Cup Lids

NON0357L	5-7 oz./147.9-207 ml, Plastic, Translucent, Cold, Straw Slot	100/bg, 2,500 ea/cs
NON0309L	9-oz./266.2 ml, Plastic, Translucent, Cold, Straw Slot	100/bg, 2,500 ea/cs
NON03010L	10-oz./295.7 ml, Plastic, Translucent, Cold, Straw Slot	50/bg, 1,000 ea/cs
NON03012L	12-14 oz./354.9-473.2 ml, Plastic, Translucent, Cold, Straw Slot	100/bg, 1,000 ea/cs
NON03016L	16 oz./473.2 ml, Plastic, Translucent, Cold, Straw Slot	100/bg, 1,000 ea/cs
NONL16BL	12-16 oz./354.9-473.2 ml, Plastic, Translucent, Cold, Straw Slot	125/bg, 2,000 ea/cs



Compostable Plastic Cup Lids

NONPLASL	2-4 oz./59.1-118.3 ml, Compostable, Clear, Cold, Flat	50/bg, 1,000 ea/cs
NONPLAFL	9-24 oz./266.2-709.8 ml, Compostable, Clear, Cold, Flat	100/bg, 1,000 ea/cs
NONPLADL	9-24 oz./266.2-709.8 ml, Compostable, Clear, Cold, Dome	100/bg, 1,000 ea/cs
NONPLAYL	9-24 oz./266.2-709.8 ml, Compostable, Clear, Cold, Parfait	50/bg, 1,000 ea/cs



Compostable Bowls, Plates & Containers

- » Compostable sugarcane/bagasse will break down in commercial composting facilities
- » Fully compostable plastic products have the look and feel of traditional plastic without damaging petro-chemicals

Item No.	Description	Pkg
Plates		
ECOPLATE9	9"/22.9 cm, Round, Sugarcane, Medline	125/bg, 500 ea/cs
BEZBEFRP6	6"/15.2 cm, Round, Bagasse	125/bg, 1,000 ea/cs
BEZBEFRP7	7"/17.8 cm, Round, Bagasse	125/bg, 1,000 ea/cs
BEZBEFRP9	9"/22.9 cm, Round, Bagasse	125/bg, 500 ea/cs
BEZBEFRP93	9"/22.9 cm, Round, Bagasse, 3 comp	125/bg, 500 ea/cs
BEZBEFRP10	10"/25.4 cm, Round, Bagasse	125/bg, 500 ea/cs
BEZBEFRP103	10"/25.4 cm, Round, Bagasse, 3 comp	125/bg, 500 ea/cs
BEZBEFSP6	6"/15.2 cm, Square, Sugarcane	50/bg, 1,000 ea/cs
BEZBEFSP8	8"/20.3 cm, Square, Sugarcane	50/bg, 500 ea/cs
BEZBEFSP10	10"/25.4 cm, Square, Sugarcane	50/bg, 500 ea/cs
NONP004	6"/15.2 cm, Round, Sugarcane, Planet	50/bg, 1,000 ea/cs
NONP013	9"/22.9 cm, Round, Sugarcane, Planet	50/bg, 500 ea/cs

Bowls/Containers

ECOBOWL12	12-oz./354.9 ml, Bowl, Sugarcane, Medline	250/bg, 1,000 ea/cs
BEZBEFRB12	12-oz./354.9 ml, Bowl, Bagasse	50/bg, 1,000 ea/cs
BEZBEFRB16	16-oz./473.2 ml, Bowl, Bagasse	50/bg, 1,000 ea/cs
BEZBESC8PLA	8-oz./236.6 ml, Container, Paper	50/bg, 1,000 ea/cs
BEZBESC12PLA	12-oz./354.9 ml, Container, Paper	25/bg, 500 ea/cs
BEZBESC16PLA	16-oz./473.2 ml, Container, Paper	25/bg, 500 ea/cs
BEZBESC24PLA	24-oz./709.8 ml, Container, Paper	25/bg, 500 ea/cs
BEZBESC32PLA	32-oz./946.4 ml, Container, Paper	25/bg, 500 ea/cs
NONB003	6 x 6 x 3"/15.2 x 15.2 x 7.6 cm, Container, Sugarcane	50/bg, 250 ea/cs
NONBCS09	8 x 8 x 3"/20.3 x 20.3 x 7.6 cm, Container, Sugarcane	50/bg, 200 ea/cs
NONBCS093	8 x 8 x 3"/20.3 x 20.3 x 7.6 cm, Container, Sugarcane, 3 comp	50/bg, 200 ea/cs
NONBCS10	9 x 9 x 3"/22.9 x 22.9 x 7.6 cm, Container, Sugarcane	50/bg, 200 ea/cs
NONPLAEC6	6 x 6"/15.2 x 15.2 cm, Container, Hinged, Plastic	240 ea/cs
NONPLAEC8	8 x 8"/20.3 x 20.3 cm, Container, Hinged, Plastic	160 ea/cs

Lids for Bowls/Containers

BEZBEFL90CPLA	8-oz./236.6 ml, Plastic, Trans, Vented	50/bg, 1,000 ea/cs
BEZBEFL115PLA	12-32-oz./354.9-946.4 ml, Plastic, Trans, Vented	50/bg, 1,000 ea/cs



Standard Bowls, Plates and Containers

- » Bowls, plates and containers are available in assorted sizes and materials including foam, paper and plastic
- » Foam food containers feature a hinged lid with two secure closure tabs that protect against accidental openings and spills during transport

Item No.	Description	Pkg
Plates, White		
GPK80600	6"/15.2 cm, Round, Foam	125/bg, 1,000 ea/cs
GPK80900	9"/22.9 cm, Round, Foam	125/bg, 500 ea/cs
GNP70600	6"/15.2 cm, Round, Plastic	125/bg, 1,000 ea/cs
GNP71300	10"/25.4 cm, Round, Plastic	125/bg, 1,000 ea/cs
AJMPP6GREWH	6"/15.2 cm, Round, Paper	100/bg, 1,000 ea/cs
AJMPP9GREWH	9"/22.9 cm, Round, Paper	100/bg, 1,000 ea/cs
Bowls/Containers		
GPK82100	12-oz./354.9 ml, Bowl, Foam, White	125/bg, 1,000 ea/cs
GNPSN240	8 x 8 x 3"/20.3 x 20.3 x 7.6 cm, Container, Foam, 1 comp	200 ea/cs
GNPSN203	9 x 9 x 3"/22.9 x 22.9 x 7.6 cm, Container, Foam, 3 comp	200 ea/cs
WNCPF6	6-oz./177.4 ml, Container, Foam, White	40/bg, 1,000 ea/cs
WNCPF8	8-oz./236.6 ml, Container, Foam, White	25/bg, 500 ea/cs
WNCPF12	12-oz./354.9 ml, Container, Foam, White	25/bg, 500 ea/cs
Lids for Bowls/Containers		
WNCPFL6V	6-oz./177.4 ml, Plastic, Trans, Vented	100/bg, 1,000 ea/cs
WNCPFL8V	8-12-oz./236.6-354.9 ml, Plastic, Trans, Vented	100/bg, 1,000 ea/cs



Compostable Disposable Cutlery

» Fully compostable plastic products have the look and feel of traditional plastic without damaging petro-chemicals

Item No.	Description	Pkg
----------	-------------	-----

Standard-Weight Utensils

BEZBEFMW	Fork, White, Compostable	1,000/cs
BEZBEFMWINV	Fork, White, Compostable, Wrapped	1,000/cs
BEZBEFMB	Fork, Black, Compostable	1,000/cs
BEZBEFMBINV	Fork, Black, Compostable, Wrapped	1,000/cs
BEZBESMW	Spoon, White, Compostable	1,000/cs
BEZBESMWINV	Spoon, White, Compostable, Wrapped	1,000/cs
BEZBESMB	Spoon, Black, Compostable	1,000/cs
BEZBESMBINV	Spoon, Black, Compostable, Wrapped	1,000/cs
BEZBESSW	Soup Spoon, White, Compostable	1,000/cs
BEZBESSWINV	Soup Spoon, White, Compostable, Wrap	1,000/cs
BEZBESSB	Soup Spoon, Black, Compostable	1,000/cs
BEZBESSBINV	Soup Spoon, Black, Compostable, Wrap	1,000/cs
BEZBEKMW	Knife, White, Compostable	1,000/cs
BEZBEKMWINV	Knife, White, Compostable, Wrapped	1,000/cs
BEZBEKMB	Knife, Black, Compostable	1,000/cs
BEZBEKMBINV	Knife, Black, Compostable, Wrapped	1,000/cs
BEZBECKMW	Fork/Spoon/Knife/Napkin Combo, White	500/cs
BEZBECKMB	Fork/Spoon/Knife/Napkin Combo, Black	500/cs

Heavy-Weight Utensils

NONCPLA-002	Fork, White, Compostable	1,000/cs
BEZBEFHW	Fork, White, Compostable	1,000/cs
BEZBEFHWINV	Fork, White, Compostable, Wrapped	1,000/cs
BEZBEFHB	Fork, Black, Compostable	1,000/cs
BEZBEFHBINV	Fork, Black, Compostable, Wrapped	1,000/cs
NONCPLA-003	Spoon, White, Compostable	1,000/cs
BEZBESHW	Spoon, White, Compostable	1,000/cs
BEZBESHWINV	Spoon, White, Compostable, Wrapped	1,000/cs
BEZBESHB	Spoon, Black, Compostable	1,000/cs
BEZBESHBINV	Spoon, Black, Compostable, Wrapped	1,000/cs
NONCPLA-001	Knife, White, Compostable	1,000/cs
BEZBEKHW	Knife, White, Compostable	1,000/cs
BEZBEKHWINV	Knife, White, Compostable, Wrapped	1,000/cs
BEZBEKHB	Knife, Black, Compostable	1,000/cs
BEZBEKHBINV	Knife, Black, Compostable, Wrapped	1,000/cs
BEZBECKHW	Fork/Spoon/Knife/Napkin Combo, White	500/cs
BEZBECKHB	Fork/Spoon/Knife/Napkin Combo, Black	500/cs



Standard Disposable Cutlery

- » Disposable utensils are perfect for diners to put in take-out bags or grab from a self-service station when dining in while saving dish washing time
- » Standard utensils are available in medium and heavy weights
- » Wrapped utensils prevent utensils from being exposed to unwanted germs

Item No.	Description	Pkg
----------	-------------	-----

Standard-Weight Utensils

NON042002	Fork, White, Plastic, Bulk	1,000/cs
NON0420020	Fork, White, Plastic, Wrapped	1,000/cs
NON042000	Spoon, White, Plastic, Bulk, 5"/12.7 cm	2,000/cs
NON042201	Spoon, White, Plastic, Inner Boxes	1,000/cs
NON042001	Spoon, White, Plastic, Bulk	1,000/cs
NON0420015	Spoon, White, Plastic, Wrapped	1,000/cs
NON042003	Knife, White, Plastic, Bulk	1,000/cs
GEN4KITMW	Fork/Knife/Spoon/Napkin Cutlery Kit	250/bx

Heavy-Weight Utensils

NON00291	Fork, White, Polystyrene	1,000/cs
NON00293	Spoon, White, Polystyrene	1,000/cs
NON00292	Knife, White, Polystyrene	1,000/cs



Straws

- » Classic white design flex straws can bend or remain straight to best fit your needs
- » Individually wrapped

Item No.	Description	Pkg
NON02315	7¼"/19.7 cm, Wrapped, Flex, Mini Pack	400/bx, 4 bx/cs, 1,600 ea/cs
NON02325	7¼"/19.7 cm, Wrapped, Flex	400/bx, 25 bx/cs, 10,000 ea/cs
NON02350	10"/25.4 cm, Wrapped, Flex	400/bx, 4 bx/cs, 1,600 ea/cs



NONZIP

NON132000206

Bags

- » Durable, convenient and resealable plastic bags offer an economical solution for a variety of applications
- » Paper bags are available in assorted sizes and weights to suit your needs
- » Bags with white write-on strip offer a convenient space to write dates and other labels for easy identification

Item No.	Description	Pkg
----------	-------------	-----

Zip, White Write-On Block

NONZIP35	3 x 5"/7.6 x 12.7 cm, 2 mil	1,000/cs
NONZIP48	4 x 8"/10.2 x 20.3 cm, 2 mil	1,000/cs
NONZIP58	5 x 8"/12.7 x 20.3 cm, 2 mil	1,000/cs
NONZIP69	6 x 9"/15.2 x 22.9 cm, 2 mil	1,000/cs
NONZIP810	8 x 10"/20.3 x 25.4 cm, 2 mil	1,000/cs
NONZIP912	9 x 12"/22.9 x 30.5 cm, 2 mil	1,000/cs
NONZIP1012	10 x 12"/25.4 x 30.5 cm, 2 mil	1,000/cs
NONZIP1215	12 x 15"/30.5 x 38.1 cm, 2 mil	1,000/cs

White Paper, Duro

BAGGW2500	#2, 4 x 2 x 7"/10.2 x 5.1 x 17.8 cm	500/cs
BAGGW4500	#4, 5 x 3 x 9.5"/12.7 x 7.6 x 24.1 cm	500/cs
NON13200211	#6, 6 x 3.5 x 11"/15.2 x 8.9 x 27.9 cm	500/cs
BAGGW8500	#8, 6 x 4 x 12.5"/15.2 x 10.2 x 31.8 cm	500/cs
NON13200209	#12, 7 x 4.5 x 13.75"/17.8 x 11.4 x 34.9 cm	500/cs

Item No.	Description	Pkg
----------	-------------	-----

Brown Paper

NON13200340	#2, Duro, 4 x 2 x 7"/10.2 x 5.1 x 17.8 cm	500/cs
BAGGX2500	#2, Duro Hvy Duty, 4 x 2 x 7"/10.2 x 5.1 x 17.8 cm	500/cs
BAGGK4500	#4, Duro, 5 x 3 x 9.5"/12.7 x 7.6 x 24.1 cm	500/cs
NON13201082	#5, Duro, 5 x 3.5 x 11"/12.7 x 8.9 x 27.9 cm	500/cs
NON13200204	#6, Duro, 6 x 3.5 x 11"/15.2 x 8.9 x 27.9 cm	500/cs
NON132000206	#8, Duro, 6 x 4 x 12.5"/15.2 x 10.2 x 31.8 cm	500/cs
NON13200985	#10, Duro, 6.5 x 4 x 13.5"/16.5 x 10.2 x 34.3 cm	500/cs
NON13200982	#12, Duro, 7 x 4.5 x 14"/17.8 x 11.4 x 35.6 cm	500/cs
NON13200178	#16, Duro, 7.5 x 5 x 16"/17.8 x 12.7 x 40.6 cm	500/cs
NON13200177	#20, Duro, 8 x 5.5 x 16"/20.3 x 14 x 40.6 cm	500/cs
NON13200978	#25, Duro, 8 x 5.5 x 18"/20.3 x 14 x 45.7 cm	500/cs
NON132000595	#57, Duro, 12 x 7 x 17"/30.5 x 17.8 x 43.2 cm	500/cs
BAGSK1466SOST	#66, BBL, 12 x 7 x 23.5"/30.5 x 17.8 x 59.7 cm	250/cs
MDH13204281	#65, Duro, 13 x 7 x 16.5"/33 x 17.8 x 41.9 cm, Straw Handle	250/cs

Plastic, Thank You

NON11460395	Plastic, Grocery, Thank You	975/cs
-------------	-----------------------------	--------



Disposable Paper Buckets

- » An economical solution to plastic buckets while saving dish washing time

Item No.	Description	Pkg
----------	-------------	-----

Buckets

HUH71846	32-oz./946.4 ml, Bucket/Lid Combo, Paper, White	250/cs
DTM3TI	53-oz./1.6 L, Unwaxed, Paper, White	300/cs
DTM5T3N0196	83-oz./2.5 L, Waxed, Paper, White	100/cs
DTM5T1N0195	83-oz./2.5 L, Unwaxed, Paper, White	100/cs
DTM5T1GR0205	83-oz./2.5 L, Grease-Resist, Paper, White	100/cs
DTM10T3N0199	165-oz./4.9 L, Waxed, Paper, White	100/cs
DTM10T1N0198	165-oz./4.9 L, Unwaxed, Paper, White	100/cs

Lids for Buckets

DTM3VT19S	53-oz./1.6 L, Lid, Paper, White, Vent, Pull	600/cs
DTM5VT19SN11	83-oz./2.5 L, Lid, Paper, White, Vent, Pull	100/cs
DTM10V19S	165-oz./4.9 L, Lid, Paper, White, Vent, Pull	100/cs



DTM3VT19S



DTM3TI

KEURIG®

K-Cup®

- » Choose from a variety of hearty, delicious K-Cups
- » For use with Keurig® brewing systems



Item No.	Description	Pkg
K-Cup Pods		
GMT6082	Bigelow® Earl Grey Tea	24/bx
GMT6085	Bigelow Green Tea	24/bx
GMT6080	Bigelow English Breakfast Tea	24/bx
GMT14734	Celestial Seasonings® Natural Antioxidant Green Tea	24/bx
GMT14731	Celestial Seasonings English Breakfast Tea	24/bx
GMT14732	Celestial Seasonings Lemon Zinger® Tea	24/bx
GMT14738	Celestial Seasonings India Spice Chai	24/bx
GMT14737	Celestial Seasonings Decaf Natural Antioxidant Green Tea	24/bx
GMT14735	Celestial Seasonings Mandarin Orange Spice® Tea	24/bx
GMT6825	Celestial Seasonings Southern Sweet Perfect Iced Tea	22/bx
GMT6520	Green Mountain Coffee® Breakfast Blend	24/bx
GMT7522	Green Mountain Coffee Breakfast Blend Decaf	24/bx
GMT4067	Green Mountain Coffee Dark Magic® Decaf	24/bx
GMT6003	Green Mountain Coffee Colombian Fair Trade Select	24/bx
GMT6694	Green Mountain Coffee French Roast	24/bx
GMT6732	Green Mountain Coffee French Vanilla	24/bx
GMT6663	Green Mountain Coffee Nantucket Blend®	24/bx
GMT6792	Green Mountain Coffee Hazelnut	24/bx
GMT4050	Newman's Own® Organics Newman's Special Blend	24/bx
DIE60052101	The Original Donut Shop® Regular	24/bx
GMT6870	Snapple® Lemon Iced Tea	22/bx
GMT6872	Snapple Peach Iced Tea	22/bx
GMT6871	Snapple Raspberry Iced Tea	22/bx
GMT6622	Snapple Diet Peach Iced Tea	22/bx
GMT6113	Snapple Diet Half 'n Half	22/bx
GMT1252	Swiss Miss® Milk Chocolate Hot Cocoa	24/bx



GMT20155



GMT23145

KEURIG®

Brewers

- » Brews a perfect cup of coffee, tea, hot cocoa, or iced beverage in under one minute at the touch of a button
- » Now everyone in your office can brew what they love
- » All brewers are certified by the National Sanitation Foundation (NSF) and are UL-listed commercial brewers

Item No.	Description	Pkg
GMT23145	K145 Office Pro Keurig K-Cup® Brewer	1 ea/cs
GMT20143*	K150P Keurig K-Cup Brewer	1 ea/cs
GMT20155	K155 Keurig K-Cup Brewer	1 ea/cs
GMT23000*	K3000SE Keurig K-Cup Brewer	1 ea/cs
GMT5572	Filter Kit for Use with Keurig Commercial Brewers	1/ea

*Plumbed brewer, free installation

Coffee Condiments and Accessories

Item No.	Description	Pkg
Creamers		
NES35110BX	Nestle® Coffee-Mate® Creamer, Non-Dairy, 0.375-oz. Per Single Serve Cup, Original Flavor	50/bx
NES35170BX	Nestle Coffee-Mate Creamer, Non-Dairy, 0.375-oz. Per Single Serve Cup, French Vanilla Flavor	50/bx
ITD827976	Mini Moo's® Half and Half Creamer, 0.05-oz. Per Single Serve Cup	192/bx
OFX00020	Office® Snax Non-Dairy Creamer: 12-oz.	1/ea
NES30212	Nestle Coffee-Mate Original Powdered Creamer: 22-oz.	1/ea
NES31803	Nestle Coffee-Mate Coffee Creamer, Non-Dairy, Pump Dispenser, French Vanilla: 1.59-qt.	1/ea

Sugar

JOJ200411	Splenda® No Calorie Sweetener Packets, Single Serve	400/bx
TRU8845	Truvia® Natural Sugar Substitute Packets, Single Serve	140/bx
OFX20015445	Equal® Sweetener Packets, Zero Calorie, Single Serve	115/bx
SMU50150	Sweet N Low® Sugar Substitute Packets, Single Serve	400/bx
OFX00019	Office Snax® Pourable Sugar Canister: 20-oz.	1/ea



NES31803



OFX00019



universal.

Binders

- » Slant rings allow smooth page turning and hold more documents
- » Two inside pockets provide storage for loose documents
- » Features back mounted rings
- » Poly material is non-stick so inserts, pockets and cover won't lift print

Item No.	Description	Item No.	Description	Pkg
----------	-------------	----------	-------------	-----

Round Ring Economy View Binders

White		Black		
UNV20962	1" Capacity	UNV20961	1" Capacity	1/ea
UNV20952	½" Capacity,	UNV20951	½" Capacity	1/ea
UNV20972	1 ½" Capacity	UNV20971	1 ½" Capacity	1/ea
UNV20982	2" Capacity	UNV20981	2" Capacity	1/ea
UNV20992	3" Capacity	UNV20991	3" Capacity	1/ea

D-Ring Economy View Binders

White		Black		
UNV20742	1" Capacity	UNV20741	1" Capacity	1/ea
UNV20748	3" Capacity	UNV20745	2" Capacity	1/ea
UNV20994	4" Capacity	UNV20995	4" Capacity	1/ea
UNV20746	2" Capacity	UNV20747	3" Capacity	1/ea
UNV20744	1 ½" Capacity	UNV20743	1 ½" Capacity	1/ea



UNV43118

universal.

Office Basics

Item No.	Description	Pkg
UNV43623	Dry Erase Board, Melamine, 36 x 24", Satin-Finished Aluminum Frame	1/ea
UNV40304	Hardboard Clipboard, 1 ¼" Capacity, Holds 8 ½ x 11", Brown	1/ea
UNV40305	Hardboard Clipboard, 1 ¼" Capacity, Holds 8 ½ x 14", Brown	1/ea
UNV75602	Correction Tape with Two-Way Dispenser, Non-Refillable, ½" x 39.3'	2 ea/pk
UNV43663	Dry Erase Eraser, Synthetic Wool Felt, 5 x 1 ¾ x 1"	1/ea
UNV92009	Economy Scissors, 8" Length, Straight Handle, Stainless Steel, Black	1/ea
UNV79000	Standard Chisel Point 210 Strip Count Staples, 5,000/Box	5,000/bx
UNV43118	Economy Full Strip Stapler, 12-Sheet Capacity, Black	1/ea
UNV10700	Wand Style Staple Remover, Black	1/ea
UNV00700	Jaw Style Staple Remover, Black	1/ea
UNV83436	Invisible Tape, ¾" x 108', 1" Core, Clear	1/ea
UNV15001	Desktop Tape Dispenser, 1" Core, Weighted Nonskid Base, Black	1/ea
UNV83410	Invisible Tape, ¾" x 83.3', 1" Core, Clear	6 ea/pk

universal.

Clips and Rubberbands



UNV10220

Item No.	Description	Pkg
UNV10200	Small Binder Clips, Steel Wire, ¾" Capacity, ¾" Wide, Black/Silver, Dozen	12 ea/pk
UNV10210	Medium Binder Clips, Steel Wire, ¾" Cap., 1 ¼" Wide, Black/Silver, Dozen	12 ea/pk
UNV10220	Large Binder Clips, Steel Wire, 1" Capacity, 2" Wide, Black/Silver, Dozen	12 ea/pk
UNV72210	Paper Clips, Smooth Finish, No. 1, Silver, 1,000/Pack	1,000 ea/pk
UNV72220	Smooth Paper Clips, Wire, Jumbo, Silver, 100/Box, 10 Boxes/Pack	1,000 ea/pk
UNV10200VP3	Small Binder Clips, Steel Wire, ¾" Capacity, ¾" Wide, Black/Silver, 36/Pack	36/pk
UNV00133	Rubber Bands, Size 33, 3 ½ x ⅛", 640 Bands/1-lb. Pack	640/pk
UNV00164	Rubber Bands, Size 64, 3 1.2 x ¼", 320 Bands/1-lb. Pack	320/pk
UNV01117	Rubber Bands, Size 117, 7 x ⅛", 210 Bands/1-lb. Pack	210/pk

universal.

Paper Products



UNV11289

Item No.	Description	Pkg
File Folders		
UNV12113	File Folders, 1/3 Cut Assorted, One-Ply Top Tab, Letter, Manila	100 ea/bx
UNV14115	Hanging File Folders, 1/3 Tab, 11 Point Stock, Letter, Standard Green	25 ea/bx
UNV15113	File Folders, 1/3 Cut Assorted, One-Ply Top Tab, Legal, Manila	100 ea/bx
UNV14113	Hanging File Folders, 1/3 Tab, 11 Point Stock, Letter, Standard Green	25 ea/bx
UNV16113	File Folders, 1/3 Cut Assorted, Two-Ply Top Tab, Letter, Manila	100 ea/bx
UNV14215	Hanging File Folders, 1/3 Tab, 11 Point Stock, Legal, Standard Green	25 ea/bx
UNV10506	File Folders, 1/3 Cut Single-Ply Top Tab, Letter, Assorted Colors	100 ea/bx
Copy Paper		
MBFTC99004	Copy Paper, 92 Brightness, 20-lb., 8 1/2 x 11", White	5,000 ea/ct
UNV11289	Copy Paper Convenience Carton, 92 Brightness, 20-lb., 8 1/2 x 11", White	2,500 ea/ct
MBFTC990019	Copy Paper, 92 Brightness, 20-lb., 11 x 17", White	2,500 ea/ct
MBFTC990020	Copy Paper, 92 Brightness, 20-lb., 8 1/2 x 14", White	5,000 ea/ct
UNV95200	Multipurpose Paper, 98 Brightness, 20-lb., 8 1/2 x 11", Bright White	5,000 ea/ct
UNV20030	30% Recycled Copy Paper, 92 Brightness, 20-lb., 8 1/2 x 11", White	5,000 ea/ct

Item No.	Description	Pkg
Notebooks and Writing Pads		
UNV96920	Steno Book, Gregg Rule, 6 x 9", White, 80 Sheets	1/ea
UNV86920	Steno Book, Gregg Rule, 6 x 9", Green, 80 Sheets	1/ea
UNV20435	Wirebound Memo Books, Narrow Rule, 3 x 5", White, 12 50-Sheet Pads/Pack	12 ea/pk
UNV66600	Wirebound Notebook, 8-1/2 x 11", College Ruled, 100 Sheets, Assorted Color Cover	1/ea
UNV20630	Perforated Edge Writing Pad, Legal Ruled, Letter, White, 50-Sheet, Dozen	12 ea/pk
UNV10630	Perforated Edge Writing Pad, Legal/Margin Rule, Letter, Canary, 50-Sheet, Dozen	12 ea/pk
UNV46300	Perforated Edge Writing Pad, Narrow Rule, 5 x 8", White, 50-Sheet, Dozen	12 ea/pk
UNV46200	Perforated Edge Writing Pad, Narrow Rule, 5 x 8", Canary, 50-Sheet, Dozen	12 ea/pk
UNV40000	Perforated Edge Writing Pad, Legal/Margin Rule, Legal, Canary, 50-Sheet, Dozen	12 ea/pk

Envelopes		
UNV35264	Kraft Clasp Envelopes, Side Seam, 28-lb., 9 x 12", Light Brown	100 ea/bx
UNV35210	Business Envelope, V-Flap, #10, 24-lb., White	500 ea/bx
UNV35267	Kraft Clasp Envelope, Side Seam, 28-lb., 10 x 13", Light Brown	100 ea/bx
UNV36301	Double Window Check Envelope, #9, White	500 ea/bx

Self-Stick Notes		
UNV35668	Standard Self-Stick Notes, 3 x 3, Yellow, 12 100-Sheet Pads/Pack	12 ea/pk
UNV35662	Standard Self-Stick Notes, 1 1/2 x 2, Yellow, 12 100-Sheet Pads/Pack	12 ea/pk
UNV35672	Standard Self-Stick Notes, 3 x 5, Yellow, 12 100-Sheet Pads/Pack	12 ea/pk

universal.
Writing Supplies



UNV15510

Item No.	Description	Pkg
UNV27410	Economy Ballpoint Stick Oil-Based Pen, Black Ink, Medium	12 ea/pk
UNV15510	Comfort Grip Ballpoint Retractable Pen, Black Ink, Medium	12 ea/pk
UNV08851	Pocket Clip Highlighter, Chisel Tip, Fluorescent Yellow Ink	12 ea/pk
UNV08850	Pocket Highlighter, Chisel Tip, Fluorescent Colors	5/st
UNV08860	Desk Highlighter, Chisel Tip, Fluorescent Colors	5/st
UNV07051	Permanent Markers, Chisel Tip, Black	12 ea/pk
UNV07071	Pen Style Permanent Markers, Fine Point, Black	12 ea/pk
UNV43651	Dry Erase Marker, Chisel Tip, Black	12 ea/pk
UNV43671	Dry Erase Marker, Bullet Tip, Black	12 ea/pk
UNV55400	Economy Woodcase Pencil, HB #2, Yellow Barrel	12 ea/pk
UNV08861	Desk Highlighter, Chisel Tip, Fluorescent Yellow	12 ea/pk

- 3M is a registered trademark of 3M Company
- 3M Littmann is a registered trademark of 3M Company
- A1CNow+ is a registered trademark of Bayer Healthcare LLC
- Aaron 940 is a trademark of Bovie Medical Corporation
- AccuCare is a registered trademark of Medline Industries, Inc.
- Accu-Therm is a registered trademark of Medline Industries, Inc.
- Accutouch is a trademark of Tillotson Healthcare Corporation
- Activation-Assist is a trademark of Becton, Dickinson and Company
- Adacel is a registered trademark of Sanofi Pasteur Limited/ Sanofi Pasteur Limitee
- Alere Acceava is a registered trademark of Inverness Medical - Biostar Inc.
- Alere is a trademark of Alere Inc.
- Aloetouch Ice is a registered trademark of Medline Industries, Inc.
- Angel Soft ps is a registered trademark of Georgia-Pacific Corporation
- ANGEL WING is a registered trademark of Becton, Dickinson and Company
- Aplisol is a registered trademark of Parke, Davis & Company
- Aquaconductive is a trademark of Medline Industries, Inc.
- Autoguard is a trademark of Becton, Dickinson and Company
- AVE is a registered trademark of Medline Industries, Inc.
- B Braun is a registered trademark of B. Braun Melsungen Aktiengesellschaft CORPORATION FED REP GERMANY
- Baxter INTERLINK is a registered trademark of Baxter International Inc.
- BD Veritor is a registered trademark of Becton, Dickinson and Company
- Bigelow is a registered trademark of R. C. Bigelow, Inc.
- BigFold Z is a registered trademark of Georgia-Pacific Consumer Products Lp
- BinaxNOW is a registered trademark of Binax, Inc.
- BioCon is a registered trademark of Mcubetechnology Co., Ltd.
- Biosyn is a trademark of Covidien LP Covidien Holding Inc.
- Bladderpoint is a trademark of Medline Industries, Inc.
- boardwalk is a registered trademark of Essendant Co.
- Bolt is a registered trademark of Tampa's Team, LLC DBA Buddy Brew Coffee
- Boostrix is a registered trademark of SmithKline Beecham Biologicals S.A.
- Bovie is a registered trademark of Bovie Medical Corporation
- Braun is a registered trademark of Braun GmbH CORPORATION FED REP GERMANY
- Bulkee II is a registered trademark of Medline Industries, Inc.
- burton is a registered trademark of Burton Medical Products Corporation
- CardinalHealth is a trademark of Cardinal Health Technologies, LLC
- CardioChek is a registered trademark of Polymer Technology Systems, Inc.
- CardioChek Plus Smart Bundle is a trademark of Polymer Technology Systems, Inc.
- CaviCide is a registered trademark of Microecology Products Corp.
- CaviWipes is a registered trademark of Metrex Research LLC
- Celestial Seasonings Lemon Zinger is a registered trademark of Celestial Seasonings, Inc.
- Celestial Seasonings Mandarin Orange Spice is a registered trademark of Celestial Seasonings, Inc.
- Celestial Seasonings is a registered trademark of Celestial Seasonings, Inc.
- Cepheid is a registered trademark of Cepheid
- Chek-Stix is a registered trademark of Bayer Healthcare LLC
- Cidex is a registered trademark of Arbrook Manufacturing Corporation
- Clearview is a registered trademark of Inverness Medical Switzerland GmbH
- Clinitek is a registered trademark of Bayer Corporation
- CLINTON is a trademark of Clinton Industries, Inc.
- CoaguChek is a registered trademark of Roche Diagnostics GmbH, LLC
- Coban is a registered trademark of 3M Company
- Coffee-Mate is a registered trademark of Societe Des Produits Nestle S.A. Corporation Switzerland
- CoFlex is a registered trademark of Andover Healthcare, Inc.
- Colpo-Master is a trademark of Bovie Medical Corporation
- ComfortEase is a trademark of Medline Industries, Inc.
- ComfortMed is a registered trademark of Medline Industries, Inc.
- Compact is a registered trademark of James River Paper Company, Inc.
- CompuGroup is a trademark of Tremco Consultants, Inc.
- Conform is a registered trademark of The Kendall Company
- Covidien is a registered trademark of Covidien AG CORPORATION SWITZERLAND
- Covidien-Medtronic is a registered trademark of Covidien LP Covidien Holding Inc.
- CryOmega is a registered trademark of CryoConcepts LP STC Consulting LLC
- CUBEScan is a registered trademark of Mcubetechnology Co., Ltd.
- CURAD is a registered trademark of Medline Industries, Inc.
- DCA Vantage is a registered trademark of Siemens Medical Solutions Diagnostics

- Defibtech is a registered trademark of Defibtech, LLC
- DermaBlade is a registered trademark of Hasting's Biomedical Corporation
- Dermabond is a registered trademark of Johnson & Johnson Corporation New Jersey
- DermAssure is a trademark of Derma Gator, LLC
- Determine is a trademark of Alere Switzerland GmbH
- Duracell EasyTab is a registered trademark of The Gillette Company Corporation
- Duracell Procell is a registered trademark of The Gillette Company Corporation
- Durapore is a registered trademark of 3M Company
- Eclipse is a trademark of Becton, Dickinson and Company
- ECO Products is a registered trademark of Eco-Products, Inc.
- eGLU is a registered trademark of Polymer Technology Systems, Inc.
- Endocell is a registered trademark of CooperSurgical, Inc.
- Engerix-B is a registered trademark of GlaxoSmithKline Biologicals S.A. CORPORATION BELGIUM
- enMotion Impulse is a registered trademark of Georgia-Pacific Consumer Products LP
- enMotion Tranquil Aloe is a registered trademark of Georgia-Pacific Consumer Products LP
- enMotion is a registered trademark of Georgia-Pacific Consumer Products LP
- enthermics is a registered trademark of Enthermics, Inc.
- Envirocide is a registered trademark of Viro Research International, Inc.
- Envision is a registered trademark of Georgia-Pacific Consumer Products LP
- Enzol is a registered trademark of Johnson & Johnson Corporation
- Equal is a registered trademark of Merisant Company Corporation
- Ethicon is a registered trademark of ETHICON, INC.
- Ethilon is a registered trademark of ETHICON, INC.
- EVENCARE is a registered trademark of Medline Industries, Inc.
- EvenCare G2 is a registered trademark of Medline Industries, Inc.
- Excel is a registered trademark of Medline Industries, Inc.
- Exuderm is a registered trademark of Medline Industries, Inc.
- E-Z-Disk is a trademark of E-Z-EM Canada Inc.
- E-Z-Paque is a registered trademark of E-Z-EM Canada Inc.
- EZ Access is a registered trademark of Adalis Corporation
- FastStart is a trademark of Masimo Corporation
- FirmFlex is a trademark of Integra LifeSciences Corporation
- First Aid Only is a registered trademark of Acme United Corporation
- FitGuard is a registered trademark of Medline Industries, Inc.
- Flex-Fabric is a trademark of Medline Industries, Inc.
- Fourflex is a registered trademark of Medline Industries, Inc.
- Freezpoint is a trademark of Bovie Medical Corporation
- Fundus Sound is a trademark of Medline Industries, Inc.
- G2 is a registered trademark of Medline Industries, Inc.
- Generation Pink is a registered trademark of Medline Industries, Inc.
- Germ Shield is a trademark of Medline Industries, Inc.
- Green Mountain Coffee Dark Magic is a registered trademark of Keurig Green Mountain, Inc.
- Green Mountain Coffee Nantucket Blend is a registered trademark of Keurig Green Mountain, Inc.
- Green Mountain Coffee is a registered trademark of Keurig Green Mountain, Inc.
- Green Tree is a registered trademark of Medline Industries, Inc.
- GreenSmart is a registered trademark of Medline Industries, Inc.
- Guardian is a registered trademark of Medline Industries, Inc.
- Halogen HPX is a registered trademark of Welch Allyn, Inc.
- Harmonic is a registered trademark of Johnson & Johnson Corporation
- Harmony is a registered trademark of Medline Industries, Inc.
- Havrix is a registered trademark of SmithKline Biologicals, S.A.
- Health o Meter is a registered trademark of Sunbeam Products, Inc.
- Hemocult is a registered trademark of Schieffelin & Co. Corporation
- HEMOCUE is a registered trademark of Hemocue AB
- HemoPoint is a registered trademark of STANBIO Laboratory LLP
- Hemosure is a registered trademark of Hemosure Inc.
- Hospira is a registered trademark of Hospira Inc. Corporation Delaware
- iCup is a registered trademark of Instant Technologies, Inc.
- Infanrix is a registered trademark of SmithKline Beecham P.L.C.
- Insyte is a trademark of Becton, Dickinson and Company
- Integra is a registered trademark of Integra LifeSciences Corporation
- Integra is a registered trademark of Integra LifeSciences Corporation
- IPOL is a registered trademark of Connaught Laboratories Limited Corporation Canada
- IQmanager is a registered trademark of Midmark Corporation

iScreen is a registered trademark of Hospital Central Services, Inc. DBA Miller Keystone Blood Center non-pr
 Isolex is a trademark of Baxter International Inc.
 ivNOW is a registered trademark of Enthermics, Inc.
 Johnson & Johnson is a registered trademark of Johnson & Johnson Corporation
 K-Cup is a registered trademark of Keurig Green Mountain, Inc.
 Kendall Conform is a registered trademark of The Kendall Company
 Kerlix is a registered trademark of Tyco Healthcare Group LP
 Keurig K-Cup is a registered trademark of Keurig Green Mountain, Inc.
 Keurig Office Pro is a registered trademark of Keurig Green Mountain, Inc.
 Kinrix is a registered trademark of GlaxoSmithKline Biologicals S.A. CORPORATION BELGIUM
 KleenSpec is a registered trademark of Welch Allyn, Inc.
 König is a registered trademark of Medline Industries, Inc.
 Laboratory solutions made simple is a trademark of Quantimetrix Corporation
 LeadCare is a registered trademark of ANDCARE, Inc.
 Lifeline View AED is a registered trademark of Defibtech, LLC
 Littmann Classic II is a trademark of 3M Company
 Littmann Classic III is a trademark of 3M Company
 Low Perfusion is a trademark of Masimo Corporation
 Luer-Lok is a registered trademark of Becton, Dickinson and Company
 Marathon is a registered trademark of Medline Industries, Inc.
 Masimo SET is a registered trademark of Masimo Corporation
 Masimo Pronto is a registered trademark of Masimo Corporation
 Matrix is a registered trademark of Medline Industries, Inc.
 maxOat+ is a registered trademark of Medline Industries, inc.
 Maxorb is a registered trademark of Medline Industries, Inc.
 Medfix is a trademark of MedFix International, LLC
 MedGel is a registered trademark of Medline Industries, Inc.
 Medigrip is a trademark of Medline Industries, Inc.
 Mediplast is a registered trademark of Medline Industries, Inc.
 Medi-Strips is a registered trademark of Medline Industries, Inc.
 Medline Compli-Mates is a registered trademark of Medline Industries, Inc.
 Medtoons is a registered trademark of Medline Industries, Inc.
 Menactra is a registered trademark of Connaught Technology Corporation

Metrex is a registered trademark of Metrex Research, LLC
 MetriCide is a registered trademark of Metrex Research Corporation
 MetriTest is a trademark of Metrex Research, LLC
 MicroCoat is a registered trademark of American Safety Razor Company
 Micro-Kill is a trademark of Medline Industries, Inc.
 Micropore is a registered trademark of 3M Company
 MicroScaffold is a trademark of Medline Industries, Inc.
 Midmark IQspiro is a registered trademark of Midmark Corporation
 MIDMARK is a registered trademark of Midmark Corporation
 Miltex is a registered trademark of MILTEX, INC.
 Mini Moo's is a registered trademark of Land O'Lakes, Inc.
 Monosof is a trademark of Covidien LP Covidien Holding Inc.
 Multistix PRO is a registered trademark of Bayer Corporation
 Napkleen is a registered trademark of NAPKLEEN LLC
 Neolon is a registered trademark of Medline Industries, Inc.
 Nestle is a registered trademark of Societe Des Produits Nestle S.A. Corporation Switzerland
 Newman's Own is a registered trademark of No Limit, LLC
 Office Snax is a registered trademark of Office Snax, Inc.
 On the Right Track is a registered trademark of On the Right Track Systems, Inc.
 OneTouch WheelBase is a registered trademark of United Metal Fabricators, Inc. DBA UMF Medical CORPORATION
 Opticell is a registered trademark of Optical Skin Services, Inc.
 Orbit is a trademark of QRS Diagnostic
 OSOM BVBLUE is a registered trademark of Gryphus Diagnostics, LLC
 OSOM is a registered trademark of Wyntek Diagnostics, Inc.
 Outpatient is a trademark of Burton Medical Products Corporation
 PassPort is a registered trademark of SPS Medical Supply Corp.
 Pediarix is a registered trademark of Smithkline Beecham Biologicals S.A. Corporation Belgium
 PedvaxHIB is a registered trademark of Merck & Co., Inc. Corporation New Jersey
 Pentacel is a registered trademark of Sanofi Pasteur Limited/Sanofi Pasteur Limitee Corporation Canada
 Performs Multistix is a registered trademark of Miles Laboratories, Inc.
 Physio-Control LIFEPAK CR is a registered trademark of Medtronic Emergency Response Systems, Inc.
 Physio-Control LIFEPAK Express is a registered trademark of Medtronic Emergency Response Systems, Inc.

Physio-Control LIFEPAK is a trademark of Medtronic Emergency Response Systems, Inc.

POCKETPEAK is a registered trademark of Ferraris Medical Inc.

Polysorb is a trademark of Covidien LP Covidien Holding Inc.

PrecisionGlide is a registered trademark of Becton, Dickinson and Company

PrecisionGlide is a trademark of Becton, Dickinson and Company

Preference is a registered trademark of Fort Howard Corporation

Prevnar 13 is a registered trademark of Wyeth LLC Limited Liability Company Delaware

Prolene is a registered trademark of ETHICON, INC.

Protection Plus is a registered trademark of Medline Industries, Inc.

Puracol is a registered trademark of Medline Industries, Inc.

Q200/HE is a trademark of QRS Diagnostic

Qualigen is a registered trademark of Qualigen, Inc.

Quantimetrix is a registered trademark of Quantimetrix Corporation

Quantum is a registered trademark of M&M Plastics, Inc.

QuickVue is a registered trademark of Quidel Corporation

QUIDEL is a registered trademark of Quidel Corporation

Rapicide is a registered trademark of Minntech Corporation

Red Dot is a registered trademark of Guardian Products, Inc.

Restore Sense is a registered trademark of Medline Industries, Inc.

Restore is a registered trademark of Medline Industries, Inc.

Rice Lake is a registered trademark of Rice Lake Weighing Systems, Inc.

Ritter is a trademark of Midmark Corporation

RotaTeq is a registered trademark of Merck & Co., Inc. Corporation New Jersey

SafetyGlide is a registered trademark of Becton, Dickinson and Company

Safety-Lok is a registered trademark of Becton, Dickinson and Company

seca is a registered trademark of Vogel & Halke Partnership Fed Rep Germany

SensiCare Ice is a registered trademark of Medline Industries, Inc.

SensiCare Silk is a registered trademark of Medline Industries, Inc.

SensiCare is a registered trademark of Medline Industries, Inc.

SilvaSorb is a registered trademark of ACRYMED, Inc. Corporation Oregon

Skintegrity is a registered trademark of Medline Industries, Inc.

SmartGuard is a registered trademark of Medline Industries, Inc.

SmartGuard is a trademark of Medline Industries, Inc.

SmartSite is a registered trademark of IVAC Holdings, Inc.

Smoke Shark is a trademark of Bovie Medical Corporation

Snapple is a registered trademark of Snapple Beverage Corp.

SoftFlex is a trademark of Integra LifeSciences Corporation

Soft-Wrap is a registered trademark of Medline Industries, Inc.

Spectrum is a trademark of Medline Industries, Inc.

Splenda is a registered trademark of Heartland Consumer Products LLC

Spring Frost is a registered trademark of Medline Industries, Inc.

STANBIO LABOATORY is a trademark of STANBIO Laboratory LLP

STAR is a trademark of Bernal Incorporated

Steri-Strip is a registered trademark of 3M Company

Stratasorb is a registered trademark of Medline Industries, Inc.

Super Exam 50 is a registered trademark of Burton Medical Products Corporation

Sureprep is a registered trademark of Medline Industries, Inc.

SureSite IV is a trademark of Medline Industries, Inc.

Suresite is a registered trademark of Medline Industries, Inc.

Sure-Wrap is a registered trademark of Medline Industries, Inc.

Surgilube is a registered trademark of Day Chemical Company Corporation

Surgipro is a trademark of United States Surgical Corporation

Sweet N Low is a registered trademark of Cumberland Packing Corp.

Swift-Wrap is a registered trademark of Medline Industries, Inc.

Swiss Miss is a registered trademark of ConAgra Foods RDM, Inc.

Synergy Classic is a registered trademark of Medline Industries, Inc.

Telfa is a registered trademark of Tyco Healthcare Group LP

Tenivac is a registered trademark of Sanofi Pasteur Limited/Sanofi Pasteur Limitee

The Grabber is a trademark of On the Right Track Systems, Inc.

The Original Donut Shop is a registered trademark of Keurig Green Mountain, Inc.

TheraHoney is a registered trademark of Medline Industries, Inc.

Thermoscan is a registered trademark of THERMOSCAN INC.

Threeflex is a registered trademark of Medline Industries, Inc.

Transpore is a registered trademark of 3M Company

Triumph LT is a registered trademark of Medline Industries, Inc.

Triumph is a registered trademark of Medline Industries, Inc.

Tru-Loc is a trademark of Draper Corporation Corporation

Truvia is a registered trademark of Cargill, Incorporated

Tubersol is a registered trademark of Sanofi Pasteur Limited/Sanofi Pasteur Limitee

Twinrix is a registered trademark of Smithkline Beecham Biologicals S.A Company Belgium

Typhim VI is a registered trademark of Connaught Technology Corporation

UltiCare is a registered trademark of UltiMed Inc.

UMFmedical is a trademark of United Metal Fabricators, Inc. DBA UMF Medical CORPORATION PENNSYLVANIA

Universal ECG is a trademark of QRS Diagnostic

Universal is a registered trademark of The Anaconda Company

Unna-Z is a registered trademark of Medline Industries, Inc.

Uristix is a registered trademark of AMES COMPANY, INC.

Urisys 1100 is a registered trademark of Roche Diagnostics GmbH LIMITED LIABILITY COMPANY FED REP GERMANY

VACUETTE is a registered trademark of Greiner Bio-One GmbH

Vacutainer is a registered trademark of Becton, Dickinson and Company

Valueklave is a registered trademark of TUTTNAUER USA CO. LTD.

Versatel is a registered trademark of Medline Industries, Inc.

Vialon is a trademark of Becton, Dickinson and Company

Vicryl is a registered trademark of Johnson & Johnson Corporation

Vivotif is a registered trademark of PAXVAX, INC.

WarmRight is a registered trademark of Enthermics Medical Systems Inc.

Welch Allyn SureTemp is a registered trademark of Diatek Instruments, Inc.

WheelBase is a registered trademark of United Metal Fabricators, Inc. DBA UMF Medical CORPORATION PENNSYLVANIA

Windows is a registered trademark of Microsoft Corporation

ZR is a registered trademark of Excelsior Medical Corporation

Item Number Index

0002-8031-01.....	126	0024-2792-10.....	126	0338-0644-04.....	130
0002-8215-01.....	127	0025-1520-31.....	132	0338-1073-02.....	126
0002-8315-01.....	127	0054-0017-25.....	134	0378-0315-93.....	134
0002-8715-01.....	127	0054-3117-63.....	132	0378-0353-01.....	134
0003-0293-05.....	127	0054-3177-63.....	132	0378-1454-01.....	133
0003-0293-20.....	127	0054-3505-47.....	135	0378-1610-01.....	133
0003-0293-28.....	127	0054-4184-25.....	133	0378-6987-89.....	123
0003-0494-20.....	127	0054-4299-25.....	133	0386-0004-04.....	135
0003-0526-20.....	131	0054-4741-25.....	134	0386-0008-03.....	136
0004-0191-09.....	124	0054-4742-25.....	134	0409-0805-01.....	125
0004-0800-85.....	134	00591-5543-01.....	132	0409-1036-30.....	125
0004-0801-85.....	134	0065-0359-05.....	131	0409-1038-50.....	125
0004-6316-01.....	133	0065-0741-14.....	132	0409-1041-30.....	125
0005-0100-05.....	137	0065-0795-15.....	131	0409-1144-02.....	130
0005-0100-10.....	137	0065-0795-50.....	131	0409-1144-05.....	130
0005-1971-02.....	137	0065-8533-02.....	135	0409-1159-01.....	124
0005-1971-05.....	137	0069-0058-01.....	126	0409-1159-02.....	124
0006-0717-54.....	133	0071-0418-13.....	134	0409-1160-01.....	124
0006-3843-71.....	127	0078-0180-01.....	129	0409-1162-01.....	124
0006-4047-20.....	137	0078-0359-34.....	133	0409-1162-02.....	124
0006-4047-41.....	137	0078-0435-61.....	129	0409-1163-01.....	124
0006-4093-02.....	137	0093-0752-01.....	132	0409-1171-01.....	126
0006-4095-02.....	138	0093-2268-01.....	132	0409-1207-03.....	126
0006-4096-02.....	138	0093-3125-01.....	133	0409-1215-01.....	129
0006-4119-03.....	137	0093-3145-01.....	132	0409-1219-01.....	129
0006-4121-02.....	137	0093-4029-01.....	133	0409-1323-05.....	128
0006-4171-00.....	137	0096-0707-35.....	135	0409-1412-04.....	124
0006-4681-00.....	137	0096-0707-37.....	135	0409-1412-10.....	124
0006-4827-00.....	138	0096-0707-60.....	135	0409-1443-04.....	127
0006-4831-41.....	138	0121-0759-08.....	134	0409-1464-01.....	128
0006-4837-03.....	137	0132-0202-20.....	131	0409-1467-01.....	128
0006-4841-00.....	138	0143-2112-50.....	133	0409-1559-10.....	128
0006-4841-41.....	138	0143-3142-05.....	133	0409-1559-30.....	128
0006-4897-00.....	137	0143-3142-50.....	133	0409-1560-10.....	128
0006-4943-00.....	137	0143-9738-01.....	134	0409-1560-29.....	128
0006-4963-00.....	138	0143-9784-10.....	126	0409-1582-10.....	128
0006-4963-41.....	138	0143-9785-10.....	130	0409-1582-29.....	128
0006-4992-00.....	137	0143-9857-25.....	125	0409-1587-50.....	128
0008-4990-20.....	130	0143-9887-01.....	132	0409-1610-50.....	128
0009-0003-02.....	130	0143-9890-01.....	129	0409-1631-10.....	124
0009-0013-05.....	130	0143-9899-91.....	126	0409-1639-10.....	126
0009-0013-06.....	130	0143-9923-90.....	125	0409-1746-10.....	128
0009-0018-20.....	130	0143-9924-90.....	125	0409-1746-30.....	128
0009-0039-28.....	130	0168-0004-80.....	136	0409-1749-10.....	128
0009-0047-22.....	130	0168-0006-15.....	136	0409-1749-29.....	128
0009-0271-01.....	125	0168-0040-15.....	135	0409-1752-50.....	128
0009-0274-01.....	125	0168-0080-31.....	135	0409-1755-50.....	128
0009-0280-02.....	125	0168-0204-37.....	135	0409-1761-02.....	128
0009-0280-03.....	125	0168-0258-15.....	135	0409-1778-05.....	128
0009-0280-51.....	125	0169-1833-11.....	129	0409-1778-35.....	128
0009-0306-02.....	125	0169-1834-11.....	129	0409-1782-69.....	129
0009-0555-02.....	128	0169-1837-11.....	129	0409-1902-01.....	129
0009-0626-01.....	125	0173-0682-20.....	123	0409-1903-01.....	129
0009-0746-30.....	125	0185-0074-01.....	133	0409-1932-02.....	123
0009-0775-26.....	125	0185-0122-01.....	134	0409-2066-05.....	128
0009-0825-01.....	130	0185-0674-01.....	133	0409-2122-01.....	126
0009-3073-01.....	125	0186-1990-04.....	123	0409-2122-02.....	126
0009-3073-03.....	125	0187-0746-33.....	133	0409-2267-20.....	127
0009-3475-01.....	125	0223-1721-01.....	134	0409-2267-54.....	127
0009-3475-03.....	125	0228-2530-10.....	134	0409-2285-05.....	128
0009-4709-13.....	125	0235-0378-82.....	135	0409-2287-31.....	127
0023-0106-05.....	131	0281-0326-30.....	136	0409-2290-31.....	126
0023-3240-15.....	131	0281-0326-60.....	136	0409-2339-34.....	127
0023-4554-30.....	131	0283-0871-31.....	135	0409-2344-02.....	126, 133
0023-9163-30.....	131	0283-1185-20.....	135	0409-2347-32.....	126

0409-2720-01.....	126	0409-7337-01.....	125	0641-0376-25.....	126
0409-2720-03.....	126	0409-7337-01H.....	125	0641-0928-25.....	129
0409-2721-01.....	126	0409-7338-01.....	125	0641-1397-35.....	125
0409-2723-01.....	127	0409-7517-16.....	126	0641-1495-35.....	128
0409-3177-01.....	127	0409-7620-03.....	127	0641-1496-35.....	129
0409-3178-01.....	127	0409-7650-62.....	127	0641-2569-44.....	126
0409-3178-02.....	127	0409-7793-62.....	127	0641-6013-10.....	126
0409-3178-03.....	127	0409-8062-01.....	123	0641-6022-25.....	126
0409-3182-01.....	128	0409-9042-01.....	124	0641-6082-25.....	129
0409-3182-02.....	128	0409-9043-01.....	124	0641-6142-01.....	129
0409-3182-03.....	128	0409-9046-01.....	124	0713-0280-31.....	135
0409-3183-01.....	128	0409-9104-20.....	126	0743-0042-11.....	136
0409-3308-03.....	123	0409-9137-05.....	127	0781-1123-01.....	134
0409-3414-01.....	128	0409-9157-01.....	130	0781-1123-05.....	134
0409-3578-01.....	130	0409-9158-01.....	130	0781-2020-01.....	132
0409-3793-01.....	127	0409-9631-04.....	126	0781-2613-01.....	132
0409-3795-01.....	127	0469-6501-89.....	127	0781-2613-05.....	132
0409-3795-01H.....	127	0472-0166-30.....	136	0781-3206-85.....	125
0409-3796-01.....	127	0487-5901-99.....	123	0781-3224-64.....	126
0409-3796-01H.....	127	0487-9301-03.....	123	0781-5061-01.....	132
0409-4011-01.....	130	0487-9501-25.....	122	0781-5710-01.....	134
0409-4276-01.....	127	0487-9501-60.....	122	0781-6041-46.....	132
0409-4276-02.....	127	0487-9801-30.....	123	0781-6135-95.....	129
0409-4277-01.....	128	0487-9904-25.....	123	0869-2452-10.....	132
0409-4277-02.....	128	0496-0882-15.....	135	0884-0449-04.....	136
0409-4278-01.....	127	0496-0882-30.....	135	0884-6297-30.....	131
0409-4279-02.....	127	0496-0883-30.....	136	0904-5309-09.....	133
0409-4282-01.....	128	0501-3203-02.....	135	0904-5854-40.....	133
0409-4282-02.....	128	0501-3712-05.....	136	0904-5986-60.....	134
0409-4283-01.....	128	0501-3798-01.....	136	0904-7915-40.....	133
0409-4332-01.....	130	0517-0031-25.....	125	10019-360-40.....	123
0409-4646-01.....	129	0517-0032-25.....	125	10019-360-60.....	123
0409-4688-02.....	126	0517-0130-054.....	125	10019-641-34.....	123
0409-4713-02.....	127	0517-0720-01.....	124	10019-651-64.....	123
0409-4713-32.....	127	0517-1010-25.....	124	10019-955-01.....	125
0409-4755-03.....	129	0517-2340-10.....	130	10223-201-03.....	135
0409-4759-01.....	129	0517-4010-01.....	129	10481-3008-1.....	136
0409-4777-02.....	125	0517-4810-25.....	129	10702-007-01.....	132
0409-4900-34.....	130	0517-5601-25.....	127	1254-7237-79.....	136
0409-4902-34.....	126	0521-2326-00.....	136	12843-158-59.....	133
0409-4903-34.....	128	0527-1050-01.....	132	13533-131-01.....	137
0409-4910-34.....	124	0527-1343-01.....	133	16477-505-01.....	133
0409-4911-34.....	124	0536-4078-10.....	134	16729-108-11.....	128
0409-4916-34.....	129	0536-4406-01.....	134	16729-116-38.....	128
0409-4921-34.....	126	0536-5150-26.....	136	16729-183-01.....	133
0409-4928-34.....	125	0536-7180-01.....	132	17478-045-32.....	132
0409-5534-34.....	130	0548-5400-00.....	128	17478-064-12.....	131
0409-5820-01.....	126	0548-5701-00.....	128	17478-070-35.....	131
0409-5921-01.....	123	0555-0572-02.....	133	17478-100-02.....	131
0409-5922-01.....	123	0555-9027-42.....	133	17478-100-12.....	131
0409-6102-02.....	126	0555-9131-67.....	134	17478-101-12.....	132
0409-6102-04.....	126	0574-0601-15.....	136	17478-102-12.....	132
0409-6102-10.....	126	0574-0850-05.....	126	17478-107-20.....	127
0409-6533-01.....	130	0574-0850-10.....	126	17478-107-30.....	127
0409-6625-02.....	130	0574-4022-35.....	131	17478-215-05.....	131
0409-6629-02.....	129	0591-2888-30.....	127	17478-223-1.....	131
0409-6637-34.....	130	0591-3128-79.....	129	17478-263-12.....	131
0409-6651-06.....	129	0591-3222-47.....	129	17478-284-35.....	131
0409-6653-05.....	129	0591-5544-01.....	132	17478-404-01.....	131
0409-6695-01.....	123	0591-5708-01.....	132	17478-503-05.....	123
0409-6695-02.....	123	0603-2116-21.....	132	17478-504-10.....	128
0409-6729-09.....	128	0603-2959-21.....	132	17478-538-02.....	129
0409-7332-01.....	125	0603-5337-21.....	134	17478-542-02.....	123
0409-7334-10.....	125	0603-5338-21.....	134	17478-701-02.....	127
0409-7335-03.....	125	0641-0367-25.....	126	17478-703-11.....	131

Item Number Index

17478-711-10.....	135	49502-501-02.....	126	52584-584-93.....	127
17478-711-30.....	135	49884-093-03.....	133	52584-637-34.....	130
17478-713-10.....	131	50111-333-01.....	133	52584-656-61.....	132
17478-716-10.....	131	50111-334-01.....	133	52584-664-01.....	126
17478-931-01.....	124	50242-041-64.....	125	52584-902-34.....	126
18657-117-04.....	127	50242-044-13.....	123	52584-921-01.....	123
23155-294-42.....	129	50383-042-24.....	134	53746-190-10.....	133
24208-002-02.....	130	50383-584-04.....	133	53746-272-01.....	134
24208-580-60.....	131	50383-667-30.....	135	54162-011-01.....	135
24208-631-10.....	131	50383-741-20.....	122	54288-103-10.....	126
24208-780-55.....	131	50383-775-04.....	135	55150-114-20.....	124
24208-920-64.....	134	50383-901-10.....	135	55150-237-01.....	126
24338-126-13.....	133	50383-917-50.....	134	5542 STYLE.....	61
24385-003-46.....	135	50383-933-55.....	135	55566-1501-1.....	129
25021-106-10.....	125	50419-188-01.....	121	55566-4100-1.....	126
25021-109-10.....	125	50419-188-05.....	121	5570 STYLE.....	61
25021-110-20.....	125	50419-188-11.....	121	5582 STYLE.....	61
25021-115-02.....	125	50419-188-15.....	121	57664-477-52.....	133
25021-115-06.....	125	50419-188-36.....	121	57894-030-01.....	129
25021-121-20.....	125	50419-188-37.....	121	58160-810-11.....	137
31722-534-01.....	133	50419-188-38.....	121	58160-810-52.....	137
31722-713-90.....	134	50419-188-58.....	121	58160-811-52.....	137
39822-0190-1.....	123	50419-320-05.....	120	58160-812-11.....	137
400 STYLE.....	61	50419-325-27.....	120	58160-812-52.....	137
42023-104-01.....	136	50419-325-28.....	120	58160-815-52.....	138
42023-104-05.....	136	50419-325-29.....	120	58160-818-11.....	137
42023-111-01.....	125	50419-342-10.....	122	58160-819-12.....	137
42023-118-01.....	130	50419-344-05.....	122	58160-820-52.....	137
42023-123-06.....	125	50419-344-10.....	122	58160-821-11.....	137
42023-164-25.....	130	50419-344-15.....	122	58160-821-52.....	137
42023-168-01.....	123	50419-344-58.....	122	58160-823-11.....	137
42515-002-01.....	137	50419-346-10.....	122	58160-825-11.....	137
42702-102-15.....	131	50419-346-58.....	122	58160-825-52.....	137
42998-302-02.....	129	5041932511CS.....	120	58160-826-52.....	137
43528-002-05.....	137	5041932512CS.....	120	58160-842-11.....	136
43598-0210-40.....	136	5041932513CS.....	120	58160-842-52.....	136
43598-210-25.....	136	50458-197-15.....	134	58160-964-12.....	137
43598-210-55.....	136	50484-010-30.....	136	58160-976-02.....	136
43825-102-01.....	129	51079-105-20.....	133	58160-976-06.....	136
44206-300-01.....	129	51079-348-20.....	134	58468-090-01.....	130
44206-300-10.....	129	51079-966-20.....	133	58468-090-03.....	130
45802-021-46.....	135	51672-1263-1.....	136	58914-080-52.....	124
45802-046-35.....	135	51672-4026-6.....	133	59365-5030-0.....	131
45802-056-35.....	135	51672-4116-6.....	136	59584-138-01.....	125
45802-425-78.....	133	52544-092-76.....	130	59676-360-01.....	129
45802-980-64.....	135	52544-188-76.....	130	59746-115-06.....	134
45963-539-30.....	134	52544-892-28.....	133	59762-3051-2.....	132
46028-208-01.....	137	52584-004-34.....	126	59762-4537-1.....	128
49281-215-10.....	137	52584-023-10.....	125	59762-5008-1.....	133
49281-215-15.....	137	52584-035-05.....	135	60758-119-05.....	131
49281-225-10.....	137	52584-039-30.....	130	60758-188-05.....	131
49281-250-51.....	137	52584-044-25.....	124	60758-880-15.....	131
49281-286-10.....	136	52584-047-25.....	130	60793-130-10.....	129
49281-400-10.....	136	52584-080-01.....	125	60793-131-10.....	129
49281-400-20.....	136	52584-102-02.....	126	60793-215-05.....	130
49281-510-05.....	137	52584-178-02.....	127	60793-600-10.....	124
49281-545-03.....	136	52584-178-03.....	127	60793-602-10.....	124
49281-562-10.....	137	52584-182-03.....	128	60793-700-10F.....	124
49281-589-05.....	137	52584-276-01.....	127	60793-701-10.....	124
49281-752-21.....	137	52584-277-01.....	128	60793-702-10.....	124
49281-7522-2.....	137	52584-277-02.....	128	6131-4646-10.....	135
49281-790-20.....	138	52584-399-10.....	125	61314-204-15.....	131
49281-790-51.....	138	52584-418-13.....	134	61314-630-06.....	131
49281-860-10.....	137	52584-485-57.....	130	61314-631-36.....	131
49502-500-02.....	126	52584-489-27.....	130	61314-637-05.....	131

61314-643-05.....	132	66794-017-10.....	123	83044RNNXL.....	59
61314-645-11.....	135	66794-017-25.....	123	83044RNNXXL.....	59
61314-647-05.....	132	66993-022-27.....	123	83044RNNXXXL.....	59
61314-647-25.....	132	66993-023-27.....	123	8307-8133-14.....	135
61314-656-05.....	131	67457-118-50.....	124	8800 STYLE.....	61
61314-701-01.....	131	67457-134-50.....	128	88018QHW10.....	60
61703-309-06.....	130	67457-146-30.....	130	88018QHW12.....	60
61748-013-10.....	133	67457-273-10.....	128	88018QHW14.....	60
61748-018-01.....	134	67457-522-37.....	129	88018QHW16.....	60
61958-701-01.....	134	67877-124-40.....	136	88018QHW18.....	60
62175-118-32.....	134	67877-198-90.....	132	88018QHW20.....	60
62584-897-01.....	133	67877-319-05.....	133	88018QHW22.....	60
63304-098-19.....	134	67877-321-01.....	133	88018QHW24.....	60
63323-010-02.....	126	68001-237-00.....	132	88018QHW26.....	60
63323-012-01.....	129	68001-238-00.....	132	88018QHW28.....	60
63323-012-10.....	129	68001-246-17.....	134	88018QHW30.....	60
63323-013-02.....	130	68084-485-01.....	134	88018QHW4.....	60
63323-044-01.....	125	68094-204-62.....	134	88018QHW6.....	60
63323-0542-07.....	127	68180-122-01.....	132	88018QHW8.....	60
63323-064-03.....	128	68180-122-02.....	132	8832 STYLE.....	61
63323-064-11.....	128	68180-160-06.....	132	89122-724-20.....	127
63323-088-63.....	130	68180-160-13.....	132	A06000005.....	33
63323-0915-01.....	127	68308-152-15.....	136	ABB079830262H.....	150
63323-117-10.....	126	68382-041-01.....	134	ABB79830939H.....	150
63323-165-05.....	126	68382-050-01.....	133	ABO03P6001.....	19
63323-165-30.....	126	68382-187-01.....	133	ABO03P8025.....	20
63323-180-01.....	129	68382-261-01.....	132	ABO03P8125.....	20
63323-185-10.....	130	68382-423-06.....	129	ABO03P8225.....	20
63323-187-30.....	129	68462-180-22.....	136	ABO03P8325.....	20
63323-229-05.....	129	68462-190-01.....	133	ABO03P8425.....	20
63323-282-04.....	125	68462-396-10.....	134	ABO03PTESTKIT.....	19
63323-282-60.....	125	68462-564-17.....	136	ABO04J6020.....	20
63323-329-31.....	124	69401-000-02.....	138	ABO07P0201.....	20
63323-330-60.....	128	70003-300-20.....	134	ABO07P0204.....	20
63323-373-02.....	129	70003-300-65.....	134	ABO07P0205.....	20
63323-424-10.....	126	76204-600-30.....	123	ABO07P0207.....	20
63323-466-17.....	129	76329-3012-5.....	135	ABO07P0208.....	20
63323-467-57.....	129	76329-3013-05.....	135	ABO07P0210.....	20
63323-477-27.....	129	76329-3015-5.....	135	ABO07P0501.....	20
63323-484-57.....	130	76329-3369-1.....	129	ABO09P3125.....	20
63323-485-57.....	130	83044GRYL.....	59	ACNUC011.....	7
63323-486-27.....	130	83044GRYM.....	59	ACNUC071.....	7
63323-486-57.....	130	83044GRYS.....	59	AHS2131090H.....	68
63323-492-27.....	130	83044GRYXL.....	59	AJMPP6GREWH.....	203
63323-494-05.....	130	83044GRYXS.....	59	AJMPP9GREWH.....	203
63323-506-01.....	126	83044GRYXXL.....	59	ALR1115015.....	17
63323-516-10.....	126	83044GRYXXXL.....	59	ALR1115178.....	17
63323-545-01.....	127	83044QHWL.....	59	ALR1116554.....	17
63323-545-05.....	127	83044QHWM.....	59	ALR427000.....	12
63323-604-01.....	129	83044QHWS.....	59	ALR4581235020.....	9
63323-614-01.....	128	83044QHWXL.....	59	ALR4581245020.....	9
63323-616-09.....	124	83044QHWXS.....	59	ALR4581255020.....	9
63323-708-00.....	124	83044QHWXXL.....	59	ALR4581505020.....	14
63323-877-15.....	123	83044QHWXXS.....	59	ALR7D2648.....	15
64253-333-35.....	127	83044QHWXXXL.....	59	ALRI734000.....	12
64380-803-07.....	134	83044RCWL.....	59	ALRNAT024.....	12
64764-119-01.....	132	83044RCWM.....	59	AMV1364Z.....	8
64980-171-01.....	133	83044RCWS.....	59	AMV2161.....	8
65162-361-10.....	133	83044RCWXL.....	59	AMV2161Z.....	8
65162-702-94.....	135	83044RCWXS.....	59	AMV2166Z.....	8
65597-103-90.....	132	83044RCWXXL.....	59	AMV2184Z.....	8
65862-077-01.....	132	83044RCWXXXL.....	59	AMV5035.....	17
65862-175-01.....	134	83044RNNL.....	59	AMV5035Z.....	17
66220-719-30.....	133	83044RNNM.....	59	AMV5068AZ.....	17
66582-313-54.....	134	83044RNNNS.....	59	ARNA1250S.....	184

ARNA1250SG.....	184	B-D305902Z.....	142	BCC044540H.....	120
ARNA804.....	184	B-D305903Z.....	142	BCC111101BX.....	122
ARNA805.....	184	B-D305904ZZ.....	142	BCC111102BX.....	122
ARNA806.....	184	B-D305906Z.....	142	BCC111103BX.....	122
ARNA807.....	184	B-D305916Z.....	142	BCC111104BX.....	122
ARNA910.....	184	B-D305945Z.....	143	BCC131525.....	121
ARNA910ST.....	184	B-D305950.....	142	BCC131530.....	121
ARNA950.....	184	B-D305950Z.....	142	BCC131535.....	121
ARNA952.....	184	B-D309570Z.....	141	BCC131541.....	121
ARNCS103TDV.....	189	B-D309571Z.....	141	BCC131547.....	121
ARNCS103TLED.....	189	B-D309572Z.....	141	BCC131550.....	121
ARNES09.....	184	B-D309581Z.....	141	BCC131598.....	121
ARNES13.....	184	B-D309589Z.....	141	BCC131604.....	121
ARNES14.....	184	B-D309626Z.....	141	BCC131635.....	121
ARNES16.....	184	B-D309643Z.....	141	BCC131637.....	121
ARNES42.....	184	B-D309646Z.....	141	BCC131641.....	121
ARNES53.....	184	B-D309657Z.....	141	BCC131698.....	121
ARNXLDPSC.....	189	B-D363083Z.....	148	BCC131702.....	121
ARNXLDSC.....	189	B-D364815Z.....	143	BCC131705.....	121
ASC720001.....	90	B-D366430Z.....	148	BCC131709.....	121
ASC720001Z.....	90	B-D366450Z.....	148	BCC141030.....	120
ATT1G1V4H.....	152	B-D367281Z.....	145	BCC141111.....	121
ATT1Q1V4H.....	152	B-D367283Z.....	145	BCC141125.....	120
ATT2G1V4H.....	152	B-D367285Z.....	144	BCC141215.....	120
AUTODISPB.....	77	B-D367290.....	144	BCC141215H.....	120
AUTODISPW.....	77	B-D367292Z.....	144	BCC270131652.....	120
AZRES600.....	26	B-D367297Z.....	144	BCC44535.....	120
AZRES609.....	26	B-D367324Z.....	146	BCC450104.....	120
AZRES921.....	26	B-D367342Z.....	146	BCC450301.....	120
B-B-D305946Z.....	143	B-D367344Z.....	146	BCC450304.....	120
B-C60151AZ.....	10	B-D367587Z.....	148	BCC450305.....	120
B-C61130AZ.....	10	B-D367815Z.....	148	BCC450307.....	120
B-C64130AZ.....	10	B-D367844Z.....	148	BCC450701.....	120
B-C64151AZ.....	10	B-D367856Z.....	148	BCC516412BX.....	121
B-D302832Z.....	141	B-D367861Z.....	148	BCC516413BX.....	121
B-D302995Z.....	141	B-D367871Z.....	148	BCC516414BX.....	121
B-D305106Z.....	142	B-D367925Z.....	148	BCC516415BX.....	121
B-D305122Z.....	141	B-D367960Z.....	148	BCC526416.....	121
B-D305124Z.....	142	B-D367961Z.....	148	BCC526417.....	121
B-D305125.....	142	B-D367962Z.....	148	BCC900001.....	120
B-D305125Z.....	142	B-D367986Z.....	148	BCC900001H.....	120
B-D305127Z.....	141	B-D367988Z.....	148	BCC900002.....	120
B-D305128Z.....	141	B-D368607Z.....	144	BCC900004.....	120
B-D305136Z.....	141	B-D368608CSZ.....	144	BCC900005.....	120
B-D305155Z.....	141	B-D368650CS.....	144	BCC900006.....	120
B-D305156Z.....	141	B-D368651CS.....	144	BCC901407.....	120
B-D305195Z.....	141	B-D369714Z.....	148	BCC901702.....	120
B-D305196Z.....	141	B-D381412Z.....	147	BCC901901.....	120
B-D305269Z.....	142	B-D381423Z.....	147	BCC902001.....	120
B-D305270Z.....	142	B-D381444Z.....	147	BCC902101.....	120
B-D305271Z.....	142	BAGGK4500.....	205	BCC902102HH.....	120
B-D305500Z.....	141	BAGGW2500.....	205	BCC902801.....	120
B-D305761Z.....	142	BAGGW4500.....	205	BCC902801H.....	120
B-D305762Z.....	142	BAGGW8500.....	205	BDD256038.....	12
B-D305763Z.....	142	BAGGX2500.....	205	BDD256040.....	12
B-D305766Z.....	142	BAGSK1466SOST.....	205	BDD256045.....	12
B-D305767Z.....	142	BAH1554Z.....	8	BDD256055.....	12
B-D305778Z.....	142	BAH1759.....	8	BEZBECKHB.....	204
B-D305780Z.....	142	BAH1760.....	8	BEZBECKHW.....	204
B-D305781Z.....	142	BAH1780.....	8	BEZBECKMB.....	204
B-D305782Z.....	142	BAH2083Z.....	8	BEZBECKMW.....	204
B-D305783Z.....	142	BAH5075US.....	17	BEZBEDWH12PLA.....	201
B-D305789Z.....	142	BCC014957.....	120	BEZBEDWH16PLA.....	201
B-D305900Z.....	142	BCC0270131630.....	121	BEZBEDWH20PLA.....	201
B-D305901Z.....	142	BCC044540.....	120	BEZBEDWHC8PLA.....	201

BEZBEFHB	204	BHL2F7113CS	150	CHLS88773	16
BEZBEFHBINV	204	BHL2F7113HH	150	CNS1101KL	179
BEZBEFHW	204	BHL2F7122HH	150	CNS201HR	182
BEZBEFHWINV	204	BHL2F7123HH	150	CNS2101KL	179
BEZBEFL15PLA	203	BHL2N3379Z	141	CNSADPT31	180
BEZBEFL90CPLA	203	BHM03000	38	CPRM1116	99
BEZBEFMB	204	BHM03000H	38	CPRM1126	99
BEZBEFMBINV	204	BHM03100	38	CPRM2216	99
BEZBEFMW	204	BHM03300	38	CPRM2226	99
BEZBEFMWINV	204	BID0326076316	8	CPRM3312	99
BEZBEFRB12	203	BID1008552	8	CPRM3322	99
BEZBEFRB16	203	BID1621173001	19	CRI1000	53
BEZBEFRP10	203	BID1895362160	8	CRI1001	53
BEZBEFRP103	203	BID1895427160	8	CRI1002	53
BEZBEFRP6	203	BID260763602	8	CRI1003	53
BEZBEFRP7	203	BID3617556001	8	CRI1004	53
BEZBEFRP9	203	BID4348150001	19	CRI1202	53
BEZBEFRP93	203	BID4625315160	19	CRI1203	53
BEZBEFSP10	203	BID4837975001	19	CRI2000	54
BEZBEFSP6	203	BID7341890190	12	CRI2001	54
BEZBEFSP8	203	BID7341911190	12	CRI2002	54
BEZBEHC10PLA	201	BID7341920190	12	CRI2003	54
BEZBEHC12PLA	201	BID8160104190	12	CRI2006	54
BEZBEHC16PLA	201	BLDGLV3	89	CRI2007	54
BEZBEHC20PLA	201	BLK7757SSHR	186	CRI2010	54
BEZBEHC4PLA	201	BLP39182	131	CRI4000	53
BEZBEHC8PLA	201	BMG573114	146	CRI4000B	53
BEZBEHL1024PL	201	BMGL8001H	150	CRI4000Z	53
BEZBEHL8PLA	200-202	BMGL8500H	150	CRI4001	53
BEZBEKHB	204	BMGL850101	150	CRI4001G	53
BEZBEKHBINV	204	BMGL8502H	150	CRI4003	53
BEZBEKHW	204	BMGR520001H	150	CRI4200	53
BEZBEKHWINV	204	BMGR520101H	150	CRI4201	53
BEZBEKMB	204	BOVCRY16	184	CRI5000	53
BEZBEKMBINV	204	BOVCRY8	184	CRI5000Z	53
BEZBEKMW	204	BOVSE02	184	CRI5001	53
BEZBEKMWINV	204	BRW3000	194	CRU911447	40
BEZBEL1024PLA	202	BRW5000	191	CT24BLK	113, 164
BEZBESC12PLA	203	BRW5001	191	CT24BLU	113, 164
BEZBESC16PLA	203	BRW6501	191	CT24PNK	113, 164
BEZBESC24PLA	203	BRW6800	191	CT24PRP	113, 164
BEZBESC32PLA	203	BRW6801	191	CT24WHT	113, 164
BEZBESC8PLA	203	CCD144X99COR	196	CT34BLK	113, 164
BEZBESHB	204	CCD144X99ENC	196	CT34BLU	113, 164
BEZBESHBINV	204	CCD144X99PP	196	CT34PNK	113, 164
BEZBESHW	204	CCD144X99TAN	196	CT34PRP	113, 164
BEZBESHWINV	204	CCD240X99COR	196	CT34WHT	113, 164
BEZBESMB	204	CCD240X99ENC	196	CT44BLK	113, 164
BEZBESMBINV	204	CCD240X99PP	196	CT44BLU	113, 164
BEZBESMW	204	CCD240X99TAN	196	CT44PNK	113, 164
BEZBESMWINV	204	CCYEMS052	182	CT44PRP	113, 164
BEZBESSB	204	CHLS10228	16	CT44WHT	113, 164
BEZBESSBINV	204	CHLS10311	16	CTRO01201	96
BEZBESSW	204	CHLS10959	16	CTRO01305H	96
BEZBESSWINV	204	CHLS10986	16	CUR00001RB	162
BHL2B0062QH	150	CHLS10987	16	CUR00002RB	162
BHL2B0088	150	CHLS10988	16	CUR001109Z	154
BHL2B0089	150	CHLS10989	16	CUR001209Z	154
BHL2B1308Z	150	CHLS10990	16	CUR001231	154
BHL2B1309Z	150	CHLS10991	16	CUR003501	154
BHL2B1322QH	150	CHLS11010	16	CUR003545Z	154
BHL2B1324XH	150	CHLS11846	16	CUR005331	154
BHL2B2322Q	150	CHLS13454	16	CUR005345Z	154
BHL2B7114X	150	CHLS26333	16	CUR01496	154
BHL2B7114XH	150	CHLS52193	16	CUR01496H	154

CUR110556.....	154	DTM12J12.....	200	DYND12582H.....	103
CUR110556H.....	135	DTM12JL.....	200, 202	DYND12584.....	103
CUR250317Z.....	159	DTM14J16.....	200	DYND12588.....	103
CUR250330Z.....	159	DTM16FTLS.....	200, 202	DYND13508.....	103
CUR250381Z.....	159	DTM16RCL.....	200, 202	DYND13510.....	103
CUR250383Z.....	159	DTM16SL.....	200, 202	DYND13512.....	103
CUR251136Z.....	159	DTM200PC.....	201	DYND13514.....	103
CUR251180Z.....	159	DTM32SL.....	200, 202	DYND13515.....	103
CUR251272Z.....	159	DTM3TI.....	205	DYND13516.....	103
CUR251318Z.....	159	DTM3VT19S.....	205	DYND13518.....	103
CUR251336Z.....	159	DTM400PC.....	201	DYND13520.....	103
CUR251390Z.....	159	DTM5T1GR0205.....	205	DYND13522.....	103
CUR251636Z.....	159	DTM5T1N0195.....	205	DYND13614.....	103
CUR253180WZ.....	159	DTM5T3N0196.....	205	DYND13714.....	103
CUR253180Z.....	159	DTM5VT195N11.....	205	DYND15205.....	104
CUR253200H.....	159	DTM8FTL.....	200, 202	DYND16900.....	104
CUR253220Z.....	159	DTM8J8.....	200	DYND30100.....	25
CUR253440Z.....	159	DTM8JL.....	200, 202	DYND30102.....	25
CUR253590WZ.....	159	DTM8SL.....	200, 202	DYND30103.....	25
CUR253590Z.....	159	DTMPL1N.....	201	DYND30260.....	24
CUR45701RB.....	162	DTMPL2N.....	201	DYND30261.....	24
CUR45702RB.....	162	DTMPL4N.....	201	DYND30271.....	24
CUR5305.....	154	DYNC1810.....	103	DYND30281.....	24
CUR8103.....	71	DYNC1816.....	104	DYND30330.....	25, 105
CUR8104.....	71	DYNC1820.....	104	DYND30331.....	25
CUR8105.....	71	DYND10300.....	104	DYND30335.....	25
CUR8106.....	71	DYND10305.....	104	DYND30350.....	25, 105
CUR8107.....	71	DYND10350.....	104	DYND30351.....	25
CUR9224.....	71	DYND10402.....	104	DYND30382.....	25
CUR9225.....	71	DYND10407.....	104	DYND30385.....	25
CUR9226.....	71	DYND10440.....	103	DYND30387.....	25
CUR9227.....	71	DYND10440H.....	103	DYND30389.....	25
CURSALINE7.....	165	DYND10500.....	104	DYND30510.....	25
DAY020536H.....	95	DYND10700.....	102	DYND34255.....	24
DAY020543Z.....	95	DYND10701.....	102	DYND34260.....	24
DEF5070A.....	40	DYND10702.....	102	DYND34265.....	24
DEFM5066HOME.....	40	DYND10703.....	102	DYND34270.....	24
DEFM5066PACK.....	40	DYND10712.....	102	DYND34275.....	24
DEFM5071A.....	40	DYND10713.....	102	DYND36500H.....	106
DEFM5072A.....	40	DYND10714.....	102	DYND36600H.....	105
DFTDCF1310EN.....	41	DYND10721.....	102	DYND36605H.....	105
DFTDCFA100EN.....	41	DYND10722.....	102	DYND371603Z.....	95
DFTDCFA120.....	41	DYND10723.....	102	DYND40540.....	150
DIE60052101.....	206	DYND10724.....	102	DYND40570.....	150
DLS137038.....	25	DYND10725.....	102	DYND4712H.....	105
DLS137135.....	25	DYND10800.....	105	DYND500MLS.....	150
DRC333661186.....	67	DYND10805.....	105	DYND500MLW.....	150
DRC333661247Z.....	67	DYND10806.....	105	DYND69STAT.....	24
DRC333661261.....	67	DYND10810.....	105	DYND70277.....	112
DRC4133301201.....	67	DYND10815.....	105	DYND70279.....	112
DRCPC1300Z.....	67	DYND10820.....	105	DYND70280.....	112
DRCPC1400Z.....	67	DYND11214.....	102	DYND70292H.....	112
DRCPC1500Z.....	67	DYND11216.....	102	DYND70292Z.....	112
DRCPC1604Z.....	67	DYND11218.....	102	DYND70293H.....	112
DRCPC2400Z.....	67	DYND11220.....	102	DYND70293Z.....	112
DRUG10DIPDXC.....	15	DYND11752.....	102	DYND70400.....	98
DRUG1127023.....	15	DYND11754.....	102	DYND70401L.....	98
DRUG1147022.....	15	DYND11756.....	102	DYND70401LZ.....	98
DRUG8CUPDX.....	15	DYND11758.....	102	DYND70401M.....	98
DTM100PC.....	201	DYND11760.....	102	DYND70401MZ.....	98
DTM10J10.....	200	DYND11762.....	102	DYND70401S.....	98
DTM10SL.....	200, 202	DYND12200.....	104	DYND70401SZ.....	98
DTM10T1N0198.....	205	DYND12550.....	104	DYND70420L.....	98
DTM10T3N0199.....	205	DYND12574.....	103	DYND70420LZ.....	98
DTM10V19S.....	205	DYND12578.....	103	DYND70420M.....	98

GPC20887	48	HCS54514S	99	JPS16332	143
GPC33587	48	HCSU4511B	99	KAWD3K223GZ	146
GPC48580HH	49	HCSU4514S	99	KAWDBMA21G	145
GPC56797A	49	HCSU4517S	99	KAWDBMA21GZ	145
GPC59447A	47	HCSU4518S	99	KAWDBMA23G	145
GPC59488A	48	HH62G02	74	KAWDBMA23GZ	145
GPC89460	47	HH62G1000	74	KAWDBMA25G	145
GPC89480	48	HH62G15	74	KAWDBMA25GZ	145
GPCGP27700	48	HH70F1000	74	KAWDBMH21G	145
GPK80600	203	HH70F18	74	KAWDBMH21GZ	145
GPK80900	203	HH70F75	74	KAWDBMH23G	145
GPK82100	203	HH70G02	74	KAWDBMH23GZ	145
GRI454052	148	HH70G04	74	KAWDBMH25G	145
GRI454081	148	HH70G08	74	KAWDBMH25GZ	145
GRI454222	148	HH70G1000	74	KDL500500	39
GRI454236	148	HH70G15	74	KDL500500Z	39
GRI454237	148	HH70GPACKET	74	KUT64024MEDA	76
GRI454238	148	HH70W160	75	KUT64041MEDA	76
GRI454302	148	HH70W50	75	KUT64078MEDA	76
GRI454320	148	HH85G02	74	KUT64341MED	77
GRI454322	148	HH85G1000	74	KUT65036	76
GRI454428CS	148	HH85G15	74	KUT68241MEDA	75
HCS4483H	100	HHABSP04	76	KUT69041MED	77
HCS4485H	100	HHABSP1000	76	LBI095305	121
HCS4486H	100	HHABSP1000FM	76	LBI095350	120
HCS4511B	99	HHABSP16	76	LBI132306	121
HCS4514	99	HHABSP75	76	LBI132316	122
HCS4516B	99	HHABSPG	76	LBI132321	122
HCS4517	99	HHSP04	76	LBI132327	122
HCS4518	99	HHSP1000	76	LBI132352	122
HCS4561	99	HHSP1000FM	77	LBI132361	122
HCS4565	99	HHSP16	76	LBI132387	122
HCS4566	99	HHSP75	76	LBI132390	122
HCS4568	99	HHSPG	76	LBI132395	122
HCS4568F	99	HLS5072	10	LBI132406	121
HCS4621B	99	HLS5073	10	LBI132411	122
HCS4630	99	HMEFCUP4144A2	15	LBI132487	122
HCS4630B	99	HMESCCUP13124	15	LBI133200	122
HCS4631B	99	HMESCDIP12104	15	LBI133206	122
HCS64184	99	HMET1CK30	10	LBI133216	122
HCS70004	100	HMET1CK30T	10	LBI133221	122
HCS70004F	100	HMET1CK50	10	LBI133306	121
HCS70004H	100	HMET1CM50	10	LBI133311	121
HCS70004RD	100	HMET1CT30	10	LBI133316	122
HCS70004RDH	100	HMET1TB50	10	LBI133321	122
HCS93250	44	HMET1TC01	10	LBI133361	122
HCS93250H	44	HMQ110706	18	LBI133387	122
HCS93500	44	HMQ110723	18	LBI133390	122
HCS93500H	44	HMQ111715	17	LBI133395	122
HCSA100	100	HMQ111716	17	LBI318305	121
HCSA113	100	HMQ120706	18	LBI318350	121
HCSA114	100	HMQ121721	17	LBI481604	121
HCSFISHNEB	100	HUD1801	44	M-B703443	16
HCSFISHNEBH	100	HUD1801H	44	M-B703447	16
HCSFSATSPORT	101	HUH71846	205	M-B706760	16
HCSM70	101	IIS2282906	135	M-B706762	16
HCSM70P	101	IMS3013SP	136	M5240	191
HCSMD400	101	IMS3015SP	135	MANDISPB	77
HCSMD400A	101	IMS6300	136	MANDISPW	77
HCSMD400P	101	ITD827976	206	MAP101400H	62
HCSPF100	44, 100	J-J20390H	62	MAP102800	62
HCSPF100H	44	J-J20392H	62	MAP102800H	62
HCSPF200	100	J-J2252HH	62	MAP102805	62
HCSPUPPY	100	J-J2920CS	62	MAP102805H	62
HCSPUPPYH	100	JOJ200411	206	MAP10303H	62

MAP10304H.....	62	MDR7211142.....	188	MDS0157107048.....	42
MAP106000.....	62	MDR7211216FL.....	188	MDS0228820.....	80
MAP106000H.....	62	MDR7211602.....	189	MDS032273H.....	95
MAP131000H.....	62	MDR721606.....	188	MDS032280Z.....	95
MAP131024CS.....	62	MDR721607.....	188	MDS032285H.....	95
MAP131100Z.....	62	MDR73102030.....	194	MDS032285Z.....	95
MAP131150.....	62	MDR731024.....	194	MDS032290H.....	95
MAP131155Z.....	62	MDR731027.....	194	MDS032290Z.....	95
MAP133300H.....	62	MDR731030.....	194	MDS0391650.....	82
MBFTC990019.....	208	MDR732024.....	194	MDS0391656.....	82
MBFTC990020.....	208	MDR732027.....	194	MDS0391657.....	82
MBFTC99004.....	208	MDR7330102.....	194	MDS0440317.....	88
MBH39802000.....	97	MDR738565.....	194	MDS046002ZZ.....	160
MBH39803000.....	97	MDR73857.....	194	MDS046003ZZ.....	160
MCHEM100.....	26	MDR741202.....	192	MDS046004ZZ.....	160
MCHEM101.....	26	MDR741700.....	192	MDS046006ZZ.....	160
MCHEM110.....	26	MDR741701.....	192	MDS055002H.....	160
MCHEM111.....	26	MDR748870.....	192	MDS055002Z.....	160
MCHEM130.....	26	MDR75454.....	180	MDS055003H.....	160
MCHEM131.....	26	MDR75457.....	180	MDS055003Z.....	160
MCHEM140.....	26	MDR75470.....	180	MDS055004H.....	160
MCHEM145.....	26	MDR75473.....	180	MDS055004Z.....	160
MCHEM150.....	26	MDR776250.....	186	MDS055006H.....	160
MCHEMPF120.....	24	MDREX21.....	194	MDS055006Z.....	160
MCHEMPF180.....	24	MDREX23.....	194	MDS0610116.....	89
MCHEMPF20.....	24	MDREX27.....	194	MDS0610300.....	89
MCHEMPF250.....	24	MDRNEWP15.....	195	MDS0610301.....	80, 89
MCHEMPF40.....	24	MDRNEWP1W.....	195	MDS0610400.....	89
MCHEMPF500.....	24	MDRNEWPB33.....	195	MDS0610700.....	89
MCHEMPF60.....	24	MDRNEWPL.....	195	MDS066002ZZ.....	113
MCHEMPF90.....	24	MDRQUS210C.....	174	MDS066003Z.....	113
MCHEMPF940.....	24	MDRQUS220C.....	174	MDS066004Z.....	113
MDH13204281.....	205	MDRQUS224C.....	174	MDS066006Z.....	113
MDR0224132.....	188	MDRQUS230C.....	174	MDS067003Z.....	113
MDR102613.....	179	MDRQUS234C.....	174	MDS067004Z.....	113
MDR107423.....	181	MDRQUS235C.....	174	MDS067006Z.....	113
MDR110611.....	179	MDRQUS239C.....	174	MDS0725964.....	80
MDR112570.....	181	MDRQUS240C.....	174	MDS0728981.....	80
MDR121304.....	180	MDRQUS245C.....	174	MDS0728986.....	80
MDR133077.....	179	MDRQUS250C.....	174	MDS0740503.....	82
MDR141447.....	179	MDRQUS260C.....	174	MDS0740504.....	82
MDR158463.....	181	MDRQUS265C.....	174	MDS0740581.....	82
MDR158464.....	181	MDRQUS270C.....	174	MDS0740852.....	82
MDR158465.....	181	MDRSE50FL.....	189	MDS0740975.....	82
MDR158466.....	181	MDRSE50FL.....	188	MDS0744091.....	82
MDR301327.....	194	MDRSE50LSC.....	188	MDS0746773.....	80
MDR3020270.....	194	MDRSE50SC.....	189	MDS077002H.....	160
MDR440FD.....	180	MDRT77.....	194	MDS077002Z.....	160
MDR505.....	185	MDS0033310.....	81	MDS077003H.....	160
MDR505B.....	185	MDS0033315.....	81	MDS077003Z.....	160
MDR505LC.....	185	MDS0033320.....	81	MDS077004H.....	160
MDR66000.....	195	MDS0033325.....	81	MDS077004Z.....	160
MDR66010.....	195	MDS0033330.....	81	MDS077006H.....	160
MDR66020.....	195	MDS0033335.....	81	MDS077006Z.....	160
MDR712102.....	185	MDS0033340.....	81	MDS0800411.....	80
MDR712130.....	185	MDS0033350.....	81	MDS0816014.....	80
MDR712130B.....	185	MDS0033360.....	81	MDS0816017.....	80
MDR712135B.....	185	MDS0033380.....	81	MDS0816114.....	80
MDR712135SH.....	185	MDS0033520.....	81	MDS0816117.....	80
MDR712155BS.....	185	MDS0033530.....	81	MDS0820909.....	80
MDR712155S.....	185	MDS0033540.....	81	MDS0828115.....	80
MDR717150.....	186	MDS0033550.....	81	MDS0834011.....	80
MDR721010.....	188	MDS0033570.....	81	MDS0838510.....	80
MDR721109.....	189	MDS0157040172.....	40, 42	MDS0839011.....	80
MDR7211122.....	188	MDS0157040173.....	40, 42	MDS086001B.....	161

MDS086002	161	MDS094188H	94	MDS15112	89
MDS086002H	161	MDS098001H	94	MDS15115	89
MDS086003	161	MDS098001Z	94	MDS15115L	89
MDS086003H	161	MDS098002	94	MDS15120	89
MDS086004	161	MDS098003Z	94	MDS15121	89
MDS086004H	161	MDS098003ZH	94	MDS15122	89
MDS087002LFH	160	MDS098011H	94	MDS15210	89
MDS087002LFZ	160	MDS098014	94	MDS15211	89
MDS087003LFH	160	MDS098014H	94	MDS15212	89
MDS087003LFZ	160	MDS098015	94	MDS15215	89
MDS087004LFH	160	MDS098015H	94	MDS15220	89
MDS087004LFZ	160	MDS098017	94	MDS15221	89
MDS087006LFH	160	MDS098017H	94	MDS15310	89
MDS087006LFZ	160	MDS1001	28	MDS15311	89
MDS087104LFZ	160	MDS1010214	80	MDS15315	89
MDS087106LFH	160	MDS1012014	80	MDS15410	89
MDS087106LFZ	160	MDS1017012	80	MDS15411	89
MDS087156LFH	160	MDS1018012	80	MDS15415	89
MDS0880111	80	MDS1018612	80	MDS160604	115
MDS0880209	80	MDS10749	83	MDS160624	115
MDS0880414	80	MDS10755	83	MDS160644	115
MDS089001	161	MDS1080209	80	MDS160664	115
MDS089001H	161	MDS1080211	80	MDS160684	115
MDS089002	161	MDS10803533	35	MDS160694	115
MDS089002H	161	MDS10803534	35	MDS160824	115
MDS089003	161	MDS108055LT	72	MDS160844	115
MDS089003CP	161	MDS108060LT	72	MDS160864	115
MDS089003H	161	MDS108065LT	72	MDS160884	115
MDS089004	161	MDS108070LT	72	MDS160894	115
MDS089004H	161	MDS108075LT	72	MDS1617025	84
MDS089006	161	MDS108080LT	72	MDS1623126	84
MDS089006CP	161	MDS108085LT	72	MDS1702ATH	114
MDS0890114	80	MDS108090LT	72	MDS1702ATSH	114
MDS0890118	80	MDS1211856	80	MDS1702BTSH	114
MDS089202	161	MDS1211951	80	MDS1702BTSH	114
MDS089204	161	MDS1222012	80	MDS1702CTH	114
MDS089207	161	MDS1222112	80	MDS1702CTSH	114
MDS0895014	80	MDS1224014	80	MDS1702DTH	114
MDS0895015	80	MDS1224114	80	MDS1702DTSH	114
MDS0895016	80	MDS1225014	80	MDS1702ETH	114
MDS0895018	80	MDS1225114	80	MDS1702ETSH	114
MDS090670Z	94	MDS129CS	190	MDS1702FTH	114
MDS090735Z	94	MDS131040	32	MDS1702FTSH	114
MDS090780	94	MDS131040H	32	MDS1702GTH	114
MDS090780H	94	MDS137000	98	MDS1702GTSH	114
MDS090855	94	MDS137000H	98	MDS1703ABH	114
MDS090855H	94	MDS138020	97	MDS1703ABSH	114
MDS092005H	96	MDS138020H	97	MDS1703BBH	114
MDS092010H	96	MDS139008	97	MDS1703BBSH	114
MDS093810Z	94	MDS139009	97	MDS1703CBH	114
MDS093901ZZ	95	MDS139010	97	MDS1703CBSH	114
MDS093902ZZ	95	MDS148065	98	MDS1703DBH	114
MDS093904H	95	MDS15010	89	MDS1703DBSH	114
MDS093906H	95	MDS15011	89	MDS1703EBH	114
MDS093917H	95	MDS15012	89	MDS1703EBSH	114
MDS093940HH	95	MDS15015	89	MDS1703FBH	114
MDS093941H	95	MDS15020	89	MDS1703FBSH	114
MDS093943H	95	MDS15020L	89	MDS1703GBH	114
MDS093944H	95	MDS15021	89	MDS1703GBSH	114
MDS093945H	95	MDS15022	89	MDS1712ATH	114
MDS093946H	95	MDS15023	89	MDS1712ATSH	114
MDS093947H	95	MDS15110	89	MDS1712BTH	114
MDS094181H	94	MDS15110L	89	MDS1712BTSH	114
MDS094184H	94	MDS15111	89	MDS1712CTH	114
MDS094186Z	94	MDS15111L	89	MDS1712CTSH	114

MDS1713ABH.....	114	MDS21463H.....	96	MDS51514-10.....	176
MDS1713ABSH.....	114	MDS2200202E.....	35	MDS52514-10.....	176
MDS1713BBH.....	114	MDS2200204E.....	35	MDS6300100.....	86
MDS1713BBSH.....	114	MDS2410013.....	80	MDS6300101.....	86
MDS1713CBH.....	114	MDS2410413.....	80	MDS6300102.....	86
MDS1713CBSH.....	114	MDS2413013.....	80	MDS6300103.....	86
MDS1821302.....	89	MDS2414413.....	80	MDS6300104.....	86
MDS1821413.....	88	MDS2418014.....	80	MDS6300105.....	86
MDS1821601.....	88	MDS2418414.....	80	MDS6300106.....	86
MDS1821602.....	88	MDS2441414.....	80	MDS6300107.....	86
MDS190505.....	105	MDS2583.....	70	MDS6300200.....	86
MDS190510.....	105	MDS2584.....	70	MDS6300201.....	86
MDS192074.....	71	MDS2585.....	70	MDS6300202.....	86
MDS192075.....	71	MDS2586.....	70	MDS6300203.....	86
MDS192076.....	71	MDS2587.....	70	MDS6300204.....	86
MDS192077.....	71	MDS2591211.....	82	MDS6300205.....	86
MDS192094B.....	70	MDS3001.....	28	MDS6300206.....	86
MDS193094.....	70	MDS3001LA.....	28	MDS6300207.....	86
MDS195174.....	70	MDS3001PLUS.....	28	MDS6300300.....	84
MDS195175.....	70	MDS3001U.....	28	MDS6300301.....	84
MDS195176.....	70	MDS3002.....	28	MDS6300302.....	84
MDS195177.....	70	MDS3002LA.....	28	MDS6300303.....	84
MDS195283.....	70	MDS36533001.....	36	MDS6300304.....	84
MDS195284.....	70	MDS36533002.....	36	MDS6300305.....	84
MDS195285.....	70	MDS36533003.....	36	MDS6300306.....	84
MDS195286.....	70	MDS36533004.....	36	MDS6300307.....	84
MDS195287.....	70	MDS36533005.....	36	MDS6300308.....	84
MDS200100.....	91	MDS36533006.....	36	MDS6300309.....	84
MDS200360.....	62	MDS36533007.....	36	MDS6300310.....	84
MDS200565.....	91	MDS36533008.....	36	MDS6300400.....	84
MDS200590.....	91	MDS36533009.....	36	MDS6300401.....	84
MDS200595.....	91	MDS36533010.....	36	MDS6300402.....	84
MDS200615.....	91	MDS36533011.....	36	MDS6300403.....	84
MDS200836H.....	92	MDS36533012.....	36	MDS6300404.....	84
MDS200846H.....	92	MDS36533013.....	36	MDS6300405.....	84
MDS200848H.....	92	MDS36533014.....	36	MDS6300406.....	84
MDS2010616.....	84	MDS36533015.....	36	MDS6300407.....	84
MDS201222.....	81	MDS36533016.....	36	MDS6300408.....	84
MDS202000Z.....	96	MDS36533017.....	36	MDS6300409.....	84
MDS202010H.....	96	MDS36533018.....	36	MDS6300410.....	84
MDS202010Z.....	96	MDS36533018B.....	36	MDS6300500.....	84
MDS202050Z.....	96	MDS36533019.....	36	MDS6300501.....	84
MDS202055H.....	96	MDS36533019B.....	36	MDS6300502.....	84
MDS202055Z.....	96	MDS3701.....	34	MDS6300503.....	84
MDS202060Z.....	97	MDS3703.....	34	MDS6300504.....	84
MDS202065H.....	97	MDS3704.....	34	MDS6300505.....	84
MDS202070HH.....	97	MDS3716.....	34	MDS6300506.....	84
MDS202073Z.....	97	MDS3725.....	35	MDS6300507.....	84
MDS202075Z.....	97	MDS3726.....	35	MDS6300600.....	85
MDS202077Z.....	97	MDS4001.....	28	MDS6300601.....	85
MDS202080Z.....	97	MDS4001LAT.....	28	MDS6300602.....	85
MDS202085Z.....	96	MDS4003.....	29	MDS6300603.....	85
MDS202090Z.....	96	MDS4406001.....	88	MDS6300604.....	85
MDS202095Z.....	96	MDS4406002.....	88	MDS6300605.....	85
MDS202105Z.....	96	MDS4407202.....	88	MDS6300606.....	85
MDS202110Z.....	97	MDS4410203.....	88	MDS6300607.....	85
MDS202999.....	81	MDS4410205.....	88	MDS6300700.....	85
MDS203444.....	81	MDS4418001.....	88	MDS6300701.....	85
MDS204111.....	81	MDS4418002.....	88	MDS6300702.....	85
MDS204888.....	81	MDS4420614.....	88	MDS6300703.....	85
MDS205777.....	81	MDS4426108.....	88	MDS6300704.....	85
MDS206333.....	81	MDS4615720.....	88	MDS6300705.....	85
MDS209555.....	81	MDS4641518.....	80	MDS6300706.....	85
MDS21460Z.....	96	MDS5001.....	28	MDS6300707.....	85
MDS21462H.....	96	MDS50514-10.....	176	MDS6300800.....	85

MDS6300801.....	85	MDS6301505.....	85	MDS6802.....	70
MDS6300802.....	85	MDS6301506.....	85	MDS6803.....	70
MDS6300803.....	85	MDS6301507.....	85	MDS6804.....	70
MDS6300804.....	85	MDS6301508.....	85	MDS7011001.....	84
MDS6300805.....	85	MDS6301509.....	85	MDS7011002.....	84
MDS6300806.....	85	MDS6301600.....	86	MDS7011003.....	84
MDS6300807.....	85	MDS6301601.....	86	MDS7012001.....	84
MDS6300808.....	85	MDS6301602.....	86	MDS7012002.....	84
MDS6300900.....	85	MDS6301603.....	86	MDS7012003.....	84
MDS6300901.....	85	MDS6301604.....	86	MDS7036024.....	84
MDS6300902.....	85	MDS6301605.....	86	MDS7050425.....	84
MDS6300903.....	85	MDS6301606.....	86	MDS705110H.....	152
MDS6300904.....	85	MDS6301607.....	86	MDS705115H.....	152
MDS6300905.....	85	MDS6301608.....	86	MDS705153H.....	152
MDS6300906.....	85	MDS6301609.....	86	MDS705154H.....	152
MDS6300907.....	85	MDS6301700.....	86	MDS705201H.....	152
MDS6300908.....	85	MDS6301701.....	86	MDS705202H.....	152
MDS6300909.....	85	MDS6301702.....	86	MDS705203H.....	152
MDS6301000.....	85	MDS6301703.....	86	MDS705208.....	152
MDS6301001.....	85	MDS6301704.....	86	MDS705208SLH.....	152
MDS6301002.....	85	MDS6301705.....	86	MDS705308.....	152
MDS6301003.....	85	MDS6301706.....	86	MDS7058132.....	84
MDS6301004.....	85	MDS6301707.....	86	MDS7059440.....	84
MDS6301005.....	85	MDS6301800.....	86	MDS7063101.....	84
MDS6301006.....	85	MDS6301801.....	86	MDS7067201.....	84
MDS6301007.....	85	MDS6301802.....	86	MDS7067211.....	84
MDS6301008.....	85	MDS6301803.....	86	MDS707011.....	152
MDS6301100.....	85	MDS6301804.....	86	MDS707930H.....	152
MDS6301101.....	85	MDS6301805.....	86	MDS707953H.....	152
MDS6301102.....	85	MDS6301806.....	86	MDS708555H.....	90
MDS6301103.....	85	MDS6301807.....	86	MDS7405812.....	82
MDS6301104.....	85	MDS6301900.....	86	MDS7405813.....	82
MDS6301105.....	85	MDS6301901.....	86	MDS7406013.....	81
MDS6301106.....	85	MDS6301902.....	86	MDS7411110.....	81
MDS6301107.....	85	MDS6301903.....	86	MDS74113420.....	82
MDS6301108.....	85	MDS6301904.....	86	MDS7411661.....	81
MDS6301109.....	85	MDS6301905.....	86	MDS7412803.....	88
MDS6301200.....	85	MDS6301906.....	86	MDS7412804.....	88
MDS6301201.....	85	MDS6301907.....	86	MDS7432711.....	81
MDS6301202.....	85	MDS6301908.....	86	MDS7433614.....	81
MDS6301203.....	85	MDS6301909.....	86	MDS7434113.....	81
MDS6301204.....	85	MDS6302000.....	85	MDS7436013.....	81
MDS6301205.....	85	MDS6302001.....	85	MDS7436213.....	81
MDS6301206.....	85	MDS6302002.....	85	MDS7436511.....	81
MDS6301207.....	85	MDS6302003.....	85	MDS7440572.....	82
MDS6301208.....	85	MDS6302004.....	85	MDS7440573.....	82
MDS6301209.....	85	MDS6302005.....	85	MDS7440574.....	82
MDS6301302.....	86	MDS6302006.....	85	MDS7444089.....	82
MDS6301303.....	86	MDS6302007.....	85	MDS7444090.....	82
MDS6301304.....	86	MDS6302008.....	85	MDS80337Z.....	177
MDS6301305.....	86	MDS6303000.....	86	MDS80430I.....	186
MDS6301306.....	86	MDS6303001.....	86	MDS80494.....	190
MDS6301307.....	86	MDS6303002.....	86	MDS80500.....	190
MDS6301308.....	86	MDS6303003.....	86	MDS80525.....	186-187
MDS6301402.....	86	MDS6303010.....	86	MDS80529.....	186
MDS6301403.....	86	MDS6303012.....	86	MDS80529P.....	186
MDS6301404.....	86	MDS6303015.....	86	MDS80530KD.....	187
MDS6301405.....	86	MDS6305001.....	86	MDS80540.....	177
MDS6301406.....	86	MDS6305002.....	86	MDS80540-8.....	177
MDS6301407.....	86	MDS6305003.....	86	MDS806250D.....	176
MDS6301408.....	86	MDS6306000.....	86	MDS806300D.....	176
MDS6301501.....	85	MDS6309001.....	86	MDS806400.....	176
MDS6301502.....	85	MDS6310000.....	86	MDS806700.....	176
MDS6301503.....	85	MDS6800.....	70	MDS806800.....	176
MDS6301504.....	85	MDS6801.....	70	MDS806900.....	176

MDS85183.....	176	MDS92160.....	32	MDS9283704.....	88
MDS85189.....	176	MDS92200.....	32	MDS9283801.....	88
MDS85190.....	176	MDS92250.....	32	MDS9283802.....	88
MDS85192.....	176	MDS92260.....	32	MDS9283803.....	88
MDS85196.....	176	MDS92270.....	32	MDS9283804.....	88
MDS851965SH.....	176	MDS92280.....	32	MDS9284001.....	88
MDS85197A.....	176	MDS92290.....	32	MDS9284002.....	88
MDS864104.....	177	MDS92500.....	32	MDS9284003.....	88
MDS864104H.....	177	MDS926101.....	32	MDS9284004.....	88
MDS86410J4.....	177	MDS926102.....	32	MDS9284005.....	88
MDS86410JW4.....	177	MDS926103.....	32	MDS9284006.....	88
MDS86410JW54.....	177	MDS926105.....	32	MDS9284007.....	88
MDS86410W4.....	177	MDS926106.....	32	MDS9284008.....	88
MDS86410W54.....	177	MDS926107.....	32	MDS9284009.....	88
MDS86410XW.....	177	MDS926108.....	32	MDS9284010.....	88
MDS86410XWW.....	177	MDS926109.....	32	MDS92842.....	88
MDS8700.....	33	MDS926110.....	32	MDS9284801.....	89
MDS8701.....	33	MDS926201.....	32	MDS9284802.....	89
MDS8701H.....	33	MDS926202.....	32	MDS9284803.....	89
MDS8701Z.....	33	MDS926203.....	32	MDS9284804.....	89
MDS88000B1H.....	91	MDS926205.....	32	MDS9284901.....	89
MDS88000B91H.....	91	MDS926206.....	32	MDS9380.....	28, 31
MDS88000T2.....	91	MDS926207.....	32	MDS9380LF.....	29, 31
MDS8800B101H.....	91	MDS926208.....	32	MDS9381.....	31
MDS8800FOAMH.....	91	MDS926209.....	32	MDS9386.....	28, 31
MDS8800T32H.....	91	MDS926210.....	32	MDS9386LF.....	29, 31
MDS8800PA28H.....	91	MDS926301.....	32	MDS9387.....	28, 31
MDS88BRUSH.....	92	MDS926302.....	32	MDS9387LF.....	29, 31
MDS88BRUSHSS.....	92	MDS926303.....	32	MDS9388.....	28, 31
MDS88GLUT28MH.....	91	MDS926305.....	32	MDS9388LF.....	29, 31
MDS88INSTWP.....	92	MDS926306.....	32	MDS9389LF.....	29
MDS9104.....	30	MDS926307.....	32	MDS9400.....	29, 31
MDS9106.....	30	MDS926311.....	32	MDS9400LF.....	29, 31
MDS9108.....	30	MDS926312.....	32	MDS9407.....	29, 31
MDS9109.....	30	MDS926501.....	32	MDS9407LF.....	29, 31
MDS9111.....	30	MDS926502.....	32	MDS9410.....	29, 31
MDS9112.....	30	MDS926503.....	32	MDS9411.....	29, 31
MDS9113.....	30	MDS926504.....	32	MDS9413.....	29, 31
MDS9114.....	30	MDS926601.....	32	MDS9414.....	31
MDS9116.....	30	MDS926801.....	32	MDS9415.....	31
MDS9117.....	30	MDS9283001.....	89	MDS9418.....	29, 31
MDS9118.....	30	MDS9283002.....	89	MDS9420.....	31
MDS9119.....	30	MDS9283003.....	89	MDS9421.....	31
MDS9120.....	30	MDS9283004.....	89	MDS9423.....	31
MDS9121.....	30	MDS9283005.....	89	MDS9472.....	29, 31
MDS9123.....	30	MDS9283006.....	89	MDS9475.....	31
MDS9124.....	30	MDS9283101.....	88	MDS9475LF.....	31
MDS9125.....	30	MDS9283102.....	88	MDS9481.....	31
MDS9128.....	30	MDS9283103.....	88	MDS9505.....	31
MDS9133.....	30	MDS9283104.....	88	MDS9505LF.....	31
MDS9136.....	30	MDS9283201.....	88	MDS9506.....	31
MDS9140.....	30	MDS9283202.....	88	MDS9506LF.....	31
MDS91420.....	31	MDS9283203.....	88	MDS9607.....	33
MDS91420LF.....	31	MDS9283204.....	88	MDS9608.....	33
MDS91421.....	31	MDS9283301.....	88	MDS9630.....	31
MDS91421LF.....	31	MDS9283302.....	88	MDS9698L.....	34
MDS91422.....	31	MDS9283401.....	88	MDS9699.....	34
MDS91422LF.....	31	MDS9283402.....	88	MDS9700.....	34
MDS91423.....	31	MDS9283501.....	88	MDS9701.....	34
MDS91423LF.....	31	MDS9283502.....	88	MDS9928.....	33
MDS91424.....	31	MDS9283601.....	88	MDS9929.....	33
MDS91481.....	31	MDS9283602.....	88	MDS9950.....	33
MDS91481LF.....	31	MDS9283701.....	88	MDS9950B.....	33
MDS91482.....	31	MDS9283702.....	88	MDS9950BZ.....	33
MDS92150.....	32	MDS9283703.....	88	MDS9950H.....	33

Item Number Index

MDS9950Z.....	33	MDSSPO2ADL.....	39	MDTCS4P80Z.....	51
MDS9952.....	33	MDSSPO2DISP.....	39	MDTCS4P84.....	51
MDS9952B.....	33	MDSSPO2EAR.....	39	MDTDS4M72.....	51
MDS9952BZ.....	33	MDSSPO2EXT.....	39	MDTDS4M72Z.....	51
MDS9952H.....	33	MDSSPO2PED.....	39	MDTDS4M81.....	51
MDS9952Z.....	33	MDSSPO2PEDDSP.....	39	MDTDS4P72.....	51
MDS9970.....	28	MDSV80534.....	177	MDTDS4P72Z.....	51
MDS9971.....	28	MDSV80535.....	177	MDTDS4P81.....	51
MDS9972.....	28	MDSV80536.....	177	MDTFS4M04.....	51
MDS9973.....	28	MDT10WHT32E.....	60	MDTFS4M04N.....	51
MDS9974.....	28	MDT10WHT34E.....	60	MDTFS4M04NZ.....	51
MDS99901.....	33	MDT10WHT36E.....	60	MDTFS4M04Z.....	51
MDS99901H.....	33	MDT10WHT38E.....	60	MDTFS4M08.....	51
MDS99901Z.....	33	MDT10WHT40E.....	60	MDTFS4M15.....	51
MDS99902.....	33	MDT10WHT42E.....	60	MDTFS4P04.....	51
MDS99902H.....	33	MDT10WHT44E.....	60	MDTFS4P04N.....	51
MDS99902Z.....	33	MDT10WHT46E.....	60	MDTFS4P04NZ.....	51
MDSB10700STND.....	42	MDT10WHT48E.....	60	MDTFS4P04Z.....	51
MDSB10750.....	42	MDT10WHT50E.....	60	MDTFS4P08.....	51
MDSB10900.....	41	MDT10WHT52E.....	60	MDTFS4P08N.....	51
MDSBATTERY42.....	39	MDT10WHT54E.....	60	MDTFS4P15.....	51
MDSBPTUBE.....	39	MDT10WHT56E.....	60	MDTFS4P15N.....	51
MDSECG2003.....	43	MDT10WHT58E.....	60	MDTHR9R04LIL.....	58
MDSECG2012.....	43	MDT11WHT10E.....	59	MDTHR9R04PINK.....	58
MDSECGDE15.....	43	MDT11WHT12E.....	59	MDTHR9R04SAG.....	58
MDSECGDP12.....	43	MDT11WHT14E.....	59	MDTNC4J10BLU.....	52
MDSECGDX12.....	43	MDT11WHT16E.....	59	MDTNC4J10WHT.....	52
MDSECGSE1201.....	42	MDT11WHT18E.....	59	MDTPC4M34.....	51
MDSECGSTAND01.....	42	MDT11WHT20E.....	59	MDTPC4M34Z.....	51
MDSK10321.....	35	MDT11WHT22E.....	59	MDTPC4P34.....	51
MDSK10331.....	35	MDT11WHT24E.....	59	MDTPC4P34Z.....	51
MDSM3A101927.....	39	MDT11WHT26E.....	59	MDTPG6IVCAB.....	57
MDSM3A7208058.....	39	MDT11WHT28E.....	59	MDTPG6IVCAG.....	57
MDSM3A78035.....	39	MDT11WHT2E.....	59	MDTPG6IVCHEV.....	57
MDSM3ABASKET.....	39	MDT11WHT30E.....	59	MDTPG6IVPOP.....	57
MDSM3ABASKET1.....	39	MDT11WHT32E.....	59	MDTPG6IVRRB.....	57
MDSM3ACUFFA.....	39	MDT11WHT34E.....	59	MDTPG6IVRSB.....	57
MDSM3ACUFFC.....	39	MDT11WHT36E.....	59	MDTPG6IVTIDI.....	57
MDSM3ACUFFI.....	39	MDT11WHT4E.....	59	MDTPG7PTSPKT.....	57
MDSM3ACUFFL.....	39	MDT11WHT6E.....	59	MDTPG7RFOFLG.....	58
MDSM3ACUFFT.....	39	MDT11WHT8E.....	59	MDTSG7ITSFLB.....	57
MDSM3ASN.....	39	MDT12WHT34E.....	60	MEN6065.....	130
MDSM3ASN1.....	39	MDT12WHT36E.....	60	MIL1202.....	82
MDSM3ASNT.....	39	MDT12WHT38E.....	60	MIL13015.....	82
MDSM3ASNTC.....	39	MDT12WHT40E.....	60	MIL181107.....	83
MDSM3ASNTCP.....	39	MDT12WHT42E.....	60	MIL181652.....	83
MDSM3ASNTP.....	39	MDT12WHT44E.....	60	MIL18868.....	83
MDSM3ASNTPIR.....	39	MDT12WHT46E.....	60	MIL18SC1652.....	83
MDSM3ASNTT.....	39	MDT12WHT48E.....	60	MIL19102.....	87
MDSM3ASNTTP.....	39	MDT12WHT50E.....	60	MIL19104.....	87
MDSM3ASNTTR.....	39	MDT12WHT52E.....	60	MIL19106.....	87
MDSM3ASTAND01.....	39	MDT12WHT54E.....	60	MIL192100.....	87
MDSM3AWALL01.....	39	MDT12WHT56E.....	60	MIL19282.....	87
MDSM3AWRRNTY.....	39	MDT12WHT58E.....	60	MIL19284.....	87
MDSM611903.....	40	MDT12WHT60E.....	60	MIL19286.....	87
MDSM611905.....	40	MDT214055.....	107	MIL19314.....	87
MDSM611950.....	40	MDT219683.....	50	MIL19316.....	87
MDSM616101.....	40	MDT219683H.....	50	MIL19320.....	87
MDSM616101Z.....	40	MDT221203.....	113	MIL19321BX.....	87
MDSM616201.....	40	MDT221204.....	113	MIL19322.....	87
MDSM616201Z.....	40	MDT221206.....	113	MIL19323.....	87
MDSMRC900.....	41	MDTCS4M80.....	51	MIL19370.....	87
MDSNOTOUCH.....	33	MDTCS4M80Z.....	51	MIL19371.....	87
MDSP3043.....	29	MDTCS4M84.....	51	MIL19385.....	87
MDSPPC107.....	177	MDTCS4P80.....	51	MIL19434.....	87

MIL19580	87	MIL7622	87	MLAB22X221	26
MIL19582	87	MIL842	82	MLAB22X2215	26
MIL19584	87	MIL842TC	83	MLAB22X222	26
MIL3010	87	MIL86	82	MLAB24X401	26
MIL301351	87	MIL86TC	83	MLAB24X4015	26
MIL301382	87	MIL88ATC	83	MLAB24X402	26
MIL301442WL	87	MIL9100	82	MLAB24X501	26
MIL301443WL	87	MIL9104	82	MLAB24X5015	26
MIL301482	87	MILV93060	87	MLAB24X502	26
MIL3015	87	MILV93061	87	MLAB24X601	26
MIL3020	87	MILV940200	82	MLAB24X6015	26
MIL303010	87	MILV940205	82	MLAB24X602	26
MIL303012	87	MILV940207	82	MLAB3501	25
MIL303015	87	MIMO210137814	193	MLAB3515	25
MIL3050	87	MIMO210137843	193	MLAB3525	25
MIL3055	87	MIMO210137846	193	MLAB3561	25
MIL3060	87	MIMO210137847	193	MLAB3571	25
MIL306000	87	MIMO210137853	193	MLABSLIDEF	26
MIL3062	87	MIMO210137855	193	MLABSLIDEFW	26
MIL30650	87	MIMO210137857	193	MLABSLIDEP	26
MIL30965	87	MIMO210137860	193	MLABT15PED	25
MIL31000	82	MIMO210138814	193	MLABT15GRAD	25
MIL33200	83	MIMO210138843	193	MLABT15LGBLB	25
MIL33201	83	MIMO210138846	193	MLABT15NGRAD	25
MIL33202	83	MIMO210138847	193	MLABT17GRAD	25
MIL3331	83	MIMO210138853	193	MMO3418	38
MIL3331HH	83	MIMO210138855	193	MMO3419	38
MIL3332	83	MIMO210138857	193	MMO3420	38
MIL3332HH	83	MIMO210138860	193	MMO3421	38
MIL3333	83	MIM204011	193	MMO3580	38
MIL3333HH	83	MIM204012	193	MMO3581	38
MIL3334	83	MIM21001206Z	44	MMO3583	38
MIL3334HH	83	MIM40000026	44	MMO9167	38
MIL3335	83	MIMM11022	183	MMO9540	38
MIL3335H	83	MIMM9022	183	MMO9711	38
MIL3336	83	MIMM9D022	183	MNK0005	171
MIL3336HH	83	MIT6M066X20V	197	MNK0015	171
MIL3352	83	MIT6M066X28V	197	MNK0077	171
MIL3353	83	MIT6M124X20V	197	MNK0082	171
MIL3354	83	MIT6M124X28V	197	MNK0087	171
MIL3355	83	MIT6M184X20V	197	MNK0089	171
MIL3357	83	MIT6M184X28V	197	MNK1344	171
MIL4110	83	MIT6M242X20V	197	MPH01WML3	175
MIL4111	83	MIT6M242X28V	197	MPH01WML4	175
MIL4115	83	MIT6M300X20V	197	MPH01WML5	175
MIL4410	82	MIT6M300X28V	197	MPH01WML6	175
MIL4411	82	MITQCSF101	196	MPH01WML9PED	174
MIL4415	82	MITQCSF102	196	MPH01WMLA1	175
MIL4510	82	MITQCSF103	196	MPH01WMLA2	175
MIL4511	82	MITQCSF104	196	MPH01WMLA5	175
MIL4515	82	MITQCSF105	196	MPH01WMLA9	175
MIL47	82	MITQCSF106	196	MPH01WMLA9P	174
MIL51000	82	MITQCSF107	196	MPH01WNCS2	174
MIL516	82	MITQCSF108	196	MPH07H499KL	180
MIL5180	83	MITQCSF109	196	MPH07HP4	180
MIL5300	83	MITQCSF110	196	MPH07SD1	181
MIL5302	83	MITQCSF111	196	MPH07SM1	181
MIL5304	82	MITQCSF112	196	MPH08015	187
MIL5306	82	MITQCSF113	196	MPH08017	187
MIL5514	82	MITQCSF114	196	MPH0807SSHB	187
MIL5516	82	MITQCSF115	196	MPH0808SSHB	187
MIL6118	82	MLAB10X75B	25	MPH08301	187
MIL6120	82	MLAB12X75B	25	MPH08741	187
MIL6197	87	MLAB13X100B	25	MPH08752	187
MIL7600	87	MLAB140020	26	MPH08MS02	187

MPH08MS03.....	187	MPHSAFETY281.....	22, 151	MSC097039A.....	74
MPH08SD01C.....	186	MPHSAFETY281Z.....	22, 151	MSC097047.....	75
MPH08SS01C.....	186	MPHSAFETY285.....	22, 151	MSC097048.....	75
MPH11025HCG.....	6, 9	MPHSAFETY285Z.....	22, 151	MSC097049.....	75
MPH12025HCG.....	6, 9	MPHSAFETY28Z.....	22, 151	MSC097062.....	74
MPH1540.....	18	MPHSAFETYBL.....	22, 151	MSC097063.....	74
MPH1550.....	18	MPHSAFETYBLZ.....	22, 151	MSC097065.....	74
MPH1550Z.....	18	MPHSS150.....	186	MSC097070.....	75
MPH1560.....	18	MPHST21.....	22, 151	MSC097071.....	75
MPH1560Z.....	18	MPHST21Z.....	22, 151	MSC097072.....	75
MPH17CD410.....	68	MPHST28.....	22, 151	MSC097618A.....	75
MPH17CD410H.....	68	MPHST28Z.....	22, 151	MSC097688A.....	75
MPH17CE210.....	68	MPHSTSIDE21.....	22, 151	MSC098103.....	77
MPH17CE210H.....	68	MPHSTSIDE21Z.....	22, 151	MSC098104.....	77
MPH17CP610.....	68	MPHSTSIDE28.....	22, 151	MSC098114.....	77
MPH22025HCG.....	6, 9	MPHSTSIDE28Z.....	22, 151	MSC098200A.....	76
MPH30025SA.....	6, 11	MPHTODDLER2.....	22, 151	MSC098203A.....	76
MPH3540.....	17, 19	MPHTODDLER2Z.....	22, 151	MSC098204A.....	76
MPH3550.....	17, 19	MPHUA10SG.....	7	MSC1044EP.....	167
MPH3550FL.....	17, 19	MPHUA120MP.....	7	MSC1065EP.....	167
MPH3550FLZ.....	17, 19	MPHUA2GP.....	7	MSC1066EP.....	167
MPH3550Z.....	17, 19	MPHUA40B.....	7	MSC1102.....	167
MPH3560.....	17, 19	MPHUA50B.....	7	MSC1104.....	167
MPH3560Z.....	17, 19	MPHUAMAC.....	7	MSC1133.....	167
MPH40020HP.....	6, 14	MPHUAURO7.....	7	MSC1133F.....	167
MPH50025M.....	6, 14	MPL1228.....	106	MSC1145.....	167
MPH6540.....	18	MPP100520GS.....	92	MSC1200EP.....	167
MPH6550.....	18	MPP100525GS.....	92	MSC1244EP.....	167
MPH6550Z.....	18	MPP100545GS.....	92	MSC1266EP.....	167
MPH6560.....	18	MPP100546GS.....	92	MSC1500Z.....	171
MPH6560Z.....	18	MPP100550GS.....	92	MSC1505Z.....	171
MPH70025CT.....	7, 10	MPP100555GS.....	92	MSC1510Z.....	171
MPH70025CTM.....	7, 10	MPP100558GS.....	92	MSC1513Z.....	171
MPH70025M.....	7, 10	MPP100565GS.....	92	MSC1523EP.....	167
MPH70025T.....	7, 10	MPP100650N.....	92	MSC1544EP.....	167
MPHADJDEPTH.....	22, 151	MPP100655N.....	92	MSC1823EPZ.....	171
MPHADJDEPTHZ.....	22, 151	MPP100660N.....	92	MSC1834EPZ.....	171
MPHB9VZ.....	67	MPP100670N.....	92	MSC1845EPZ.....	171
MPHBAAAZ.....	18, 67	MPP100680N.....	92	MSC1847EPZ.....	171
MPHBAAZ.....	67	MPP100685N.....	92	MSC18812EPZ.....	171
MPHBCZ.....	67	MPP200215.....	92	MSC20109EP.....	166
MPHBDZ.....	67	MPP200225.....	92	MSC20162EP.....	166
MPHCP610H.....	68	MPP200245.....	92	MSC2033EP.....	166
MPHM191M.....	188	MPP200246.....	92	MSC2044EP.....	166
MPHMMPREEMIE2.....	22, 151	MPP200250.....	92	MSC2066EP.....	166
MPHMMPREEMIE2Z.....	22, 151	MPP200255.....	92	MSC2077EP.....	166
MPHNEWBORN2.....	22, 151	MPP200258.....	92	MSC2244EP.....	166
MPHNEWBORN2Z.....	22, 151	MPP200265.....	92	MSC2266EP.....	166
MPHPREEMIE2.....	22, 151	MPP200270.....	92	MSC2288EP.....	166
MPHPREEMIE2Z.....	22, 151	MPP200275.....	92	MSC2302.....	172
MPHPRESS21.....	22, 151	MPP200410R.....	92	MSC2304.....	172
MPHPRESS21Z.....	22, 151	MPP200410S.....	92	MSC263040H.....	107
MPHPRESS23.....	22, 151	MPP200415R.....	92	MSC2701Z.....	172
MPHPRESS23Z.....	22, 151	MPP200415S.....	92	MSC2703Z.....	172
MPHPRESS26.....	22, 151	MPP200420R.....	92	MSC2705Z.....	172
MPHPRESS26Z.....	22, 151	MPP200420S.....	92	MSC2706ZZ.....	172
MPHPRESS28.....	22, 151	MPP200422R.....	92	MSC2710ZZ.....	172
MPHPRESS28Z.....	22, 151	MPP200422S.....	92	MSC2712ZZ.....	172
MPHPRESS30.....	22, 151	MPP200425R.....	92	MSC281224C.....	106
MPHPRESS30Z.....	22, 151	MPP200425S.....	92	MSC281226C.....	106
MPHSAFETY21.....	22, 151	MRXMX1405H.....	62	MSC281230.....	106
MPHSAFETY21Z.....	22, 151	MSC0002.....	22	MSC281236P.....	106
MPHSAFETY23.....	22, 151	MSC093001.....	171	MSC30410Z.....	172
MPHSAFETY23Z.....	22, 151	MSC093005.....	171	MSC30414Z.....	172
MPHSAFETY28.....	22, 151	MSC097034NA.....	74	MSC3044Z.....	172

MSC3066Z	172	MSC8744EPZ	170	MSG1070	72
MSC3068Z	172	MSC9301EP	170	MSG1075	72
MSC3222Z	172	MSC93025EP	170	MSG1080	72
MSC322505	107	MSC9303EP	170	MSG1085	72
MSC323000C	107	MSC9308EP	170	MSG1090	72
MSC3236Z	172	MSC9405Z	170	MSG2255	72
MSC323950C	107	MSC9410Z	170	MSG2260	72
MSC32410Z	172	MSC9500	170	MSG2265	72
MSC32414Z	172	MSC9501	170	MSG2270	72
MSC3244Z	172	MSC9502	170	MSG2275	72
MSC3245Z	172	MSC9503	170	MSG2280	72
MSC3248Z	172	MSC9504	170	MSG2285	72
MSC3266Z	172	MSC9504YDH	170	MSG2290	72
MSC350320	75	MSC9505	170	MSG6055	72
MSC351000	65	MSC9505YDH	170	MSG6060	72
MSC351010	65	MSC9506	170	MSG6065	72
MSC351128	67	MSC9507	170	MSG6070	72
MSC351200	66	MSC9508	170	MSG6075	72
MSC351210	66	MSC9544EP	167	MSG6080	72
MSC351230	66	MSC9566EP	167	MSG6085	72
MSC351240	66	MSC9588EP	167	MSG6090	72
MSC351300	66	MSC9604EP	168	MSG6555	72
MSC351310	66	MSC9606EP	168	MSG6560	72
MSC351400ANZ	63	MSC9614EP	168	MSG6565	72
MSC351410AN	63	MSC96310	168	MSG6570	72
MSC351420ANZ	63	MSC96314	168	MSG6575	72
MSC351440AN	63	MSC9636	168	MSG6580	72
MSC4002	171	MSC9644EP	167	MSG6585	72
MSC4004	171	MSC9666EP	167	MSG6590	72
MSC4006	171	MSC97310	168	MSG9055	72
MSC4300	169	MSC97314	168	MSG9060	72
MSC4400	169	MSC9736	168	MSG9065	72
MSC4400XL	169	MSC97410	168	MSG9070	72
MSC5422	170	MSC97412	168	MSG9075	72
MSC5444	170	MSC97414	168	MSG9080	72
MSC5466	170	MSC9744EP	168	MSG9085	72
MSC5470	170	MSC9746	168	MSG9090	72
MSC5475	170	MSC9748	168	MSS4011	166
MSC5488	170	MSC9766EP	168	MSS4011H	166
MSC5822	168	MSC9818RZ	169	MSS4022	166
MSC5844	168	MSC9822EPZ	169	MSS4022H	166
MSC5868	168	MSC9845EPZ	169	MSS4044	166
MSC6008H	172	MSC9866EPZ	169	MSS4048	166
MSC6016H	172	MSC98812EPZ	169	MSUSL5627GMDL	107
MSC6433EP	168	MSC9912EPZ	169	MSUSL5628GMDL	107
MSC6444EP	168	MSC9922EPZ	169	MSUSL5687GMDL	107
MSC6455EP	168	MSC9941BLK	77	MSUSL5691GMDL	107
MSC64610EP	168	MSC9941GRA	77	MSUSM5637GMDL	107
MSC64812EP	168	MSC9945EPZ	169	MSUSM5687GMDL	107
MSC68817A	75	MSC9948EPZ	169	MSUSN5667GMDL	107
MSC68824BA	75	MSC9950	77	MSUSN5696GMDL	107
MSC68841A	75	MSC9951	77	MSUSN5698GMDL	107
MSC69024MED	77	MSC9962BLK	77	MSUSN5699GMDL	107
MSC7312EPZ	169	MSC9962GRA	77	MSUSP5695GMDL	107
MSC7318EPZ	169	MSC9966EPZ	169	MSUSP5698GMDL	107
MSC7322EPZ	169	MSC99812EPZ	169	NES30212	206
MSC7344EPZ	169	MSCB281245LB	106	NES31803	206
MSC7348EPZ	169	MSCEX33EP	166	NES35110BX	206
MSC7366EPZ	169	MSCEX44EP	166	NES35170BX	206
MSC8488	170	MSCEX66EP	166	NON00291	204
MSC861X8EPZ	170	MSCEX77EP	166	NON00292	204
MSC8622EPZ	170	MSCEX99EP	166	NON00293	204
MSC8644EPZ	170	MSG1055	72	NON02315	204
MSC871X8EPZ	170	MSG1060	72	NON02325	204
MSC8722EPZ	170	MSG1065	72	NON02350	204

NON024214.....	201	NON21450H.....	157	NON24337.....	52
NON024215.....	201	NON21451.....	157	NON24338.....	52
NON024220.....	201	NON21452.....	157	NON24339A.....	52
NON024230.....	201	NON21453H.....	157	NON24339B.....	52
NON030035.....	200	NON21454H.....	157	NON24339C.....	52
NON03005.....	200	NON21456.....	157	NON24339M.....	52
NON03007.....	200	NON21457.....	157	NON24340.....	52
NON03009.....	200	NON21842Z.....	155	NON24345.....	51
NON03010BAR.....	200	NON23319.....	50	NON24346.....	51
NON03010L.....	200, 202	NON23321.....	50	NON24354.....	56
NON03012.....	200	NON23322.....	50	NON24355.....	56
NON03012L.....	200, 202	NON23323.....	50	NON24356.....	56
NON03014.....	200	NON23324.....	50	NON24356B.....	51
NON03016.....	200	NON23325.....	50	NON24356W.....	51
NON03016L.....	200, 202	NON23326.....	50	NON24357B.....	51
NON0309L.....	200, 202	NON23328.....	50	NON24357W.....	51
NON034215.....	201	NON23339.....	52	NON24358B.....	51
NON0357L.....	200, 202	NON241278.....	106	NON24358W.....	51
NON042000.....	204	NON241278H.....	106	NON24359.....	51
NON042001.....	204	NON241279.....	106	NON24361.....	51
NON0420015.....	204	NON241279H.....	106	NON24398.....	52
NON042002.....	204	NON241280.....	107	NON24398H.....	52
NON0420020.....	204	NON241281.....	107	NON24398NV.....	52
NON042003.....	204	NON241283.....	107	NON24399.....	52
NON042201.....	204	NON241286.....	106	NON245277.....	46
NON05003.....	200	NON241286Z.....	106	NON245277Z.....	46
NON05003BAR.....	200	NON241288.....	106	NON245278.....	46
NON05005.....	200	NON241288H.....	106	NON245288.....	46
NON05007.....	200	NON241289.....	106	NON24770.....	54
NON05009.....	200	NON241289Z.....	106	NON24774.....	54
NON05012.....	200	NON24247.....	56	NON24775.....	54
NON05016.....	200	NON24248.....	56	NON24776.....	54
NON05016BAR.....	200	NON24285.....	165	NON24777V.....	54
NON05020.....	200	NON24286.....	165	NON249005.....	56
NON11460395.....	205	NON24288.....	165	NON250114Z.....	164
NON122424.....	49	NON24288H.....	165	NON250314Z.....	164
NON122433.....	49	NON24289.....	165	NON250318Z.....	164
NON123143.....	49	NON24289H.....	165	NON250412Z.....	164
NON132000206.....	205	NON24290.....	165	NON250414Z.....	164
NON132000595.....	205	NON24291.....	165	NON250501Z.....	164
NON13200177.....	205	NON24291H.....	165	NON25208H.....	155
NON13200178.....	205	NON24292.....	165	NON25212H.....	155
NON13200204.....	205	NON24293.....	165	NON25223Z.....	156
NON13200209.....	205	NON24293H.....	165	NON25224Z.....	156
NON13200211.....	205	NON24294.....	165	NON25249.....	56
NON13200340.....	205	NON24297.....	165	NON25300.....	51
NON13200978.....	205	NON24297H.....	165	NON25312H.....	155
NON13200982.....	205	NON243200.....	50	NON253200CS.....	50
NON13200985.....	205	NON24322.....	50	NON25334Z.....	156
NON13201082.....	205	NON24323.....	50	NON25408H.....	155
NON151717h.....	49	NON24324.....	50	NON25412H.....	155
NON21001Z.....	106	NON24325.....	50	NON25416Z.....	155
NON21002Z.....	106	NON24326.....	50	NON25443Z.....	156
NON21224H.....	156	NON243275.....	46	NON25444H.....	156
NON21334Z.....	156	NON243276.....	46	NON25449.....	156
NON21420H.....	155	NON243277.....	46	NON254572.....	56
NON21422H.....	155	NON24329.....	50	NON25491H.....	158
NON21424H.....	155	NON24330A.....	52	NON25492H.....	158
NON21425Z.....	155	NON24330B.....	52	NON25493H.....	158
NON21426H.....	155	NON24330BLU.....	52	NON25494H.....	158
NON21428H.....	155	NON24333.....	52	NON254955H.....	158
NON21444H.....	156	NON24335.....	52	NON25495H.....	158
NON21446Z.....	156	NON24336.....	52	NON25496H.....	158
NON21448H.....	155	NON24336A.....	52	NON25497H.....	158
NON21449H.....	156	NON24336C.....	52	NON25498H.....	158

NON25499H.....	158	NON270003Z.....	163	NON30233B.....	53
NON25500Z.....	162	NON270012Z.....	163	NON32100.....	52
NON255015H.....	159	NON270101Z.....	163	NON32200.....	52
NON25501H.....	162	NON270102Z.....	163	NON32300.....	51
NON255025H.....	159	NON270103Z.....	163	NON32500.....	51
NON25502Z.....	162	NON270112Z.....	163	NON33025.....	52
NON25504Z.....	162	NON270201SZ.....	163	NON33100.....	52
NON25509Z.....	162	NON270201Z.....	163	NON33100H.....	52
NON25510Z.....	162	NON270202SZ.....	163	NON34000.....	52
NON255125H.....	159	NON270202Z.....	163	NON34600.....	52
NON25513Z.....	162	NON270203Z.....	163	NON34700.....	52
NON255145H.....	159	NON270212Z.....	163	NON37250.....	52
NON25524Z.....	162	NON27036.....	56	NON494109.....	143
NON256000H.....	156	NON27038.....	56	NON494109H.....	143
NON256001Z.....	156	NON27046.....	56	NON7008HDL.....	200
NON25600Z.....	162	NON27048.....	56	NON70600.....	165
NON256015H.....	159	NON27120.....	55	NON70600H.....	165
NON256025H.....	159	NON27122.....	55	NON80328.....	165
NON256125H.....	159	NON27146.....	56	NON80328Z.....	165
NON256130Z.....	162	NON27146SL.....	56	NON80329.....	165
NON256131Z.....	162	NON27146SLXL.....	56	NON80329Z.....	165
NON256145H.....	159	NON27148.....	56	NONASMBLD.....	54
NON25650Z.....	162	NON27148XL.....	56	NONB003.....	203
NON25660Z.....	162	NON27209.....	56	NONBCS09.....	203
NON25700A.....	157	NON27346SL.....	56	NONBCS093.....	203
NON25700H.....	157	NON27375.....	55	NONBCS10.....	203
NON25710A.....	157	NON27378.....	55	NONCPLA-001.....	204
NON25710Z.....	157	NON27382FR.....	55	NONCPLA-002.....	204
NON25720Z.....	157	NON27405.....	55	NONCPLA-003.....	204
NON25744.....	162	NON27408.....	55	NONCXB36H.....	49
NON25750.....	162	NON27408EL.....	55	NONDB4080.....	51
NON25760.....	162	NON27410.....	55	NONDWHC8PE.....	200
NON25812ZZ.....	155	NON27410EL.....	55	NONEMBGDISP.....	165
NON25815.....	47	NON27411AG.....	55	NONEMBGDISPH.....	165
NON25825.....	46	NON27412.....	55	NONFRAME.....	54
NON25850H.....	158	NON27412EL.....	55	NONFS100.....	54
NON25852Z.....	158	NON27710.....	55	NONFS300.....	54
NON25853H.....	158	NON27710EL.....	55	NONFS400.....	54
NON25854H.....	158	NON27712.....	55	NONHCR24LB.....	50
NON25855.....	158	NON27712EL.....	55	NONHCR24MC.....	50
NON25855H.....	158	NON27800.....	46	NONHCR33MC.....	50
NON25861H.....	158	NON27810.....	55	NONHCR40HC.....	50
NON25865H.....	158	NON27810EL.....	55	NONL16BL.....	200, 202
NON260001DZ.....	163	NON27812EL.....	55	NONLENS.....	54
NON260002DZ.....	163	NON27SMS2.....	53	NONP004.....	203
NON260301H.....	164	NON27SMS2XL.....	53	NONP013.....	203
NON260302H.....	164	NON27SMS2XXXL.....	53	NONPBM350.....	46
NON260302Z.....	164	NON28234P.....	53	NONPBM800B.....	46
NON260303HH.....	164	NON28533B.....	53	NONPBM800N.....	46
NON260315H.....	164	NON28625Z.....	54	NONPLA12.....	201
NON26800.....	46	NON28626Z.....	54	NONPLA16.....	201
NON26800H.....	46	NON28628Z.....	54	NONPLA2.....	201
NON26801.....	46	NON28752.....	54	NONPLA20.....	201
NON26805.....	46	NON28758.....	54	NONPLA24.....	201
NON26810.....	47	NON28759.....	54	NONPLA4.....	201
NON26810H.....	47	NON28800.....	46	NONPLA9.....	201
NON26813.....	47	NON28805.....	46	NONPLADL.....	201-202
NON26813H.....	47	NON28830.....	46	NONPLAEC6.....	203
NON26818.....	47	NON28852.....	54	NONPLAEC8.....	203
NON26820.....	47	NON28858.....	54	NONPLAFL.....	201-202
NON26830.....	46	NON28859.....	54	NONPLASL.....	201-202
NON26835.....	47	NON29858.....	54	NONPLAYL.....	201-202
NON26871.....	46	NON29859.....	54	NONPLC12.....	201
NON270001Z.....	163	NON30223B.....	53	NONPLC12DW.....	201
NON270002Z.....	163	NON30232B.....	53	NONPLC16.....	201

NONPLC16DW.....	201	NYOV560.....	122	ORT19400RXS.....	116
NONPLC20.....	201	NYOV562.....	122	ORT27200.....	116
NONPLC20DW.....	201	NYOV564.....	122	ORT27200S.....	116
NONPLC8.....	201	NYOV566.....	122	ORT27200T.....	116
NONPLC8DW.....	201	NYOV568B.....	122	ORT27220.....	116
NONPLCLID.....	201-202	NYOY101.....	120	ORT27220S.....	116
NONPLCLID20DW.....	201-202	NYOY102.....	120	ORT27220T.....	116
NONPLCLID8.....	200-202	NYOY203.....	120	ORT28100L.....	117
NONTGG39X.....	49	NYOY220.....	120	ORT28100M.....	117
NONTGG58X.....	49	NYOY306.....	120	ORT28100S.....	117
NONTH150.....	53	NYOY420.....	120	ORT28100XL.....	117
NONTH180.....	53	NYOY503.....	121	ORT28100XS.....	117
NONTH200.....	53	NYOY510.....	121	ORT28200L.....	117
NONUNNAS130.....	169	NYOY522.....	121	ORT28200M.....	117
NONUNNAS140.....	169	NYOY524.....	120	ORT28200S.....	117
NONUNNAS30.....	169	NYOY526.....	121	ORT28200XL.....	117
NONUNNAS40.....	169	NYOY530.....	120	ORT28200XS.....	117
NONZIP1012.....	205	NYOY532.....	121	ORT28350L.....	116
NONZIP1215.....	96, 205	NYOY534.....	121	ORT28350M.....	116
NONZIP1215Z.....	96	NYOY538B.....	121	ORT28350S.....	116
NONZIP35.....	96, 205	NYOY540.....	120	ORT28350XL.....	116
NONZIP35Z.....	96	NYOY542.....	120	ORT28400L.....	116
NONZIP48.....	96, 205	NYOY546.....	121	ORT28400S.....	116
NONZIP58.....	96, 205	NYOY548B.....	120	ORT28610L.....	117
NONZIP58Z.....	96	OAT2583.....	70	ORT28610M.....	117
NONZIP69.....	96, 205	OAT2584.....	70	ORT28610S.....	117
NONZIP69Z.....	96	OAT2585.....	70	ORT28610XL.....	117
NONZIP810.....	96, 205	OAT2586.....	70	ORT28610XS.....	117
NONZIP810Z.....	96	OAT2587.....	70	ORT28620L.....	117
NONZIP912.....	96, 205	OAT6800.....	70	ORT28620M.....	117
NONZIP912Z.....	96	OAT6801.....	70	ORT28620S.....	117
NPMCI1450.....	147	OAT6802.....	70	ORT28620XL.....	117
NPMCI1632.....	147	OAT6803.....	70	ORT28620XS.....	117
NPMCI1650.....	147	OAT6804.....	70	ORT28630L.....	117
NPMCI1832.....	147	OFX00019.....	206	ORT28630M.....	117
NPMCI1850.....	147	OFX00020.....	206	ORT28630S.....	117
NPMCI2025.....	147	OFX20015445.....	206	ORT28630XL.....	117
NPMCI2032.....	147	ORA10010292Z.....	111	ORT28630XS.....	117
NPMCI2050.....	147	ORA36M1C.....	111	ORT28640L.....	117
NPMCI2225.....	147	ORT11100L.....	115	ORT28640M.....	117
NPMCI2232.....	147	ORT11100M.....	115	ORT28640S.....	117
NPMCI2420.....	147	ORT11100S.....	115	ORT28640XL.....	117
NPMCI2425.....	147	ORT11100XL.....	115	ORT28640XS.....	117
NPMNM20G25Z.....	143	ORT11100XS.....	115	ORT30300ML.....	116
NPMNM20G38Z.....	143	ORT11300L.....	115	ORT30300MM.....	116
NPMNM21G25Z.....	143	ORT11300M.....	115	ORT30300MS.....	116
NPMNM21G38Z.....	143	ORT11300S.....	115	ORT30300MXL.....	116
NPMNM22G25Z.....	143	ORT11300XL.....	115	ORT30300WL.....	116
NPMNM22G38Z.....	143	ORT11300XS.....	115	ORT30300WM.....	116
NPMSPR21G19.....	146	ORT13300L.....	115	ORT30300WS.....	116
NPMSPR23G19.....	146	ORT13300M.....	115	ORT30300WXS.....	116
NPMSPR25G19.....	146	ORT13300MXL.....	115	ORT32200L.....	112
NYO140J.....	122	ORT13300S.....	115	ORT32200M.....	112
NYO420Y.....	121	ORT13300SXL.....	115	ORT32200S.....	112
NYO520Y.....	121	ORT18000L.....	116	ORT32700.....	112
NYO532Y.....	120	ORT18000R.....	116	ORT33100LM.....	112
NYO540Y.....	121	ORT19400LL.....	116	ORT33100RM.....	112
NYO544Y.....	120	ORT19400LM.....	116	ORT33200LM.....	112
NYO548Y.....	120	ORT19400LS.....	116	ORT33200LS.....	112
NYOJ100.....	121	ORT19400LXL.....	116	ORT33200RM.....	112
NYOJ120.....	122	ORT19400LXS.....	116	ORT33200RS.....	112
NYOJ140.....	122	ORT19400RL.....	116	OSOM101.....	9
NYOJ640.....	122	ORT19400RM.....	116	OSOM102.....	9
NYOV552.....	122	ORT19400RS.....	116	OSOM134.....	9
NYOV554.....	122	ORT19400RXL.....	116	OSOM145.....	14

OSOM147.....	11	PAIN9012S.....	110	POD14208SM.....	118
OSOM149.....	11	PAIN9013S.....	109	POD14209LXL.....	118
OSOM181.....	15	PAIN9014S.....	109	POD14209SM.....	118
OSOM183.....	15	PAIN9016.....	109	POD14210L.....	118
OTC000252.....	131	PAIN9019S.....	110	POD14210M.....	118
OTC006009.....	136	PAIN9020S.....	110	POD14210S.....	118
OTC010002.....	135	PAINDRAPE1.....	108	POD14211N.....	118
OTC010015.....	131	PAINDRAPE2.....	108	POD14211W.....	118
OTC052547.....	135	PAINDRAPE3.....	108	POD142125M.....	118
OTC071710.....	131	PAINDRAPE4.....	108	POD14212LXL.....	118
OTC073230.....	131	PAINLOR10LL25.....	108	POD14213L.....	118
OTC105656.....	135	PAINLOR10LS25.....	108	POD14213S.....	118
OTC110231.....	136	PAINLOR5LL25.....	108	PRM086001.....	161
OTC211402.....	131	PAINLOR5LS25.....	108	PRM086002.....	161
OTC211602.....	131	PBOBSP510.....	11	PRM086003.....	161
OTC27716.....	132	PCS1550.....	150	PRM086004.....	161
OTC430122.....	14	PCS1650.....	150	PRM086006.....	161
OTC571135.....	136	PEDP36.....	185	PRM088001.....	161
OTC589430.....	134	PGL050.....	166	PRM088001CP.....	161
OTC664524CS.....	132	PGL050H.....	166	PRM088002.....	161
OTC92001A.....	11	PGL050TUBE.....	166	PRM088003.....	161
PAIN0609.....	109	PGL050TUBEH.....	166	PRM088004.....	161
PAIN0611.....	110	PGL400.....	166	PRM088006.....	161
PAIN8001.....	111	PGL400H.....	166	PRM08802CP.....	161
PAIN8002.....	111	PHY0403000149.....	41	PRM08803CP.....	161
PAIN8003.....	111	PHY0427000134.....	41	PRM21208CZ.....	155
PAIN8004.....	111	PHY80403148.....	41	PRM21212CZ.....	155
PAIN8005.....	111	PHY9425000023.....	41	PRM21224.....	156
PAIN8006.....	111	PINK2583.....	71	PRM21312CZ.....	155
PAIN8007.....	111	PINK2584.....	71	PRM21408CZ.....	155
PAIN8008.....	111	PINK2585.....	71	PRM21412CZ.....	155
PAIN8009.....	111	PINK2586.....	71	PRM21416CZ.....	155
PAIN8010.....	111	PINK2587.....	71	PRM21419.....	155
PAIN8011.....	111	PINK5083.....	71	PRM21423.....	155
PAIN8012.....	111	PINK5084.....	71	PRM21426.....	155
PAIN8013.....	111	PINK5085.....	71	PRM21428.....	155
PAIN8014.....	111	PINK5086.....	71	PRM21444.....	156
PAIN8015.....	111	PINK5087.....	71	PRM21444H.....	156
PAIN8016.....	111	PINK6073.....	71	PRM21450H.....	157
PAIN8017.....	111	PINK6074.....	71	PRM21450Z.....	157
PAIN8018.....	111	PINK6075.....	71	PRM21451.....	157
PAIN8019.....	111	PINK6076.....	71	PRM21453H.....	157
PAIN8021.....	111	PINK6077.....	71	PRM21453Z.....	157
PAIN8022.....	111	PINK6800.....	71	PRM21454Z.....	157
PAIN8023.....	111	PINK6801.....	71	PRM21456.....	157
PAIN8024.....	111	PINK6802.....	71	PRM2208Z.....	155
PAIN8025.....	111	PINK6803.....	71	PRM22DS.....	155
PAIN8028.....	111	PINK6804.....	71	PRM25224.....	156
PAIN8029.....	111	PLI3801Z.....	40	PRM25224H.....	156
PAIN8030.....	111	PLI3803Z.....	40	PRM254192.....	158
PAIN8031.....	111	PM1622-22.....	50	PRM254193.....	158
PAIN8032.....	111	PM1824-15.....	50	PRM25444.....	156
PAIN8033.....	111	PM1824-15H.....	50	PRM25444H.....	156
PAIN8034.....	111	PM2127-18.....	50	PRM25500H.....	163
PAIN8036.....	111	PM2127-22.....	50	PRM256000H.....	156
PAIN9001S.....	110	PM2127-22H.....	50	PRM256001H.....	156
PAIN9002S.....	110	PMU9034000.....	111	PRM25600H.....	163
PAIN9003S.....	110	PMU9034002.....	111	PRM25855.....	158
PAIN9004S.....	110	POD14201.....	118	PRM25865.....	158
PAIN9005S.....	110	POD14202.....	118	PRM25865H.....	158
PAIN9006S.....	110	POD14203.....	118	PRM25865P.....	158
PAIN9007S.....	110	POD14204.....	118	PRM260001H.....	163
PAIN9009S.....	110	POD14206.....	118	PRM260001HH.....	163
PAIN9010S.....	110	POD14207LXL.....	118	PRM260002H.....	163
PAIN9011S.....	110	POD14207SM.....	118	PRM260003H.....	163

PRM260101H.....	164	SMU50150.....	206	SYR105010.....	140
PRM260102H.....	164	SNAP66PSFBLM.....	197	SYR105020.....	140
PRM260103.....	164	SNAP66PSFOAT.....	197	SYR105205.....	140
PRM260201H.....	163	SNAP66PSFSPA.....	197	SYR105207.....	140
PRM260202H.....	163	SNAP66PUTEDEN.....	197	SYR105217.....	140
PRM3312Z.....	155	SNT91000300.....	39	SYR110010.....	140
PRM354194.....	158	STAPLER35RB.....	90	SYR110020.....	140
PRM354195.....	158	STAPLER35W.....	90	SYR110205.....	140
PRM354196.....	158	STAPLER35WA.....	90	SYR110207.....	140
PRM354197.....	158	STAPLER35WPB.....	90	SYR110215.....	140
PRM354198.....	158	SUP10FTSTR.....	197	SYR110217.....	140
PRM354199.....	158	SUP2FTSTR.....	197	SYR110225.....	140
PRM4408Z.....	155	SUP2X2FTW45.....	197	SYR110227.....	140
PRM4412H.....	155	SUP2X2FTW90.....	197	SYR110235.....	140
PTH1710.....	16	SUP3FTSTR.....	197	SYR120010.....	140
PTH2700.....	16	SUP4FTSTR.....	197	SYR120020.....	140
PTH2713.....	16	SUP5FTSTR.....	197	SYR130010.....	140
PTH2729.....	16	SUP6FTSTR.....	197	SYR130020.....	140
PTH2750.....	16	SUP7FTSTR.....	197	SYR160010.....	140
PTH3021.....	18	SUS1020B.....	107	SYR160020.....	140
PTH3024.....	18	SUS1085B.....	107	SYRA101261.....	141
QDL0272.....	9	SUS1282B.....	107	SYRA101272.....	141
QDL0343.....	11	SUS1291B.....	107	SYRS100185.....	140
QDL0345.....	11	SUS493B.....	107	SYRS100187.....	140
QDL0354.....	11	SUS494B.....	107	SYRS100197.....	140
QDL20108.....	11	SUS750B.....	107	SYRS100205.....	140
QDL20109.....	9	SUS928B.....	107	SYRS100207.....	140
QDL20110A.....	9	SUS931B.....	107	SYRS100215.....	140
QDL20183.....	11	SUSSP104BX.....	107	SYRS100217.....	140
QDL20193.....	14	SWD216033.....	143	SYRS100225.....	140
QDL20218.....	12	SWD216066.....	143	SYRS100227.....	140
QDL20227.....	11	SWD225190.....	145	SYRS100235.....	140
QDL20260.....	12	SWD225257Z.....	144	SYRS100237.....	140
QDL20274.....	12	SWD225281.....	145	SYRS100253.....	140
QDL20299.....	12	SWD225299.....	145	SYRS100255.....	140
QDL281.....	9	SWD610102.....	143	SYRS100257.....	140
QMX144001.....	8	SYR100165.....	140	SYRS100272.....	140
QRS70000500.....	42	SYR100167.....	140	SYRS100273.....	140
QRS70100000.....	42	SYR100185.....	140	SYRS100276.....	140
QRSZ50002608.....	44	SYR100187.....	140	SYRS100302.....	140
QRSZ70000101.....	44	SYR100195.....	140	SYRS101235F.....	140
QRSZ70000301.....	42	SYR100197.....	140	SYRS101255F.....	140
RAYCR2032.....	17, 19	SYR100205.....	140	SYRS101272F.....	140
RAYCR2032.....	18	SYR100207.....	140	SYRS101392F.....	140
RDI30295.....	150	SYR100215.....	140	SYRS103205.....	140
RDI30296.....	150	SYR100225.....	140	SYRS103207.....	140
S-D297000Z.....	44	SYR100227.....	140	SYRS103215.....	140
S-D298020.....	44	SYR100235.....	140	SYRS103217.....	140
S-D298040050.....	44	SYR100252.....	140	SYRS103225.....	140
S-03015100.....	17	SYR100253.....	140	SYRS103227.....	140
S-03015200.....	17	SYR100255.....	140	SYRS103235.....	140
S-03025050.....	17	SYR100257.....	140	SYRS103253.....	140
S-OG3000CONTC.....	17	SYR100302.....	140	SYRS103255.....	140
SAL5415047.....	150	SYR101010.....	140	SYRS105205.....	140
SCA201.....	181	SYR101020.....	140	SYRS105207.....	140
SCA206.....	181	SYR101272.....	140	SYRS105215.....	140
SCA213.....	182	SYR101282.....	140	SYRS105227.....	140
SCA217.....	182	SYR103010.....	140	SYRS110215.....	140
SDCEMS012BX.....	182	SYR103020.....	140	SYRS110227.....	140
SDCML020127.....	62	SYR103205.....	140	SYRS1101292.....	140
SDCML020137.....	62	SYR103215.....	140	SYRS1193292.....	140
SDCPP012.....	91, 182	SYR103217.....	140	SYRTCAPZ.....	140
SDCPP052.....	182	SYR103225.....	140	TRI66140.....	83
SDCPS012BX.....	182	SYR103227.....	140	TRI66145.....	83
SDCPS052BX.....	182	SYR103253.....	140	TRI66155.....	83

TRI66195.....	83	UNV35267.....	208	W-A52434UZ.....	38
TRI66590.....	83	UNV35662.....	208	W-A59000.....	98
TRI66790.....	83	UNV35668.....	208	W-A59000Z.....	98
TRI66855.....	83	UNV35672.....	208	W-A59001.....	98
TRU8845.....	206	UNV36301.....	208	W-A59001Z.....	98
TUN1730.....	183	UNV40000.....	208	W-A59004CSZ.....	98
ULT101253.....	140	UNV40304.....	207	W-A59010.....	98
UME4040650100.....	190	UNV40305.....	207	W-A590XS.....	98
UME4040650200.....	190	UNV43118.....	207	W-A71000A.....	37
UME4040650300.....	190	UNV43623.....	207	W-A71900.....	37
UME4070650100.....	190	UNV43651.....	208	W-A72200.....	38
UME4070650200.....	190	UNV43663.....	207	W-A72300.....	38
UME4070650300.....	190	UNV43671.....	208	W-A77710.....	37
UNV00133.....	207	UNV46200.....	208	W-A777912MPX.....	37
UNV00164.....	207	UNV46300.....	208	W-A77791MX.....	37
UNV00700.....	207	UNV55400.....	208	W-A777921MP.....	37
UNV01117.....	207	UNV66600.....	208	W-A77792M.....	37
UNV07051.....	208	UNV72210.....	207	W-A791MXPROBP.....	37
UNV07071.....	208	UNV72220.....	207	W-A80010.....	98
UNV08850.....	208	UNV75602.....	207	W-ACP1501ENB.....	43
UNV08851.....	208	UNV79000.....	207	W-ACP150A1ENB.....	43
UNV08860.....	208	UNV83410.....	207	W-ACPRUIUBD.....	43
UNV08861.....	208	UNV83436.....	207	W-ACPRUNUBD.....	43
UNV10200.....	207	UNV86920.....	208	WCA6100.....	176
UNV10200VP3.....	207	UNV92009.....	207	WCA6100.....	176
UNV10210.....	207	UNV95200.....	208	WCA806991.....	176
UNV10220.....	207	UNV96920.....	208	WNCPC24C18.....	200
UNV10506.....	208	URO3PVC1016C.....	102	WNCPC1618.....	200
UNV10630.....	208	URO3PVC1016L.....	102	WNCPC1618CC.....	200
UNV10700.....	207	URO3PVC1216C.....	102	WNCPC20J16.....	200
UNV11289.....	208	URO3PVC1216L.....	102	WNCPC32TJ32.....	200
UNV12113.....	208	URO3PVC1416C.....	102	WNCPC6.....	200
UNV14113.....	208	URO3PVC1416L.....	102	WNCPF12.....	203
UNV14115.....	208	URO3PVC1616C.....	102	WNCPF6.....	203
UNV14215.....	208	URO3PVC1616L.....	102	WNCPF8.....	203
UNV15001.....	207	URO3PVC1816C.....	102	WNCPFL6V.....	203
UNV15113.....	208	URO3PVC1816L.....	102	WNCPFL8V.....	203
UNV15510.....	208	URO4PVCPC12MC.....	103	WSD908015ZZ.....	98
UNV16113.....	208	URO4PVCPC12MH.....	103	ZOL21110201CC.....	40
UNV20030.....	208	URO4PVCPC16MC.....	103	ZOL21110201PK.....	40
UNV20435.....	208	URO4PVCPC16MH.....	103	ZOL8000080701.....	40
UNV20630.....	208	URO5454.....	103	ZOL8900080001.....	40
UNV20741.....	207	URO70283.....	112	ZOL8900081001.....	40
UNV20742.....	207	W-A01690200.....	34		
UNV20743.....	207	W-A02892000.....	34		
UNV20744.....	207	W-A02893000.....	34		
UNV20745.....	207	W-A03000U.....	38		
UNV20746.....	207	W-A03100U.....	38		
UNV20747.....	207	W-A05031750.....	34		
UNV20748.....	207	W-A05031750Z.....	34		
UNV20951.....	207	W-A06000005ZZ.....	33		
UNV20952.....	207	W-A06000200.....	33		
UNV20961.....	207	W-A06000300.....	33		
UNV20962.....	207	W-A06500U.....	38		
UNV20971.....	207	W-A105660.....	43		
UNV20972.....	207	W-A11710.....	37		
UNV20981.....	207	W-A11720.....	37		
UNV20982.....	207	W-A150AS1ENB.....	43		
UNV20991.....	207	W-A20200.....	37		
UNV20992.....	207	W-A21501.....	37		
UNV20994.....	207	W-A21504.....	38		
UNV20995.....	207	W-A23820.....	38		
UNV27410.....	208	W-A25020.....	38		
UNV35210.....	208	W-A52401H.....	38		
UNV35264.....	208	W-A52432UZ.....	38		

A

- Abbott, 20, 150
 - Abdominal Pads, 157
 - Absorbent Pads, 157
 - Non-Sterile, 157
 - Sterile, 157
 - Acme United, 68
 - Adhesive Bandages, 162, 163
 - CURAD, 162
 - Fabric, 162
 - Kid-Friendly, 162
 - Plastic, 162, 163
 - Adhesive Island Wound Dressing, 172
 - Adhesive Remover, 94
 - Adhesive Tapes, 163, 164
 - Cloth, 163, 164
 - CURAD, 163
 - Paper, 163
 - Sport Tape, 164
 - Transparent, 163
 - Wound Closure Strips, 164
 - Advanced Wound Care, 166–172
 - Alginates, 169
 - Antimicrobial Silver Dressings, 170
 - Bordered Gauze, 172
 - Burn Dressings, 166
 - Collagen, 170
 - Composite Dressings, 172
 - Compression Systems, 169, 170
 - Dressing Retention Tape, 171
 - Foam Dressings, 166–168
 - Gelling Fibers, 169
 - Hydrocolloids, 170
 - Post-Op Dressings, 168
 - Silicone Contact Layer, 171
 - Skin Protectants, 171
 - Super Absorbent Polymer, 168
 - TheraHoney Dressings, 171
 - Transparent Film Dressings, 172
 - Waterproof Composite Dressings, 172
 - Wound Cleansers, 172
 - Wound Device, 166
 - Wound Dressings, 166, 172
 - Wound Gel, 170
 - AJM Packaging, 203, 205
 - Alcohol, 94
 - Prep Pads, 94
 - Swabsticks, 94
 - Alere, 9, 11–12, 14–17
 - Alginates, 169
 - AloeTouch Exam Gloves, 70
 - AloeTouch Ice Exam Gloves, 70
 - Aneroids, 29, 31
 - Handheld, 29, 31
 - Mobile, 29
 - Wall-Mounted, 29, 31
 - Anesthesia, 108–111
 - Ankle, 116, 117
 - Air/Gel, 116
 - Air Splints, 116
 - Brace, 116
 - Dorsal Night Splint, 116
 - Foam Splints, 116
 - Lacer, 116
 - Stirrup, 116
 - Ankle (continued)
 - Walkers, 117
 - Antibacterial Soap, 76
 - Anti-Embolic Stockings, 115
 - Antimicrobial Silver Dressings, 169, 170
 - Argales Powder, 170
 - Gelling Fibers, 169
 - Hydrogel, 170
 - Antimicrobial Wound Gel, 154
 - Apparel, Patient, 56
 - Disposable, 56
 - Capes, 56
 - Exam Shorts, 56
 - Gowns, 56
 - Robes, 56
 - Apparel, Protective, 53–55
 - Disposable Bouffant Caps, 53
 - Disposable Boundary Bouffant Caps, 53
 - Disposable Cover Gowns, 53
 - Disposable Isolation Gowns, 53
 - Eyewear, 54
 - Face Mask, 55
 - Face Shields, 54
 - Lab Goggles, 54
 - Safety Glasses, 54
 - Shoe Covers, 54
 - Shoe Covers, Disposable, 54
 - Surgeon's Cap, 54
 - Applicators, 96, 97
 - Cotton-Tipped, 96
 - Jars, 97
 - OB/GYN Swab, 96
 - Procto Swab, 96
 - Wooden Spoons, 97
 - Wooden Stick, 97
 - Applicators & Tongue Blades, 96, 97
 - Argales Powder, 170
 - Arm Boards, 165
 - MRI-Safe, 165
 - Arm Slings, 115
 - ave, 61
- ## B
- Bacitracin, 154
 - Bags, 24, 49–50, 96, 205
 - Biohazard, 24, 49
 - High Density, 50
 - Low Density, 49
 - Paper, 205
 - Patient Belonging, 24
 - Plastic, 96, 205
 - Red Biohazard, 24
 - Specimen, 24
 - STAT, 24
 - Valuables, 24
 - Write-On, 96
 - Zip Closure, 96, 205
 - Bandages, 158, 160–161
 - CoFlex, 161
 - Conforming, 158
 - Elastic, 160
 - Gauze, 100% Cotton, 158
 - Self-Adherent, 161
 - Wrap, 161
 - Batteries, 67
 - Batteries (continued)
 - Alkaline, 67
 - Duracell Procell, 67
 - Hearing Aid, 67
 - MedCell Advantage, 67
 - Baxter, 141, 150
 - BD Medical, 12, 141, 142–148
 - Beckman Coulter, 10
 - Benzalkonium Chloride Towelettes, 94
 - better earth, 201–204
 - Billeau Ear Loops, 87
 - Binders, 207
 - Bins, 174
 - Hanging, 174
 - Stacking, 174
 - Biohazard Bags, 24
 - Biopsy Punches, 81, 83–84
 - Dermal Curettes, 83
 - Bladder Scanners, 41, 42
 - BioCon 900, 42
 - Blade Glove Cartridges, 89
 - Blankets, 51
 - Disposable, 51
 - Blanket Warmers, 174
 - Blickman, 187
 - Blood Collection, 143–145
 - Low Volume, 148
 - Needles, 143, 144
 - Tube Holder, 143, 144
 - Sets, 144, 145
 - Tubes, 148
 - Blood Draw, 143–145
 - Blood Collection Needles, 143, 144
 - Blood Collection Sets, 144, 145
 - Infusion Sets, 144, 145
 - Luer Adapters, 144
 - Blood Glucose Testing, 16–19
 - Blood Pressure, 29–31
 - Aneroids, 29–31
 - Digital Monitors, 28, 29
 - Monitors, 28
 - Semi-Automatic, 28
 - Sphygmomanometers, 30
 - Wrist Monitor, 29
 - Blood Testing, 19, 20
 - Bordered Gauze, 172
 - Bovie, 184, 189
 - Bowls, 203
 - Disposable, 203
 - Bracket, Wall Mount, 67
 - Braun, 146, 150
 - Breakroom Supplies, 206
 - Coffee Creamers, 206
 - Coffee Maker, 206
 - Coffee Pods, 206
 - K-Cups, 206
 - Keurig Brewer, 206
 - Stirrers, 206
 - Sugar, 206
 - Breast Cancer Awareness, 71
 - Generation Pink, 71
 - Generation Pink 3G, 71
 - Generation Pink Pearl, 71
 - Generation Pink Sense, 71
 - Brewer, 188, 191, 194, 206
 - Brushes, 92, 95

- Brushes (continued)
 Instrument Cleaning, 92
 Surgical Scrub, 95
 Bunion Guard, Gel, 118
 Burn Dressings, 166
 Burton, 188, 189
- C**
- Cabinets, 174
 Narcotics, 174
 Calculi Strainer, 105
 Cannulas, 87, 99
 Endo Biopsy, 87
 Oxygen Therapy, 99
 Cantel Medical, 62, 91, 182
 Carts, 174, 175
 Medical, 174, 175
 Pediatric, 174
 Casting Tape, 113
 Cast Padding, 113
 Catheters & Drain Bags, 102–105
 Catheter Plugs, 104
 Catheter Straps, 104
 Clear Vinyl, 102
 Coude-Tip, 102
 Drain Bags, 104
 Drain Tube Protector, 104
 Foley, 102
 Hydrophilic Straight, 102
 Insertion Tray, 103
 Intermittent, 102, 103
 Intermittent Catheterization Trays, 104
 Latex, 102
 Leg Bag Extension Tubing, 104
 Leg Bags, 103
 Pre-Connected Intermittent Catheter Trays, 104
 Red Rubber All-Purpose Intermittent, 103
 Silicone-Elastomer Coated, 102
 Speci-Cath Kits, 105
 Touch-Free Self Catheter System, 103
 Vinyl Coude, 102
 Cervical Collars, 115
 Serpentine, 115
 Chairs, 195
 Lobby and Reception, 195
 Phlebotomy, 195
 Charts, 181
 Boys and Girls Growth, 181
 Chemicals, Lab, 26
 Cholestech LDX Systems, 16
 Accessories, 16
 Capillary Plungers, 16
 Capillary Tubes, 16
 Optics Cassette w/Case, 16
 Pipette Tips, 16
 Thermal Printer Labels, 16
 Cholesterol Testing, 16
 Test Cassettes, 16
 Cleaners and Disinfectants, 62, 63–66
 Detergent, 62
 Disinfectant, 62
 Micro-Kill AF² Wipes, 65
 Micro-Kill AF+ Wipes, 66
 Cleaners and Disinfectants (continued)
 Micro-Kill Bleach Solutions, 64
 Micro-Kill Bleach Wipes, 63
 Micro-Kill One, 66
 Sanitizer, 62
 Test Strips, 62
 Cleaning and Sterilization, 91
 Detergent, 91
 Presoaks, 91
 Cleansers, 172
 Wound Cleansers, 172
 Cleansing Towelettes, 94
 Clear Jars, 97
 Clinton Industries, 185, 186, 188, 192, 194–195
 Clips, 207
 Coag-Sense, 19
 Coffee, 206
 Creamers, 206
 Maker, 206
 Pods, 206
 Stirrers, 206
 Sugar, 206
 Co-Flex, 161
 AFD, 161
 Cold Packs, 97
 Collagen, 170
 Composite Dressings, 172
 Compression Stockings, 114
 Compression Systems, 169, 170
 Containers, 24, 25, 112, 152
 Biohazard, 152
 Denture, 112
 General Use, 25
 Multipurpose, 152
 Pathology, 24
 Pre-Filled Formalin, 24
 Sharps, 152
 Mail Back, 152
 Patient Room, 152
 PG-II, 152
 Phlebotomy, 152
 Specimen, 25
 Copy Paper, 208
 Cotton Balls, 96
 Cotton Gauze Sponges, 155
 Cotton-Tipped Applicators, 96
 Cover Glass, 26
 Microscope, 26
 Covers, Equipment, 52
 Disposable, 52
 Covidien, 107
 Cream, 206
 Crutches, 176, 177
 3-in-1 Quick Fit, 177
 Basic, 176
 Standard, 177
 Crutch Pillows, 177
 Cryosurgery Device, 184
 Cryosurgical Systems, 111
 Disposable, 111
 Portable, 111
 Cubicle Curtains, 196, 197
 Mesh Panels, 197
 Snap Panels, 197
 Cubicle Track System, 197
 Cups, 200–202
 Disposable, 200–202
 Lids, 202
 Portion/Souffle, 201
 CURAD, 71, 154, 157, 159, 162–165
 Curettes, 87
 Curettes-Excavators, 82
 Curtains, 196, 197
 Cubicle, 196, 197
 Cubicle Track System, 197
 Disposable, 196
 Mesh Panels, 197
 Snap Panels, 197
 Cuticle Nippers, 81
- D**
- Dart Solo, 205
 Dazor Lighting, 188
 Defibrillators, 40, 41
 Defibtech, 41
 Denture Containers, 112
 Dermal Curettes, 81, 83
 Desiccators, Electrosurgery, 184
 Diagnostic Integrated Wall System, 36, 37
 Accessories, 36, 37
 Aneroid, 36, 37
 Ophthalmoscopes, 36, 37
 Otoscopes, 36, 37
 Specula Dispenser, 36, 37
 Thermometers, 36, 37
 Wall Transformer, 36, 37
 Digital, 28, 29, 33–34, 44, 179–181
 Baby Scale, 181
 Blood Pressure Monitors, 28, 29
 Diaper Scale, 181
 Dietary Scales, 180
 Eye- and Waist-Level Scales, 180
 Floor Level Scales, 180
 Physician Scales, 180
 Platform Scales, 179
 Scales, 179
 Spirometers, 44
 Step-On Scales, 180
 Thermometers, 33, 34
 Wrist Monitor, 29
 Disinfectant, 91
 OPA Disinfectant, 91
 Disinfectant, Wipes & Cleaners, 63–66
 Micro-Kill AF² Wipes, 65
 Micro-Kill AF+ Wipes, 66
 Micro-Kill Bleach Solutions, 64
 Micro-Kill Bleach Wipes, 63
 Micro-Kill One, 66
 Dispensers, 47–49, 70, 77
 Automated Touchless Towel, 47, 48
 C-Fold/Multi-Fold Towel, 47
 Glove, 70
 Manual, 47
 Sanitizer, 77
 Soap, 77
 Toilet Paper, 47, 49
 Disposable, 32, 33, 35, 50–56, 81–83, 87, 89–90, 98, 106, 143, 196, 200–205
 Biopsy Punches, 81
 Blankets, 51

Disposable (continued)
 Bouffant Caps, 53
 Bowls, 203
 Buckets, 205
 Cover Gowns, 53
 Cubicle Curtains, 196
 Cup Lids, 202
 Cups, 200, 201
 Cutlery, 204
 Dermal Currettes, 81
 Drapes, 106
 Draw Sheets, 52
 Ear Currettes, 81, 87
 Equipment Covers, 52
 Eyewear, 54
 Face Mask, 55
 Face Shields, 54
 Flat Sheets, 52
 Forceps, 83
 Isolation Gowns, 53
 Lab Goggles, 54
 Patient Capes, 56
 Patient Exam Shorts, 56
 Patient Gowns, 56
 Patient Robes, 56
 Pen Lights, 32
 Pillows, 50
 Plates, 203
 Portion/Souffle Cups, 201
 Probe Covers, Thermometer, 33
 Professional Towels, 51
 Safety Glasses, 54
 Scalpels, 82, 89
 Serveware, 200–202
 Shave Biopsy, 90
 Shoe Covers, 54
 Skin Stapler, 90
 Speculas, 35
 Sterile Drapes, 106
 Sterile, Fenestrated, Drapes, 106
 Stretcher Sheets, 52
 Surgeon's Cap, 54
 Underpads, 106
 Utensils, 204
 Vaginal Speculums, 98
 Drain Bags, 104
 Urinary, 104
 Drain Sponges, 156
 Drapes, 106, 108
 Procedural, 108
 Sterile, Disposable, 106
 Sterile, Fenestrated, Disposable, 106
 Dressing Jars, 97
 Dressing Retention Tape, 171
 Dressings, 159, 166–172
 Adhesive Island, 172
 Alginates, 169
 Antimicrobial Silver, 169, 170
 Burn and Wound, 166
 Composite Dressings, 172
 Emulsion, 159
 Foam Dressings, 166–168
 Gelling Fibers, 169
 Hydrocolloids, 170
 Petrolatum, 159
 Post-Op Dressings, 168

Dressings (continued)
 Super Absorbent Polymer, 168
 TheraHoney, 171
 Transparent Film, 172
 Waterproof Composite Dressings, 172
 Wound Dressings, 172
 Drug Screen Testing, 15
 Drug Testing, 15
 Alere Toxicology Drug Screening
 Devices, 15
 Dry Erase, 207, 208
 Board, 207
 Eraser, 207
 Marker, 208
 Duracell, 67
 Dura Hooks, 88
 Frazier, 88
 Duro, 205

E
 Ear Currettes, 81, 88
 Ear Currettes, Disposable, 87
 Ear Syringes, 87, 112
 Ear Ulcer Syringes, 112
 ECG, 42, 43
 12 Channel Interpretive, 42
 Combination Event Recorder, 42
 Holter Software Kit, 42
 PC Based, 43
 Universal Diagnostic, 42
 ECG Adapter Clips, 40
 Edan, 39, 43
 EKG Caliper, 82
 Elastic Bandages, 160
 Electrodes, 40
 Electrosurgery Accessories, 184
 Elevators, 80
 Freer, 80
 Periosteal, 80
 Emesis Bags, 165
 Emesis Basins, 112
 Endometrial Cell Sample, 98
 Endometrial Sampling Set, 87
 Endospeculum, 84
 Enthermics, 174
 Envelopes, 208
 Epidural Needles, 111
 Epidural Trays, 110
 Ethicon, 108
 Evacuator, 184
 Smoke, 184
 Exam, 185, 188–194
 LED, 189
 Lights, 188, 189
 Patient Lift, 194
 Stools, 185
 Tables, 190, 191–193
 Exam Gloves, 70, 71
 Generation Pink, 71
 Generation Pink 3G, 71
 Generation Pink Pearl, 71
 Generation Pink Sense, 71
 Exam Table Paper, 50
 Excavators, 82

Extractor, 81
 Extractors, 83
 Eyewash, 68
 Eyewash Stations, 68
 Eyewear, Disposable, 54
 E-Z Lubricating Jelly, 95

F
 Face Mask, Disposable, 55
 Face Shields, Disposable, 54
 Facial Tissue, 46, 49
 Facility Supplies, 200–206
 Coffee Creamers, 206
 Coffee Maker, 206
 Coffee Pods, 206
 Disposable, 204, 205
 K-Cups, 206
 Keurig Brewer, 206
 Paper Bags, 205
 Plastic Bags, 205
 Serveware, 200–205
 Stirrers, 206
 Sugar, 206
 Zip Bags, 205
 Feminine Hygiene, 106, 107
 Maternity Pads, 106, 107
 Maxi Pads, 106
 Sanitary Pads, 106
 Tampons, 106
 Fiber Casting Tape, 164
 File Folders, 208
 Finger Splints, 112
 Frog, 112
 Finger Splints and Accessories, 112
 Stackies Kits, 112
 First Aid, 97, 98, 154
 Antibiotic Ointment, 154
 Antimicrobial Wound Gel, 154
 Bacitracin, 154
 Cold Packs, 97
 CURAD, 154
 Hot/Cold Packs, 97
 Hot Packs, 98
 Ointments and Creams, 154
 Petroleum Jelly, 154
 Vitamin A & D, 154
 Wart Pads, 154
 First Aid Only, 68
 FitGuard Touch Exam Gloves, 70
 Flow Meters, 44
 Fluid Clean-Up Kit, 68
 Foam Dressings, 166–168
 Foot & Ankle Supports, 116
 Forceps, 80, 81–84, 87–88
 Adson, 80
 Disposable, 83
 Dressing, Serrated, 80
 Platypus Nail-Pulling, 81
 Single-Use, 83
 Splinter, Serrated, 80
 Tissue, Standard, 80
 Forearm Splints, 112
 Formalin Containers, 24

Furnishings, 195
 Chairs, 195
 Lobby and Reception, 195

G

Gauze, 155, 156
 Sponges, 155, 156
 Dental, 155, 156
 Non-Sterile, 155, 156
 Sterile, 155, 156
 Woven, 155, 156
 Gauze Bandages, 158
 100% Cotton, 158
 Conforming Bandages, 158
 Sponges, 158
 Gauze Dressing, 158
 100% Cotton Bandages, 158
 Bulkee II Sponges, 158
 Conforming Bandages, 158
 Gauze Sponges, 156
 Drain Sponges, 156
 Gel, 96, 118, 154
 Antimicrobial Wound, 154
 Bunion Guard, 118
 CURAD, 154
 Toe Cap, 118
 Toe Crest, 118
 Toe Pads, 118
 Toe Separators, 118
 Toe Spacers, 118
 Toe Straightener, 118
 Toe Tubes, 118
 Ultrasound, 96
 Gelling Fibers, 169
 Generators, Electrosurgery, 184
 Genpak, 203, 204
 Georgia Pacific, 47, 48–49
 Glove Dispensers, 70
 Gloves, 70–72
 AloeTouch, 70
 AloeTouch Ice, 70
 Aloe vera, 70–72
 CURAD, 71
 DermAssure, 72
 Exam, 70, 71
 FitGuard Touch, 70
 Generation Pink, 71
 Generation Pink 3G, 71
 Generation Pink Pearl, 71
 Generation Pink Sense, 71
 Latex, 71, 72
 Latex-free, 70
 Nitrile, 70, 71
 Nitrile Exam, 70
 Powder-free, 70–72
 Restore, 70
 Restore Sense, 70
 SensiCare, 70, 72
 SensiCare Ice, 70
 SmartGuard, 70
 Sterile, 72
 Stretch Vinyl, 71
 Surgical, 72
 Synthetic, 72

Gloves (continued)
 Triumph LT, 72
 Vinyl, 70, 71
 Glucose Meter, 17
 Gowns, 53, 56–58
 Disposable Cover Gowns, 53
 Disposable Isolation Gowns, 53
 IV, 57
 Nursing, 57
 Patient Gowns, 56
 Unisex Teen, 57
 Women's Front Open, 58
 Women's Wellness Gowns, 58
 Greiner bio-one, 148
 Growth Charts, 181
 Boys and Girls, 181

H

Hammers, 80
 Percussion, 80
 Reflex, 82
 Hamper Stand, 187
 Stainless Steel, 187
 X-Frame, 187
 Hamper Stands, 186
 Handheld Oximeters, 101
 Hand Hygiene, 74–77
 Aerosol Foam, 75
 Antibacterial Soap, 76
 Foam, 74
 Foaming Soap, 77
 Gel, 74
 Hand Soap, 76, 77
 Lotion Soap, 76, 77
 Sanitizers, 74, 75
 Spectrum, 74, 75
 Sterillium, 74
 Sterillium Med Liquid, 75
 Sterillium Surgical Scrub, 75
 Hand Soap, 76, 77
 Antibacterial, 76
 Foaming, 76, 77
 HealthGuard, 76, 77
 Lotion, 76, 77
 Skintegrity, 76, 77
 Spectrum, 76
 Spectrum Antibacterial, 76
 Spectrum Foaming, 77
 hCG Pregnancy Testing, 6, 9
 Health o meter, 179, 180, 182
 Heart Rate Watch, 29
 Helena, 10
 HemoCue, 17, 18
 Hemoglobin Sensors, 38
 Hemoglobin Testers, 38
 Hemoglobin Testing, 17–19
 Hemostat, 80
 Crile, 80
 Dressing/Tissue, 80
 Halsted, Mosquito, 80
 Kelly, 80
 Swiss, Cilia, 80
 HemoSure, 10, 15
 Highlighter, 208

HIV/Drug Tests, 15
 Honey Wound Gel, 171
 Hot/Cold Therapy, 97, 98
 Hot Packs, 97, 98
H. Pylori Test, 6, 14
 H-R Pharmaceuticals, 95
 Huhtamaki, 205
 Hydrocolloids, 170
 Hydrogel, 170
 Silver, 170
 Hydrogen Peroxide, 94

I

iFOB Test Kits, 7, 10
 Impregnated Gauze, 159
 CURAD, 159
 Emulsion Dressing, 159
 Packing Strips, 159
 Petrolatum Dressing, 159
 Incision and Drainage Trays, 90
 Influenza Testing, 11, 12
 Infusion Sets, 144, 146
 Butterfly, 146
 Push-Button Safety, 146
 Safety Needle, 146
 Small Vein, 146
 Instrument Brush, 82
 Instruments, 80, 81–89, 91–92
 Bandage Scissors, 80
 Baron Suction Tubes, 87
 Billeau Ear Loops, 87
 Biopsy Punches, 81, 83–84
 Blade Glove Cartridges, 89
 Cannulas, Endo Biopsy, 87
 Cleaning Brushes, 92
 Curettes, 87
 Curettes-Excavators, 82
 Cuticle Nippers, 81
 Dermal Curettes, 81, 83
 Detergent, 91
 Disinfectant, 91
 Disposable Forceps, 83
 Disposable Scalpels, 82
 Dura Hooks, 88
 Ear Curettes, 81, 88
 Ear Curettes, Disposable, 87
 Ear Syringe, 87
 EKG Caliper, 82
 Elevators, 80
 Endometrial Sampling Set, 87
 Endospeculum, 84
 Excavators, 82
 Extractor, 81
 Extractors, 83
 Forceps, 80, 82–84, 87–88
 Hemostat, 80
 Instrument Brush, 82
 Iris, Knapp Scissors, 80
 Knife Handle, 80, 82
 Mayo Scissors, 80
 Nail Nippers, 81, 82
 Nail Splitters, 81
 Nasal Tenacula Hooks, 88
 Needle Holders, 80, 82–83
 Nucleus Knives, 82

Instruments (continued)

- Operating Scissors, 80
- Percussion Hammer, 80
- Pessaries, 84–86
- Phlebectomy Extractors, Varady, 89
- Phlebectomy Hooks, 88
- Pin-Strippers, Oesch, 89
- Platypus Nail-Pulling Forceps, 81
- Presoaks, 91
- Reflex Hammer, 82
- Scalpel Knife Handle, 89
- Scalpels, 82
- Scalpels, Disposable, 89
- Scissors, 82, 83
- Single-Use Forceps, 83
- Skin Hooks, Kleinert-Kutz, 89
- Speculums, 87
- Stain Remover Wipes, 92
- Stitch, Littauer Scissors, 80
- Stitch, Spencer Scissors, 80
- Suction Tubes, 88
- Surgical Blades, 83, 89
- Tuning Forks, 87, 88
- Uterine Sounds, 84, 87
- Vaginal Specula, 84
- Wipes, 92
- Integra, 82, 83, 87
- Iodine, 95
- Isopropyl Rubbing Alcohol, 94
- IV Catheters, 147
 - Auto Safety, 147
 - Shielded, 147
 - Slide Safety, 147
 - Standard, 147
- IV Extension Sets, 149
- IV Fluid Warmers, 174
- IV Poles, 190
 - Aluminum, 190
- IV Solutions, 150
 - Sterile, 150
- IV Stands, 190
 - Deluxe, 190
 - Heavy Duty, 190
- IV Supplies, 149
 - Accessories, 149
 - Administration Sets, 149

J

- Jars, 97, 194
 - Clear, 97
 - Dressing, 97
 - Sundry, 97, 194
- Johnson & Johnson, 62

K

- KAWASUMI, 145, 146
- K-Cups, 206
- Keurig, 206
 - Pods, 206
- Knife Handle, 82
- KUTOL, 74–77

L

- Lab Chemicals, 26
 - Acetone, 26
 - Alcohol, Isopropyl, 26
 - Alcohol, Reagent Grade, 26
 - Deionized Water, 26
 - Formalin, 26
 - Methanol, 26
 - Wright Giemsa Stain, 26
 - Xylene, 26
- Lab Coats, 59, 60
 - Men's Staff-Length, 60
 - Unisex Consultation, 60
 - Unisex Knee-Length, 59
 - Women's Consultation, 60
- Labeling Marker, 26
- Lab Goggles, Disposable, 54
- Laceration Trays, 91
- Lancets, 22, 151
 - Adjustable-Depth, 151
 - Contact Activated, 151
 - Disposable, 22, 151
 - Heel, 22, 151
 - Premie, 151
 - Pressure Activated, 151
 - Push Button, 151
 - Safety, 22, 151
 - Side Button, 151
- Latex Exam Gloves, 71
- Latex Surgical Gloves, 72
- Leg Bags, 103
- Leg Walkers, 117
 - Short, 117
- Lesro, 195
- LIFEPAK, 41
- LifeSign, 11
- Lights, 188, 189
 - Coloscopes, 189
 - Exam, 188, 189
 - Floorstand, 188, 189
 - Gooseneck, 188
 - Handheld, 189
 - Incandescent, 188
 - LED, 188, 189
 - Procedure, 188
 - Swing Arm, 189
 - Swivel Reflector, 188
 - Ultraviolet, 189
- Linens, 51, 52
 - Blankets, Disposable, 51
 - Drape Sheets, 52
 - Draw Sheets, Disposable, 52
 - Flat Sheets, Disposable, 52
 - Pillow Cases, 51
 - Professional Towels, Disposable, 51
 - Sheets, 51
 - Stretcher Sheets, Disposable, 52
 - Stretcher Sheets, Reusable, 52
- Liquids & Wet Pack Products, 94, 95–96
- Lobby and Reception, 195
 - Chairs, 195
- Loss of Resistance Syringe, 108
- Lubricating Jelly, 95
 - E-Z Jelly, 95
 - Surgilube, 95

M

- Mabis, 97
- Marker, Labeling, 26
- Masimo, 38
- Masks, 99
 - Aerosol, 99
 - O₂, 99
- Mat, 181
 - Babies Measuring, 181
- Maternity Pads, 107
- Mayo Stands, 187
 - Knob-Operated, 187
 - Stainless Steel, 187
- Measuring Mat, 181
 - Babies, 181
- Measuring Rod, 182
 - Wall-Mount, 182
- Measuring Tape, 181
 - Circumference, 181
 - Mechanical, 181
- Medical Carts, 174, 175
 - Medical, 175
 - Pediatric, 174
- Medical Spoons, 97
- Medigrip, 170
- Medtronic, 143–145
- meridian, 16
- Metacarpal Splints, 112
- Metrex, 62
- Micro-Kill Bleach Solutions, 64
- Micro-Kill Germicidal Wipes, 66
- Micro-Kill Wipes, 63, 65–66
- Microscope, 26
 - Cover Glass, 26
- Microscope Slides, 26
 - Frosted, 26
 - Plain, 26
- Mid Central Medical, 187
- Midmark, 44, 183, 193
- Monitor, 43
 - ECG, 43
- Monitors, 28, 29
 - Blood Pressure, 28
 - Blood Pressure, Wrist, 29
- Monofilaments, 22
- Mononucleosis Testing, 6, 14
- MSD, 152

N

- Nail Brushes, 95
 - Surgical, 95
- Nail Nippers, 81, 82
- Nail Polish Remover, 94
- Nail Splitters, 81
- Narcotics Cabinets, 174
- Nasal Tenacula Hooks, 88
 - Joseph, 88
- Nebulizers, 100
 - Compact, 100
 - Compressor, 100
 - Fish, 100
 - Handheld, 100
 - Kits, 100
 - MDI Spacers, 100

Nebulizers (continued)
 Peak Flow Meters, 100
 Puppy, 100
 Needle Holders, 80, 82–83
 Halsey, Serrated, 80
 Halsey, Smooth, 80
 Halsted-Mosquito, 80
 Mayo-Hegar, 80
 Olsen Hegar, 80
 Splinter, Serrated, 80
 Webster, Smooth, 80
 Needles, 111, 143–144
 Blood Collection, 143, 144
 Epidural, 111
 Spinal, 111
 Needles/Syringes, 140–144
 Allergist Tray w/Needle, 142
 Allergy Syringe, 141
 Allergy Syringe Trays, 141
 Beveled Needles, 141
 Blood Collection Needles, 143, 144
 Hypodermic, 140
 Luer Lock, 140, 142
 Luer Lock Adapter, 141
 PrecisionGlide Needles, 142
 Prefilled Saline Flush Syringes, 141
 Retracting Syringes, 142, 143
 Safety Needles, 140
 Safety Syringes, 140, 142
 Standard Syringes, 141
 Tuberculin Syringes, 141, 143
 Nipro, 143, 146–147
 Nitrile Exam Gloves, 70, 71
 Non-Adherent Pads, 157
 CURAD, 157
 Non-Woven Gauze Sponges, 156
 Non-Sterile, 156
 Sterile, 156
 Notebooks, 208
 Nucleus Knives, 82

O

Obstetrical Towelettes, 94
 Office Supplies, 200, 201–208
 Binders, 207
 Clips & Rubberbands, 207
 Coffee Creamers, 206
 Coffee Maker, 206
 Coffee Pods, 206
 Copy Paper, 208
 Dry Erase, 207
 Envelopes, 208
 File Folders, 208
 K-Cups, 206
 Keurig Brewer, 206
 Notebooks & Writing Pads, 208
 Paper Bags, 205
 Paper Products, 208
 Plastic Bags, 205
 Scissors, 207
 Self-Stick Note Pads, 208
 Serveware, 200, 201–205
 Stapler & Staples, 207
 Stirrers, 206
 Sugar, 206

Office Supplies (continued)
 Tape, 207
 Writing Supplies, 208
 Zip Closure Bags, 205
 Ointments, 154
 Antimicrobial Wound Gel, 154
 CURAD, 154
 Petroleum Jelly, 154
 Triple Antibiotic, 154
 Vitamin A & D, 154
 Ophthalmoscopes, 35, 37
 OraSure Technologies, 111
 Otoscopes, 34, 35, 37–38
 Oximeters, 101
 Easy-Grip, 101
 Handheld, 101
 High-Impact, 101
 Pediatric, 101
 Oxygen Delivery, 99
 Oxygen Therapy, 99
 Cannulas, 99
 Masks, 99
 Oxygen Tubing, 99
 Resuscitation, Manual, 99

P

Packing Strips, 159
 Packs, 97, 98
 Cold, 97
 Hot/Cold, 97, 98
 Pain Management, 108–111
 Block Trays, 109
 Epidural Needles, 111
 Single Shot Epidural Trays, 110
 Spinal Needles, 111
 Paper, 50
 Exam Table, 50
 Paper Products, 208
 Copy Paper, 208
 Envelopes, 208
 File Folders, 208
 Notebooks & Writing Pads, 208
 Self-Stick Note Pads, 208
 Paper Products and Supplies, 46, 47–49
 Automated Dispenser, 47
 Automated Touchless Dispenser, 48
 Dispenser, Toilet Paper, 47, 49
 Facial Tissue, 46, 49
 Manual Dispenser, 47
 Toilet Tissue, 46, 49
 Toilet Tissue, Coreless, 49
 Towels, C-Fold, 47, 48
 Towels, Multi-Fold, 47
 Towels, Roll, 46, 47–48
 Pap Smear Tray, 98
 Parker Labs, 40
 Patient Apparel, 56–58
 Disposable, 56
 Capes, 56
 Exam Shorts, 56
 Gowns, 56
 Robes, 56
 IV Gowns, 57
 Nursing Gowns, 57
 Unisex Teen Gowns, 57

Patient Apparel (continued)
 Women's Front Open Gowns, 58
 Women's Wellness Gowns, 58
 Patient Care, 95
 Brushes, 95
 Scrub Brush, 95
 Surgical Scrub Brush, 95
 Patient Lift, 194
 Peak Flow Meters, 100
 Pediatric Carts, 174
 Pedometer, 29
 Pencil, 208
 Pen Lights, 32
 Disposable, 32
 Pens, 208
 Peroxide, Hydrogen, 94
 Personal Care, 106, 107
 Pessaries, 84–86
 Accessories, 86
 Cube, 84
 Cup, 86
 Dish, 85
 Donut, 84
 Gehrung, 85
 Gellhorn, 85
 Hodge, 85
 Knobbed Ring, 86
 Marland, 86
 Oval, 85
 Ring, 86
 Shaatz, 85
 Petri Dishes, 25
 Culture, 25
 Petroleum Jelly, 95, 154
 E-Z Jelly, 95
 Surgilube, 95
 Pharmaceuticals, 120–138
 Contrast, 120–122
 Inhalation, 122, 123
 Injectables, 123–130
 Miscellaneous, 131
 Ophthalmics, 131, 132
 Oral Tablets, Capsules & Liquids, 132, 133–134
 Otics, 135
 RX/OTC Topicals & Ointments, 135, 136
 Vaccines, 136–138
 Philips, 40
 Phlebotomy Extractors, 89
 Varady, 89
 Phlebotomy Hooks, 88
 Crochet, 88
 Mueller-Style, 88
 Oesch-Style, 88
 Ramelet, 88
 Phlebotomy Chairs, 195
 Pillows, 50
 Disposable, 50
 Pin-Strippers, Oesch, 89
 Pipettes, 25
 Graduated, 25
 Non-Graduated, 25
 Transfer, 25
 Plates, 203
 Disposable, 203
 Podiatry Products, 118

Post-Op Dressings, 168
 Post-Op Shoes, 116
 Povidone Iodine, 95
 Prep Pads, 95
 Prep Solution, 95
 Scrub Solution, 95
 Swabsticks, 95
 Premier, 111
 Prep Pads, 94
 Alcohol, 94
 Prep Solution, 95
 Povidone Iodine, 95
 Prep Pads, 95
 Swabsticks, 95
 Procedural Drapes, 108
 Protective Apparel, Disposable, 53–55
 Bouffant Caps, 53
 Boundary Bouffant Caps, 53
 Cover Gowns, 53
 Eyewear, 54
 Face Mask, 55
 Face Shields, 54
 Isolation Gowns, 53
 Lab Goggles, 54
 Safety Glasses, 54
 Shoe Covers, 54
 Surgeon's Cap, 54
 PT/INR System, 19
 Hemoglobin, 19
 pts Diagnostics, 16, 18
 Pulse Oximeters, 38, 101

Q

QRS, 42, 44
 Quantimetrix, 8
 Quantum Storage Systems, 174
 Quidel, 9, 11–12, 14

R

Rechargeable Handles, 37
 Reflex Hammers, 82
 Restore Exam Gloves, 70
 Resuscitation, Manual, 99
 Oxygen Therapy, 99
 Rice Lake, 179–181
 Roche, 8, 12, 19
 Room Signs, 194
 Rubberbands, 207
 Rubbing Alcohol, 94

S

Safety Glasses, Disposable, 54
 SafetyGlide Needles, 142
 Safety Lancets, 22, 151
 Adjustable-Depth, 151
 Contact Activated, 151
 Disposable, 151
 Heel, 151
 Premie, 151
 Pressure Activated, 151
 Push Button, 151
 Safety Lancets (continued)
 Side Button, 151
 Safety Needles, 146
 Safety Needles/Syringes, 140
 Saline Wound Wash, 165
 CURAD, 165
 Sanitary Pads, 106
 Maternity Pads, 106
 Maxi Pads, 106
 Sanitizers, Hand, 74, 75, 77
 Aerosol Foam, 75
 Alcohol-free, 75
 Dispensers, 77
 Foaming, 74, 75
 Gel, 74
 HealthGuard, 74, 75
 Instant, 74
 KUTOL, 74, 75
 Spectrum, 74, 75
 Sterillium, 74
 Sterillium Med Liquid, 75
 Sterillium Surgical Scrub, 75
 Scales, 179–181
 Baby, 181
 Diaper, 181
 Dietary, 180
 Digital, 179–181
 Digital Platform, 179
 Digital Wheelchair, 179
 Eye- and Waist-Level Scales, 180
 Floor Level, 180
 iPC Digital Dietary, 180
 Mechanical, 179
 Physician, 179, 180
 Platform, 179
 Step-On, 180
 Wheelchair, 179
 Scalpel Knife Handle, 89
 Scalpels, 82
 Scissors, 80, 82–83, 207
 Bandage, 80
 Curved, 80
 Iris, 80
 Mayo, 80
 Operating, 80
 Stitch, 80
 Straight, 80
 Scrub Brush, 95
 Scrubs, 61
 Unisex Jacket, 61
 Unisex Pants, 61
 Unisex Top, 61
 Women's Pants, 61
 Women's Top, 61
 Scrub Solution, 95
 Povidone Iodine, 95
 sdi Diagnostics, 44
 seca, 181, 182
 Sekisui, 9, 11, 14–15
 Self-Adherent Bandages, 161
 Co-Flex, 161
 Wrap, 161
 Self-Adherent Wrap, 161
 SensiCare Ice Exam Gloves, 70

Serveware, 200, 201–205
 Bowls & Plates, 203
 Buckets, 205
 Cup Lids, 202
 Cups, 200, 201
 Cutlery, 204
 Disposable, 203
 Straws, 204
 Utensils, 204
 Sharps Containers, 152
 Biohazard, 152
 Mail Back, 152
 Multipurpose, 152
 Patient Room, 152
 PG-II, 152
 Phlebotomy, 152
 Shave Biopsy, 90
 Sheets, 51, 52
 Contour, 51
 Drape Sheets, 52
 Draw, 51
 Draw, Disposable, 52
 Flat, Disposable, 52
 Full, 51
 Pillow Cases, 51
 Stretcher, Disposable, 52
 Stretcher Sheets, Reusable, 52
 Shoe Covers, Disposable, 54
 Shoes, 116
 Post-Op, 116
 Short & Long Leg Walkers, 117
 Walkers, 117
 Siemens, 8, 17
 Signs, 194
 Exit, 194
 Room, 194
 Silicone Contact Layer, 171
 Silver Dressings, 169, 170
 Argales Powder, 170
 Gelling Fibers, 169
 Hydrogel, 170
 Skin Hooks, 89
 Kleinert-Kutz, 89
 Skin Marker, 90
 Sterile, 90
 Surgical, 90
 Skin Protectants, 171
 Skin Stapler, 90
 Slides, Microscope, 26
 Frosted, 26
 Plain, 26
 Slings, 115
 Deep Pocket, Arm, 115
 Standard, Arm, 115
 Smoke Evacuator, 184
 Soap, 76, 77
 Antibacterial, 76
 Automatic Dispensers, 77
 Dispensers, 77
 Foaming, 76, 77
 HealthGuard, 76, 77
 Lotion, 76, 77
 Skintegrity, 76, 77
 Spectrum, 76, 77

- Soap (continued)
 Spectrum Antibacterial, 76
 Solution Stands, 187
 Stainless Steel, 187
 Specimen Bags, 24
 Specimen Collection, 105, 106
 Bottle, 105, 106
 Calculi Strainer, 105
 Containers, 105
 Infant, 105
 Speci-Pan, 105
 Stool, 105, 106
 Two Piece for Commode, 106
 Urine, 105, 106
 Specimen Containers, 25, 105–106
 Speculas, 35
 Disposable, 35
 Speculums, 87, 98
 Disposable, 98
 Vaginal, 98
 Sphygmomanometers, 30
 Spinal Needles, 111
 Spirometers, 44
 Digital, 44
 Splints, 112, 116
 Ankle, 116
 Dorsal Night, 116
 Forearm, 112
 Frog Finger Splints, 112
 Metacarpal, 112
 Stackies Kits, 112
 Wrist, 116
 Wrist and Forearm, 116
 Sponges, 155, 156, 158
 Bulkee II, 158
 C-Fold, 155
 Cotton Gauze, 155
 Dental, 155
 Drain, 156
 Non-Woven Gauze, 156
 Woven Gauze, 155
 Stack and Hang Bins, 174
 Stackies Kits, 112
 Stadiometer, 182
 Freestanding, 182
 Portable, 182
 Stainless Steel, 186, 187, 190
 Foot Stool, 186
 Hamper Stands, 186, 187
 IV Poles, 190
 IV Stands, 190
 Mayo Stands, 187
 Solution Stands, 187
 Stacking Foot Stool, 186
 Stalk Market, 201–204
 STANBIO, 17
 Stapler & Staples, 207
 STAT Bags, 24
 Steam Autoclave Tape, 92
 Step Stools, 186
 Sterile Solutions, 150
 Saline, 150
 Sodium Chloride, 150
 Water, 150
 Sterilization, 91, 92
 Detergent, 91
 Disinfectant, 91
 PassPort Plus Monitoring Service, 91
 Pouches, 92
 Presoaks, 91
 Stain Remover Wipes, 92
 Steam Autoclave Tape, 92
 Steam Indicator Strips, 91
 Wipes, 92
 Sterilizers, 182, 183
 EMS, 182
 Monitoring Services, 182
 PassPort and PassPortPlus, 182
 UltraClave Automatic, 183
 Valueklave 1730, 183
 Stethoscopes, 32
 Sticky Notes, 208
 Stirrers, Coffee, 206
 Stirrup, 116
 Air/Gel Ankle, 116
 Stockinettes, 113
 Stockings, 115
 Anti-Embolism, 115
 Stools, 185, 186
 Exam, 185
 Foot Stool, 186
 Foot Stool, Stacking, 186
 Height-Adjustable, 185
 Pneumatic, 185
 Step, 186
 Traditional, 185
 Strep A Testing, 6, 11
 Stretcher Sheets, 52
 Disposable, 52
 Reusable, 52
 Suction Tubes, 88
 Suction Tubes, Baron, 87
 Sugar, 206
 Sundry Jars, 97, 194
 Super Absorbent Dressings, 168
 Surgeon's Cap, Disposable, 54
 Surgery, 95
 Scrub Brush, 95
 Surgical Blades, 83, 89
 Carbon, 89
 Disposable, 89
 Safety, 89
 Scalpels, 89
 Stainless Steel, 89
 Surgical Gloves, 72
 Surgical Scrub, 75, 95
 Brushes, 95
 Fragrance Free, 75
 Povidone Iodine, 95
 Sterillium, 75
 Surgical Scrub Brush, 95
 Surgical Specialties, 107
 Surgilube Lubricating Jelly, 95
 Suture Removal Trays, 90
 Sutures, 107, 108
 Swabs, 96
 OB/GYN, 96
 Procto, 96
 Swabsticks, 94, 95
 Alcohol, 94
 Povidone Iodine, 95
 Swab Sticks, 96
 Cotton-Tipped, 96
 Synthetic Surgical Gloves, 72
 Syringes, 108, 112
 Ear Ulcer, 112
 Loss of Resistance, 108
 Syringes/Needles, 140–144
 Allergist Tray w/Needle, 142
 Allergy Syringe, 141
 Allergy Syringe Trays, 141
 Beveled Needles, 141
 Blood Collection Needles, 143, 144
 Hypodermic, 140
 Luer Lock, 140, 142
 Luer Lock Adapter, 141
 PrecisionGlide Needles, 142
 Prefilled Saline Flush Syringes, 141
 Retracting Syringes, 142, 143
 SafetyGlide Needles, 142
 Safety Needles, 140
 Safety Syringes, 140, 142
 Standard Syringes, 141
 Tuberculin Syringes, 141, 143
T
 Table Paper, Exam, 50
 Tables, 190, 191–194
 Exam, 190, 191–194
 Pediatric, 192
 Power, 190, 191
 Treatment, 194
 Tampax, 106
 Tape, 92, 207
 Steam Autoclave, 92
 Tape Remover, 94
 Teleflex, 44
 Test Cassettes, 6, 9, 14
 Alere hCG Pregnancy, 9
 Alere *H. Pylori* Test, 14
 Medline hCG Pregnancy, 6, 9
 Medline *H. Pylori* Test, 6, 14
 Medline Mononucleosis, 6, 14
 Quidel hCG Pregnancy, 9
 Sekisui hCG Pregnancy, 9
 Testing, 6, 9–12, 14–20
 Blood Analysis, 20
 Blood Chemistry Testing, 20
 Blood Glucose, 16–19
 BVBLUE, 15
 Cholesterol, 16
 Drug Screen, 15
 Fecal Occult Blood, 10
 hCG Pregnancy, 6, 9
 Hemoglobin, 17–19
 HIV/Drug Tests, 15
 Influenza, 11, 12
 Influenza A & B, 12
 i-STAT, 20
 Lead Analyzer, 16
 Medline Strep A, 11

Testing (continued)
 Mononucleosis, 6, 14
 RSV, 14
 Strep A, 6, 11
 Trichomonas, 15
 Test Strips, 11, 62
 Influenza A & B, 11
 Test Tubes, 25
 TheraHoney Dressings, 171
 Thermometers, 33, 34
 Digital, 33, 34
 Forehead, 33
 Probe Covers, 33
 Temple, 34
 Tissue, 46, 49
 Dispenser, 49
 Facial, 46, 49
 Toilet, 46, 49
 Toilet, Coreless, 49
 Toe, 118
 Caps, 118
 Gel, 118
 Ribbed Gel, 118
 Crest, Gel, 118
 Pads, Gel, 118
 Separators, Gel, 118
 Spacers, Gel, 118
 Straightener, Gel, 118
 Tube, Ribbed Gel, 118
 Toilet Paper, 46, 49
 Coreless, 49
 Tongue Depressors, 97
 Tourniquets, 143
 Disposable, 143
 Towelettes, 94
 Benzalkonium Chloride, 94
 Cleansing, 94
 Obstetrical, 94
 Wash-Up, 94
 Towels, 46, 47–48, 51
 Automated Dispenser, 47
 Automated Touchless Dispenser, 48
 Manual Dispenser, 47
 Paper, C-Fold, 47, 48
 Paper, Multi-Fold, 47
 Paper, Roll, 46, 47–48
 Professional, Disposable, 51
 Traditional Wound Care, 154, 155–165
 Absorbent Pads, 157
 Adhesive Bandages, 162, 163
 Adhesive Tapes, 163, 164
 Arm Boards, 165
 Elastic Bandages, 160
 Emesis Bags, 165
 Fiber Casting Tape, 164
 First Aid Ointments & Creams, 154
 Gauze Bandage Rolls, 158
 Gauze Sponges, 155, 156
 Impregnated Gauze, 159
 Non-Adherent Pads, 157
 Saline Wound Wash, 165
 Self-Adherent Bandages, 161
 Transparent Film Dressings, 172
 Trash Liners, 49, 50
 Biohazard, 49
 High Density, 50

Trash Liners (continued)
 Low Density, 49
 Trays, 90, 91, 98, 103–104, 109–110,
 141–142
 Allergist Tray w/Needle, 142
 Allergy Syringe, 141
 Block Trays, 109
 Foley Catheter Insertion Tray, 103
 Incision and Drainage, 90
 Intermittent Catheterization, 104
 Laceration, 91
 Pain Management, 109, 110
 Pap Smear Tray, 98
 Pre-Connected Intermittent Catheter
 Trays, 104
 Single Shot Epidural, 110
 Suture Removal, 90
 Treatment Tables, 194
 Triple Antibiotic Ointment, 154
 Tubes, 25
 Round Bottom Test, 25
 Tubing, Oxygen, 99
 Tuning Forks, 87, 88
 Tuttnauer, 183

U
 Ultrasound Gel, 96
 Ultrasound Probe Covers, 40
 UMF Medical, 186, 190, 191
 Undercast Padding, 113
 Underpads, 106
 Incontinence, Disposable, 106
 Universal, 207, 208
 Unna-Boots, 169
 Urinal Holders, 106
 Urinalysis, 7, 8
 Reagent Strips, 7, 8
 Testing, 7, 8
 Urinary Drain Bags, 104
 Urine, 7, 8
 Analyzer, 7, 8
 Dipstick Control, 8
 Urine Collection, 105, 106
 Bottle, 105, 106
 Calculi Strainer, 105
 Infant, 105
 Speci-Pan, 105
 Uterine Sounds, 84, 87

V
 Vacutainers, 148
 Low Volume, 148
 Vaginal Specula, 84
 Vaginal Speculums, 98
 Vinyl Exam Gloves, 70, 71
 Vital Signs Monitor, 39
 Vitamin A & D Ointment, 154

W
 Walkers, 117, 177
 Ankle, 117
 Bariatric, 177
 Folding, 177
 Short Leg, 117
 Short & Long Leg, 117
 Walkers & Post-Op Shoes, 116, 117
 Wallach, 98
 Warmers, 174
 Blanket, 174
 IV Fluid, 174
 Wart Pads, 154
 Wash Basins, 112
 Wash-Up Towelette, 94
 Watch, Heart Rate, 29
 Waterloo Healthcare, 174
 Welch Allyn, 33, 34, 37–38, 43, 98
 Wheelchairs, 176
 Excel 2000, 176
 Extra-Wide, 176
 Scales, 179
 Winco, 194
 WinCup, 200, 202, 203
 Wipes, 63, 65–66, 92, 107
 Alcohol, 66
 Alcohol-Free, 107
 Baby, 107
 Bleach, 63
 Fragrance-Free, 107
 Instrument, 92
 Micro-Kill, 63, 65–66
 Presoaks, 92
 Wire Glove Dispensers, 70
 Wooden Applicator Sticks, 97
 Wooden Spoons, 97
 Wound Cleansers, 172
 Wound Device, 166
 Hyalomatrix Hyaluronic Acid, 166
 Wound Dressings, 166
 Woven Gauze Sponges, 155
 Dental, 155
 Non-Sterile, 155
 Sterile, 155
 Wrist & Forearm Splints, 116
 Writing Pads, 208
 Writing Supplies, 208

Z
 Zoll, 40

With a broad and deep product line, immense inventory and distribution capacities and 24/7 online ordering, Medline can help you provide the best possible service for your Physician Office.

STRENGTH IN NUMBERS

Our more than 20 million square feet of warehouse space, stocked to deliver within 24 hours, puts us in a unique position to quickly deliver what you need, where and when you need it. Here's how we do it:

100+

YEARS OF EXPERTISE

\$11 BILLION
GLOBAL ENTERPRISE

1,600+

DEDICATED SALES REPRESENTATIVES SERVING THE CONTINUUM

DEDICATED PHYSICIAN OFFICE SALES REPRESENTATIVES

18

MANUFACTURING FACILITIES WORLDWIDE

SERVE CUSTOMERS IN MORE THAN **90 COUNTRIES**

18,000+

EMPLOYEES WORLDWIDE


750+

MEDTRANS VEHICLES

40+ DISTRIBUTION CENTERS

IN NORTH AMERICA, ALLOWING US TO SERVE 99 PERCENT OF OUR CUSTOMERS IN ONE DAY





Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

FOLLOW US

Medline United States
1-855-294-9618
medline.com/physician | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com