

Exclusive soft formulation forms to hands for a second skin feel

- Designed with ERGOFORM™ Design Technology proven to support musculoskeletal health
- Conforms to the shape of your hands through the soft nitrile formulation
- Increases productivity and performance compared to using bare hands
- Packaged as 300 gloves per box



Industries

- Dental
- Emergency Medical Services
- Healthcare
- Industrial Safety
- Laboratory
- Law Enforcement
- Manufacturing

Applications

- Analytical testing/ Measurements
- Assembly and inspection
- Biotechnologies
- Crime scene evaluation
- Emergency services
- Handling small and sensitive parts and materials
- Laboratory analysis and testing
- Light duty maintenance and clean-up
- Manufacturing & preparing pharmaceutical products
- Patient examination during transport on EMS vehicles

TECHNICAL DATA SHEET:

PRODUCT INFORMATION:

ULTRAFORM®	
Material	Nitrile
Color	Blue
Glove Design	Ambidextrous, Powder-Free, Chlorinated, Textured Fingers
Cuff	Beaded Cuff
Manufacturing/ QMS Audit Standards	ISO 13485
Regulatory/Standards Compliance	ASTM D6319, US QSR/FDA 510K MEDICAL (EXAMINATION) GRADE, COMPLIANT WITH U.S. FOOD CONTACT REGULATIONS
Packaging	300 gloves per dispenser, 10 dispensers per shipper case, 1,000 gloves per shipper case
Storage	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Country of Origin	Malaysia
User Needs Segment	High Touch Glove - Ultra-thin gloves that provide the highest level of tactile sensitivity, especially at the fingertip.

PHYSICAL PROPERTIES:

	TYPICAL VALUES					TESTING METHOD
	XS	S	M	L	XL	
Length (mm/in.)	230 / 9.1	230 / 9.1	230 / 9.1	230 / 9.1	230 / 9.1	ASTM D3767
Average Palm Width (mm/in.)	75 / 2.95	85 / 3.34	95 / 3.74	105 / 4.13	115 / 4.52	
Freedom from Holes (Inspection Level I)	1.5 AQL					ASTM D5151
Palm Thickness Single Wall	(mm : 0.06) / (mil : 2.4)					ASTM D3767
	BEFORE AGING		AFTER AGING			
Ultimate Tensile Strength	23 Mpa		14 Mpa		ASTM D 412-06a	
Elongation at Break (%)	< 500		> 400		ASTM D 412-06a	

ORDERING INFORMATION:

	SIZE	XS	S	M	L	XL
ULTRAFORM®	REORDER NO.	UF-524-XS	UF-524-S	UF-524-M	UF-524-L	UF-524-XL

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application

For additional information visit us at www.ansell.com, or call us at 1-877-750-6110

Ansell
111 Wood Avenue South
Suite 210
Iselin, NJ 08830
USA

Ansell
105 Lauder
Cowansville, QC J2K 2K8
Canada

Technology

