

MPM Wound & Skin Cleanser

Non-toxic – Contains Poloxamer 188



Used to:

- Remove particulate matter & other debris from dermal wounds without harming living tissue

Directions:

Spray wound & surrounding area. Allow time for debris to soften. Adjust nozzle to "STREAM" to loosen debris. Adjust nozzle to "SPRAY" to gently flush debris from wound bed. No rinsing necessary.

	Product #	Size	Pkg.
MPM Wound & Skin Cleanser	MP00030	16 oz. Spray Bottle	6/Cs.
MPM Wound & Skin Cleanser	MP00031	8 oz. Spray Bottle	9/Cs.

Specifications:

Ph Range 6.8 to 7.5

PSI 5.8 to 6.2 (Depends on distance)

Spray can be adjusted



16 oz

8 oz