

Bathing Cushion

The Bathing Cushion is used to help someone who has difficulty getting in and out of the bath



"Warm water reduces musculoskeletal pain."
Arthritis Foundation



Technical Specifications

Max user weight	150 kg / 24 stone
Max height	40 cm / 15 ¾ "
Min height	Flat
Seat depth	44 cm / 17 ½ "
Width	58 cm / 23 "
Total Product weight	2 kg / 4 ½ lbs
Back/headrest length	33 cm / 13 ¾ "
Airflo Compressor	Mk3

Product code: MPCA030200

Easy to fit

- Fitted and removed within seconds
- Some air can be left in the backrest for greater levels of comfort
- The design adapts to the size and shape of the bath, so no expensive bathroom alterations are required

Safe to use

- Only lowers a person into the bath if there is enough battery power to raise them back out again
- Battery operated – no electric cables
- Fully waterproof hand control (IP67)

Simple to maintain

- Lightweight and portable, it can be removed and cleaned with ease
- Rinse with clean water after each bath and leave it to dry inflated
- Compact and easy to store

Improves wellbeing

- Perfect for a person who has difficulty getting in and out of the bath
- The cushion deflates allowing the user to lie flat in the bath, fully covered in water so benefiting from the therapeutic properties of bathing
- The user can take the bathing cushion with them when travelling away
- Can be removed from the bath when not in use so enabling a more flexible environment for shared bathrooms

