

3 & 4-WHEEL SCOOTER

MAXIMA

HEAVY-DUTY
PERFORMANCE
AT AN
EXCEPTIONAL
VALUE

The **Maxima** Scooter features an upgraded electronics package to deliver heavy-duty performance at an exceptional value. Engineered with the bariatric market in mind, the Maxima features a weight capacity of 500 lbs., deep, comfortable seating and an ergonomically enhanced design.



Maxima 4-wheel



Maxima 3-wheel

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

www.pridemobility.com

SPECIFICATIONS

MAXIMA 3W: SC900 4W: SC940

Weight capacity:	500 lbs.
Maximum speed ¹ :	Up to 5.25 mph
Ground clearance ³ :	3.25"
Turning radius ³ :	3W: 49.25" 4W: 54"
Overall length ³ :	3W: 46" 4W: 46.9"
Overall width ³ :	3W: 25.5" 4W: 25.5"
Seat-to-deck height:	16" - 18"
Seat-to-floor height:	22.5" - 24.5"
Front tires:	3" x 10" solid
Rear tires:	4" x 10" solid
Suspension:	3W: Rear 4W: Front and rear
Range per charge ^{1,2} :	3W: Up to 8.7 miles with U-1 batteries 3W: Up to 15.4 miles with NF-22 batteries 4W: Up to 10.2 miles with U-1 batteries 4W: Up to 16 miles with NF-22 batteries
Total weight without batteries ^{3,7} :	3W: 186 lbs. 4W: 192 lbs.
Heaviest piece when disassembled ³ :	76 lbs. (rear section)
Standard seating:	Type: Deluxe Contoured Medium-Back Weight: 40 lbs. ³ Material: Black Vinyl Dimensions: 22" x 18"
Drive system:	Rear-wheel drive, sealed transaxle with 24-volt DC motor
Dual braking system:	Regenerative and electromechanical
Battery requirements ^{3,4,5} :	(2) U-1: 38 lbs. each; (2) NF-22 (included): 24.5 lbs. each
Battery charger:	5-amp, off-board

Color:



Candy Apple Red



Viper Blue

FEATURES

- 500 lbs. weight capacity
- Upgraded electronics package
- Black, flat-free, non-marking tires
- Upgraded front basket
- Ultra heavy-duty drivetrain with 110 amp controller
- Maximum speed up to 5.25 mph
- Redesigned seat for improved comfort
- Easily adjustable tiller
- Durable and stylish rear polyurethane bumper
- Convenient tie-down points
- User-friendly console
- NF-22 batteries included

OPTIONS

- Power elevating seat
- Forearm crutch holder
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover
- Saddlebag
- Quad cane holder

Pride® Lifts & Ramps

Compatible Lifts:

- Outlander
- Outlander DE
- Outlander XL
- Commander 400 (4-wheel only)
- Zeus

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

² Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ AGM or gel cell type required.

⁵ Battery weight may vary based on battery manufacturer.

⁶ U-1 batteries required with power elevating seat.

⁷ Includes standard seat weight.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.



PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

Beamsville, Ontario • 888-570-1113
www.pridemobility.com