

# SILIPOS®

## ORTHOPAEDIC

### QUICK REFERENCE GUIDE

- SYMBOL KEY**
- Black Material Available.
  - Retail Packaging Available.
  - Antibacterial Available.
  - Bulk Packaging Available.

#### SILIPOS MINERAL OIL GEL

VS

#### SILICONE GEL

- Moisturises and softens the skin by releasing mineral oil over time.
- Great for scar management.
- Reduces friction and shear.
- Gel moves with the skin, and does not pull on hair or fragile tissues.
- Does not bottom out like foams or felts.
- Washable and hypoallergenic.

- May dry out the skin without the use of additional creams.
- Often causes rashes and follicular skin disorders.
- Will not glide along skin without additional creams which feels hot due to friction.
- No benefit to scar, burn, and graft tissues, and may cause injury due to friction.
- Occlusive environment is dry unless creams are used.

#### FINGERS & TOES

- Digital caps, tubes, and pads ease pressure and friction.
- They moisturise the skin making them great for corns, calluses, and blisters.



##### DIGITAL CAP

- Fully coated, covers tip of finger or toe.
- Available in 4 sizes.
- Mesh = more stretch.
- Ribbed = more durable.



##### TOE TUBE

- Fully coated.



##### CORN PAD

- Half coated.
- Narrow & Medium available.



##### DIGITAL PAD

- Gel dot for targeted use.
- Mini, S/M & L/XL available.

#### BULK TUBES & STRIPS

- Tubes can be cut to desired size.
- Cost-effective alternative to pre-cut solutions.
- Medical grade mineral oil gel to provide all around digital protection.



##### GEL TUBING

- Fully lined with gel.
- Mesh = more stretch.
- Ribbed = more durable.



##### DIGI-STRIP

- Half coated.
- Narrow, medium, and wide available.
- Creates 20 3cm long corn pads.



##### DIGITAL PADS ON A STRIP

- Gel dot for targeted use.
- Mini - 17 pads/strip.
- S/M - 16 pads/strip.
- L/XL - 14 pads/strip.

#### METATARSAL PADS



##### GEL METATARSAL PAD

- Stays in place with toe loop.
- Absorbs shock, impact, and provides even distribution of pressure to the forefoot.



##### GEL FOOT COVER

- Fully lined with gel.
- Gel is covered in fabric.
- Ideal for activities that place pressure on toes and metatarsals.



##### GEL METATARSAL PAD with SPREADER

- Silicone construction with built-in toe spreader helps reduce discomfort between toes whilst securing product in place.
- Spreader helps to align the big toe and prevent friction.



##### UNIVERSAL STRAP

- Stitching can be cut to fit instep, bunions, elbow, hand or wrist.
- Available in 2 sizes with gel covered or uncovered.
- Reduces pain and discomfort in the forefoot.



##### BALL OF FOOT GEL CUSHION

- Stays in place with toe loop.
- Low profile design fits most shoes.
- Gel is covered in fabric.



##### SLIM GEL-FIT METATARSAL PAD

- Gel is covered in fabric.
- Helps to relieve shoe pressure, friction, and shear forces on the bunion.

#### TOE ALIGNMENT

- Spreaders and separators relieve pressure from toes rubbing together.
- Toe crests and gel splints keep toes straight and are good for hammer, claw, and mallet toes.



##### SEPARATORS

- Fits between small toes.
- Gently divides, cushions, and aligns toes.



##### SPREADERS

- Fits between 1st and 2nd toe.
- Used as an inter-digital wedge for bunion deformity overlapping toes, or toe drift.



##### TOE CREST

- Helps to relieve pressure and pain caused by hammer, claw, and mallet toes.
- Reinforced toe loop to ensure proper fit.



##### GEL SPLINTS

- Can be used for buddy wrapping.
- Helps to stabilise injured digits.
- Relieves rubbing and irritation.

#### BUNION PADS

- Bunion products provide protection from rubbing and irritation.



##### ALL GEL BUNION CUSHION

- Hallux Bunion = big toe
- Tailor's Bunion = small toe



##### GEL BUNION

- Medical grade mineral oil gel pad cushions, protects, and reduces pressure on affected area.
- Gel dot can be placed on foam shell for exact positioning.



##### SLIM GEL-FIT BUNION SLEEVE

- Gel covered by elastic sleeve.



##### BUNION CARE GEL SLEEVE

- Thick uncovered gel pad directly touches skin.



##### DELUXE GEL BUNION SLEEVE

- Gel ring sits around the bunion to disperse pressure.

# ORTHOPAEDIC QUICK REFERENCE GUIDE

- SYMBOL KEY**
- Black Material Available.
  - Retail Packaging Available.
  - Antibacterial Available.
  - Bulk Packaging Available.



Email Orders to:  
orders@algeos.com

T: (+44) 151 448 1228  
W: www.algeos.com

## HEEL CUPS & INSOLES

- Made with medical grade silicone for maximum durability.
- Gel will not release oil onto the skin or into shoe.
- Help to reduce shock loads on the feet, knees, and back by providing even pressure distribution throughout the orthotic.



### WONDERSPUR

- For heel spurs.
- Absorbs shock, cushions, and protects.
- Dual density, coloured area is soft.
- Thick, high-profile cup.



### WONDERCUP

- For the entire heel.
- Absorbs shock, cushions, and protects.
- Dual density, coloured area is soft.
- Thick, high-profile cup.



### THERASTEP HEEL CUP

- Thin, low-profile cup.
- Absorbs shock, cushions, and protects.



### SILIPOS ACTIVE HEEL CUSHION

- Thin, low-profile flat cushion.
- Tiny waves help distribute pressure.
- Helps reduce heel pain.



### SOFTZONE HEEL SPUR WITH REMOVABLE PLUG

- Firm cushion for heel spurs.
- Plug can be removed for pressure relief.
- Wear with plug for general heel shock absorption.



### SOFTZONE 3/4 LENGTH INSOLE

- Soft multi-density, supports heel and arch.
- Anatomically designed to support longitudinal and transverse arches, the deep heel cup ensures superior comfort and control.



### SOFTZONE FULL LENGTH INSOLE

- Soft multi-density, supports heel, arch, and metatarsals. Relief for heel, arch, and forefoot problems.
- Multiple vents enhance air circulation.

## SLEEVES

- Sleeves are made with a stretchy fabric and offer gentle compression to the affected area.



### CARPAL GEL SLEEVE

- Can be used to reduce scars after surgery.
- Ambidextrous sleeve can be worn on either hand.



### ACHILLES HEEL PAD

- Reduces pressure on Achilles.
- Reduces pain caused by footwear pressure on irregular heel bumps.



### MALLEOLAR SLEEVE

- Best for inside sports wear.
- Protects sides of ankle from pressure.



### BOOT BUMPER

- Best for inside sports wear.
- Protects top of foot and Achilles.
- Provides light compression to the foot and ankle.



### PLANTAR FASCIITIS ARCH SLEEVE

- Prevents foot from rolling over.
- Helps with plantar fasciitis, tendonitis, heel spurs, and heel neuromas.




### HALF-COATED BODY SLEEVE

- Gentle compression is great for skin grafts, burns, and scars.
- Available in 8cm, 10cm, and 13cm dia.
- Cut to desired length.


## SHEETING & PADDING

- Can be used to custom make and pad orthotics, braces, and prosthetics.
- Gel conforms to almost any area of the body.



### SOFT SHEAR SHEETING

- 3mm sheet of gel available lined or unlined with stretch fabric.
- Can be used as orthotic interface.



### PRESSURE RELIEF PADDING

- 2mm thick, 2 way stretch top cover.
- Adhesive backing.
- Can use to line braces, and splints.



### GEL-E ROLL

- Elastic bandage coated on one side with gel.
- Can be used for compression, edema, cold therapy, and burn care.



### GEL DOTS

- Can be used for spot pressure relief.
- 3mm thick with adhesive backing.



### BODY DISCS

- Multipurpose gel discs.
- 6cm or 10cm dia. disc can be cut to size.



### GEL SQUARES

- 10cm x 10cm squares come with or without adhesive backing.
- Gel Squares are medical grade mineral oil gel cushions designed to soften, moisturise, protect, and cushion the skin against shear and abrasive forces.

## DIABETIC GEL SOCK



- Approved for people with poor circulation, diabetes, edema, and neuropathy.
- 0.16cm gel pad on the insole reduces shear and abrasion.

## KITS DIGITAL CARE / POST BURN



- An assortment of popular medical grade mineral oil gel products in a portable kit.
- Digital Care Kit.
- Post Burn Kit.

