




RESPIRATORY CARE



**Respiratory Supplies
Across the Continuum of Care**



Innovating to inspire.

From neonates to adults and across the continuum of care, Medline is continuously innovating to help improve outcomes, enhance the patient experience, and increase caregiver and operational efficiencies.

Table Of Contents

Airway Management

Endotracheal Tubes, Subglottic	5
Endotracheal Tube Holders, mBrace	5
Cuff Pressure Monitor	6
Inner Cannulas and Phoneris/Aeroflex System	7
Trach Ties	8

Bronchial Hygiene

Volume Incentive Spirometers	10
------------------------------	----

Gas Delivery & Suction Devices

Oxygen & Air Flowmeters	12-13
Air & Oxygen Blenders	14
Medical Gas Hose Assemblies	15
O2 Regulators	16
O2 Conservers	17
Suction Regulators	18-19
Suction Aspirator	20

Humidification & Filtration

Heat & Moisture Exchangers	22
Heat & Moisture Exchanger Filters	23
TrachVox HME	24
Bacterial/Viral Filters	24
Ventilator Circuits, Unheated	24
Pulmonary Function Filters	25

Medication Delivery & Asthma Management

Medication Nebulizers, Small Volume	27
Aerosol Masks	28
Valved Tee Adapters	28
Nebulizer Compressors, Compact	29
Nebulizer Compressors, Characters	30
Nebulizer Compressors, Portable	31
MDI Spacers	32
Peak Flowmeters	32

Non-Invasive Interfaces

NIV Masks	34
NIV Circuits	34
CPAP Masks AirGel	35-36
CPAP Masks Silicone	37
CPAP Accessories	38

O2 & CO2 Disposables

Cannulas, O2	40
Cannulas, CO2	41
Tubing, O2	42
Tubing, Aerosol	42
Masks, O2	43
Masks, Aerosol	44
Bubble Humidifiers	45
O2 Accessories	45

Pulse Oximeters

Fingertip Pulse Oximeters	47
Handheld Oximeters	48
Pulse Oximeter Probes & Accessories	49

Resuscitation

Resuscitation Bags, Self Inflating	51-53
Mouth-to-Mask Resuscitator	54
Resuscitation Accessories	54
Hyperinflation Bags	55

Suction Catheters & Trach Care

Closed Suction	57-58
Open Suction Catheters	59
Open Suction Catheters, Sleeved	59
Open Suction Kits	60
Open Suction Kits with Solutions	61
Sterile Solutions	61
Trach Care & Cleaning Trays	62
Mucus/Specimen Traps	63

Top Respiratory Codes for Seasonal Flu Demand

64-65

Airway Management

Endotracheal Tubes, Subglottic	5
Endotracheal Tube Holders, mBrace	5
Cuff Pressure Monitor	6
Inner cannulas & Phoneris/Aeroflex System	7
Trach Ties	8

SubG™ Endotracheal Tube

Our SubG subglottic endotracheal tubes are designed to reduce the risk of nosocomial infections due to aspiration of oral and upper airway secretions.

- » Polyurethane cuff helps to prevent channel formation
- » High volume, low pressure barrel-shaped cuff for optimal sealing
- » Kink-resistant tubing with easy-to-read insertion markings
- » Subglottic suction port for removing secretions above the cuff
- » Suction line cap keeps contaminants out of lumen

Item No.	Description	Pkg
DYNJ18860	SubG ET Tube, Size 6.0	10/bx
DYNJ18865	SubG ET Tube, Size 6.5	10/bx
DYNJ18870	SubG ET Tube, Size 7.0	10/bx
DYNJ18875	SubG ET Tube, Size 7.5	10/bx
DYNJ18880	SubG ET Tube, Size 8.0	10/bx
DYNJ18885	SubG ET Tube, Size 8.5	10/bx



mBrace

The new mBrace from Medline provides the enhancements hospitals have been looking for in ET tube holder technology. Designed with patient comfort and caregiver efficiency in mind, the mBrace can help customers achieve their skin breakdown initiatives.

- » Clear silicone adhesive pads provide complete visibility to the skin and are gentle upon removal
- » Silicone pads can be peeled back, adjusted and reapplied, giving the caregiver access to clean sensitive skin without removing the device
- » Additional strap connection points allow caregiver flexibility for difficult patients and added security to help prevent lip ulcers
- » Sliding ET tube track allows the ET tube to be moved side to side in accordance with hospital protocols
- » ET tube strap securely locks the ET tube in place
- » Adjustable neck strap provides security
- » A 4 point, full-head strap available separately
- » Available with and without a bite block

Item No.	Description	Pkg
DYNJMBRACE	ET Tube Holder	Coming Soon! Call for Availability
DYNJMBRACEBB	ET Tube Holder with Bite Block	Coming Soon! Call for Availability
DYNJ4PTSTRAP	4 Point Harness	Coming Soon! Call for Availability



mBrace shown with 4 point harness.
*Sold separately

IntelliCuff® Cuff Pressure Controller

Continuous and automatic cuff management system.

- » IntelliCuff cuff pressure controller continuously measures and automatically maintains set cuff pressure of ET or tracheal tubes during mechanical ventilation
- » Secures airway management in intensive care units, operating rooms and air transport
- » Large display keeps important data visible at all times, while the convenient buttons make adjusting and verifying settings easy
- » Automatic regulation of cuff pressure helps keep pressure consistently within range, which could help reduce the risk of ventilator associated pneumonia (VAP) and tracheal injuries

Item No.	Description	Pkg
HCS951001	IntelliCuff Pressure Controller	1 ea
HCS282016	IntelliCuff Disposable Tubing	10/cs
HCS951012	IntelliCuff Power Cord	1 ea
HCS951019	IntelliCuff Mount, Click	1 ea
HCS951055	IntelliCuff Mount, Rail	1 ea
HCS951057	IntelliCuff Mount, Multi-angle	1 ea



Optimal cuff pressure with IntelliCuff



Over-inflated cuff




Under-inflated cuff

For more information
or to view an In-Service
Video to see how
Intellcuff works, visit
www.Intellcuff.com

HCS951001





Phoneris® Inner Cannula and Aeroflex® Ventilator Circuit Connector

The Phoneris/Aeroflex system uses an innovative design to comfortably and securely connect tracheostomy patients to ventilator circuits.

- » **Secure** – A tight connection helps avoid accidental disconnects
- » **Gentle** – Swivel enables comfortable movement and reduces torque on the stoma
- » **Easy** – Puts an end to improvised, unapproved securing devices

Item No.	Description	Pkg
HCSTIC4A	Phoneris Inner Cannula; Jackson 4	10/cs
HCSTIC6A	Phoneris Inner Cannula; Jackson 6	10/cs
HCSTIC8A	Phoneris Inner Cannula; Jackson 8	10/cs
HCSTIC10A	Phoneris Inner Cannula; Jackson 10	10/cs
HCSTIC4681A	Phoneris AeroFlex	50/cs
HCS3222	15 mm x 22 mm Flex Connector	50/cs
HCS3215	15 mm x 15 mm Flex Connector	50/cs





HCS541



HCS597S

Tracheostomy Ties

These soft, comfortable straps have adjustable hook-and-loop closures to keep tracheostomy tubes securely in place, ensuring uninterrupted airflow.

Premium Tracheostomy Ties with Technology from Pepper Medical

- » Soft, hypoallergenic, no stretch cotton ensures patient safety and comfort while reducing skin irritation
- » Offered in wide range of sizes

Item No.	Description	Pkg
HCS540	Tracheostomy Tube Holder, Adult, Up to a 20" neck, 2 piece	10/cs
HCS541	Tracheostomy Tube Holder, Up to an 18" neck, 2 piece	10/cs
HCS542	Tracheostomy Tube Holder, Up to an 9" neck, 2 piece	10/cs
HCS597S	Tracheostomy Tie, Small/Pediatric, 1 piece	10/cs
HCS597M	Tracheostomy Tie, Medium/Adult, 1 piece	10/cs
HCS597LA	Tracheostomy Tie, Large/Bariatric, 1 piece	10/cs

Standard Medline Tracheostomy Ties

- » Available in two sizes to fit a broad range of sizes

Item No.	Description	Pkg
ZPP501	Tracheostomy Tube Holder, 2-Piece Adult	20/cs
ZPP501P	Tracheostomy Tube Holder, 2-Piece Pediatric	20/cs
ZPP301	Tracheostomy Tube Holder, 1-Piece Adult	20/cs
ZPP301P	Tracheostomy Tube Holder, 1-Piece Pediatric	20/cs
ZPP401	Anti-Disconnect Tracheostomy Tube Holder	20/cs



Bronchial Hygiene

Volume Incentive Spirometers 10



HCS93250

Incentive Spirometers

Our Volume Incentive Spirometer uses an intuitive design to help patients correctly perform and monitor their post-surgical breathing exercises.

- » One way valve ensures patients inhale rather than exhale through the device
- » Easy to use Floating ball indicator and guidelines help patients inhale at the proper rate
- » Sliding goal indicator can be easily adjusted allowing patients to monitor their progress
- » Flexible tubing and dual sided calibrations allow use in various positions
- » Available in 2,500ml and 5,000ml capacities

Item No.	Description	Pkg
HCS93250	Incentive Spirometer, 2500 ml	12/cs
HCS93500	Incentive Spirometer, 5000 ml	12/cs



Gas Delivery & Suction Devices















Oxygen & Air Flowmeters	12-13
Air & Oxygen Blenders	14
Medical Gas Hose Assemblies	15
O ₂ Regulators	16
O ₂ Conservers	17
Suction Regulators	18-19
Suction Aspirator	20

Oxygen and Air Flowmeters

Antimicrobial-treated knobs*, environmentally friendly**, and exclusively distributed nationwide by Medline.

- » Antimicrobial and environmentally friendly additive in the rotary knobs*
- » For both oxygen and medical air
- » Thoroughly tested for accuracy
- » Durable polycarbonate hood and metering tube
- » Available in nickel plated solid brass or aluminum
- » Stainless steel ball float
- » Warranty: ten years (manufacturer)

Flowmeter Ordering Matrix

FLOWMETER TYPE ADAPTER						
	Air Flowmeter, 0-15 LPM	Air Flowmeter with Power Take-Off, 0-15 LPM	Air Dual Flowmeter, 0-15 LPM	Oxygen Flowmeter, 0-15 LPM	Oxygen with Power Take-Off, 0-15 LPM	Oxygen Dual Flowmeter, 0-15 LPM
 Standard Inlet 1/8" NPT Female	GTFM187A15L00	GTFM187A1500P	GFMD187A15L00	GTFM197A15L00	GTFM197A1500P	GFMD197A15L00
 Ohmeda® Adapter	GTFM187A15LOH	GTFM187A150HP	GFMD187A15LOH	GTFM197A15LOH	GTFM197A150HP	GFMD197A15LOH
 Chemetron® Adapter	GTFM187A15LCH	GTFM187A15CHP	GFMD187A15LCH	GTFM197A15LCH	GTFM197A15CHP	GFMD197A15LCH
 DISS Hand-Tight Nut and Nipple Adapter	GTFM187A15LDH	GTFM187A15DHP	GFMD187A15LDH	GTFM197A15LDH	GTFM197A15DHP	GFMD197A15LDH
 DISS Nut and Nipple Adapter	GTFM187A15LDS	GTFM187A15DSP	GFMD187A15LDS	GTFM197A15LDS	GTFM197A15DSP	GFMD197A15LDS
 Puritan-Bennett Adapter	GTFM187A15LPB	GTFM187A15PBP	GFMD187A15LPB	GTFM197A15LPB	GTFM197A15PBP	GFMD197A15LPB
 Oxequip® Medstar Adapter	GTFM187A15LOE	GTFM187A150EP	GFMD187A15LOE	GTFM197A15LOE	GTFM197A150EP	GFMD197A15LOE
 Oxequip Twist Adapter	GTFM187A15LOQ	GTFM187A150QP	GFMD187A15LOQ	GTFM197A15LOQ	GTFM197A150QP	GFMD197A15LOQ

Options

- » Flowmeters listed above are aluminum. For brass, change the A to a B (example: GTFM187A15LOH to GTFM187B15LOH)
- » Dual flowmeter with 2 power take-offs, also available

Specifications

Weight (No Adapter):	Aluminum Body: 4.24 oz/122 gm Brass Body: 6.64 oz/188 gm
Increment:	0-6 LPM in ½ LPM increment, 6-15 LPM in 1 LPM Increment
Flush Flow Rate:	65 LPM (depending on supply flow)
Max Inlet Pressure:	100 psig
PTO Flow Rate:	Over 100 LPM (depending on supply flow)

* This product is treated with an antimicrobial to protect the product itself.

**Agion's silver antimicrobial has been certified for its material content, recyclability, and manufacturing characteristics. EPA 71227-4-88165

Neonatal, Perinatal and Specialty Flowmeters

The performance and accuracy you need for critical care applications for neonate through adult patients.

- » 200 cc and 1 LPM models are designed for perinatal and neonatal applications where extreme accuracy under very low flow is required
- » Available for oxygen and air
- » Acrylic block provides a strong flow tube
- » Maximum inlet pressure: 100 PSI
- » Outlet: DISS male (9/16"-18UNF)
- » Calibrated at 50 PSI
- » Warranty: five years (manufacturer)



FLOWMETER ADAPTER		200 CC Oxygen	1 LPM Oxygen	3 LPM Oxygen	16 LPM Oxygen	75 Oxygen LPM	16 LPM Air	75 LPM Air
	1/8" NPT Female	GTFM293200CXX	GTFM29301XX	GTFM29303XX	GTFM29316LXX	GTFM29375LXX	GTFM28316LXX	GTFM28375LXX
	Chemetron® Adapter	GTFM293200CCH	GTFM29301CH	GTFM29303CH	GTFM29316LCH	GTFM29375LCH	GTFM28316LCH	GTFM28375LCH
	DISS Hex	GTFM293200CDS	GTFM29301DS	GTFM29303DS	GTFM29316LDS	GTFM29375LDS	GTFM28316LDS	GTFM28375LDS
	DISS Hand-Tight	GTFM293200CDH	GTFM29301DH	GTFM29303DH	GTFM29316LDH	GTFM29375LDH	GTFM28316LDH	GTFM28375LDH
	Ohmeda® Adapter	GTFM293200COH	GTFM29301OH	GTFM29303OH	GTFM29316LOH	GTFM29375LOH	GTFM28316LOH	GTFM28375LOH
	Puritan-Bennett Adapter	GTFM293200CPB	GTFM29301PB	GTFM29303PB	GTFM29316LPB	GTFM29375LPB	GTFM28316LPB	GTFM28375LPB
	Oxequip® Medstar Adapter	GTFM293200COE	GTFM29301OE	GTFM29303OE	GTFM29316LOE	GTFM29375LOE	GTFM28316LOE	GTFM28375LOE

Durable Housing

- » The curved surface is smooth and easy to clean

Control Knob

- » Ergonomic design provides smooth and easy control
- » Easy adjustment of oxygen concentration from 21% to 100%

Primary Outlet

- » CGA V-5 1240 (DISS) connector
- » With check valve
- » Compatible with GENTEC flowmeter

Oxygen Inlet

- » USA standard: CGA V-5 1240 (DISS) oxygen connector
- » Compatible with GENTEC medical hose

Side Mount

- » Easily mounted with the mounting bracket




Auxiliary Outlet

- » CGA V-5 1240 (DISS) connector
- » With check valve
- » Compatible with GENTEC flowmeter

Air Inlet

- » USA standard: CGA V-5 1160 (DISS)
- » Compatible with GENTEC medical hose



Primary Flowmeter

Auxiliary Flowmeter



Blenders

Our blenders require very little maintenance and last for years.

- » Suitable for use in Emergency Room and ICU where accurate oxygen concentration is required
- » Double-ended proportioning valve enables easy adjustment of oxygen concentration from 21% to 100%
- » Two-stage pressure-balancing mechanism is adopted to deliver steady oxygen concentration
- » Two outlet ports and bleed valve allow a wide range of applications
- » Internal bypass and alarm mechanism ensure safe use
- » Accuracy: $\pm 3\%$ FS

Step 1. Choose your Blender

- GTBH - High-Low Air/Oxygen Blender
- GTBL - Low Flow Air/Oxygen Blender

Step 2. Choose your (Left) Primary Flowmeter

- NN - None
- O3 - 0-6 LPM
- O4 - 0-15LPM

Step 3. Choose your (Right) Auxiliary Flowmeter

- NN - None
- O1 - 0-1LPM
- O2 - 0-3LPM
- O3 - 0-6LPM
- O4 - 0-15LPM

Example: GTBH0405 - High-Low Flow Air/O2 Blender
15-120LPM,0-15LPM L Flowmeter,0-50LPM R Flowmeter

Medical Gas Hose Assemblies

Medline offers custom-configured medical gas hose assemblies in four easy steps:

» 1/4" ID Non-Conductive

1. Select hose type
2. Select the inlet fitting (gas service and type)
3. Select the outlet fitting (gas service and type)
4. Specify the hose length

GT = Hose Assembly

1. Hose Type

OXY = Oxygen

VAC = Vacuum

N2O = Nitrous Oxide

AIR = Medical Air

N2 = Nitrogen

CO2 = Carbon Dioxide

WAGD = Waste Anesthetic Gas Disposal

2. Inlet Fitting (Connects to the Source)

OH: OHIO

DS: DISS Male

DH: DISS Hand-tight

CH: Chemetron® Male

PB: Puritan-Bennett Male

FOH: OHIO Coupler Female

FDS: DISS Adapter Female

FDH: DISS Hand-tight Female

FCH: Chemetron Female Coupler

FPB: Puritan-Bennett Female Coupler

3. Outlet Fitting (Connects to the Medical Device)

OH: OHIO

DS: DISS Male

DH: DISS Hand-tight

CH: Chemetron Male

PB: Puritan-Bennett Male

FOH: OHIO Coupler Female

FDS: DISS Adapter Female

FDH: DISS Hand-tight Female

FCH: Chemetron Female Coupler

FPB: Puritan-Bennett Female Coupler

4. Hose Length, in feet (e.g. 2= 2ft, 4 =4ft, 7=7ft, etc)

EXAMPLE: GTOXYFDHDS5 1/4" ID non-conductive oxygen hose assembly with a DISS Hand-tight female inlet and DISS Male outlet, 5' in length.

For custom configurations, Contact your Medline Representative.





CPRB87251B2C

HCS8715M

HCS5415M

HCS8715MA

Medline Regulators

For all medical and emergency applications.

- » Fitted with brass sleeves around key parts
- » Laser-etched knob and body

Economy

Item No.	Connection	Flow Range	Flow Increments	Pkg
HCS8708M	870 CGA Connection	0-8 LPM	0.12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8	ea
HCS8715M	870 CGA Connection	0-15 LPM	0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15	ea
HCS8725M	870 CGA Connection	0-25 LPM	0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25	ea
HCS5408M	540 CGA Connection	0-8 LPM	0.12, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8	ea
HCS5415M	540 CGA Connection	0-15 LPM	0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15	ea

Super Mini

HCS8715MA	870 connection	0-15 LPM	0.12, 0.5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15	ea
-----------	----------------	----------	--	----

Accessories

Item No.	Description	Pkg
HCS166840	Brass Seal Washer for O ₂ regulators	50/bx
HCS0700	Cylinder Wrench	50/bx

Medline Regulators for EMS

Easily transport ventilator patients.

- » Maximum inlet pressure: 3,000 PSI
- » Chrome plated brass body
- » Solid brass core


Item No.	50 PSI Power Take-Off	Hose Barb	Liters Per Minute	Cylinder Style	Pin Index	Pkg
CPRB34501C	1			CGA346	Air	ea
CPRB87251B	0	●	0-5-1-2-3-4-5-6-8-10-15-25	CGA870	Oxygen	ea
CPRB87251B2C	2	●	0-5-1-2-3-4-5-6-8-10-15-25	CGA870	Oxygen	ea
CPRB54251B2C	2	●		CGA540	Oxygen	ea
CPRB54402C	2			CGA540	Oxygen	ea
CPRB54501C	1			CGA540	Oxygen	ea

EconO₂mizer II Single-Lumen Pneumatic Conservor

Less noise for better compliance

Among the quietest pneumatic conserving devices available.

- » Very quiet with no hissing during pulse delivery
- » Aids in compliance with ambulatory patients
- » 6 conserving flow settings (1, 2, 2.5, 3, 4, & 5)
- » 3 continuous settings (2, 3, & 4)
- » Features all-brass construction in high-pressure zones
- » Weight: 13-oz.
- » Size: 4.8" x 2.4" x 1.5" (12.2 cm x 6.1 cm x 3.8 cm)
- » 5.7:1 conserving ratio



HCS2CON1L

Item No.	Description	Pkg
HCS2CON1L	Single-Lumen Pneumatic Conservor-870 connection	ea
HCSTHANDLE	Replacement T Handle, Black	ea
HCS166970	Replacement Tubing Barb Port Outlet to both HCSO2CON1 and HCSO2CONS	ea

EconO₂mizer Dual-Lumen Conserving Device

A breakthrough in conserving-device technology!

- » All-brass diaphragm and 4:1 conserving ratio for fewer service calls
- » Light and compact
- » 0-6 LPM
- » Dual mode for conserving and continuous flows
- » 11 settings for each mode (0,1, 1.5, 2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5 & 6) (no additional regulator required)
- » Solid brass core



Dual Mode Switch

HCSO2CONS

Item No.	Description	Pkg
HCSO2CONS	Conservor, 0-6 Settings, 870 Connection	20/cs
HCSTHANDLE	Replacement T Handle, Black	20/cs

Suction Regulators

Exclusively distributed nationwide by Medline.

- » Antimicrobial* and environmentally friendly** additive in the rotary knobs and selector switches
- » Extremely durable
- » Continuous-only and continuous/intermittent designs
- » Available with a variety of inlet and outlet fittings and suction trap assemblies
- » Warranty: ten years (manufacturer)

User-friendly mode selector

Large, user-friendly mode selector switch and a 1.25" (3.2 cm) diameter adjustment knob

Large diameter gauge and Easy-to-read numbers

Large 2.5" (6.4 cm) diameter gauge with glow-in-the dark scale and arrow


Color-coded sticker

Color-coded to help facilities manage inventory and allocation across different departments




Antimicrobial additive

Rotary knobs and selector switches are made with an antimicrobial additive*



Continuous/Intermittent (881 Series)

- » Available with ranges of 0-160 mmHg and 0-300 mmHg
- » Settings include continuous, off and intermittent
- » Factory default settings for on/off cycles are 16 seconds (± 2) on, 8 seconds (± 2) off; on/off settings can be adjusted without opening the case (controls are not accessible without tools)



Continuous Only (882 Series)

- » Available with ranges of 0-300 mmHg and 0-760 mmHg
- » Settings include regulated, off and full



























*This product is treated with an antimicrobial to protect the product itself.

**Agion's silver antimicrobial has been certified for its material content, recyclability, and manufacturing characteristics. EPA 71227-4-88165

Suction Regulator Ordering Matrix

Instructions: Select suction regulator, adapter then fitting (steps 1-3).
Combine codes to create model number, then write chosen sticker color. See example below.

STEPS:	1 SUCTION REGULATOR	2 ADAPTER	3 FITTING	MODEL NO.	COLOR				
Example:	GT881VR160	+	OH	+	T	=	GT881VR160OHT	in	YL
Please Write Your Selection:	_____	+	_____	+	_____	=	_____	in	_____

1 SUCTION REGULATORS	2 ADAPTERS		3 FITTINGS		STICKER COLOR OPTIONS
 <p>Continuous / Intermittent 0-160 GT881VR160</p>	<p>NO ADAPTER</p> <p>1/8" NPT Female</p> <p>OO</p>	 <p>Chematron Adapter GT3103VACNT2 CH</p>	<p>NO FITTING</p> <p>1/8" NPT Female</p> <p>A</p>	 <p>Hose Barb GT880VR304 D</p>	<ul style="list-style-type: none">  Blank (White)  YL (Yellow)  BL (Blue)  RD (Red)  PL (Purple)  PC (Peach)  PK (Pink)  GY (Gray)  GN (Green)  BB (Baby Blue) <p>* Color must be specified in sales notes.</p>
 <p>Continuous / Intermittent 0-300 GT881VR300</p>	 <p>Ohio* Adapter GT3111UVACNT2 OH</p>	 <p>DISS Hand-Tight Adapter GT3222UVACNT2 DH</p>	 <p>Trap Fitting GT881VR14K B</p>	 <p>DISS Male GTDISSV18NPT C</p>	
 <p>Continuous Only 0-300 GT882VR300</p>	 <p>Puritan-Bennett Adapter GT3104VACNT2 PB</p>	 <p>DISS Nut-Nipple Adapter GT3102VACNT2 DS</p>	 <p>Complete Suction Trap with Trap Fitting GT880VT5 T</p>		
 <p>Continuous Only 0-760 GT882VR760</p>	 <p>Oxequip Twist Adapter GT3311VACNT2 OQ</p>	 <p>Oxequip Medstar Adapter GT3301VACNT2 DE</p>	 <p>Complete Suction Trap with DISS GT880VT6 S</p>	 <p>Complete Suction Trap with DISS Hand-Tight GT880VT7 H</p>	

ACCESSORIES	Suction Traps				
	 <p>Suction Trap with No Fitting on Top GT880VT0</p>	 <p>Suction Trap with Male Trap Fitting Only GT880VT1</p>	 <p>Suction Trap with DISS Female Nut and Nipple GT880VT2</p>	 <p>Suction Trap with DISS Female Hand-Tight GT880VT3</p>	 <p>Suction Trap with Tubing Nipple GT880VT4</p>
 <p>Vacuum Adapter DISS to Hose Barb GTDISSMVHCHK</p>	 <p>Extension Adapter 2-1/2" Long GTEXTADPT18</p>	 <p>Security Cable GTFMS2</p>			



HCS7000

Suction Aspirator and Disposable Canister Kits

Our easy-to-use aspirators are ultra-quiet.

- » Weight: 12.1-lb./5.49 kg
- » Maximum suction: -21 inHg/520 mmHg
- » Rugged construction is ideal for bedside carts
- » Disposable 850 cc or 1200 cc suction canister sets with float lids.
- » HCS7851 is a kit that also includes suction tubing (1 each of 6' and 18" length).
- » Replacement lids are not available

Item No.	Description	Pkg
HCS7000	Portable Suction Aspirator	1/ea

Canister with Float Lid

HCS7850	Canister with Float Lid: 850 cc	12/cs
HCS7851	Canister with Float Lid and Tubing: 850 cc	12/cs
HCS7120	Canister with Float Lid: 1,200 cc	12/cs
HCS7121	Canister with Float Lid and Tubing: 1,200 cc	12/cs

Canister with Filter Lid

HCS7852	Canister with Filter Lid and Tubing: 850 cc	50/cs
HCS7853	Canister with Filter Lid: 850 cc	100/cs

Accessories

HCS4051	Hydrophobic Filter Set, Barb Both Ends	1/pk
HCS4053	Hydrophobic Filter Set, Barb Both Ends	3/pk



HCS7850

HCS7120



Tubing included with HCS7851, HCS7852, and HCS7121

Humidification & Filtration

Heat & Moisture Exchangers	22
Heat & Moisture Exchanger Filters	23
TrachVox HME	24
Bacterial/Viral Filters	24
Ventilator Circuits, Unheated	24
Pulmonary Function Filters	25

Heat and Moisture Exchangers and Filtered Heat and Moisture Exchangers

Medline's broad portfolio of HMEs and HMEFs meet your patient's passive humidification needs throughout the care continuum.

- » Quality product for confidence in delivering passive humidification
- » Wide range of tidal volumes
- » Includes non-filtered, electrostatic-filtered and HEPA-filtered options



	DYNJAAHME10	DYNJAPHME18	DYNJAPHME20
DESCRIPTION	HME, Small with Port	HME, Neonatal	HME, Pediatric with Port
Type	Straight	Straight	Straight
Connection	22M/15F x 22F/15M	15M X 15F	22M/15F x 22F/15M
Vt (ml)	120-1500	10-50	90-1500
Vd (ml)	32.1	2	20.9
Weight (g)	14.7	2	13.3
CO2 Port	Yes	No	Yes
Resistance (cm H2O)	0.6 @ 30 lpm	145pa @ 10 lpm	0.4 @ 30 lpm
Moisture Output (Mg H2O/L)*	33.7 @ 500 mL	32	34.0 @ 250 mL
Media	Foam	Foam	Foam
Pkg	50/cs	50/cs	50/cs



	DYNJAAHME12	DYNJAAHME40
DESCRIPTION	HME, Small, No Port	HME, Trach with Oxygen Port
Type	Straight	N/A
Connection	22M/15F x 22F/15M	15MM ID
Vt (ml)	120-1500	>50
Vd (ml)	32.1	13.4
Weight (g)	14.7	7.1
CO2 Port	No	No
Resistance (cm H2O)	0.6 @ 30 lpm	0.42 @ 30 lpm
Moisture Output (Mg H2O/L)*	33.7 @ 500 mL	26
Media	Foam	
Pkg	50/cs	50/cs

Our Humidification and Filtration offering is expanding, for the latest options visit our website or call your Medline Representative

Filter HMEs

Medline's wide range of choices ensures compatibility with your equipment.

- » Choose from both angled and straight designs
- » Electrostatic and HEPA configurations available
- » Translucent housing makes it easy to see possible obstructions
- » HME media foam composition



	DYNJAAHME1	DYNJAAHME5	DYNJAAHME9	DYNJAAHME11
Description	HMEF, Large with Port	HMEF, Large, No Port	HMEF, Small with Port	HMEF, Small, No Port
Type	Straight	Straight	Straight	Straight
Connection	22M/15F x 22F/15M	22M/15F x 22F/15M	22M/15F x 22F/15M	22M/15F x 22F/15M
Vt (ml)	150-1500	150-1500	120-1500	120-1500
Vd (ml)	40.7	40.7	31.3	31.3
Weight (g)	18.5	18.5	15	15
CO2 Port	Yes	No	Yes	Yes
Resistance (cm H2O)	0.8 @ 30 lpm	0.8 @ 30 lpm	2.1 @ 30 lpm	2.1 @ 30 lpm
Moisture Output (Mg H2O/L)*	34.1 @ 500 mL	34.1 @ 500 mL	33.6 @ 500 mL	33.6 @ 500 mL
Filter Media	ElectroStatic	ElectroStatic	Electrostatic	Electrostatic
BFE(%)	99.999	99.999	99.999	99.999
VFE(%)	99.999	99.999	99.999	99.999
Pkg	50/cs	50/cs	50/cs	50/cs



	DYNJAAHME13	DYNJAPHME17	DYNJAAHMEHEPA	DYNJ603011
Description	HMEF, Small, Angled with Port	HMEF, Pediatric with Port	HEPA Filtered HME with CO2 Port	Filter,HME,Flat/Rec, Large, Sample Port
Type	Angled	Straight	8	Straight
Connection	22M/15F x 22F/15M	22M/15F x 22F/15M	22M/15F x 22F/15M	22ID x 22OD/15ID
Vt (ml)	120-1500	90-1500	200-1500	200 - 1500
Vd (ml)	33.7	20.5	52	62
Weight (g)	17.2	13.6	40	37.2
CO2 Port	Yes	Yes	Yes	Yes
Resistance (cm H2O)	1.9 @ 30 lpm	2.4 @ 30 lpm	1.5 @ 30 lpm	<1 @ 30 lpm
Moisture Output (Mg H2O/L)*	33.4 @ 500 mL	34.0 @ 250 mL	30	32.5
Filter Media	Electro-Static	Electro-Static	HEPA	
BFE(%)	99.999	99.999	99.9999	>99.99%
VFE(%)	99.999	99.999	99.9999	>99.99%
Pkg	50/cs	50/cs	50/cs	50/cs

TrachVox® Speaking HME

TrachVox HME helps spontaneously breathing tracheostomy patients verbally communicate. When depressed, the HME changes the airway flow to make speaking possible. Additionally, it heats and humidifies inspired air and delivers supplemental oxygen through an integrated port.

- » Allows for force inspiration, expiration and suctioning
- » Low flow oxygen can be attached directly to the TrachVox oxygen port
- » Uses standard 15 mm connector

Item No.	Description	Pkg
EMS3131010	Trachvox Speaking valve with HME	30/cs



EMS3131010

Filters

Exceptional bacterial and viral filtration efficiency—up to 99.9999 percent!

- » Translucent housing makes it easy to see vomitus or other possible obstructions in the system

	DYNJAAHEPA	DYNJAABV1
Orientation	Straight	Straight
Vt (ml)	200-1500	90-1500
Vd (ml)	52	30
Weight (g)	40	23
CO2 Port	YES	NO
Resistance @ 15LPM (cm H ₂ O)	—	—
Resistance @ 30LPM (cm H ₂ O)	1.48	0.77
Resistance @ 60LPM (cm H ₂ O)	3.33	1.78
Filter	HEPA	Electro-Static
BFE(%)	99.9999	99.999
VFE(%)	99.9999	99.999
Pkg	50/cs	50/cs



DYNJAABV1

DYNJAAHEPA



Non-Heated Breathing Circuits

Medline's family of non-heated breathing circuits support your mechanically ventilated patients with several configurations, which includes "Y" and parallel "Y" options.

Item No.	Description	Pkg
HCS1571215	Dual limb circuit, with temperature sensing port, 5 ft.	10/cs
HCS1830	Dual limb circuit, 6 ft.	10/cs
HCS1831	Dual limb circuit with wye and temperature sensing port, 6 ft.	10/cs
HCS1831P	Dual limb circuit with parallel wye and temperature sensing port, 6 ft.	10/cs
HCS72R7002	Dual limb circuit with wye, 6 ft.	10/cs
HCS72R7401	Dual limb circuit with parallel wye and HME	10/cs

Pulmonary Function Filters

Medline offers a variety of pulmonary function filters that are compatible with all spirometry machines. Please contact your Medline rep or dial 1-800-Medline(633-5463), for ordering information.

- » Independent testing shows >99.99% microbiological filtration effectiveness
- » Low resistance performance ensures efficacy of results
- » Low functional volume
- » Flexible, comfortable, disposable mouth piece helps with patient compliance and improves effectiveness
- » Comfortable single-patient nose clip prevents patient-to-patient touch contamination



Item No.	Description	OD	ID	PFT Equipment	Pkg
DYNJ756ELB	Elliptical w/Lip	48.6	45.5	Koko, KokoLegend, Keystone, Koko-Trek, Cybermedic Full Gas Systems, Collins 4G and 3G Full Gas Systems	100/cs
DYNJ756NB	Round	48.6	45.5		100/cs
DYNJ756SB	Children and Adolescent	48.6	45.5		100/cs
DYNJ757BMM	Kit w/mouthpiece & nose-clip	34.3	30.1	SensorMedics, 2200, 2400, 2450, 2800, Vmax, Keystone, I&II, Pail, Cybermedics, Moose, Gould 5004	100/cs
DYNJ757ELB	Elliptical w/Lip	34.3	30.1		100/cs
DYNJ757NB	Round	34.3	30.1		100/cs
DYNJ757SB	Children and Adolescent	34.3	30.1		100/cs
DYNJ758B	Round	30.1	25.4	Incheck, InCheck DIAL, Multispiro, Creative Biomedics, Vigtalograph, Schiller, MicroDirect, SDI (PulmogardII), Most Cosmed Models, MIR (All Models), P.K.Morgan 232, Future Med(Spirobank), PB-400, PB-600, CDX Spiro 850, Spirometrics (Flowmate, Flowmate+), SMII, Vacumed (MicroPlus, Micro GP-GL, Micro Loop, Pony)	100/cs
DYNJ758BMM	Kit w/mouthpiece & nose-clip	30.1	25.4		100/cs
DYNJ758ELB	Elliptical w/Lip	30.1	25.4		100/cs
DYNJ760B	Round	38.8	34.5	Brentwood 4000, Medical Graphics	100/cs
DYNJ760BMM	Kit w/mouthpiece & nose-clip	38.8	34.5		100/cs
DYNJ760ELB	Elliptical w/Lip	38.8	34.5		100/cs
DYNJ780B	Round	34.3	30.8	Morgan	100/cs
DYNJ780BMM	Kit w/mouthpiece & nose-clip	34.3	30.8		100/cs
DYNJ780ELB	Elliptical w/Lip	34.3	30.8		100/cs

Accessories

- » Flexible, comfortable, disposable mouth piece helps with patient compliance and improves effectiveness
- » Comfortable single-patient nose clip prevents patient-to-patient touch contamination

Item No.	Description	Pkg
DYNJ750	Standard Nose Clip	100/bag
DYNJ751	Marshmallow Nose Clip	100/bag
DYNJ752	Large Nose Clip	50/bag
DYNJ754	Disposable Rubber Mouthpiece	50/bag




Medication Delivery & Asthma Management

Medication Nebulizers, Small Volume	27
Aerosol Masks	28
Valved Tee Adapters	28
Nebulizer Compressors, Compact	29
Nebulizer Compressors, Characters	30
Nebulizer Compressors, Portable	31
MDI Spacers	32
Peak Flowmeters	32

Handheld Nebulizer Kits

Medline offers a wide selection of reusable and disposable nebulizer kits for medication delivery.

- » Ideal particle size of less than five microns, increasing the amount of medication delivered deep within lungs
- » Standard fittings accommodate most clinical applications
- » High respirable dose, leading to more efficient treatments and less wasted medication
- » Not Made with Natural Rubber Latex
- » Medicine can be added without removing the cap
- » Kits include a medicine cup, 7-ft. (2.1 m) of kink-resistant extension tubing, and either a mask or a mouthpiece with a T adaptor
- » Packaged 50/cs



Upstream Nebulizers

Item #	Tee Adaptor	Valved Tee	Mouthpiece	Reservoir Tube	Tubing Length	Tubing Connector	Mask Type	Disposable/Reusable**
HCS4485					7 foot	Standard	Adult	Disposable
HCSU4485					7 foot	Universal	Adult	Disposable
HCS4483	22mm OD/ID		●	●	7 foot	Standard		Disposable
HCS448314	22mm OD/ID		●	●	14 foot	Standard		Disposable
HCSU4483	22mm OD/ID		●	●	7 foot	Universal		Disposable
HCSU448314	22mm OD/ID				14 foot	Universal		Disposable
HCS4482	22mm OD/ID		●		7 foot	Standard		Disposable
HCSU4482	22mm OD/ID		●		7 foot	Universal		Disposable
HCS4486					7 foot	Standard	Pediatric	Disposable
HCSU4486					7 foot	Universal	Pediatric	Disposable
HCS4480*					None	None		Disposable
HCS4483VT		22mm OD/OD		●	7 foot	Standard		Disposable
HCS4481VT		22mm OD/ID			7 foot	Standard		Disposable
HCS4482R	22mm OD/ID		●		7 foot	Standard		Reusable
HCS4483R	22mm OD/ID		●	●	7 foot	Standard		Reusable
HCS4482RR	45° Angled		●		7 foot	Standard		Reusable



Sidestream Nebulizers

Item #	Tee Adaptor	Valved Tee	Mouthpiece	Reservoir Tube	Tubing Length	Tubing Connector	Mask Type	Disposable/Reusable**
HCS4493	22mm OD/OD		●	●	7 foot	Standard		Disposable
HCS4492	22mm OD/OD		●		7 foot	Standard		Disposable
HCS4495	22mm OD/OD				7 foot	Standard	Adult	Disposable



*Neb Cup only, no additional components **Disposable up to 2 weeks of use **Reusable up to 6 months of use



HCS4631B



HCS4630

Aerosol Masks

Extra patient comfort comes from an anatomical form and a soft, light and comfortable vinyl material.

- » A secure fit with elastic straps
- » Standard 22 mm connector
- » Not made with natural rubber latex

Adult

Item No.	Description	Pkg
HCS4630B	Adult Elongated Aerosol Mask with Straight Connector	50/cs
HCS4630	Adult Elongated Aerosol Mask with Swivel Connector	50/cs



HCS4630B

Pediatric

HCS4631B	Pediatric Elongated Aerosol Mask with Swivel Connector	50/cs
----------	--	-------

Valved Tee Adapters

- » For inline medication delivery
- » Available in 22 mm OD/OD and 22 mm OD/ID

Item No.	Description	Pkg
HCSA1008	Valved Tee Adapter, 22 mm OD/OD	30/cs
HCSA1009	Valved Tee Adapter, 22 mm OD/ID	30/cs



Aeromist Compact™ Nebulizer Compressors

Our smallest nebulizer compressor is less bulky and takes up less of your valuable space. Choose from the standard white model to a more fun and stylish alternative with our color options.

- » Kits include a disposable nebulizer medicine cup, 7-ft. (2.13 m) of tubing, 7' power cord and a mouthpiece
- » Color choices: electric blue, neon green, neon pink, charcoal, or white
- » Durable hard plastic shell
- » Rugged compressor rated for continuous use
- » Compact and stable, with a small footprint, low center of gravity and four non-skid rubber feet
- » Warranty: 5 years
- » Replacement Kits: HCS4483 (Disposable) and HCS4483R (Reusable)



Specifications

Method of nebulization: dry piston compressor
 Electrical Requirements: 120VAC/60Hz
 Maximum jet nebulizer capacity: 10 ml
 Sound noise level: ≤ 55 dBA
 Pressure: 30 PSI, Flow rate: 8 Lpm



Item No.	Description	Pkg
HCS70004G	Neon Green Compressor and Nebulizer Kit	8/cs
HCS70004BL	Electric Blue Compressor and Nebulizer Kit	8/cs
HCS70004P	Neon Pink Compressor and Nebulizer Kit	8/cs
HCS70004C	Charcoal Compressor and Nebulizer Kit	8/cs
HCS70004	White Compressor and Nebulizer Kit	8/cs
HCS70004RD	White Compressor with Reusable and Disposable Neb Kit	8/cs
HCS70004DP	Sample Pack, Two of each color	8/cs
HCS70004F	Replacement Filters	5/pk
HCS70004CAP	Replacement Cap for HCS70004	8/pk



Components included with each compressor

Aeromist Buddies™ Nebulizer Compressors

Kid-friendly nebulizers help put young patients at ease, while still providing safe and effective aerosol treatments. Their designs help to calm fears and make treatments more fun.

- » Kit includes a disposable nebulizer medicine cup, 7-ft. (2.13 m) of tubing, mouthpiece, and pediatric mask
- » Durable hard plastic shell
- » Rugged compressor rated for continuous use
- » Compact and stable, with a small footprint, low center of gravity and four non-skid rubber feet
- » Warranty: 5 years

Item No.	Description	Pkg
HCSMONKEY	MonkeyNeb Nebulizer Character Compressor	8/cs
HCSMONKEYCAP	Replacement Filter Cap for HCSMONKEY	8/pk
HCSFISHNEB*	FishNeb Nebulizer Character Compressor	8/cs
HCSFISHNEBCAP	Replacement Filter Cap for HCSFISHNEB	8/pk
HCSPUDDY*	PuppyNeb Nebulizer Character Compressor	8/cs
HCSPUDDYCAP	Replacement Filter Cap for HCSPUDDY	8/pk
HCSTURTLE	Turtle Nebulizer Character Compressor	8/cs
HCSTURTLECAP	Replacement Filter Cap for HCSTURTLE	8/pk
HCS70004F	Replacement Filters	5/pk

*Comes with Grab-N-Go backpack for convenient storage.

All compressors include a complementary nebulizer kit.



Specifications

Method of nebulization: dry piston compressor
 Electrical Requirements: 120 VAC/ 60Hz
 Maximum nebulizer capacity: 10 ml
 Sound noise level: ≤ 60 dBA
 Pressure: 30 PSI, Flow rate: 8 Lpm



Aeromist Trek Portable Nebulizer Compressor

Better compliance with easy treatments, anytime and anywhere.

- » Suitable for travel with a sleek design, extra small size, light weight and battery operation
- » Durable hard plastic shell
- » Compact and stable, with a small footprint
- » Kit includes a disposable nebulizer medicine cup, 7-ft. (2.13 m) of tubing, mouthpiece, rechargeable battery, carry case, AC and DC adapters, adult and pediatric masks, and 5 replacement filters
- » Warranty: 3 years for the compressor, 6 months for the battery



Item No.	Description	Pkg
HCS2GO	Portable Battery-Operated Nebulizer Compressor Kit	each
HCS2GONB	Portable Nebulizer Compressor, Battery not included	each
HCS2GOBATT	Replacement Battery for Portable Nebulizer Compressors (Items HCS2GO and HCS2GONB)	each
HCS2GOF	Replacement filters	5/pk
HCS2GOAC	Replacement AC Adapter	each
HCS2GODC	Replacement Car Charger	each

Medline MDI Spacers

Medline spacers are designed for patients who have difficulty with metered dose inhalers (MDIs). With spacers, more medication is deposited in the lungs and less is deposited in the mouth and throat.

- » A simple and economical way to improve effectiveness of inhaler
- » Easy to use and helps avoid some of the unpleasant taste of inhaled medication
- » Anti-static and compatible with any standard mask with a 22 mm fitting
- » Masks are made of a soft silicone and fit snugly to ensure optimal drug delivery
- » The collapsible model is portable and convenient

Item No.	Description	Pkg
HCSA100	MDI Spacer Only	10/bx
HCSA200	MDI Spacer Only, Collapsible	10/bx
HCSA112	MDI Spacer with Infant Mask	10/bx
HCSA113	MDI Spacer with Pediatric Mask	10/bx
HCSA114	MDI Spacer with Adult Mask	10/bx



Peak Flow Meters

Peak flow meters help patients manage their asthma and gauge the effectiveness of their treatment program, which may help them predict their next attack.

- » Individually calibrated to ensure accuracy
- » Durable for years of reliable performance
- » The pocket model offers superb convenience and portability
- » Available in low range and full range versions – low range models are better for geriatric, pediatric or other patients with diminished flow rates



Peak Flow Meter Color Zones²

Green Zone

80- 100% of normal peak flow rate. Asthma is reasonably under control. Continue prescribed management program.

Yellow Zone

50- 80% of normal peak flow rate. Airways are narrowing and may require extra treatment. Alert physician.

Red Zone

Less than 50% of normal peak flow rate. Severe airway narrowing signal. Immediate decisions and actions need to be taken. Alert physician.

Item No.	Description	Pkg
HCSPF100	Medline Portable Peak Flow Meter, Full Range, 60-720 lpm	12/bx
HCSPF200	Standard Peak Flow Meter, Full Range, 60-850 lpm	12/bx
HCSFPEDIH	Mini Wright Pediatric Peak Flow Meter	each
HCSPF220	Disposable Pediatric Mouthpiece, Adapter	20/bx
HCSPF262	Disposable One Way Valve Mouthpiece, Adult	100/cs

Non-Invasive Interfaces

NIV Masks	34
NIV Circuits	34
CPAP Masks AirGel	35-36
CPAP Masks Silicone	37
CPAP Accessories	38



Veraseal 2 Non-Invasive Ventilation Masks

Veraseal 2 NIV Masks feature the advanced, Air^ogel cushion technology, which promotes advanced skin protection, superior performance, and comfort. It is available in multiple sizes and configurations to meet your non-invasive ventilation needs.

Item No.	Description	Pkg
CPAP50827BX	Veraseal [®] 2 NIV Mask, Vented w/ AAV, Small	10/bx
CPAP50828BX	Veraseal [®] 2 NIV Mask, Vented w/ AAV, Medium	10/bx
CPAP50829BX	Veraseal [®] 2 NIV Mask, Vented w/ AAV, Large	10/bx
CPAP50830BX	Veraseal [®] 2 NIV Mask, Vented w/ AAV, X-Large	10/bx
CPAP50937BX	Veraseal [®] 2 NIV Mask, Non-Vented w/ AAV, Small	10/bx
CPAP50938BX	Veraseal [®] 2 NIV Mask, Non-Vented w/ AAV, Medium	10/bx
CPAP50939BX	Veraseal [®] 2 NIV Mask, Non-Vented w/ AAV, Large	10/bx
CPAP50940BX	Veraseal [®] 2 NIV Mask, Non-Vented w/ AAV, X-Large	10/bx
CPAP50933BX	Veraseal [®] 2 NIV Mask, Non-Vented, Small	10/bx
CPAP50934BX	Veraseal [®] 2 NIV Mask, Non-Vented, Medium	10/bx
CPAP50935BX	Veraseal [®] 2 NIV Mask, Non-Vented, Large	10/bx
CPAP50936BX	Veraseal [®] 2 NIV Mask, Non-Vented, X-Large	10/bx



Non-Invasive Circuits

Medline offers smooth bore single limb circuits for use with non-invasive ventilation. These 72" circuits include a disposable exhalation port and pressure monitoring tubing as well as an option with a bacterial/viral filter.

Item No.	Description	Pkg
DYNJCIRC10*	Universal NIV Circuit, 72"	10/cs
DYNJCIRC10F*	Universal NIV Circuit with Filter, 72"	10/cs

*Coming soon. Call for availability.

AIRgel™ CPAP Masks

Amazing comfort and customized fit technology. The patented AIRgel replaceable cushioning minimizes discomfort and skin irritation while effectively sealing airflow. It is easily detached and reattached for easy cleaning, saving time and money.

- » Proprietary AIRgel cushioning minimizes discomfort and skin irritation while effectively sealing airflow
- » The flexible and formable mask provides a custom fit—flex, fit and forget!

Aura Nasal Mask

- » The shapeable shell of each mask can be easily molded to fit individual needs
- » Flexible spacebar is designed to provide superior stability without touching the skin
- » Active headgear connector moves with the patient to help maintain an effective seal

Aura-Nasal

Item No.	Description	Mask Size	Pkg
CPAP50176	Nasal Mask With Adjustable Gel Cushion	Small	ea
CPAP50177	Nasal Mask With Adjustable Gel Cushion	Medium	ea
CPAP50178	Nasal Mask With Adjustable Gel Cushion	Large	ea

Mojo® Full-Face Gel Mask

- » Full-face design in four sizes
- » Adjustable tension wheel minimizes discomfort while ergonomic cushions adapt to unique facial features for a comfortable fit
- » Magnetic headgear clips hold firmly and disconnects with one hand easy on-and-off

Item No.	Description	Headgear Style	Mask Size	Pkg
CPAP50835	Vented	One Size	Small	ea
CPAP50836	Vented	One Size	Medium	ea
CPAP50837	Vented	One Size	Large	ea
CPAP50838	Vented	One Size	XL	ea
CPAP55045	Headgear only	One Size	One Size	ea



AIRgel™ CPAP Masks

iQ Blue® Nasal Gel Mask

- » One-size nasal mask design fits a wide variety of patients
- » Clearest field of vision possible with no forehead pad required – patients can keep their glasses on while reading or watching TV
- » Easy on-and-off thanks to quick release, breathable headgear



CPAP50655

iQ Blue Nasal Gel

Item No.	Mask Style	Headgear Style	Mask Size	Pkg
CPAP50655	Vented	Stable Fit	One Size	ea
CPAP55044	Headgear Only	Stable Fit	One Size	ea

MiniMe® 2 Pediatric Gel Mask

- » Nasal mask in two sizes, designed for children ages 2 to 12
- » Touchless spacebar eliminates the need for any forehead contact
- » Soft, breathable headgear for all-night comfort
- » Available in vented and non vented designs
- » Adjustable cushion with memory wire for a true custom fit.



CPAP60217

Item No.	Description	Headgear Style	Mask Size	Pkg
CPAP60217	Vented w/Leak Port	Small	Small/Medium	ea
CPAP60218	Vented w/Leak Port	Large	Medium/Large	ea
CPAP60257	Non-Vented	Small	Small/Medium	ea
CPAP60258	Non-Vented	Large	Medium/Large	ea

Ascend® Masks

Patented AIRgel replaceable cushion is easily detached and reattached for easy cleaning and replacement. Leak ports are built-in, so a separate exhalation port is not needed.

- » Minimizes pressure points and skin irritation
- » Ergonomic cushion technology delivers an effective seal for the most comfortable fit
- » Flexible spacebar with a unique bendable design provides superior stability, holding a seal without contacting the skin
- » Lightweight, breathable headgear slides off easily with one hand
- » 360° rotating swivel for optimal tube positioning and easier mobility



CPAP50815

Ascend Full Face Mask

Item No.	Description	Pkg
CPAP50825	CPAP Full-Face Mask Starter kit with Small, Medium, and Large Cushions	1 ea
CPAP50813	CPAP Full-Face Mask, Small	1 ea
CPAP50814	CPAP Full-Face Mask, Medium	1 ea
CPAP50815	CPAP Full-Face Mask, Large	1 ea
CPAP50207	CPAP Full-Face Mask Cushion, Small	1 ea
CPAP50208	CPAP Full-Face Mask Cushion, Medium	1 ea
CPAP50209	CPAP Full-Face Mask Cushion, Large	1 ea
CPAP55051	CPAP Headgear Replacement, One Size	1 ea

Ascend Nasal Mask

Item No.	Description	Pkg
CPAP50825	CPAP Nasal Mask Starter Kit with Small, Medium, and Large Cushion	1 ea
CPAP50813	CPAP Nasal Mask, Small	1 ea
CPAP50814	CPAP Nasal Mask, Medium	1 ea
CPAP50815	CPAP Nasal Mask, Large	1 ea
CPAP50207	CPAP Nasal Mask Cushion, Small	1 ea
CPAP50208	CPAP Nasal Mask Cushion, Medium	1 ea
CPAP50209	CPAP Nasal Mask Cushion, Large	1 ea

Silicone CPAP Masks

A consistent seal without excessive tightness, delivering the prescribed therapy with enough comfort to help improve compliance.

- » OmniClip and silicone forehead pad for a custom fit at the bridge of the nose for patients of all head shapes and sizes
- » SilentVent elbow flushes out exhaled air quietly, without an additional inline exhalation port
- » Headgear attaches with convenient clips, includes a handy wash bag to keep it clean and fresh
- » Cushions, forehead pads, elbows, cones and clips are individually replaceable
- » Available in four sizes, with and without headgear



Item No.	Size	Description	Pkg
CPAP9400	Small	Full Face Mask with Head Strap » SilentVent Anti-Asphyxia Elbow, 22mm Male » OmniClip with Silicone Forehead Pad	1/ea
CPAP9401	Medium		
CPAP9402	Large		
CPAP9403	X-Large		
CPAP9404	Small	Full Face Mask, No Head Strap » SilentVent Anti-Asphyxia Elbow, 22mm Male » OmniClip with Silicone Forehead Pad	1/ea
CPAP9405	Medium		
CPAP9406	Large		
CPAP9407	X-Large		
CPAP9408	Small	Nasal Mask with Head Strap » SilentVent Elbow, 22mm Male » Silicone Forehead Pad » OmniClip with Silicone Forehead Pad	1/ea
CPAP9409	Medium		
CPAP9410	Large		
CPAP9411	X-Large		
CPAP9412	Small	Nasal Mask, No Head Gear Strap » SilentVent Elbow, 22mm Male » Silicone Forehead Pad » OmniClip with Silicone Forehead Pad	1/ea
CPAP9413	Medium		
CPAP9414	Large		
CPAP9415	X-Large		

Aloha Nasal Pillow Systems

Medline offers the Aloha nasal pillow system for CPAP therapy patients. This lightweight, nasal pillows system is quiet and offers a stable and comfortable seal.

Item No.	Description	Pkg
SUHCM014	Aloha Nasal Pillow System, Includes: One small, medium and large pillow Headgear	ea
SUHCMP014SP	Nasal Pillows for Aloha, Small	ea
SUHCMP014MP	Nasal Pillows for Aloha, Medium	ea
SUHCMP014LP	Nasal Pillows for Aloha, Large	ea
SUHHGO14	Headgear for Aloha	ea



CPAP Accessories

Tubing

Item No.	Description	Pkg
SUH0610	6' Tubing with Flex Connector	10/bx
SUHTUB010	10' Tubing with Flex Connectors	ea
SUHTUB008	8' Tubing with Flex Connectors	ea
SUHTUB0018	18" Tubing with Flex Connectors	ea



Chin Straps

Item No.	Description	Pkg
SUHCS001	Standard Chin Strap	ea
SUHCS003	Adjustable Chin Strap	ea
SUHCS004	Deluxe Chin Strap	ea
SUHCS006	Premium Chin Strap	ea



Citrus II® Cleaning Supplies

Item No.	Description	Pkg
SUH1001	Citrus II® Spray Cleaner, 1.5 oz.	24/cs
SUH1002	Citrus II® Spray Cleaner, 8 oz.	12/cs
SUH1003	Citrus II® Mask Wipes	62 Wipes/ Dispenser, 12/cs
SUH1004	Citrus II® Cleaner Concentrate, 32 oz.	ea
SUH1005	Citrus II® CPAP Wipes	12 bx/cs
SUHCAP1003S	CPAP Cleaning Wipes, Unscented	62 Wipes/ Dispenser, 12/cs



O₂ & CO₂ Disposables

Cannulas, O ₂	40
Cannulas, CO ₂	41
Tubing, O ₂	42
Tubing, Aerosol	42
Masks, O ₂	43
Masks, Aerosol	44
Bubble Humidifiers	45
O ₂ Accessories	45

SuperSoft Nasal Cannulas

Medical device-related pressure ulcers (MDRPU) have been receiving greater attention. MDRPUs tend to progress rapidly, as they often occur over areas with minimal fatty tissue. Findings indicate that patients that use a medical device are up to four times more likely to develop a pressure ulcer of any kind.¹

Even short term use of an oxygen cannula can lead to the development of pressure ulcers on the ears and cheeks. The Medline SuperSoft nasal cannula is made of an advanced material that is soft, flexible and non-tacky. Also features curved, non-flared nasal prongs. You will notice a difference right away:

- » Maintain skin integrity – Help reduce pressure and friction on the ears and cheeks.
- » Promote patient compliance with added comfort – Patients are less apprehensive about wearing a cannula that feels soft and comfortable.
- » Reduce costs – No need for the additional time and expense of padding cannulas with foam ear cushions.
- » Not made with natural rubber latex.

Adult

Item No.	Description	Connector	Pkg
HCSU4514S	7' (2.1 m) Tubing	Universal	50/cs
HCS4514S	7' (2.1 m) Tubing	Standard	50/cs
HCS4510S	No Tubing, (Cannula Harness only)	Male	50/cs
HCSU4504S	4' (1.2 m) Tubing	Universal	50/cs
HCSU4516S	14' (4.3 m) Tubing	Universal	50/cs
HCSU4515S	25' (7.6 m) Tubing	Universal	25/cs

Pediatric

Item No.	Description	Connector	Pkg
HCSU4518S	Pediatric, 7' (2.1 m) Tubing	Universal	50/cs
HCSU4517S	Infant, 7' (2.1 m) Tubing	Universal	50/cs
HCSU4521S	Preemie, 7' (2.1 m) Tubing	Universal	50/cs
HCSU4520S	Neonatal, 7' (2.1 m) Tubing	Universal	50/cs

Soft-Touch Cannulas

- » Curved, non-flared nasal prongs
- » Available in a wide variety of sizes and lengths
- » Crush-resistant tubing
- » Not made with natural rubber latex

Adult

Item No.	Description	Connector	Pkg
HCS4514	7' (2.1 m) Tubing	Standard	50/cs
HCS4510	Cannula Harness w/o Tubing	Male	50/cs
HCS4504B	4' (1.2 m) Tubing	Standard	50/cs
HCS4516B	14' (4.3 m) Tubing	Standard	50/cs
HCS4515B	25' (7.6 m) Tubing	Standard	25/cs
HCSU4514	7' (2.1 m) Tubing	Universal	50/cs
HCSU4504B	4' (1.2 m) Tubing	Universal	50/cs
HCSU4516B	14' (4.3 m) Tubing	Universal	50/cs
HCSU4515B	25' (7.6 m) Tubing	Universal	25/cs
HCSO2KIT1	Oxygen Setup Kit*	N/A	10/cs

*Contains two adult soft-touch cannulas (HCS4514), Christmas tree adapter (HCS4555), empty bubble humidifier (HCS64377), oxygen tubing connector (HCS65100) and 25' violet oxygen tubing (HCS4525V)



Crush-Resistant Tubing
This multi-channel interior feature helps to maintain a strong airflow



Pediatric

Item No.	Description	Connector	Pkg
HCS4518	7' (2.1 m) Tubing	Standard	50/cs
HCSU4518	7' (2.1 m) Tubing	Universal	50/cs

Infant

Item No.	Description	Connector	Pkg
HCS4517	7' (2.1 m) Tubing	Standard	50/cs
HCSU4517	7' (2.1 m) Tubing	Universal	50/cs

Standard Cannulas

- » Crush-resistant tubing with standard and universal oxygen connectors
- » Not made with natural rubber latex

Item No.	Description	Connector	Pkg
HCS4511B	Straight Prongs, 7' (2.1 m) Tubing	Standard	50/cs
HCSU4511B	Straight Prongs, 7' (2.1 m) Tubing	Universal	50/cs



High-Flow SuperSoft Cannulas

- » Made of an exceptionally soft material for increased patient comfort
- » Adult liter flow from 6LPM to 15LPM
- » For use with high-flow bubble humidifiers
- » Not made with natural rubber latex

Item No.	Description	Pkg
HCS4514HF	Adult, 7' (2.1 m) Tubing, Standard Connector	25/cs
HCSS4516HF	Adult, 14' (4.3 m) Tubing, Standard Connector	25/cs
HCS4515HF	Adult, 25' (7.6 m) Tubing, Standard Connector	10/cs
HCS4514HFK	High Flow Kit: 1 HCS4514HF, 1 VTB6900	7/cs



CO₂ Nasal Sampling Cannulas

Our comprehensive line of CO₂ nasal sampling cannulas features a "Pull N' Peel" design to resist tangling, while also providing accurate, uninterrupted CO₂ sampling and oxygen delivery. Universal O₂ connectors save you from searching for adaptors.



- » Universal O₂ Connector is designed to fit any flowmeter without the need for adaptors
- » Male and female luer connectors are available to meet caregiver needs

CO₂ Nasal Sampling Cannulas with Male Luer Connectors

Item No.	Description	Gas Sampling Line	Pkg
HCS4560	Adult, 7 ft.	4" pigtail	25 ea/cs
HCS4568	Adult, 7 ft.	7 ft.	25 ea/cs
HCS4568B	Adult, 7 ft.	7 ft.	100 ea/cs
HCS4571*	Adult, 7 ft.	7 ft.	25 ea/cs
HCS4565	Adult, 10 ft.	10 ft.	25 ea/cs
HCS4565B	Adult, 10 ft.	10 ft.	100 ea/cs
HCS4566	Adult, 14 ft.	14 ft.	25 ea/cs
HCS4566B	Adult, 14 ft.	14 ft.	100 ea/cs
HCS4562	Pediatric, 7 ft.	4" pigtail	25 ea/cs

CO₂ Nasal Sampling Cannulas with Female Luer Connectors

Item No.	Description	Gas Sampling Line	Pkg
HCS4561	Adult, 7 ft.	4" pigtail	25 ea/cs
HCS4568F	Adult, 7 ft.	7 ft.	25 ea/cs
HCS4568FB	Adult, 7 ft.	7 ft.	100 ea/cs
HCS4572*	Adult, 7 ft.	7 ft.	25 ea/cs
HCS4565F	Adult, 10 ft.	10 ft.	25 ea/cs
HCS4563	Pediatric, 7 ft.	4" pigtail	25 ea/cs



Oxygen Tubing

- » Unique, relaxed, crush resistant, low-memory material helps prevent kinking
- » Selected products in green color to provide greater visual contrast to assist in patient safety
- » For standard and high flow applications
- » Not made with natural rubber latex

Clear

Item No.	Description	Connector	Pkg
HCS4507	7' (2.1 m) Oxygen Tubing, Clear	Standard	50/cs
HCS4524	14' (4.3 m) Oxygen Tubing, Clear	Standard	50/cs
HCS4525	25' (7.6 m) Oxygen Tubing, Clear	Standard	25/cs
HCS4540FT	40' (1.0 m) Oxygen Tubing, Clear	Standard	10/cs
HCS455010	50' (15.2 m) Oxygen Tubing, Clear	Standard	10/cs
HCSU4507	7' (2.1 m) Oxygen Tubing, Clear	Universal	50/cs
HCSU4524	14' (4.3 m) Oxygen Tubing, Clear	Universal	50/cs
HCSU4525	25' (7.6 m) Oxygen Tubing, Clear	Universal	25/cs
HCSU455010	50' (15.2 m) Oxygen Tubing, Clear	Universal	50/cs



Green

Item No.	Description	Connector	Pkg
HCS4507G	7' (2.1 m) Oxygen Tubing, Green	Standard	50/cs
HCS4524G	14' (4.3 m) Oxygen Tubing, Green	Standard	50/cs
HCS4525G	25' (7.6 m) Oxygen Tubing, Green	Standard	25/cs
HCS4540FTG	40' (1.0 m) Oxygen Tubing, Green	Standard	25/cs
HCS4550GB	50' (15.2 m) Oxygen Tubing, Green	Standard	10/cs
HCSU4507G	7' (2.1 m) Oxygen Tubing, Green	Universal	50/cs
HCSU4524G	14' (4.3 m) Oxygen Tubing, Green	Universal	50/cs
HCSU4525G	25' (7.6 m) Oxygen Tubing, Green	Universal	25/cs
HCSU4550GB	50' (15.2 m) Oxygen Tubing, Green	Universal	10/cs



Universal Connector




Standard Connector

Corrugated Aerosol Tubing

Item No.	Description	Pkg
HSK600	Corrugated, Segmented Tubing, 100' (30.5 m), Opaque with cuttable cuff every 6" (15.2 cm)	1/cs
HSK626	Corrugated, Segmented Tubing, 6' (1.8 m), Opaque	50/cs
HCS4541	6" (15.2 cm) Corrugated Tubing, Blue	50/cs



HCS4541





Oxygen Masks

Choose from an industry-leading selection of masks for all your oxygen-delivery needs.

- » Soft, anatomical form for patient comfort
- » Crush-resistant tubing
- » A secure fit with elastic straps
- » Convenient choice of universal or standard connectors
- » Not made with natural rubber latex

Adult

Item No.	Description	Connector	Pkg
HCS4600B	Adult, Medium-Concentration Mask with 7' (2.1 m) Tubing	Standard	50/cs
HCSU4600B	Adult, Medium-Concentration Mask with 7' (2.1 m) Tubing	Universal	50/cs
HCS4610B	Adult, High-Concentration Mask with 7' (2.1 m) Tubing	Standard	50/cs
HCSU4610B	Adult, High-Concentration Mask with 7' (2.1 m) Tubing	Universal	50/cs
HCS4640B	Adult, Partial Non-Rebreather Mask with 7' (2.1 m) Tubing	Standard	50/cs
HCSU4640B	Adult, Partial Non-Rebreather Mask with 7' (2.1 m) Tubing	Universal	50/cs
HCS120E	Adult, Venturi Mask with 6 Jets and 7' (2.1 m) Tubing	Standard	50/cs
HCSU120E	Adult, Venturi Mask with 6 Jets and 7' (2.1 m) Tubing	Universal	50/cs
HCS121E	Adult, Venturi Mask with Single Dial and 7' (2.1 m) Tubing	Standard	50/cs
HCSU121E	Adult, Venturi Mask with Single Dial and 7' (2.1 m) Tubing	Universal	50/cs
HCS123E	Venturi Single Dial with Corrugated Tubing and 7" Tubing	Standard	50/cs
HCS4670B	Adult, Total Non-Rebreather Mask with 7' Tubing	Standard	*
HCSU4670B	Adult, Total Non-Rebreather Mask with 7' Tubing	Universal	*
HCS4680B	Adult, 3-in-1 Mask with 7' Tubing	Standard	*
HCSU4680B	Adult, 3-in-1 Mask with 7' Tubing	Universal	*



HCS4600B

Pediatric

Item No.	Description	Connector	Pkg
HCS4601B	Pediatric, Medium-Concentration Mask with 7' (2.1 m) Tubing	Standard	50/cs
HCSU4601B	Pediatric, Medium-Concentration Mask with 7' (2.1 m) Tubing	Universal	50/cs
HCS4642B	Pediatric, Partial Non-Rebreather Mask with 7' (2.1 m) Tubing	Standard	50/cs
HCSU4642B	Pediatric, Partial Non-Rebreather Mask with 7' (2.1 m) Tubing	Universal	50/cs
HCS4675B	Pediatric, Total Non-Rebreather Mask with 7' Tubing	Standard	*
HCSU4675B	Pediatric, Total Non-Rebreather Mask with 7' Tubing	Universal	*
HCS4685B	Pediatric, 3-in-1 Mask with 7' Tubing	Standard	*
HCSU4685B	Pediatric, 3-in-1 Mask with 7' Tubing	Universal	*



HCS4640B



Universal Connector



Standard Connector

*Coming soon. Call for availability.



HCS4630B



HCS4631B



HCS4630

Aerosol Masks

Extra patient comfort comes from an anatomical form and a soft, light and comfortable vinyl material.

- » A secure fit with elastic straps
- » Standard 22 mm connector
- » Not made with natural rubber latex

Adult

Item No.	Description	Pkg
HCS4630B	Adult Elongated Aerosol Mask with Straight Connector	50/cs
HCS4630	Adult Elongated Aerosol Mask with Swivel Connector	50/cs

Pediatric

HCS4631B	Pediatric Elongated Aerosol Mask with Swivel Connector	50/cs
----------	--	-------

Tent and Trach Masks

Adult

HCS4632B	Adult Face Tent with 22 mm Connector	50/cs
HCS245E	Adult Trach Mask	50/cs

Pediatric

HCS4621B	Pediatric Trach Mask	50/cs
----------	----------------------	-------



HCS4632B



HCS245E

Bubble Humidifiers

Choose from a wide variety of sizes and configurations.

Prefilled

Item No.	Description	Pkg
HCS00350	Prefilled Bubble Humidifier, 350 ml, with Adaptor	20/cs
HCS00550	Prefilled Bubble Humidifier, 550 ml, with Adaptor	12/cs
HCS00640	Prefilled Bubble Humidifier, 640 ml, with Adaptor	10/cs



Empty Bubble Humidifier

HCS64377	300 ml Empty, Disposable, Single Patient Use, 6 PSI Bubble Humidifier	50/cs
VTB6900	350 ml Empty, Disposable, Single Patient Use, Hi-Flow Bubble Humidifier	10/cs
HCS02KIT1	Oxygen Setup Kit	10/cs
HCS4514HFK	High Flow Kit: 1 HCS4514HF, 1 VTB6900	7/cs



Large Volume Nebulizer

Item No.	Description	Pkg
HCS03128	O28 Nebulizer Adaptor, F102 range 28-98%	50/cs
HCS04000	1,070 ml Sterile Water Bottle	10/cs
HCS04128	Prefilled Aerosol Humidifier, 1070 ml, with O28 Nebulizer Adaptor	10/cs



Accessories

Medline's complete selection includes a wide variety of everyday connectors and adaptors.

Item No.	Description	Pkg
HCS0440	EZ Wrap Cushions for Cannula	50/cs
HCS0590	Swivel Connector for Oxygen Tubing	25/cs
HCS1642	Pressure Line Adapter, 5-7 mm ID	50/cs
HCS2503	Disposable Water Trap	25/cs
HCS4555	Green Christmas Tree Adaptor with Swivel	50/cs
GT3102OXYHB4	Metal Nut and Nipple Adapter for oxygen	1 ea
HCS65100	Rigid Connector for Oxygen Tubing	200/cs
HSK962A	Circuit Connector, 22 OD/15 ID 22 OD	50/cs
HCS974	Aerosol Drainage Bag with "T" Adapter Hanger	50/cs
HCS974Y	Aerosol Drainage Bag with "Y" Adapter Hanger	50/cs
WST219	Aerosol Tee Adaptor	50/cs
DYNJ750	Oxygen Nose Clip	100/cs
AGISH749	Y Connector	5/ea
HCS676	15" Humidifier Connector, Tubing Adapter	50/cs
HCS1790	21" Humidifier Connector, Tubing Adapter	50/cs
HCS8986	Humidifier Connector with Rigid Black Elbow Connector	1 ea



Pulse Oximeters

Fingertip Pulse Oximeters	47
Handheld Oximeters	48
Pulse Oximeter Probes & Accessories	49

Finger Pulse Oximeters

Get quick and precise SpO₂ and pulse readings in seconds with our high-quality fingertip pulse oximeters. Small and lightweight, these clip sensors automatically take readings when a finger is inserted. They fit a wide range of finger sizes, have a long battery life and are solidly constructed for worry-free use. All models include AAA batteries, lanyard and warranty.



pulSTAT™ Finger Pulse Oximeter

Patented multi-directional display offers 6 viewing modes.

- » SpO₂, pulse rate, pulse bar, plethysmograph
- » Low power consumption provides 30 hours of use
- » Bright, digital OLED with adjustable readout sizes
- » Soft padding for patient comfort
- » Hinged battery cover
- » 2 year warranty

HCSM70D, 1/ea



Easy-Grip Finger Pulse Oximeter

Hinges adjust to comfortably accommodate larger fingers.

- » SpO₂, pulse rate, pulse bar, plethysmograph
- » Backlit LED display with 4 viewing modes
- » Slip-resistant grip
- » Hinged battery cover
- » 2 year warranty

HCSM70, 1/ea



Soft-Touch Finger Pulse Oximeter

Soft padding enhances patient comfort.

- » SpO₂, pulse rate, pulse bar, plethysmograph
- » Backlit LED display with 4 viewing modes
- » 1 year warranty

HCSM70C, 1/ea



Basic Finger Pulse Oximeter

Digital display lights up for easy reading in all lighting conditions.

- » SpO₂, pulse rate, pulse bar
- » Bright LCD screen
- » 1 year warranty

HCSM70B, 1/ea



High-Impact Finger Pulse Oximeter

Rugged casing withstands accidental drops.

- » SpO₂, pulse rate and pulse bar readings
- » Multi-directional display with full-color OLED screen
- » Durable silicone boot construction resists cracking and chipping
- » Hinged battery cover
- » Non Rx
- » 2 year warranty

HCSFSATSPORT, 1/ea



Pediatric Finger Pulse Oximeter

Child-friendly frog design makes kids smile.

- » SpO₂, pulse rate, pulse bar, plethysmograph
- » Patented multi-directional display has 6 viewing modes
- » Bright, digital OLED with adjustable readout sizes
- » Soft padding for patient comfort
- » Hinged battery cover
- » 2 year warranty

HCSM70P, 1/ea

Handheld Oximeters

Dependable and easy-to-use pulse oximeters with alarms for spot checks and continuous monitoring.

Medline MD300

- » SpO₂ and pulse rate
- » Adjustable display: 2 readout sizes and 10 brightness levels
- » Includes 1 adult finger probe
- » Automatic power off
- » Includes protective silicone cover



Item No.	Description	Pkg
HCSMD300	Handheld Spot Check Oximeter with 30" (70 cm) Adult Finger Probe	ea

Medline MD400

- » SpO₂, pulse rate, pulse bar, waveform
- » Adjustable audible and visual alarms
- » Two display modes
- » Comfortable finger clip probe
- » Data transfer to computer via USB cable
- » Low battery indicator



Item No.	Description	Pkg
HCSMD400	Handheld Continuous Economy Pulse Oximeter with Adult Finger Probe	ea

Replacement Probes

Item No.	Description	Pkg
HCSMD400A	Replacement Probe, Adult, for HCSMD400 and HCSMD300	ea
HCSMD400P	Replacement Probe, Pediatric, for HCSMD400 and HCSMD300	ea
HCSMD400I	Replacement Probe, Infant, for HCSMD400 and HCSMD300	ea





7 - Pin Nellcor



9 - Pin BCI

Pulse Oximeter Probes

Medline has probes for a wide variety of pulse oximeters, including handheld and tabletop models, such as the popular 7-pin Nellcor® and 9-pin BCI®.

- » Probe choices for unique requirements, materials and applications
- » Disposable probes in neonatal through adult sizes
- » Reusable probes in pediatric and adult sizes
- » Micro-Foam for easy sensor application and repositioning, Medplast for maximum adhesive for long-term use in one location, or Transpore for transparency with repositioning



Finger Clip
PCL0603120

Adult

Item No.	Category	Type	Compatible With	Pkg
PCL0603104	Reusable	Finger Clip	BCI 9-PIN	ea
PCL0202301	Disposable	Micro-Foam Style	BCI 9-PIN	24/cs
PCL0603120	Reusable	Finger Clip	Nellcor 7-PIN	ea
PCL1003120	Reusable	Soft Boot	Nellcor 7-PIN	ea
PCL0342121	Disposable	Medplast Style	Nellcor 7-PIN	24/cs
PCL0322101	Disposable	Transpore Style	Nellcor 7-PIN	24/cs
BCI3078	Reusable	Ear	BCI 9-PIN	ea



Medplast
PCL0342121

Pediatric

Item No.	Description	Type	Compatible With	Pkg
PCL0623104	Reusable	Finger Clip	BCI 9-PIN	ea
PCL1023104	Reusable	Soft Boot	BCI 9-PIN	ea
PCL0202302	Disposable	Micro-Foam Style	BCI 9-PIN	24/cs
PCL0322302	Disposable	Transpore Style	BCI 9-PIN	24/cs
PCL0623120	Reusable	Finger Clip	Nellcor 7-PIN	ea
PCL1023120	Reusable	Soft Boot	Nellcor 7-PIN	ea
PCL0202122	Disposable	Micro-Foam Style	Nellcor 7-PIN	24/cs
PCL0322122	Disposable	Transpore® Style	Nellcor 7-PIN	24/cs



Micro-Foam
PCL0202301

Infant

Item No.	Description	Type	Compatible With	Pkg
PCL0202303	Disposable	Micro-Foam Style	BCI 9-PIN	24/cs
PCL0202123	Disposable	Micro-Foam Style	Nellcor 7-PIN	24/cs



Transpore
PCL0322101

Neonatal

Item No.	Description	Type	Compatible With	Pkg
PCL0202304	Disposable	Micro-Foam Style	BCI 9-PIN	24/cs
PCL0342124	Disposable	Medplast Style	Nellcor 7-PIN	24/cs
PCL0202124	Disposable	Micro-Foam Style	Nellcor 7-PIN	24/cs
PCL0703120	Disposable	Multi Site Sensor	Nellcor 7-PIN	ea



Soft Boot
PCL1003120

Universal

Item No.	Description	Type	Compatible With	Pkg
PCL0703104	Reusable	Multi Site Sensor	BCI 9-PIN	ea

Resuscitation

Resuscitation Bags, Self Inflating	51-53
Mouth-to-Mask Resuscitator	54
Resuscitation Accessories	54
Hyperinflation Bags	55



Also available with manometer.

Manual Resuscitators

Innovations provide clinicians with improved ease of use. Available in infant, pediatric and adult sizes.

A Better Bag

- » Delivers up to 99% FIO₂
- » Maneuverable – The body of the bag swivels 360 degrees for multiple clinician hand positions, helping to reduce clinician fatigue
- » Translucent – The translucent, lightweight bag body and components allow for visualization of vomitus or occlusions in the system
- » Textured – Provides an enhanced, non-slip grip
- » Repositionable – The patient port swivels 360 degrees, allowing the clinician freedom to move around the patient without placing excessive torque or stress on the tracheal tube or mask
- » Manometer allows for a clear visual reading of the centimeters of water pressure being delivered to the patient

An Improved Mask

- » Firmly Sealing–Premium material provides a soft cushion and a firm seal

	Adult (Mask size #5)	Small Adult (Mask size #5)	Pediatric (Mask size #4)	Infant (Mask size #2)
Bag Volume	1,900 ml	1,500 ml	1,000 ml	300 ml
Stroke Volume	950-960 ml	640-660 ml	650-670 ml	140-170 ml
Dead Space with Mask	10 ml	10 ml	10 ml	10 ml
Inspiratory Resistance	<5 cm H ₂ O	<5 cm H ₂ O	<5 cm H ₂ O	<5 cm H ₂ O
Exhalation Resistance	<5 cm H ₂ O 50 LPM Flow	<5 cm H ₂ O 50 LPM Flow	<5 cm H ₂ O 50 LPM Flow	<5 cm H ₂ O 50 LPM Flow
Attainable Delivery Pressure	>80 cm H ₂ O	>80 cm H ₂ O	>80 cm H ₂ O	>80 cm H ₂ O
Ventilatory Frequency	>60 BPM (No Load)	>80 BPM (No Load)	>60 BPM (No Load)	>90 BPM (No Load)
Body Mass Use Range	>40 kg (88-lbs.)	>40 kg (88-lbs.)	<40 kg (88-lbs.)	<10 kg (22-lbs.)
Minimal Tidal Volume	600 ml	600 ml	500 ml	150 ml

Adult Manual Resuscitation Bags

Item No.	Description	Pkg
CPRM1116	Adult, Bag Reservoir, Size 5 Mask, 7' Oxygen Tubing	6
CPRM1116C	Adult, Bag Reservoir, Size 5 Mask, CO2 Detector, 7' Oxygen Tubing	6
CPRM1116F	Adult, Bag Reservoir, Size 5 Mask, Filter, 7' Oxygen Tubing	6
CPRM1116K	Adult, Bag Reservoir, Size 5 Mask, Flocap CO2 Detector, 7' Oxygen Tubing	6
CPRM1116M	Adult, Bag Reservoir, Size 5 Mask, Manometer, 7' Oxygen Tubing	6
CPRM1116MC	Adult, Bag Reservoir, Size 5 Mask, Manometer, CO2 Detector, 7' Oxygen Tubing	6
CPRM1116P	Adult, Bag Reservoir, Size 5 Mask, Peep, 7' Oxygen Tubing	6
CPRM1116PC	Adult, Bag Reservoir, Size 5 Mask, Peep, CO2 Detector, 7' Oxygen Tubing	6
CPRM1116PM	Adult, Bag Reservoir, Size 5 Mask, Peep, Manometer 7' Oxygen Tubing	6
CPRM1126	Adult, Tube Reservoir, Size 5 Mask, 7' Oxygen Tubing	6
CPRM1126C	Adult, Tube Reservoir, Size 5 Mask, CO2 Detector, 7' Oxygen Tubing	6
CPRM1126F	Adult, Tube Reservoir, Size 5 Mask, Filter, 7' Oxygen Tubing	6
CPRM1126K	Adult, Tube Reservoir, Size 5 Mask, Flocap CO2 Detector, 7' Oxygen Tubing	6
CPRM1126M	Adult, Tube Reservoir, Size 5 Mask, Manometer, 7' Oxygen Tubing	6
CPRM1126MC	Adult, Tube Reservoir, Size 5 Mask, Manometer, CO2 Detector, 7' Oxygen Tubing	6
CPRM1126P	Adult, Tube Reservoir, Size 5 Mask, Peep, 7' Oxygen Tubing	6

Small Adult Manual Resuscitation Bags

Item No.	Description	Pkg
CPRM4416	Small Adult, Bag Reservoir, Size 5 Mask, 7' Oxygen Tubing	6
CPRM4416M	Small Adult, Bag Reservoir, Size 5 Mask, Manometer, 7' Oxygen Tubing	6
CPRM4426	Small Adult, Tube Reservoir, Size 5 Mask, 7' Oxygen Tubing	6
CPRM4426M	Small Adult, Tube Reservoir, Size 5 Mask, Manometer, 7' Oxygen Tubing	6

Pediatric Manual Resuscitation Bags

Item No.	Description	Pkg
CPRM2216	Pediatric, Bag Reservoir, Size 4 Mask, 7' Oxygen Tubing	6
CPRM2216K	Pediatric, Bag Reservoir, Size 4 Mask, Flocap CO2 Detector, 7' Oxygen Tubing	6
CPRM2216P	Pediatric, Bag Reservoir, Size 4 Mask, Peep, 7' Oxygen Tubing	6
CPRM2216PC	Pediatric, Bag Reservoir, Size 4 Mask, Peep, CO2 Detector, 7' Oxygen Tubing	6
CPRM2216PC	Pediatric, Bag Reservoir, Size 4 Mask, Peep, CO2 Detector, 7' Oxygen Tubing	6
CPRM2216PM	Pediatric, Bag Reservoir, Size 4 Mask, Peep, Manometer, 7' Oxygen Tubing	6
CPRM2226	Pediatric, Tube Reservoir, Size 4 Mask, 7' Oxygen Tubing	6
CPRM2226C	Pediatric, Tube Reservoir, Size 4 Mask, CO2 Detector, 7' Oxygen Tubing	6
CPRM2226F	Pediatric, Tube Reservoir, Size 4 Mask, Filter, 7' Oxygen Tubing	6
CPRM2226K	Pediatric, Bag Reservoir, Size 4 Mask, Flocap CO2 Detector, 7' Oxygen Tubing	6
CPRM2226M	Pediatric, Tube Reservoir, Size 4 Mask, Manometer, 7' Oxygen Tubing	6
CPRM2226MC	Pediatric, Tube Reservoir, Size 4 Mask, Manometer, CO2 detector, 7' Oxygen Tubing	6
CPRM2226P	Pediatric, Tube Reservoir, Size 4 Mask, Peep, 7' Oxygen Tubing	6

Infant Manual Resuscitation Bags

Item No.	Description	Pkg
CPRM3312	Infant, Bag Reservoir, Size 2 Mask, 7' Oxygen Tubing	12
CPRM3312K	Infant, Bag Reservoir, Size 2 Mask, Flocap CO2 Detector, 7' Oxygen Tubing	12
CPRM3312P	Infant, Bag Reservoir, Size 2 Mask, Peep, 7' Oxygen Tubing	12
CPRM3312P	Infant, Bag Reservoir, Size 2 Mask, Peep, 7' Oxygen Tubing	12
CPRM3322	Infant, Tube Reservoir, Size 2 Mask, 7' Oxygen Tubing	12
CPRM3322C	Infant, Tube Reservoir, Size 2 Mask, CO2 Detector, 7' Oxygen Tubing	12
CPRM3322M	Infant, Tube Reservoir, Size 2 Mask, Manometer, 7' Oxygen Tubing	12
CPRM3322MC	Infant, Tube Reservoir, Size 2 Mask, Manometer, CO2 Detector, 7' Oxygen Tubing	12
CPRM3322P	Infant, Tube Reservoir, Size 2 Mask, Peep, 7' Oxygen Tubing	12
CPRM3322PM	Infant, Tube Reservoir, Size 2 Mask, Peep, Manometer, 7' Oxygen Tubing	12

If you do not see your preferred configuration, Medline can customize a resuscitator bag to your preference.

Manual Resuscitator Ordering Guide

Item No.	Bag Size				Mask Size	Reservoir Style		Tubing Length		Pressure Relief Valve	CO2 Detector	On the Exhalation Port		Mano-meter	Pkg
	Infant	Pediatric	SM Adult	Adult		Bag	40" Tubing	7'	14'			Peep Valve	Filter		
CPRM1116				•	5	•		•							6/cs
CPRM1116C				•	5	•		•			•				6/cs
CPRM1116F				•	5	•		•					•		6/cs
CPRM1116FC				•	5	•		•			•		•		6/cs
CPRM1116FM				•	5	•		•					•	•	6/cs
CPRM1116FPC				•	5	•		•			•	•	•		6/cs
CPRM1116FPM				•	5	•		•				•	•	•	6/cs
CPRM1116L				•	5	•		•							6/cs
CPRM1116M				•	5	•		•						•	6/cs
CPRM1116MC				•	5	•		•			•			•	6/cs
CPRM1116P				•	5	•		•				•			6/cs
CPRM1116PC				•	5	•		•			•	•			6/cs
CPRM1116PM				•	5	•		•				•		•	6/cs
CPRM1116PMT				•	5	•			•			•		•	6/cs
CPRM1116T				•	5	•			•						6/cs
CPRM1126				•	5		•	•							6/cs
CPRM1126C				•	5		•	•			•				6/cs
CPRM1126F				•	5		•	•					•		6/cs
CPRM1126FL				•	5		•	•					•		6/cs
CPRM1126FM				•	5		•	•					•	•	6/cs
CPRM1126FP				•	5		•	•				•	•		6/cs
CPRM1126FPC				•	5		•	•			•	•	•		6/cs
CPRM1126M				•	5		•	•						•	6/cs
CPRM1126MC				•	5		•	•			•			•	6/cs
CPRM1126P				•	5		•	•				•			6/cs
CPRM1126PC				•	5		•	•			•	•			6/cs
CPRM2216		•			4	•		•		40 CM					6/cs
CPRM2216C		•			4	•		•		40 CM	•				6/cs
CPRM2216F		•			4	•		•		40 CM			•		6/cs
CPRM2216FM		•			4	•		•		40 CM			•	•	6/cs
CPRM2216FP		•			4	•		•		40 CM		•	•		6/cs
CPRM2216M		•			4	•		•		40 CM				•	6/cs
CPRM2216P		•			4	•		•		40 CM		•			6/cs
CPRM2216PC		•			4	•		•		40 CM	•	•			6/cs
CPRM2216PCM		•			4	•		•		40 CM		•		•	6/cs
CPRM2216PM		•			4	•		•		40 CM		•		•	6/cs
CPRM2216PMT		•			4	•			•	40 CM		•		•	6/cs
CPRM2226		•			4		•	•		40 CM					6/cs
CPRM2226C		•			4		•	•		40 CM	•				6/cs
CPRM2226F		•			4		•	•		40 CM			•		6/cs
CPRM2226M		•			4		•	•		40 CM				•	12/cs
CPRM2226MC		•			4		•	•		40 CM	•			•	12/cs
CPRM2226P		•			4		•	•		40 CM		•			12/cs
CPRM3312	•				2	•		•		40 CM					12/cs
CPRM3312C	•				2	•		•		40 CM	•				12/cs
CPRM3312F	•				2	•		•		40 CM			•		12/cs
CPRM3312M	•				2	•		•		40 CM				•	12/cs
CPRM3312M3	•				2	•		•		40 CM				•	12/cs
CPRM3312P	•				2	•		•		40 CM		•			12/cs
CPRM3312PM	•				2	•		•		40 CM		•		•	12/cs
CPRM3312PMC2	•				1&2	•		•		40 CM	•	•		•	12/cs
CPRM3322	•				2		•	•		40 CM					12/cs
CPRM3322C	•				2		•	•		40 CM	•				12/cs
CPRM3322F	•				2		•	•		40 CM			•		12/cs
CPRM3322FPC	•				2		•	•		40 CM	•	•	•		12/cs
CPRM3322M	•				2		•	•		40 CM				•	12/cs
CPRM3322MC	•				2		•	•		40 CM	•			•	12/cs
CPRM3322P	•				2		•	•		40 CM		•			12/cs
CPRM3322PM	•				2		•	•		40 CM		•		•	12/cs
CPRM4416			•		5	•		•							12/cs
CPRM4416M			•		5	•		•						•	6/cs
CPRM4426			•		5		•	•							6/cs
CPRM4426M			•		5		•	•						•	6/cs

If you do not see your preferred configuration, Medline can customize a resuscitator bag to your preference.

Resuscitation Masks and Accessories

Mouth-To-Mask Resuscitator Mask

This ready-to-use CPR mask is a highly effective one-way filtering system that shields the rescuer and victim from cross-contamination.

- » A clear plastic design allows for easy placement, and the soft inflated cuff improves the fit and patient comfort
- » Standard 22 mm fittings are compatible with any resuscitator



HCS64184

Item No.	Description	Pkg
HCS64184	Mouth-to-mask Resuscitator with Valve	12/cs

Medline offers a wide variety of masks which are not shown here. Please contact us for more information.

Resuscitation Masks

Medline's resuscitation masks provide an ergonomic cushion design for an optimal seal and clear cone to visualize patient breathing.

- » Removable hook ring for hand-held or head-strap application (hook rings included on sizes 3-6)
- » No valve prevents mask deflation
- » Available in all sizes, including preemie



Item No.	Description	Pkg
DYNJAAMASK11	Cushion Mask, Neonatal, Size 1	40/cs
DYNJAAMASK12	Cushion Mask, Infant, Size 2	40/cs
DYNJAAMASK13	Cushion Mask, Child, Size 3	40/cs
DYNJAAMASK14	Cushion Mask, Small Adult, Size 4	40/cs
DYNJAAMASK15	Cushion Mask, Medium Adult, Size 5	20/cs
DYNJAAMASK16	Cushion Mask, Large Adult, Size 6	20/cs

Medline offers a wide variety of masks which are not shown here. Please contact us for more information.

Resuscitation Accessories

Medline's high quality resuscitation accessories are designed to work with all Medline resuscitation devices and allow you to customize to patient need.

Item No.	Description	Pkg
CPRMP19	19 MM Adjustable Peep Valve	20/cs
CPRMP19MRI	19 MM Adjustable Peep Valve, Mri Compatible	100/cs
CPRMP30	19 MM Adjustable Peep Valve With 30 mm Adapter	20/cs
CPRMFILTER	Expiratory Filter	50/cs
CERFLOCAP	Flowcap Co2 Detector With Flow Indicator	24/cs
CPRMCO2	CO2 Detector	10/cs

CPRM19



CPRMFILTER



CERFLOCAP

Hyperinflation Systems

Help prevent or treat lung collapse with a hyperinflation system.


- » Bag design uses a “dipped” process to enhance tactile sensitivity and improve control
- » Uniform inflation and the highest level of performance
- » 100% compliant with ISO 5362 standards
- » Variety of sizes to more accurately fit every patient



Hyperinflation Systems

7' (2.13 m) of O₂ tubing included*

Item No.	Breathing Bag	O ₂ Connector	Mask	Manometer	Pkg.
DYNJAAHY005	0.5 L	Standard	N/A		20/cs
DYNJAAHY600	0.5 L w/Tail-End Bleed Valve	Universal	N/A		50/cs
DYNJAAHY30	0.5 L	Standard	1 Neonatal		20/cs
DYNJAAHY31	0.5 L	Standard	2 Infant		20/cs
DYNJAAHY29	0.5 L	Universal	2 Infant (Pre-Inflated)		20/cs
DYNJAAHY01	1 L	Standard	N/A		20/cs
DYNJAAHY601	1 L w/Tail-End Bleed Valve	Universal	N/A		50/cs
DYNJAAHY01M2	1 L	Universal	2 Infant		20/cs
DYNJAAHY02	2 L	Standard	N/A		20/cs
DYNJAAHY02M3	2 L	Universal	3 Child		20/cs
DYNJAAHY03	3 L	Standard	N/A		20/cs
DYNJAAHY03M5	3 L	Universal	5 Adult		20/cs
DYNJAAHY5005	0.5 L	Standard	N/A	●	50/cs
DYNJAAHY5053P	0.5 L	Universal	N/A	●	50/cs
DYNJAAHY32	0.5 L	Standard	1 Neonatal	●	50/cs
DYNJAAHY33	0.5 L	Standard	2 Infant	●	50/cs
DYNJAAHY35	0.5 L	Universal	2 Infant	●	50/cs
DYNJAAHY34	0.5 L	Universal	2 Infant (Pre-Inflated)	●	50/cs
DYNJAAHY501	1 L	Universal	N/A	●	50/cs
DYNJAAHY501M2	1 L	Universal	2 Infant	●	50/cs
DYNJAAHY502	2 L	Universal	N/A	●	50/cs
DYNJAAHY502M3	2 L	Universal	3 Child	●	50/cs
DYNJAAHY503	3 L	Universal	N/A	●	50/cs
DYNJAAHY503M5	3 L	Universal	5 Adult	●	50/cs



*Custom options are available

Suction Catheters and Trach Care

Closed Suction	57-58
Open Suction Catheters	59
Open Suction Catheters, Sleeved	59
Open Suction Kits	60
Open Suction Kits with Solutions	61
Sterile Solutions	61
Trach Care & Cleaning Trays	62
Mucus/Specimen Traps	63

ClearPro™ Closed Suction Catheter System

The latest addition to our line of closed suction catheters, ClearPro features a more advanced cleaning chamber and stronger sheath than competitors for greater performance.

- » Labeled for facilities to determine change out policy
- » DEHP-free
- » Isolated cleaning chamber for rapid catheter irrigation
- » Pressure wiper seal for PEEP maintenance and catheter decontamination
- » Color-coded connectors by French size
- » Kitted with all necessary adapters and connectors to make SKU selection simple
- » Catheter depth markings on all sizes
- » High-quality translucent sleeve to reduce cross contamination and improve ability to advance the catheter
- » One-way needlessly irrigation port perfect for saline vials

Endotracheal

Item Number	Description	Pkg.	Thumb Valve Color
DYNCPDS10	ClearPro, Double Swivel Elbow, 10 French	20/cs	Black
DYNCPDS12	ClearPro, Double Swivel Elbow, 12 French	20/cs	White
DYNCPDS14	ClearPro, Double Swivel Elbow, 14 French	20/cs	Green
DYNCPDS16	ClearPro, Double Swivel Elbow, 16 French	20/cs	Orange
DYNCPTP10	ClearPro, T-Piece Adapter, 10 French	20/cs	Black
DYNCPTP12	ClearPro, T-Piece Adapter, 12 French	20/cs	White
DYNCPTP14	ClearPro, T-Piece Adapter, 14 French	20/cs	Green
DYNCPTP16	ClearPro, T-Piece Adapter, 16 French	20/cs	Orange

Tracheostomy

Item Number	Description	Pkg.	Thumb Valve Color
DYNCPDS10T	ClearPro, Double Swivel Elbow, 10 French	20/cs	Black
DYNCPDS12T	ClearPro, Double Swivel Elbow, 12 French	20/cs	White
DYNCPDS14T	ClearPro, Double Swivel Elbow, 14 French	20/cs	Green
DYNCPDS16T	ClearPro, Double Swivel Elbow, 16 French	20/cs	Orange
DYNCPTP10T	ClearPro, T-Piece Adapter, 10 French	20/cs	Black
DYNCPTP12T	ClearPro, T-Piece Adapter, 12 French	20/cs	White
DYNCPTP14T	ClearPro, T-Piece Adapter, 14 French	20/cs	Green
DYNCPTP16T	ClearPro, T-Piece Adapter, 16 French	20/cs	Orange



DYNCPTP14

Also available with Double Swivel Elbow

Medline Closed Suction Catheter Systems

Adult

- » Labeled for facilities to determine change out policy
- » On/off valve allows suction catheter tip to be isolated and creates an enclosed cleaning chamber
- » Soft tip helps reduce damage to mucous membranes
- » Transparent connector for easy observation
- » High-quality, translucent EVA sheath helps reduce cross-contamination and facilitates operation when advancing the catheter

Endotracheal

Item Number	Description	Pkg.	Thumb Valve Color
DYNCSDS12	Adult, Double Swivel Elbow, 12 Fr	10/CS	White
DYNCSDS14	Adult, Double Swivel Elbow, 14 Fr	10/CS	Green
DYNCSTP12	Adult, T-Piece Adapter, 12 Fr	10/CS	White
DYNCSTP14	Adult, T-Piece Adapter, 14 French	10/CS	Green

Tracheostomy

Item Number	Description	Pkg.	Thumb Valve Color
DYNCSDS12T	Adult, Double Swivel Elbow, 12 Fr	10/CS	White
DYNCSDS14T	Adult, Double Swivel Elbow, 14 Fr	10/CS	Green
DYNCSTP12T	Adult, T-Piece Adapter, 12 Fr	10/CS	White
DYNCSTP14T	Adult, T-Piece Adapter, 14 Fr	10/CS	Green

Pediatric

- » Available with an elbow or Y-style modular patient adapter
- » Three adapter sizes to ensure compatibility
- » Labeled for facilities to determine change out policy
- » Compatible with all modes of ventilation
- » On/off valve allows suction catheter tip to be isolated and creates an enclosed cleaning chamber

Elbow Adapter

Item Number	Description	Pkg.	Thumb Valve Color
DYNCPE5	Pediatric CSC, Elbow Adapter, 5 Fr	10/CS	Black
DYNCPE6	Pediatric CSC, Elbow Adapter, 6 Fr	10/CS	Light Green
DYNCPE7	Pediatric CSC, Elbow Adapter, 7 Fr	10/CS	Purple
DYNCPE8	Pediatric CSC, Elbow Adapter, 8 Fr	10/CS	Light Blue

Y-Connector

Item Number	Description	Pkg.	Thumb Valve Color
DYNCPY5	Pediatric CSC, Y-Connector, 5 Fr	10/CS	Black
DYNCPY6	Pediatric CSC, Y-Connector, 6 Fr	10/CS	Light Green
DYNCPY7	Pediatric CSC, Y-Connector, 7 Fr	10/CS	Purple
DYNCPY8	Pediatric CSC, Y-Connector, 8 Fr	10/CS	Light Blue



Open Suction Catheters

Medline's D419 series catheters offers a full range of straight packaging.

- » Available in a complete range of continuous adult and pediatric sizes
- » All Medline suction catheters feature easy-to-read depth markings
- » DeLee tip is flat, soft and smooth to minimize airway trauma for pediatric patients
- » Angled, soft and smooth whistle tip maximizes suctioning capability for adult patients, while minimizing trauma potential

Single Sterile-Straight Pack

Item No.	Description	Pkg
DYND41906	Suction Catheter, Single Sterile, 6 Fr	100/cs
DYND41908	Suction Catheter, Single Sterile, 8 Fr	100/cs
DYND41900	Suction Catheter, Single Sterile, 10 Fr	100/cs
DYND41901	Suction Catheter, Single Sterile, 12 Fr	100/cs
DYND41902	Suction Catheter, Single Sterile, 14 Fr	100/cs
DYND41903	Suction Catheter, Single Sterile, 16 Fr	100/cs
DYND41904	Suction Catheter, Single Sterile, 18 Fr	100/cs

Coiled

Item No.	Description	Pkg
DYND41906C	Suction Catheter, DeLee, 6 Fr	100/cs
DYND41908C	Suction Catheter, DeLee, 8 Fr	100/cs
DYND41900C	Suction Catheter, Whistle, 10 Fr	100/cs
DYND41901C	Suction Catheter, Whistle, 12 Fr	100/cs
DYND41902C	Suction Catheter, Whistle, 14 Fr	100/cs
DYND41903C	Suction Catheter, Whistle, 16 Fr	100/cs
DYND41904C	Suction Catheter, Whistle, 18 Fr	100/cs

Poly-Cath

DYND40990	Suction Catheter, Poly-Cath, 10 Fr	100/cs
DYND40992	Suction Catheter, Poly-Cath, 14 Fr	100/cs

Continuous Valve


DYND41960	Suction Catheter, Continuous Valve, 10 Fr	100/cs
DYND41962	Suction Catheter, Continuous Valve, 14 Fr	100/cs
DYND41964	Suction Catheter, Continuous Valve, 18 Fr	100/cs

Sleeved Catheters

Our D407 series open-suction catheters are contained in a sleeve to help avoid cross-contamination.

Item No.	Description	Pkg
DYND40708F	Sleeved Suction Catheter: 8 Fr	50/cs
DYND40700F	Sleeved Suction Catheter: 10 Fr	50/cs
DYND40701F	Sleeved Suction Catheter: 12 Fr	50/cs
DYND40702F	Sleeved Suction Catheter: 14 Fr	50/cs

DYND41901



DYND41962



DYND40708F





DYND40982



DYND40973



DYND48982



Kits and Trays

Medline's D409 series kits and trays provide all the necessities for airway clearance.

- » Irrigation solutions can be poured into the tray or the pop-up cup
- » All kits and trays contain Medline Aloetouch® 3G gloves
- » Catheters have easy-to-read depth markings
- » DeLee tip is flat, soft and smooth to minimize airway trauma for pediatric patients
- » Angled, soft and smooth whistle tip maximizes suctioning capability for adult patients, while minimizing trauma potential

Suction Catheter Trays

Item No.	Description	Pkg
DYND40986	Suction Catheter Mini Tray with Gloves: 6 Fr	100/cs
DYND40988	Suction Catheter Mini Tray with Gloves: 8 Fr	100/cs
DYND40980	Suction Catheter Mini Tray with Gloves: 10 Fr	100/cs
DYND40981	Suction Catheter Mini Tray with Gloves: 12 Fr	100/cs
DYND40982	Suction Catheter Mini Tray with Gloves: 14 Fr	100/cs
DYND40983	Suction Catheter Mini Tray with Gloves: 16 Fr	100/cs
DYND40984	Suction Catheter Mini Tray with Gloves: 18 Fr	100/cs

Suction Catheter Kits

Item No.	Description	Pkg
DYND40976	Suction Catheter Pack with gloves and Pop-up cup, 6 Fr	50/cs
DYND40978	Suction Catheter Pack with gloves and Pop-up cup, 8 Fr	50/cs
DYND40970	Suction Catheter Pack with gloves and Pop-up cup, 10 Fr	50/cs
DYND40971	Suction Catheter Pack with gloves and Pop-up cup, 12 Fr	50/cs
DYND40972	Suction Catheter Pack with gloves and Pop-up cup, 14 Fr	50/cs
DYND40973	Suction Catheter Pack with gloves and Pop-up cup, 16 Fr	50/cs
DYND40974	Suction Catheter Pack with gloves and Pop-up cup, 18 Fr	50/cs

Suction Catheter Trays, Poly-cath

Item No.	Description	Pkg
DYND48980A	Suction Catheter Mini Tray with Gloves, Poly-Cath: 10 Fr	100/cs
DYND48982	Suction Catheter Mini Tray with Gloves, Poly-Cath: 14 Fr	100/cs

MEDLINE MANUFACTURES
OVER 125,000
 sterile procedure kits *every day.*

Suction Catheter Kits with Sterile Solutions

Our D414 series suction catheters come with everything needed for airway clearance, including a choice of irrigation fluids.

- » All kits and trays contain Medline Aloetouch 3G gloves
- » Catheters have easy-to-read depth markings
- » DeLee tip is flat, soft and smooth to minimize airway trauma for pediatric patients
- » Angled, soft and smooth whistle tip maximizes suctioning capability for adult patients, while minimizing trauma potential



DYND41472



DYND40442



Water

Item No.	Description	Pkg
DYND41476	Suction Catheter Tray with Sterile Water, 100 ml; 6 Fr	36/cs
DYND41478	Suction Catheter Tray with Sterile Water, 100 ml; 8 Fr	36/cs
DYND41470	Suction Catheter Tray with Sterile Water, 100 ml; 10 Fr	36/cs
DYND41471	Suction Catheter Tray with Sterile Water, 100 ml; 12 Fr	36/cs
DYND41472	Suction Catheter Tray with Sterile Water, 100 ml; 14 Fr	36/cs
DYND41474	Suction Catheter Tray with Sterile Water, 100 ml; 18 Fr	36/cs

Saline

Item No.	Description	Pkg
DYND41478SA	Suction Catheter Tray with Sterile Saline, 100 ml; 8 Fr	36/cs
DYND41470SA	Suction Catheter Tray with Sterile Saline, 100 ml; 10 Fr	36/cs
DYND41442	Suction Catheter Tray with Sterile Saline, 100 ml; 14 Fr	36/cs
DYND40442	Suction Catheter Tray with Saline and Drape: 14 Fr	36/cs

Sterile Solutions

These convenient, ready-to-use containers reduce the contamination risks and potential waste associated with bulk containers of sterile water and saline.

- » Ideal for device irrigation, suction, dressing moistening and wound debridement
- » Cup size provides the proper amount for suctioning

Water

Item No.	Description	Pkg
RD130295	Sterile Water, 100 ml	48/cs
DYND40570	Sterile Water, 110 ml	48/cs
PCS1550	Sterile Water, 250 ml	24/cs
DYND500MLW	Sterile Water, 500 ml	18/cs
DYND1000MLW	Sterile Water, 1000 ml	6/cs

Saline

Item No.	Description	Pkg
RD130296	Sterile Saline, 0.9%, 100 ml	48/cs
DYND40540	Sterile Saline, 0.9%, 110 ml	48/cs
PCS1650	Sterile Saline, 0.9%, 250 ml	24/cs
DYND500MLS	Sterile Saline, 0.9%, 500 ml	18/cs
DYND1000MLS	Sterile Saline, 0.9%, 1000 ml	6/cs



DYND1000MLS

PCS1650

RD130296

DYND500MLS

Tracheostomy Care and Cleaning Trays

Medline offers a wide variety of tracheostomy care and cleaning trays in standard and customized configurations to meet diverse clinical needs.



- » Available in economy, basic and full-compartment kits
- » May be customized



	DYND40582	DYND40622	DYND40610	DYND40620	DYND40589	DYND40615	DYND40580	DYND40650	DYND40511	DYND40594	CC3T3030A	CC3T4691A	CC4682A	CC4681A
Tray Compartments	3	1	1	1	3	3	3	3	3	3	1	3	3	3
AloeTouch Gloves	2		2	2	2	2	2	2	2	2		2		2
Drape	1		1	1	1	1	1		1	1		1	1	1
Pipe Cleaners		2	2						3	3	3	3	3	3
Non-woven Gauze 4 x 4	4	2	2	2	2	2	2	2	1	3	4	1	3	3
Cotton Tip Applicator	2	2	5	5	2	5	2	2	4	2	2	2	3	3
Pliable Cleaning Brush	1	1	1		1	1	1		1	1	1	1	1	1
Twill Tape 34"	1	1	1	1	1	1	1	1	1		1	1	1	
6-ply trach Dressing	1		1	1	1	1	1	1	1	1	1	1	1	1
Hydrogen Peroxide					2	2		1		1				
Saline 100ml					1				1	1				
Water 100ml							1							
Suction Catheter 14 fr							1		1					
Forceps													1	1
Pkg	20/cs	30/cs	20/cs	20/cs	20/cs	20/cs	20/cs	20/cs	20/cs	20/cs	30/cs	20/cs	20/cs	20/cs

Did You Know?

We can build a custom kit to your exact specifications. Medline is a leading manufacturer of custom procedure kits.



Mucus/Specimen Traps

Our mucus and specimen traps have leak-proof transport caps for more secure specimen collection.

- » Sterile plastic vials, graduated in 1cc increments
- » Includes threaded cap, tubing not made from natural rubber latex and patient ID label

Sterile DeLee Mucus Trap

Item No.	Description	Pkg
DYND44108	Trap with catheter and cap, capacity 20cc: 8 Fr	50/cs
DYND44110	Trap with catheter and cap, capacity 20cc: 10 Fr	50/cs

Sterile Mucus Specimen Trap

Item No.	Description	Pkg
DYND44140	Trap With Transport Cap, Capacity: 40 cc	50/cs
DYND44180	Trap With Transport Cap, Capacity: 80 cc	50/cs

Top Respiratory Codes for Seasonal Flu Demand



Item No.	Description
Cannula	
HCSU4518	Cannula, Pediatric, Soft Touch, 7' Tbg with Universal Connector
HCS4518	Cannula, Pediatric, Soft Touch, 7' Tbg with Standard Connector
HCS4514	Cannula, Adult, Soft Touch, 7' Tbg with Standard Connector
HCSU4514	Cannula, Adult, Soft Touch, 7' Tbg with Universal Connector
HCS4516B	Cannula, Adult, Soft Touch, 14' Tbg with Standard Connector
HCSU4516B	Cannula, Adult, Soft Touch, 14' Tbg with Universal Connector
HCSS4514S	Cannula, Adult, Super Soft, 7' Tbg with Standard Connector
HCSU4514S	Cannula, Adult, Super Soft, 7' Tbg with Universal Connector
HCSU4516S	Cannula, Adult, Super Soft, 14' Tbg with Universal Connector
HCSU4518S	Cannula, Pediatric, Super Soft, 7' Tbg with Universal Connector
HCS4514HF	Cannula, Adult, Soft Touch, High Flow, 7' Tbg with Standard Connector



Item No.	Description
Mask	
HCS4600B	Mask, Adult, Medium Concentration, 7' Tbg with Standard Connector
HCSU4600B	Mask, Adult, Medium Concentration, 7' Tbg with Universal Connector
HCS4640B	Mask, Adult, Non-Rebreather Standard, 7' Tbg with Standard Connector
HCSU4640B	Mask, Adult, Non-Rebreather Standard, 7' Tbg with Universal Connector
HCS4601B	Mask, Pediatric, Medium Concentration, 7' Tbg with Standard Connector
HCSU4601B	Mask, Pediatric, Medium Concentration, 7' Tbg with Universal Connector
HCS4642B	Mask, Pediatric, Non-Rebreather Standard, 7' Tbg with Standard Connector
HCSU4642B	Mask, Pediatric, Non-Rebreather Standard, 7' Tbg with Universal Connector
HCS4630B	Mask, Adult, Elongated Aerosol Mask
HCS4631B	Mask, Pediatric, Elongated Aerosol Mask



Item No.	Description
Nebulizer	
HCS4485	Nebulizer Kit, Mask, Adult, 7' Tbg with Standard Connector
HCSU4485	Nebulizer Kit, Mask, Adult, 7' Tbg with Universal Connector
HCS4483	Nebulizer Kit, Tee, Mouthpiece, 6" Reservoir, 7' Tbg with Standard Connector
HCSU4483	Nebulizer Kit, Tee, Mouthpiece, 6" Reservoir, 7' Tbg with Universal Connector
HCS4482	Nebulizer Kit, Tee, Mouthpiece, 7' Tbg with Standard Connector
HCSU4482	Nebulizer Kit, Tee, Mouthpiece, 7' Tbg with Universal Connector
HCS4486	Nebulizer Kit, Mask, Pediatric, 7' Tbg with Standard Connector
HCSU4486	Nebulizer Kit, Mask, Pediatric, 7' Tbg with Universal Connector
HCS70004	Nebulizer Compressor



Item No.	Description
Incentive Spirometer	
HCS93250	Incentive Spirometer, 2500ml
HCS93500	Incentive Spirometer, 5000ml



MDI Spacer	
HCSA100	MDI Spacer
HCSA112	MDI Spacer w Infant Mask
HCSA113	MDI Spacer w Pediatric Mask
HCSA114	MDI Spacer w Adult Mask



Bacterial/Viral Filter	
DYNJAABV1	Bacterial/Viral Filter
DYNJAAHME5	Large HME, filtered, no port
DYNJAAHME12	Large HME, non-filtered, no port



Prefilled Bubble Humidifier	
HCS00350	Prefilled Bubble Humidifier, 350ml
HCS00550	Prefilled Bubble Humidifier, 550ml

Retail packaging available on a variety of products, ask your Medline Representative for more details.





Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US    

REFERENCES 1. Black JM, Cuddigan JE, Walko MA, Diddier LA, Lander MJ, Kelpel MR. Medical device related pressure ulcers in hospitalized patients. *International Wound Journal*. 2010; 7(5):358-365. 2. American Lung Association website. Measuring Your Peak Flow Rate. Available at: <http://www.lung.org/lung-health-and-diseases/lung-disease-lookup/asthma/living-with-asthma/managing-asthma/measuring-your-peak-flow-rate.html>. Accessed September 27, 2016.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure. © 2016 Medline Industries, Inc. Aeromist Compact, Aeromist Buddies, and Aeromist Colors are trademarks, and Aloetouch, DoggyNeb, and Medline are registered trademarks of Medline Industries, Inc. Aeroflex and Phoneris are registered trademarks of Lazarus Medical, LLC, Airgel, iQ Blue, MiniMe and Mojo are registered trademarks of Sleepnet Corporation. Oxequip is a registered trademark of Allied Healthcare Products, Inc. Citrus II is a registered trademark of Beaumont Products, Inc. Eclipse is a registered trademark of SeQual Technologies Inc. Gentec is a registered trademark of Genstar Technologies Co., Inc. Nellcor is a registered trademark of Covidien LP Corporation. Ohio is a registered trademark of Ohio Medical Corporation. Ohmeda is a registered trademark of GE Healthcare. QuietRide is a registered trademark of The Mail Store, Inc. Trachvox is a registered trademark of Engineered Medical Systems, Inc. Transpore is a registered trademark of 3M Corporation. MKT1792025 / LIT085 / 2.7M / SG / 15