

Safe n Simple™
Clean. Comfortable. Confident.



Ostomy Product Guide

Introduction...

Safe n Simple™ is a medical care innovator, specializing in the promotion of healthy skin, complete security and confident ostomy and wound care management. We view our product development through the eyes of our customers, who want active, productive lives unimpeded by their conditions. Our goal is to make our customer's ostomy and wound care safe and simple.

We follow two overriding principles: create innovative, high-performance products, and deliver the greatest possible value. We continually push our technicians to develop better ideas, more advanced formulations and superior products, while at the same time making those products convenient and affordable.

We've spent years in **Ostomy care** learning from our customers what they really want. The result is a complete line of ostomate care products designed to *Cleanse, Protect, Secure, and Assure.*

And, a complete line of wound care products created to promote healing and promote cost savings.

Above all, our customers want safety and simplicity. With Safe n Simple, that's our promise.

Table of Contents

Ostomy Care

Page 1	Introduction
Page 2-3	Ostomy & Wound Care - Table of Contents
Page 4	Peri-Stoma Cleanser & Adhesive Remover Adhesive Remover
Page 5	No-Sting Skin Barrier Film No-Sting Skin Barrier Wand No-Sting Skin Barrier Spray
Page 6	Skin Barrier Film Ostomy Skin Barrier Powder
Page 7	Skin Barrier Ring Conforming Seal Integrity Skin Barrier Seal
Page 8	Skin Barrier Arcs 1/2 Arc Skin Barrier Strips Skin Barrier Sheet
Page 9	No-Sting Skin Barrier Paste Skin Barrier Paste Skin Barrier Paste-Fast Drying
Page 10.....	Assure C Odor Eliminator Ostomy Pouch Deodorant
Page 11.....	PouchesStoma Caps
Page 12-14	Security Support Belts
Page 15.....	DermaPro Waterproof Silicone Tape DermaPro Waterproof Plastic Tape

Special Information

Page 16..... Ostomy Accessory Allowables

Page 17-18 Stay Connected with SNS



Safe n Simple provides focused care that you can trust in. Our customer service staff are dedicated specialists who are trained to support you with product recommendations to help you improve your quality of life.

We are proud to offer our SimpleStart Ostomy Care Kits to give you an opportunity to trial our ostomy accessory product line.

Call our toll free customer service line to request a kit:

844-767-6334

Email Us: simplestart@sns-medical.com

or fax us at: 248-575-4172

Peri-Stoma Cleanser & Adhesive Remover

**ALCOHOL
FREE!**



Use these soft wipes to breakdown and remove adhesives, dirt and oils from peri-stomal skin. Does not leave an oily residue. No-sting, alcohol-free, and no-rinse formulation. Gentle, convenient and scent-free. Convenient re-sealable pouch.



Item	Product#	Size	HCPC
Large Wipes, 50 Pack	SNS00525	5"x7"	A4456
Large Individual Wipes 5 Box	SNS00505	5"x7"	A4456
Large Individual Wipes 50 Box	SNS00550	5"x7"	A4456
Large Individual Wipes 75 Box	SNS00575	5"x7"	A4456

Customer Testimonial:

"I had some red and irritated skin using other wipes and when I received the Safe n Simple wipes from my medical supply co., they gave me instant relief and the skin healed. I like also, that they are big and so soft. I have nothing but praise for them and I will always use them, nothing compares to them." Sincerely, Hildegard.

Adhesive Remover



If a more aggressive adhesive remover is needed, then this is the wipe to use. The alcohol based wipe has a very gentle formula to clean and remove debris from the skin. Rinse is required.

Item	Product#	Size	HCPC
Individual Wipes 50 Box	SNS00650	2"x2"	A4456

Protect

No-Sting Skin Barrier Film

ALCOHOL FREE!



This wipe forms a barrier to protect the skin from adhesives, friction and body fluids. Increases adhesion of tapes and wafers. Gentle, convenient and scent-free.

Item	Product#	Size	HCPC
Large Wipes, 25 Box	SNS00807	5"x7"	A5120
Individual Wipes, 25 Box	SNS80725	2"x2"	A5120
Individual Wipes, 75 Box	SNS80775	2"x2"	A5120
Individual Wipes, 100 Box	SNS80744	2"x2"	A5120

No-Sting Skin Barrier Wand

ALCOHOL FREE!



Applicator wand provides a no-sting, alcohol-free transparent protective coating to the area applied. Increases adhesion of tapes and wafers. Conveniently packaged applicator wands are gentle on the skin, sterile and scent-free.

Item	Product#	Size	HCPC
Wands, 25 Box	SNS80711	1mL	A4369

No-Sting Skin Barrier Spray

ALCOHOL FREE!



No-sting, alcohol-free spray that provides a transparent protective coating onto the area sprayed. Increase adhesion of tapes and wafers. Gentle convenient and scent-free. Ideal for use over powder to create a crusting technique to protect traumatized skin.

Item	Product#	Size	HCPC
Spray Bottle	SNS80792	2 oz.	A4369

Skin Barrier Film



Alcohol based wipe has a very gentle formula that dries quickly. Helps to protect the skin against irritation and adhesive buildup.

Item	Product#	Size	HCPC
Individual Wipes, 50 Box	SNS81850	2"x2"	A5120

Ostomy Skin Barrier Powder



Protects and promotes healing of irritated or broken skin. Dusted onto the skin it absorbs moisture and enhances comfort and wear time of the pouching system. Perfect for crusting.

Item	Product#	Size	HCPC
Applicator Bottle	SNS92305	5 oz.	A4371
Applicator Bottle	SNS92301	1 oz.	A4371



Skin Barrier Ring Conforming Seal



The conforming seals are ideal for any type of ostomy. It will create a secure seal at the base of the stoma and is highly absorbent. These rings “Hold Up Better”, promoting longer pouch wear time and do not leave a mess when removed from the skin. They are very moldable and stretchable with a strong tack and require no cutting. Great for replacing paste.

Item	Product#	Size	HCPC
Individual Seals, 10 Box	SNS684U2	2"	A4385
Individual Seals, 10 Box	SNS684U4	4"	A4385

Integrity Skin Barrier Ring



These durable skin barrier rings maintain their integrity and do not melt down from excessive moisture or body heat. Integrity skin barrier rings are ideal for those with a Colostomy that do not need absorption in a ring. These rings provide a comfortable protective seal to prevent leakage. Sticky, thin, low profile seals stretch and mold easily to fit around the stoma.

Item	Product#	Size	HCPC
Individual Seals, 20 Box	SNS68002	2"	A4385

Customer Testimonial:

"I am using three of your products now, Integrity Seals for fill-in areas, Conforming Seals for around my stoma and also your deodorant and have been getting three days longer wear time now over one day wear time from the other products. You've kinda saved my life."
David B.



Skin Barrier Arcs



Under the pouching system: The skin barrier arcs can be placed directly on the irritated skin around the stoma to create a smooth platform to adhere the pouch or wafer. The arcs protect the skin from trauma caused by repetitive removal of tapes and adhesives.

Over the pouching system: The skin barrier arcs can be placed around the wafer to secure the edges from lifting up, creasing or loosening. Use 1-3 arcs around the edges. Skin barrier arcs provide an excellent alternative to waterproof tapes and give you a waterproof seal around the wafer when bathing, showering or swimming.

Item	Product#	Size	HCPC
Tray - Regular, 30 Arcs	SNS20630	1"Width	A4362
Tray - X-tra Wide, 30 Arcs	SNS21130	2"Width	A4362

1/2 ARC Skin Barrier Strips

NEW SIZE NOW AVAILABLE!



A hydrocolloid 1/2 circle arc that is water resistant and perfect for picture framing around the ostomy pouch. Thin profile for comfortable wear. Non-sterile.

Item	Product#	Size	HCPC
Tray-1/2 Arc, 20 Strips	SNS22222	1.25" Width	A4362

Skin Barrier Sheet



A thin hydrocolloid sheet that protects the skin from stoma output. Highly absorbent. Cut to fit - no starter hole. Non-sterile.

Item	Product#	Size	HCPC
Individual Packaged			
5 Box	SNS21605	4"x4"	A4362

No-Sting Skin Barrier Paste



An alcohol-free formulation to reduce the sting and pain caused by conventional alcohol based pastes. Easy to squeeze tube helps you get the paste out smooth and quick. Creates an optimum seal for a better fitting pouch and longer wear time.

Item	Product#	Size	HCPC
Screw Cap Tube	SNS92802	2 oz.	A4406

Skin Barrier Paste



A protective pectin-based paste that helps to prevent leakage and skin irritation around the stoma providing an optimum seal for a better fit and longer wear time. Easy to squeeze tube helps you get the paste out smooth and quick.

Item	Product#	Size	HCPC
Screw cap Tube	SNS90502	2 oz.	A4406

Skin Barrier Paste - Fast Drying



A protective pectin-based paste that hardens upon application. Contains alcohol for faster drying. Seals and prevents skin irritation around the stoma and is great as a filler for those problem areas.

Item	Product#	Size	HCPC
Screw cap Tube	SNS90516	2 oz.	A4406

Ostomy Pouch Deodorant



Fantastic odor control with no perfumes or artificial scents added for discreet protection. Blue formula helps to see the liquid inside the pouch.

Item	Product#	Size	HCPC
Push Top Dispenser Bottle	SNS40202	2 oz.	A4394
Push Top Dispenser Bottle	SNS40208	8 oz.	A4394

Assure C Odor Eliminator



Assure C is an odor eliminator that also contains a lubricant to help empty pouch contents. This clear formulation eliminates the fear or concern of staining associated with the blue colored deodorants. No perfumes or artificial scents added.


Item	Product#	Size	HCPC
Single Serve Packet	SNS41425	1/4 oz.	A4394
Tilt Tip Dispenser Bottle	SNS41404	4 oz.	A4394
Tilt Tip Dispenser Bottle	SNS41408	8 oz.	A4394



**Be Confident
and Secure!**

Pouches

One-Piece Drainable Pouch



Standard wear drainable skin barrier with a comfortable fit panel and dual adhesive tapes. Cut-to-fit, clamp closure. Latex free.

Item	Product#	Size	HCPC
Transparent Pouch 10 Box	SNS12T100	12"	A5061
Opaque Pouch 10 Box	SNS12B100	12"	A5061

Two (2) clips per box!

One-Piece Drainable Pouch with Filter



Standard wear drainable skin barrier with a comfortable fit panel and dual adhesive tapes. Odor controlling filter to reduce gas buildup and ballooning. Cut-to-fit, clamp closure. Latex free.

Item	Product#	Size	HCPC
Transparent Pouch 10 Box	SNS12T10F	12"	A4424
Opaque Pouch 10 Box	SNS12B10F	12"	A4424

Stoma Cap



A one-piece, lightweight flexible system. Easy to manage. A temporary option to assist with showering, swimming or moments of intimacy. Latex free.

Item	Product#	Size	HCPC
Stoma Cap w/Hydrocolloid Collar 30 Box	SNS14502	3/4-2"	A5055
Stoma Cap w/Tape Collar 30 Box	SNS14506	3/4-2"	A5055



Security™

Hernia/Ostomy
Support Belt

- ★ Revolutionary 2-way stretch fabric is supportive yet comfortable
- ★ The fabric conforms to fit an individual's shape during daily activity
- ★ A unique glove design allows easy control of the belt



"As an ostomate and ET nurse, I've been a long time user and advocate of support belts. In my personal experience, both versions of the Safe n Simple Security support belts are amazingly comfortable while providing the support I need in my physically active lifestyle. The conforming fit glides with body movement and isn't visible when wearing snug fitting clothing. Whether the goal is prevention or hernia support, I'd highly recommend this product."

- K. Wickham, MI

Hernia/Ostomy Support Belts

Made from a unique fabric that gives equal and even stretch. This ensures pressure is evenly distributed throughout the garment. The fabric provides complete support to peristomal hernias while still having the flexibility to ensure the normal function of the stoma.

- Security belts were found to be the most comfortable while providing excellent support
- Security belts can be worn for longer periods of time without feeling restricted
- Security belts are fully adjustable so can be used at all times
- Security belt has a soft Velcro fastener which is ideal for easy fitting and can be adjusted as required

Sizes	Measurement Width
Extra Small	up to 29.5 inches up to 75cms
Small	29.5 to 36.5 inches 75cms to 90cms
Medium	36.5 to 41.25 inches 90cms to 105cms
Large	41.25 to 47.25 inches 105cms to 120cms
Extra Large	47.25 to 53.25 inches 120cms to 135cms
XX Large	53.25 to 59 inches 135cms to 150 cms

Larger sizes available by special order

Visit our website to find out where to order:

www.sns-medical.com

Questions? Please call: 844-767-6334

Safe n Simple SECUREWEAR
Clean. Comfortable. Confident.

Security Support Belts


Security
 Hernia/Ostomy
 Support Belt



Support Belt without Hole

Size	Depth	Color	Product#	Size	Depth	Color	Product#
XS	6"/16cm	Beige	SNS36BG1	XS	6"/16cm	Black	SNS36BK1
S	6"/16cm	Beige	SNS36BG2	S	6"/16cm	Black	SNS36BK2
M	6"/16cm	Beige	SNS36BG3	M	6"/16cm	Black	SNS36BK3
L	6"/16cm	Beige	SNS36BG4	L	6"/16cm	Black	SNS36BK4
XL	6"/16cm	Beige	SNS36BG5	XL	6"/16cm	Black	SNS36BK5
XXL	6"/16cm	Beige	SNS36BG6	XXL	6"/16cm	Black	SNS36BK6
XS	8"/20cm	Beige	SNS3BG21	XS	8"/20cm	Black	SNS3BK21
S	8"/20cm	Beige	SNS3BG22	S	8"/20cm	Black	SNS3BK22
M	8"/20cm	Beige	SNS3BG23	M	8"/20cm	Black	SNS3BK23
L	8"/20cm	Beige	SNS3BG24	L	8"/20cm	Black	SNS3BK24
XL	8"/20cm	Beige	SNS3BG25	XL	8"/20cm	Black	SNS3BK25
XXL	8"/20cm	Beige	SNS3BG26	XXL	8"/20cm	Black	SNS3BK26
XS	10.25"/26cm	Beige	SNS3BG61	XS	10.25"/26cm	Black	SNS3BK61
S	10.25"/26cm	Beige	SNS3BG62	S	10.25"/26cm	Black	SNS3BK62
M	10.25"/26cm	Beige	SNS3BG63	M	10.25"/26cm	Black	SNS3BK63
L	10.25"/26cm	Beige	SNS3BG64	L	10.25"/26cm	Black	SNS3BK64
XL	10.25"/26cm	Beige	SNS3BG65	XL	10.25"/26cm	Black	SNS3BK65
XXL	10.25"/26cm	Beige	SNS3BG66	XXL	10.25"/26cm	Black	SNS3BK66

Need a larger or smaller size? Two holes? Special measurements? We can make it custom for you!

HCPC: A4396

Also Available on Amazon

	Size	Depth	Color	Hole	Product#	Size	Depth	Color	Hole	Product#
Support Belt with Hole - 6"/16cm	XS	6"/16cm	Beige	Right	SNS36GR1	XS	6"/16cm	Beige	Left	SNS36GL1
	S	6"/16cm	Beige	Right	SNS36GR2	S	6"/16cm	Beige	Left	SNS36GL2
	M	6"/16cm	Beige	Right	SNS36GR3	M	6"/16cm	Beige	Left	SNS36GL3
	L	6"/16cm	Beige	Right	SNS36GR4	L	6"/16cm	Beige	Left	SNS36GL4
	XL	6"/16cm	Beige	Right	SNS36GR5	XL	6"/16cm	Beige	Left	SNS36GL5
	XXL	6"/16cm	Beige	Right	SNS36GR6	XXL	6"/16cm	Beige	Left	SNS36GL6
	XS	6"/16cm	Black	Right	SNS36KR1	XS	6"/16cm	Black	Left	SNS36KL1
	S	6"/16cm	Black	Right	SNS36KR2	S	6"/16cm	Black	Left	SNS36KL2
	M	6"/16cm	Black	Right	SNS36KR3	M	6"/16cm	Black	Left	SNS36KL3
	L	6"/16cm	Black	Right	SNS36KR4	L	6"/16cm	Black	Left	SNS36KL4
	XL	6"/16cm	Black	Right	SNS36KR5	XL	6"/16cm	Black	Left	SNS36KL5
	XXL	6"/16cm	Black	Right	SNS36KR6	XXL	6"/16cm	Black	Left	SNS36KL6
Support Belt with Hole - 8"/20cm	XS	8"/20cm	Beige	Right	SNS3BGR1	XS	8"/20cm	Beige	Left	SNS3BGL1
	S	8"/20cm	Beige	Right	SNS3BGR2	S	8"/20cm	Beige	Left	SNS3BGL2
	M	8"/20cm	Beige	Right	SNS3BGR3	M	8"/20cm	Beige	Left	SNS3BGL3
	L	8"/20cm	Beige	Right	SNS3BGR4	L	8"/20cm	Beige	Left	SNS3BGL4
	XL	8"/20cm	Beige	Right	SNS3BGR5	XL	8"/20cm	Beige	Left	SNS3BGL5
	XXL	8"/20cm	Beige	Right	SNS3BGR6	XXL	8"/20cm	Beige	Left	SNS3BGL6
	XS	8"/20cm	Black	Right	SNS3BKR1	XS	8"/20cm	Black	Left	SNS3BKL1
	S	8"/20cm	Black	Right	SNS3BKR2	S	8"/20cm	Black	Left	SNS3BKL2
	M	8"/20cm	Black	Right	SNS3BKR3	M	8"/20cm	Black	Left	SNS3BKL3
	L	8"/20cm	Black	Right	SNS3BKR4	L	8"/20cm	Black	Left	SNS3BKL4
	XL	8"/20cm	Black	Right	SNS3BKR5	XL	8"/20cm	Black	Left	SNS3BKL5
	XXL	8"/20cm	Black	Right	SNS3BKR6	XXL	8"/20cm	Black	Left	SNS3BKL6
Support Belt with Hole - 10.25"/26cm	XS	10.25"/26cm	Beige	Right	SNS3GR61	XS	10.25"/26cm	Beige	Left	SNS3GL61
	S	10.25"/26cm	Beige	Right	SNS3GR62	S	10.25"/26cm	Beige	Left	SNS3GL62
	M	10.25"/26cm	Beige	Right	SNS3GR63	M	10.25"/26cm	Beige	Left	SNS3GL63
	L	10.25"/26cm	Beige	Right	SNS3GR64	L	10.25"/26cm	Beige	Left	SNS3GL64
	XL	10.25"/26cm	Beige	Right	SNS3GR65	XL	10.25"/26cm	Beige	Left	SNS3GL65
	XXL	10.25"/26cm	Beige	Right	SNS3GR66	XXL	10.25"/26cm	Beige	Left	SNS3GL66
	XS	10.25"/26cm	Black	Right	SNS3KR61	XS	10.25"/26cm	Black	Left	SNS3KL61
	S	10.25"/26cm	Black	Right	SNS3KR62	S	10.25"/26cm	Black	Left	SNS3KL62
	M	10.25"/26cm	Black	Right	SNS3KR63	M	10.25"/26cm	Black	Left	SNS3KL63
	L	10.25"/26cm	Black	Right	SNS3KR64	L	10.25"/26cm	Black	Left	SNS3KL64
	XL	10.25"/26cm	Black	Right	SNS3KR65	XL	10.25"/26cm	Black	Left	SNS3KL65
	XXL	10.25"/26cm	Black	Right	SNS3KR66	XXL	10.25"/26cm	Black	Left	SNS3KL66

DermaPro Tapes

DermaPro Waterproof Plastic Tape



Sempurity DermaPro Waterproof Plastic Tape is superior for securing wound dressings and repelling moisture from outside contaminants. DermaPro has a latex free adhesive and is gentle to remove from sensitive skin areas. The roll contains a start strip for easy dispensing.

Item	Product#	Size	HCPC
Plastic Tape, 1 Roll	SNS57220	2"x15'	A4452

DermaPro Waterproof Silicone Tape



Sempurity DermaPro Waterproof Silicone Tape is latex free and has a silicone adhesive that is gentle to remove from sensitive skin areas. Superior for securing wound dressings and repelling moisture from outside contaminants. Contains 12-1" rolls per box, and 6-2" rolls per box.

Item	Product#	Size	HCPC
Silicone Tape, 12 Rolls	SNS57230	1"x15'	A4452
Silicone Tape, 6 Rolls	SNS57232	2"x15'	A4452



Ostomy Allowables

Medicare Allowables...30 day supply

Peri-Stoma Cleanser & Adhesive Remover	50
Adhesive Remover	50
No-Sting Skin Barrier Film	25
No-Sting Skin Barrier Spray.....	2 oz.
Skin Barrier Film	25
Ostomy Skin Barrier Powder.....	5 oz.
Skin Barrier Ring Conforming Seal	20
Skin Barrier Arcs	20
Skin Barrier Sheet	20
Skin Barrier Paste	4 oz.
No-Sting Skin Barrier Paste	4 oz.
Ostomy Appliance Deodorant.....	as necessary
Assure C Odor Eliminator.....	as necessary
Pouches - Drainable.....	20
Stoma Cap.....	31
Tape 1" x 10 yds.....	2
Tape 2" x 10 yds.....	1

Safe n Simple has provided the above information as a guideline of Medicare policies and assumes no responsibility for inconsistencies. All medical and claims determinations will be made by Medicare and your insurance plan.



Stay Connected....

Contact Our Team...

1. Call Safe n Simple Customer Service
Toll Free 844-767-6334 or
Amy Ham 248-214-4877
Mon-Fri 9am-5pm EST
2. Find a Supplier or Order from our E-Store at:
www.sns-medical.com
3. Mailing Address: 8112 Rosebud Lane
Clarkston, MI 48348
Shipping Address: 5827 Terex
Clarkston, MI 48346

Inside Sales Team...

Holly Loos	248-496-2508
Michelle Pitylak	248-212-2243

Outside Sales Team...

Midwest:	Joseph LeRoy	248-974-3245
Northeast:	Michele Endreszl	248-977-8324
Southeast:	Beth Peter	561-526-6718
Ind. Rep. CA:	Dan Penrose	909-730-0707
Ind. Rep. IN/OH:	Joseph DeHeer	317-771-4319
Ind. Rep. MD/VA:	Jim Spradlin	804-752-2897
Ind. Rep. Canada:	Peter West	519-726-5959

Samples Team...

To request samples call our customer service team
or email us at: samples@sns-medical.com

Clinicians/New Ostomy Patients:
Request a **SimpleStart Ostomy Sample Kit**
email us at: simplestart@sns-medical.com
or fax us at: 248-575-4172

Stay Connected....

Supplier Orders:


1. Call Safe n Simple Customer Service
Toll Free 844-767-6334 or
Amy Ham 248-214-4877
Mon-Fri 9am-5pm EST
2. Email Us: purchaseorder@sns-medical.com
3. FAX Us: 248-698-3454
4. Mailing Address: Shipping Address:
8112 Rosebud Lane 5827 Terex
Clarkston, MI 48348 Clarkston, MI 48346

Marketing Team:

Have a special event or sponsorship opportunity you would like to share with Safe n Simple?

Please feel free to email us: marketing@sns-medical.com

Call Us:
Michelle Bliszack - Marketing Manager
248-408-9501 Direct
844-767-6334 Toll Free



Our Mission:

Safe n Simple is an innovator, developer and marketer of new medical products. We enhance the patient's quality of life with superior products at an affordable price bringing value to the communities in which we live.

Safe n Simple Ostomy Care Products

Adhesive Remover
Peri-Stoma Cleanser & Adhesive Remover
No-Sting Skin Barrier Film
Skin Barrier Film
Ostomy Skin Barrier Powder
Skin Barrier Ring - Conforming Seal
Skin Barrier Ring - Integrity Seal
Skin Barrier Arc
X-Tra Wide Skin Barrier Arc
1/2 Arc Skin Barrier Strips
Skin Barrier 4"x4" Sheet
Skin Barrier Paste
No-Sting Skin Barrier Paste
Skin Barrier Paste - Fast Drying
Ostomy Pouch Deodorant
Assure C Odor Eliminator
One-Piece, CTF Pouches
One-Piece, CTF Pouches with Filters
Stoma Caps
Security Support Belts

Simpurity Wound Care Products

Alginate - Dressing Pads & Rope
Alginate Silver - Dressing Pads & Rope
Collagen - Dressing Pads & Powder
Collagen Silver - Dressing Pads
DermaPro Waterproof Plastic Tape
DermaPro Waterproof Silicone Tape
Fibergel - Dressing Pads
Fibergel Silver - Dressing Pads
Foam - Dressing Pads
Foam - Adhesive Bordered Dressing Pads
Foam - Silicone Adhesive Dressing Pads
Foam - Silicone Bordered Dressing Pads
Foam - Silicone Silver Foam Dressing Pads
Foam - Silver Dressing Pads
Hydrocolloid Dressing Pads
Hydrogel Dressing Sheets
Hydrogel Gauze Dressing Pads
Polymer Dressing Pads
Wound Care Kits



Safe n Simple[™]
Clean. Comfortable. Confident.

Clarkston, MI 48346 • P: 844.767.6334 • F: 248.698.3454

www.sns-medical.com

Email: samples@sns-medical.com

Like us on Facebook:

<https://www.facebook.com/Safe-n-Simple>