

In-Room™ Sharps Containers with Tortuous Path Lid



In-Room™ System with Tortuous Path Lid

- One-hand, one-step, no-touch sharps depositing.
- Horizontal loading, non-mechanical sharps opening.
- Tortuous-path design impedes hand entry.
- Various mounting brackets and glove box holders are available.

Code	Description	Dimensions	Ship Case
31144010	5 Quart, Transparent Red	12.25"H x 4.25"D x 11"W	30
31143897	5 Quart, Clear	12.25"H x 4.25"D x 11"W	30

GatorGuard™ Safety In-Room™ Sharps Containers



GatorGuard™ In-Room™ Sharps Containers

- Counter-balanced door design allows for disposal of various sized items, from a lightweight butterfly to a 60mL syringe.
- Easy visibility of recommended fill line.
- Signals when container's capacity is reached with a FULL indicator.
- Horizontal drop system maximizes full use of container.
- Thumb tabs lock the container closed for final disposal when full.
- Complete line of brackets and accessories available.

Code	Description	Dimensions	Ship Case
GatorGuard™ Jr. Sharps Containers			
31353603	5 Quart, Transparent Red	14"H x 6"D x 13"W	14
GatorGuard™ Sharps Containers			
31323333	2 Gallon, Transparent Red	17.5"H x 6"D x 14"W	12
31314886	3 Gallon, Transparent Red	20.5"H x 6"D x 14"W	12