

Comfort**Foam**™ Border Self-Adherent Soft Silicone Foam Dressing with Border

ComfortFoam Border is a self-adherent silicone foam island dressing consisting of a soft silicone contact layer, a flexible, absorbent polyurethane foam pad and a vapor-permeable, moisture-proof outer film. The silicone adhesive sticks to surrounding skin but not to the wound bed, minimizing pain and the risk of damage to the wound area.

Benefits:

- Silicone adhesive is gentle to both wound bed and skin for
 Supports moist wound healing safe, easy removal
- Low profile border for added adhesion and flexibility
- Absorbs exudate while minimizing the risk of maceration
- Border is showerproof and prevents exudate leakage
- Provides thermal insulation and protection
- Can be safely worn for several days, depending on wound condition and exudate levels
- Latex free

Indications:

For the management of partial and full thickness wounds, including:

- Pressure ulcers
- Arterial ulcers
- Traumatic wounds

- Diabetic ulcers
- Venous ulcers

	Product	Size	Reorder #	Units
NEW	ComfortFoam Border	2x2	43220	10/Box
NEW	ComfortFoam Border	2x5	43250	10/Box
NEW	ComfortFoam Border	3x3	43330	10/Box
	ComfortFoam Border	4x4	00317E	10/Box
NEW	ComfortFoam Border	4x8	43480	5/Box
NEW	ComfortFoam Border	4x12	43412	5/Box
	ComfortFoam Border	6x6	00318E	10/Box
NEW	ComfortFoam Border	6x8	43680	5/Box

	Product	Size	Reorder #	Units
NEW	ComfortFoam Border	7x7	43770	5/Box
NEW	ComfortFoam Border	Heel/Elbow	43580	5/Box
NEW	ComfortFoam Border	Sacral (7x7)	43880	5/Box
NEW	ComfortFoam Border	Large Sacral (9x9)	43990	5/Box

