

Dressing Change Tray – Sterile

Convenience in a Complete and Sterile Kit

Dynarex Dressing Change Tray is a convenient sterile kit for aseptic dressing changes. All necessary items are packed in a single, disposable tray, eliminating the need to maintain inventory on multiple items and assembling them for a procedure.

- Sterile
- Latex free
- Components arranged in procedural sequence for convenience
- Cost-effective and disposable

Contents:

- (1) Instructions for use
- (1) Wound care tray
- (1) Non-woven face mask(1) Non-woven split sponge,
- 2" x 2", 4 ply
- (1) Non-adherent sponge, 2" x 2"
- (1) Gauze sponge, 2" x 2", 8 ply
- (1) Gauze sponge, 4" x 4", 8 ply
- (1) Waterproof drape, 45 x 34cm
- (1) Dressing change label
- (1) PVC gloves, size M,
- powder free, white, with wallet
- (1) Tape measure, 36"
- (3) PVP swabsticks
- (3) Alcohol swabsticks
- (1) Roll Transparent Dressing Tape
- (1) Transparent Dressing, 4" x 4 $\frac{3}{4}$ "

Reorder No.	Description	Quantity
4705	Dressing Change Tray – Sterile	20/cs

Depending on market demand and Dynarex stock, we will be happy to address any orders for items/sizes currently not shown.

Dressing Change Tray

4 5 6 7

Visit our website to see our Minor Laceration Tray (#4523) and our Staple/Suture Removal Kits (#4521-4522)



For more information contact your Dynarex sales representative 888-396-2739