



Enduro XL3⁺

The Enduro XL3⁺ is a heavy duty three-wheel scooter with a 500 lb. weight capacity. It has a delta style tiller which houses the controls for the lighting system, turn signals, wig wag and horn. It has a dual purpose wig wag which allows the user to push forward or pull back on the drive lever depending on dexterity for ease of driving. The unit comes standard with a 20in. seat that can accommodate most riders, and has swivel and recline features for added comfort and ease of use.

- + 500lb Weight Capacity
- + Slide & Swivel Seating
- + Full Suspension System
- + Full Lighting Package
- + Adjustable Tiller

● Burgundy

Model Number: 778XLSBN
Medicare Billing Code: K0899

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

Technical Specifications

Overall Dimensions (L x W x H)	in	53" x 28" x 52"
Number, Size of Tires	Front, Rear	1, 12" (2, 12")
Suggested User Weight ¹	lb	Max. 500
Weight of Heaviest Part ²	lb	104
Weight ³	lb	222
Each Battery Module Weight	lb	50
Battery Size	Type x PCS	12V50-75AH x 2
Charger	Amp	Off board (8A)
Output Power of Motor	HP	2
Maximum Speed	MPH	5
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(°)	14(8)
Range (Per Charge) ⁶	Mile	25-39
Turning Radius	in	49"
Ground Clearance	in	4"
Handle Bar	Type	Delta
Model #		778XLSBN
Standard Seat	W x D x BH	SV201820B

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.