

## How does **MEDIHONEY®** / Active *Leptospermum* Honey help to support the healing of routine and challenging wounds?

### *Gel* vs *Paste* Dressings

- Helps to maintain a moist environment for optimal healing<sup>1-6</sup>
- High osmotic potential pulls fluid from deeper tissues to the wound surface, aiding in autolytic debridement<sup>1-6</sup>
- Low pH helps to lower the pH within the wound environment<sup>7,8</sup>, which has been shown to have wound healing benefits.<sup>9</sup>
- Is easy-to-use, non-toxic, natural and safe with the potential for extended wear times

#### *Gel*



#### *Paste*




- 80% Active *Leptospermum* Honey and 20% Natural Gelling agents
- More viscous formulation with increased stability at the site of the wound
- Demonstrated positive clinical outcomes<sup>1, 5, 6</sup>
- Optimal for superficial to full thickness wounds

- 100% Active *Leptospermum* Honey
- Demonstrated positive clinical outcomes<sup>2-4</sup>
- Optimal for deep tunneled or undermined wounds

# MEDIHONEY® Ordering Information

Order Code	Description	Packaging unit/Case	HCPCS*
<b>Gel</b>			
31805	0.5 oz tube	10/box, 4 boxes/case	A4649
31815	1.5 oz tube	1/box, 12 boxes/case	A4649
<b>HCS</b>			
<b>Non-adhesive</b>			
31622	2.4" x 2.4"	10/box, 5 boxes/case	A4649
31644	4.33" x 4.33"	10/box, 5 boxes/case	A4649
31688	8" x 8"	5/box, 4 boxes/case	N/A
31612	8" x 12"	2/box, 5 boxes/case	N/A
<b>Fenestrated - Non-adhesive</b>			
31618	1.8" x 1.8"	10/box, 5 boxes/case	A4649
<b>Adhesive</b>			
31722	2.8" x 2.8" (4.3" x 4.3" with adhesive border)	10/box, 5 boxes/case	A4649
31744	4 ½" x 4 ½" (6" x 6" with adhesive border)	10/box, 5 boxes/case	A4649
<b>Calcium Alginate</b>			
31012	¾" x 12"	5/box, 4 boxes/case	A4649
31022	2" x 2"	10/box, 10 boxes/case	A4649
31045	4" x 5"	10/box, 5 boxes/case	A4649
<b>Honeycolloid™</b>			
<b>Non-adhesive</b>			
31222	2" x 2"	10/box, 10 boxes/case	A4649
31245	4" x 5"	10/box, 5 boxes/case	A4649
<b>Adhesive</b>			
31422	2" x 2" (3 ½" x 3 ½" with adhesive border)	10/box, 10 boxes/case	A4649
31445	4 ½" x 4 ½" (6" x 6" with adhesive border)	10/box, 5 boxes/case	A4649
<b>Paste</b>			
31505	0.5 oz tube	10/box, 4 boxes/case	A4649
31515	1.5 oz tube	1/box, 12 boxes/case	A4649
31535	3.5 oz tube	1/box, 12 boxes/case	A4649



\*Refer to [www.dmeptac.com](http://www.dmeptac.com) for the most current HCPCS coding of MEDIHONEY® surgical dressings.

Pair MEDIHONEY® with our super absorbent cover dressing, XTRASORB®. It's osmotic gradient pulls exudate to the back of the dressing and converts it into a gel, locking it away - even under compression!

CATMHGELPASTE

 SUPER ABSORBENT DRESSING  
**XTRASORB® HCS**



 SUPER ABSORBENT DRESSING  
**XTRASORB® Foam**



 SUPER ABSORBENT DRESSING  
**XTRASORB® Classic**



[Click here to view detailed product descriptions and prices.](#)

**Call 800-397-5899 to order today!**