

3. PRODUCT-SPECIFIC CHECKS ON REDUCED SAMPLE SIZE (as applicable for specific product / order) :
For Rubber Gloves such as Examination Gloves, Surgical Gloves, Household Gloves, Nitrile Gloves, s etc.

3.1 Watertight Testing, Sample Size = 500 pieces
 Requested AQL : 1.5 Defect max. allowed : 14 Found : 7 Within Beyond AQL
 Defects breakdown : S: 1 pc M: 3pcs L: 3pcs XL: 0pc
 Remarks (if any): -

3.2 Dimension check on: 13 pieces sub-samples:
 Requested AQL 4.0 Defect max. Allowed : 1 Found : 0 Within AQL Beyond AQL

Size	S	M	M	M	M	M	L	L	L	L	L	L	XL
Width (± 10mm) S: 80 / M:95 / L:110 / XL: ≥110	85	95	97	96	95	93	108	103	104	106	107	105	113
Length Min. 240mm	242	250	249	245	241	245	248	242	251	245	251	250	245
Thickness at finger Min: 0.05mm	0.11	0.09	0.09	0.10	0.10	0.10	0.10	0.11	0.13	0.11	0.11	0.11	0.11
Thickness at palm Min : 0.05mm	0.07	0.07	0.08	0.07	0.07	0.08	0.07	0.08	0.07	0.07	0.08	0.07	0.08

3.3 Physical Properties

Not tested

Testing conducted at WRP observed by SGS on 03 Agu, 2020. Refer to factory test result as below. Sample size = 13 pieces. .

3.4 Dusting Powder

Not tested

Testing in progress
Sample size =

3.5 Protein Test

Not tested

Testing in progress
Sample size =

3.6 Quantity per bags 100 pcs		Sample size = 8 dispenser	
Size:	Small	Size:	Large
Qty:	102, 101	Qty:	98, 100
Size:	Medium	Size:	X-Large
Qty:	98, 99	Qty:	101, 100

Test equipment used:

Equipment owner	Equipment description	Equipment No.	Brand	Model No.	Calibration date	Due date	Calibration certificate issued by accredited lab
SGS	Thickness Gauge	CRS/TH/3	Mitutoyo	-	29.01.20	29.01.21	Sendi Mahir
SGS	Steel ruler	CRS/SR/5	-	-	30.01.20	20.01.20	Sendi Mahir
Factory	Watertight Machine	WL-39	-	-	15.07.20	15.09.20	Factory

"This document is issued by the Company under its General Conditions of Service accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein.

Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law."