



SilverDerm 7™ Antimicrobial Wound Contact Silver Dressing

Antimicrobial contact dressing for wound bed protection. Non-adherent and porous - allows wound exudate to pass through to absorbent dressing materials.

Benefits:

- Provides broad spectrum antimicrobial activity.
- Silver ions are released within 30 minutes.
- 7-day wear time - sustained delivery of silver ions.
- Assists with autolysis of necrotic wounds.
- Reduces pain and trauma associated with dressing change.
- Hypoallergenic

Indications:

- Partial-thickness wounds
- Full-thickness wounds
- Skin grafts
- Donor sites

Product	Size	Reorder #	Units	HCPC
SilverDerm 7	4x4	00550E	10/Box	A2970

DermaRite®