

Juzo Dynamic Compression Stump Shrinkers

Art. 3511, 3512 · FG

Lengths in cm

Standard sizes · Circumferences in cm

I / XS

II / S

III / M


IV / L

V / XL

VI / XXL

PFG

short	18.0 – 22.0
medium	23.0 – 27.0
long	28.0 – 32.0
knee exart.	33.0 – 37.0



Silicone border - Top border

cG

42,0 – 45,0

45,0 – 49,0

49,0 – 53,0

53,0 – 57,0

57,0 – 61,0

61,0 – 65,0

cF*

30,0 – 32,0

32,0 – 35,0

35,0 – 38,0

38,0 – 41,0

41,0 – 44,0

44,0 – 47,0

* cF: 5 cm from the distal end

Juzo Dynamic Compression Stump Shrinkers

Art. 3511, 3512 · CF

Lengths in cm

Standard sizes · Circumferences in cm

I / XS

II / S

III / M

IV / L

V / XL

VI / XXL

Silicone border · Top border

PEF

24.0 – 26.0

cF

40,0 – 42,0

42,0 – 45,0

45,0 – 48,0

48,0 – 51,0

51,0 – 54,0

54,0 – 57,0

PCE

4.0 – 6.0

cE

30,0 – 32,0

32,0 – 35,0

35,0 – 38,0

38,0 – 41,0

41,0 – 44,0

44,0 – 47,0

short

12.0 – 14.0

medium

20.0 – 22.0

long

cC*

28,0 – 30,0

30,0 – 33,0

33,0 – 36,0

36,0 – 39,0

39,0 – 42,0

42,0 – 45

* cC: 5 cm from the distal end

