

ottobock.



Size 3 pictured


Everyday Activity Seat

Everyday seating for everyday experiences

Children and young adults with special needs will be better supported for independent activities of daily living with the redesigned Leckey Everyday Activity Seat. With new positioning options, updated fabrics, and a sleeker design, your client or child can be more successful at daily activities rather than being limited by their adaptive equipment. The Everyday Activity Seat is designed specifically for children or young adults with moderate needs who still require pelvic stability but do not need the same advanced level of support as those with complex needs. The system includes a range of postural options, padding, and adjustments to help ensure that the child is free to use their arms and hands more effectively.

Quality for life

Everyday Activity Seat



	Pal	Everyday	Squiggles	Mygo	Kit
Chailey	4-7	3-6	1-4	1-3	1-3
GMFCS	II-III	III-IV	III-IV	IV-V	IV-V

Technical data

Size	1	2	3
Max load	66 lbs/30 kg	110 lbs/50 kg	220 lbs/100 kg
Weight seat and basic chassis	50 lbs/22.5 kg	60 lbs/27 kg	79 lbs/35.9 kg
Weight seat and hi-low chassis	65 lbs/29.5 kg	75 lbs/34 kg	108 lbs/49 kg
Seat height basic chassis	17.7 - 23.6 in / 450 - 600 mm	17.7 - 23.6 in / 450 - 600 mm	17.7 - 23.6 in / 450 - 600 mm
Seat height hi-low chassis	9.8 - 23.6 in / 250 - 600 mm	12.6 - 23.6 in / 320-600 mm	17.7 - 23.6 in / 450 - 600 mm
Seat width	6.3 - 11 in / 160 - 280 mm	10.2 - 15 in / 260 - 380 mm	13.4 - 19.7 in / 340 - 500 mm
Seat depth (standard)	5.9 - 11 in / 150 - 280 mm	9 - 15 in / 230 - 380 mm	13 - 18.9 in / 330 - 480 mm
Seat depth (extended)	n/a	11.8 - 17.7 in / 300 - 500 mm	n/a
Lateral distance - chest	5.9 - 8.7 in / 150 - 220 mm	6.7 - 9.4 in / 170 - 240 mm	9.4 - 14.6 in / 240 - 370 mm
Back height	11.8 - 15.7 in / 300 - 400 mm	15.7 - 20.5 in / 400 - 520 mm	17.7 - 24.8 in / 450 - 630 mm
Back angle	30° supine 10° prone	30° supine 10° prone	30° supine 10° prone
Armrest height	5.9 - 9.8 in / 150 - 250 mm	7.9 - 11.8 in / 200 - 300 mm	7 - 12.7 in / 180 - 325 mm
Lower leg length	5.9 - 9.8 in / 150 - 250 mm	9.8 - 17.7 in / 250 - 450 mm	9.8 - 17.7 in / 250 - 450 mm
Lower leg angle	-10° to +50°	-10° to +50°	-10° to +50°
Footplate angle	-15° to +10°	-15° to +10°	-15° to +10°
Tilt-in-space	30° supine 10° prone	30° supine 10° prone	30° supine 10° prone

Color swatches

Fabric and Vinyl



Accessories



• Trunk harness (black)



• Hip guides and adduction pad with covers



• Chest harness (cover included - size 1 and 2 only)



• Contoured headrest hardware (black vinyl, size 3 only)



• Basic chest laterals (cover included)



• Rigid PU chest laterals (cover included)



• Armrests



• Four-point pelvic harness



• Sandals, including straps. Sandal raiser (single, not shown)



• Activity tray (grey)