

Bite Blocks

McKesson Bite Blocks are used to stabilize the mouth and provide supplemental oxygenation during gastrointestinal endoscopy. They help provide a safe passage through the patient's mouth, while helping protect the teeth as well as valuable endoscope equipment during the procedure. Their supportive, rigid design helps avert collapse of airway and reduce hypoxia risk in sedated patients. The unique molding process creates a smooth finish to help prevent sharp edges, while strict quality control procedures and careful selection of raw materials help assure our bite blocks are manufactured without natural rubber latex.

Features and Benefits

- Supportive, rigid designed helps avert collapse of airway
- Aids in delivery of supplemental oxygen; helps avoid risk of hypoxia
- Eliminates need for separate nasal prongs/oxygen masks
- Precision manufacturing provides smooth edges for greater patient comfort
- Single-use design helps prevent cross-contamination
- Includes strap with size adjustments



ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Bite Blocks, Green	16-429	50ea/bx, 10bx/cs

Satisfaction Guaranteed
 If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support, clinical training and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.