

## SPECIFICATION

Overall Length	41"
Overall Width	27.4"
Overall Height	40"
Front Wheels	7.9"
Rear Wheels	13"
Weight W/ Batteries	143.5 lbs
Max. Speed	5 mph
Weight Capacity	350 lbs
Ground Clearance	3.2"
Grade Climable	8 degree
Curb Climable	1"
Turning Radius	33.5"
Brake	Electro-Mechanical
Seat Type	Padded Breathable Ballistic Nylon Back Rest W/ Adjustable Back Angle & Seat Belt
Seat Depth	15.6"
Armrest Height	7.5"
Seat Back Angle	100 degree
Seat Back Height	18.5"
Seat Width	20"
Drive Train	2 Motor, Rear-Wheel Drive
Motor Size	420W, 4600 r.p.m
Battery Size	(2) 12V. 36Ah
Weight of Batteries	55 lbs
Travel Range	12 Miles
Battery Charger	5A Off Board

\*Subject to change without notice. ( Rev. 0, 2020/02/07 )