

GO GO[®]

FOLDING SCOOTER

4-WHEEL SCOOTER

THE **ULTIMATE**
COMPACT, **FOLDABLE**
TRAVEL **SCOOTER**

Ready for transport and easily folded, the 4-wheel **Go-Go[®] Folding Scooter** is perfect for the active lifestyle. Offering great performance indoors and out at a maximum speed of 4 mph and a range per charge of up to 12.5 miles, the Go-Go Folding Scooter provides a comfortable and convenient ride. It is a super portable space saver that adjusts to fit perfectly inside any closet, corner, or vehicle trunk. The Go-Go Folding Scooter is sleek, stylish, and extremely easy to operate, making it the perfect travel mobility scooter choice.



Scooter Folded



Easy to Transport



Optional Armrests




PRIDE[®] **MOBILITY SCOOTERS**
LIVE YOUR BEST[®]

pridemobility.com |  

SPECIFICATIONS

GO-GO® FOLDING SCOOTER 4W: S19

Weight capacity:	250 lbs.
Maximum speed ¹ :	Up to 4 mph
Ground clearance ² :	2.3"
Turning radius ² :	32"
Overall length ² :	37"
Overall width ² :	19"
Drive wheels ² :	2" x 8" solid
Front wheels ² :	1.6" x 7" solid
Range per charge ^{1,3} :	Up to 12.5 miles (lithium-ion); Up to 9.3 miles (sealed lead acid)
Total weight without batteries:	51.4 lbs.
Heaviest piece when disassembled:	51.4 lbs. (unit does not disassemble - unit weight without battery pack)
Standard seating:	Type: Foldable, compact, vinyl-covered molded plastic Material: Gray Vinyl Dimensions: 16" x 13"
Battery requirements ⁴ :	Size: (2) 12 AH (sealed lead acid) standard or (1) 11.5 AH (lithium-ion) optional Weight: 12AH sealed lead acid - 9 lbs. each (2 required), 11.5 AH lithium-ion - 5 lbs. each (1 required)
Battery charger:	2-amp, off-board (sealed lead acid) or 2.5-amp, off-board (lithium-ion)
Color:	 White

FEATURES

- Unique folding design
- Folds in just three easy steps
- Height adjustable tiller
- Up to 12.5 miles range per charge
- Easy to transport at home, hotels, or airports
- Can be stowed upright or laid down
- Articulated front end adds greater stability
- Front and rear LED lights
- Optional armrests
- Folded: 30"H, 19"W, 18"D

TO FOLD:



TO UNFOLD:



- ¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ³ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.
- ⁴ Battery weight may vary based on manufacturer.

All specifications provided by Heartway® Medical Products.
Product distributed by Pride Mobility®.

NOTE: All specifications subject to change without notice.
The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.



Note: Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go-Go® Folding Scooter and its lithium-ion battery is permitted.

Pride
Mobility Products Corp.
— LIVE YOUR BEST —

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

(US) 800-800-8586 • (Canada) 888-570-1113
pridemobility.com