



Life Begins at the End of the Ramp !!

Don't be stopped by the edges of ramps, or beach mats! DeBug Beach Wheelchairs give you the freedom to enjoy the beach, or any outdoor terrain the way you deserve.

Without Limitations!

The DeBugs are designed by the Deming's who know first hand the obstacles wheelchair users face on a daily basis. They have used their personal experiences and incorporated them into the design and quality of their products.

So come on, roll off that edge and enjoy life

LONG LASTING #316L STAINLESS STEEL FRAME / ADJUSTABLE TELESCOPIC LEG REST / FITS THROUGH A STANDARD ADA DOORWAY / LOW MAINTENANCE DELRIN BUSHINGS / POLYURETHANE PNEUMATIC WHEELS ROLL ACROSS SAND, SNOW AND GRAVEL / DUAL STITCHED SLING WITH TENERA 130 THREAD / PHIFERTEX PLUS SLING MATERIAL RESISTS STRETCHING, MILDEW AND BACTERIA, COLD CRACK AND ULTRAVIOLET DAMAGE AND IS MACHINE WASHABLE.



SPECIFICATIONS

SEAT SIZES	ENGLISH (INCH)	METRIC (CM)
Widths	14, 18, 22	36, 46, 56
Average depth	17	43
Seat-to-floor height	17 – 21	43 – 53
Backrest height	18	46
Telescope footrest knee to heel	15 – 22	38 – 56

WEIGHTS

Capacity	250 – 400 lbs	114 – 182 kg
Chair weight	40 – 65 lbs	18 – 30 kg

ARMREST STYLES

Removable swing away armrest
Adjustable height or padded <i>optional</i>
"Contoured" or specialized <i>optional</i>

Manufactured in the USA,
DeBug has a 10 year warranty on stainless steel frames and a 30 day money back guarantee.

© Copyright Deming Designs Inc., 2001–2016. All rights reserved.

ORDER ONLINE WWW.DEBUGPRODUCTS.COM