



EW-M39

Features include a maximum Speed of 5mph, a 300 lb. capacity, a range of 9 miles, wireless design for trouble free dismantling, an on board battery charger, and delta handlebar.

EW-M39 Specifications

Overall Length

40 “

Overall Width

20”

Overall Height

34”

Wheels:

Front 8”

Rear 8”

Weight w/ Batteries

93 lbs

Weight of Heaviest Piece

39 lbs

Max Speed
5mph

Weight Capacity
300 lbs

Ground Clearance
1.6"

Grade Climbable
8 degree

Turning Radius
44 "

Brake
Electromagnetic

Seat Type
Swivel Padded Foldable

Seat Width
17"

Motor Size
200 W

Battery Size
12AH or 18AH

Battery Weight
20lbs

Battery Charger
On Board

Electronics
On/Off Key Switch, Battery Indicator, Speed Control Knob

