




Sprinter XL4

The Sprinter XL4 is a heavy duty four-wheel scooter with a 350 lb. weight capacity. The unit comes standard with a deluxe captain's seat that can accommodate most users. It has an adjustable tiller for added user comfort and a front head light for added convenience for evening operation. It is available optionally with a full lighting package adding turn indicators and hazard lights. The 12in. solid wheels give the unit a 3in. ground clearance for operation on multiple terrains.

- + 350lb Weight Capacity
- + Slide & Swivel Seating
- + Adjustable Tiller
- + Optional Rear Light Package
- + Adjustable Seat Height

 Burgundy

Model Number: 889B-4
Medicare Billing Code: K0899

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

Technical Specifications

Overall Dimensions (L x W x H)	in	52" x 23" x 30-42"
Number, Size of Tires	Front, Rear	2, 12" (2, 12")
Suggested User Weight ¹	lb	Max. 350
Weight of Heaviest Part ²	lb	108
Weight ³	lb	181
Each Battery Module Weight	lb	32
Battery Size	Type x PCS	12V33-40AH x 2
Charger	Amp	Off board (5A)
Output Power of Motor	HP	1
Maximum Speed	MPH	5
Maximum Climbing Angle/ Safe Climbing Angle ⁴	%(°)	14(8)
Range (Per Charge) ⁶	Mile	25
Turning Radius	in	54"
Ground Clearance	in	3"
Handle Bar	Type	T-Bar
Model #		889B-4
Standard Seat	W x D x BH	SV191821B

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.