

Suction Connecting Tubing

McKesson Suction Connecting Tubing is a clear, non-conductive tubing that connects the suction catheter to the suction canister and is designed to evacuate fluid and debris. Each tubing end has two soft universal female connectors for convenience and provides a secure attachment to any standard male connection. Tubing has minimal coil memory and comes in a wide range of sizes. Selection also includes an economy option in both 3/16" and 1/4" diameter, in lengths ranging up to 20'.

Features and Benefits

- Strong material provides high collapse resistance
- Wide range of sizes and configurations allows for efficient, cost-effective setup
- Clear PVC-formulation with striated walls provides clear visualization, superior strength and secure suction
- Two female connectors; One straight male connector
- 100-foot bulk roll features 16 male connectors
- Single use
- Sterile and non-sterile configurations
- Not made with natural rubber latex



Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them — making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

Suction Connecting Tubing



ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Suction Connecting Tubing, 3/16" x 6', Non-Conductive, Sterile	16-66301	50/cs
Suction Connecting Tubing, 3/16" x 10', Non-Conductive, Sterile	16-66302	50/cs
Suction Connecting Tubing, 3/16" x 12', Non-Conductive, Sterile	16-66303	20/cs
Suction Connecting Tubing, 3/16" x 1 1/2'	16-66304	50/cs
Suction Connecting Tubing, 1/4" x 6'	16-66305	50/cs
Suction Connecting Tubing, 1/4" x 10'	16-66306	50/cs
Suction Connecting Tubing, 1/4" x 12'	16-66307	20/cs
Suction Connecting Tubing, 1/4" x 20'	16-66308	20/cs
Suction Connecting Tubing, Bulk, Non-Conductive, Sterile, 100 Feet, 16 Male Connectors	16-N6100	20/cs