



中国认可
国际互认
检测
TESTING
CNAS L6760

171020340088

Test Report

Report No.: HAPEC2009066302

Page 1 of 5

Applicant Niujian technology co., ltd
Address Room 401, Building 1, Security City, No.919 Qunxian East Road, Yuecheng District, Shaoxing City, Zhejiang Province
Manufacturer Niujian technology co., ltd
Production unit Room 401, Building 1, Security City, No.919 Qunxian East Road, Yuecheng District, Shaoxing City, Zhejiang Province

Sample Information

Sample Name Disposable Vinyl Examination Glove

Sample Brand WTS PPE SUPPLIES

*The information above is provided and confirmed by the applicant.

Sample Received Date Sep.23,2020

Testing Period Sep.23,2020 to Sep.25,2020

Testing Requested As per client's request, to determine the RoHS 6 (Pb, Cd, Hg, Cr⁶⁺, PBBs, PBDEs) in the submitted sample according to RoHS Directive 2011/65/EU Annex II.

Testing Results Please refer to next page(s)

Signer:

Paul. Chen

Authorized signatory



Date: Sep.27,2020

[Website of verification report: CNCA verification platform yz.cnca.cn](http://cnca.cn)

江苏环谱检测技术服务有限公司

JIANGSU HAPEC2009066302 TESTING SERVICE CO.,LTD

地址: 扬州市经济技术开发区吴州东路 198 号

Address: NO.198 Wuzhou East Road, economic and technological development zone, Yangzhou

URL: www.HAPEC2009066302-test.com

E-mail: HAPEC2009066302@HAPEC2009066302-test.com

400-6600-776



中国认可
国际互认
检测
TESTING
CNAS L6760

171020340088

Test Report

Report No.: HAPEC2009066302

Page 2 of 5

Determination of RoHS 6 (unit: mg/kg)

Testing method

- (1) With reference to IEC 62321-3-1:2013. By XRF for measuring;
- (2) With reference to IEC 62321-5:2013. By ICP-OES for measuring;
- (3) With reference to IEC 62321-4: 2013+Amd1: 2017. By ICP-OES for measuring;
- (4) With reference to IEC 62321-7-2:2017. By UV-VIS for measuring;
- (5) With reference to IEC 62321-6:2015. By GC-MS for measuring.

Testing Item(s)	MDL	Limit	Result
Lead (Pb)	2	1000	ND
Cadmium (Cd)	2	100	ND
Mercury (Hg)	2	1000	ND
Chromium(VI) (Cr ⁶⁺)	8	1000	ND
Polybrominated Biphenyls (PBBs)	—	1000	ND
Polybrominated Diphenyl Ethers (PBDEs)	—	1000	ND

Note:

- (1) 1 mg/kg=1 ppm=0.0001%
- (2) MDL=Method Detection Limit
- (3) ND=Not Detected (<MDL)
- (4) "—" =Not Regulated
- (5) Polybrominated Biphenyls, Polybrominated Diphenyl Ethers list, and detection limit (MDL)

Polybrominated Biphenyls (PBBs)	MDL	Polybrominated Diphenyl Ethers (PBDEs)	MDL
Bromobiphenyl	5	Bromobiphenyl ether	5
Dibromobiphenyl	5	Dibromobiphenyl ether	5
Tribromobiphenyl	5	Tribromobiphenyl ether	5
Tetrabromobiphenyl	5	Tetrabromodiphenyl ether	5
Pentabromobiphenyl	5	Pentabromodiphenyl ether	5
Hexabromobiphenyl	5	Hexabromodiphenyl ether	5
Heptabromobiphenyl	5	Heptabromodiphenyl ether	5
Octabromobiphenyl	5	Octabromobiphenyl ether	5
Nonabromobiphenyl	5	Nonabromobiphenyl ether	5
Decabromodiphenyl	5	Decabromobiphenyl ether	5

江苏环谱检测技术服务有限公司

JIANGSU HAPEC2009066302 TESTING SERVICE CO.,LTD

地址: 扬州市经济技术开发区吴州东路 198 号

Address: NO.198 Wuzhou East Road, economic and technological development zone, Yangzhou

URL: www.HAPEC2009066302-test.com

E-mail: HAPEC2009066302@HAPEC2009066302-test.com

400-6600-776


TEST
江苏环谱检测
检验

Test Report

Report No.: HAPEC2009066302

Page 3 of 5

RoHS Testing Flow Chart





中国认可
国际互认
检测
TESTING
CNAS L6760

171020340088

Test Report

Report No.: HAPEC2009066302

Page 4 of 5

Sample photo:



HAPEC2009066302

End of report





中国认可
国际互认
检测
TESTING
CNAS L6760

171020340088

Test Report

Report No.: HAPEC2009066302

Page 5 of 5

注 意 事 项

Notice Items

1、本报告无批准人签字、报告专用章及骑缝章无效。

This report shall be invalid without the signature of the approver, the special seal for the report and the cross-page seal.

2、本报告不得擅自修改、增加或删除。

This report shall not be modified, added or deleted without authorization.

3、报告结果只对本次受检样品负责。若对检测结果有异议，请在报告签发日期后十五天内书面提出，逾期不予受理。

The results of the report are only responsible for the samples tested this time. If there is any objection to the test result, please submit it in writing within fifteen days after the date of issuance of the report. Overdue will not be accepted.

4、未经环谱检测书面同意，不得部分复制本报告，亦不可作为宣传品使用。

Without the written consent of HAPEC2009066302, this report shall not be partially copied or used as publicity materials.

5、当需要结果符合性评价时，若客户无要求，一般不考虑结果的测量不确定度的影响。

When the result conformity evaluation is required, if the customer does not require it, the influence of the measurement uncertainty of the result is generally not considered.

6、相关项目（△）未取得资质认定时，检测数据仅限科研、教学、内部质量控制或研发等活动使用。

When the relevant Items (△) has not been qualified, the test data is limited to scientific research, teaching, internal quality control or research and development activities.

7、本公司出具纸质正版报告与电子数字签名版具备相同效力，当两者内容有差异时以电子数字签名版为准。

The company's paper genuine report has the same legal effect as the electronic digital signature version. When the contents are different, subject to electronic digital signature version.

（具体通用条款详见 <http://www.HAPEC2009066302-test.com/company/tytk>）

(For specific general terms, please see <http://www.HAPEC2009066302-test.com/company/tytk>)