

Prevent Falls & Wandering

New
Low Cost!

Low Cost Monitor Signals Pager For Quiet Wireless Fall Prevention



TL-2020

Key Features:

- Range: up to 300 feet
- Call Button on front of monitor
- Optional AC adapter available (AC-05)
- Works with 433-CMU - 40/60
- Pager Included in Systems
- Batteries not included (3 x AA)



TL-2020P

CareGiver
PAGER™

Signals Pager So They Don't Walk Alone!

Remove Alarm Noise - Volume Adjusts From Silent To The Level You Need

The Wireless Economy Quiet Fall Prevention Monitor to Pager can be set to silent to reduce in-room alert noise. When a resident attempts to get up, a silent wireless signal is sent to a caregiver pager or Economy Central Monitor.

Flexible & Versatile - Works With All Corded Sensors

The Wireless Economy Quiet Fall Prevention Monitor to Pager is designed to work with Corded Bed & Chair Sensor Pads, Floor Mats, and Seat Belts. Also works with Wireless Pager and Central Monitoring Units. This monitor allows caregivers the flexibility to monitor residents based on their unique needs.

Smart Wireless Caregiver Personal Pager™

Caregiver can attend to other duties and receive an alert with use of a pager when someone is getting up from a bed or chair. Can also be used with a floor mat to notify when resident is on their feet.

Anti-Wandering Solutions
SMART®
Quieter Fall Prevention™

Call Us Today! (800) 650-3637 | www.smartcaregiver.com