

Medical Leg Wear Below Knee Open Toe

Below Knee Open Toe Soft Top

- 0865 20-30 mmHg
- 0865S (short) 20-30 mmHg
- 0845 30-40 mmHg
- 0845S (short) 30-40 mmHg

Below Knee Open Toe Dot Top

- 0864 20-30 mmHg
- 0844 30-40 mmHg

Silicone
dot top



0864 BEIGE

Material Information

0865(S) | 0845(S) | 0864 | 0844
Manufactured using 80% nylon and 20% spandex

Product Description

Designed to help relieve moderate and pronounced conditions associated with poor circulation. Truform medical stockings feature classic styling and a smooth knit fabric that helps conceal unsightly conditions. The open toe feature allows for easy toe inspection and permits one standard size to fit a variety of foot dimensions. Color options are basic neutrals for all-occasion wear. For men and women.

Product Features

- Therapeutic graduated compression
- Firm or Extra Firm support
- Classic, concealing smooth knit texture fabric
- Available in soft top and non-constricting, silicone dot top
- Standard and larger sizes

LATEX FREE

Available Colors



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

	CIRCUMFERENCE		LENGTH
	ANKLE	CALF	TO KNEE
SMALL	7" - 8 1/4"	11" - 14"	 Regular: Greater than 15" Short: 15" or less
MEDIUM	8 3/8" - 9 5/8"	13 1/2" - 16"	
LARGE	9 3/4" - 11"	15 1/2" - 18"	
X-LARGE	11 1/8" - 12 3/8"	17 1/2" - 20"	
2X-LARGE	12 1/2" - 13 3/4"	19 1/2" - 22"	
3X-LARGE	12 1/2" - 13 3/4"	21 1/2" - 24"	up to 18"

Classic Medical Knee High Open Toe Compression Socks

Model	Compression	Description	Latex Free
0864	20 - 30 mmHg	Below Knee Dot Top	✓
0865	20 - 30 mmHg	Below Knee Soft Top	✓
0865S (Short)	20 - 30 mmHg	Below Knee Soft Top	✓
0844	30 - 40 mmHg	Below Knee Dot Top	✓
0845	30 - 40 mmHg	Below Knee Soft Top	✓
0845S (Short)	30 - 40 mmHg	Below Knee Soft Top	✓