

A&D Medical ...*Clearly a Better Value*

Multi-User Upper Arm Blood Pressure Monitor UA-767F/FAC



Clinically
Validated

The UA-767 Blood Pressure Monitor has been A&D Medical's most popular blood pressure monitor for 20 years. A&D Medical has taken a great product and made it even better for users with the introduction of the UA-767F series (UA-767F and UA-767FAC).

The UA-767F series has the core features users have come to know, such as clinically validated accuracy and Irregular Heartbeat (IHB) Detection, with the UA-767 line. They additionally feature the ability for up to four users to use the monitor while storing up to 60 results for each user. The AccuFit™ Plus Cuff System that provides a customizable cuff size for users that offers the width of both a medium and large cuff. The cuff also features a flexible inflation technology to provide a more comfortable inflation.

The UA-767FAC has all the features of the UA-767F, with the additional benefit of an AC adapter being included with the monitor for those seeking an alternative to batteries.

New useful features are added to the legacy product. The UA-767F series are 4 users for family use. Each memory bank stores 60 readings. The cuff is an advanced system called AccuFit™ Plus cuff system. The arm circumference almost covers both existing medium and large cuff and a gentle inflation mechanism gives users a comfort measurement.

The performance of the UA-767F series is driven by:

- Advanced Measuring & Reporting Capabilities:
 - o Family-Friendly for Sharing Data – Stores 60 results each with time and date stamp for up to four users
 - o Easy retrieval of your results to provide reporting to your physician
 - o Clinically Validated – Utilizes accuracy that is clinically validated by the latest protocol from the European Society of Hypertension

For questions about A&D products, call us toll free at 1-888-726-9966 and a specially trained representative will assist you, or visit us at www.andmedical.com.

Blood Pressure Monitor UA-767F/FAC

- Risk Factor Indicators:
 - Pressure Rating Indicator™ – Gives a snap shot of your BP classification in accordance with US JNC7 guidelines
 - IHB Detection – Alerts the user to the presence of an irregular heartbeat risk
- Advanced Cuff System:
 - AccuFit™ Plus Cuff – Designed to fit the width of most peoples' arm with 8.6" – 16.5" (22 – 42cm) circumference range
 - Gentle inflation – Minimizes the likelihood of re-inflation with a flexible inflation technology
 - Cuff Wrap Error – Indicates the cuff is not properly wrapped to ensure users get accurate results



PART NUMBER	UPC NUMBER	DESCRIPTION	BATTERIES	AC ADAPTER
UA-767F	0-93764-60336-5	8.6" – 16.5" (22cm – 42cm)	Included	Optional
UA-767FAC	0-93764-60337-2	8.6" – 16.5" (22cm – 42cm)	Optional	Included

Accessories:

- Carrying case included
- Bilingual instruction manual with blood pressure log sheets (English and Spanish)

Technical Specifications:

Dimensions	5.5" (140mm) x 2.3" (60mm) x 4.1" (105mm)
Weight	9.0 oz (255 g) without batteries
Display	1.7" (42mm) x 1.9" (49mm)
Memory	4 users, 60 readings each
Pressurization	Automatic, using micropump
Deflation	Constant air release valve system
Type	Oscillometric
Accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Measurement Range	Pressure: 0 – 299 mmHg Pulse: 40 – 180 beats/minute
Power Source	4 x 1.5V batteries (R6P, LR6 or AA) or AC adaptor TB-233
Operating Condition	+50°F to +104°F (10°C to 40°C) 30% RH to 85% RH
Storage Condition	-4°F to +140°F (-20°C to 60°C) 10% RH to 95% RH

Packaging and Shipping Specifications

Individual Unit (Box)	W x D x H (inches)	WT (lb.)	Individual Item Cube (cubic feet)
UA-767F	5.1 x 3.9 x 7.7	1.57	0.05
UA-767FAC	5.1 x 3.9 x 7.7	1.67	0.05
Unit Case (10 units)	W x D x H (inches)	WT (lb.)	Case Cube (cubic feet)
UA-767F	20 x 11 x 9	16.54	1.08
UA-767FAC	20 x 11 x 9	17.6	1.08

Ask your A&D Medical Sales Representative about other A&D Medical products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Activity Monitors

For questions about A&D products, call us toll free at 1-888-726-9966 and a specially trained representative will assist you, or visit us at