



Urine Collection Bags & Bottles

Uro-Safe® Disposable Vinyl Leg Bags



THE FIRST SIGNIFICANT IMPROVEMENT IN URINARY LEG BAGS IN OVER 50 YEARS!



Uro-Safe® bags have the first ever, built-in, anti-reflux valve that is virtually **self-cleaning** and **interference free**... ensuring uninterrupted flow and problem free use.


Because our anti-reflux valve is enclosed in a hard valve housing, it can't be held shut even when in sitting or reclining positions unlike other inferior products.

Uro-Safe® Why settle for anything less?

LATEX-FREE

STERILE (EO)


Item #77180 Shown




Exclusive "HIGH-FLOW" Anti-Reflux Valve



Improved, Reinforced Strap Eyelets



Easy, Twist-Type Drain Valve Closure



Two-Position Thumb Clamp Closure

CONTENTS:

- 1 Uro-Safe® Disposable Vinyl Urinary Leg Bag† with
- 1 Pair of "Latex-Free" Leg Straps and Buttons (#60062)
- Illustrated Use and Care Instructions

† Please refer to the chart below for specific box contents.

Sold: 1 bag with straps per package/12 packages per case.*
This is a stocked item.
Sterile (EO), Latex-Free.



LEARN MORE



URO-SAFE® BAGS WITH SEMI-TRANSPARENT FRONT AND OPAQUE-WHITE BACK

Item No.	Description	Dimensions†	Capacities‡	Case Item No.
77321*	Large with Thumb Clamp	16.5" x 5.25"	32 fl. oz. (1,000 mL)	7732112 (12/CS)
77320*	Large with Twist Drain	15" x 5.25"	32 fl. oz. (1,000 mL)	7732012 (12/CS)
77181*	Medium with Thumb Clamp	14.38" x 4.75"	18 fl. oz. (550 mL)	7718112 (12/CS)
77180*	Medium with Twist Drain	13" x 4.75"	18 fl. oz. (550 mL)	7718012 (12/CS)

URO-SAFE® BAGS WITH SEMI-TRANSPARENT FRONT AND BACK

Item No.	Description	Dimensions†	Capacities‡	Case Item No.
76321*	Large with Thumb Clamp	16.5" x 5.25"	32 fl. oz. (1,000 mL)	7632112 (12/CS)
76320*	Large with Twist Drain	15" x 5.25"	32 fl. oz. (1,000 mL)	7632012 (12/CS)
76181*	Medium with Thumb Clamp	14.38" x 4.75"	18 fl. oz. (550 mL)	7618112 (12/CS)
76180*	Medium with Twist Drain	13" x 4.75"	18 fl. oz. (550 mL)	7618012 (12/CS)

* These items are not sold individually, sold only by the case.

† Approximate dimensions in inches(") and centimeters(cm) including adaptors and drain valves.

‡ Approximate volumes in fluid ounces (fl. oz.) and milliliters (mL)

* SOLD ONLY BY THE CASE

* SOLD ONLY BY THE CASE