

Shown in Badlands Walnut Fabric



Elevate Your Lifestyle!



Shown in Badlands Walnut Fabric

When looking for true elegance, look no further than the **Elegance** line from **VivaLift!™** Power Recliners. Experience comfort in a whole new way with power lumbar, which offers personalized ergonomic support. The full-width headrest lifts your head, neck and shoulders and infinite lay flat positioning ensures maximum relaxation. Customize your chair with your choice of fashion forward fabrics and make your home truly elegant with the Elegance line!

(US) 800-800-8586 | (CAN) 888-570-1113 pridemobility.com | prideliftchairs.com









Position	Relaxed Reclined Infinite Lay Flat Upward Lift 5'4" - 6'
Weight Capacity	375 lbs.
Height	43"
Seat to Floor	20″
Seat Depth	20″
Seat Width	22″
Top of Back to Seat	27″
Back Style	Blown fiber, Pillow back
Fully Padded Chaise	Yes
Footrest Extension	Standard
Distance from the Wall	14"
Overall Chair Width	35 ¹ / ₂ "
Chair Weight	152 lbs.
Warranty	7-Year Pro-rated



STANDARD FABRIC SELECTION GUIDE











UPWARD LIFT

STANDARD FEATURES

POWER HEADREST POWER LUMBAR USB REMOTE LITHIUM BATTERY BACKUP INFINITE LAY FLAT

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

ICON KEY = Indicates shortest available lead time

LC-ELEGANCE-VIVA-SS-11-03-17

(US) 800-800-8586 | (CAN) 888-570-1113 pridemobility.com | prideliftchairs.com

