

Winnie Lite Supreme/Go Lite 3 Wheel Aluminum Rollator

HCPCS:




FIGURE A




FIGURE B



FIGURE C



FIGURE D

- Special loop lock made of internal aluminum casting operates easily and ensures safety (Figure A)
- 7.5" casters with soft-grip tires are ideal for indoor and outdoor use
- Comes standard with carry pouch
- Handles easily adjust in height with self-threading knob (Figure B)
- Easy, one-hand folding

PRODUCTS

	Description
RTL199	Winnie Lite Supr.w/TanPlaidPou,1/cs

SPECIFICATIONS

	Description
Carton Shipping Weight	15 lbs.
Handle Height	32" - 38"
Unit Dimensions	26" (W) x 22.75" (L)
Warranty	Limited Lifetime
Weight	11 lbs
Weight Capacity	300 lbs

PARTS

	Description
10208KN	Knob,Ht.Adj.Rollator,w/Attch.1/ea, EA
808	BLACK TOTE,1/ea, EA
9500001232	Brake Cable/Adjust,Lower,1/ea, EA
9500004101	Cable,Fastner,1/ea, EA
9500002211	Axle Bearing,Rollator,1ea, EA
9503HBRAKE	BrakeAssm.,ThreeWheelRollator,1/ea, EA
9500H30001	Rollator Wheel,8"6 Spoke, 1/ea, EA
9500002051	Lower Pivot Bolt,197,1/ea, EA
950330049901	Hand Grip, Right, 3-wheel Rollator, 1/ea, EA
950330039901	Hand Grip, Left, 3-wheel Rollator, 1/ea, EA
950000199901	Generic Brake Cable Only for Rollat,1/ea, EA