



Published on *WoundSource* (<https://www.woundsource.com>)

[Home](#) > Simpurity™ Silver Alginate Dressing Pads

Simpurity™ Silver Alginate Dressing Pads [1]



Simpurity™ Silver Alginate Dressing Pads have gel forming properties, silver particles, conform to wound contours and have strong wet tensile strength.

Safe n' Simple [2]

Safe n' Simple is a developer and marketer of ostomy and wound care products designed to enhance patients' quality life. Products are safe, effective, reimbursement-friendly and easy to use.

Website: www.sns-medical.com

Email: marketing@sns-medical.com

Phone: (248) 214-4877

Toll-free: (844) 767-6334

Fax: (248) 698-3454

Benefits:

- Strong wet tensile strength
- Easy to remove from the wound bed
- High absorption capacity
- Conforms to wound contours
- Reduces bioburden
- Gel forming properties

Indications:

Simpurity™ Silver Alginate Wound Dressings are indicated for moderately to highly exuding wounds, such as pressure ulcers, leg ulcers, cavity wounds and post-debridement hemostasis.

Contraindications:

Contraindicated for use on third-degree burns, dry wounds and surgical implantations.

Storage Requirements:

Store at room temperature.

Active Ingredients:

Calcium alginate with silver particles

How Supplied/Sizing:

2"x2", 4"x5", 10 per box.

HCPCS Code (text):

A6196, A6197

Recommended Use:

Acute Wounds
Cavity Wounds
Chronic Wounds
Deep Wounds
Diabetic Foot
Infected Wounds
Moderate/Highly Exudating Wounds
Pressure Ulcers
Surgical Wounds
Venous Ulcers

Mode of Use/Application:

Wet Simpurity™ Silver Alginate Wound Dressing with saline or water. Apply directly to the whole wound bed. If dry, wet with saline or water and apply into wound bed. Apply secondary dressing.

Additional Recommended Dressings:

Secondary dressing needed.

Product features:

Absorptive dressing format

Alginate dressing format

Antimicrobial effects up to 48 hours

Antimicrobial effects up to 3 days

Antimicrobial effects up to 7 days

Bacterial binding

Bacterial trapping

Bactericidal

Bactericidal broad-spectrum

Bacteriostatic •

Collagen dressing format

Composite dressing format

Contact layer dressing format

Foam dressing format

Fungicidal

Gauze or non-woven dressing format

Gelling fiber dressing format

Hydrogel (amorphous) dressing format

Hydrogel (impregnated) dressing format

Hydrogel (sheet) dressing format

Impermeable to outside contaminants

Impregnated dressing format

Microbicidal

Moisture vapor permeable

Secondary dressing required •

Silicone-based adhesive

Sporicidal

Super absorbent polymers

Sustained release formula

Transparent film dressing format

Virucidal

Wound filler dressing format

Other features:

Educational Material Available

Free Samples/Trials Available •

Published Clinical Article Available

Source URL (modified on 2020-06-07 17:23): <https://www.woundsource.com/product/simpurity-silver-alginate-dressing-pads>

Links

[1] <https://www.woundsource.com/product/simpurity-silver-alginate-dressing-pads>

[2] <https://www.woundsource.com/company/safe-n-simple>