**Product Name:** CaviWipes HP

**Application:** All fill sizes, see www.metrex.com for catalogue numbers and safety data sheets.

**Company:** METREX™ RESEARCH, LLC
28210 Wick Rd.
Romulus, MI 48174
U.S.A.
1-800-841-1428

**Date of Disclosure**
January 24, 2022

**Ingredients**

| Ingredient Name | CAS | Present on Designated List | Function (HCPA) |
|---|---|---|---|
| Purified Water | 7732-18-5 | No | Diluent |
| Proprietary Aromatic Alcohol | **Withheld for CBI** | EU Fragrance Allergen | Solvent |
| Organic Acid | **Withheld for CBI** | No | **Solvent** |
| Hydrogen Peroxide | 7722-84-1 | No | Active ingredient |
| Amines, coco alkyldimethyl, N-oxide | **Withheld for CBI** | No | Surfactant |
| Alcohols ethoxylated | **Withheld for CBI** | No | Surfactant |
| Phosphoric Acid | 7664-38-2 | OEHHA REL ; CA Non-Cancer Hazards | pH adjustor |

**Contact for Additional Ingredient Information:** 1-800-841-1428

**Links**

OEHHA REL/CA Non-Cancer Hazards https://oehha.ca.gov/air/general-info/oehha-acute-8-hour-and-chronic-reference-exposure-level-rel-summary

EU Fragrance Allergens https://ec.europa.eu/health/scientific_committees/opinions_layman/perfume-allergies/en/l-3/1-introduction.htm